Edwards

Edwards eSheath Introducer Set

Edwards eSheath -införingsset

Edwards eSheath indføringssæt

DIRECTORY

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Instructions For Use

1.0 Device Description

The Edwards eSheath Introducer Set contains:

- a) an expandable sheath (eSheath) (Fig 1) that provides access into the target vessel while maintaining hemostasis and temporarily enlarges its diameter to allow for passage of a device.


Figure 1

Model	eSheath I.D. (unexpanded)	eSheath O.D. (unexpanded)
9610ES14	14F (4.6 mm)	6.0 mm
9610ES16	16F (5.3 mm)	6.7 mm

Edwards, Edwards Lifesciences, the stylized E logo, Edwards eSheath, and eSheath are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

THV252

Figure 2

- b) two dilators (Fig 2) with hydrophilic coating that can either be used to dilate the vessel to accommodate the sheath and/or facilitate entry and trackability of the sheath into the vessel.

2.0 Intended Use

The Edwards eSheath Introducer Set is intended for introduction of interventional devices into the vascular system. The product is intended for use by physicians trained and experienced in interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.

3.0 Contraindications

This product is contraindicated for patients with tortuous or calcified vessels that would prevent safe entry of the dilators and eSheath.

4.0 Warnings

The devices are designed, intended, and distributed for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.

The Edwards eSheath Introducer Set must be used with a compatible 0.035" (0.89 mm) guidewire to prevent vessel injury.

Do not mishandle the device or use it if the packaging or any components are not sterile, have been opened or are damaged (i.e., kinked or stretched, etc.), or the expiration date has elapsed.

5.0 Precautions

- The sheath temporarily enlarges to allow the passage of devices; ensure that the vasculature can accommodate the maximum diameter of the expanded sheath.
- When inserting, manipulating or withdrawing a device through the sheath, always maintain orientation of the sheath position.
- When puncturing, suturing or incising the tissue near the sheath, use caution to avoid damage to the sheath.

6.0 Potential Adverse Events

Complications associated with standard catheterization and use of angiography include, but are not limited to, allergic reaction to anesthesia or to contrast media; injury including perforation or dissection of vessels; injury at the site of access that might require vessel repair; thrombosis and/or plaque dislodgment which may result in emboli formation; distal vessel obstruction; stroke; ischemia and/or death.

7.0 Directions for Use

1. Visually inspect device components for damage.
2. Flush the dilators using heparinized saline through the guidewire lumen.
3. Flush the sheath using heparinized saline through the flush port; close the flush port.
4. Hydrate the length of the introducer/dilators and sheath with heparinized saline to activate the hydrophilic coating.
5. Insert one dilator completely into the sheath.
6. Using standard catheterization techniques, gain access to the vessel and dilate as necessary with the other dilator to accommodate the sheath.
7. Orient the sheath with the Edwards logo facing up and maintain orientation throughout the procedure. Fully insert the sheath assembly using standard technique and advance into the vessel while following its progression under fluoroscopy.

NOTE: The proximal tapered end of the sheath working length is larger in diameter.

8. Suture the eSheath into place using the suture rings and remove the dilator from the sheath.
9. Insert the device into the eSheath.

NOTE: The eSheath should be intermittently flushed with heparinized saline throughout the procedure, per standard interventional technique.

10. After the completion of the procedure and removal of the device, remove the suture, and then remove the sheath entirely without torquing and do not reinsert.

8.0 How Supplied

The Edwards eSheath Introducer Set is supplied in a pouch and sterilized with ethylene oxide.

9.0 Storage

The Edwards eSheath Introducer Set should be stored in a cool, dry place.

10.0 Device Disposal

Used sheath sets may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

These products are manufactured and sold under one or more of the following US patents: US Patent No. 8,690,936; 8,790,387; 9,301,840; 9,301,841; and corresponding foreign patents. Additional patents are pending.

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© Copyright 2019, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC