Edwards

Ascendra Balloon Aortic Valvuloplasty Catheter

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English

Instructions for Use – Transapical/Transaortic Approach

For single use only

The maximum guidewire diameter which may be used with the 9100BAVC is 0.035" (0.89 mm).

Minimum sheath compatibility is 14F (4.62 mm).

Minimum guidewire length is 100 cm.

STERILE EO

1. Device Description

The Ascendra balloon aortic valvuloplasty catheter (Figure 1) is a coaxial designed catheter with a distal inflatable balloon. Two radiopaque marker bands indicate the dilating section of the balloon and aid in balloon placement. At the proximal end of the catheter, there is a standard "Y" connector for balloon inflation and a guidewire lumen. An additional balloon extension tubing is optional for user preference. The balloon is inflated by injecting a diluted contrast medium solution through the luer port (marked "BALLOON") on the "Y" connector.

Contents of unopened, undamaged package are:

- Sterile
- Nonpyrogenic

**Table 1 Compliance Table
Pressure and/or Volume vs. Diameter**

Volume (ml)	Volume w/ Balloon Extension (ml)	Applied Pressure ATM (kPa)	Diameter (mm) Ø ±10%
10.0	11.0	2.0 (202.7)	18.3
10.5	11.5	2.5 (253.3)	18.6
11.0	12.0	3.0 (304.0)	18.8
11.5	12.5	3.5 (354.6)	19.0
12.0	13.0	4.0 (405.3)	19.2
12.5	13.5	4.5 (456.0)	19.4
13.0	14.0	5.0 (506.6)	19.5
13.3	14.3	5.5 (557.3)	19.6
13.5	14.5	6.0 (608.0)	19.8
13.8	14.8	6.5 (658.6)	19.9
14.0	15.0	7.0 (709.3)	20.0
14.3	15.3	7.5 (759.9)	20.1
14.5	15.5	8.0 (810.6)	20.2

Note: Nominal sizing must be based on volume only.

Note: For proper volume sizing, the balloon aortic valvuloplasty catheter should be used with the inflation device provided by Edwards Lifesciences.

Edwards, Edwards Lifesciences, the stylized E logo, and Ascendra are trademarks of Edwards Lifesciences Corporation. All other trademarks are property of their respective owners.

The upper figures in bold face represent the balloon diameter at Nominal Inflation Volume, Pressure and Diameter.

The lower figures in bold face represent the balloon diameter at Rated Burst Volume, Pressure and Diameter.

2. Indications

The 9100BAVC is indicated for valvuloplasty of a stenotic aortic valve prior to implantation of an Edwards transcatheter heart valve.

3. Contraindications

The device is contraindicated for patients with:

- Non-valvular aortic stenosis
- Congenital aortic stenosis, unicuspid or bicuspid aortic valve
- Non-calcific acquired aortic stenosis
- Presence of mitral bioprosthesis
- Evidence of intracardiac mass, thrombus, or vegetation
- Severe deformation of the chest
- Severe coagulation problems
- Active bacterial endocarditis or other active infections
- Patients unable to tolerate anticoagulation therapy
- Recent (within 6 months) cerebrovascular accident (CVA)
- Hypertrophic cardiomyopathy with or without obstruction (HOCM)

The device is not to be used if the implanting physician believes its use would be contrary to the best interest of the patient.

There are no known contraindications for the device other than standard risks associated with balloon aortic valvuloplasty.

4. Warnings

- This device is intended for single use only. DO NOT re-sterilize and/or reuse it. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.
- Use the catheter prior to the "Use By" date specified on the package.
- Catheter balloon inflation diameter must be carefully considered in selecting a particular size for any patient. The inflated balloon diameter should not be significantly greater than the annulus diameter being pre-dilated.
- Catheter is not intended for redilation of deployed transcatheter heart valves.
- While the catheter is exposed within the body, advancement and retrieval should not be done without the aid of fluoroscopic equipment. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum.
- Do not use if package is opened or damaged.

5. Precautions

- Prior to use, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.
- Dilation procedure should be conducted under fluoroscopic guidance with appropriate fluoroscopic equipment.
- Careful attention must be paid to the maintenance of tight catheter connections and aspiration before proceeding to avoid air introduction into the system.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Minimum acceptable sheath French size is printed on the package label. Do not attempt to pass the Ascendra balloon aortic valvuloplasty catheter through a smaller sheath introducer than indicated on the label.
- Before withdrawing the catheter from the sheath, the balloon must be fully deflated.

6. Potential Adverse Events

- Perforation
- Conduction System Injury
- Thromboembolic Events
- Hematoma
- Cardiovascular Injury and Infundibulum
- Arrhythmia Development
- Valvular Tearing or Trauma
- Inflammation
- Infection
- Chordae Damage or Trauma

7. Directions for Use

7.1 Ascendra Balloon Aortic Valvuloplasty Catheter Preparation

Step	Procedure
1	Remove catheter from package.
2	Remove balloon cover and inspect the catheter for damage.
3	Flush the guidewire lumen (located at the "Y" connector) with saline.
4	Wipe outside of catheter and balloon extension using saline solution.
5	Prepare an inflation device with a mixture of contrast medium and saline (15:85 medium to saline dilution).
6	If using the balloon extension tubing, flush tubing with saline.
7	Attach balloon extension tubing to "Y" connector.
8	Attach the inflation device to the inflation port of the catheter.
9	Induce a negative pressure to remove any air from the balloon and inflation lumen. Repeat until all air is expelled. Close stopcock.
10	Remove inflation device from catheter. Fill inflation device with the appropriate volume for desired balloon diameter (Table 1).
11	Attach filled inflation device to stopcock and open stopcock. Allow the inflation lumen to fill with the diluted contrast medium and maintain at neutral pressure.

7.2 Catheter Use

Step	Procedure
1	Wipe the exposed guidewire with normal saline to remove residual tissue or contrast medium.
2	Advance the Ascendra balloon aortic valvuloplasty catheter over the guidewire to the intended site. An introducer should be utilized to facilitate catheter insertion.
3	Under fluoroscopic guidance, continue to advance the balloon aortic valvuloplasty catheter to the intended site. Position balloon markers at the intended site.
4	Once desired position is achieved, begin rapid pacing.
5	Perform dilation using inflation device. Patient monitoring is required during dilatation. Balloon can be either partially or fully inflated to achieve dilatation. DO NOT EXCEED THE RATED BURST VOLUME (TABLE 1).
6	Fully inflate and rapidly deflate balloon by drawing a vacuum with an inflation device.
7	When the balloon is fully deflated, the pacemaker may be turned off.

Step	Procedure
8	Gently withdraw the deflated catheter into the sheath, then remove the sheath and balloon catheter. Note: The greater the vacuum applied and held during withdrawal, the lower the deflated balloon profile. Note: If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit under fluoroscopic guidance, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
9	Discard the Ascendra balloon aortic valvuloplasty catheter after use.

8. Storage

The model 9100BAVC must be stored in a cool, dry place.

9. Device Disposal

The device may be disposed of in the same manner that the hospital waste and biohazardous materials are handled. There are no special or unusual risks related to the disposal of the device.

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© Copyright 2019, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC