Edwards

Edwards SAPIEN 3 Ultra System
Edwards SAPIEN 3 Ultra Transcatheter Heart Valve
Edwards Commander Delivery System
Transfemoral, Subclavian / Axillary

Edwards SAPIEN 3 Ultra system
Edwards SAPIEN 3 Ultra transkateter hjerteklaff
Edwards Commander leveringssystem
transfemoralt, subklavisk/aksillært

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English

Instructions for Use

Implantation of transcatheter heart valves should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in balloon aortic valvuloplasty. It is at the physician's discretion to choose the appropriate access route to implant the THV based on the patient anatomy and associated risks.

1.0 Device Description

• **Edwards SAPIEN 3 Ultra System**

The Edwards SAPIEN 3 Ultra system consists of the Edwards SAPIEN 3 Ultra transcatheter heart valves and delivery systems.

• **Edwards SAPIEN 3 Ultra Transcatheter Heart Valve (Figure 1)**

The Edwards SAPIEN 3 Ultra transcatheter heart valve (THV) is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) inner and outer fabric skirts. The leaflets are treated according to the Carpentier-Edwards ThermoFix process.

The THV is intended to be implanted in a native annulus size range associated with the three-dimensional area of the aortic annulus measured at the basal ring during systole:

Table 1

Native Valve Annulus Size (TEE)*	Native Valve Annulus Size (CT)		THV Size
	Native Annulus Area (mm ²)	Area-derived diameter (mm)	
16-19 mm	273 - 345	18.6 - 21.0	20 mm
18-22 mm	338 - 430	20.7 - 23.4	23 mm
21-25 mm	430 - 546	23.4 - 26.4	26 mm

THV size recommendations are based on native valve annulus size, as measured by transesophageal echocardiography (TEE) or computed tomography (CT). Patient anatomical factors and multiple imaging modalities should be considered during THV size selection.

Note: Risks associated with undersizing and oversizing should be considered to minimize the risk of paravalvular leak, migration, and/or annular rupture.

*Due to limitations in two-dimensional images, 2-D TEE imaging should be supplemented with 3-D area measurements.

Edwards, Edwards Lifesciences, the stylized E logo, Carpentier-Edwards, Edwards Commander, Edwards SAPIEN, Edwards SAPIEN 3, Edwards SAPIEN 3 Ultra, eSheath, Qualcrimp, SAPIEN, SAPIEN 3, SAPIEN 3 Ultra, and ThermoFix are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

• Edwards Commander Delivery System (Figure 2)

The Edwards Commander delivery system (Figure 2) facilitates the placement of the bioprosthesis. It consists of a Flex Catheter to aid in valve alignment to the balloon, tracking, and positioning of the THV. The delivery system includes a tapered tip to facilitate crossing of the native valve. The handle contains a Flex Wheel to control flexing of the Flex Catheter, and a Balloon Lock and Fine Adjustment Wheel to facilitate valve alignment and positioning of the valve within the native annulus. A stylet is included within the guidewire lumen of the delivery system. The Balloon Catheter has radiopaque Valve Alignment Markers defining the working length of the balloon. A radiopaque Center Marker in the balloon is provided to help with valve positioning. A radiopaque Triple Marker proximal to the balloon indicates the Flex Catheter position during deployment.

The inflation parameters for THV deployment are:

Table 2

Model	Nominal Balloon Diameter	Nominal Inflation Volume	Rated Burst Pressure (RBP)
9610TF20	20 mm	11 ml	7 atm (709 kPa)
9610TF23	23 mm	17 ml	7 atm (709 kPa)
9610TF26	26 mm	23 ml	7 atm (709 kPa)

• Qualcrimp Crimping Accessory (Figure 3)

The Qualcrimp crimping accessory is used during crimping of the THV.

• Loader (Figure 4)

The loader allows for the delivery of the crimped valve through the hemostasis valves of the sheath.

• Crimper and Crimp Stopper (Figure 5)

The crimper reduces the diameter of the THV to mount it to the delivery system. The crimper is comprised of a compression mechanism that is closed with a handle located on the housing. The crimper is used with a 2-piece crimp stopper to correctly crimp the THV.

• Edwards Sheath

Refer to the Edwards sheath instructions for use for device description.

• Inflation Devices

An Inflation device with locking mechanism is used during native valve predilation and THV deployment.

Note: For proper volume sizing, the Edwards Commander delivery system and the Edwards transfemoral balloon catheter should be used with the inflation device provided by Edwards Lifesciences.

2.0 Indications

The Edwards SAPIEN 3 Ultra system is indicated for use in patients with heart disease due to native calcific aortic stenosis at any or all levels of surgical risk for open heart surgery.

3.0 Contraindications

Use of the Edwards SAPIEN 3 Ultra system is contraindicated in patients with:

- Evidence of intracardiac mass, thrombus, vegetation, active infection or endocarditis.
- Inability to tolerate anticoagulation/antiplatelet therapy.

4.0 Warnings

- The devices are designed, intended, and distributed STERILE for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- Correct sizing of the THV is essential to minimize the risk of paravalvular leak, migration, and/or annular rupture.
- The physician must verify correct orientation of the THV prior to its implantation; the inflow (outer skirt end) of the THV should be oriented distally towards the tapered tip to prevent the risk of severe patient harm.
- Accelerated deterioration of the THV may occur in patients with an altered calcium metabolism.
- Observation of the pacing lead throughout the procedure is essential to avoid the potential risk of pacing lead perforation.
- The THV must remain hydrated at all times and cannot be exposed to solutions, antibiotics, chemicals, etc. other than its shipping storage solution and sterile physiologic saline solution to prevent leaflet damage that may impact valve functionality. THV leaflets mishandled or damaged during any part of the procedure will require replacement of the THV.
- Patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, and/or polymeric materials may have an allergic reaction to these materials.
- Do not use the THV if the tamper evident seal is broken, as sterility may be compromised.
- Do not use the THV if the temperature indicator has been activated, as valve function may be compromised.
- Do not use the THV if the expiration date has elapsed, as either sterility or valve function may be compromised.
- Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging sterile barriers and any components have been opened or damaged, cannot be flushed, or the expiration date has elapsed.
- Access characteristics such as severe obstructive or circumferential calcification, severe tortuosity, vessel diameters less than 5.5 mm (for size 20, 23 and 26 mm SAPIEN 3 Ultra transcatheter heart valve) may preclude safe placement of the sheath and should be carefully assessed prior to the procedure.

5.0 Precautions

- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Material Safety Data Sheet available from Edwards Lifesciences.
- The safety and effectiveness of the THV implantation has not been established in patients who have:
 - Congenital unicuspid aortic valve
 - Pre-existing prosthetic heart valve or prosthetic ring in any position
 - Severe ventricular dysfunction with ejection fraction < 20%
 - Hypertrophic cardiomyopathy with or without obstruction
 - Aortic stenosis characterized by a combination of AV low flow, low gradient
- If a significant increase in resistance occurs when advancing the catheter through the vasculature, stop advancement and investigate the cause of resistance before proceeding. Do not force passage, as this could increase the risk of vascular complications.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- THV recipients should be maintained on anticoagulant/antiplatelet therapy to minimize the risk of valve thrombosis or thromboembolic events, as determined by their physicians.
- Long-term durability has not been established for the THV. Regular medical follow-up is advised to evaluate valve performance.
- Based on the treating physician's consideration of risks and benefits, the THV may be implanted in relatively young patients, although the longer-term durability is still the subject of ongoing clinical research.
- Do not overinflate the deployment balloon, as this may prevent proper valve leaflet coaptation and thus impact valve functionality.
- The risks of subclavian/axillary access are low and acceptable, but subclavian/axillary access should be considered when the physician determines there is an increased risk associated with transfemoral access.
- For left axillary approach, a left subclavian takeoff angle $\sim \geq 90^\circ$ from the aortic arch causes sharp angles, which may be responsible for potential sheath kinking, subclavian/axillary dissection and aortic arch damage.
- For left / right axillary approach, ensure there is flow in the Left Internal Mammary Artery (LIMA) / Right Internal Mammary Artery (RIMA), respectively, during procedure and monitor pressure in homolateral radial artery.

6.0 Potential Adverse Events

Potential risks associated with the overall procedure including access, cardiac catheterization, local and/or general anesthesia:

- Allergic reaction to antithrombotic therapy or contrast medium or anesthesia
- Anemia
- Aneurysm
- Angina
- Arrhythmias including ventricular fibrillation (VF) and ventricular tachycardia (VT)
- AV fistula or pseudoaneurysm
- Cardiogenic shock
- Compartment syndrome
- Death
- Dissection: aortic or other vessels
- Emboli, distal (air, tissue or thrombotic emboli)
- Hematoma
- Hypertension or hypotension
- Inflammation
- Myocardial ischemia or infarction
- Pain or changes at the access site
- Perforation or rupture of cardiac structures
- Perforation or rupture of vessels
- Pericardial effusion or cardiac tamponade
- Peripheral ischemia or nerve injury or brachial plexus injury
- Pulmonary edema
- Renal insufficiency or renal failure
- Respiratory insufficiency or respiratory failure
- Syncope
- Thoracic bleeding
- Vasovagal response
- Vessel spasm
- Vessel thrombosis/occlusion
- Vessel trauma requiring surgical repair or intervention

Additional potential risks associated with the TAVR procedure, the bioprosthesis, and the use of its associated devices and accessories include:

- Allergic/immunologic reaction to the implant
- Atrial fibrillation/Atrial flutter
- Bleeding requiring transfusion or intervention
- Cardiac arrest

- Cardiac failure or low cardiac output
- Cardiogenic shock
- Conduction system injury (defect) including AV block, which may require a permanent pacemaker
- Coronary occlusion
- Dissection, rupture, trauma of the aortic annulus and surrounding structures including ascending aorta, coronary ostia and ventricular septum
- Emergency cardiac surgery
- Hemolysis
- Infection, fever, septicemia, abscess, endocarditis
- Injury to mitral valve
- Mechanical failure of delivery system, and/or accessories, including balloon rupture and tip separation
- Silent cerebral ischemia, stroke, transient ischemic attack, cognitive impairment
- Structural valve deterioration (wear, fracture, calcification, stenosis)
- Valve deployment in unintended location
- Valve explants
- Valve migration, malposition or embolization requiring intervention
- Valve regurgitation, paravalvular or transvalvular
- Valve thrombosis

7.0 Directions for Use

7.1 System Compatibility

	20 mm System	23 mm System	26 mm System
Product Name	Model		
Edwards SAPIEN 3 Ultra Transcatheter Heart Valve	9750TFX (20 mm)	9750TFX (23 mm)	9750TFX (26 mm)
Edwards Commander Delivery System	9610TF20	9610TF23	9610TF26
Sheath provided by Edwards Lifesciences			
Inflation device, Qualcrimp Crimping Accessory, Crimp Stopper and Loader provided by Edwards Lifesciences			
Edwards Crimper	9600CR		

Additional Equipment

- Standard cardiac catheterization lab equipment
- Fluoroscopy (fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal or transthoracic echocardiography capabilities
- Exchange length 0.035 inch (0.89 mm) extra-stiff guidewire
- Pacemaker (PM) and pacing lead

- Edwards Transfemoral Balloon catheter or equivalent
- Sterile rinsing bowls; sterile physiological saline solution; sterile heparinized saline solution, and diluted radiopaque contrast medium (15:85 medium to saline dilution)
- Sterile table for THV and device preparation
- 20 cc syringe or larger
- 50 cc syringe or larger
- High-pressure 3-way stopcock (x2)

7.2 THV Handling and Preparation

Follow sterile technique during device preparation and implantation.

7.2.1 THV Rinsing Procedure

Before opening the valve jar, carefully examine for evidence of damage (e.g. a cracked jar or lid, leakage, or broken or missing seals).

CAUTION: If the container is found to be damaged, leaking, without adequate sterilant, or missing intact seals, the THV must not be used for implantation, as sterility may be compromised.

Step	Procedure
1	Set up two (2) sterile bowls with at least 500 ml of sterile physiological saline to thoroughly rinse the glutaraldehyde sterilant from the THV.
2	Carefully remove the valve/holder assembly from the jar without touching the tissue. Verify the valve serial identification number with the number on the jar lid and record in the patient information documents. Inspect the valve for any signs of damage to the frame or tissue.
3	Rinse the THV as follows: - Place the THV in the first bowl of sterile, physiological saline. Be sure the saline solution completely covers the THV and holder. - With the valve and holder submerged, slowly agitate (to gently swirl the valve and holder) back and forth for a minimum of 1 minute. - Transfer the THV and holder to the second rinsing bowl of sterile physiological saline and gently agitate for at least one more minute. Ensure the rinse solution in the first bowl is not used. - The valve should be left in the final rinse solution until needed to prevent the tissue from drying. CAUTION: Do not allow the valve to come into contact with the bottom or sides of the rinse bowl during agitation or swirling in the rinse solution. Direct contact between the identification tag and valve is also to be avoided during the rinse procedure. No other objects should be placed in the rinse bowls. The valve should be kept hydrated to prevent the tissue from drying.

7.2.2 Prepare the System

Step	Procedure
1	Visually inspect all the components for damage. Ensure the delivery system is fully unflexed and the balloon catheter is fully advanced in the flex catheter. WARNING: To prevent possible damage to the balloon shaft, ensure that the proximal end of the balloon shaft is not subjected to bending.
2	Flush the delivery system with heparinized saline through the flush port.
3	Remove the distal balloon cover from the delivery system. Remove the stylet from the distal end of the guidewire lumen and set aside.
4	Flush the guidewire lumen with heparinized saline. Insert the stylet back into the guidewire lumen. Note: Failure to replace the stylet in the guidewire lumen may result in damage to the lumen during the THV crimping process.
5	Place the delivery system into the Default Position (end of strain relief is aligned between the two white markers on the balloon shaft) and make sure that the flex catheter tip is covered by the proximal balloon cover.
6	Unscrew the loader cap from the loader and flush the loader cap with heparinized saline.
7	Place the loader cap onto the delivery system with the inside of the cap oriented towards the distal tip. Fully advance the balloon catheter in the flex catheter. Peel off the proximal balloon cover over the blue section of the balloon shaft.
8	Attach a 3-way stopcock to the balloon inflation port. Fill a 50 cc or larger syringe with 15-20 ml of diluted contrast medium and attach to the 3-way stopcock.
9	Fill the inflation device with excess volume of diluted contrast medium relative to the indicated inflation volume. Lock and attach to the 3-way stopcock. Close stopcock to the inflation device.
10	Pull vacuum with the syringe to remove air. Slowly release the plunger to ensure that the contrast medium enters the lumen of the delivery system. Repeat until all air bubbles are removed from the system. Leave zero-pressure in the system. WARNING: Ensure there is no residual fluid left in the balloon to avoid potential difficulty with valve alignment during the procedure. Close stopcock to the delivery system.
11	Rotate the knob of the inflation device to remove the contrast medium into the syringe and achieve the appropriate volume required to deploy the THV. Close the stopcock to the syringe and remove syringe.

Step	Procedure
12	Verify that the inflation volume in the inflation device is correct. CAUTION: Maintain the inflation device in the locked position until THV deployment to minimize the risk of premature balloon inflation and subsequent improper THV deployment.

7.2.3 Mount and Crimp the THV on the Delivery System

Step	Procedure
1	Completely submerge the Qualcrimp crimping accessory in a bowl of 100 ml physiological saline. Gently compress until fully saturated. Swirl for a minimum of 1 minute. Repeat this process in a second bowl.
2	Remove the THV from the holder and remove the ID tag.
3	Rotate the crimper handle until the aperture is fully open. Attach the 2-piece Crimp Stopper to the base of the crimper and click into place.
4	If necessary, partially crimp the THV in the crimper until it snugly fits inside the Qualcrimp crimping accessory. Note: Partial crimping is not necessary for the 20 mm valve.
5	Place the Qualcrimp crimping accessory over the THV aligning the edge of the Qualcrimp crimping accessory with the outflow of the THV.
6	Place the THV and Qualcrimp crimping accessory in crimper aperture. Insert the delivery system coaxially within the THV 2-3 mm distal to the blue balloon shaft (in the Valve Crimp Section) of the delivery system with the inflow of the THV towards the distal end of the delivery system.
7	Center the balloon shaft coaxially within the THV. Crimp the THV until it reaches the Qualcrimp stop.
8	Remove the Qualcrimp crimping accessory from the THV and Qualcrimp stop from the Crimp Stopper, leaving the Final Stop in place.
9	Center the THV within the crimper aperture. Fully crimp the THV until it reaches the Final Stop and hold for 5 seconds. Repeat this crimp step two (2) more times for a total of 3 crimps. Note: Ensure that the Valve Crimp Section is coaxial within the THV.
10	Pull the balloon shaft and engage the Balloon Lock so the delivery system is in Default Position.
11	Flush the loader with heparinized saline. Immediately advance the THV into the loader until it is completely inside the loader. CAUTION: The THV should not remain fully crimped and/or in the loader for over 15 minutes, as leaflet damage may result and impact valve functionality.

Step	Procedure
12	Attach the loader cap to the loader, re-flush the Flex Catheter and close the stopcock to the delivery system. Remove the stylet and flush the guidewire lumen of the delivery system. CAUTION: Keep THV hydrated until ready for implantation to prevent damage to the leaflets which may impact valve functionality. WARNING: The physician must verify correct orientation of the THV prior to its implantation; the inflow (outer skirt end) of the THV should be oriented distally towards the tapered tip to prevent the risk of severe patient harm.

7.3 Native Valve Predilation and THV Delivery

Native valve predilation and THV delivery should be performed under local and/or general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

Administer heparin to maintain the ACT at ≥ 250 sec.

CAUTION: Contrast media usage should be monitored to reduce the risk of renal injury.

7.3.1 Baseline Parameters

Step	Procedure
1	Perform a supra-aortic angiogram with the projection of the native aortic valve perpendicular to the view.
2	Evaluate the distance of the left and right coronary ostia from the aortic annulus in relation to the THV frame height.
3	Introduce a pacemaker (PM) lead until its distal end is positioned in the right ventricle.
4	Set the stimulation parameters to obtain 1:1 capture, and test pacing.

7.3.2 Native Valve Predilation

Refer to Edwards Transfemoral Balloon Catheter Instructions for Use.

7.3.3 THV Delivery

Step	Procedure
1	Prepare the Edwards sheath per its instructions for use.
2	If necessary, predilate the vessel.
3	Introduce the sheath per its instructions for use.
4	Insert the loader assembly into the sheath until the loader stops.

Step	Procedure
5	Advance the delivery system until the THV exits the sheath. CAUTION: The THV should not be advanced through the sheath if the sheath tip is not past the aortic bifurcation to minimize the risk of damage to the iliac vessel(s). CAUTION: The THV should not remain in the sheath for over 5 minutes as leaflet damage may result and impact valve functionality.
6	In a straight section of the vasculature, initiate valve alignment by disengaging the Balloon Lock and pulling the balloon catheter straight back until part of the Warning Marker is visible. Do not pull past the Warning Marker. WARNING: To prevent possible damage to the balloon shaft, ensure that the proximal end of the balloon shaft is not subjected to bending. WARNING: If valve alignment is not performed in a straight section, there may be difficulties performing this step which may lead to delivery system damage and inability to inflate the balloon. Utilizing alternate fluoroscopic views may help with assessing curvature of the anatomy. If excessive tension is experienced during valve alignment, repositioning the delivery system to a different straight section of the vasculature and relieving compression (or tension) in the system will be necessary. Engage the Balloon Lock. Utilize the Fine Adjustment Wheel to position the THV between the Valve Alignment Markers. Note: Do not turn the Fine Adjustment Wheel if the Balloon Lock is not engaged. WARNING: Do not position the THV past the distal Valve Alignment Marker to minimize the risk of improper THV deployment or THV embolization. CAUTION: Maintain guidewire position in the left ventricle during valve alignment to prevent loss of guidewire position.
7	Utilize the Flex wheel to traverse the aortic arch and cross the native valve. Note: Verify the Edwards logo is facing up. Note: The delivery system articulates in a direction opposite from the flush port.
8	Disengage the Balloon Lock and retract the tip of the Flex Catheter to the center of the Triple Marker. Engage the Balloon Lock.
9	Position the THV with respect to the native valve.
10	As necessary, utilize the Flex wheel to adjust the co-axiality of the THV and the Fine Adjustment Wheel to adjust the position of the THV.

Step	Procedure
11	Before deployment, ensure that the THV is correctly positioned between the Valve Alignment Markers and the Flex Catheter tip is over the Triple Marker.
12	Begin THV deployment: a) Unlock the inflation device. b) Ensure hemodynamic stability is established and begin rapid pacing; once arterial blood pressure has decreased to 50 mmHg or below, balloon inflation can commence. c) Using slow controlled inflation, deploy the THV with the entire volume in the inflation device, hold for 3 seconds and confirm that the barrel of the inflation device is empty to ensure complete inflation of the balloon. d) Deflate the balloon. When the balloon catheter has been completely deflated turn off the pacemaker.

7.3.4 System Removal

Step	Procedure
1	Unflex the delivery system while traversing the aortic arch. Verify that the Flex Catheter tip is locked over the Triple Marker. Retract the loader to the proximal end of the delivery system. Remove the delivery system from the sheath. Note: For subclavian-axillary approach, keep delivery system inside sheath until ready to remove all devices as one unit. CAUTION: Completely unflex the delivery system prior to removal to minimize the risk of vascular injury.

7.4 Verification of Prosthetic Valve Position and Measurements

Measure and record hemodynamic parameters.

Step	Procedure
1	Perform a supra-aortic angiogram to evaluate device performance and coronary patency.
2	Measure and record the transvalvular pressure gradients.
3	Remove all devices when the ACT level is appropriate (e.g., reaches < 150 sec). Refer to the Edwards sheath instructions for use for device removal.
4	Close the access site.

8.0 How Supplied

The THV is supplied sterile and nonpyrogenic packaged in buffered glutaraldehyde, in a plastic jar to which a tamper evident seal has been applied. Each jar is shipped in a shelf box containing a temperature indicator to detect exposure of the THV to extreme temperature. The shelf box is enclosed in Styrofoam prior to shipping.

The delivery system and accessories are supplied sterilized by ethylene oxide.

8.1 Storage

The THV must be stored at 10 °C to 25 °C (50 °F to 77 °F). Each jar is shipped in an enclosure containing a temperature indicator to detect exposure of the THV to extreme temperature.

The delivery system and accessories should be stored in a cool, dry place.

9.0 MR Safety

MR Conditional

Non-clinical testing has demonstrated that the Edwards SAPIEN 3 Ultra transcatheter heart valve is MR Conditional. A patient with this device can be scanned safely, immediately after placement of this device under the following conditions:

- Static magnetic field of 1.5 tesla (T) or 3.0 tesla (T)
- Maximum spatial gradient field of 2500 Gauss/cm (25 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the transcatheter heart valve is expected to produce a maximum temperature rise of 3.0 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 14.5 mm from the implant for spin echo images and 30 mm for gradient echo images when scanned in a 3.0T MRI system. The artifact obscures the device lumen in gradient echo images.

The implant has not been evaluated in MR systems other than 1.5 T or 3.0 T.

10.0 Patient Information

A patient registration form is provided with each THV. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the THV. Return the original form to the Edwards Lifesciences address indicated on the form and provide the temporary identification card to the patient prior to discharge.

11.0 Clinical Studies

SAPIEN 3 THV FOR PATIENTS AT LOW SURGICAL RISK- STS/ACC TRANSCATHETER VALVE THERAPY REGISTRY (TVTR) ANALYSIS

A database extract was performed on May 2020, which yielded 9,479 patients treated with an Edwards SAPIEN 3 transcatheter heart valve placed in a native aortic valve (via subclavian/transaxillary or transfemoral access). The patients were treated between July 16, 2015 and April 29, 2020. The procedure was performed in 592 participating hospitals.

There were 120 subclavian/transaxillary (SC/TAx) and 9,238 transfemoral (TF) within the low risk patient population. The mortality rate for SC/TAx was 2.6% and 0.7% for TF. The all strokes rate was SC/TAx 0.8% and 1.3% for TF. Device embolization for SC/TAx was 0% and 0% for TF. Migration for SC/TAx was 0% and 0% for TF. Life threatening bleeding for SC/TAx was 0% and 0% for TF. Major bleeding for SC/TAx was 0% and 0.3% for TF. Major vascular complication for SC/TAx was 0% and 0.6% for TF.

12.0 Recovered THV and Device Disposal

The explanted THV should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit.

Used devices may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

These products are manufactured and sold under one or more of the following US patent(s): US Patent No. 7,530,253; 7,780,723; 7,895,876; 8,382,826; 8,591,575; 8,690,936; 8,790,387; 9,301,840; 9,301,841; and 9,393,110; and corresponding foreign patents.

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• Edwards Commander leveringssystem (figur 2)

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© Copyright 2021, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC