Edwards

Edwards SAPIEN 3 System
Edwards SAPIEN 3 Transcatheter Heart Valve
Edwards Certitude Delivery System
Transapical and Transaortic
Sistema Edwards SAPIEN 3
Válvula cardíaca transcáteter Edwards SAPIEN 3
Sistema de suministro Edwards Certitude
Transapical y transaórtico
Sistema Edwards SAPIEN 3
Válvula cardíaca transcáteter Edwards SAPIEN 3
Sistema de colocação Edwards Certitude
Transapical e transaórtico

DIRECTORY

English	1
Español	9
Português	17
References / Referencias / Referências	25
Figures / Figuras	26-27
Symbol Legend / Significado de los símbolos / Legenda de símbolos	30-31

English

Instructions for Use

Implantation of the transcatheter heart valve should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in standard catheterization techniques.

1.0 Device Description

Edwards SAPIEN 3 System

The Edwards SAPIEN 3 system consists of the Edwards SAPIEN 3 transcatheter heart valve and delivery systems.

• **Edwards SAPIEN 3 Transcatheter Heart Valve (Figure 1)**

The Edwards SAPIEN 3 transcatheter heart valve (THV) is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) inner and outer fabric skirts. The leaflets are treated according to the Carpentier-Edwards TheraFix process.

The THV is recommended to be implanted in a native annulus size range associated with the three-dimensional area of the aortic annulus measured at the basal ring during systole:

Table 1

Native Valve Annulus Size (TEE)*	Native Valve Annulus Size (CT)		THV Size
	Area	Area Derived Diameter	
16 – 19 mm	273 – 345 mm ²	18.6 – 21.0 mm	20 mm
18 – 22 mm	338 – 430 mm ²	20.7 – 23.4 mm	23 mm
21 – 25 mm	430 – 546 mm ²	23.4 – 26.4 mm	26 mm
24 – 28 mm	540 – 683 mm ²	26.2 – 29.5 mm	29 mm

THV size recommendations are based on native valve annulus size, as measured by transesophageal echocardiography (TEE) or computed tomography (CT). Patient anatomical factors and multiple imaging modalities should be considered during THV size selection.

NOTE: Risks associated with undersizing and oversizing should be considered to minimize the risk of paravalvular leak, migration, and/or annular rupture.

* Due to limitations in two-dimensional imaging, 2-D TEE imaging should be supplemented with 3-D area measures.

Sizing recommendations for implanting the Edwards SAPIEN 3 transcatheter heart valve in a failing bioprosthesis are provided in the table below:

Surgical Valve True Inner Diameter (ID) ^[1]	THV-in-THV (Native Valve Annulus Size)	THV Size
16.5-19.0 mm	18.6-21.0 mm	20 mm
18.5-22.0 mm	20.7-23.4 mm	23 mm
22.0-25.0 mm	23.4-26.4 mm	26 mm
25.0-28.5 mm	26.2-29.5 mm	29 mm

NOTE: Surgical valve 'True ID' may be smaller than the labeled valve size. For THV-in-THV, the native valve annulus size should be considered to determine the appropriate THV size to implant. For a failing stentless bioprosthesis, consider sizing recommendations for a native annulus. The dimensions of the failed bioprosthesis should be determined so that the appropriate THV size can be implanted; and is best determined by using computed tomography, magnetic resonance imaging, and/or transesophageal echocardiography.

NOTE: Exact volume required to deploy the THV may vary depending on the bioprosthesis inner diameter. Factors such as calcification and pannus tissue growth may not be accurately visualized in imaging and may reduce the effective inner diameter of the failing bioprosthesis to a size smaller than the 'True ID'. These factors should be considered and assessed in order to determine the most appropriate THV size to achieve nominal THV deployment and sufficient anchoring. Do not exceed the rated burst pressure. See Table 2 for inflation parameters.

• Edwards Certitude Delivery System (Figure 2)

The Edwards Certitude delivery system facilitates the placement of the bioprosthesis. The delivery system consists of a flex catheter to aid in tracking and valve positioning. The delivery system includes a tapered tip to facilitate crossing of the valve. The handle contains a Flex Wheel to control flexing of the Balloon Catheter. A stylet is included within the guidewire lumen of the delivery system. A radiopaque Center Marker in the balloon is provided to help with valve positioning. The extension tubing is used during THV deployment. The inflation parameters for the valve deployment are:

Table 2

Model	Nominal Balloon Diameter	Nominal Inflation Volume	Rated Burst Pressure (RBP)
9620TA20	20 mm	12 mL	7 atm
9620TA23	23 mm	17 mL	7 atm
9620TA26	26 mm	23 mL	7 atm
9620TA29	29 mm	30 mL	7 atm

• Edwards Certitude Introducer Sheath Set (Figure 3)

The Edwards Certitude introducer sheath set facilitates the introduction and removal of devices utilized with the SAPIEN 3 transcatheter heart valve. The sheath has a radiopaque marker for visualization of the sheath tip and non-radiopaque depth markings on the distal end of the body of the sheath. The proximal end of the sheath includes a flush tube and three hemostasis valves. An introducer is supplied with the sheath. The entire introducer is radiopaque.

Introducer Sheath Set Information

Model	9620IS18	9620IS21
Sheath Inside Diameter	18F (6.1 mm)	21F (6.9 mm)
Sheath Effective Length	21 cm	21 cm
Introducer Size	OD: 6.3 mm	OD: 7.0 mm
Introducer Effective Length	33 cm	
Diameter of the largest guidewire that can be used	0.035" (0.89 mm)	

• Qualcrimp Crimping Accessory (Figure 4a)

The Qualcrimp crimping accessory (Figure 4a) is used during THV crimping.

• Loader (Figure 4b)

The loader is used to aid insertion of the delivery system into the sheath.

• Crimper and Crimp Stopper (Figure 4c)

The crimper reduces the diameter of the valve to mount it onto the delivery system. The crimper is comprised of a housing and a compression mechanism that is closed with a handle located on the housing. A 2-piece crimp stopper is used to crimp the valve to its intended diameter.

• Inflation Device

An inflation device with locking mechanism is used during valve deployment.

NOTE: For proper volume sizing, the delivery system and the Ascendra balloon aortic valvuloplasty catheter must be used with the inflation device provided by Edwards Lifesciences.

2.0 Indications

1. The Edwards SAPIEN 3 system is indicated for use in patients with heart disease due to native calcific aortic stenosis at any or all levels of surgical risk for open heart surgery.
2. The Edwards SAPIEN 3 system is indicated for use in patients with symptomatic heart disease due to a failing aortic bioprosthetic valve or a failing mitral surgical bioprosthetic valve (stenosed, insufficient, or combined) who are judged by a heart team to be at high or greater risk for open surgical therapy (i.e., predicted risk of surgical mortality $\geq 8\%$ at 30 days, based on the Society of Thoracic Surgeons (STS) risk score and other clinical co-morbidities unmeasured by the STS risk calculator).

3.0 Contraindications

Use of the Edwards SAPIEN 3 system is contraindicated in patients with:

- Evidence of intracardiac mass, thrombus, vegetation, active infection, or endocarditis;
- Inability to tolerate anticoagulation/antiplatelet therapy.

4.0 Warnings

- The devices are designed, intended, and distributed **STERILE** for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- Correct sizing of the THV is essential to minimize the risk of paravalvular leak, migration, and/or annular rupture.
- The physician must verify correct orientations of the THV prior to its implantation.
- Accelerated deterioration of the THV may occur in patients with altered calcium metabolism.
- When using venous pacing, observation of the pacing lead throughout the procedure is essential to avoid the potential risk of the pacing lead causing a cardiovascular perforation.
- The THV must remain hydrated at all times and cannot be exposed to any solutions, chemicals, antibiotics, etc. other than its shipping storage solution and sterile physiologic saline solution to prevent leaflet damage that may impact valve functionality. THV leaflets mishandled or damaged during any part of the procedure will require replacement of the THV.
- Patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, and/or polymeric materials may have an allergic reaction to these materials.
- Do not use the THV if the tamper evident seal is broken, as sterility may be compromised.
- Do not use the THV if the temperature indicator has been activated, as valve function may be compromised.
- Do not use the THV if the expiration date has elapsed, as either sterility or valve function may be compromised.

- Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging sterile barriers and any components have been opened or damaged, cannot be flushed, or the expiration date has elapsed.

- Balloon valvuloplasty should be avoided in the treatment of failing bioprostheses as this may result in embolization of bioprosthesis material and mechanical disruption of the valve leaflets.

5.0 Precautions

- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Material Safety Data Sheet available from Edwards Lifesciences.
- The safety of the THV implantation has not been established in patients who have:
 - Congenital unicuspid aortic valve
 - Pre-existing prosthetic ring in any position.
 - Severe ventricular dysfunction with ejection fraction $< 20\%$
 - Hypertrophic cardiomyopathy with or without obstruction
 - Aortic stenosis characterized by a combination of AV low flow, low gradient
- Residual mean gradient may be higher in a “THV-in-failing bioprosthesis” configuration than that observed following implantation of the valve inside a native aortic annulus using the same size device. Patients with elevated mean gradient post procedure should be carefully followed. It is important that the manufacturer, model and size of the preexisting bioprosthetic valve be determined, so that the appropriate valve can be implanted and a prosthesis-patient mismatch be avoided. Additionally, pre-procedure imaging modalities must be employed to make as accurate a determination of the inner diameter as possible.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Special care must be exercised in mitral valve replacement if chordal preservation techniques were used in the primary implantation to avoid entrapment of the subvalvular apparatus.
- THV recipients should be maintained on anticoagulant/antiplatelet therapy to minimize risk of valve thrombosis or thromboembolic events, as determined by their physicians.
- Long-term durability has not been established for the THV. Regular medical follow-up is advised to evaluate valve performance.
- Based on the treating physician's consideration of risks and benefits, the SAPIEN 3 valve may be implanted in relatively young patients, although the longer-term durability is still the subject of ongoing clinical research.
- Do not overinflate the deployment balloon, as this may prevent proper valve leaflet coaptation and thus impact valve functionality.
- Patients with pre-existing bioprostheses should be carefully evaluated before implantation of the THV to ensure proper THV positioning and deployment.

6.0 Potential Adverse Events

Potential risks associated with the overall procedure including access, cardiac catheterization, local and/or general anesthesia:

- Allergic reaction to antithrombotic therapy or contrast medium or anesthesia
- Anemia
- Aneurysm
- Angina
- Arrhythmias including ventricular fibrillation (VF) and ventricular tachycardia (VT)
- AV fistula or pseudoaneurysm
- Cardiogenic shock
- Compartment syndrome
- Death
- Dissection: aortic or other vessels
- Emboli, distal (air, tissue or thrombotic emboli)
- Hematoma
- Hypertension or hypotension
- Inflammation
- Myocardial ischemia or infarction
- Pain or changes at the access site
- Perforation or rupture of cardiac structures
- Perforation or rupture of vessels
- Pericardial effusion or cardiac tamponade
- Peripheral ischemia or nerve injury
- Pulmonary edema
- Renal insufficiency or renal failure
- Respiratory insufficiency or respiratory failure
- Syncope
- Vasovagal response
- Vessel spasm
- Vessel thrombosis/occlusion
- Vessel trauma requiring surgical repair or intervention

Additional potential risks associated with the transcatheter heart valve implantation procedure, the bioprosthesis, and the use of its associated devices and accessories include:

- Allergic/immunologic reaction to the implant
- Atrial fibrillation/Atrial flutter
- Bleeding requiring transfusion or intervention
- Cardiac arrest
- Cardiac failure or low cardiac output

- Cardiogenic shock
- Conduction system injury (defect) including AV block, which may require a permanent pacemaker
- Coronary occlusion
- Dissection, rupture, trauma of the aortic annulus and surrounding structures including ascending aorta, coronary ostia and ventricular septum
- Emergency cardiac surgery
- Hemolysis
- Infection, fever, septicemia, abscess, endocarditis
- Injury to mitral valve
- Left ventricular outflow tract obstruction
- Mechanical failure of delivery system, and/or accessories, including balloon rupture and tip separation
- Mediastinitis
- Mediastinal bleeding
- Silent cerebral ischemia, stroke, transient ischemic attack, cognitive impairment
- Structural valve deterioration (wear, fracture, calcification, stenosis)
- Valve deployment in unintended location
- Valve explants
- Valve migration, malposition or embolization requiring intervention
- Valve regurgitation, paravalvular or transvalvular
- Valve thrombosis

7.0 Directions for Use

7.1 System Compatibility

Product Name	20 mm	23 mm	26 mm	29 mm
	Model/REF			
Edwards SAPIEN 3 Transcatheter Heart Valve	9600TFX (20 mm)	9600TFX (23 mm)	9600TFX (26 mm)	9600TFX (29 mm)
Edwards Certitude Delivery System	9620TA20	9620TA23	9620TA26	9620TA29
Edwards Certitude Introducer Sheath Set	9620IS18 (18F)			9620IS21 (21F)
Crimper	9600CR			
Inflation device, Qualcrimp crimping accessory, 2-piece crimp stopper, loader and extension tubing provided by Edwards Lifesciences				

Additional Equipment

- Cardiac catheterization/hybrid OR suite
- Standard cardiac catheterization lab equipment and supplies, and access to standard heart valve operating room equipment and supplies

- Fluoroscopy (fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal or transthoracic echocardiography system
- 18 gauge Seldinger needle (for transaortic)
- 145 cm x 0.035" (0.89 mm) soft guidewire
- 180 cm or 260 cm x 0.035" (0.89 mm) & Exchange length 0.035" (0.89 mm) extra-stiff guidewires
- Temporary pacemaker (PM) and pacing leads
- 20 mm Ascendra balloon aortic valvuloplasty catheter (BAVC) or equivalent for 23 mm, 26 mm, and 29 mm valves
- 16 mm commercially available valvuloplasty balloon catheter for 20 mm valve
- Sterile rinsing basins; physiological saline; heparinized saline; radiopaque 15% diluted radiopaque contrast medium
- Sterile table for THV and accessories preparation
- 20 cc syringe or larger
- 50 cc syringe or larger
- High-pressure 3-way stopcock

7.2 Valve Handling and Preparation

Follow sterile technique during device preparation and implantation.

7.2.1 THV Rinsing Procedure

Before opening, carefully examine the jar for evidence of damage (e.g., a cracked jar or lid, leakage, or broken or missing seals).

CAUTION: If the container is found to be damaged, leaking, without adequate sterilant, or missing intact seals, the THV must not be used for implantation, as sterility may be compromised.



Step	Procedure
1	Set up two (2) sterile bowls with at least 500 mL of sterile physiological saline to thoroughly rinse the THV.
2	Carefully remove the valve/holder assembly from the jar without touching the tissue. Verify the valve serial identification number with the number on the jar lid and record in the patient information documents. Inspect the valve for any signs of damage to the frame or tissue.

Step	Procedure
3	Rinse the THV as follows: • Place the THV in the first bowl of sterile, physiological saline. Be sure the saline solution completely covers the THV and holder. • With the valve and holder submerged, slowly agitate (to gently swirl the valve and holder) back and forth for a minimum of 1 minute. • Transfer the THV and holder to the second rinsing bowl of sterile physiological saline and gently agitate for at least one more minute. Ensure the rinse solution in the first bowl is not used. • The valve should be left in the final rinse solution until needed to prevent the tissue from drying. CAUTION: Do not allow the valve to come into contact with the bottom or sides of the rinse bowl during agitation or swirling in the rinse solution. Direct contact between the identification tag and valve is also to be avoided during the rinse procedure. No other objects should be placed in the rinse bowls. The valve should be kept hydrated to prevent the tissue from drying.

7.2.2 Prepare the System

Step	Procedure
1	Visually inspect all components for damage. Ensure the system is fully unflexed.
2	Prime and flush the introducer and sheath with heparinized saline. Hydrate the length of the introducer and sheath.
3	Advance the introducer fully into the sheath housing.
4	Unscrew the loader cap from the loader and flush the loader cap with heparinized saline.
5	Place the loader cap onto the delivery system with the inside of the cap oriented towards the tapered tip.
6	Flush the extension tubing and connect to the delivery system.
7	Partially fill a 50 mL or larger syringe with diluted contrast medium, and connect to the extension tubing.
8	Fill the inflation device with 20 mL of diluted contrast medium, lock the inflation device, and connect to the extension tubing. Close 3-way stopcock to inflation device.
9	De-air the delivery system using the luer lock syringe. Leave zero-pressure in the system. Close the 3-way stopcock to the luer lock syringe.
10	Remove 3 mL fluid from the delivery system by turning the knob of the locked inflation device. Keep the inflation device locked for THV crimping steps.

7.2.3 Mount and Crimp the THV onto the Delivery System

Step	Procedure
1	Completely submerge the Qualcrimp crimping accessory in a bowl of 100 mL physiological saline solution. Gently compress until fully saturated. Swirl for a minimum of 1 minute. Repeat this process in a second bowl.
2	Rotate the crimper until the aperture is fully opened. Attach the 2-piece Crimp Stopper to the crimper.
3	Remove the THV from the holder and remove the ID tag.
4	If necessary, partially crimp the THV in the crimper until it fits snugly inside the Qualcrimp crimping accessory. NOTE: Partial crimping is not necessary for the 20 mm valve.
5	Place the Qualcrimp crimping accessory over the THV.
6	The orientation of the THV on the delivery system is described below: Antegrade Approach: Inflow (outer skirt) end of the THV towards the proximal end of the delivery system.  Retrograde Approach: Inflow (outer skirt) end of the THV towards the distal end of the delivery system. 
7	Place the THV and Qualcrimp crimping accessory in the crimper. Insert the delivery system coaxially into the THV.
8	Crimp the THV between the two internal shoulders of the delivery system until it reaches the Qualcrimp stop.
9	Remove the Qualcrimp crimping accessory from the THV/balloon assembly and Qualcrimp stop from the Crimp Stopper, leaving the Final Stop in place. NOTE: Ensure that the THV remains centered and coaxial within the two internal shoulders.
10	Place the THV/balloon assembly back in the crimper aperture, fully crimp the THV until it reaches the Final Stop and hold for 5 seconds.

Step	Procedure
11	Repeat the full crimp of the THV two times for a total of 3 crimps.
12	Flush the loader with heparinized saline. Immediately advance the loader over the THV until the tapered tip of the delivery system is exposed and the THV is within the distal end of the loader tube. CAUTION: The THV should not remain fully crimped and/or in the loader over 15 minutes, as leaflet damage may result and impact valve functionality.
13	Attach the loader cap to the loader and flush through the flush port on the loader. Remove the stylet and flush the guidewire lumen of the delivery system. CAUTION: Keep the THV hydrated until ready for implantation to prevent damage to the leaflets which may impact valve functionality. WARNING: The physician must verify correct orientation of the THV prior to its implantation to prevent the risk of severe patient harm.
14	With 3-way stopcock still closed to the luer lock syringe, unlock the inflation device. Allow the delivery system to reach zero-pressure.
15	Close the 3-way stopcock to the delivery system. Use the luer lock syringe to de-air the inflation device if necessary.
16	Adjust the inflation device to the inflation volume required to deploy the THV, per table 2. Re-lock the inflation device. Close the 3-way stopcock to the luer lock syringe and remove syringe. CAUTION: Maintain the inflation device in a locked position until THV deployment to prevent premature balloon inflation and subsequent improper THV deployment.

7.3 Native Valve Predilation and THV Delivery

Native valve predilation and THV delivery should be performed under general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

The following table shows the minimum required distances from the valvular plane to the distal tip of the Edwards Certitude sheath to allow the Edwards Certitude delivery system balloon to inflate properly during THV deployment. **These distances do not include sheath insertion depth**, which should be considered during the transaortic approach when selecting the access site on the ascending aorta.

Delivery System	THV	Minimum Required Distance From Sheath Tip to Valvular Plane
Model 9620TA20	20 mm	3.5 cm
Model 9620TA23	23 mm	3.5 cm
Model 9620TA26	26 mm	3.5 cm
Model 9620TA29	29 mm	4.0 cm

Administer heparin to maintain the ACT at ≥ 250 sec.

Balloon valvuloplasty should be avoided in the treatment of failing bioprostheses as this may result in embolization of bioprosthesis material and mechanical disruption of the valve leaflets.

CAUTION: Contrast media use should be monitored to reduce the risk of renal injury.

7.3.1 Baseline Parameters

Step	Procedure
1	Advance a 5F (1.67 mm) or 6F (2.0 mm) pigtail catheter and perform an angiogram with the projection of the valve perpendicular to the view.
2	For aortic implantation, evaluate the distances of the right and left coronary ostia from the aortic annulus in relation to the THV frame height.
3	Introduce a pacemaker (PM) lead and position appropriately.
4	Set the stimulation parameters to obtain 1:1 capture, and test pacing.

7.3.2 Access

CAUTION: Care should be taken to avoid damage to soft tissue, chordae, aorta, native leaflet or ventricular wall during insertion, positioning and removal of devices.

Transapical Access	
Step	Procedure
1	Access the apex through an anterior mini thoracotomy at the 5th or 6th intercostal space. Incise the pericardium to expose the apex of the left ventricle (LV).
2	Attach epicardial pacing leads to left ventricle or insert transvenous pacing leads and secure proximal ends of leads into pacemaker. Set the stimulation parameters, test rapid pacing.
3	Place a reinforced double purse string on the LV apex to access the left ventricle.
4	Gain access through standard transapical techniques.
5	Insert the tip of the Edwards Certitude introducer sheath set or desired introducer sheath for BAV through the apex of the LV to approximately 4 cm and locate the sheath tip in the LV immediately below the valve; withdraw the introducer slowly, keeping the sheath in place. Maintain guidewire position across the valve.

Transaortic Access	
Step	Procedure
1	Access the ascending aorta using standard surgical technique (e.g., a partial J-sternotomy or right parasternal mini thoracotomy).
2	Place two reinforced purse string sutures at the intended access site in the ascending aorta. NOTE: The selected access site should be soft by digital palpation.

Transaortic Access	
Step	Procedure
3	Introduce a pacemaker lead until its distal end is positioned in the right ventricle. Set the stimulation parameters and test pacing.
4	Gain aortic valve access through standard transaortic techniques.
5	Insert the Edwards Certitude introducer sheath set, or desired introducer sheath for BAV, into the aorta to approximately 2 cm. Withdraw the introducer slowly, keeping the sheath in place. Maintain guidewire position across the aortic valve.

7.3.3 Native Valve Predilation

CAUTION: Care should be taken to avoid damage to soft tissue, chordae, aorta, native leaflet or ventricular wall during insertion, positioning and removal of devices.

Step	Procedure
1	Prepare the valvuloplasty balloon catheter per its instructions for use.
2	Advance the prepared valvuloplasty balloon catheter through the sheath over the guidewire, cross the aortic valve, and position the balloon.
3	Begin predilation: - Begin rapid pacing. Once arterial blood pressure has decreased to 50 mmHg or below, balloon inflation can commence. - Inflate the valvuloplasty balloon catheter as per its instructions for use. - Completely deflate the balloon. Stop rapid pacing.
4	Remove the valvuloplasty balloon catheter, leaving the guidewire in place. NOTE: If not using the Edwards Certitude sheath for native valve predilation, remove the sheath used for the valvuloplasty and advance the Edwards Certitude introducer sheath set over the guidewire.

7.3.4 THV Delivery

CAUTION: Care should be taken to avoid damage to soft tissue, chordae, aorta, native leaflet or ventricular wall during insertion, positioning and removal of devices.

Step	Procedure
1	Confirm that the THV is oriented properly and the volume in the inflation device matches the indicated volume.
2	Advance the THV/balloon assembly with the loader over the guidewire.
3	Engage loader into the sheath housing while maintaining a firm grip.

Step	Procedure
4	Advance the valve out of the loader into the large section of the sheath. Tap on the sheath housing to release air bubbles to the proximal end of the loader. Depress button valve on loader to aspirate the loader.
5	Advance the THV/balloon assembly through the sheath and position within the target valve. If needed, rotate the flex wheel on the handle to articulate the THV/balloon assembly into position. CAUTION: To prevent possible leaflet damage that may impact valve functionality, the THV should not remain in the sheath for over 5 minutes.
6	Ensure that the THV is correctly positioned between the two internal shoulders of the delivery system.
7	Begin THV deployment: – Unlock the inflation device. – Ensure hemodynamic stability is established and begin rapid pacing. Once arterial blood pressure has decreased to 50 mmHg or below, balloon inflation can commence. – Using a slow, controlled inflation, deploy the THV by inflating the balloon with the entire volume in the inflation device, hold for 3 seconds and confirm that the barrel of the inflation device is empty to ensure complete inflation of the balloon. – Once the THV has been deployed, rapidly deflate the balloon catheter. When the delivery system balloon has been completely deflated, turn off the pacemaker.
8	If articulation was used, return the delivery system to the straight position prior to removal. Retract the delivery system and guidewire into the sheath. Remove the loader and delivery system from the sheath. CAUTION: Properly deflate the balloon and straighten the delivery system prior to removal to prevent patient injury.

7.4 Verification of THV Position and Measurements

Step	Procedure
1	Perform an angiogram to evaluate device performance and coronary patency (if applicable).
2	Measure and record the transvalvular pressure gradients and assess valve competency.
3	Upon satisfactory deployment, remove all devices when the ACT level is appropriate (e.g., reaches < 150 sec).
4	Tie the purse string sutures in place and confirm hemostasis.

8.0 How Supplied

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system and accessories are supplied sterilized by ethylene oxide gas.

The THV is supplied sterile and non-pyrogenic, packaged in buffered glutaraldehyde, in a plastic jar to which a tamper evident seal has been applied. Each jar is shipped in a shelf box containing a temperature indicator to detect exposure of the THV to extreme temperature. The shelf box is enclosed in Styrofoam prior to shipping.

8.1 Storage

The THV must be stored at 10 °C to 25 °C (50 °F to 77 °F). Each jar is shipped in an enclosure containing a temperature indicator to detect exposure of the THV to extreme temperature.

The delivery system and accessories should be stored in a cool, dry place.

9.0 MR Safety

MR Conditional

Non-clinical testing has demonstrated that the SAPIEN 3 transcatheter heart valve is MR Conditional. A patient with this device can be scanned safely, immediately after placement of this device under the following conditions:

- Static magnetic field of 1.5 tesla (T) or 3 tesla
- Maximum spatial gradient field of 2500 Gauss/cm (25 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the transcatheter heart valve is expected to produce a maximum temperature rise of 3.0 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 14.5 mm from the implant for spin echo images and 30 mm for gradient echo images when scanned a 3.0 T MRI system. The artifact obscures the device lumen in gradient echo images.

The implant has not been evaluated in MR systems other than 1.5 T or 3.0 T.

For valve-in-valve implantation or in the presence of other implants, please refer to the MRI safety information for the surgical valve or other devices prior to MR imaging.

10.0 Patient Information

A patient registration form is included with each THV. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the THV. Return the original form to the Edwards Lifesciences address indicated on the form and provide the temporary identification card to the patient prior to discharge.

11.0 Recovered THV and Device Disposal

The explanted THV should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit.

Used delivery system may be disposed of in the same manner that the hospital waste and biohazardous materials are handled. There are no special risks related to the disposal of these devices.

This product is manufactured and sold under one or more of the following US patents: 7,530,253; 7,895,876; 8,439,970; 8,475,522; 8,764,820; 8,945,208; and 9,393,110; and corresponding foreign patents.

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
1. Bapat V, Attia R, Thomas M. Effect of Valve Design on the Stent Internal Diameter of a Bioprosthetic Valve: A Concept of True Internal Diameter and Its Implications for the Valve-in-Valve Procedure. *JACC: Cardiovascular Interventions*. Vol. 7, No. 2 2014: 115–127.

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This page intentionally left blank.

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INTL_ESP_140x182_4

Edwards

Edwards Lifesciences Services GmbH

Edisonstr. 6
85716 Unterschleißheim
Germany

2021-08
10037390003 A
© Copyright 2021, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525