Edwards

Edwards SAPIEN 3 Ultra Transcatheter Heart Valve System

Edwards SAPIEN 3 and Edwards SAPIEN 3 Ultra Transcatheter Heart Valve Edwards Certitude Delivery System

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English

Instructions for Use

Implantation of the transcatheter heart valve should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in standard catheterization techniques.

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system, Edwards sheath, and crimper are supplied sterilized with ethylene oxide gas.

1.0 Device Description

Edwards SAPIEN 3 Ultra Transcatheter Heart Valve (THV) System

The Edwards SAPIEN 3 Ultra transcatheter heart valve system consists of the Edwards SAPIEN 3 Ultra and SAPIEN 3 transcatheter heart valves and delivery systems.


• Edwards SAPIEN 3 Ultra Transcatheter Heart Valve - (Figure 1)

The Edwards SAPIEN 3 Ultra transcatheter heart valve is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) inner and outer fabric skirts. The leaflets are treated according to the Carpentier-Edwards ThermaFix process.

• Edwards SAPIEN 3 Transcatheter Heart Valve - (Figure 2)

The Edwards SAPIEN 3 transcatheter heart valve is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) fabric skirt. The leaflets are treated according to the Carpentier-Edwards ThermaFix process.

Edwards, Edwards Lifesciences, the stylized E logo, Carpentier-Edwards, Certitude, Edwards SAPIEN, Edwards SAPIEN 3, Edwards SAPIEN 3 Ultra, eSheath, INSPIRIS, INSPIRIS RESILIA, Qualcrimp, RESILIA, SAPIEN, SAPIEN 3, SAPIEN 3 Ultra, ThermaFix, and VFit are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.


9750TFX

Figure 1. Edwards SAPIEN 3 Ultra Transcatheter Heart Valve

Table 1

Valve Size	Valve Height
20 mm	15.5 mm
23 mm	18 mm
26 mm	20 mm

9600TFX

Figure 2. Edwards SAPIEN 3 Transcatheter Heart Valve

Table 2

Valve Size	Valve Height
20 mm	15.5 mm
23 mm	18 mm
26 mm	20 mm
29 mm	22.5 mm

Sizing recommendations for implanting the Edwards SAPIEN 3 Ultra transcatheter heart valve and the Edwards SAPIEN 3 transcatheter heart valve in a native annulus are provided in the table below:

Table 3

Native Valve Annulus Size (TEE)	Native Valve Annulus Size (CT)		THV Size
	Area	Area Derived Diameter	
16 – 19 mm	273 – 345 mm ²	18.6 – 21.0 mm	20 mm
18 – 22 mm	338 – 430 mm ²	20.7 – 23.4 mm	23 mm
21 – 25 mm	430 – 546 mm ²	23.4 – 26.4 mm	26 mm
24 – 28 mm	540 – 683 mm ²	26.2 – 29.5 mm	29 mm

Valve size recommendations are based on native valve annulus size, as measured by transesophageal echocardiography (TEE) or computed tomography (CT). Patient anatomical factors and multiple imaging modalities should be considered during valve size selection.

NOTE: Risks associated with undersizing and oversizing should be considered.

Sizing recommendations for implanting the Edwards SAPIEN 3 Ultra transcatheter heart valve and the Edwards SAPIEN 3 transcatheter heart valve in a failing surgical and transcatheter bioprosthesis, except for the INSPIRIS RESILIA aortic valve sizes 19-25 mm, are provided in the table below:

Table 4

Surgical Valve True Inner Diameter (ID) ^[1]	THV-in-THV (Native Valve Annulus Size)	SAPIEN 3 Ultra / SAPIEN 3 THV Size
16.5 – 19.0 mm	18.6-21.0 mm	20 mm
18.5 – 22.0 mm	20.7-23.4 mm	23 mm
22.0 – 25.0 mm	23.4-26.4 mm	26 mm
25.0 – 28.5 mm	26.2-29.5 mm	29 mm

NOTE: Surgical valve 'True ID' may be smaller than the labeled valve size. For a failing stentless bioprosthesis, consider sizing recommendations for a native annulus. The dimensions of the failed bioprosthesis should be determined so that the appropriate THV size can be implanted; and is best determined by using computed tomography, magnetic resonance imaging, and/or transesophageal echocardiography.

NOTE: Exact volume required to deploy the THV may vary depending on the bioprosthesis inner diameter. Factors such as calcification and pannus tissue growth may not be accurately visualized in imaging and may reduce the effective inner diameter of the failing bioprosthesis to a size smaller than the 'True ID'. These factors should be considered and assessed in order to determine the most appropriate THV size to achieve nominal THV deployment and sufficient anchoring. Do not exceed the rated burst pressure. See Table 6 for inflation parameters.

Sizing recommendations for implanting the Edwards SAPIEN 3 Ultra transcatheter heart valve and the SAPIEN 3 transcatheter heart valve in a failing INSPIRIS RESILIA aortic surgical bioprosthesis in sizes 19 – 25 mm are provided in the table below:

Table 5

INSPIRIS RESILIA Aortic Valve (model 11500A)* Labeled Size	THV Size
19 mm	20 mm or 23 mm
21 mm	23 mm or 26 mm
23 mm	26 mm
25 mm	29 mm

*INSPIRIS RESILIA aortic valve model 11500A sizes 19 – 25 mm incorporate VFit technology which consists of expandable bands and fluoroscopically visible size markers designed for potential future valve-in-valve procedures. Clinical data are not currently available on the INSPIRIS RESILIA aortic valve Model 11500A valve-in-valve procedure or expansion feature. The impact of tissue ingrowth on the expansion feature of the INSPIRIS RESILIA aortic valve has not been assessed.

WARNING: Do not perform stand-alone balloon aortic valvuloplasty procedures in the INSPIRIS RESILIA aortic valve for the sizes 19-25 mm. This may expand the valve causing aortic incompetence, coronary embolism or annular rupture.

NOTE: INSPIRIS RESILIA aortic valve model 11500A sizes 27 – 29 mm do not incorporate VFit technology and therefore follow the surgical valve True ID sizing provided in Table 4.

• **Edwards Certitude Delivery System (Figure 3)**

The Edwards Certitude delivery system facilitates the placement of the bioprosthesis. The delivery system consists of a Flex Catheter to aid in tracking and valve positioning. The delivery system includes a tapered tip to facilitate crossing of the valve. The handle contains a Flex Wheel to control flexing of the balloon catheter. A stylet is included within the guidewire lumen of the delivery system. A radiopaque Center Marker in the balloon is provided to assist with valve positioning. The inflation parameters for valve deployment are:

Table 6

Model	Nominal Balloon Diameter	Nominal Inflation Volume	Rated Burst Pressure (RBP)
9620TA20	20 mm	12 mL	7 atm
9620TA23 9630TA23	23 mm	17 mL	7 atm
9620TA26 9630TA26	26 mm	23 mL	7 atm
9620TA29	29 mm	30 mL	7 atm


Figure 3. Edwards Certitude Delivery System

• **Edwards Certitude Introducer Sheath Set (Figure 4)**

The Edwards Certitude introducer sheath set facilitates the introduction and removal of devices utilized with the SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valves. The sheath has a radiopaque marker for visualization of the sheath tip and non-radiopaque depth markings on the distal end of the body of the sheath. The proximal end of the sheath includes a flush tube and three hemostasis valves. An introducer is supplied with the sheath. The entire introducer is radiopaque.


Figure 4. Edwards Certitude Introducer Sheath Set

Additional Accessories


Figure 6. Qualcrimp Crimping Accessory


Figure 5. Loader


Figure 7. Extension Tubing


Figure 8. Edwards Crimper and 2-piece Crimp Stopper

• Loader (Figure 5)

The loader allows for the delivery of the crimped valve through the hemostasis valves of the sheath

• Qualcrimp Crimping Accessory (Figure 6)

The Qualcrimp crimping accessory (packaged with the Edwards Certitude System) is used during THV crimping of the bioprosthesis.

• Extension Tubing (Figure 7)

The extension tubing is used during clinical use of the delivery system.

• Edwards Crimper (Figure 8)

The Edwards crimper reduces the diameter of the valve to mount it onto the delivery system. The crimper is comprised of a housing and a compression mechanism that is closed with a handle located on the housing. A 2-piece crimp stopper is used to crimp the valve to its intended diameter.

2.0 Indications

The Edwards SAPIEN 3 Ultra Transcatheter Heart Valve System is indicated for relief of aortic stenosis in patients with symptomatic heart disease due to severe native calcific aortic stenosis or in patients with symptomatic heart disease due to a failing aortic bioprosthetic valve or a failing mitral surgical bioprosthetic valve (stenosed, insufficient, or combined). Patients must be deemed by a Heart Team to be 'non-operable', defined as the probability of death or serious, irreversible morbidity exceeding 50%.

3.0 Contraindications

The valve and delivery systems are contraindicated in patients who cannot tolerate an anticoagulation/antiplatelet regimen or who have active bacterial endocarditis or other active infections.

4.0 Warnings

- Observation of the pacing lead throughout the procedure is essential to avoid the potential risk of pacing lead perforation.
- There may be an increased risk of stroke in transcatheter aortic valve replacement procedures, as compared to balloon aortic valvuloplasty or other standard treatments in high or greater risk patients.
- The devices are designed, intended, and distributed for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.

- Incorrect sizing of the valve may lead to paravalvular leak, migration, embolization, residual gradient (patient-prosthesis mismatch) and/or annular rupture.
- The physician must verify correct orientation of the THV prior to its implantation to prevent the risk of severe patient harm.
- Accelerated deterioration of the valve may occur in patients with an altered calcium metabolism.
- Prior to delivery, the valve must remain hydrated at all times and cannot be exposed to solutions other than its shipping storage solution and sterile physiologic rinsing solution. Valve leaflets mishandled or damaged during any part of the procedure will require replacement of the valve.
- Caution should be exercised in implanting a valve in patients with clinically significant coronary artery disease.
- Patients with pre-existing bioprostheses should be carefully assessed prior to implantation of the valve to ensure proper valve positioning and deployment.
- Do not use the valve if the tamper evident seal is broken, the storage solution does not completely cover the valve, the temperature indicator has been activated, the valve is damaged, or the expiration date has elapsed.
- Do not mishandle the delivery system or use it if the packaging or any components are not sterile, have been opened or are damaged (e.g. kinked or stretched), cannot be flushed, or the expiration date has elapsed.
- Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.
- Patient injury could occur if the delivery system is not un-flexed prior to removal.
- Care should be exercised in patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, and/or polymeric materials.
- The procedure should be conducted under fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.
- Valve recipients should be maintained on anticoagulant/antiplatelet therapy, except when contraindicated, as determined by their physician. This device has not been tested for use without anticoagulation.
- Do not add or apply antibiotics to the storage solution, rinse solutions, or to the valve.
- Balloon valvuloplasty should be avoided in the treatment of failing bioprostheses as this may result in embolization of bioprosthesis material and mechanical disruption of the valve leaflets.

5.0 Precautions

- Long-term durability has not been established for the valve. Regular medical follow-up is advised to evaluate valve performance.
- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Material Safety Data Sheet available from Edwards Lifesciences.
- To maintain proper valve leaflet coaptation, do not overinflate the deployment balloon.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Special care must be exercised in mitral valve replacement if chordal preservation techniques were used in the primary implantation to avoid entrapment of the subvalvular apparatus.
- Safety and effectiveness have not been established for patients with the following characteristics/comorbidities:
 - Non-calcified aortic annulus
 - Severe ventricular dysfunction with ejection fraction < 20%
 - Congenital unicuspid aortic valve
 - Pre-existing prosthetic ring in any position
 - Severe mitral annular calcification (MAC), severe (> 3+) mitral insufficiency, or Gorlin syndrome
 - Blood dyscrasias defined as: leukopenia (WBC < 3000 cells/mL), acute anemia (Hb < 9 g/dL), thrombocytopenia (platelet count < 50,000 cells/mL), or history of bleeding diathesis or coagulopathy
 - Hypertrophic cardiomyopathy with or without obstruction (HOCM)
 - Echocardiographic evidence of intracardiac mass, thrombus, or vegetation
 - A known hypersensitivity or contraindication to aspirin, heparin, ticlopidine (Ticlid™), or clopidogrel (Plavix™), or sensitivity to contrast media, which cannot be adequately premedicated
 - Bulky calcified aortic valve leaflets in close proximity to coronary ostia

-
- A concomitant paravalvular leak where the failing bioprosthesis is not securely fixed in the native annulus or is not structurally intact (e.g. wireform frame fracture)
 - A partially detached leaflet of the failing bioprosthesis that in the aortic position may obstruct a coronary ostium
 - Residual mean gradient may be higher in a “THV-in-failing bioprosthesis” configuration than that observed following implantation of the valve inside a native aortic annulus using the same size device. Patients with elevated mean gradient post procedure should be carefully followed. It is important that the manufacturer, model and size of the preexisting bioprosthetic valve be determined, so that the appropriate valve can be implanted and a prosthesis-patient mismatch be avoided. Additionally, pre-procedure imaging modalities must be employed to make as accurate a determination of the inner diameter as possible.

6.0 Potential Adverse Events

Potential risks associated with the overall procedure including potential access complications associated with standard cardiac catheterization, balloon valvuloplasty, the potential risks of conscious sedation and/or general anesthesia, and the use of angiography:

- Death
- Stroke/transient ischemic attack, silent cerebral ischemia, clusters or neurological deficit
- Paralysis
- Permanent disability
- Aneurysm
- Respiratory insufficiency or respiratory failure
- Hemorrhage requiring transfusion or intervention
- Cardiovascular injury including perforation or dissection of vessels, ventricle, atrium, septum, myocardium, aortic annulus, coronary ostia or valvular structures that may require intervention
- Pericardial effusion or cardiac tamponade
- Thoracic bleeding
- Embolization including air, calcific valve material or thrombus
- Infection including abscess, septicemia and endocarditis
- Heart failure
- Myocardial ischemia or infarction
- Renal insufficiency or renal failure
- Conduction system defect which may require a permanent pacemaker
- Arrhythmia
- Retroperitoneal bleed
- Arteriovenous (AV) fistula or pseudoaneurysm
- Reoperation
- Ischemia or nerve injury
- Restenosis
- Pulmonary edema
- Bleeding
- Anemia
- Abnormal lab values (including electrolyte imbalance)
- Hypertension or hypotension
- Allergic reaction to antithrombotic therapy, anesthesia, contrast media, or device materials
- Hematoma
- Syncope
- Pain or changes at the access site
- Exercise intolerance or weakness

- Inflammation
- Angina
- Heart murmur
- Fever

Additional potential risks associated with the use of the valve, delivery system, and/or accessories include:

- Cardiac arrest
- Cardiogenic shock
- Emergency cardiac surgery
- Cardiac failure or low cardiac output
- Coronary flow obstruction/transvalvular flow disturbance
- Device thrombosis requiring intervention
- Valve thrombosis
- Device embolization
- Device migration or malposition requiring intervention
- Left ventricular outflow tract obstruction
- Valve deployment in unintended location
- Valve stenosis
- Injury to mitral valve
- Structural valve deterioration (wear, fracture, calcification, leaflet tear/tearing from the stent posts, leaflet retraction, suture line disruption of components of a prosthetic valve, thickening, stenosis)
- Device degeneration
- Paravalvular or transvalvular leak
- Valve regurgitation
- Hemolysis
- Vasovagal response
- Mediastinitis
- Mediastinal bleeding
- Device explants
- Nonstructural dysfunction
- Mechanical failure of delivery system, and/or accessories including balloon rupture and tip separation
- Non-emergent reoperation

7.0 Directions for Use

7.1 System Compatibility

Table 7

Product Name	20 mm System	23 mm System	26 mm System
	Model		
Edwards SAPIEN 3 Ultra Transcatheter Heart Valve	9750TFX20	9750TFX23	9750TFX26
Edwards Certitude Delivery System	9620TA20	9630TA23	9630TA26
Edwards Certitude Introducer Sheath Set	9620IS18	9620IS21	
Inflation device, Qualcrimp crimping accessory, 2-piece crimp stopper, loader and extension tubing provided by Edwards Lifesciences			
Edwards Crimper	9600CR		

Table 8

Product Name	20 mm System	23 mm System	26 mm System	29 mm System
	Model			
Edwards SAPIEN 3 Transcatheter Heart Valve	9600TFX20	9600TFX23	9600TFX26	9600TFX29
Edwards Certitude Delivery System	9620TA20	9620TA23	9620TA26	9620TA29
Edwards Certitude Introducer Sheath Set	9620IS18			9620IS21
Inflation device, Qualcrimp crimping accessory, 2-piece crimp stopper, loader and extension tubing provided by Edwards Lifesciences				
Edwards Crimper	9600CR			

Additional Equipment:

- Balloon catheter, per the discretion of the physician
- 20 cc syringe or larger (x2)
- 50 cc syringe or larger
- Standard cardiac catheterization lab equipment
- Fluoroscopy (fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal or transthoracic echocardiography capabilities
- 18 gauge Seldinger needle (for transaortic)
- 145 cm x 0.035 inch (0.89 mm) soft guidewire
- 180 cm or 260 cm x 0.035 inch (0.89 mm) & Exchange length 0.035 inch (0.89 mm) extra-stiff guidewire
- Temporary pacemaker (PM) and pacing lead
- Sterile rinsing basins, physiological saline, heparinized saline, 15% diluted radiopaque contrast medium
- Sterile table for valve and device preparation

7.2 Valve Handling and Preparation

Follow sterile technique during device preparation and implantation.

7.2.1 Valve Rinsing Procedure

Before opening the valve jar, carefully examine for evidence of damage (e.g. a cracked jar or lid, leakage, or broken or missing seals).

CAUTION: Valves from containers found to be damaged, leaking, without adequate sterilant, or missing intact seals must not be used for implantation.

Step	Procedure
1	Set up two (2) sterile bowls with at least 500 mL of sterile physiological saline to thoroughly rinse the glutaraldehyde sterilant from the valve.
2	Carefully remove the valve/holder assembly from the jar without touching the tissue. Verify the valve serial identification number with the number on the jar lid and record in the patient information documents. Inspect the valve for any signs of damage to the frame or tissue.
3	Rinse the valve as follows: Place the valve in the first bowl of sterile, physiological saline. Be sure the saline solution completely covers the valve and holder. With the valve and holder submerged, slowly agitate (to gently swirl the valve and holder) back and forth for a minimum of 1 minute. Transfer the valve and holder to the second rinsing bowl of sterile, physiological saline and gently agitate for at least one more minute. Ensure the rinse solution in the first bowl is not used. The valve should be left in the final rinse solution until needed to prevent the tissue from drying. CAUTION: Do not allow the valve to come into contact with the bottom or sides of the rinse bowl during agitation or swirling in the rinse solution. Direct contact between the identification tag and valve is also to be avoided during the rinse procedure. No other objects should be placed in the rinse bowls. The valve should be kept hydrated to prevent the tissue from drying.

7.2.2 Prepare the Components

Step	Procedure
1	Visually inspect all components for damage. Ensure the Edwards Certitude delivery system is fully unflexed.
2	Prime and flush the introducer and sheath with heparinized saline. Hydrate the length of the introducer and sheath.
3	Advance the introducer fully into the sheath housing.
4	Unscrew the loader cap from the loader and flush the loader cap with heparinized saline.
5	Place the loader cap onto the delivery system with the inside of the cap oriented towards the tapered tip.
6	Flush the extension tubing and connect to the delivery system.
7	Partially fill a 50 mL or larger syringe with diluted contrast medium, and connect to the extension tubing.

Step	Procedure
8	Fill the inflation device with 20 mL of diluted contrast medium, lock the inflation device, and connect to the extension tubing. Close 3-way stopcock to inflation device.
9	De-air the delivery system using the luer lock syringe. Leave zero-pressure in the system. Close the 3-way stopcock to the luer lock syringe.
10	Remove 3 mL fluid from the delivery system by turning the knob of the locked inflation device. Keep the inflation device locked for valve crimping steps.

7.2.3 Mount and Crimp the Valve on the Delivery System

Step	Procedure
1	Set up two (2) additional sterile bowls with at least 100 mL of sterile, physiological saline to thoroughly rinse the Qualcrimp crimping accessory.
2	Completely submerge the Qualcrimp crimping accessory in the first bowl and gently compress it to ensure complete saline absorption. Slowly swirl the Qualcrimp crimping accessory for a minimum of 1 minute. Repeat this process in the second bowl.
3	Remove crimper from packaging.
4	Rotate the crimper handle until the aperture is fully open.
5	Remove the valve from the holder and remove the ID tag.
6	Attach the 2-piece crimp stopper to the base of the crimper and click into place.
7	With the crimper in the open position, gently place the valve into the crimper aperture. Gradually crimp the valve (if necessary) until it fits into the Qualcrimp crimping accessory.
8	Place the Qualcrimp crimping accessory over the valve making sure the valve is parallel to the edge of the Qualcrimp crimping accessory.
9	Place the valve and Qualcrimp crimping accessory in the crimper aperture. Insert the delivery system coaxially within the valve. The orientation of the valve on the delivery system is described below: Antegrade approach: Inflow (outer skirt end) of the valve towards the proximal end of the delivery system. Retrograde approach: Inflow (outer skirt end) of the valve towards the distal end of the delivery system.
10	Crimp the valve between the two internal shoulders of the delivery system until it reaches the Qualcrimp Stop located on the 2-piece Crimp Stopper.
11	Gently remove the Qualcrimp crimping accessory from the valve. Remove the Qualcrimp Stop from the Final Stop, leaving the Final Stop in place.
12	Fully crimp the valve until it reaches the Final Stop. NOTE: Ensure that the delivery system remains coaxial within the valve.
13	Repeat the full crimp of the valve two more times for a total of three full crimps.
14	Flush the loader with heparinized saline. Immediately advance the valve into the loader until the tapered tip of the delivery system is exposed and the valve is within the distal end of the loader tube. CAUTION: To prevent possible leaflet damage, the valve should not remain fully crimped and/or in the loader for over 15 minutes.
15	Attach the loader cap to the loader and flush through the port on the loader. CAUTION: Keep the valve hydrated until ready for implantation. CAUTION: The physician must verify correct orientation of the valve prior to its implantation.
16	With 3-way stopcock still closed to the luer lock syringe, unlock the inflation device. Allow the delivery system to reach zero-pressure.
17	Close the 3-way stopcock to the delivery system. Use the luer lock syringe to de-air the inflation device if necessary.

Step	Procedure															
18	<p>Adjust the inflation device to the inflation volume required to deploy the valve, per the following:</p> <table border="1"> <thead> <tr> <th>Delivery System</th> <th>Valve</th> <th>Inflation Volume</th> </tr> </thead> <tbody> <tr> <td>Model 9620TA20</td> <td>20 mm</td> <td>12 mL</td> </tr> <tr> <td>Model 9620TA23 9630TA23</td> <td>23 mm</td> <td>17 mL</td> </tr> <tr> <td>Model 9620TA26 9630TA26</td> <td>26 mm</td> <td>23 mL</td> </tr> <tr> <td>Model 9620TA29</td> <td>29 mm</td> <td>30 mL</td> </tr> </tbody> </table> <p>Re-lock the inflation device. Close the 3-way stopcock to the luer lock syringe and remove syringe.</p> <p>CAUTION: Maintain the inflation device in a locked position until valve deployment.</p>	Delivery System	Valve	Inflation Volume	Model 9620TA20	20 mm	12 mL	Model 9620TA23 9630TA23	23 mm	17 mL	Model 9620TA26 9630TA26	26 mm	23 mL	Model 9620TA29	29 mm	30 mL
Delivery System	Valve	Inflation Volume														
Model 9620TA20	20 mm	12 mL														
Model 9620TA23 9630TA23	23 mm	17 mL														
Model 9620TA26 9630TA26	26 mm	23 mL														
Model 9620TA29	29 mm	30 mL														
19	Remove the stylet and flush the guidewire lumen of the delivery system.															

7.3 Valvuloplasty and Valve Delivery

Valvuloplasty and valve delivery should be performed under conscious sedation and/or general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

The following table shows the minimum required distances from the valvular plane to the distal tip of the sheath to allow the Edwards Certitude delivery system balloon to inflate properly during valve deployment. These distances do not include sheath insertion depth, which should be considered during the transaortic approach when selecting the access site on the ascending aorta.

Delivery System	Valve	Minimum Required Distance From Sheath Tip to Valvular Plane
9620TA20	20 mm	3.5 cm
9620TA23 9630TA23	23 mm	3.5 cm
9620TA26 9630TA26	26 mm	3.5 cm
9620TA29	29 mm	4.0 cm

Administer heparin to maintain the ACT at ≥ 250 sec during the procedure.

Balloon valvuloplasty should be avoided in the treatment of failing bioprostheses as this may result in embolization of bioprosthesis material and mechanical disruption of the valve leaflets.

CAUTION: Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

CAUTION: Care should be taken to avoid damage to soft tissue, chordae, aorta, native leaflet or ventricular wall during insertion, positioning and removal of devices.

7.3.1 Baseline Parameters

Step	Procedure
1	Perform an angiogram with fluoroscopic view perpendicular to the valve.
2	Evaluate the distance of the left and right coronary ostia from the aortic annulus in relation to the valve frame height.
3	Introduce a pacemaker (PM) lead and position appropriately.
4	Set the stimulation parameters to obtain 1:1 capture, and test pacing.

7.3.2 Access

Step	Procedure
1	Gain access using standard catheterization techniques.
2	Using the sheath depth markers, advance the introducer and sheath over the guidewire to the desired depth while following its progression on fluoroscopy.
3	Withdraw the introducer slowly, keeping the sheath in place. Maintain guidewire position across the valve.

7.3.3 Valvuloplasty

Pre-dilate the native aortic valve, per the discretion of the physician, according to the instructions for use for the selected balloon aortic valvuloplasty catheter.

CAUTION: Valve implantation should not be carried out if the balloon cannot be fully inflated during valvuloplasty.

7.3.4 Valve Delivery

Step	Procedure
1	Confirm that the valve is oriented properly and the volume in the inflation device matches the indicated volume.
2	Advance the valve/balloon assembly with the loader over the guidewire.
3	Engage loader into sheath housing while maintaining a firm grip.
4	Advance the valve out of the loader into the large section of the sheath. Tap on the sheath housing to release air bubbles to the proximal end of the loader. Depress button valve on loader to de-air.
5	Advance the valve/balloon assembly through the sheath and position within the target valve. If needed, rotate the flex wheel on the handle to articulate the valve/balloon assembly into position. CAUTION: To prevent possible leaflet damage, the valve should not remain in the sheath for over 5 minutes.
6	Ensure that the valve is correctly positioned between the two internal shoulders of the delivery system.
7	Begin valve deployment: <ul style="list-style-type: none">• Unlock the Inflation device provided by Edwards Lifesciences.• Begin rapid pacing; once systolic blood pressure has decreased to 50 mmHg or below, balloon inflation can commence.• Deploy the valve by inflating the balloon with the entire volume in the Inflation device provided by Edwards Lifesciences, hold for 3 seconds and confirm that the barrel of the inflation device is empty to ensure complete inflation of the balloon.• Deflate the balloon. When the balloon has been completely deflated, turn off the pacemaker.

7.3.5 System Removal

Step	Procedure
1	If articulation was used, completely unflex the delivery system. Retract the delivery system and guidewire into the sheath. Remove the loader and delivery system from the sheath. CAUTION: Properly deflate the balloon and unflex the delivery system prior to removal.
2	Remove all devices when the ACT level is appropriate.
3	Remove the sheath from the access site, close the access site and confirm hemostasis.

8.0 How Supplied

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system is supplied sterilized with ethylene oxide gas.

8.1 Storage

The valve must be stored at 10 °C to 25 °C (50 °F to 77 °F). Each jar is shipped in an enclosure containing a temperature indicator to detect exposure of the valve to extreme temperature.

The delivery system should be stored in a cool, dry place.

9.0 MR Safety

MR Conditional

Non-clinical testing has demonstrated that the Edwards SAPIEN 3 Ultra transcatheter heart valve and the Edwards SAPIEN 3 transcatheter heart valve are MR Conditional. A patient with this device can be scanned safely, immediately after placement of this device under the following conditions:

- Static magnetic field of 1.5 tesla or 3 tesla.
- Maximum spatial gradient field of 2500 gauss/cm (25 T/m) or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the transcatheter heart valve is expected to produce a maximum temperature rise of 3.0 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 14.5 mm from the implant for spin echo images and 30 mm for gradient echo images when scanned in a 3.0T MRI system. The artifact obscures the device lumen in gradient echo images.

The implant has not been evaluated in MR systems other than 1.5 or 3.0 T.

For valve-in-valve implantation or in the presence of other implants, please refer to the MRI safety information for the surgical valve or other devices prior to MR imaging.

10.0 Patient Information

A patient registration form is provided with each THV. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the THV. Return the original form to the Edwards Lifesciences address indicated on the form and provide the temporary identification card to the patient prior to discharge.

11.0 Recovered Valve and Device Disposal

The explanted valve should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit.

Used delivery system may be disposed of in the same manner that hospital waste and biohazardous materials are handled. There are no special risks related to the disposal of these devices.

12.0 References

[1] Bapat V, Attia R, Thomas M. Effect of Valve Design on the Stent Internal Diameter of a Bioprosthetic Valve: A Concept of True Internal Diameter and Its Implications for the Valve-in-Valve Procedure. *JACC: Cardiovascular Interventions*. Vol. 7, No. 2 2014: 115-127

These products are manufactured and sold under one or more of the following US patent(s): US Patent No. 7,530,253; 7,780,723; 7,895,876; 8,382,826; 8,591,575; 8,690,936; 8,790,387; 9,061,119; 9,301,840; 9,301,841; 9,339,384; 9,393,110; and corresponding foreign patents.

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© Copyright 2022, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC