Edwards

Edwards eSheath+ Introducer Set Edwards eSheath+ införingssats Edwards eSheath+ indföringssæt

Directory ■ Register ■ Register	
English (en).....	1
Svenska (sv).....	2
Dansk (da).....	4
Symbol Legend ■ Symbolförklaring ■ Symbolforklaring.....	6

English

Instructions For Use

The product is intended for use by physicians trained and experienced in interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.

1.0 Device Description

The Edwards eSheath+ introducer set contains:

1. an expandable sheath (eSheath+) (Fig 1) that provides access into the target vessel while maintaining hemostasis and temporarily enlarges its diameter to allow for passage of a device.
2. an introducer (Fig 2) with hydrophilic coating that is used to facilitate entry and trackability of the sheath into the vessel.
3. a dilator (Fig 3) with hydrophilic coating that is used to dilate the vessel to accommodate the sheath.
4. an expansion tool (Fig 4) that is used to pre-expand the sheath during device preparation.

Figure 1: Sheath


Figure 2: Introducer


Figure 3: Dilator


Figure 4: Expansion Tool

	914ESP	916ESP
Sheath I.D. (unexpanded)	14F (4.6 mm)	16F (5.3 mm)
Sheath O.D. (unexpanded)	6.0 mm	6.7 mm

Edwards, Edwards Lifesciences, the stylized E logo, Edwards eSheath+, eSheath+, SAPIEN, SAPIEN 3, and SAPIEN 3 Ultra are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

	914ESP	916ESP
Compatible THV	20 mm 23 mm 26 mm	29 mm
Introducer O.D.	14F	16F
Dilator O.D.	16F	18F

2.0 Intended Use

The product is intended to be used to gain access to the vasculature.

3.0 Indications

The Edwards eSheath+ introducer set is indicated for the introduction and removal of SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valve systems into the vascular system.

4.0 Contraindications

There are no known contraindications.

5.0 Warnings

The devices are designed, intended, and distributed for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.

The Edwards eSheath+ introducer set must be used with a compatible 0.035 in (0.89 mm) guidewire to prevent vessel injury.

Do not mishandle the device or use it if the packaging or any components are not sterile, have been opened or are damaged (i.e., kinked or stretched, etc.), or the expiration date has elapsed.

6.0 Precautions

- Expansion tool does not contain a hydrophilic coating. Do not use as a dilator.
- The sheath temporarily enlarges to allow the passage of devices; ensure that the vasculature can accommodate the maximum diameter of the expanded sheath.
- When inserting, manipulating or withdrawing a device through the sheath, always maintain orientation of the sheath position.
- When puncturing, suturing or incising the tissue near the sheath, use caution to avoid damage to the sheath.
- Caution should be used in vessels that have diameters less than 5.5 mm or 6 mm as it may preclude safe placement of the 14F and 16F Edwards eSheath+ introducer set respectively.
- Use caution in tortuous or calcified vessels that would prevent safe entry of the introducer set.

7.0 Potential Adverse Events

Complications associated with standard catheterization and use of angiography include, but are not limited to, allergic reaction to anesthesia or to contrast media; injury including perforation or dissection of vessels; injury at the site of access that might require vessel repair; thrombosis and/or plaque dislodgment which may result in emboli formation; distal vessel obstruction; stroke; ischemia and/or death.

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident

has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

8.0 Directions for Use

1. Visually inspect the introducer, dilator, expansion tool and sheath for surface defects and damage.
2. Flush the introducer and dilator using heparinized saline through the guidewire lumen.
3. Hydrate the length of the introducer, dilator, and sheath with heparinized saline to activate the hydrophilic coating.
4. Wet the surface of the expansion tool.
5. Flush the sheath using heparinized saline through the flush port; close the flush port.
6. Use the expansion tool to pre-expand the partially expandable portion of the sheath prior to procedural use.

Note: After pre-expanding the sheath, inspect the length of the expandable portion for damage prior to use.

7. After removing the expansion tool, flush the sheath a second time using the heparinized saline through the flush port; close the flush port.
8. Insert the introducer completely into the sheath and turn clockwise to lock the introducer hub to the sheath hub.
9. Using standard catheterization techniques, gain access to the vessel and dilate as necessary with the dilator to accommodate the sheath.
10. Orient the sheath appropriately and maintain orientation throughout the procedure. Insert the sheath assembly using standard technique while following its progression under fluoroscopy.

Note: The proximal tapered end of the sheath working length is larger in diameter.

11. If possible, suture the sheath into place using the suture ring(s).
12. Remove the introducer from the sheath by turning counterclockwise to unlock the introducer hub from the sheath.
13. Insert the device into the sheath.

Note: The sheath should be intermittently flushed with heparinized saline throughout the procedure, per standard interventional technique.

14. After the completion of the procedure and removal of the device, remove the suture, and then remove the sheath entirely without torquing and do not reinsert.

9.0 How Supplied

The Edwards eSheath+ introducer set is supplied in a pouch and sterilized with ethylene oxide.

10.0 Storage

The Edwards eSheath+ introducer set should be stored in a cool, dry place.

11.0 Device Disposal

Used sheath sets may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

12.0 Hazardous Substances

This medical device does not contain hazardous substances.

13.0 Summary of Safety and Clinical Performance (SSCP)

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information.

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/ Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

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
	English	Svenska	Dansk
	Contains biological material of animal origin	Innehåller biologiskt material av animaliskt ursprung	Indeholder biologisk materiale af animalsk oprindelse
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Time & Temperature Sensitive	Tids- och temperaturkänslig	Følsom overfor tid og temperatur
	Contains hazardous substances	Innehåller farliga ämnen	Indeholder farlige stoffer
SZ	Size	Storlek	Størrelse

Note: Not all symbols may be included in the labeling of this product. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne.

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2023-07
10056055001 A
© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU

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Edwards Lifesciences GmbH