

# Edwards SAPIEN 3 Transcatheter Pulmonary Valve System with Alterra Adaptive Prestent

## **Instructions for Use**

Implantation of the Edwards Alterra adaptive prestent should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in balloon valvuloplasty.

## **1.0 Device Description**

### **Edwards Alterra Adaptive Prestent System**

The Edwards Alterra adaptive prestent system consists of an Alterra adaptive prestent that is fully loaded in an Alterra delivery system and supplied together in one package.

#### • Edwards Alterra Adaptive Prestent (Figure 1)

The Edwards Alterra adaptive prestent is used as a docking adaptor for the 29 mm Edwards SAPIEN 3 transcatheter heart valve (THV). It is comprised of a self-expanding, radiopaque, nitinol frame assembly with polyethylene terephthalate (PET) fabric covering. The prestent has designated inflow and outflow sides. The proximal inflow section is identifiable by the presence of two triangular tabs (prestent connector) that are attached to the catheter of the delivery system. The distal outflow section is distinguished by the open cells for blood flow. The PET fabric is attached by sutures to the inside surface of the frame to create sealing at the inflow section and opening for the outflow. Sutures are also used in the center to support the middle section when an Edwards SAPIEN 3 transcatheter heart valve is implanted. Three (3) radiopaque fluoroscopically visible markers at the prestent waist help with positioning.



#### Figure 1: Edwards Alterra Adaptive Prestent

#### Table 1

Inflow Sealing OD	Outflow Sealing OD	Height	
39 mm	41 mm	49 mm	

Sizing recommendation for implanting the prestent in the right ventricular outflow tract/pulmonary valve (RVOT/PV) landing zone are provided in Table 2:

#### **Table 2: Prestent Sizing in RVOT landing zone**

Perimeter	Perimeter Derived Diameter <sup>1</sup>	Prestent Size Diameter <sup>2</sup> x Length	Valve Size
84.9 mm - 119.3 mm 27 mm - 38 mm		40 mm x 49 mm	29 mm

<sup>1</sup> Diameter range during systole

<sup>2</sup> Diameter is average of inflow and outflow diameters

Note: For Edwards SAPIEN 3 transcatheter heart valve implantation, refer to the Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestent instructions for use.

Edwards Alterra Delivery System (Figure 2)

Edwards, Edwards Lifesciences, the stylized E logo, Alterra, Edwards SAPIEN, Edwards SAPIEN 3, Qualcrimp, SAPIEN, and SAPIEN 3 are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

The delivery system includes a handle which consists of a wheel that allows for deployment, two primary shafts with a flush port to flush the delivery system, and a long tapered tip at the distal end to facilitate tracking through the vasculature. A fluoroscopically visible radiopaque marker band shows the location of the tip of the outer shaft. The prestent is fully loaded in the delivery system. A stylet is included within the guidewire lumen.



#### **Additional Accessories**

#### Edwards Sheath

Refer to the Edwards sheath instructions for use for device description.

#### 2.0 Intended Use

The bioprosthesis with prestent is intended for use in patients requiring pulmonary heart valve replacement. The delivery systems and accessories are intended to facilitate placement of the bioprosthesis and prestent via the transfemoral access approach.

#### 3.0 Indications

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestent is indicated for use in the management of patients with pulmonary regurgitation who have a native or surgically-repaired right ventricular outflow tract and are clinically indicated for pulmonary valve replacement.

#### 4.0 Contraindications

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestent is contraindicated in patients who:

· Cannot tolerate an anticoagulation/antiplatelet regimen or who have active bacterial endocarditis or other active infections.

#### 5.0 Warnings

- The devices are designed, intended, and distributed STERILE for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- Assessment for coronary compression risk prior to implantation is essential to prevent the risk of severe patient harm.
- Patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, and/or polymeric materials may have an allergic reaction to these materials.
- Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging sterile barriers and any components have been opened or damaged, cannot be flushed, or the expiration date has elapsed.
- The procedure should be conducted under fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.

#### 6.0 Precautions

- Correct sizing of the prestent into the RVOT is essential to minimize risks such as paravalvular leak, migration, embolization, and/or RVOT rupture.
- Long-term durability has not been established for the prestent. Medical follow-up is advised so that device related complications can be diagnosed and properly managed.
- Patient venous anatomy should be evaluated to prevent the risk of access that would preclude the delivery and deployment of the device.
- Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.
- Patient should be heparinized to maintain the ACT at ≥ 250 sec prior to introduction of the delivery system in order to prevent thrombosis.
- Prestent recipients should be maintained on anticoagulant/antiplatelet therapy, except when contraindicated, as determined by their physician. This device has not been tested for use without antiplatelet therapy.
- It is recommended that all prestent recipients be prophylactically treated for endocarditis to minimize the risk of infection.
- If a prestent fracture is detected with significant loss in valve functionality, reintervention should be considered.
- Safety and effectiveness have not been established for patients with the following characteristics/comorbidities:
- Blood dyscrasias defined as: leukopenia, acute anemia, thrombocytopenia, or history of bleeding diathesis or coagulopathy
- A known hypersensitivity or contraindication to aspirin, heparin, ticlopidine (Ticlid<sup>™</sup>), or clopidogrel (Plavix<sup>™</sup>), or sensitivity to contrast media, which cannot be adequately premedicated
- · Positive urine or serum pregnancy test in female patients of child-bearing potential

## 7.0 Potential Adverse Events

Potential risks associated with the anesthesia, interventional procedure and imaging include but are not limited to:

- Death
- Stroke/transient ischemic attack
- Respiratory insufficiency or respiratory failure
- Cardiovascular or vascular injury, such as perforation or damage (dissection) of vessels, myocardium or valvular structures including rupture of the RVOT that may require intervention
- Pericardial effusion/cardiac tamponade
- Cardiac failure
- Embolic event: air, calcific material, thrombus
- Infection including incisional site infection, septicemia and endocarditis
- Myocardial infarction
- · Renal insufficiency or renal failure
- Conduction system injury
- Arrhythmia
- Deep vein thrombosis
- Arteriovenous (AV) fistula
- Systemic or peripheral nerve injury
- Systemic or peripheral ischemia
- Pulmonary edema
- Pneumothorax
- Pleural effusion
- Dyspnea
- Atelectasis
- · Dislodgement of previously implanted devices (i.e., pacing lead)
- Blood loss requiring transfusion
- Anemia
- Radiation injury
- Electrolyte imbalance
- Hypertension or hypotension
- · Allergic reaction to anesthesia, contrast media, antithrombotic therapy, device materials
- Hematoma or ecchymosis
- Syncope
- Pain
- Exercise intolerance or weakness
- Inflammation
- Angina
- Fever

Potential risks that may or may not require intervention associated with the prestent, delivery system and/or accessories include, but may not be limited to, the following:

- Cardiac arrest
- Cardiogenic shock
- Coronary flow obstruction/transvalvular flow disturbance
- Device thrombosis
- · Injury to tricuspid valve
- Device fracture
- Device embolization
- Device migration or malposition
- Endocarditis
- Chest pain/discomfort
- Device penetration/perforation into surrounding vasculature
- Device dysfunction (regurgitation and/or stenosis)
- Aortic root distortion
- Embolic event: device fragments
- · Mechanical failure of delivery system, and/or accessories

## 8.0 Directions for Use

## 8.1 System Compatibility

Table 3

Product Name	Model
Edwards Alterra adaptive prestent system <sup>[1]</sup>	29AP4045

Product Name	Model
Edwards eSheath+ introducer set <sup>[2]</sup>	916ESP
or	or
Edwards eSheath introducer set <sup>[2]</sup>	9610ES16

<sup>[1]</sup> Includes an Alterra adaptive prestent that is fully loaded in an Alterra delivery system

<sup>[2]</sup> Sheath provided by Edwards Lifesciences or equivalent

Additional Equipment:

- Balloon tip catheter
- Sizing balloons
- 20 cc syringe or larger
- 50 cc syringe or larger
- · Standard cardiac catheterization lab equipment and supplies, and access to standard heart valve operating room equipment and supplies
- Fluoroscopy (bi-plane, fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Exchange length 0.035 in (0.89 mm) stiff guidewire
- Physiological saline
- Sterile table for device preparation

Refer to the Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestent instructions for use for additional materials required to prepare the Edwards SAPIEN 3 transcatheter heart valve system.

## 8.2 Device Handling and Preparation

Follow sterile technique during device preparation and implantation.

#### 8.2.1 Prepare the System

Refer to the Edwards sheath instructions for use for device preparation.

- 1. Remove the delivery catheter with the preloaded prestent from packaging. Visually inspect all components for damage.
- 2. Ensure there is a small gap between the outer shaft and the tapered tip to facilitate flushing of the inner lumen. If needed retract the outer shaft using the deployment wheel.

#### Note: Do not allow the end of the prestent to begin to exit the delivery system.

- 3. With the stylet still in place, flush the guidewire lumen with heparinized saline.
- 4. Using the deployment wheel, re-advance the outer shaft until it is even with the tapered tip.

## Note: Do not overdrive the outer shaft onto the tapered tip.

5. Remove the stylet and repeat flushing the guidewire lumen.

## 8.3 Prestent Delivery

Prestent delivery should be performed under local and/or general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

Administer heparin to maintain the ACT at  $\geq$  250 sec.

## CAUTION: Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

- 1. Gain access using standard catheterization techniques.
- 2. If necessary, predilate the vessel.
- 3. Introduce the sheath per its instructions for use.
- 4. Insert and advance the delivery system to the RVOT landing zone.
- Note: Advance the delivery system from the shaft. Do not push the delivery system in by using the handle. Do not rotate the deployment wheel during advancement of the delivery system.
- 5. Position the delivery system fluoroscopically visible marker band distal to the intended landing zone.
- 6. Begin deployment by rotating the deployment wheel to retract the outer shaft.

Note: The delivery system marker band is located slightly proximal to the distal edge of the outer shaft.

## Note: Waist markers on the prestent indicate the middle of the prestent.

Note: The prestent can be recaptured into the outer shaft and repositioned if deployed approximately 65%.

CAUTION: Once deployment has begun, do not reposition the device more distally. Advancement of the device with the prestent exposed may increase the risk for vascular damage.

- 7. Continue deploying the Alterra pausing at approximately 30%, 50%, and 65% to assess for Alterra positioning, coaxiality, and deployment angulation.
- 8. If needed, recapture and reposition the prestent by rotating the deployment wheel in the reverse direction from deployment until the outer shaft fully covers the prestent as shown by the delivery system marker band.

Note: The prestent can be recaptured and redeployed one time. If a second recapture of the partially deployed prestent is performed, remove and replace the device.

Note: Several rotations of the deployment wheel may be necessary before the Alterra begins to be recaptured.

CAUTION: Do not overdrive the outer shaft onto the tapered tip when recapturing the prestent. This may cause the delivery system to cinch down on the guidewire preventing independent movement of the delivery system and guidewire.

CAUTION: Recapturing and redeploying a prestent more than one time may impact implant integrity.

CAUTION: Recapturing a prestent that has been deployed more than 65% may cause system damage.

- 9. After achieving an acceptable position, completely deploy the prestent by continuing to rotate the deployment wheel until the delivery system marker band is beyond the prestent connector.
- 10. Confirm release of prestent.

CAUTION: Failure to identify release of the prestent connector tabs from the prestent connector may lead to prestent embolization during removal of the Alterra delivery system.

#### 8.4 System Removal

- Slowly retract the system through the prestent. Remove the delivery system.
  CAUTION: Ensure that the tapered tip and delivery system do not interfere with the prestent upon removal to prevent movement of the prestent.
- Assess Alterra prestent stability by evaluating apices engagement in surrounding tissue, wall apposition, and/or motion of prestent within the anatomy. If adequate stability is not noted, consider staging valve deployment after allowing sufficient time for prestent endothelialization.
   CAUTION: Failure to identify prestent instability may lead to prestent migration/embolization when tracking interventional devices through the prestent.

## 9.0 How Supplied

STERILE: The Edwards Alterra adaptive prestent system is supplied pouched and sterilized by e-beam sterilization. The sheath is supplied sterilized with ethylene oxide gas.

#### 9.1 Storage

The prestent and delivery system must be stored in a cool, dry place.

## 10.0 MR Safety



## MR Conditional

Non-clinical testing has demonstrated that the Edwards Alterra adaptive prestent, alone or with a deployed SAPIEN 3 transcatheter heart valve, is MR Conditional. A patient can be scanned safely immediately after placement of this implant in an MR system meeting the following conditions:

- Static magnetic fields of 1.5 Tesla and 3.0 Tesla
- Spatial magnetic gradient field of 3000 Gauss/cm (30 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) scanning per sequence
- · Gradient system is in normal operating mode

Under the scan conditions defined above, the Edwards Alterra adaptive prestent, alone or with a deployed SAPIEN 3 transcatheter heart valve, is expected to produce a maximum temperature rise of 4.0 °C or less after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 15 mm for gradient echo images when scanned using a 3.0 T MRI system. The artifact obscures the device lumen in spin and gradient echo images.

The delivery system has not been evaluated for MR compatibility and is considered MR unsafe.

## **11.0 Patient Information**

Patient education brochures are provided to each site and should be given to the patient to inform them of the risks and benefits of the procedure and alternatives in adequate time before the procedure to be read and discussed with their physician. A copy of this brochure may also be obtained from Edwards Lifesciences by calling 1.800.822.9837.

A patient implant card is provided with each prestent. After implantation, please complete all requested information and provide the implant card to the patient. The serial number is found on the package. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

## **12.0 Recovered Prestent and Device Disposal**

The explanted prestent should be placed into a suitable container and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an explant kit.

Used devices may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

## Symbol Legend

	English		English		English
REF	Reorder Number		Sterilized using irradiation		Importer
#	Model Number	STERNIZE	Do not resterilize		Balloon diameter
	Usable length	eSheath	eSheath compatibility		Balloon working length
$\otimes$	Do not re-use	eSheath™	eSheath compatibility	20 mm	For use with size 20 mm Edwards transcatheter heart valve
LOT	Lot Number	$\bigcup$	Single sterile barrier system	23 mm	For use with size 23 mm Edwards transcatheter heart valve
			Single sterile barrier system with protective packaging inside	26 mm	For use with size 26 mm Edwards transcatheter heart valve
	Caution	ΟΤΥ	Quantity	29 mm	For use with size 29 mm Edwards transcatheter heart valve
eifu.edwards.com	Consult instructions for use		Use-by date	MR	[Implant only] The implant device has been determined to be MR Conditional when used under the conditions listed in the instruc-
+ 1 888 570 4016	website				tions for use.
	Do not use if package is damaged and consult instructions for use	SN	Serial Number		Contents
	Exterior diameter		Manufacturer	X	Non-pyrogenic
$\bigcirc$	Inner diameter		Date of manufacture	MD	Medical device
**	Store in a cool, dry place	EC REP	Authorized representative in the European Community/European Union	BIO	Contains biological material of an- imal origin
Ť	Keep dry	GWC	Guidewire compatibility	Rx only	Caution: Federal (USA) law re- stricts this device to sale by or on the order of a physician.
*	Keep away from sunlight	NP	Nominal Pressure		Time & Temperature Sensitive
UDI	Unique Device Identifier	RBP	Rated burst pressure		Contains hazardous substances
	Temperature limit		Recommended guidewire length	SZ	Size
		Sheath 🖉	Minimum sheath size		
	Sterile	Catheter 🖯	Catheter shaft size		
STERILEEO	Sterilized using ethylene oxide	<u> </u>	I	l	

Note: Not all symbols may be included in the labeling of this product.



2024-01 10059844001 A <sup>©</sup> Copyright 2024, Edwards Lifesciences LLC All rights reserved.

Telephone +1.949.250.2500 +1.800.424.3278 FAX +1.949.250.2525

Web IFU