Edwards

Edwards eSheath+ Introducer Set  
Conjunto de introductor Edwards eSheath+  
Conjunto introdutor Edwards eSheath+

| Directory ■ Directorio ■ Diretório                                     |   |
|------------------------------------------------------------------------|---|
| English (en).....                                                      | 1 |
| Español (es).....                                                      | 4 |
| Português (pt).....                                                    | 6 |
| Symbol Legend ■ Significado de los símbolos ■ Legenda de símbolos..... | 8 |

English

Instructions For Use

The product is intended for use by physicians trained and experienced in interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.

1.0 Device Description

The Edwards eSheath+ introducer set contains:

- 1. an expandable sheath (eSheath+) (Fig 1) that provides access into the target vessel while maintaining hemostasis and temporarily enlarges its diameter to allow for passage of a device.
- 2. an introducer (Fig 2) with hydrophilic coating that is used to facilitate entry and trackability of the sheath into the vessel.
- 3. a dilator (Fig 3) with hydrophilic coating that is used to dilate the vessel to accommodate the sheath.
- 4. an expansion tool (Fig 4) that is used to pre-expand the sheath during device preparation.


Figure 1: Sheath


Figure 2: Introducer

Figure 3: Dilator

Figure 4: Expansion Tool

|                          | 914ESP                  | 916ESP       |
|--------------------------|-------------------------|--------------|
| Sheath I.D. (unexpanded) | 14F (4.6 mm)            | 16F (5.3 mm) |
| Sheath O.D. (unexpanded) | 6.0 mm                  | 6.7 mm       |
| Compatible THV           | 20 mm<br>23 mm<br>26 mm | 29 mm        |
| Introducer O.D.          | 14F                     | 16F          |

|              | 914ESP | 916ESP |
|--------------|--------|--------|
| Dilator O.D. | 16F    | 18F    |

## 2.0 Intended Use

The product is intended to be used to gain access to the vasculature.

## 3.0 Indications

The Edwards eSheath+ introducer set is indicated for the introduction and removal of SAPIEN 3, SAPIEN 3 Ultra, and SAPIEN 3 Ultra RESILIA transcatheter heart valve systems into the vascular system.

## 4.0 Contraindications

There are no known contraindications.

## 5.0 Warnings

The devices are designed, intended, and distributed for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.

The Edwards eSheath+ introducer set must be used with a compatible 0.035 in (0.89 mm) guidewire to prevent vessel injury.

Do not mishandle the device or use it if the packaging or any components are not sterile, have been opened or are damaged (i.e., kinked or stretched, etc.), or the expiration date has elapsed.

## 6.0 Precautions

- Expansion tool does not contain a hydrophilic coating. Do not use as a dilator.
- The sheath temporarily enlarges to allow the passage of devices; ensure that the vasculature can accommodate the maximum diameter of the expanded sheath.
- When inserting, manipulating or withdrawing a device through the sheath, always maintain orientation of the sheath position.
- When puncturing, suturing or incising the tissue near the sheath, use caution to avoid damage to the sheath.
- Caution should be used in vessels that have diameters less than 5.5 mm or 6 mm as it may preclude safe placement of the 14F and 16F Edwards eSheath+ introducer set respectively.
- Use caution in tortuous or calcified vessels that would prevent safe entry of the introducer set.

## 7.0 Potential Adverse Events

Complications associated with standard catheterization and use of angiography include, but are not limited to, allergic reaction to anesthesia or to contrast media; injury including perforation or dissection of vessels; injury at the site of access that might require vessel repair; thrombosis and/or plaque dislodgment which may result in emboli formation; distal vessel obstruction; stroke; ischemia and/or death.

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en).

## 8.0 Directions for Use

1. Visually inspect the introducer, dilator, expansion tool and sheath for surface defects and damage.
2. Flush the introducer and dilator using heparinized saline through the guidewire lumen.
3. Hydrate the length of the introducer, dilator, and sheath with heparinized saline to activate the hydrophilic coating.
4. Wet the surface of the expansion tool.
5. Flush the sheath using heparinized saline through the flush port; close the flush port.
6. Use the expansion tool to pre-expand the partially expandable portion of the sheath prior to procedural use.

**Note: After pre-expanding the sheath, inspect the length of the expandable portion for damage prior to use.**

7. After removing the expansion tool, flush the sheath a second time using the heparinized saline through the flush port; close the flush port.
8. Insert the introducer completely into the sheath and turn clockwise to lock the introducer hub to the sheath hub.
9. Using standard catheterization techniques, gain access to the vessel and dilate as necessary with the dilator to accommodate the sheath.
10. Orient the sheath appropriately and maintain orientation throughout the procedure. Insert the sheath assembly using standard technique while following its progression under fluoroscopy.

**Note: The proximal tapered end of the sheath working length is larger in diameter.**

11. If possible, suture the sheath into place using the suture ring(s).
12. Remove the introducer from the sheath by turning counterclockwise to unlock the introducer hub from the sheath.
13. Insert the device into the sheath.

**Note: The sheath should be intermittently flushed with heparinized saline throughout the procedure, per standard interventional technique.**

14. After the completion of the procedure and removal of the device, remove the suture, and then remove the sheath entirely without torquing and do not reinsert.

## 9.0 How Supplied

The Edwards eSheath+ introducer set is supplied in a pouch and sterilized with ethylene oxide.

## 10.0 Storage

The Edwards eSheath+ introducer set should be stored in a cool, dry place.

---

## **11.0 Device Disposal**

Used sheath sets may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

## **12.0 Hazardous Substances**

This medical device does not contain hazardous substances.

## **13.0 Summary of Safety and Clinical Performance (SSCP)**

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information.

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

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
|                                                                                   | English                                                                                     | Español                                                                                     | Português                                                                                   |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|  | Medical device                                                                              | Producto sanitario                                                                          | Dispositivo médico                                                                          |
|  | Contains biological material of animal origin                                               | Contiene materiales biológicos de origen animal                                             | Contém material biológico de origem animal                                                  |
| <b>Rx only</b>                                                                    | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|  | Time & Temperature Sensitive                                                                | Sensible al tiempo y a la temperatura                                                       | Sensível à passagem do tempo e à temperatura                                                |
|  | Contains hazardous substances                                                               | Contiene sustancias peligrosas                                                              | Contém substâncias perigosas                                                                |
|  | Size                                                                                        | Tamaño                                                                                      | Tamanho                                                                                     |

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2023-10  
10057702001 A  
© Copyright 2023, Edwards Lifesciences LLC  
All rights reserved.

**Edwards Lifesciences LLC**  
One Edwards Way  
Irvine, CA 92614 USA

**Edwards Lifesciences B.V.**  
Verlengde Poolseweg 16  
4818 CL Breda, Netherlands

Telephone +1.949.250.2500  
+1.800.424.3278  
FAX +1.949.250.2525

Web IFU

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**Edwards Lifesciences GmbH**