Edwards

Edwards SAPIEN 3 Transcatheter Pulmonary Valve System with Alterra Adaptive PreStent
Sistema de válvula pulmonar transcatéter Edwards SAPIEN 3 con endoprótesis adaptativa previa Alterra
Sistema de válvula pulmonar transcateter Edwards SAPIEN 3 com pré-stent adaptativo Alterra

Directory ■ Directorio ■ Diretório	
English (en).....	1
Español (es).....	11
Português (pt).....	22
Symbol Legend ■ Significado de los símbolos ■ Legenda de símbolos.....	33

English

Instructions for Use - Pulmonic

Implantation of the transcatheter heart valve and the adaptive preStent should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in balloon valvuloplasty.


1.0 Device Description

Edwards SAPIEN 3 Transcatheter Pulmonary Valve (TPV) System with Alterra Adaptive PreStent

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive preStent consists of the Edwards 29 mm SAPIEN 3 transcatheter heart valve (THV), the Edwards Commander delivery system, the Edwards Alterra adaptive preStent system and accessories.

• **Edwards SAPIEN 3 Transcatheter Heart Valve (Figure 1)**

The Edwards SAPIEN 3 transcatheter heart valve is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) fabric skirt. The leaflets are treated according to the Carpentier-Edwards TheraFix process.

9600TFX

Table 1

Valve Size	Valve Height
29 mm	22.5 mm

Figure 1: Edwards SAPIEN 3 Transcatheter Heart Valve

• **Edwards Commander Delivery System (Figures 2, 3, 4)**


Edwards, Edwards Lifesciences, the stylized E logo, Alterra, Carpentier-Edwards, Commander, Edwards Commander, Edwards eSheath, Edwards eSheath+, Edwards SAPIEN, Edwards SAPIEN 3, eSheath, eSheath+, Qualcrimp, SAPIEN, SAPIEN 3, and TheraFix are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

The Edwards Commander delivery system facilitates the placement of the bioprosthesis. It consists of a flex catheter to aid in valve alignment to the balloon, tracking, and positioning of the valve. The delivery system includes a tapered tip to facilitate crossing of the valve. The handle contains a flex wheel to control flexing of the flex catheter, and a balloon lock and fine adjustment wheel to facilitate valve alignment and positioning of the valve within the target location. A stylet is included within the guidewire lumen of the delivery system. The balloon catheter has radiopaque valve alignment markers defining the working length of the balloon. A radiopaque center marker in the balloon is provided to help with valve positioning. A radiopaque triple marker proximal to the balloon indicates the flex catheter position during deployment.

The inflation parameters for valve deployment are:

Table 2

Model	Nominal Balloon Diameter	Nominal Inflation Volume	Rated Burst Pressure (RBP)
9610TF29	29 mm	33 ml	7 atm

1. Center Marker
2. Valve Crimp Section
3. Flex Catheter
4. Flex Indicator
5. Flush Port
6. Strain Relief
7. Volume Indicator
8. Tapered Tip
9. Valve Alignment Markers
10. Triple Marker
11. Flex Wheel
12. Fine Adjustment Wheel
13. Balloon Lock
14. Balloon Catheter
15. Balloon Inflation Port
16. Guidewire Lumen

Figure 2: Edwards Commander Delivery System

• Edwards Alterra Adaptive Prestant System

Refer to the Edwards Alterra adaptive prestant system instructions for use.

Additional Devices and Accessories


Figure 3: Loader


Figure 4: Qualcrimp Crimping Accessory


Figure 5: Edwards Crimper and 2-piece Crimp Stopper

- **Loader (Figure 3)**

The loader is used to aid insertion of the delivery system into the sheath.

- **Edwards Sheath**

Refer to the Edwards sheath instructions for use for device description.

- **Qualcrimp Crimping Accessory (Figure 4)**

The Qualcrimp crimping accessory is used during THV crimping. It is packaged with the Edwards Commander delivery system.

- **Edwards Crimper and Crimp Stopper (Figure 5)**

The Edwards crimper reduces the diameter of the valve to mount it onto the delivery system. The crimper is comprised of a housing and a compression mechanism that is closed with a handle located on the housing. A 2-piece crimp stopper is used to crimp the valve to its intended diameter. The 2-piece crimp stopper is packaged with the Edwards Commander delivery system.

2.0 Intended Use

The bioprosthesis with prestent is intended for use in patients requiring pulmonary heart valve replacement. The delivery systems and accessories are intended to facilitate placement of the bioprosthesis and prestent via the transfemoral access approach.

3.0 Indications

The 29 mm Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestent is indicated for use in the management of patients with pulmonary regurgitation who have a native or surgically-repaired right ventricular outflow tract and are clinically indicated for pulmonary valve replacement.

4.0 Contraindications

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestent is contraindicated in patients who:

- Cannot tolerate an anticoagulation/antiplatelet regimen or who have active bacterial endocarditis or other active infections.

5.0 Warnings

- The devices are designed, intended, and distributed STERILE for single use only. Do not resterilize or reuse the devices. There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- The physician must verify correct orientation of the valve prior to its implantation; the inflow (outer skirt end) of the valve should be oriented towards the proximal end (handle) of the delivery system to prevent the risk of severe patient harm.
- Accelerated deterioration of the THV may occur in patients with an altered calcium metabolism.
- Assessment for coronary compression risk prior to valve implantation is essential to prevent the risk of severe patient harm.
- The THV must remain hydrated at all times and cannot be exposed to solutions, antibiotics, chemicals, etc. other than its shipping storage solution and sterile physiologic saline solution to prevent leaflet damage that may impact valve functionality. THV leaflets mishandled or damaged during any part of the procedure will require replacement of the THV.
- Patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, bovine tissue, and/or polymeric materials may have an allergic reaction to these materials.
- Do not use the THV if the tamper evident seal is broken, as sterility may be compromised.
- Do not use the THV if the temperature indicator has been activated, as valve function may be compromised.
- Do not use the THV if the expiration date has elapsed, as either sterility or valve function may be compromised.
- Do not use the THV if the storage solution does not completely cover the THV or the THV is damaged.
- Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging sterile barriers and any components have been opened or damaged, cannot be flushed, or the expiration date has elapsed.
- Patient injury could occur if the Edwards Commander delivery system is not un-flexed prior to removal.
- Valve recipients should be maintained on anticoagulant/antiplatelet therapy, except when contraindicated, to minimize the risk of valve thrombosis or thromboembolic events, as determined by their physicians. This device has not been tested for use without anticoagulation.
- The procedure should be conducted under fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.

6.0 Precautions

- Long-term durability has not been established for the THV. Regular medical follow-up is advised to evaluate valve performance.
- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Material Safety Data Sheet available from Edwards Lifesciences.

- The safety and effectiveness of the THV implantation has not been established for patients who have:
 - Blood dyscrasias defined as: leukopenia, acute anemia, thrombocytopenia, or history of bleeding diathesis or coagulopathy
 - A known hypersensitivity or contraindication to aspirin, heparin, ticlopidine (Ticlid™), or clopidogrel (Plavix™), or sensitivity to contrast media, which cannot be adequately premedicated
 - Positive urine or serum pregnancy test in female subjects of child-bearing potential
 - A concomitant paravalvular leak where the failing bioprosthesis is not securely fixed in the native annulus or is not structurally intact (e.g., wireframe fracture)
- If a significant increase in resistance occurs when advancing the catheter through the vasculature, stop advancement and investigate the cause of resistance before proceeding. Do not force passage, as this could increase the risk of vascular complications.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Do not overinflate the deployment balloon, as this may prevent proper valve leaflet coaptation and thus impact valve functionality.
- Patient venous anatomy should be evaluated to prevent the risk of access that would preclude the delivery and deployment of the device.
- Patient should be heparinized to maintain the ACT at ≥ 250 sec prior to introduction of the delivery system in order to prevent thrombosis.
- Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

7.0 Potential Adverse Events

Potential risks associated with the anesthesia, interventional procedure and imaging include but are not limited to:

- Death
- Stroke/transient ischemic attack
- Respiratory insufficiency or respiratory failure
- Cardiovascular or vascular injury, such as perforation or damage (dissection) of vessels, myocardium or valvular structures including rupture of the RVOT that may require intervention
- Pericardial effusion/cardiac tamponade
- Cardiac failure
- Embolic event: air, calcific material, thrombus, device fragments
- Infection including incisional site infection, septicemia and endocarditis
- Myocardial infarction
- Renal insufficiency or renal failure
- Conduction system injury
- Arrhythmia
- Deep vein thrombosis
- Arteriovenous (AV) fistula
- Systemic or peripheral nerve injury
- Systemic or peripheral ischemia
- Pulmonary edema
- Pneumothorax
- Pleural effusion
- Dyspnea
- Atelectasis
- Dislodgement of previously implanted devices (i.e., pacing lead)
- Blood loss requiring transfusion
- Anemia
- Radiation injury
- Electrolyte imbalance
- Hypertension or hypotension
- Allergic reaction to anesthesia, contrast media, antithrombotic therapy, device materials
- Hematoma or ecchymosis
- Syncope
- Pain
- Exercise intolerance or weakness
- Inflammation
- Angina
- Fever

Potential risks that may or may not require intervention associated with the valve, delivery system and/or accessories include, but may not be limited to, the following:

- Cardiac arrest
- Cardiogenic shock
- Coronary flow obstruction/transvalvular flow disturbance
- Device thrombosis
- Injury to tricuspid valve
- Device embolization
- Device acute migration or malposition
- Endocarditis
- Chest pain/discomfort
- Hemolysis / hemolytic anemia
- Device dysfunction (regurgitation and/or stenosis)
- Aortic root distortion
- Embolic event: device fragments

- Mechanical failure of delivery system, and/or accessories
- See Section 12 for adverse events that occurred in the clinical study.

8.0 Directions for Use

8.1 System Compatibility

Table 3

Product Name	Model
Edwards SAPIEN 3 transcatheter heart valve (29 mm)	9600TFX29
Edwards Commander delivery system	9610TF29
Edwards eSheath+ introducer set ^[1] or Edwards eSheath introducer set ^[1]	916ESP or 9610ES16
Edwards crimper	9600CR
Edwards Alterra adaptive prestant system ^[2]	29AP4045
Inflation device	96406
Qualcrimp crimping accessory, crimp stopper and loader provided by Edwards Lifesciences	

^[1] Sheath provided by Edwards Lifesciences or equivalent

^[2] Includes an Alterra adaptive prestant that is fully loaded in an Alterra delivery system

Additional Equipment:

- Other compatible sheath:
 - Valve size: 29 mm, GORE DrySeal Flex Introducer Sheath (26 F, 65 cm)
- Balloon tip catheter
- Sizing balloons
- 20 cc syringe or larger (x2)
- 50 cc syringe or larger
- High-pressure 3-way stopcock
- Standard cardiac catheterization lab equipment and supplies, and access to standard heart valve operating room equipment and supplies
- Fluoroscopy (bi-plane, fixed, mobile, or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Exchange length 0.035 in (0.89 mm) stiff guidewire
- Sterile rinsing basins, physiological saline, heparinized saline, 15% diluted radiopaque contrast medium
- Sterile table for valve and accessories preparation

8.2 Edwards Alterra Adaptive Prestent System Procedure

See Edwards Alterra adaptive prestant system instructions for use for device preparation and implantation prior to transcatheter heart valve preparation and deployment.

Prior to valve implantation, assess Alterra prestant stability by evaluating apices engagement in surrounding tissue, wall apposition, and/or motion of prestant within the anatomy. Three radiopaque markers are provided at the waist of the prestant to help with positioning during fluoroscopy. If adequate stability is not noted, consider staging valve deployment after allowing sufficient time for prestant endothelialization.

CAUTION: Failure to identify prestant instability may lead to prestant migration/embolization when tracking interventional devices through the prestant.

8.3 Valve Handling and Preparation

Follow sterile technique during device preparation and implantation.

8.3.1 Valve Rinsing Procedure

Before opening the valve jar, carefully examine for evidence of damage (e.g., a cracked jar or lid, leakage, or broken or missing seals).

CAUTION: Valves from containers found to be damaged, leaking, without adequate sterilant, or missing intact seals must not be used for implantation.

1. Set up two (2) sterile bowls with at least 500 ml of sterile physiological saline to thoroughly rinse the glutaraldehyde sterilant from the valve.
2. Carefully remove the valve/holder assembly from the jar without touching the tissue. Verify the valve serial identification number with the number on the jar lid and record in the patient information documents. Inspect the valve for any signs of damage to the frame or tissue.
3. Rinse the valve as follows: Place the valve in the first bowl of sterile, physiological saline. Be sure the saline solution completely covers the valve and holder. With the valve and holder submerged, slowly agitate (to gently swirl the valve and holder) back and forth for a minimum of 1 minute. Transfer the valve and holder to the second rinsing bowl of physiological saline and gently agitate for at least one more minute. Ensure the rinse solution in the first bowl is not used. The valve should be left in the final rinse solution until needed to prevent the tissue from drying.

CAUTION: Do not allow the valve to come into contact with the bottom or sides of the rinse bowl during agitation or swirling in the rinse solution. Direct contact between the identification tag and valve is also to be avoided during the rinse procedure. No other objects should be placed in the rinse bowls. The valve should be kept hydrated to prevent the tissue from drying.

8.3.2 Prepare the System

Refer to the Edwards sheath, GORE DrySeal Flex Introducer Sheath, and Edwards crimper instructions for use for device preparation.

1. Visually inspect all the components for damage. Ensure the Edwards Commander delivery system is fully unflexed and the balloon catheter is fully advanced in the flex catheter.

WARNING: To prevent possible damage to the balloon shaft, ensure that the proximal end of the balloon shaft is not subjected to bending.

2. Flush the flex catheter.
3. Carefully remove the distal balloon cover from the delivery system.
4. Remove the stylet from the distal end of the guidewire lumen and set aside. Flush the guidewire lumen with heparinized saline and insert the stylet back into the distal end of the guidewire lumen.

CAUTION: Failure to replace the stylet in the guidewire lumen may result in damage to the lumen during the transcatheter heart valve crimping process.

5. Place the delivery system into the default position and make sure that the flex catheter tip is covered by the proximal balloon cover.
6. If using the Edwards provided sheath, unscrew the loader cap from the loader tube and flush the loader cap. Place the loader cap over the proximal balloon cover and onto the flex catheter with the inside of the cap oriented towards the distal tip.
If using the GORE DrySeal Flex Introducer Sheath, proceed to step 7.
7. Fully advance the balloon catheter in the flex catheter.
Peel off the proximal balloon cover over the blue section of the balloon shaft.
8. Attach a 3-way stopcock to the balloon inflation port. Fill a 50 cc or larger syringe with 15-20 ml of diluted contrast medium and attach to the 3-way stopcock.
9. Fill the inflation device provided by Edwards Lifesciences with excess volume of diluted contrast medium relative to the indicated inflation volume. Lock and attach to the 3-way stopcock.
10. Close 3-way stopcock to the inflation device provided by Edwards Lifesciences and de-air the system using the 50 cc or larger syringe. Slowly release the plunger and leave zero-pressure in the system.

WARNING: Ensure there is no residual fluid left in the balloon to avoid potential difficulty with valve alignment during the procedure.

11. Close the stopcock to the delivery system. By rotating the knob of the inflation device provided by Edwards Lifesciences, transfer the contrast medium into the syringe to achieve the appropriate volume required to deploy the valve, per the inflation parameters.
12. Close the stopcock to the 50 cc or larger syringe. Remove the syringe. Verify that the inflation volume is correct and lock the inflation device provided by Edwards Lifesciences.

CAUTION: Maintain the inflation device provided by Edwards Lifesciences in the locked position until valve deployment.

8.3.3 Mount and Crimp the Valve on the Delivery System

8.3.3.1 Procedure with the Edwards Provided Sheath

1. Set up two (2) additional sterile bowls with at least 100 ml of sterile physiological saline to thoroughly rinse the Qualcrimp crimping accessory.
2. Completely submerge the Qualcrimp crimping accessory in the first bowl and gently compress it to ensure complete saline absorption. Slowly swirl the Qualcrimp crimping accessory for a minimum of 1 minute. Repeat this process in the second bowl.
3. Remove the valve from the holder and remove the ID tag.
4. Attach the 2-piece crimp stopper to the base of the crimper and click into place.
5. With the crimper in the open position, gently place the valve into the crimper aperture. Gradually crimp the valve until it fits into the Qualcrimp crimping accessory.
6. Place the Qualcrimp crimping accessory over the valve making sure the valve is parallel to the edge of the Qualcrimp crimping accessory.
7. Place the valve and Qualcrimp crimping accessory in crimper aperture. Insert the delivery system coaxially within the valve on the valve crimp section (2-3 mm distal to the balloon shaft) with the orientation of the valve on the delivery system with the Inflow (outer skirt end) of the valve towards the proximal end of the delivery system.

Note: Verify correct valve orientation with the inflow (outer sealing skirt) oriented towards the handle.

8. Crimp the valve until it reaches the Qualcrimp stop located on the 2-piece crimp stopper.
9. Gently remove the Qualcrimp crimping accessory from the valve. Remove the Qualcrimp stop from the final stop, leaving the final stop in place.
10. Fully crimp the valve until it reaches the final stop for 5 seconds.

Note: Ensure that the valve crimp section remains coaxial within the valve.

11. Repeat the full crimp of the valve two more times (5 seconds each) for a total of three full crimps.
12. Pull the balloon shaft and lock in default position.
13. Flush the loader with heparinized saline. Immediately advance the valve into the loader until the tapered tip of the delivery system is exposed.

CAUTION: To prevent possible leaflet damage, the valve should not remain fully crimped and/or in the loader for over 15 minutes.

14. Attach the loader cap to the loader, re-flush the delivery system through the flush port and close the stopcock to the delivery system.
Remove the stylet and flush the guidewire lumen of the delivery system.

CAUTION: Keep valve hydrated until ready for implantation.

WARNING: The physician must verify correct orientation of the valve prior to its implantation.

8.3.3.2 Procedure with GORE DrySeal Flex Introducer Sheath

1. Set up two (2) additional sterile bowls with at least 100 ml of sterile physiological saline to thoroughly rinse the Qualcrimp crimping accessory.
2. Completely submerge the Qualcrimp crimping accessory in the first bowl and gently compress it to ensure complete saline absorption. Slowly swirl the Qualcrimp crimping accessory for a minimum of 1 minute. Repeat this process in the second bowl.
3. Remove the valve from the holder and remove the ID tag.
4. Attach the 2-piece crimp stopper to the base of the crimper and click into place.
5. With the crimper in the open position, gently place the valve into the crimper aperture. Gradually crimp the valve until it fits into the Qualcrimp crimping accessory.
6. Place the Qualcrimp crimping accessory over the valve making sure the valve is parallel to the edge of the Qualcrimp crimping accessory.
7. Place the valve and Qualcrimp crimping accessory in crimper aperture. Insert the delivery system coaxially within the valve on the valve crimp section (2-3 mm distal to the balloon shaft) with the orientation of the valve on the delivery system with the Inflow (outer skirt end) of the valve towards the proximal end of the delivery system.

Note: Verify correct valve orientation with the inflow (outer sealing skirt) oriented towards the handle.

8. Crimp the valve until it reaches the Qualcrimp stop located on the 2-piece crimp stopper.
9. Gently remove the Qualcrimp crimping accessory from the valve. Remove the Qualcrimp stop from the final stop, leaving the final stop in place.
10. Fully crimp the valve until it reaches the final stop for 5 seconds.

Note: Ensure that the valve crimp section remains coaxial within the valve.

11. Repeat the full crimp of the valve two more times (5 seconds each) for a total of three full crimps.
12. Pull the balloon shaft and lock in default position.
13. Flush the catheter with heparinized saline.

CAUTION: To prevent possible leaflet damage, the valve should not remain fully crimped for over 15 minutes.

14. Close the stopcock to the delivery system.

CAUTION: Keep valve hydrated until ready for implantation.

WARNING: The physician must verify correct orientation of the valve prior to its implantation.

15. Initiate valve alignment by disengaging the balloon lock and pulling the balloon catheter straight back until part of the warning marker is visible. Do not pull past the warning marker.

WARNING: To prevent possible damage to the balloon shaft, ensure that the proximal end of the balloon shaft is not subjected to bending.

16. Open the stopcock and flush the flex catheter using heparinized saline. Close the stopcock.
17. Engage the balloon lock.
18. Under fluoroscopy, utilize the fine adjustment wheel to position the valve between the valve alignment markers.

CAUTION: Do not turn the fine adjustment wheel if the balloon lock is not engaged.

WARNING: Do not position the valve past the distal valve alignment marker. This will prevent proper valve deployment.

19. Remove the stylet and flush the guidewire lumen of the delivery system.

8.4 Valve Delivery

Valve delivery should be performed under local and/or general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

CAUTION: Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

8.4.1 Procedure with the Edwards Sheath

1. If necessary, gain access using standard catheterization techniques, prepare and insert the sheath per its instructions for use.
2. If not present, insert the guidewire into the vasculature. Advance the guidewire into the intended landing zone per standard technique.
3. Insert the loader assembly into the sheath until the loader stops.
4. Advance the delivery system, with the Edwards logo in the proper orientation (the delivery system articulates in a direction opposite from the flush port), through the sheath until the valve exits the sheath. Retract the loader to the proximal end of the delivery system.

Note: The delivery system articulates in a direction opposite from the flush port.

CAUTION: The valve should not be advanced through the sheath if the sheath tip is not past the IVC bifurcation to minimize the risk of damage to the iliac vessel(s).

CAUTION: To prevent possible leaflet damage, the valve should not remain in the sheath for over 5 minutes.

5. In the vena cava, initiate valve alignment by disengaging the balloon lock and pulling the balloon catheter straight back until part of the warning marker is visible. Do not pull past the warning marker.

WARNING: To prevent possible damage to the balloon shaft, ensure that the proximal end of the balloon shaft is not subjected to bending.

Engage the balloon lock.

Utilize the fine adjustment wheel to position the valve between the valve alignment markers.

CAUTION: Do not turn the fine adjustment wheel if the balloon lock is not engaged.

WARNING: Do not position the valve past the distal valve alignment marker. This will prevent proper valve deployment.

CAUTION: Maintain guidewire position during valve alignment.

WARNING: If valve alignment is not performed in a straight section, there may be difficulties performing this step which may lead to delivery system damage and inability to inflate the balloon. Utilizing alternate fluoroscopic views may help with assessing curvature of the anatomy. If excessive tension is experienced during valve alignment, repositioning the delivery system to a different straight section of the vena cava and relieving compression (or tension) in the system will be necessary.

6. Advance the catheter and use the flex wheel, if needed, and cross the landing zone.

Note: Verify the orientation of the Edwards logo to ensure proper articulation. The delivery system articulates in a direction opposite from the flush port.

CAUTION: Use caution when advancing devices/delivery systems into the implanted Alterra adaptive prenent to avoid engagement with the inflow apices.

7. If additional working length is needed, remove the loader by unscrewing the loader cap and peeling the loader tubing from the delivery system.
8. Disengage the balloon lock and retract the tip of the flex catheter to the center of the triple marker. Engage the balloon lock.
9. Verify the correct position of the valve with respect to the landing zone.
10. As necessary, utilize the flex wheel to adjust the coaxial orientation of the valve and the fine adjustment wheel to adjust the position of the valve.
11. Before deployment, ensure that the valve is correctly positioned between the valve alignment markers and the flex catheter tip is locked over the triple marker.
12. Verify final position and begin valve deployment:
 - Unlock the inflation device provided by Edwards Lifesciences.
 - Using slow controlled inflation, deploy the valve by inflating the balloon with the entire volume in the inflation device, hold for 3 seconds and confirm that the barrel of the inflation device is empty to ensure complete inflation of the balloon.

- Deflate the balloon.

8.4.2 Procedure with GORE DrySeal Flex Introducer Sheath

1. If necessary, gain access using standard catheterization techniques, prepare and insert the sheath per its instructions for use.
2. If not present, insert the guidewire into the vasculature. Advance the guidewire into the intended landing zone per standard technique.
3. Insert the delivery system into the sheath.
4. Advance the delivery system, with the Edwards logo in the proper orientation (the delivery system articulates in a direction opposite from the flush port), through the sheath.

Note: The delivery system articulates in a direction opposite from the flush port.

CAUTION: The valve should not be advanced through the sheath if the sheath tip is not past the IVC bifurcation to minimize the risk of damage to the iliac vessel(s).

CAUTION: To prevent possible leaflet damage, the valve should not remain in the sheath for over 5 minutes.

5. Advance the catheter to the landing zone.

CAUTION: Use caution when advancing devices/delivery systems into the implanted Alterra adaptive prestant to avoid engagement with the inflow apices.

6. Expose the valve by retracting the GORE DrySeal Flex Introducer Sheath tip beyond the triple marker.
7. Disengage the balloon lock and retract the tip of the flex catheter to the center of the triple marker. Engage the balloon lock.
8. Verify the correct position of the valve with respect to the landing zone.
9. As necessary, utilize the flex wheel to adjust the coaxial orientation of the valve and the fine adjustment wheel to adjust the position of the valve.
10. Before deployment, ensure that the valve is correctly positioned between the valve alignment markers and the flex catheter tip is locked over the triple marker.
11. Begin valve deployment:
 - Unlock the inflation device provided by Edwards Lifesciences.
 - Using slow controlled inflation, deploy the valve by inflating the balloon with the entire volume in the inflation device, hold for 3 seconds and confirm that the barrel of the inflation device is empty to ensure complete inflation of the balloon.
 - Deflate the balloon.

8.5 System Removal

1. Once the balloon is fully deflated, unflex the delivery system. Verify that the flex catheter tip is locked over the triple marker.
If using the Edwards provided sheath, remove the delivery system from the sheath.
If using the GORE DrySeal Flex Introducer Sheath, retract the sheath and delivery system into the vena cava, then remove the delivery system from the sheath.

CAUTION: Patient injury could occur if the delivery system is not unflexed prior to removal.

2. Remove all devices when the ACT level is appropriate.

Refer to the Edwards sheath or the GORE DrySeal Flex Introducer Sheath instructions for use for device removal.

Note: A sheath or other device may need to be inserted per standard of care.

3. Close the access site.

9.0 How Supplied

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system, sheath, and crimper are supplied sterilized with ethylene oxide gas. The Edwards Alterra adaptive prestant system is supplied pouched and sterilized by e-beam sterilization.

The THV is supplied nonpyrogenic packaged in buffered glutaraldehyde, in a plastic jar to which a tamper evident seal has been applied. Each jar is shipped in a shelf box containing a temperature indicator to detect exposure of the THV to extreme temperature. The shelf box is enclosed in Styrofoam prior to shipping.

9.1 Storage

The transcatheter heart valve must be stored at 10 °C to 25 °C (50 °F to 77 °F). Each jar is shipped in an enclosure containing a temperature indicator to detect exposure of the THV to extreme temperature. The delivery system and accessories should be stored in a cool, dry place. The prestant and delivery system must be stored in a cool, dry place.

10.0 MR Safety

MR Conditional

Non-clinical testing has demonstrated that the Edwards Alterra adaptive prestant, alone or with a deployed SAPIEN 3 transcatheter heart valve, is MR Conditional. A patient can be scanned safely immediately after placement of this implant in an MR system meeting the following conditions:

- Static magnetic fields of 1.5 Tesla and 3.0 Tesla
- Spatial magnetic gradient field of 3000 Gauss/cm (30 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) scanning per sequence
- Gradient system is in normal operating mode

Under the scan conditions defined above, the Edwards Alterra adaptive prestant, alone or with a deployed SAPIEN 3 transcatheter heart valve, is expected to produce a maximum temperature rise of 4.0 °C or less after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 15 mm for gradient echo images when scanned using a 3.0 T MRI system. The artifact obscures the device lumen in spin and gradient echo images.

The delivery system has not been evaluated for MR compatibility and is considered MR unsafe.

11.0 Qualitative and Quantitative Information related to the Alterra Adaptive Presept System and SAPIEN 3 Transcatheter Heart Valve

The Alterra adaptive presept system and SAPIEN 3 transcatheter heart valve contain the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

For the Alterra adaptive presept, the following table shows the qualitative and quantitative information on the materials and substances:

Table 4

Substance	CAS	Model Mass Range (mg)
Nickel	7440-02-0	430 - 450
Titanium	7440-32-6	337 - 359
Polyethylene terephthalate	25038-59-9	146
Polyethylene	9002-88-4	27.5
Tantalum	7440-25-7	9.68 - 9.70
Titanium dioxide	13463-67-7	0.319 - 0.613
Iron	7439-89-6	0 - 0.396
Cobalt	7440-48-4	0 - 0.395
Oxygen	7782-44-7	0 - 0.317
Carbon	7440-44-0	0 - 0.317
Niobium	3/1/7440	0 - 0.207
Antimony trioxide	1309-64-4	0.176
Chromium	7440-47-3	0 - 0.0789
Copper	7440-50-8	0 - 0.0789
Nitrogen	7727-37-9	0 - 0.0404
Hydrogen	1333-74-0	0 - 0.0396
Tungsten	7440-33-7	0 - 0.00485
Molybdenum	7439-98-7	0 - 0.00194
Erucamide	112-84-5	0.00149 - 0.00152
Silicon	7440-21-3	0 - 0.000485
4-Dodecylbenzenesulfonic acid	121-65-3	0.000160

For the SAPIEN 3 transcatheter heart valve, the following table shows the qualitative and quantitative information on the materials and substances:

Table 5

Substance	CAS	Model Mass Range (mg)
Cobalt	7440-48-4	131 - 427
Nickel	7440-02-0	148 - 405
Chromium	7440-47-3	85.2 - 230
Polyethylene terephthalate	25038-59-9	102 - 170
Collagens, bovine, polymers with glutaraldehyde	2370819-60-4	58.3 - 141
Molybdenum	7439-98-7	40.3 - 115
Polytetrafluoroethylene	9002-84-0	17.5 - 25.5
Polyethylene	9002-88-4	14.2 - 19.7
Iron	7439-89-6	0 - 10.9
Titanium	7440-32-6	0 - 10.9
Manganese	7439-96-5	0 - 1.64
Silicon	7440-21-3	0 - 1.64
Titanium dioxide	13463-67-7	0.219 - 0.752
Polybutylate	24936-97-8	0.273 - 0.383
Carbon	7440-44-0	0 - 0.274
Antimony trioxide	1309-64-4	0.112 - 0.190
Boron	7440-42-8	0 - 0.164

Substance	CAS	Model Mass Range (mg)
Phosphorus	7723-14-0	0 - 0.164
Sulfur	7704-34-9	0 - 0.109
D&C Green No. 6	128-80-3	0.0394 - 0.0578
Silicon dioxide	7631-86-9	0.00422 - 0.00592
Erucamide	112-84-5	0.000683 - 0.00128
4-Dodecylbenzenesulfonic acid	121-65-3	0.000286 - 0.000430

12.0 Summary of Safety and Clinical Performance (SSCP)

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information.

The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the Alterra platform and the SAPIEN 3 platform.

Conformity with the Alterra platform and the SAPIEN 3 platform of the performance requirements (GSPR) for safety (MDR GSPR1), performance (MDR GSPR1), acceptability of side-effects (MDR GSPR8), usability (MDR GSPR5), device lifetime (MDR GSPR6), and acceptable benefit-risk profile (MDR GSPR8) has been established for the labelled indications.

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

13.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

Table 6

Product Name	Model	Basic UDI-DI
Edwards SAPIEN 3 transcatheter heart valve (29 mm)	9600TFX29	0690103D003SAP000VP
Edwards SAPIEN 3 Commander delivery system	9610TF29	0690103D00 3COM000TC
Edwards eSheath+ introducer set or Edwards eSheath introducer set	916ESP or 9610ES16	0690103D003S3E000NT
Edwards crimper	9600CR	0690103D003CRI000TH
Edwards Alterra adaptive prestant system	29AP4045	0690103D003AAP000ND
Inflation device	96406	0690103D003IND000TG

14.0 Expected Lifetime of the Device

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant has been subjected to rigorous pre-clinical durability testing per the testing requirements and in clinical studies and post market studies. The valve with prestant has been successfully tested to 5 years of simulated wear. In addition, clinical data show durability with follow-up to 2 years. The actual lifetime performance is continuing to be studied and varies from patient to patient.

15.0 Patient Information

Patient education brochures are provided to each site and should be given to the patient to inform them of the risks and benefits of the procedure and alternatives in adequate time before the procedure to be read and discussed with their physician. A copy of this brochure may also be obtained from Edwards Lifesciences by calling 1.800.822.9837. A patient implant card is provided with each prestant and THV. After implantation, please complete all requested information and provide the implant card to the patient. The serial number is found on the package. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

16.0 Recovered Valve, Prestent and Device Disposal

The explanted valve should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to the company. The explanted prestant should be placed into a suitable container and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit.

Used devices may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

17.0 Clinical Studies

Refer to the SSCP for clinical benefits.

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
	English	Español	Português
	Medical device	Producto sanitario	Dispositivo médico
	Contains biological material of animal origin	Contiene materiales biológicos de origen animal	Contém material biológico de origem animal
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Time & Temperature Sensitive	Sensible al tiempo y a la temperatura	Sensível à passagem do tempo e à temperatura
	Contains hazardous substances	Contiene sustancias peligrosas	Contém substâncias perigosas
	Size	Tamaño	Tamanho

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2024-04
10057933001 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU

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Edwards Lifesciences GmbH