Edwards

Edwards Alterra Adaptive Prestent System
Edwards Alterra adaptivt förstentsystem
Edwards Alterra adaptivt præ-stentsystem

Directory ■ Register ■ Register	
English (en).....	1
Svenska (sv).....	8
Dansk (da).....	14
Symbol Legend ■ Symbolförklaring ■ Symbolforklaring.....	20

English

Instructions for Use

Implantation of the Edwards Alterra adaptive prestent should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in balloon valvuloplasty.

1.0 Device Description

Edwards Alterra Adaptive Prestent System

The Edwards Alterra adaptive prestent system consists of an Alterra adaptive prestent that is fully loaded in an Alterra delivery system and supplied together in one package.

• Edwards Alterra Adaptive Prestent (Figure 1)

The Edwards Alterra adaptive prestent is used as a docking adaptor for the 29 mm Edwards SAPIEN 3 transcatheter heart valve (THV). It is comprised of a self-expanding, radiopaque, nitinol frame assembly with polyethylene terephthalate (PET) fabric covering. The prestent has designated inflow and outflow sides. The proximal inflow section is identifiable by the presence of two triangular tabs (prestent connector) that are attached to the catheter of the delivery system. The distal outflow section is distinguished by the open cells for blood flow. The PET fabric is attached by sutures to the inside surface of the frame to create sealing at the inflow section and opening for the outflow. Sutures are also used in the center to support the middle section when an Edwards SAPIEN 3 transcatheter heart valve is implanted. Three (3) fluoroscopically visible radiopaque markers at the prestent waist help with positioning.


Figure 1: Edwards Alterra Adaptive Prestent

Table 1

Inflow Sealing OD	Outflow Sealing OD	Height
39 mm	41 mm	49 mm

Sizing recommendation for implanting the prestant in the right ventricular outflow tract/pulmonary valve (RVOT/PV) landing zone are provided in Table 2:

Table 2: Prestent Sizing in RVOT landing zone

Perimeter	Perimeter Derived Diameter ¹	Prestent Size Diameter ² x Length	Valve Size
84.9 mm - 119.3 mm	27 mm - 38 mm	40 mm x 49 mm	29 mm

¹ Diameter range during systole

² Diameter is average of inflow and outflow diameters

Note: For Edwards SAPIEN 3 transcatheter heart valve implantation, refer to the Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant instructions for use.

• Edwards Alterra Delivery System (Figure 2)

The delivery system includes a handle which consists of a wheel that allows for deployment, two primary shafts with a flush port to flush the delivery system, and a long tapered tip at the distal end to facilitate tracking through the vasculature. A fluoroscopically visible radiopaque marker band shows the location of the tip of the outer shaft. The prestant is fully loaded in the delivery system. A stylet is included within the guidewire lumen.


Figure 2: Edwards Alterra Delivery System

Additional Accessories

• Edwards Sheath

Refer to the Edwards sheath instructions for use for device description.

2.0 Intended Use

The bioprosthesis with prestant is intended for use in patients requiring pulmonary heart valve replacement. The delivery systems and accessories are intended to facilitate placement of the bioprosthesis and prestant via the transfemoral access approach.

3.0 Indications

The 29 mm Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant is indicated for use in the management of patients with pulmonary regurgitation who have a native or surgically-repaired right ventricular outflow tract and are clinically indicated for pulmonary valve replacement.

4.0 Contraindications

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant is contraindicated in patients who:

- Cannot tolerate an anticoagulation/antiplatelet regimen or who have active bacterial endocarditis or other active infections.

5.0 Warnings

- The devices are designed, intended, and distributed STERILE for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- Assessment for coronary compression risk prior to implantation is essential to prevent the risk of severe patient harm.
- Patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, and/or polymeric materials may have an allergic reaction to these materials.
- Do not mishandle the delivery system or use the delivery system and accessory devices if the packaging sterile barriers and any components have been opened or damaged, cannot be flushed, or the expiration date has elapsed.
- The procedure should be conducted under fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.

6.0 Precautions

- Correct sizing of the prestant into the RVOT is essential to minimize risks such as paravalvular leak, migration, embolization, and/or RVOT rupture.
- Long-term durability has not been established for the prestant. Medical follow-up is advised so that device related complications can be diagnosed and properly managed.
- Patient venous anatomy should be evaluated to prevent the risk of access that would preclude the delivery and deployment of the device.
- Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.
- Patient should be heparinized to maintain the ACT at ≥ 250 sec prior to introduction of the delivery system in order to prevent thrombosis.
- Prestent recipients should be maintained on anticoagulant/antiplatelet therapy, except when contraindicated, as determined by their physician. This device has not been tested for use without antiplatelet therapy.
- It is recommended that all prestant recipients be prophylactically treated for endocarditis to minimize the risk of infection.
- If a prestant fracture is detected with significant loss in valve functionality, reintervention should be considered.
- Safety and effectiveness have not been established for patients with the following characteristics/comorbidities:
 - Blood dyscrasias defined as: leukopenia, acute anemia, thrombocytopenia, or history of bleeding diathesis or coagulopathy
 - A known hypersensitivity or contraindication to aspirin, heparin, ticlopidine (Ticlid[™]), or clopidogrel (Plavix[™]), or sensitivity to contrast media, which cannot be adequately premedicated
 - Positive urine or serum pregnancy test in female patients of child-bearing potential

7.0 Potential Adverse Events

Potential risks associated with the anesthesia, interventional procedure and imaging include but are not limited to:

- Death
- Stroke/transient ischemic attack
- Respiratory insufficiency or respiratory failure
- Cardiovascular or vascular injury, such as perforation or damage (dissection) of vessels, myocardium or valvular structures including rupture of the RVOT that may require intervention
- Pericardial effusion/cardiac tamponade
- Cardiac failure
- Embolic event: air, calcific material, thrombus
- Infection including incisional site infection, septicemia and endocarditis
- Myocardial infarction
- Renal insufficiency or renal failure
- Conduction system injury
- Arrhythmia
- Deep vein thrombosis
- Arteriovenous (AV) fistula
- Systemic or peripheral nerve injury
- Systemic or peripheral ischemia
- Pulmonary edema
- Pneumothorax
- Pleural effusion
- Dyspnea
- Atelectasis
- Dislodgement of previously implanted devices (i.e., pacing lead)
- Blood loss requiring transfusion
- Anemia
- Radiation injury
- Electrolyte imbalance
- Hypertension or hypotension
- Allergic reaction to anesthesia, contrast media, antithrombotic therapy, device materials
- Hematoma or ecchymosis
- Syncope
- Pain
- Exercise intolerance or weakness
- Inflammation
- Angina
- Fever

Potential risks that may or may not require intervention associated with the prestant, delivery system and/or accessories include, but may not be limited to, the following:

- Cardiac arrest
- Cardiogenic shock
- Coronary flow obstruction/transvalvular flow disturbance
- Device thrombosis
- Injury to tricuspid valve
- Device fracture
- Device embolization
- Device migration or malposition
- Endocarditis
- Chest pain/discomfort

- Device penetration/perforation into surrounding vasculature
- Device dysfunction (regurgitation and/or stenosis)
- Aortic root distortion
- Embolic event: device fragments
- Mechanical failure of delivery system, and/or accessories

For a patient/user/third party in the European Economic area; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and your national competent authority, which can be found at https://ec.europa.eu/growth/sectors/medical-devices/contacts_en.

8.0 Directions for Use

8.1 System Compatibility

Table 3

Product Name	Model
Edwards Alterra adaptive prestant system ^[1]	29AP4045
Edwards eSheath+ introducer set ^[2] or Edwards eSheath introducer set ^[2]	916ESP or 9610ES16

^[1] Includes an Alterra adaptive prestant that is fully loaded in an Alterra delivery system

^[2] Sheath provided by Edwards Lifesciences or equivalent

Additional Equipment:

- Balloon tip catheter
- Sizing balloons
- 20 cc syringe or larger
- 50 cc syringe or larger
- Standard cardiac catheterization lab equipment and supplies, and access to standard heart valve operating room equipment and supplies
- Fluoroscopy (bi-plane, fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Exchange length 0.035 in (0.89 mm) stiff guidewire
- Physiological saline
- Sterile table for device preparation

Refer to the Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant instructions for use for additional materials required to prepare the Edwards SAPIEN 3 transcatheter heart valve system.

8.2 Device Handling and Preparation

Follow sterile technique during device preparation and implantation.

8.2.1 Prepare the System

Refer to the Edwards sheath instructions for use for device preparation.

1. Remove the delivery catheter with the preloaded prestant from packaging. Visually inspect all components for damage.
2. Ensure there is a small gap between the outer shaft and the tapered tip to facilitate flushing of the inner lumen. If needed retract the outer shaft using the deployment wheel.

Note: Do not allow the end of the prestant to begin to exit the delivery system.

3. With the stylet still in place, flush the guidewire lumen with heparinized saline.
4. Using the deployment wheel, re-advance the outer shaft until it is even with the tapered tip.

Note: Do not overdrive the outer shaft onto the tapered tip.

5. Remove the stylet and repeat flushing the guidewire lumen.

8.3 Prestent Delivery

Prestent delivery should be performed under local and/or general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

Administer heparin to maintain the ACT at ≥ 250 sec.

CAUTION: Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

1. Gain access using standard catheterization techniques.
2. If necessary, predilate the vessel.
3. Introduce the sheath per its instructions for use.
4. Insert and advance the delivery system to the RVOT landing zone.

Note: Advance the delivery system from the shaft. Do not push the delivery system in by using the handle. Do not rotate the deployment wheel during advancement of the delivery system.

5. Position the delivery system fluoroscopically visible marker band distal to the intended landing zone.
6. Begin deployment by rotating the deployment wheel to retract the outer shaft.

Note: The delivery system marker band is located slightly proximal to the distal edge of the outer shaft.

Note: Waist markers on the prestant indicate the middle of the prestant.

Note: The prestant can be recaptured into the outer shaft and repositioned if deployed approximately 65%.

CAUTION: Once deployment has begun, do not reposition the device more distally. Advancement of the device with the prestant exposed may increase the risk for vascular damage.

- Continue deploying the Alterra pausing at approximately 30%, 50%, and 65% to assess for Alterra adaptive prestant positioning, coaxiality, and deployment angulation.
- If needed, recapture and reposition the prestant by rotating the deployment wheel in the reverse direction from deployment until the outer shaft fully covers the prestant as shown by the delivery system marker band.

Note: The prestant can be recaptured and redeployed one time. If a second recapture of the partially deployed prestant is performed, remove and replace the device.

Note: Several rotations of the deployment wheel may be necessary before the Alterra adaptive prestant begins to be recaptured.

CAUTION: Do not overdrive the outer shaft onto the tapered tip when recapturing the prestant. This may cause the delivery system to cinch down on the guidewire preventing independent movement of the delivery system and guidewire.

CAUTION: Recapturing and redeploying a prestant more than one time may impact implant integrity.

CAUTION: Recapturing a prestant that has been deployed more than 65% may cause system damage.

- After achieving an acceptable position, completely deploy the prestant by continuing to rotate the deployment wheel until the delivery system marker band is beyond the prestant connector.
- Confirm release of prestant.

CAUTION: Failure to identify release of the prestant connector tabs from the prestant connector may lead to prestant embolization during removal of the Alterra delivery system.

8.4 System Removal

- Slowly retract the system through the prestant. Remove the delivery system.

CAUTION: Ensure that the tapered tip and delivery system do not interfere with the prestant upon removal to prevent movement of the prestant.

- Assess Alterra prestant stability by evaluating apices engagement in surrounding tissue, wall apposition, and/or motion of prestant within the anatomy. If adequate stability is not noted, consider staging valve deployment after allowing sufficient time for prestant endothelialization.

CAUTION: Failure to identify prestant instability may lead to prestant migration/embolization when tracking interventional devices through the prestant.

9.0 How Supplied

STERILE: The Edwards Alterra adaptive prestant system is supplied pouched and sterilized by e-beam sterilization. The sheath is supplied sterilized with ethylene oxide gas.

9.1 Storage

The prestant and delivery system must be stored in a cool, dry place.

10.0 MR Safety

MR Conditional

Non-clinical testing has demonstrated that the Edwards Alterra adaptive prestant, alone or with a deployed SAPIEN 3 transcatheter heart valve, is MR Conditional. A patient can be scanned safely immediately after placement of this implant in an MR system meeting the following conditions:

- Static magnetic fields of 1.5 Tesla and 3.0 Tesla
- Spatial magnetic gradient field of 3000 Gauss/cm (30 T/m) or less
- Maximum MR system-reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) scanning per sequence
- Gradient system is in normal operating mode

Under the scan conditions defined above, the Edwards Alterra adaptive prestant, alone or with a deployed SAPIEN 3 transcatheter heart valve, is expected to produce a maximum temperature rise of 4.0 °C or less after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends as far as 15 mm for gradient echo images when scanned using a 3.0 T MRI system. The artifact obscures the device lumen in spin and gradient echo images.

The delivery system has not been evaluated for MR compatibility and is considered MR unsafe.

11.0 Qualitative and Quantitative Information related to the Alterra Prestent System

This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

For the Alterra adaptive prestant, the following table shows the qualitative and quantitative information on the materials and substances:

Table 4

Substance	CAS	Model Mass Range (mg)
Nickel	7440-02-0	430 - 450
Titanium	7440-32-6	337 - 359
Polyethylene terephthalate	25038-59-9	146
Polyethylene	9002-88-4	27.5

Substance	CAS	Model Mass Range (mg)
Tantalum	7440-25-7	9.68 - 9.70
Titanium dioxide	13463-67-7	0.319 - 0.613
Iron	7439-89-6	0 - 0.396
Cobalt	7440-48-4	0 - 0.395
Oxygen	7782-44-7	0 - 0.317
Carbon	7440-44-0	0 - 0.317
Niobium	3/1/7440	0 - 0.207
Antimony trioxide	1309-64-4	0.176
Chromium	7440-47-3	0 - 0.0789
Copper	7440-50-8	0 - 0.0789
Nitrogen	7727-37-9	0 - 0.0404
Hydrogen	1333-74-0	0 - 0.0396
Tungsten	7440-33-7	0 - 0.00485
Molybdenum	7439-98-7	0 - 0.00194
Erucamide	112-84-5	0.00149 - 0.00152
Silicon	7440-21-3	0 - 0.000485
4-Dodecylbenzenesulfonic acid	121-65-3	0.000160

12.0 Summary of Safety and Clinical Performance (SSCP)

The SSCP has been adapted in accordance with the clinical evaluation assessment by the Notified Body on which CE certification has been granted. The SSCP contains a relevant summary of the same information.

The Notified Body has taken notice of and agreed with the benefit-risk rationales for the short- and long-term safety and effectiveness of the Alterra platform.

Conformity with the Alterra platform of the performance requirements (GSPR) for safety (MDR GSPR1), performance (MDR GSPR1), acceptability of side-effects (MDR GSPR8), usability (MDR GSPR5), device lifetime (MDR GSPR6), and acceptable benefit-risk profile (MDR GSPR8) has been established for the labelled indications.

Refer to <https://meddeviceinfo.edwards.com/> for a SSCP for this medical device.

After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

13.0 Basic Unique Device Identification-Device Identifier (UDI-DI)

The Basic UDI-DI is the access key for device-related information entered in the Eudamed. The Basic UDI-DI for the prestant and sheath can be used to locate the SSCP.

The following tables contain the Basic UDI-DIs:

Table 5

Product	Model	Basic UDI-DI
Edwards Alterra Adaptive Prestent	29AP4045	0690103D003AAP000ND
Edwards eSheath+ introducer set or Edwards eSheath introducer set	916ESP or 9610ES16	0690103D003S3E000NT

14.0 Expected Lifetime of the Device

The Edwards SAPIEN 3 transcatheter pulmonary valve system with Alterra adaptive prestant has been subjected to rigorous pre-clinical durability testing per the testing requirements and in clinical studies and post market studies. The valve with prestant has been successfully tested to 5 years of simulated wear. In addition, clinical data show durability with follow-up to 2 years. The actual lifetime performance is continuing to be studied and varies from patient to patient.

15.0 Patient Information

Patient education brochures are provided to each site and should be given to the patient to inform them of the risks and benefits of the procedure and alternatives in adequate time before the procedure to be read and discussed with their physician. A copy of this brochure may also be obtained from Edwards Lifesciences by calling 1.800.822.9837.

A patient implant card is provided with each prestant. After implantation, please complete all requested information and provide the implant card to the patient. The serial number is found on the package. This implant card allows patients to inform healthcare providers what type of implant they have when they seek care.

16.0 Recovered Prestent and Device Disposal

The explanted prestant should be placed into a suitable container and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an explant kit.

Used devices may be handled and disposed of in the same manner as hospital waste and biohazardous materials. There are no special risks related to the disposal of these devices.

17.0 Clinical Studies

Refer to the SSCP for clinical benefits.

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
	English	Svenska	Dansk
	Contains biological material of animal origin	Innehåller biologiskt material av animaliskt ursprung	Indeholder biologisk materiale af animalsk oprindelse
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Time & Temperature Sensitive	Tids- och temperaturkänslig	Følsom overfor tid og temperatur
	Contains hazardous substances	Innehåller farliga ämnen	Indeholder farlige stoffer
SZ	Size	Storlek	Størrelse

Note: Not all symbols may be included in the labeling of this product. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne.

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2024-05
10057977001 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1 949.250.2525

Web IFU

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Edwards Lifesciences GmbH