



Edwards

Edwards SAPIEN 3 and SAPIEN 3 Ultra Transcatheter Heart Valve System

SAPIEN 3 and SAPIEN 3 Ultra Transcatheter Heart Valve Edwards Certitude Delivery System

Instructions for Use

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

Implantation of the transcatheter heart valve should be performed only by physicians who have received Edwards Lifesciences training. The implanting physician should be experienced in balloon aortic valvuloplasty.

Please verify that you have the latest version of the instructions for use prior to using the device by visiting <http://eIFU.edwards.com> or by calling 1.800.822.9837. In order to access the instructions for use, an IFU Code will be required.

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system, introducer sheath set, and crimper are supplied sterilized with ethylene oxide gas.

1.0 Device Description

Edwards SAPIEN 3 and SAPIEN 3 Ultra Transcatheter Heart Valve (THV) System

The Edwards SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valve system consists of the Edwards SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valves and delivery systems.

- **Edwards SAPIEN 3 Ultra Transcatheter Heart Valve - (Figure 1)**

The Edwards SAPIEN 3 Ultra transcatheter heart valve is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) inner and outer fabric skirts. The leaflets are treated according to the Carpentier-Edwards ThermaFix process.

- **Edwards SAPIEN 3 Transcatheter Heart Valve - (Figure 2)**

The Edwards SAPIEN 3 transcatheter heart valve is comprised of a balloon-expandable, radiopaque, cobalt-chromium frame, trileaflet bovine pericardial tissue valve, and polyethylene terephthalate (PET) fabric skirt. The leaflets are treated according to the Carpentier-Edwards ThermaFix process.




Figure 1: Edwards SAPIEN 3 Ultra Transcatheter Heart Valve

9750TFX

Table 1

Valve Size	Valve Height
20 mm	15.5 mm
23 mm	18 mm
26 mm	20 mm




Figure 2: Edwards SAPIEN 3 Transcatheter Heart Valve

9600TFX

Table 2

Valve Size	Valve Height
20 mm	15.5 mm
23 mm	18 mm
26 mm	20 mm
29 mm	22.5 mm

Sizing recommendations for implanting the Edwards SAPIEN 3 Ultra transcatheter heart valve and the Edwards SAPIEN 3 transcatheter heart valve in a native annulus are provided in the table below:

Table 3

Native Valve Annulus Size (TEE)	Native Valve Annulus Size (CT)		THV Size
	Area	Area Derived Diameter	
16 - 19 mm	273 - 345 mm ²	18.6 - 21 mm	20 mm
18 - 22 mm	338 - 430 mm ²	20.7 - 23.4 mm	23 mm
21 - 25 mm	430 - 546 mm ²	23.4 - 26.4 mm	26 mm
24 - 28 mm	540 - 683 mm ²	26.2 - 29.5 mm	29 mm

Valve size recommendations are based on native valve annulus size, as measured by transesophageal echocardiography (TEE) or computed tomography (CT). Patient anatomical factors and multiple imaging modalities should be considered during valve size selection.

Note: Risks associated with undersizing and oversizing should be considered.

Sizing recommendations for implanting the Edwards SAPIEN 3 Ultra transcatheter heart valve and the Edwards SAPIEN 3 transcatheter heart valve in a failing surgical bioprosthesis, except for the INSPIRIS RESILIA aortic valve sizes 19 - 25 mm, are provided in the table below:

Table 4

Surgical Valve True Inner Diameter (ID) ^[1]	THV Size
16.5 - 19.0 mm	20 mm
18.5 - 22.0 mm	23 mm
22.0 - 25.0 mm	26 mm
25.0 - 28.5 mm	29 mm

Note: Surgical valve 'True ID' may be smaller than the labeled valve size. For THV-in-THV, the native valve annulus size should be considered to determine the appropriate THV size to implant. For a failing stentless bioprosthesis, consider sizing recommendations for a native annulus. The dimensions of the failed bioprosthesis should be determined so that the appropriate THV size can be implanted; and is best determined by using computed tomography, magnetic resonance imaging, and/or transesophageal echocardiography.

Sizing recommendations for implanting the Edwards SAPIEN 3 Ultra transcatheter heart valve and the SAPIEN 3 transcatheter heart valve in a failing INSPIRIS RESILIA aortic surgical bioprosthesis in sizes 19 - 25 mm are provided in the table below:

Table 5

INSPIRIS RESILIA Aortic Valve (model 11500A)* Labeled Size	THV Size
19 mm	20 mm or 23 mm
21 mm	23 mm or 26 mm
23 mm	23 mm or 26 mm
25 mm	26 mm or 29 mm

*INSPIRIS RESILIA aortic valve model 11500A sizes 19-25 mm incorporate VFit technology which consists of expandable bands and fluoroscopically visible size markers designed for potential future valve-in-valve procedures. Clinical data are not currently available on the INSPIRIS RESILIA aortic valve Model 11500A valve-in-valve procedure or expansion feature. The impact of tissue ingrowth on the expansion feature of the INSPIRIS RESILIA aortic valve has not been assessed.

WARNING: Do not perform stand-alone balloon aortic valvuloplasty procedures in the INSPIRIS RESILIA aortic valve for the sizes 19 - 25 mm. This may expand the valve causing aortic incompetence, coronary embolism or annular rupture.

Note: INSPIRIS RESILIA aortic valve model 11500A sizes 27 - 29 mm do not incorporate VFit technology and therefore follow the surgical valve True ID sizing provided in Table 4.

Sizing recommendations for implanting the Edwards SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valves in a failing native mitral valve with an annuloplasty ring are provided in the table below:

Table 6

CT Orifice Area		SAPIEN 3 and SAPIEN 3 Ultra valve
Min Value	Max Value	
280 mm ²	350 mm ²	23 mm
350 mm ²	450 mm ²	26 mm
450 mm ²	600 mm ²	29 mm

Note that the size 20 mm SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valve is not included due to its size being too small for the mitral position.

WARNING: Transcatheter valve replacement in mitral annuloplasty rings is not recommended in cases of partial annuloplasty ring dehiscence due to high risk of PVL.

WARNING: Transcatheter valve replacement in mitral annuloplasty rings is not recommended in cases of partial (incomplete) annuloplasty rings in the absence of annular calcium due to increased risk of valve embolization.

WARNING: Transcatheter valve replacement in mitral annuloplasty rings is not recommended in cases of rigid annuloplasty rings due to increased risk of PVL or THV deformation.

Note: Exact volume required to deploy the THV may vary depending on the prosthesis inner diameter. Factors such as calcification and pannus tissue growth may not be accurately visualized in imaging and may reduce the effective inner diameter of the failing prosthesis to a size smaller than the 'True ID'. These factors should be considered and assessed in order to determine the most appropriate THV size to achieve nominal THV deployment and sufficient anchoring. Do not exceed the rated burst pressure. See Table 7 for inflation parameters.

• Edwards Certitude Delivery System (Figure 3)

The Edwards Certitude delivery system facilitates the placement of the bioprosthesis. The delivery system consists of a flex catheter to aid in tracking and valve positioning. The delivery system includes a tapered tip to facilitate crossing of the valve. The handle contains a flex wheel to control flexing of the balloon catheter. A stylet is included within the guidewire lumen of the

delivery system. A radiopaque center marker in the balloon is provided to assist with valve positioning. The inflation parameters for valve deployment are:

Table 7

Model	Nominal Balloon Diameter	Nominal Inflation Volume	Rated Burst Pressure (RBP)
9630TA20 9600SDS20	20 mm	12 ml	7 atm
9630TA23 9600SDS23	23 mm	17 ml	7 atm
9630TA26 9600SDS26	26 mm	23 ml	7 atm
9630TA29 9600SDS29	29 mm	30 ml	7 atm




Figure 3: Edwards Certitude Delivery System

Additional Accessories




Figure 4: Qualcrimp Crimping Accessory




Figure 5: Loader






Figure 6: Extension Tubing

• Edwards Certitude Introducer Sheath Set (Figure 7)

The Edwards Certitude introducer sheath set facilitates the introduction and removal of devices utilized with the SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valves. The sheath has a radiopaque marker for visualization of the sheath tip and non-radiopaque depth markings on the distal end of the body of the sheath. The proximal end of the sheath includes a flush tube and three hemostasis valves. An introducer is supplied with the sheath. The entire introducer is radiopaque.



1. Housing
2. Flush Tube with Stopcock
3. Sheath
4. Radiopaque Marker
5. Non-Radiopaque Depth Markers
6. Introducer

Figure 7: Edwards Certitude Introducer Sheath Set

• Edwards Crimper

Refer to the Edwards Crimper instructions for use for device description.

2.0 Indications

1. The Edwards SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valve system is indicated for relief of aortic stenosis in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a Heart Team, including a cardiac surgeon, to be appropriate for the transcatheter heart valve replacement therapy.
2. The Edwards SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valve system is indicated for patients with symptomatic heart disease due to failing (stenosed, insufficient, or combined) of a surgical or transcatheter bioprosthetic aortic valve, a surgical bioprosthetic mitral valve, or a native mitral valve with an annuloplasty ring who are judged by a heart team, including a cardiac surgeon, to be at high or greater risk for open surgical therapy (i.e., predicted risk of surgical mortality $\geq 8\%$ at 30 days, based on the Society of Thoracic Surgeons (STS) risk score and other clinical co-morbidities unmeasured by the STS risk calculator).

3.0 Contraindications

The Edwards SAPIEN 3 and SAPIEN 3 Ultra transcatheter heart valve system are contraindicated in patients who cannot tolerate an anticoagulation/antiplatelet regimen, who have active bacterial endocarditis or other active infections, or who has significant annuloplasty ring dehiscence.

4.0 Warnings

- Observation of the pacing lead throughout the procedure is essential to avoid the potential risk of pacing lead perforation.
- There may be an increased risk of stroke in transcatheter aortic valve replacement procedures, as compared to balloon aortic valvuloplasty or other standard treatments in high or greater risk patients.
- The devices are designed, intended, and distributed for single use only. **Do not resterilize or reuse the devices.** There are no data to support the sterility, nonpyrogenicity, and functionality of the devices after reprocessing.
- Incorrect sizing of the valve may lead to paravalvular leak, migration, embolization, residual gradient (patient-prosthesis mismatch) and/or annular rupture.
- Accelerated deterioration of the valve may occur in patients with an altered calcium metabolism.
- Prior to delivery, the valve must remain hydrated at all times and cannot be exposed to solutions other than its shipping storage solution and sterile physiologic rinsing solution. Valve leaflets mishandled or damaged during any part of the procedure will require replacement of the valve.
- Caution should be exercised in implanting a valve in patients with clinically significant coronary artery disease.
- Patients with pre-existing prostheses should be carefully assessed prior to implantation of the valve to ensure proper valve positioning and deployment.
- Do not use the valve if the tamper evident seal is broken, the storage solution does not completely cover the valve, the temperature indicator has been activated, the valve is damaged, or the expiration date has elapsed.
- Do not mishandle the delivery system or use it if the packaging or any components are not sterile, have been opened or are damaged (e.g., kinked or stretched), or the expiration date has elapsed.
- Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.
- Care should be exercised in patients with hypersensitivities to cobalt, nickel, chromium, molybdenum, titanium, manganese, silicon, and/or polymeric materials.
- The procedure should be conducted under fluoroscopic guidance. Some fluoroscopically guided procedures are associated with a risk of radiation injury to the skin. These injuries may be painful, disfiguring, and long-lasting.
- Valve recipients should be maintained on anticoagulant/antiplatelet therapy, except when contraindicated, as determined by their physician. This device has not been tested for use without anticoagulation.
- Do not add or apply antibiotics to the storage solution, rinse solutions, or to the valve.
- Balloon valvuloplasty should be avoided in the treatment of failing bioprostheses as this may result in embolization of bioprosthesis material and mechanical disruption of the valve leaflets.

5.0 Precautions

- Long-term durability has not been established for the valve. Regular medical follow-up is advised to evaluate valve performance.
- Limited clinical data are available for transcatheter aortic valve replacement in patients with a congenital bicuspid valve who are deemed to be low surgical risk. Anatomical characteristics should be considered when using this valve in this population. In addition, patient age should be considered as long-term durability for this valve has not been established.
- Glutaraldehyde may cause irritation of the skin, eyes, nose and throat. Avoid prolonged or repeated exposure to, or breathing of, the solution. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with water; in the event of contact with eyes, seek immediate medical attention. For more information about glutaraldehyde exposure, refer to the Material Safety Data Sheet available from Edwards Lifesciences.
- To maintain proper valve leaflet coaptation, do not overinflate the deployment balloon.
- Appropriate antibiotic prophylaxis is recommended post-procedure in patients at risk for prosthetic valve infection and endocarditis.
- Special care must be exercised in mitral valve replacement to avoid entrapment of the subvalvular apparatus.
- Safety and effectiveness have not been established for patients with the following characteristics/comorbidities:
 - Non-calcified aortic annulus
 - Severe ventricular dysfunction with ejection fraction < 20%
 - Congenital unicuspid aortic valve
 - Pre-existing prosthetic ring in the tricuspid position
 - Severe mitral annular calcification (MAC), severe (> 3+) mitral insufficiency, or Gorlin syndrome
 - Blood dyscrasias defined as: leukopenia (WBC < 3000 cells/ml), acute anemia (Hb < 9 g/dL), thrombocytopenia (platelet count < 50,000 cells/ml), or history of bleeding diathesis or coagulopathy
 - Hypertrophic cardiomyopathy with or without obstruction (HOCM)

- Echocardiographic evidence of intracardiac mass, thrombus, or vegetation
- A known hypersensitivity or contraindication to aspirin, heparin, ticlopidine (Ticlid™), or clopidogrel (Plavix™), or sensitivity to contrast media, which cannot be adequately premedicated
- Excessive calcification at access site
- Bulky calcified aortic valve leaflets in close proximity to coronary ostia
- A concomitant paravalvular leak where the failing prosthesis is not securely fixed in the native annulus or is not structurally intact (e.g., wireform frame fracture, annuloplasty ring dehiscence)
- A partially detached leaflet of the failing bioprosthesis that in the aortic position may obstruct a coronary ostium
- Residual mean gradient may be higher in a "THV-in-failing prosthesis" configuration than that observed following implantation of the valve inside a native aortic annulus using the same size device. Patients with elevated mean gradient post procedure should be carefully followed. It is important that the manufacturer, model and size of the preexisting prosthesis be determined, so that the appropriate valve can be implanted and a prosthesis-patient mismatch be avoided. Additionally, pre-procedure imaging modalities must be employed to make as accurate a determination of the inner diameter as possible.
- Post-procedure and follow-up assessment of TAVR device performance by Doppler echocardiography may be impacted by inherent limitations in the Bernoulli equation used to determine measurements such as mean gradient, EOA, and prosthesis-patient mismatch. These limitations may lead to an overstating or understating of valve performance measurements after TAVR implantation. Therefore, a post-TAVR echocardiogram should be used to establish a baseline from which future follow-up visits are compared to. Confirmatory direct pressure measurement via cardiac catheterization may be considered, when indicated, prior to reintervention.

6.0 Potential Adverse Events

Potential risks associated with the overall procedure including potential access complications associated with standard cardiac catheterization, balloon valvuloplasty, the potential risks of conscious sedation and/or general anesthesia, and the use of angiography:

- Death
- Stroke/transient ischemic attack, clusters or neurological deficit
- Paralysis
- Permanent disability
- Respiratory insufficiency or respiratory failure
- Hemorrhage requiring transfusion or intervention
- Cardiovascular injury including perforation or dissection of vessels, ventricle, atrium, myocardium or valvular structures that may require intervention
- Pericardial effusion or cardiac tamponade
- Embolization including air, calcific valve material or thrombus
- Infection including septicemia and endocarditis
- Heart failure
- Myocardial infarction
- Renal insufficiency or renal failure
- Conduction system defect which may require a permanent pacemaker
- Arrhythmia
- Retroperitoneal bleed
- Arteriovenous (AV) fistula or pseudoaneurysm
- Reoperation
- Ischemia or nerve injury
- Restenosis
- Pulmonary edema
- Pleural effusion
- Bleeding
- Anemia
- Abnormal lab values (including electrolyte imbalance)
- Hypertension or hypotension
- Allergic reaction to anesthesia, contrast media, or device materials
- Hematoma
- Syncope

- Pain or changes (e.g., wound infection, hematoma, and other wound care complications) at the access site
- Exercise intolerance or weakness
- Inflammation
- Angina
- Heart murmur
- Fever

Additional potential risks associated with the use of the valve, delivery system, and/or accessories include:

- Cardiac arrest
- Cardiogenic shock
- Emergency cardiac surgery
- Cardiac failure or low cardiac output
- Coronary flow obstruction/transvalvular flow disturbance
- Device thrombosis requiring intervention
- Valve thrombosis
- Device embolization
- Device migration or malposition requiring intervention
- Left ventricular outflow tract obstruction
- Valve deployment in unintended location
- Valve stenosis
- Structural valve deterioration (wear, fracture, calcification, leaflet tear/tearing from the stent posts, leaflet retraction, suture line disruption of components of a prosthetic valve, thickening, stenosis)
- Device degeneration
- Paravalvular or transvalvular leak
- Valve regurgitation
- Hemolysis
- Injury to the mitral valve
- Device explants
- Mediastinitis
- Mediastinal bleeding
- Nonstructural dysfunction
- Mechanical failure of delivery system, and/or accessories
- Non-emergent reoperation

7.0 Directions for Use

7.1 System Compatibility

Table 8

Product Name	20 mm System	23 mm System	26 mm System
	Model		
Edwards SAPIEN 3 Ultra Transcatheter Heart Valve	9750TFX20	9750TFX23	9750TFX26
Edwards Certitude delivery system	9630TA20	9630TA23	9630TA26
Edwards Certitude introducer sheath set	9600IS18	9600IS21	
Inflation device, Qualcrimp crimping accessory, 2-piece crimp stopper, loader and extension tubing provided by Edwards Lifesciences			
Edwards Crimper	9600CR		

Table 9

Product Name	20 mm System	23 mm System	26 mm System	29 mm System
	Model			
Edwards SAPIEN 3 Transcatheter Heart Valve	9600TFX20	9600TFX23	9600TFX26	9600TFX29
Edwards Certitude delivery system	9600SDS20	9600SDS23	9600SDS26	9600SDS29
Edwards Certitude introducer sheath set	9600IS18			9600IS21
Inflation device, Qualcrimp crimping accessory, 2-piece crimp stopper, loader and extension tubing provided by Edwards Lifesciences				
Edwards Crimper	9600CR			

Additional Equipment:

- Balloon catheter, per the discretion of the physician
- 20 cc syringe or larger (x2)
- 50 cc syringe or larger
- Standard cardiac catheterization lab equipment
- Fluoroscopy (fixed, mobile or semi-mobile fluoroscopy systems appropriate for use in percutaneous coronary interventions)
- Transesophageal or transthoracic echocardiography capabilities
- 18 gauge Seldinger needle (for transaortic)
- 145 cm x 0.035 inch (0.89 mm) soft guidewire
- 180 cm or 260 cm x 0.035 inch (0.89 mm) & Exchange length 0.035 inch (0.89 mm) extra-stiff guidewire
- Temporary pacemaker (PM) and pacing lead
- Sterile rinsing basins, physiological saline, heparinized saline, 15% diluted radiopaque contrast medium
- Sterile table for valve and device preparation

7.2 Valve Handling and Preparation

Follow sterile technique during device preparation and implantation.

7.2.1 Valve Rinsing Procedure

Before opening the valve jar, carefully examine for evidence of damage (e.g., a cracked jar or lid, leakage, or broken or missing seals).

CAUTION: Valves from containers found to be damaged, leaking, without adequate sterilant, or missing intact seals must not be used for implantation.

1. Set up two (2) sterile bowls with at least 500 ml of sterile physiological saline to thoroughly rinse the glutaraldehyde sterilant from the valve.
2. Carefully remove the valve/holder assembly from the jar without touching the tissue. Verify the valve serial identification number with the number on the jar lid and record in the patient information documents. Inspect the valve for any signs of damage to the frame or tissue.
3. Rinse the valve as follows: Place the valve in the first bowl of sterile, physiological saline. Be sure the saline solution completely covers the valve and holder. With the valve and holder submerged, slowly agitate (to gently swirl the valve and holder) back and forth for a minimum of 1 minute. Transfer the valve and holder to the second rinsing bowl of sterile, physiological saline and gently agitate for at least one more minute. Ensure the rinse solution in the first bowl is not used. The valve should be left in the final rinse solution until needed to prevent the tissue from drying.

CAUTION: Do not allow the valve to come into contact with the bottom or sides of the rinse bowl during agitation or swirling in the rinse solution. Direct contact between the identification tag and valve is also to be avoided during the rinse procedure. No other objects should be placed in the rinse bowls. The valve should be kept hydrated to prevent the tissue from drying.

7.2.2 Prepare the Components

Refer to the Edwards Crimper instructions for use for device preparation.

1. Visually inspect all components for damage. Ensure the Edwards Certitude delivery system is fully unflexed.
2. Prime and flush the introducer and sheath with heparinized saline. Hydrate the length of the introducer and sheath.

3. Advance the introducer fully into the sheath housing.
4. Unscrew the loader cap from the loader and flush the loader cap with heparinized saline.
5. Place the loader cap onto the delivery system with the inside of the cap oriented towards the tapered tip.
6. Flush the extension tubing and connect to the delivery system.
7. Partially fill a 50 ml or larger syringe with diluted contrast medium, and connect to the extension tubing.
8. Fill the inflation device with 20 ml of diluted contrast medium, lock the inflation device, and connect to the extension tubing. Close 3-way stopcock to inflation device.
9. De-air the delivery system using the luer lock syringe. Leave zero-pressure in the system. Close the 3-way stopcock to the luer lock syringe.
10. Remove 3 ml fluid from the delivery system by turning the knob of the locked inflation device. Keep the inflation device locked for valve crimping steps.

7.2.3 Mount and Crimp the Valve on the Delivery System

1. Set up two (2) additional sterile bowls with at least 100 ml of sterile, physiological saline to thoroughly rinse the Qualcrimp crimping accessory.
2. Completely submerge the Qualcrimp crimping accessory in the first bowl and gently compress it to ensure complete saline absorption. Slowly swirl the Qualcrimp crimping accessory for a minimum of 1 minute. Repeat this process in the second bowl.
3. Remove the valve from the holder and remove the ID tag.
4. Attach the 2-piece crimp stopper to the base of the crimper and click into place.
5. With the crimper in the open position, gently place the valve into the crimper aperture. Gradually crimp the valve (if necessary) until it fits into the Qualcrimp crimping accessory.
6. Place the Qualcrimp crimping accessory over the valve making sure the valve is parallel to the edge of the Qualcrimp crimping accessory.
7. Place the valve and Qualcrimp crimping accessory in the crimper aperture. Insert the delivery system coaxially within the valve. The orientation of the valve on the delivery system is described below:

Antegrade approach: Inflow (outer skirt end) of the valve towards the proximal end of the delivery system.

Retrograde approach: Inflow (outer skirt end) of the valve towards the distal end of the delivery system.

8. Crimp the valve between the two internal shoulders of the delivery system until it reaches the Qualcrimp crimper Stop located on the 2-piece Crimp Stopper.
9. Gently remove the Qualcrimp crimping accessory from the valve. Remove the Qualcrimp crimper Stop from the Final Stop, leaving the Final Stop in place.
10. Fully crimp the valve until it reaches the Final Stop.

Note: Ensure that the delivery system remains coaxial within the valve.

11. Repeat the full crimp of the valve two more times for a total of three full crimps.
12. Flush the loader with heparinized saline. Immediately advance the valve into the loader until the tapered tip of the delivery system is exposed and the valve is within the distal end of the loader tube.

CAUTION: To prevent possible leaflet damage, the valve should not remain fully crimped and/or in the loader for over 15 minutes.

13. Attach the loader cap to the loader and flush through the port on the loader.

CAUTION: Keep the valve hydrated until ready for implantation.

CAUTION: The physician must verify correct orientation of the valve prior to its implantation.

14. With 3-way stopcock still closed to the luer lock syringe, unlock the inflation device. Allow the delivery system to reach zero-pressure.
15. Close the 3-way stopcock to the delivery system. Use the luer lock syringe to de-air the inflation device if necessary.
16. Adjust the inflation device to the inflation volume required to deploy the valve, per the following:

Table 10

Delivery System	Valve	Inflation Volume
Model 9630TA20 9600SDS20	20 mm	12 ml

Delivery System	Valve	Inflation Volume
Model 9630TA23 9600SDS23	23 mm	17 ml
Model 9630TA26 9600SDS26	26 mm	23 ml
Model 9630TA29 9600SDS29	29 mm	30 ml

Re-lock the inflation device. Close the 3-way stopcock to the luer lock syringe and remove syringe.

CAUTION: Maintain the inflation device in a locked position until valve deployment.

17. Remove the stylet and flush the guidewire lumen of the delivery system.

7.3 Valvuloplasty and Valve Delivery

Valvuloplasty and valve delivery should be performed under conscious sedation and/or general anesthesia with hemodynamic monitoring in a catheterization lab/hybrid operating room with fluoroscopic and echocardiographic imaging capabilities.

The following table shows the minimum required distances from the valvular plane to the distal tip of the sheath to allow the Edwards Certitude delivery system balloon to inflate properly during valve deployment. These distances do not include sheath insertion depth, which should be considered during the transaortic approach when selecting the access site on the ascending aorta.

Table 11

Delivery System	Valve	Minimum Required Distance From Sheath Tip to Valvular Plane
9630TA20 9600SDS20	20 mm	3.5 cm
9630TA23 9600SDS23	23 mm	3.5 cm
9630TA26 9600SDS26	26 mm	3.5 cm
9630TA29 9600SDS29	29 mm	4.0 cm

Administer heparin to maintain the ACT at ≥ 250 sec during the procedure.

Balloon valvuloplasty should be avoided in the treatment of failing bioprostheses as this may result in embolization of bioprosthesis material and mechanical disruption of the valve leaflets.

CAUTION: Use of excessive contrast media may lead to renal failure. Measure the patient's creatinine level prior to the procedure. Contrast media usage should be monitored.

CAUTION: Care should be taken to avoid damage to soft tissue, chordae, aorta, native leaflet or ventricular wall during insertion, positioning and removal of devices.

7.3.1 Baseline Parameters

1. Perform an angiogram with fluoroscopic view perpendicular to the valve.
2. Evaluate the distance of the left and right coronary ostia from the aortic annulus in relation to the valve frame height.
3. Introduce a pacemaker (PM) lead and position appropriately.
4. Set the stimulation parameters to obtain 1:1 capture, and test pacing.

7.3.2 Access

1. Gain access using standard catheterization techniques.
2. Using the sheath depth markers, advance the introducer and sheath over the guidewire to the desired depth while following its progression on fluoroscopy.
3. Withdraw the introducer slowly, keeping the sheath in place. Maintain guidewire position across the valve.

7.3.3 Valvuloplasty

Pre-dilate the native aortic valve, per the discretion of the physician, according to the instructions for use for the selected balloon aortic valvuloplasty catheter.

CAUTION: Valve implantation should not be carried out if the balloon cannot be fully inflated during valvuloplasty.

7.3.4 Valve Delivery

1. Confirm that the valve is oriented properly and the volume in the inflation device matches the indicated volume.
2. Advance the valve/balloon assembly with the loader over the guidewire.
3. Engage loader into sheath housing while maintaining a firm grip.
4. Advance the valve out of the loader into the large section of the sheath. Tap on the sheath housing to release air bubbles to the proximal end of the loader. Depress button valve on loader to de-air.
5. Advance the valve/balloon assembly through the sheath and position within the target valve or annuloplasty ring.

If needed, rotate the flex wheel on the handle to articulate the valve/balloon assembly into position.

CAUTION: To prevent possible leaflet damage, the valve should not remain in the sheath for over 5 minutes.

6. Ensure that the valve is correctly positioned between the two internal shoulders of the delivery system.
7. Begin valve deployment:
 - a) Unlock the Inflation device provided by Edwards Lifesciences.
 - b) Begin rapid pacing; once systolic blood pressure has decreased to 50 mmHg or below, balloon inflation can commence.
 - c) Deploy the valve by inflating the balloon with the entire volume in the Inflation device provided by Edwards Lifesciences, hold for 3 seconds and confirm that the barrel of the inflation device is empty to ensure complete inflation of the balloon.
 - d) Deflate the balloon. When the balloon has been completely deflated, turn off the pacemaker.

7.3.5 System Removal

1. If articulation was used, completely unflex the delivery system.

Retract the delivery system and guidewire into the sheath. Remove the loader and delivery system from the sheath.

CAUTION: Properly deflate the balloon and unflex the delivery system prior to removal.

2. Remove all devices when the ACT level is appropriate.
3. Remove the sheath from the access site, close the access site and confirm hemostasis.

8.0 How Supplied

STERILE: The valve is supplied sterilized with glutaraldehyde solution. The delivery system is supplied sterilized with ethylene oxide gas.

8.1 Storage

The valve must be stored at 10 °C to 25 °C (50 °F to 77 °F). Each jar is shipped in an enclosure containing a temperature indicator to detect exposure of the valve to extreme temperature.

The delivery system should be stored in a cool, dry place.

9.0 Magnetic Resonance (MR) Safety Information



MR Conditional

A person with the SAPIEN 3 or SAPIEN 3 Ultra transcatheter heart valve implant may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
Device Name	Edwards SAPIEN 3 transcatheter heart valve Edwards SAPIEN 3 Ultra transcatheter heart valve
Static Magnetic Field Strength (B0)	1.5 tesla (T) or 3.0 tesla (T)
Maximum Spatial Field Gradient	3000 gauss/cm (30 T/m) or less
RF Excitation	Circularly Polarized (CP) / Multichannel-2 (MC-2)
RF Transmit Coil Type	There are no Transmit Coil restrictions

Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2.0 W/kg
Scan Duration	2.0 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back to back series/scan without breaks)
MR Image Artifact	The presence of this implant may produce an image artifact.
For valve-in-prosthesis implantation or in the presence of other implants, please refer to the MRI safety information for the surgical valve or other devices prior to MR imaging.	

10.0 Patient Information

Patient education brochures are provided to each site and should be given to the patient to inform them of the risks and benefits of the procedure and alternatives in adequate time before the procedure to be read and discussed with their physician. A copy of this brochure may also be obtained from Edwards Lifesciences by calling 1.800.822.9837. A patient implant card request form is provided with each transcatheter heart valve. After implantation, all requested information should be completed on this form. The serial number may be found on the package and on the identification tag attached to the transcatheter heart valve. The original form should be returned to the Edwards Lifesciences address indicated on the form and upon receipt, Edwards Lifesciences will provide an identification card to the patient.

11.0 Recovered Valve and Device Disposal

The explanted valve should be placed into a suitable histological fixative such as 10% formalin or 2% glutaraldehyde and returned to the company. Refrigeration is not necessary under these circumstances. Contact Edwards Lifesciences to request an Explant Kit.

Used delivery system may be disposed of in the same manner that hospital waste and biohazardous materials are handled. There are no special risks related to the disposal of these devices.

12.0 Clinical Studies

SUMMARY OF PRIMARY CLINICAL STUDY

The PARTNER II Trial Overview, SAPIEN 3 Valve

SAPIEN 3 High Risk and Inoperable Cohort: The SAPIEN 3 High Risk and Inoperable Cohort of the PARTNER II trial (PIIS3HR) was a single arm, non-randomized, historical-controlled study to compare the third generation Edwards SAPIEN 3 system with the first generation Edwards SAPIEN system in patients who either have high risk for surgery or cannot undergo surgery (inoperable). The valve sizes used in the PIIS3HR trial included only the 23, 26 and 29 mm sizes. The 20 mm valve size was introduced into the trial after enrollment was completed with the three larger sizes, thus a separate nested registry, NR7, with identical inclusion/exclusion criteria as the PIIS3HR Cohort except for the aortic annulus diameter, was created to collect data for the 20 mm valve. Data from the PIIS3HR Cohort and NR7 are pooled for the statistical analyses. For convenience, this combined cohort is referred to as "PIIS3HR" hereafter.

The database included 583 eligible patients enrolled at 29 investigational sites in the U.S.

The PIIS3HR study used an independent Data Safety Monitoring Board (DSMB) that was instructed to notify Edwards Lifesciences of any safety or compliance issues, a Clinical Events Committee (CEC) that was responsible for adjudicating endpoint related events reported during the trial per *a priori* established VARC 2 definitions^[2], an ECG core laboratory for independent analysis of rhythm, and an echocardiographic core laboratory for independently analyzing all echocardiograms.

SAPIEN 3 Post Approval Study, High Risk and Inoperable Cohort: The post approval study consisted of continued follow-up of the patients enrolled in the PARTNER II High Risk and Inoperable Cohort and the nested registry for the 20 mm valve size (NR7). A database extract was performed on February 01, 2021, which yielded 664 operable high risk and inoperable patients that had been treated with an Edwards SAPIEN 3 transcatheter heart valve. The patients were treated between October 1, 2013 and August 11, 2015. The procedure was performed in 40 sites, 39 in the US and 1 in Canada.

Clinical Inclusion and Exclusion Criteria

Patients in the trial received a SAPIEN 3 transcatheter heart valve for symptomatic heart disease due to severe native calcific aortic stenosis who were deemed to be at high risk or inoperable for surgical aortic valve replacement.

Follow-up Schedule

Patients in the post approval study were followed-up annually through five years.

Clinical Endpoints

The endpoints analyzed in this application included death rate, adjudicated adverse events (stroke, aortic valve reinterventions, prosthetic valve dysfunction, rehospitalization from symptoms of aortic stenosis and/or complications of the valve procedure, and endocarditis), key site reported adverse events, valve performance based on echocardiographic data, New York Heart Association (NYHA) classification, length of hospitalization, and the Kansas City Cardiomyopathy Questionnaire (KCCQ) score.

SAPIEN 3 Intermediate Risk Cohort: The PIIS3i Cohort of the PARTNER II trial was a single arm, non-randomized, historical-controlled study to compare TAVR with the Edwards SAPIEN 3 system to the surgical aortic valve replacement (SAVR) arm from the previous PARTNER II trial Cohort A (PIIA-SAVR) in patients who were judged by a heart team to be at intermediate risk for open surgical therapy. The valve sizes used in the PIIS3i study included the 20, 23, 26, and 29 mm sizes.

Patients in PIIS3i were treated between February 2014 and September 2014. Patients in PIIA-SAVR were treated between January 2012 and November 2013. The database reflected data collected through December 10, 2015 and included 1,078 patients in PIIS3i enrolled at 51 investigational sites in the U.S. and 1,021 patients in PIIA-SAVR enrolled at 57 investigational sites in the U.S.

The PIIS3i study used an independent Data Safety Monitoring Board (DSMB) that was instructed to notify Edwards Lifesciences of any safety or compliance issues and a Clinical Events Committee (CEC) that was responsible for adjudicating endpoint related events reported during the trial. The CEC adjudicated the events per pre-established definitions, which were primarily Valve Academic Research Consortium-1 VARC-2 definitions^[2], with the following exceptions:

- Prosthetic valve dysfunction was adjudicated per VARC-1
- Aortic valve reintervention was adjudicated per protocol definition
- Rehospitalization for symptoms of aortic stenosis and/or complications of the valve procedure were adjudicated using the protocol and VARC-2 definitions as guidelines

The events in the PIIA-SAVR cohort were adjudicated by the CEC in accordance with the pre-specified, primarily VARC-1 definitions, with the following exceptions:

- Acute Kidney Injury (AKI) was adjudicated with a modified VARC-1 definition in which the CEC applied the 72-hour staging window to any AKI event that occurred within 30-days
- Aortic valve reintervention was adjudicated per the protocol definition
- Rehospitalization for symptoms of AS and/or complications of the valve procedure were adjudicated using the protocol and VARC-1 as guidelines
- Bleeding events were adjudicated irrespective of whether there was an identifiable, overt source of bleeding

An electrocardiogram (ECG) core laboratory was used for independent analysis of rhythm, an echocardiographic core laboratory for echocardiograms, and a computerized tomography (CT) core laboratory for baseline CTs for annulus dimensions.

The PARTNER 3 Trial Overview, SAPIEN 3 Valve

Patients were enrolled between March 2016 and June 2018. The database reflected data collected through December 21, 2018 and included 1000 patients. There were 71 investigational sites in the U.S, Australia, Canada, New Zealand, and Japan.

The PARTNER 3 trial was a prospective, randomized (1:1), controlled, multicenter study to compare TAVR with the Edwards SAPIEN 3 THV to SAVR. A subset of patients were enrolled in a computed tomography (CT) substudy to investigate the prevalence of Hypoattenuated Leaflet Thickening (HALT) and reduced leaflet mobility.

The PARTNER 3 trial used an independent Data Safety Monitoring Board (DSMB) that was instructed to notify the applicant of any safety or compliance issues and a Clinical Events Committee (CEC) that was responsible for adjudicating endpoint-related events reported during the trial. The CEC adjudicated the events per Valve Academic Research Consortium-2 (VARC-2) definitions^[2]. A CT core laboratory was used for assessment of baseline CTs for annulus dimensions and the CT images acquired in the CT substudy.

Clinical Inclusion and Exclusion Criteria

Patients in the database extract received a commercially available SAPIEN 3 transcatheter heart valve and surgical valves for symptomatic heart disease due to severe native calcific aortic stenosis who were deemed to be at low risk for surgical aortic valve replacement.

Clinical Endpoints

The endpoints analyzed in this application included: death rate, adjudicated adverse events (stroke, TIA, and aortic valve reinterventions, rehospitalization), key site reported adverse events, atrial fibrillation, length of index hospitalization, valve performance based on echocardiographic data, New York Heart Association (NYHA) classification, 6-minute walk test, and the Kansas City Cardiomyopathy Questionnaire (KCCQ) score. The analyses in the application focused on the 30-day and/or one-year time points.

SAPIEN 3 THV IN BICUSPID AORTIC VALVE FOR PATIENTS AT INTERMEDIATE OR GREATER SURGICAL RISK - STS/ACC TRANSCATHETER VALVE THERAPY REGISTRY (TVTR) ANALYSIS

A database extract was performed on November 15, 2017, which yielded 545 patients with bicuspid aortic valves that had been treated with an Edwards SAPIEN 3 transcatheter heart valve. The patients were treated between July 14, 2015 and August 15, 2016. The procedure was performed in 225 participating hospitals.

Adjudications were completed per the TVT Registry Coder's Data Dictionary by the Duke Clinical Research Institute (DCRI) for three adverse events: stroke, transient ischemic attack (TIA), and aortic valve reinterventions.

Clinical Inclusion and Exclusion Criteria

Patients in the database extract received a commercially available SAPIEN 3 transcatheter heart valve for symptomatic heart disease associated with a bicuspid aortic valve. The patients were treated based on clinical judgement of their treating physicians.

Follow-up Schedule

All patients were followed post implantation according to their local standards of care. The TVT Registry collects follow-up data at discharge, 30 days, and 1 year.

Clinical Endpoints

Data entered into the TVT Registry were collected through standardized data collection forms. The endpoints analyzed in this application included: death rate, adjudicated adverse events (stroke, TIA, and aortic valve reinterventions), key site reported adverse events, valve performance based on echocardiographic data, New York Heart Association (NYHA) classification, 5-meter walk test, and the Kansas City Cardiomyopathy Questionnaire (KCCQ) score. The analyses in the application focused on the 30-day and one-year time points.

SAPIEN 3 THV IN BICUSPID AORTIC VALVE FOR PATIENTS AT LOW SURGICAL RISK - PARTNER 3 BICUSPID REGISTRY ANALYSIS

Under the PARTNER 3 trial, a prospective, single-arm, registry was enrolled to establish a reasonable assurance of safety and effectiveness of TAVR with the Edwards SAPIEN 3 THV in patients with severe, native, calcific, aortic stenosis of a bicuspid aortic valve who are judged by a heart team to be at low risk for open surgical therapy. Patients were enrolled between October 2017 and May 2018, and data was collected through August 28, 2019 for 75 patients.

The data safety monitoring board and clinical events committee described in the PARTNER 3 trial above were also used in this registry.

Clinical Inclusion and Exclusion Criteria

Patients in the trial received a commercially available SAPIEN 3 transcatheter heart valve for symptomatic heart disease due to severe native calcific aortic stenosis of a bicuspid aortic valve who were deemed to be at low risk for surgical aortic valve replacement.

Follow-up Schedule

Patients were followed-up at discharge, 30 days, and 1 year, and will continue to be followed up to 10 years.

Clinical Endpoints

The clinical endpoints were the same as those listed above in the PARTNER 3 trial.

SAPIEN 3 THV Valve-in-Valve - STS/ACC Transcatheter Valve Therapy Registry (TVTR) Analysis

A database extract was performed on August 4, 2016, which yielded 314 patients that had been treated with an Edwards SAPIEN 3 transcatheter heart valve placed in a failed surgical aortic bioprosthesis (i.e., aortic valve-in-valve) and 311 patients that had been treated with an Edwards SAPIEN XT transcatheter heart valve (N = 241) or SAPIEN 3 transcatheter heart valve (N = 70) placed in a failed surgical mitral bioprosthesis (i.e., mitral valve-in-valve). Patients who presented with an existing valve-in-valve that was failing were excluded from the database extract.

The SAPIEN XT transcatheter heart valve was included in the database extract for the mitral valve-in-valve uses because there were fewer SAPIEN 3 transcatheter heart valve cases in the registry due to its relatively shorter commercial use history and the SAPIEN XT transcatheter heart valve data were considered to be generally applicable to the SAPIEN 3 transcatheter heart valve due to their similarities in design. The aortic valve-in-valve patients were treated between July 23rd, 2015 and June 29th, 2016 at 130 participating hospitals; the mitral valve-in-valve patients were treated at 112 participating hospitals between July 10th, 2014 and June 27th, 2016 for the SAPIEN XT transcatheter heart valve and between June 23rd, 2015 and June 15th, 2016 for the SAPIEN 3 transcatheter heart valve.

Adjudications were completed per the TVT Registry Coder's Data Dictionary by the Duke Clinical Research Institute (DCRI) for three adverse events: readmission for heart failure, stroke/transient ischemic attack (TIA), and aortic and mitral valve reinterventions.

Clinical Inclusion and Exclusion Criteria

Patients in the database extract received a commercially available SAPIEN 3 transcatheter heart valve (for both aortic and mitral valve-in-valve) or SAPIEN XT transcatheter heart valve (for mitral valve-in-valve only) for symptomatic heart disease associated with a failed (stenosed, insufficient, or combined) surgical bioprosthetic aortic or mitral valve. They were deemed to be at high or greater risk for open surgical therapy and were treated off-label based on the clinical judgement of their treating physicians.

Follow-up Schedule

All patients were followed post implantation according to their local standards of care. The TVT Registry collects follow-up data at discharge, 30 days, and 1 year.

Clinical Endpoints

Data entered into the TVT Registry were collected through standardized data collection forms. The endpoints analyzed in this application included: death rate, adjudicated adverse events (readmission for heart failure, stroke/TIA, and valve reinterventions), key site reported adverse events, valve performance based on echocardiographic data, New York Heart Association (NYHA) classification, 6-minute or 5-meter walk test, and the Kansas City Cardiomyopathy Questionnaire (KCCQ) score. The analyses in the application focused on the discharge and 30-day time points.

SAPIEN 3 and SAPIEN 3 Ultra THV THV-in-THV - STS/ACC Transcatheter Valve Therapy Registry (TVTR) Analysis

A database extract was performed on August 9, 2019. Patients were excluded if their previous transcatheter aortic valve replacement (TAVR) procedure was performed prior to April 23, 2007 (first implant in the US under the PARTNER trial). The database extract yielded 404 patients that had undergone a THV-in-THV procedure with an Edwards SAPIEN 3 (N = 402) or Edwards SAPIEN 3 Ultra (N = 2) transcatheter heart valve. These patients were treated between August 4, 2015 and July 11, 2019 at 188 participating hospitals.

To obtain more complete 1-year follow-up data, a treatment cutoff date of June 9, 2018 was then applied to the data set obtained above, which yielded 263 patients (SAPIEN 3 THV only) treated at 138 participating hospitals. The cutoff date was 14 months before the database extract date, which included a +60-day window for the 1-year visit. These 263 patients constituted the clinical data set used to support this application.

Clinical Inclusion and Exclusion Criteria

The initial database extract included all patients who received a commercially available Edwards SAPIEN 3 or Edwards SAPIEN 3 Ultra transcatheter heart valve in an aortic THV-in-THV procedure. The final data set was a subset of the initial database extract as described above.

Follow-up Schedule

All patients were followed post-implantation according to their local standards of care. The TVT Registry collects follow-up data at 30 days and 1 year.

Clinical Endpoints

Data entered into the TVT Registry were collected through standardized data collection forms. The endpoints analyzed in this application included: death, stroke/transient ischemic attack (TIA), valve reinterventions, key site reported adverse events, valve performance based on echocardiographic data, New York Heart Association (NYHA) classification, and the Kansas City Cardiomyopathy Questionnaire (KCCQ) score. The analyses in the application focused on 30-day and 1-year time points.

SAPIEN 3 Valve in Ring - STS/ACC Transcatheter Valve Therapy Registry (TVTR) Analysis & Mitral Implantation of Transcatheter Valves (MITRAL) Study Overview

A database extract from the TVT Registry was performed of all eligible patients in the registry on May 29, 2020, with a treatment cutoff date of May 28, 2019. This yielded 206 patients (SAPIEN 3 THV only) treated at 90 participating hospitals, which was pooled with 30 additional patients treated at 11 participating hospitals between February 2016 and October 2017 in the "valve in ring" arm of the MITRAL study, a sponsor-investigator IDE. While a formal poolability analysis was not conducted, the two datasets were deemed poolable because both included patients with comparable baseline characteristics, functional status, and clinical comorbidities, with similar risk for open surgical therapy. These 236 patients constituted the clinical dataset used to support this application.

Clinical Inclusion and Exclusion Criteria

The analysis population consists of all patients with a failing native mitral valve with an annuloplasty ring who had a procedure to implant the Edwards SAPIEN 3 THV by May 28, 2019 in the TVTR mitral module and patients from the "valve in ring" arm of the MITRAL study (NCT02370511: <https://clinicaltrials.gov/ct2/show/NCT02370511?term=mitral+mac&draw=3&rank=11>). The baseline characteristics of these patients show a morbid population at high surgical risk with no other good options, as evidenced by an STS score of 9.4 ± 6.4 and 81% of the patients being in New York Heart Association (NYHA) class of III/IV.

Follow-up Schedule

All patients were followed post-implantation according to their local standards of care. The TVT Registry collects follow-up data at 30 days and 1 year. All patients enrolled in the MITRAL study were scheduled for follow-up examinations at discharge, 30 days, 6 months, 1 year, and annually thereafter to a minimum of 5 years post-procedure.

Preoperative and post-operative assessments included physical assessment and patient interview, laboratory measurements, imaging tests, and health status/quality of life (QoL) questionnaire. Adverse events and complications were recorded at all visits.

Clinical Endpoints

The endpoints analyzed included: death, stroke/TIA, valve reinterventions, key site reported adverse events, valve performance based on echocardiographic data, NYHA classification, and the Kansas City Cardiomyopathy Questionnaire (KCCQ) score. The analyses focused on the 30-day and 1-year time points and were carried out descriptively.

SAPIEN 3 Ultra Confirmatory Study Overview

A prospective, single-arm, multicenter clinical study was conducted to confirm the procedural safety and effectiveness of the SAPIEN 3 Ultra system in patients with severe, calcific aortic stenosis (AS) who are at intermediate operative risk for surgical aortic valve replacement (SAVR).

The study enrolled 40 patients in Canada and United Kingdom.

It utilized an echocardiographic core laboratory for echocardiograms.

Follow-up Schedule

All patients were followed post implantation at discharge and 30 days, and will continue to be followed at 6 months, 1 year and annually thereafter for a minimum of 5 years.

Clinical Endpoints

The endpoints analyzed were procedure success (defined as freedom from mortality, conversion to surgery and moderate or severe paravalvular regurgitation at exit from the procedure room), key site reported adverse events, valve performance, and New York Heart Association (NYHA) classification.

PARTNER II SAPIEN 3 HIGH-RISK/INOPERABLE COHORT

Patient Accountability

All 583 eligible patients were successfully implanted with a SAPIEN 3 valve, which constitutes the Valve Implant (VI) population. Among the VI population, 491 patients were implanted via the transfemoral (TF) access route, and 92 patients via the transapical (TA) or transaortic (TAo) access route.

Table 12: Patient Accountability

	SAPIEN 3 Valve Overall	SAPIEN 3 Valve Transfemoral Access	SAPIEN 3 Valve Non-Transfemoral Access
Eligible Patient Population (EPP)	583	491	92
Valve Implant (VI) Population	583	491	92

Eligible Patient Population (EPP) consists of all enrolled patients who received treatment assignment from the database and entered into the catheterization laboratory/hybrid suite and who remained eligible to receive the implant.

Valve Implant (VI) Population consists of all enrolled patients who received a SAPIEN 3 valve, and retained the valve upon leaving the catheterization laboratory/hybrid suite.

Study Population Demographics and Baseline Parameters

The demographics of the study population are summarized in Table 13, which are typical of a TAVR study performed in the U.S.

Table 13: Patient Demographics and Baseline Characteristics – PIIS3HR VI Population

Characteristic	SAPIEN 3 Valve Overall (N = 583)	SAPIEN 3 Valve Transfemoral Access (N = 491)	SAPIEN 3 Valve Non-Transfemoral Access (N = 92)
Age, yr	82.6 ± 8.1	82.8 ± 8.2	81.7 ± 7.5
Male sex, no. (%)	338 (58.0%)	277 (56.4%)	61 (66.3%)
STS score	8.6 ± 3.7	8.4 ± 3.5	10.0 ± 4.3

Characteristic	SAPIEN 3 Valve Overall (N = 583)	SAPIEN 3 Valve Transfemoral Access (N = 491)	SAPIEN 3 Valve Non-Transfemoral Access (N = 92)
New York Heart Association (NYHA) class, no. (%):			
I/II	58 (9.9%)	51 (10.4%)	7 (7.6%)
III/IV	525 (90.1%)	440 (89.6%)	85 (92.4%)
Coronary artery disease, no. (%)	444 (76.2%)	360 (73.3%)	84 (91.3%)
Previous myocardial infarction, no. (%)	117 (20.1%)	87 (17.7%)	30 (32.6%)
Previous intervention, no. (%)			
Coronary-artery bypass grafting (CABG)	193 (33.1%)	145 (29.5%)	48 (52.2%)
Percutaneous coronary intervention (PCI)	199 (34.1%)	163 (33.2%)	36 (39.1%)
Prior aortic valvuloplasty	62 (10.6%)	49 (10.0%)	13 (14.1%)
Cerebral vascular accident (CVA), no. (%)	64 (11.0%)	53 (10.8%)	11 (12.0%)
Peripheral vascular disease, no. (%)	205 (35.2%)	155 (31.6%)	50 (54.3%)
Chronic obstructive pulmonary disease (COPD), no. (%):			
Any	259 (44.6%)	216 (44.1%)	43 (47.3%)
Oxygen-dependent	68 (11.8%)	58 (11.9%)	10 (11.0%)
Atrial fibrillation, no. (%)	255 (43.7%)	212 (43.2%)	43 (46.7%)
Permanent pacemaker, no. (%)	95 (16.3%)	78 (15.9%)	17 (18.5%)
Severe pulmonary hypertension, no. (%)	30 (5.1%)	24 (4.9%)	6 (6.5%)
Frailty, no. (%)	180 (30.9%)	162 (33.0%)	18 (19.6%)
Chest deformities that preclude an open chest procedure, no. (%)	4 (0.7%)	3 (0.6%)	1 (1.1%)
Cirrhosis, no. (%)	11 (1.9%)	9 (1.8%)	2 (2.2%)
Echocardiographic findings			
Effective Orifice Area (EOA), cm ²	0.7 ± 0.2	0.7 ± 0.2	0.7 ± 0.1
Mean aortic-valve gradient, mmHg	45.5 ± 14.3	45.7 ± 14.4	44.0 ± 13.2
Mean left ventricular ejection fraction (LVEF), %	56.4 ± 14.8	57.0 ± 14.5	53.2 ± 15.9
Moderate or severe mitral regurgitation, no./total no. (%)	69/541 (12.8%)	63/461 (13.7%)	6/80 (7.5%)

Safety and Effectiveness Results

Primary Endpoint

The composite rate of all-cause mortality, all stroke, and AI ≥ moderate at 30 days was 6.7% in the SAPIEN 3 cohort and 15.6% in the SAPIEN cohort, as shown in Table 14. The resulting proportion difference in the average treatment effect on the treated (ATT; ^[3]) was -6.9% (90% CI: [-13.3%, -0.5%]). Since the upper limit of the CI was < 7.5%, the non-inferiority was met.

Table 14: Primary Endpoint Analysis – Non-Inferiority Test SAPIEN 3 Valve (PIIS3HR VI Population) vs. SAPIEN Valve

Event at 30 days	SAPIEN 3 Valve (N = 583)	SAPIEN Valve (N = 326)	Weighted Proportion Difference in Average Treatment Effect on the Treated (ATT)
Composite of Death, Stroke and AI ≥ Moderate	6.7% [5.1%, 8.6%] ¹	15.6% [12.6%, 19.5%] ¹	-6.9% [-13.3%, -0.5%] ²
¹ For each individual study, the two-sided 90% stratified Wilson confidence interval was provided. ² The Wald-type two-sided 90% confidence interval using weighted mean and SD is provided.			

The Kaplan-Meier (K-M) estimates for all-cause mortality, cardiac mortality, and all stroke at 30 days for the SAPIEN 3 cohort and the SAPIEN cohort are provided in Table 15.

Table 15: Death and Stroke at 30 Days – SAPIEN 3 Valve vs. SAPIEN Valve (VI Population)

SAPIEN 3 Valve (N = 583)				SAPIEN Valve (N = 326)		
Event at 30 Days	No. Events	No. Pts with Events	K-M Estimated Event Rate ¹ (95% CI)	No. Events	No. Pts with Events	K-M Estimated Event Rate (95% CI)
Death	13	13	2.2% ([1.3%, 3.8%])	15	15	4.6% ([2.8%, 7.5%])
Cardiac Death	8	8	1.4% ([0.7%, 2.7%])	10	10	3.1% ([1.7%, 5.7%])
All Stroke	9	9	1.6% ([0.8%, 3.0%])	14	14	4.3% ([2.6%, 7.2%])

¹ Kaplan-Meier (K-M) estimates at 30 days used time to first event for each patient. Events occurring after 30 days were not included in this analysis.

Secondary Endpoints

Aortic insufficiency by visit is provided in Figure 8.




Figure 8: Aortic Insufficiency by Visit – SAPIEN 3 Valve (PIIS3HR VI Population) vs. SAPIEN Valve

The proportion of patients with AI ≥ moderate at 30 days was 3.0% in the SAPIEN 3 cohort and 14.3% in the SAPIEN cohort, which were found to be statistically significantly different (p=0.0051; Table 16).

Table 16: Aortic Insufficiency at 30 Days (SAPIEN 3 Valve vs. SAPIEN Valve VI Population)

Event at 30 Days	SAPIEN 3 Valve (N = 583)	SAPIEN Valve (N = 326)	Weighted Proportion Difference in Average Treatment Effect on the Treated (ATT)	P-value
AI ≥ Moderate, n/ Total no. (%) [95% CI]	16/532 (3.0%) [1.7%, 4.8%] ¹	40/280 (14.3%) [10.4%, 18.9%] ¹	-13.1% [-22.2%, -3.9%] ²	0.0051

¹ 95% Clopper-Pearson Exact confidence interval.

² The Wald-type two-sided 90% confidence interval using weighted mean and SD is provided.

The rate of major vascular complications at 30 days post implantation is shown in Figure 9. The rate was 5.0% for the SAPIEN 3 cohort and 10.1% for the SAPIEN cohort, which were found to be not statistically significantly different (p=0.0578; Table 17).




Figure 9: Major Vascular Complications at 30 Days – SAPIEN 3 Valve vs. SAPIEN Valve (VI Population)

Table 17: Major Vascular Complications at 30 Days – SAPIEN 3 Valve vs. SAPIEN Valve (VI Population)

Event at 30 Days	SAPIEN 3 Valve (N = 583)	SAPIEN Valve (N = 326)	Weighted Proportion Difference in Average Treatment Effect on the Treated (ATT)	P-value
Major Vascular Complications, n/ Total no. (%) [95% CI]	29/583 (5.0%) [3.4%, 7.1%]	33/326 (10.1%) [7.1%, 13.9%] ¹	-8.0% [-16.2%, 0.3%] ²	0.0578

¹ 95% Clopper-Pearson Exact confidence interval.

² The Wald-type two-sided 90% confidence interval using weighted mean and SD is provided.

Table 18 lists the hypothesis testing of the two secondary endpoints conducted with p-values in descending order for the Hochberg multiplicity adjustment steps. The largest p-value (p=0.0578 from major vascular complications) was greater than 0.05. As such, the null hypothesis was not rejected for the testing of major vascular complications at 30 days. The subsequent testing of AI ≥ moderate at 30 days had a p-value of 0.0051, which was less than 0.025. As such, the null hypothesis was rejected for AI ≥ moderate at 30 days, indicating that the SAPIEN 3 cohort was superior over the SAPIEN cohort in regards to AI ≥ moderate at 30 days.

Table 18: Secondary Endpoints for Labeling – SAPIEN 3 Valve vs. SAPIEN Valve (VI Population)

Endpoints	Original p-value	Inference
Major Vascular Complications at 30 Days	0.0578	> 0.05; reject the alternative hypothesis. Proceed to the rest of testing
AI at 30 Days	0.0051	< 0.025; claim superiority

Adverse Events

The key CEC adjudicated adverse events at 30 days are presented in Table 19.

Table 19: CEC Adjudicated Adverse Events at 30 Days (PIIS3HR VI Population)

30-Day Adverse Events	SAPIEN 3 Valve Overall	SAPIEN 3 Valve Transfemoral Access TF	SAPIEN 3 Valve Non-Transfemoral Access
Composite Event Rate of Death, All Stroke and AI \geq Moderate, n/N (%)	37/545 (6.8%)	27/463 (5.8%)	10/82 (12.2%)
Death			
From any cause, n/N (%)	13/583 (2.2%)	8/491 (1.6%)	5/92 (5.4%)
From cardiovascular cause, n/N (%)	8/583 (1.4%)	5/491 (1.0%)	3/92 (3.3%)
Stroke, n/N (%)	9/583 (1.5%)	8/491 (1.6%)	1/92 (1.1%)
AI \geq moderate, n/N (%)	16/532 (3.0%)	12/455 (2.6%)	4/77 (5.2%)
Myocardial Infarction, n/N (%)	3/583 (0.5%)	2/491 (0.4%)	1/92 (1.1%)
Major Vascular Complications, n/N (%)	29/583 (5.0%)	26/491 (5.3%)	3/92 (3.3%)
Acute Kidney Injury, Stage III, n/N (%)	6/583 (1.0%)	4/491 (0.8%)	2/92 (2.2%)
Disabling Bleeding Event, n/N (%)	37/583 (6.3%)	27/491 (5.5%)	10/92 (10.9%)
Aortic Valve Re-Intervention, n/N (%)	6/583 (1.0%)	4/491 (0.8%)	2/92 (2.2%)
Endocarditis, n/N (%)	1/583 (0.2%)	1/491 (0.2%)	0/92 (0.0%)
Conduction Disturbance Requiring Permanent Pacemaker, n/N (%)	76/583 (13.0%)	65/491 (13.2%)	11/92 (12.0%)

Other Results

Procedural Information

Overall, the mean duration in the catheterization laboratory/hybrid suite was 192.8 ± 59.3 min, the mean total procedure time was 86.3 ± 44.2 min, and the mean total anesthesia time was 193.7 ± 62.9 min. These duration times were slightly shorter in the TF patients. General anesthesia was used in the vast majority of cases; 15.9% of the TF patients had conscious sedation. Correct positioning of the valve was achieved in 99.1% of the patients. Five patients (0.9%; including 3 TF patients) were implanted with a second valve. One patient (0.2%) experienced valve embolization following rupture of the delivery balloon on annular calcium. This patient was converted to surgical aortic valve replacement and later died from aortic dissection.

Valve Performance

The mean EOA increased from 0.7 ± 0.2 cm² at baseline to 1.6 ± 0.4 cm² at 30 days, as shown in Figure 10.




Figure 10: Effective Orifice Area (PIIS3HR VI Population)

The average mean gradient decreased from 45.5 ± 14.3 mmHg at baseline to 11.1 ± 4.5 mmHg at 30 days, as shown in Figure 11.




Figure 11: Mean Gradient (PIIS3HR VI Population)

The mean peak gradient decreased from 75.8 ± 22.6 mmHg at baseline to 21.2 ± 8.5 mmHg at 30 days, as shown in Figure 12.




Figure 12: Peak Gradient (PIIS3HR VI Population)

The proportion of patients with AI \geq moderate was 7.3% at baseline and 3.0% at 30 days, as shown in Figure 13.





Figure 13: Aortic Insufficiency (PIIS3HR VI Population)

The proportion of patients with aortic paravalvular leak (PVL) \geq moderate was 2.9% at 30 days, as shown in Figure 14.





Figure 14: Aortic Paravalvular Leak (PIIS3HR VI Population)

NYHA

The NYHA class by visit is shown in Figure 15. For all patients, the mean NYHA class was 3.2 ± 0.6 at baseline and 1.7 ± 0.7 at 30 days.




Figure 15: NYHA Class by Visit (PIIS3HR VI Population)

Six Minute Walk Test (6MWT)

The improvement in mean 6MWT distance was 38.5 ± 110.2 meters from baseline to 30 days for all patients, 42.6 ± 107.8 meters for all TF patients, and 15.9 ± 121.2 meters for all TA/TAo patients.

Length of Stay (LoS)

The overall mean LoS was 6.8 ± 4.8 days, which included 3.0 ± 2.7 days in the ICU. The mean LoS was 6.1 ± 4.3 days (including 2.7 ± 2.3 days in the ICU) for the TF patients and 10.4 ± 5.4 days (including 4.8 ± 3.9 days in the ICU) for the TA/TAo patients.

Quality of Life (QoL)

QoL was measured using the visual analog scale (VAS) of the EuroQoL (EQ-5D) measure. The VAS is a self-assessment in which patients rate their well-being on a scale from 0 to 100 where 0 is the worst state they can imagine and 100 is the best state. During the trial, the mean improvement in VAS scale from baseline to 30 days was 14.6 ± 22.2 for all patients, 15.1 ± 21.5 for the TF patients, and 11.5 ± 25.7 for the TA/TAo patients.

Additional QoL instruments

The mean overall Kansas City Cardiomyopathy Questionnaire (KCCQ) summary score was 46.9 ± 22.6 at baseline, and 67.5 ± 22.6 at 30 days for the entire VI population. Except for self-efficacy which showed a small improvement, moderate to large improvements were observed in all other subscores at 30 days. In general, improvements in the TF patients were slightly larger compared to those observed in the TA/TAo patients.

Using the SF-36 norm based questionnaire, the physical component score for all patients improved from 32.0 ± 8.9 at baseline to 37.1 ± 9.7 at 30 days, and the mental component score improved from 46.9 ± 12.8 at baseline to 50.0 ± 12.5 at 30 days. In the TF patients, the physical component score improved from 31.8 ± 8.7 at baseline to 37.3 ± 9.8 at 30 days, and the mental component score improved from 46.8 ± 13.1 at baseline to 50.5 ± 12.2 at 30 days. In the TA/TAo patients, the physical component score improved from 32.9 ± 10.0 at baseline to 35.9 ± 9.4 at 30 days, and the mental component scores were 47.2 ± 11.1 at baseline and 47.2 ± 14.0 at 30 days.

SUMMARY OF SUPPLEMENTAL CLINICAL INFORMATION

Supplemental Clinical Study Design

Supplemental clinical data came from a study (referred to as "S3OUS" hereafter) conducted in Europe and Canada.

The S3OUS study was a non-randomized, prospective, multi-center study in inoperable, high surgical risk, and intermediate surgical risk patients who underwent implantation of the 23, 26, or 29 mm SAPIEN 3 valve.

Except the intermediate surgical risk patients, the inclusion/exclusion criteria of the S3OUS trial were largely similar to those of the PIIS3HR trial. The S3OUS study had a minimum age requirement (≥ 75 years) and the upper limit for AVA was higher ($< 1 \text{ cm}^2$ instead of $\leq 0.80 \text{ cm}^2$). Additionally, the S3OUS study included BAV within 30 days of the procedure (unless BAV was a

bridge to procedure), patients with planned concomitant surgical or transcatheter ablation for atrial fibrillation, hemodynamic or respiratory instability requiring inotropic support, mechanical ventilation or mechanical heart assistance within 30 days of screening; and the need for emergency surgery for any reason. Furthermore, the exclusion criteria in the S3OUS study excluded senile dementia and any neurologic disease which severely affected the ability to walk or perform everyday activities, and shortened the time interval regarding confirmed stroke or TIA (within 3 months instead of 6 months of the procedure). The follow-up periods were discharge or 7 days, whichever comes first, 30 days, 1 year, and annually thereafter to a minimum of 5 years post procedure.

Patient Accountability

Patients were treated at 14 investigational sites. Note that the intermediate risk patients enrolled in the S3OUS study were excluded from the analysis presented herein. The database included 102 “all treated” (AT) inoperable and high surgical risk patients. “All treated” population is defined to include all patients who were enrolled in the trial and for whom the study valve implantation procedures were started (i.e., the anesthesia was started).

One patient was excluded from the VI population. This patient experienced an aortic root rupture caused by displacement of a large lump of calcium with sharp edges through the native aortic annulus following balloon expansion of the SAPIEN 3 valve. The patient was subsequently converted to SAVR. After the patient was weaned off cardio-pulmonary bypass, bleeding in the region of the dorsal root occurred, and the patient died on the operating table.

A total of 56 patients were successfully implanted with a SAPIEN 3 valve via the transfemoral access route, and 45 via the transapical/transaortic access route, as shown in Table 20.

Table 20: Patient Accountability (S3OUS)

SAPIEN 3 Valve Overall		SAPIEN 3 Valve Transfemoral Access		SAPIEN 3 Valve Non-Transfemoral Access	
All Treated (AT) Population	Valve Implant (VI) Population	All Treated (AT) Population	Valve Implant (VI) Population	All Treated (AT) Population	Valve Implant (VI) Population
102	101	57	56	45	45

All Treated (AT) Population consists of all patients who were enrolled in the trial and for whom the study valve implantation procedures were started (i.e., anesthesia was started).

Valve Implant (VI) Population consists of all enrolled patients who received a SAPIEN 3 valve, and retained the valve upon leaving the catheterization laboratory/hybrid suite.

Study Population Demographics and Baseline Parameters

The demographics of the S3OUS study population are shown in Table 21.

Table 21: Patient Demographics and Baseline Characteristics (S3OUS AT Population)

Demographics and Baseline Characteristics	SAPIEN 3 Valve Overall (N = 102)	SAPIEN 3 Valve Transfemoral Access (N = 57)	SAPIEN 3 Valve Non-Transfemoral Access (N = 45)
Age, yr	84.1 ± 5.0	85.1 ± 4.6	83.0 ± 5.3
Male sex, no.(%)	40 (39.2%)	23 (40.4%)	17 (37.8%)
STS score	8.0 ± 4.7	8.2 ± 4.2	7.9 ± 5.2
Logistic EuroSCORE	24.1 ± 13.0	22.3 ± 11.3	26.4 ± 14.7
New York Heart Association (NYHA) class, no.(%):			
I/II	11 (10.8%)	6 (10.5%)	5 (11.1%)
III/IV	91 (89.2%)	51 (89.5%)	40 (88.9%)
Coronary artery disease, no.(%)	68 (66.7%)	36 (63.2%)	32 (71.1%)
Previous myocardial infarction, no.(%)	20 (19.6%)	7 (12.3%)	13 (28.9%)
Previous intervention, no.(%)			
Coronary-artery bypass grafting (CABG)	24 (23.5%)	10 (17.5%)	14 (31.1%)
Percutaneous coronary intervention (PCI)	34 (33.3%)	16 (28.1%)	18 (40.0%)
Prior aortic valvuloplasty	10 (9.8%)	8 (14.0%)	2 (4.4%)

Demographics and Baseline Characteristics	SAPIEN 3 Valve Overall (N = 102)	SAPIEN 3 Valve Transfemoral Access (N = 57)	SAPIEN 3 Valve Non-Transfemoral Access (N = 45)
Stroke, no.(%)	7 (6.9%)	4 (7.0%)	3 (6.7%)
Peripheral vascular disease, no.(%)	27 (26.5%)	10 (17.5%)	17 (37.8%)
Chronic obstructive pulmonary disease (COPD), no.(%):			
Any	25 (24.5%)	13 (22.8%)	12 (26.7%)
Oxygen-dependent	1 (1.0%)	1 (1.8%)	0 (0%)
Atrial fibrillation, no.(%)	48 (47.1%)	22 (38.6%)	26 (57.8%)
Permanent pacemaker, no.(%)	15 (14.7%)	7 (12.3%)	8 (17.8%)
Severe pulmonary hypertension, no.(%)	10 (9.8%)	6 (10.5%)	4 (8.9%)
Severe liver disease / Cirrhosis, no.(%)	1 (1.0%)	1 (1.8%)	0 (0%)
Echocardiographic findings			
Effective Orifice Area (EOA), cm ²	0.6 ± 0.2	0.6 ± 0.2	0.6 ± 0.1
Mean aortic-valve gradient, mmHg	44.8 ± 15.3	45.2 ± 14.7	44.2 ± 16.1
Mean left ventricular ejection fraction (LVEF), %	56.7 ± 9.1	57.7 ± 9.3	55.3 ± 8.7
Moderate or severe mitral regurgitation, no./total no. (%)	23/85 (27.1%)	9/48 (18.8%)	14/37 (37.8%)
Plus-minus values are means ± SD.			

Safety and Effectiveness Results

Key Adverse Events

Key adverse events as adjudicated by the CEC are presented in Table 22.

Table 22: CEC Adjudicated Adverse Events at 1 Year (S3OUS AT Population)

Outcomes	30 Days			1 Year		
	SAPIEN 3 Valve Overall	SAPIEN 3 Valve Transfemoral Access	SAPIEN 3 Valve Non-Transfemoral Access	SAPIEN 3 Valve Overall	SAPIEN 3 Valve Transfemoral Access	SAPIEN 3 Valve Non-Transfemoral Access
Composite Event Rate of Death, All Stroke and AI ≥ Moderate, n/N (%)	13/88 (14.8%)	3/50 (6.0%)	10/38 (26.3%)	25/82 (30.5%)	9/47 (19.1%)	16/35 (45.7%)
Death						
From any death, n/N (%)	8/102 (7.8%)	2/57 (3.5%)	6/45 (13.3%)	20/102 (19.6%)	7/57 (12.3%)	13/45 (28.9%)
From cardiovascular cause, n/N (%)	7/102 (6.9%)	2/57 (3.5%)	5/45 (11.1%)	9/102 (8.8%)	2/57 (3.5%)	7/45 (15.6%)
Stroke, n/N (%)	3/102 (2.9%)	1/57 (1.8%)	2/45 (4.4%)	5/102 (4.9%)	2/57 (3.5%)	3/45 (6.7%)
Aortic Insufficiency (AI) ≥ Moderate, n/N (%)	3/81 (3.7%)	1/49 (2.0%)	2/32 (6.3%)	1/62 (1.6%)	1/40 (2.5%)	0/22 (0.0%)
Disabling Stroke, n/N (%)	0/102 (0.0%)	0/57 (0.0%)	0/45 (0.0%)	1/102 (1.0%)	1/57 (1.8%)	0/45 (0.0%)
Myocardial Infarction, n/N (%)	2/102 (2.0%)	2/57 (3.5%)	0/45 (0.0%)	3/102 (2.9%)	2/57 (3.5%)	1/45 (2.2%)
Major Vascular Complications, n/N (%)	5/102 (4.9%)	1/57 (1.8%)	4/45 (8.9%)	N/A	N/A	N/A
Acute Kidney Injury - Stage III, n/N (%)	0/102 (0.0%)	0/57 (0.0%)	0/45 (0.0%)	N/A	N/A	N/A
Disabling Bleeding Event, n/N (%)	6/102 (5.9%)	3/57 (5.3%)	3/45 (6.7%)	N/A	N/A	N/A

Outcomes	30 Days			1 Year		
	SAPIEN 3 Valve Overall	SAPIEN 3 Valve Transfemoral Access	SAPIEN 3 Valve Non-Transfemoral Access	SAPIEN 3 Valve Overall	SAPIEN 3 Valve Transfemoral Access	SAPIEN 3 Valve Non-Transfemoral Access
Valve Dysfunction Requiring Intervention, n/N (%)	0/102 (0.0%)	0/57 (0.0%)	0/45 (0.0%)	N/A	N/A	N/A
Prosthetic Valve Endocarditis, n/N (%)	0/102 (0.0%)	0/57 (0.0%)	0/45 (0.0%)	1/102 (1.0%)	0/57 (0.0%)	1/45 (2.2%)
Conduction Abnormality Requiring Pacemaker, n/N (%)	14/102 (13.7%)	7/57 (12.3%)	7/45 (15.6%)	14/102 (13.7%)	7/57 (12.3%)	7/45 (15.6%)

The composite adverse event rate involving all-cause mortality, all stroke, and AI \geq moderate at 30 days for all patients is higher in the S3OUS cohort than PIIS3HR cohort (14.8% vs. 6.8%). This disparity is due to the composition of the study populations, specifically the S3OUS cohort comprises 44.1% TA/TAo patients vs. 15.8% TA/TAo patients in the PIIS3HR cohort. Note, the composite adverse event rate at 30 days for TF patients was similar, specifically, 6.0% in the S3OUS cohort and 5.8% in the PIIS3HR cohort.

The K-M estimates for all-cause mortality for all patients, the TF patients, and the TA/TAo patients are shown in Figure 16.




Figure 16: All-Cause Mortality at 1 Year (S3OUS AT Population)

Note: The confidence intervals are calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

The K-M estimates for the stroke rate for all patients, the TF patients, and the TA/TAo patients are shown in Figure 17.




Figure 17: All Stroke at 1 Year (S3OUS AT Population)

Note: The confidence intervals are calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

Valve Performance

The mean EOA increased from $0.6 \pm 0.2 \text{ cm}^2$ at baseline to $1.5 \pm 0.4 \text{ cm}^2$ at 30 days and $1.4 \pm 0.4 \text{ cm}^2$ at 1 year, as shown in Figure 18.




Figure 18: Effective Orifice Area (S3OUS VI Population)

The average mean gradient decreased from $44.8 \pm 15.4 \text{ mmHg}$ at baseline to $10.4 \pm 4.1 \text{ mmHg}$ at 30 days and maintained at $10.7 \pm 4.1 \text{ mmHg}$ at 1 year, as shown in Figure 19.




Figure 19: Mean Gradient (S3OUS VI Population)

The mean peak gradient decreased from 77.5 ± 24.9 mmHg at baseline to 21.0 ± 7.7 mmHg at 30 days, and maintained at 21.5 ± 8.2 mmHg at 1 year, as shown in Figure 20.




Figure 20: Peak Gradient (S3OUS VI Population)

The proportion of patients with aortic insufficiency \geq moderate was 9.8% at baseline, 3.7% at 30 days, and 1.6% at 1 year, as shown in Figure 21.




Figure 21: Aortic Insufficiency (S3OUS VI Population)

The proportion of patients with aortic PVL \geq moderate was 3.7% at 30 days, and 1.6% at 1 year, as shown in Figure 22.




Figure 22: Aortic Paravalvular Leak (S3OUS VI Population)

NYHA

The NYHA class by visit is shown in Figure 23. For all patients, the mean NYHA class decreased from 3.0 ± 0.5 at baseline to 1.6 ± 0.7 at 30 days and 1.8 ± 0.6 at 1 year.




Figure 23: NYHA Class by Visit (S3OUS AT Population)

PARTNER II POST APPROVAL STUDY OVERVIEW SAPIEN 3 VALVE

Patient Accountability

At the time of database extract, 216 of 664 high risk and inoperable patients were eligible for the 5-year visit and 186 (86.1%) completed the visit within the 30-day follow-up window. A detailed summary of the patient accountability at 30 days to 5 years is shown in Table 23.

Table 23: Patient Accountability (Valve Implant Population N= 664)

	30 Days	6 Months	1 Year	2 Years	3 Years	4 Years	5 Years
The patients with known visit accountability*	664 (100.0%)	664 (100.0%)	664 (100.0%)	664 (100.0%)	664 (100.0%)	664 (100.0%)	664 (100.0%)
On study at follow-up visit	647	601	555	465	389	303	216
Follow-up visit completed	644 (99.5%)	574 (95.5%)	540 (97.3%)	446 (95.9%)	363 (93.3%)	280 (92.4%)	186 (86.1%)
Missed Visit	3 (0.5%)	27 (4.5%)	15 (2.7%)	19 (4.1%)	26 (6.7%)	23 (7.6%)	30 (13.9%)
Discontinued prior to follow-up visit [†]	17	63	109	199	275	361	448
Death	17	59	99	180	241	310	376
Withdrawal	0	4	7	12	22	34	43
Lost to follow-up	0	0	1	3	7	12	24
Exit with other reason	0	0	2	4	5	5	5
Visit pending [‡]	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)

*The patients with known accountability is calculated as follows: (Number with FU visit + those with missed visit + those who have exited the study before the upper limit of the visit window) / the total population

[†]This includes all patients who exited the study prior to the end of the follow-up visit window and who have not had the visit.

[‡]Patients are not yet due or are within the follow-up visit window, but have not completed the follow-up visit yet.

The "Attempted Implant (AI)" population consisted of all patients in the dataset. The "Valve Implant (VI)" population consisted of those patients for whom the valve implant procedure succeeded. The numbers of patients in these two analysis populations are shown in Table 24.

Table 24: Analysis Populations

Analysis Population	Number of Patients
Attempted implant population	671
Valve implant population	664

Study Population Demographics and Baseline Characteristics

Patient demographics and baseline characteristics, as shown in Table 25, present an elderly, multimorbid cohort of patients, consistent with the high operative risk of the populations.

Table 25: Demographics and Baseline Characteristics (VI Population)

Demographics and Baseline Characteristics	Summary Statistics*
Age - years	82.9 ± 8.07 (664)
Male sex	51.1% (339/664)
Society of Thoracic Surgeons (STS) score	8.7 ± 3.75 (664)
New York Heart Association (NYHA) class	
Class I/II	12.0% (80/664)
Class III/IV	88.0% (584/664)
Previous myocardial infarction	18.2% (121/664)
Previous intervention	
Coronary artery bypass grafting (CABG)	31.2% (207/664)
Percutaneous coronary intervention (PCI)	33.0% (219/664)
Cerebrovascular accident (CVA)	11.4% (76/664)
Peripheral vascular disease (PVD)	35.2% (234/664)
Atrial fibrillation	42.2% (280/664)
Atrial flutter	5.6% (37/664)
Permanent pacemaker	15.5% (103/664)
Echocardiographic findings (Valve Implant Population)	
Valve area (cm ²)	0.66 ± 0.166 (622)
Mean gradient (mmHg)	45.9 ± 14.22 (649)
Peak gradient	76.4 ± 22.48 (649)

*Continuous measures - Mean ± SD (Total no.); Categorical measures % (no./Total no.)

Safety and Effectiveness Results

Safety Endpoints

The composite rate of death, all-stroke, and moderate or greater aortic insufficiency (AI) is presented in Table 26.

Table 26: Composite Event Rate of Death, All Stroke and AI ≥ Moderate (VI Population N = 664)

Outcome	30 Days	1 Year	2 Years	3 Years	4 Years	5 Years
Death, All Stroke and Aortic Insufficiency (AI) ≥ Moderate	42/627 (6.7%)	126/551 (22.9%)	208/553 (37.6%)	275/535 (51.4%)	344/507 (67.9%)	400/532 (75.2%)

The Kaplan-Meier estimates of CEC-adjudicated adverse events through five years are presented in Table 27. The all-cause mortality rate was 61.9% at 5 years, including a cardiovascular mortality rate of 42.3%.

Table 27: CEC-Adjudicated Adverse Events (VI Population N = 664)

CEC-Adjudicated Adverse Events	30 Days	1 Year	2 Years	3 Years	4 Years	5 Years
All-Cause Death	15, 15 (2.3%)	97, 97 (14.7%)	177, 177 (27.2%)	243, 243 (37.8%)	308, 308 (48.8%)	375, 375 (61.9%)
Cardiovascular Death	9, 9 (1.4%)	58, 58 (9.1%)	102, 102 (16.7%)	141, 141 (24.1%)	177, 177 (31.9%)	216, 216 (42.3%)
All Stroke	11, 11 (1.7%)	26, 26 (4.2%)	47, 43 (7.4%)	62, 55 (10.2%)	68, 60 (11.6%)	77, 69 (15.0%)
Major Stroke	5, 5 (0.8%)	13, 13 (2.1%)	30, 27 (4.8%)	42, 36 (6.8%)	45, 38 (7.4%)	54, 47 (10.8%)
Minor Stroke	6, 6 (0.9%)	13, 13 (2.1%)	47, 43 (7.4%)	62, 55 (10.2%)	68, 60 (11.6%)	77, 69 (15.0%)
Rehospitalization from Symptoms of Aortic Stenosis and/or Complications of the Valve Procedure	55, 52 (7.9%)	165, 126 (19.8%)	214, 151 (24.3%)	250, 172 (28.8%)	280, 186 (32.4%)	301, 197 (36.3%)
Kaplan-Meier estimate - no. of events, no. of subjects with the event (%)						

The percentage of subjects with moderate or greater aortic insufficiency (AI) through five years is presented in Table 28. At five years, moderate or greater AI was reported in 2.1% of subjects.

Table 28: Aortic Insufficiency through 5 Years (VI Population N = 664)

Outcome	1 Year	2 Years	3 Years	4 Years	5 Years
Aortic Insufficiency (AI) ≥ Moderate	10/449 (2.2%)	4/369 (1.1%)	5/286 (1.7%)	10/192 (5.2%)	3/146 (2.1%)

Kaplan-Meier estimates of site-reported adverse events through five years are presented in Table 29.

Table 29: Site Reported Adverse Events (VI Population N = 664)

Outcomes	1 Year	2 Years	3 Years	4 Years	5 Years
Bleeding	35.1% (297, 229)	38.5% (336, 247)	41.9% (372, 263)	46.5% (411, 280)	48.6% (433, 286)
Vascular Complication	11.3% (81, 75)	11.3% (81, 75)	11.5% (82, 76)	11.8% (84, 77)	11.8% (84, 77)
Myocardial Infarction	2.7% (19, 17)	4.5% (31, 26)	6.1% (38, 33)	7.8% (45, 39)	9.0% (49, 42)
New Permanent Pacemaker	16.2% (106, 106)	17.2% (111, 111)	18.4% (116, 116)	20.9% (126, 125)	21.8% (128, 127)
New onset Atrial Fibrillation	7.4% (53, 48)	8.6% (61, 54)	9.5% (65, 58)	11.0% (70, 63)	12.8% (74, 67)
Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)					

Effectiveness Endpoints

Valve Performance

The echocardiographic valve performance results are shown in Figure 24 through Figure 28. The effective orifice area increased from 0.66 cm² at baseline to 1.60 cm² at 1 year, which was maintained through 5 years (1.54 cm²). The mean gradient decreased from 45.8 mmHg at baseline to 12.1 mmHg at 1 year, which was maintained through 5 years (11.8 mmHg). Peak gradient decreased from 76.4 mmHg at baseline to 22.6 mmHg at 1 year, which was maintained through 5 years (21.4 mmHg). At five years, 92.5% of subjects had mild or less total AR. No subjects had severe total AR. At five years, 94.5% of subjects had mild or less paravalvular regurgitation. No subjects had severe paravalvular AR.




Figure 24: Effective Orifice Area by Visit (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 25: Mean Gradient by Visit (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 26: Peak Gradient by Visit (VI Population)

Note: Line plot with mean and standard deviation. The total number of patients at each visit time point only counted the patients with valid values.




Figure 27: Total Aortic Regurgitation (VI Population)

Note: The total number of patients at each visit time point only counted the patients with valid values.




Figure 28: Paravalvular Aortic Regurgitation (VI Population)

Note: The total number of patients at each visit time point only counted the patients with valid values.

Kaplan-Meier estimates for structural valve deterioration (SVD) through five years are presented in Table 30. At five years, the rate of SVD was 1.2%. Five events of SVD were reported for 5 patients.

Table 30: Structural Valve Deterioration through 5 Years (VI Population N = 664)

Outcomes	1 Year	2 Years	3 Years	4 Years	5 Years
Structural Valve Deterioration	0.3% (2, 2)	0.5% (3, 3)	0.7% (4, 4)	0.7% (4, 4)	1.2% (5, 5)
Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)					

NYHA Functional Class

The NYHA functional class distributions by visit are presented in Figure 29. At baseline, 88.1% of patients were in NYHA III/IV. At 5 years, the majority (88.3%) of patients were in NYHA I/II.




Figure 29: NYHA Class by Visit (VI Population)

Length of Stay

The mean index hospitalization stay was 6.8 days, which included an average of 3.0 days spent in the intensive care unit (ICU), as shown in Table 31.

Table 31: Index Hospitalization (VI Population)

Index Hospitalization	Length in Days*
Index hospitalization duration	6.8 ± 4.75 (664)
Intensive care stay	3.0 ± 2.84 (664)
* Mean ± SE (Total no.)	

Quality of Life

The results for the KCCQ overall summary score are presented in Figure 30. The mean score increased from 47.5 at baseline to 72.3 at 1 year, which was maintained through 5 years (69.0).




Figure 30: KCCQ Overall Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

Other Study Observations

Procedural Information

The procedural information is summarized in Table 32. General anesthesia was used in the majority (86.1%) of patients. Conversion to open heart surgery occurred in one patient due to device embolization.

Table 32: General Procedure Data (VI Population)

Procedural Data Summary	Summary Statistics (N=664)
Total procedure time (minute)	86.0 ± 44.09 (649)
Implanted valve size (mm)	
20	13.9% (92/664)
23	30.1% (200/664)
26	34.2% (227/664)
29	21.8% (145/664)
Total anesthesia time (minutes)	192.2 ± 63.06 (571)
Type of anesthesia used	
General	86.1% (572/664)
Conscious Sedation	13.9% (92/664)
Arterial site access	
Percutaneous	81.5% (539/661)
Surgical Cut-Down	14.2% (94/661)
Conduit	0.2% (1/661)
None	4.1% (27/661)
Fluoroscopy time (minutes)	18.6 ± 9.88 (655)

Procedural Data Summary	Summary Statistics (N=664)
Final position of the functioning transcatheter heart valve	
Acceptable - correct at intended site	99.2% (659/664)
Too ventricular	0.8% (5/664)
Too aortic	0.0% (0/664)
Valve not deployed	0/664 (0.0%)
Post-dilatation	16.3% (108/664)
Number of post-dilatation performed	
1	90.7% (98/108)
≥ 2	9.3% (10/108)
Valve in valve procedure performed (THV-in-THV)	1.1% (7/664)
Concomitant procedure(s) performed	2.4% (16/664)
Conversion to open heart surgery	
Valve dislodged	(0.2% (1/664)
*Continuous measures - mean ± SD (n); categorical measures - % (no./Total no.)	

Study Strength and Weaknesses

Data collected from the single-arm post approval study demonstrates acceptable clinical, functional, and hemodynamic long-term outcomes, and supports the safety and effectiveness of the SAPIEN 3 for the treatment of high risk subjects with severe symptomatic aortic stenosis. While efforts were made to ensure all treated subjects were followed through study completion, compliance with all follow-up requirements was difficult to obtain due to the existing comorbidities in this elderly, high surgical risk population, and their inability to attend in-office visits.

PARTNER II SAPIEN 3 INTERMEDIATE RISK COHORT

Patient Accountability

At the time of database lock, of the 1078 patients enrolled in the PMA study (PIIS3i), 99.2% (1069) patients are available for analysis at the completion of the study, the 1-year post-operative visit. Table 33 presents patient accountability in the PIIS3i and PIIA-SAVR cohorts. Of the 1,074 eligible patients (Eligible Patient or EP Population) in PIIS3i, 1,069 were successfully implanted with a SAPIEN 3 valve and constitute the PIIS3i Valve Implant (VI) population. Among the VI population, 943 patients were implanted via the transfemoral (TF) access route, and 126 patients via a non-transfemoral (non-TF; mainly transapical and transaortic) access route. Of the 938 eligible patients in the PIIA-SAVR cohort, 936 were successfully implanted with a surgical valve and constitute the PIIA-SAVR VI population.

Table 33: Patient Accountability

	All Enrolled Patients	Eligible Patient (EP) Population*	Valve Implant (VI) Population†
SAPIEN 3 Cohort	1078	1074	1069
TF	952	948	943
Non-TF	126	126	126
PIIA SAVR	1021	938	936
* Eligible Patient (EP) Population consists of all enrolled patients who were determined eligible after screening, entered into the catheterization laboratory and remained eligible to receive the assigned implant.			
† Valve Implant (VI) Population is a subset of the EP Population who received the assigned valve, and retained the valve upon leaving the catheterization laboratory.			

Study Population Demographics and Baseline Parameters

The demographics of the study population are typical for an aortic stenosis valve replacement study performed in the US, as summarized in Table 34 for the PIIS3i and PIIA-SAVR EP populations.

Table 34: Patient Demographics and Baseline Characteristics of the EP Population

Demographics & Characteristics*	SAPIEN 3 Valve			PIIA-SAVR (N = 938)
	Overall (N = 1074)	TF Only (N = 948)	Non-TF Only (N = 126)	
Age – years	81.9 ± 6.60	82.1 ± 6.57	80.7 ± 6.69	81.6 ± 6.73
Male sex	662/1074 (61.6%)	577/948 (60.9%)	85/126 (67.5%)	514/938 (54.8%)
Society of Thoracic Surgeons (STS) score	5.3 ± 1.29	5.3 ± 1.29	5.6 ± 1.28	5.8 ± 1.92
New York Heart Association (NYHA) class				
I/II	294/1074 (27.4%)	262/948 (27.6%)	32/126 (25.4%)	225/937 (24.0%)
III/IV	780/1074 (72.6%)	686/948 (72.4%)	94/126 (74.6%)	712/937 (76.0%)
Coronary artery disease	748/1074 (69.6%)	652/948 (68.8%)	96/126 (76.2%)	623/938 (66.4%)
Previous myocardial infarction	172/1074 (16.0%)	133/948 (14.0%)	39/126 (31.0%)	166/938 (17.7%)
Previous intervention				
Coronary artery bypass grafting (CABG)	301/1074 (28.0%)	248/948 (26.2%)	53/126 (42.1%)	241/938 (25.7%)
Percutaneous coronary intervention (PCI)	344/1074 (32.0%)	299/948 (31.5%)	45/126 (35.7%)	254/938 (27.1%)
Prior aortic valvuloplasty	55/1074 (5.1%)	51/948 (5.4%)	4/126 (3.2%)	45/938 (4.8%)
Cerebral vascular accident (CVA)	97/1074 (9.0%)	81/948 (8.5%)	16/126 (12.7%)	96/938 (10.2%)
Peripheral vascular disease	304/1074 (28.3%)	231/948 (24.4%)	73/126 (57.9%)	301/938 (32.1%)
Chronic obstructive pulmonary disease (COPD)				
Any	321/1072 (29.9%)	270/946 (28.5%)	51/126 (40.5%)	279/932 (29.9%)
Oxygen-dependent	53/1067 (5.0%)	46/942 (4.9%)	7/125 (5.6%)	26/925 (2.8%)
Atrial fibrillation	385/1074 (35.8%)	342/948 (36.1%)	43/126 (34.1%)	326/938 (34.8%)
Permanent pacemaker	142/1074 (13.2%)	121/948 (12.8%)	21/126 (16.7%)	113/938 (12.0%)
Severe pulmonary hypertension	25/1074 (2.3%)	19/948 (2.0%)	6/126 (4.8%)	N/A
Frailty	92/1074 (8.6%)	86/948 (9.1%)	6/126 (4.8%)	15/938 (1.6%)
Porcelain aorta	1/1074 (0.1%)	1/948 (0.1%)	0/126 (0.0%)	0/938 (0.0%)
Chest deformities that preclude an open chest procedure	1/1074 (0.1%)	1/948 (0.1%)	0/126 (0.0%)	0/938 (0.0%)
Cirrhosis	4/1074 (0.4%)	4/948 (0.4%)	0/126 (0.0%)	4/938 (0.4%)
Echocardiographic findings (Valve Implant Population)				
Effective orifice area (EOA) - cm ²	0.7 ± 0.17	0.7 ± 0.16	0.7 ± 0.18	0.7 ± 0.20
Mean aortic-valve gradient - mmHg	46.1 ± 12.63	46.1 ± 12.66	45.8 ± 12.47	44.7 ± 12.55
Mean left ventricular ejection fraction (LVEF) %	58.5 ± 13.36	58.8 ± 13.24	56.0 ± 14.05	55.4 ± 11.75
Moderate or severe mitral regurgitation	91/1033 (8.8%)	87/909 (9.6%)	4/124 (3.2%)	153/841 (18.2%)

*Continuous measures - Mean ± SD; Categorical measures – n/total no. (%)

Safety and Effectiveness Results

Primary Endpoints

The primary endpoint was a composite of all-cause death, stroke, and AI ≥ moderate at 1 year. The weighted proportion difference of the primary endpoint was -9.2% (90% CI: [-12.4%, -6.0%]) using the average treatment effect on the treated (ATT) method^[3], as shown in Table 35 and Figure 31. Since the upper limit of the CI was < 7.5%, non-inferiority was met.

Table 35: Primary Endpoint Non-Inferiority Test (VI Population)

	Observed Event rate		Propensity Score Quintile Pooled Proportion Difference (ATT Method*) [90% CI] [†]	Margin	Conclusion for Non-Inferiority Test
	SAPIEN 3 Valve (N = 1069)	PIIA-SAVR (N = 936)			
Composite of all-cause death, all stroke, and aortic insufficiency (AI) ≥ moderate at 1 year	13.0%	23.2%	-9.2% [-12.4%, -6.0%]	7.5%	Pass
* ATT: average treatment effect on the treated † Two-sided 90% Wald-type confidence interval					





Figure 31: Forest Plot – Composite of All Death, All Stroke and AI ≥ Moderate (VI Population)

The Kaplan-Meier (KM) estimates for all-cause death and all stroke at 1 year for the PIIS3i cohort and the PIIA-SAVR cohort are provided in Table 36, as well as Figure 32 and Figure 33, respectively.

Table 36: All-Cause Death and All Stroke at 1 Year (VI Population)

Endpoints	SAPIEN 3 Valve (N = 1069)			PIIA-SAVR (N = 936)			Propensity Score Quintile Pooled Proportion Difference (ATT Method [†])
	Observed Event Rate	Kaplan-Meier Event Rate*		Observed Event Rate	Kaplan-Meier Event Rate*		
		Point Estimate	Standard Error		Point Estimate	Standard Error	
All-cause death at 1 year	7.0%	7.1%	0.79%	12.4%	12.6%	1.09%	-5.2%
All stroke at 1 year	4.5%	4.6%	0.65%	7.9%	8.1%	0.91%	-3.5%
* Kaplan-Meier estimates were calculated at 365 days and included only the first event for each patient. Events occurring after 365 days were not included in this analysis. † ATT: average treatment effect on the treated							




No. at Risk

SAPIEN 3	1069	1057	1039	1014	988	960
SAVR	936	898	859	836	808	793

Figure 32: All-Cause Death through 1 Year (VI Population)

Note: The confidence intervals were calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, the confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.



No. at Risk

SAPIEN 3	1069	1030	1008	985	959	927
SAVR	936	848	805	786	758	741

Figure 33: All Stroke through 1 Year (VI Population)

Note: The confidence intervals were calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, the confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

The proportion of patients with AI \geq moderate at 1 year was 1.6% for the PIIS3i cohort and 0.3% for the PIIA-SAVR cohort, as shown in Table 37.

Table 37: Aortic Insufficiency (AI) \geq Moderate at 1 Year (VI Population)

	Observed Event Rate		Propensity Score Quintile Pooled Proportion Difference (ATT Method*)
	SAPIEN 3 Valve (N = 1069)	SAVR (N = 936)	
Aortic insufficiency (AI) \geq moderate	1.6%	0.3%	1.2%

* ATT: average treatment effect on the treated

Secondary Endpoints

The secondary endpoints were examined in a pre-specified order adjusted for the propensity quintiles using the ATT method. Table 38 summarizes the statistical conclusions on the non-inferiority hypothesis testing of the five secondary endpoints for labeling that were evaluated using a gatekeeping/hierarchical multiplicity adjustment procedure to control the overall type I error to 0.05. For each secondary endpoint, the upper limit of the confidence interval was less than the respective non-inferiority margin. Therefore, for each of the secondary endpoints for labeling, the SAPIEN 3 valve was non-inferior to SAVR.

Table 38: Secondary Endpoints for Labeling – Gatekeeping/Hierarchical Method (VI Population)

Pre-Specified Order for Gatekeeping/Hierarchical Method	Endpoints	Observed Event Rate		Weighted Proportion Difference in Average Treatment Effect on the Treated [90% CI] [†]	Margin	Conclusion for Non-Inferiority Test
		SAPIEN 3 Valve (N = 1069)	PIIA-SAVR (N = 936)			
No. 1	Composite of all death, all strokes, life threatening (disabling)/ major bleeding and major vascular complication at 30 days	18.3%	79.4%	-60.5% [-63.5%, -57.4%]	7.5%	Pass
No. 2	Major vascular and access complications through 30 days	5.8%	5.3%	0.3% [-1.5%, 2.0%]	5.0%	Pass
No. 3	Life threatening (disabling)/ major bleeding through 30 days	14.6%	78.2%	-63.2% [-66.2%, -60.2%]	5.0%	Pass
No. 4	All-cause death through 30 days	0.9%	3.7%	-2.7% [-3.9%, -1.5%]	2.5%	Pass
No. 5	All stroke through 30 days	2.6%	6.1%	-3.2% [-4.7%, -1.6%]	2.5%	Pass

[†] Two-sided 90% Wald-type confidence interval.

The forest plots for all-cause death and all stroke at 30 days are provided in Figure 34 and Figure 35, respectively.




Figure 34: Forest Plot – All-Cause Death at 30 Days (VI Population)

Note: As part of a pre-specified hierarchy, the hypothesis for this endpoint was tested using a hierarchical gatekeeping approach. The confidence interval shown here was not adjusted for multiplicity per the gatekeeping approach.




Figure 35: Forest Plot – All Stroke at 30 Days (VI Population)

Note: As part of a pre-specified hierarchy, the hypothesis for this endpoint was tested using a hierarchical gatekeeping approach. The confidence interval shown here was not adjusted for multiplicity per the gatekeeping approach.

Adverse Events

The key CEC-adjudicated adverse events through 1 year for the EP population are presented in Table 39.

Table 39: CEC-Adjudicated Adverse Events through 1 Year (EP Population)

Event*	SAPIEN 3 Valve			PIIA-SAVR
	Overall	TF Only	Non-TF Only	
7 Days				
Acute kidney injury: Stage III	5/1074 (0.5%)	3/948 (0.3%)	2/126 (1.6%)	N/A
30 Days				
Death	12/1074 (1.1%)	10/948 (1.1%)	2/126 (1.6%)	35/938 (3.7%)
Cardiac death	10/1074 (0.9%)	9/948 (0.9%)	1/126 (0.8%)	26/938 (2.8%)
Non-cardiac death	2/1074 (0.2%)	1/948 (0.1%)	1/126 (0.8%)	9/938 (1.0%)
Stroke	29/1074 (2.7%)	24/948 (2.5%)	5/126 (4.0%)	57/938 (6.1%)
Major (disabling) stroke	11/1074 (1.0%)	7/948 (0.7%)	4/126 (3.2%)	41/938 (4.4%)
Minor (non-disabling) stroke	18/1074 (1.7%)	17/948 (1.8%)	1/126 (0.8%)	16/938 (1.7%)
Myocardial infarction	3/1074 (0.3%)	3/948 (0.3%)	0/126 (0.0%)	17/938 (1.8%)
Major vascular complication	65/1074 (6.1%)	60/948 (6.3%)	5/126 (4.0%)	50/938 (5.3%)
Life threatening (disabling) or major bleeding	159/1074 (14.8%)	112/948 (11.8%)	47/126 (37.3%)	733/938 (78.1%)
Aortic valve re-intervention	1/1074 (0.1%)	1/948 (0.1%)	0/126 (0.0%)	0/938 (0.0%)
Any endocarditis	2/1074 (0.2%)	2/948 (0.2%)	0/126 (0.0%)	0/938 (0.0%)
Rhythm disturbance requiring permanent pacemaker	108/1074 (10.1%)	99/948 (10.4%)	9/126 (7.1%)	68/938 (7.2%)
1 Year				
Death	79/1074 (7.4%)	61/948 (6.4%)	18/126 (14.3%)	117/938 (12.5%)
Cardiac death	47/1074 (4.4%)	37/948 (3.9%)	10/126 (7.9%)	70/938 (7.5%)
Non-cardiac death	32/1074 (3.0%)	24/948 (2.5%)	8/126 (6.3%)	47/938 (5.0%)
Stroke	49/1074 (4.6%)	40/948 (4.2%)	9/126 (7.1%)	74/938 (7.9%)
Major (disabling) stroke	24/1074 (2.2%)	16/948 (1.7%)	8/126 (6.3%)	53/938 (5.7%)
Minor (non-disabling) stroke	25/1074 (2.3%)	24/948 (2.5%)	1/126 (0.8%)	22/938 (2.3%)
Aortic valve re-intervention	6/1074 (0.6%)	6/948 (0.6%)	0/126 (0.0%)	4/938 (0.4%)
Any endocarditis	8/1074 (0.7%)	7/948 (0.7%)	1/126 (0.8%)	6/938 (0.6%)

*Categorical measures - n. / total no. (%).

In addition, site-reported new-onset atrial fibrillation was 5.9% in the PIIS3i EP population and 29.2% in the PIIA-SAVR EP population.

Bleeding Rate

The bleeding rates utilizing the number of units transfused are presented in Table 40.

Table 40: Bleeding Rate Using Site-Reported Units Transfused (EP Population)

Event*	SAPIEN 3 Valve (N = 1074)	PIIA-SAVR (N = 938)
Transfusion units ≥ 2 and < 4	47/1074 (4.4%)	184/938 (19.6%)
Transfusion units ≥ 4	18/1074 (1.7%)	218/938 (23.2%)

*Site-reported Transfusion at Day 0 or Day 1; Categorical measures - n. / total no. (%)

Other Results

Procedural Information

In the PIIS3i EP population the mean duration in the catheterization laboratory was 187.3 ± 53.2 minutes, the mean total procedure time was 84.2 ± 40.7 minutes, and the mean total anesthesia time was 186.9 ± 61.1 minutes, all of which were slightly shorter in the TF group. General anesthesia was used in the vast majority of cases; 18.9% of the TF patients had conscious sedation. Correct positioning of the valve was achieved in 99.3% of the patients. Four (4) patients (0.4%, all TF patients) were implanted with a second valve. One (1) patient (0.1%) experienced valve embolization and two (2) patients (0.2%) experienced annular rupture.

In the PIIA-SAVR EP population, the mean duration in the operating room was 333.2 ± 96.4 min, the mean total procedure time was 237.5 ± 86.58 min, and the mean anesthesia time was 333.5 ± 108.42 min. General anesthesia was used in all patients.

Valve Performance

The measurements of EOA, mean gradient, peak gradient, total aortic regurgitation (AR), and aortic paravalvular leak (PVL) are presented in Figure 36 to Figure 40. The increase in EOA and decrease in gradient were sustained at 1 year. In PIIS3i, the proportion of patients with total AR ≥ moderate was 6.2% at baseline, 3.9% at 30 days, and 1.6% at 1 year, while in PIIA-SAVR, the proportion of patients with total AR ≥ moderate was 12.0% at baseline, 0.7% at 30 days, and 0.3% at 1 year. The proportion of patients with aortic PVL ≥ moderate was 3.8% at 30 days and 1.5% at 1 year in PIIS3i, as compared to 0.5% at 30 days and 0.3% at 1 year in PIIA-SAVR.




Figure 36: Effective Orifice Area (VI Population)




Figure 37: Mean Gradient (VI Population)




Figure 38: Peak Gradient (VI Population)




Figure 39: Total Aortic Regurgitation (VI Population)




Figure 40: Aortic Paravalvular Leak (VI Population)

NYHA

The NYHA classifications by visit are presented in Figure 41. In PIIS3i, 72.6% of the patients were in NYHA Class III or IV at baseline, which reduced to 6.3% at 30 days and 6.7% at 1 year, while in PIIA-SAVR, the percentage of patients in NYHA Class III or IV was 76.0% at baseline, 13.6% at 30 days, and 6.7% at 1 year. A side-by-side comparison of the results by access approach is presented in Figure 42.




Figure 41: NYHA Class by Visit (EP Population)




Figure 42: NYHA Class by Visit – TF versus non-TF Access (EP Population)

Six-Minute Walk Test (6MWT)

The improvements in mean 6MWT distance are presented in Table 41. The SAPIEN 3 valve patients had a similar increase in mean 6MWT distance from baseline to 1 year as the PIIA-SAVR patients.

Table 41: 6MWT Distance (EP Population)

6MWT Distance (m)*	SAPIEN 3 Valve			PIIA-SAVR
	All	TF	Non-TF	
Baseline	193.9 ± 118.1	194.1 ± 117.2	192.5 ± 125.5	179.3 ± 123.2
30 days	230.6 ± 126.1	234.6 ± 123.6	199.0 ± 140.6	166.7 ± 126.4
1 year	227.7 ± 134.7	230.6 ± 133.6	202.8 ± 142.1	219.2 ± 133.8

*Plus-minus values are means ± SD.

Length of Stay (LoS)

The results for LoS are presented in Table 42. Overall, the SAPIEN 3 valve patients had shorter LoS than the PIIA-SAVR patients.

Table 42: Length of Stay (EP Population)

Length of Stay (days)*	SAPIEN 3 Valve			PIIA-SAVR
	All	TF	Non-TF	
Overall	5.5 ± 5.7	5.0 ± 5.2	9.3 ± 7.7	11.9 ± 7.6
ICU	2.7 ± 3.0	2.5 ± 2.6	4.2 ± 4.9	5.6 ± 6.1

*Plus-minus values are means ± SD.

QoL

The QoL measurements using the Kansas City Cardiomyopathy Questionnaire (KCCQ) clinical summary score are presented in Figure 43. Except for self-efficacy which showed a small improvement, moderate to large improvements were observed in all other subscores at 30 days and were sustained at 1 year in the PIIA-SAVR EP population. A side-by-side comparison of the results by access approach is presented in Figure 44. In general, improvements in the TF group were slightly larger as compared to those observed in the Non-TF group.




Figure 43: KCCQ Clinical Summary Score (EP Population)




Figure 44: KCCQ Clinical Summary Score - TF versus non-TF Access (EP Population)

Additional QoL instruments

QoL was also measured using the visual analog scale (VAS) of the EuroQoL (EQ-5D) measure and the SF-36 Health Status Questionnaire. The VAS is a self-assessment in which patients rate their well-being on a scale from 0 to 100 where 0 is the worst state they can imagine and 100 is the best state. SF-36 uses 36 questions to measure functional health and well-being from the patient’s point of view and is generally reported in two (2) summary scores on a scale from 0 to 100 which evaluate physical (the

Physical Summary Score) and mental (the Mental Summary Score) health, with higher scores representing better functional health and well-being. The results of the VAS and SF-36 measures are presented in Table 43 and Table 44, respectively.

Table 43: EQ-5D Visual Analog Scale (EP Population)

EQ-5D Visual Analog Scale*	SAPIEN 3 Valve			PIIA-SAVR
	All	TF	Non-TF	
Baseline	60.3 ± 20.0	61.0 ± 19.8	55.1 ± 20.7	59.5 ± 20.5
30 days	74.0 ± 16.6	74.8 ± 16.6	68.5 ± 16.2	67.2 ± 19.5
1 year	74.4 ± 17.2	74.7 ± 17.1	71.8 ± 17.8	74.3 ± 16.7

*Plus-minus values are means ± SD.

Table 44: SF-36 Health Status Questionnaire Score (EP Population)

SF-36 Health Status Questionnaire Score*	SAPIEN 3 Valve			PIIA-SAVR
	All	TF	Non-TF	
Physical Component Score				
Baseline	34.7 ± 9.1	35.0 ± 9.1	33.1 ± 8.5	34.3 ± 9.0
30 days	39.7 ± 9.8	40.3 ± 9.7	34.8 ± 9.2	34.5 ± 8.4
1 year	40.0 ± 10.3	40.4 ± 10.2	37.0 ± 10.8	39.5 ± 10.4
Mental Component Score				
Baseline	48.0 ± 11.8	48.1 ± 11.8	47.0 ± 12.3	48.0 ± 12.3
30 days	51.8 ± 10.6	52.3 ± 10.4	47.8 ± 11.3	45.5 ± 13.3
1 year	52.5 ± 10.7	52.7 ± 10.8	50.7 ± 10.1	52.0 ± 11.3

*Plus-minus values are means ± SD.

PARTNER 3 SAPIEN 3 Low Risk Cohort

A. Accountability of the PMA Cohort

At the time of database lock, a total of 1000 subjects were randomized in the study, including 503 TAVR patients and 497 SAVR patients.

There were three different analysis populations defined in the protocol: Intention-to-Treat (ITT), As Treated (AT), and Valve Implant (VI), as summarized in Table 45 and Figure 45. The primary analysis was the AT analysis.

Table 45: Analysis Populations

Analysis Population	Definition	Number of Patients	
		TAVR	SAVR
Intention-To-Treat (ITT)	All randomized patients.	503	497
As Treated (AT)	All ITT patients for whom the index procedure was begun, whether or not the index procedure was completed.	496	454
Valve Implant (VI)	All AT patients who received and retained the intended valve during the index procedure.	495	453




Figure 45: Patient Population Flowchart

The overall follow-up compliance of the trial is summarized in Table 46.

Table 46: Overall Study Compliance (AT Population)

Patient Accountability	30-day Visit		1 Year Visit	
	TAVR (N=496)	SAVR (N=454)	TAVR (N=496)	SAVR (N=454)
Total patients	496	454	496	454
Non-eligible	2	11	6	30
Death	2	6	5	11
Withdrawal	0	3	0	12
Lost to follow-up	0	0	0	1
Exit with other reason	0	2	1	6
Visit not yet due	0	0	0	0
Eligible	494	443	490	424
Follow-up visit completed	96.5% (493)	96.5% (438)	97.8% (485)	91.2% (414)
Missed visit	0.2% (1)	1.1% (5)	1.0% (5)	2.2% (10)

B. Study Population Demographics and Baseline Characteristics

The demographics and baseline characteristics of the study population are typical for a TAVR study performed in the U.S., as shown in Table 47. The treatment cohorts were generally well balanced with respect to age, gender, and STS risk score.

Table 47: Patient Demographics and Baseline Characteristics (AT Population)

Demographics and Baseline Characteristics	Summary Statistics*	
	TAVR (N = 496)	SAVR (N = 454)
Age - years	73.3 ± 5.8	73.6 ± 6.1
Male sex	67.5% (335/496)	71.1% (323/454)
Society of Thoracic Surgeons (STS) score	1.9 ± 0.7	1.9 ± 0.6
New York Heart Association (NYHA) class		
I/II	68.8% (341/496)	76.2% (346/454)

Demographics and Baseline Characteristics	Summary Statistics*	
	TAVR (N = 496)	SAVR (N = 454)
III/IV	31.1% (155/496)	23.8% (108/454)
Previous myocardial infarction	5.7% (28/495)	5.8% (26/452)
Previous intervention		
Coronary artery bypass grafting (CABG)	3.0% (15/494)	1.8% (8/451)
Percutaneous coronary intervention (PCI)	18.8% (93/494)	16.2% (73/452)
Stroke or cerebrovascular accident (CVA)	3.4% (17/496)	5.1% (23/453)
Peripheral vascular disease (PVD)	6.9% (34/494)	7.3% (33/453)
Atrial fibrillation	15.7% (78/496)	18.8% (85/453)
Atrial flutter	3.0% (15/496)	2.4% (11/452)
Permanent pacemaker or defibrillator	2.4% (12/496)	2.9% (13/454)
Hostile chest	0.0% (0/496)	0.0% (0/454)
Echocardiographic findings (Valve Implant Population)		
Valve area (cm ²)	0.8 ± 0.2 (459)	0.8 ± 0.2 (424)
Mean gradient (mmHg)	49.4 ± 12.8 (484)	48.3 ± 11.8 (442)
Mean left ventricular ejection fraction (LVEF) %	65.7 ± 9.0 (472)	66.2 ± 8.6 (436)
Moderate or severe aortic regurgitation	3.9% (19/484)	2.5% (11/446)
Moderate or severe mitral regurgitation	1.3% (6/477)	3.2% (14/437)

* Continuous measures - Mean ± SD (Total no.); Categorical measures % (no./Total no.)

C. Safety and Effectiveness Results

1. Primary Endpoint

The primary endpoint results are presented in Table 48 and Figure 46. The rate of all-cause death, all stroke, or rehospitalization (valve-related or procedure-related and including heart failure) at 1-year was 8.5% in the TAVR group and 15.1% in the SAVR group. Since the upper limit of the 95% confidence interval for the difference in the primary endpoint event rate was < 6.0%, non-inferiority was achieved.

Table 48: Primary Endpoint Analysis (AT Population)

Event	Kaplan-Meier Rate*		Difference of KM Estimate (TAVR – SAVR)	95% CI* for the Difference	Non-inferiority Criterion
	TAVR (N=496)	SAVR (N=454)			
All-cause death, all stroke, or rehospitalization	8.5% (42)	15.1% (68)	-6.65%	[-10.77%, -2.52%]	Pass
All-cause death	1.0% (5)	2.5% (11)	-1.44%	[-3.13%, 0.24%]	
All stroke	1.2% (6)	3.1% (14)	-1.90%	[-3.77%, -0.02%]	
Rehospitalization	7.3% (36)	11.0% (49)	-3.74%	[-7.45%, -0.02%]	

*Kaplan-Meier estimate - % (no. of subjects with the event)




Figure 46: All-Cause Death, All Stroke, and Rehospitalization through 1 Year (AT Population)

Note: The confidence intervals were calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

2. Secondary Endpoints

Hypothesis testing

Since the primary endpoint passed the non-inferiority testing, the prespecified superiority testing was carried out on the six select secondary endpoints sequentially. TAVR with the SAPIEN 3 valve was found to be superior to SAVR in all six secondary endpoints, as shown in Table 49.

Table 49: Superiority Testing of Select Secondary Endpoints (AT Population)

No.	Endpoint	Summary Statistics*		Difference (TAVR – SAVR)	95% CI for the Difference	p-value (Superiority Test Result)
		TAVR (N=496)	SAVR (N=454)			
1	New onset atrial fibrillation at 30 days†	5.0% (21/417)	39.3% (145/369)	-34.3%	[-39.7%, -28.9%]	<.0001 (pass)
2	Length of index hospitalization (days)	2.9 ± 0.1 (496)	7.4 ± 0.2 (454)	-4.5	[-4.8, -4.1]	<.0001 (pass)
3	All-cause death, all stroke, or rehospitalization at 1 year	8.5% (42)	15.1% (68)	-6.6%	[-10.8%, -2.5%]	0.0016 (pass)
4	Death, KCCQ < 45 or KCCQ decrease from baseline ≥ 10 points at 30 days	3.9% (19/492)	30.6% (133/435)	-26.7%	[-31.4%, -22.1%]	<.0001 (pass)
5	Death or all stroke at 30 days	1.0% (5/496)	3.3% (15/454)	-2.3%	[-4.2%, -0.4%]	0.0214 (pass)
6	All stroke at 30 days	0.6% (3/496)	2.4% (11/454)	-1.8%	[-3.4%, -0.2%]	0.0284 (pass)

No.	Endpoint	Summary Statistics*		Difference (TAVR – SAVR)	95% CI for the Difference	p-value (Superiority Test Result)
		TAVR (N=496)	SAVR (N=454)			
*Continuous measures - Mean ± SE (Total no.); Categorical measures – observed rate, % (no./Total no.), except No. 3 - Kaplan-Meier rate, % (Total no.).						
† Patients with pre-procedural atrial fibrillation were excluded from the analysis.						

Valve Performance

The effective orifice area (EOA), mean aortic gradient, total aortic regurgitation (AR), and paravalvular regurgitation values obtained over time for the TAVR and SAVR patients are shown in Figure 47 through Figure 50, respectively. The increase in EOA and decrease in gradient were sustained through 1 year in both cohorts. In the TAVR cohort, the proportion of patients with total AR ≥ moderate was 0.8% at 30 days and 1.1% at 1 year, while in the SAVR cohort, the corresponding proportion was 0.4% at 30 days and 0.6% at 1 year. The proportion of patients with paravalvular regurgitation ≥ moderate was 0.8% at 30 days and 0.6% at 1 year in the TAVR cohort, as compared to 0.0% at 30 days and 0.8% at 1 year in the SAVR cohort.




Figure 47: Effective Orifice Area (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 48: Mean Aortic Gradient (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 49: Total Aortic Regurgitation (VI Population)




Figure 50: Paravalvular Regurgitation (VI Population)

New York Heart Association (NYHA) Functional Class

The NYHA classifications by visit are presented in Figure 51. At baseline, 31.3% of TAVR patients and 23.6% of SAVR patients were in NYHA III/IV. At 1 year the majority (~99%) of TAVR and SAVR patients were in NYHA Class I/II.




Figure 51: NYHA Class by Visit (VI Population)

Six-Minute Walk Test (6MWT)

The results for the 6MWT distance are presented in Figure 52. The TAVR patients showed an increase in mean 6MWT distance from 331.0 m at baseline to 349.1 m at 30 days, while SAVR patients showed a decrease from 329.4 m at baseline to 314.4 m at 30 days. The two cohorts had similar values at 1 year (347.6 m for TAVR and 351.7 m for SAVR).




Figure 52: 6MWT Distance (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

Quality of Life (QoL)

KCCQ

The results for the KCCQ overall summary score are presented in Figure 53. The mean score increased from 70.3 at baseline to 88.9 at 30 days and 89.9 at 1 year in TAVR patients and from 70.1 at baseline to 72.8 at 30 days and 88.1 at 1 year in SAVR patients.




Figure 53: KCCQ Overall Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

EuroQol (EQ-5D)

The results for the EQ-5D visual analog score (VAS) are presented in Figure 54. The mean score was 74.2 at baseline, 85.2 at 30 days, and 84.4 at 1 year in TAVR patients as compared to 75.2 at baseline, 76.5 at 30 days, and 84.7 at 1 year in SAVR patients.




Figure 54: EQ-5D Visual Analog Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

Short Form (SF)-36

The results for the SF-36 physical component summary score and mental component summary score are presented in Figure 55 and Figure 56, respectively.




Figure 55: SF-36 Physical Component Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 56: SF-36 Mental Component Summary (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

3. Adverse Events

The Kaplan-Meier estimates of the CEC-adjudicated adverse events through 1 year are presented in Table 50.

Table 50: CEC-Adjudicated Adverse Events through 1 Year (AT Population)

Event	Kaplan-Meier Rate*			
	30 Days		1 Year	
	TAVR (N = 496)	SAVR (N = 454)	TAVR (N = 496)	SAVR (N = 454)
All cause death	0.4% (2, 2)	1.1% (5, 5)	1.0% (5, 5)	2.5% (11, 11)
Cardiovascular death	0.4% (2, 2)	0.9% (4, 4)	0.8% (4, 4)	2.0% (9, 9)
All stroke	0.6% (3, 3)	2.4% (11, 11)	1.2% (6, 6)	3.1% (14, 14)
Disabling stroke	0.0% (0, 0)	0.4% (2, 2)	0.2% (1, 1)	0.9% (4, 4)
Non-disabling stroke	0.6% (3, 3)	2.0% (9, 9)	1.0% (5, 5)	2.2% (10, 10)
Death or stroke	1.0% (5, 5)	3.3% (16, 15)	1.8% (11, 9)	4.9% (25, 22)
Death or disabling stroke	0.4% (2, 2)	1.3% (7, 6)	1.0% (6, 5)	2.9% (15, 13)
Major vascular complications	2.2% (12, 11)	1.5% (8, 7)	2.8% (15, 14)	1.5% (8, 7)
Life-threatening / disabling, or major bleeding	3.6% (22, 18)	24.5% (123, 111)	7.7% (45, 38)	25.9% (132, 117)
Life-threatening / disabling bleeding	1.2% (9, 6)	11.9% (58, 54)	2.8% (17, 14)	12.8% (63, 58)
Major bleeding	2.6% (13, 13)	13.5% (65, 61)	5.3% (28, 26)	14.2% (69, 64)
Myocardial infarction	1.0% (5, 5)	1.3% (6, 6)	1.2% (6, 6)	2.2% (10, 10)
Requirement for renal replacement [†]	0.2% (1, 1)	0.7% (3, 3)	0.2% (1, 1)	0.7% (3, 3)
New permanent pacemaker implantation resulting from new or worsened conduction disturbances [‡]	6.5% (32, 32)	4.0% (18, 18)	7.3% (36, 36)	5.4% (24, 24)

Event	Kaplan-Meier Rate*			
	30 Days		1 Year	
	TAVR (N = 496)	SAVR (N = 454)	TAVR (N = 496)	SAVR (N = 454)
Coronary obstruction requiring intervention	0.2% (1, 1)	0.7% (3, 3)	0.2% (1, 1)	0.7% (3, 3)
New onset atrial fibrillation	5.0% (21, 21)	39.5% (145, 145)	7.0% (29, 29)	40.9% (150, 150)
Rehospitalization	3.4% (18, 17)	6.5% (30, 29)	7.3% (39, 36)	11.0% (59, 49)

*Kaplan-Meier rate (no. of events, no. of patients with the event).
†Requirement for renal replacement was based on the site-reported event. All the other events were based on the CEC-adjudicated results.
‡Patients with pacemaker or ICD at baseline were not counted as new events.
^{||}Rehospitalization (valve-related or procedure-related and including heart failure).

4. Subgroup Analysis

Gender Analysis

The protocol specified a subgroup analysis on gender. The primary endpoint result stratified by gender is presented in Figure 57.




Figure 57: All-Cause Death, All Stroke, and Rehospitalization through 1 Year Stratified by Gender (AT Population)

Note: The confidence intervals were calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

5. Other Study Observations

Procedural Information

The general procedural data are summarized in Table 51. Conscious sedation was used in the majority of TAVR patients (65.1%). The mean procedure time was significantly lower for TAVR compared to SAVR (58.6 minutes vs. 208.3 minutes). There were less concomitant (planned) procedures performed for TAVR patients compared to SAVR patients (6.9% vs. 26.4%). Additional TAVR and SAVR specific procedural data are presented in Table 52 and Table 53, respectively.

Table 51: General Procedural Data (AT Population)

Variable	Summary Statistics*	
	TAVR (N = 496)	SAVR (N = 454)
Subject treated according to their treatment assignment	99.8% (495/496)	99.8% (453/454)
Procedure aborted	0	1
Subject was assigned to TAVR but received SAVR	1	0
Procedure time (min)	58.6 ± 1.6 (496)	208.3 ± 2.9 (454)
Anesthesia type		
General	33.3% (165/496)	100.0% (454/454)
Conscious sedation	65.1% (323/496)	NA
Conversion from conscious sedation to general anesthesia during the procedure	1.6% (8/496)	NA
Anesthesia time (min)	138.7 ± 2.20 (496)	309.7 ± 3.7 (454)
Concomitant procedures	6.9% (34/496)	26.4% (120/454)
Annular area (mm ²)	473.5 ± 83.3 (486)	479.6 ± 87.6 (441)
*Continuous measures – mean ± SE (n) for procedure and anesthesia time, mean ± SD (n) for annular area; Categorical measures - % (no./Total no.)		

Table 52: TAVR Procedure Data (AT Population)

Variable	Summary Statistics*
	TAVR (N = 496)
Valve size	
20 mm	2.2% (11/496)
23 mm	29.2% (145/496)
26 mm	47.6% (236/496)
29 mm	21.0% (104/496)
Successful access, delivery and retrieval of the device delivery system	99.8% (494/495)
Arterial access method	
Left percutaneous	22.2% (109/490)
Right percutaneous	76.7% (376/490)
Left surgical cutdown	0.0% (0/490)
Right surgical cutdown	1.0% (5/490)
Total fluoroscopy time (min)	13.9 ± 0.3 (487)
BAV performed	57.8% (286/495)
Post dilatation performed	20.9% (103/494)
Number of post dilatations	
1	89.3% (92/103)
2	8.7% (9/103)
3	1.9% (2/103)
More than one SAPIEN 3 THV implanted	0.2% (1/495)
*Continuous measures - mean ± SE (n); categorical measures - % (no./Total no.). For patients in whom the procedure was aborted or who were converted to surgery, the rest of the procedure data except valve size were not collected.	

Table 53: SAVR Procedure Data (AT Population)

Variable	Summary Statistics*
	SAVR (N = 454)
Procedure aborted†	0.2% (1/454)
Valve size	
19 mm	2.9% (13/453)
21 mm	17.2% (78/453)
23 mm	36.6% (166/453)
25 mm	35.5% (161/453)
27 mm	6.8% (31/453)
29 mm	0.9% (4/453)
Total aortic cross clamp time (min)	74.3 ± 1.3 (453)
Total pump time (min)	97.7 ± 1.6 (453)
SAVR approach	
Sternotomy	95.4% (432/453)
Thoracotomy	0.9% (4/453)
Mini right upper thoracotomy	2.9% (13/453)
Port access	0.2% (1/453)
Other	0.7% (3/453)
Successful implantation of the surgical valve	100.0% (453/453)
*Continuous measures - mean ± SE (n); categorical measures - % (no./Total no.).	
†For patients in whom the procedure was aborted, the rest of the procedure data were not collected.	

Computed Tomography (CT) Sub-Study

There were 184 TAVR and 162 SAVR patients at 30 days and 160 and 134 patients at 1 year, respectively, who had at least one adequate CT for leaflet assessments. The HALT and leaflet mobility imaging findings are summarized in Table 54, along with the associated mean aortic pressure gradients. The mean aortic pressure gradients at 1 year stratified by HALT and leaflet mobility at 30 days are summarized in Table 55 and Table 56, respectively. The rate of death, stroke or TIA at 1 year stratified by HALT and leaflet mobility at 30 days are summarized in Table 57 and Table 58, respectively. The CT substudy was not powered to compare the relative incidence or the severity of HALT or reduced leaflet mobility between the TAVR and SAVR cohorts, or to determine whether late clinical outcomes were affected by the presence of HALT or reduced leaflet mobility.

Table 54: HALT and Leaflet Mobility Findings and Associated Mean Gradients

Findings	Summary Statistics*			
	30 Days		1 Year	
	TAVR (N = 184)	SAVR (N = 162)	TAVR (N = 160)	SAVR (N = 134)
Proportion of patients on oral anticoagulants at time of scan	6.0% (11/184)	21.0% (34/162)	8.1% (13/160)	13.4% (18/134)
HALT [†]				
No thickening	84.8% (156/184)	95.7% (155/162)	74.4% (119/160)	82.1% (110/134)
Mean gradient (mmHg)	12.5 ± 0.3 (156)	10.8 ± 0.3 (155)	13.7 ± 0.4 (115)	11.7 ± 0.4 (106)
< 25% leaflet length thickened	4.9% (9/184)	1.2% (2/162)	11.3% (18/160)	7.5% (10/134)
Mean gradient (mmHg)	11.4 ± 0.9 (9)	16.5 ± NA (1)	12.9 ± 0.7 (18)	9.3 ± 1.8 (8)
25%-50% leaflet length thickened	3.3% (6/184)	1.9% (3/162)	6.3% (10/160)	5.2% (7/134)

Findings	Summary Statistics*			
	30 Days		1 Year	
	TAVR (N = 184)	SAVR (N = 162)	TAVR (N = 160)	SAVR (N = 134)
Mean gradient (mmHg)	13.7 ± 1.7 (6)	9.4 ± 1.4 (3)	13.2 ± 1.8 (10)	15.1 ± 2.4 (7)
50%-75% leaflet length thickened	6.5% (12/184)	0.6% (1/162)	5.0% (8/160)	3.7% (5/134)
Mean gradient (mmHg)	15.2 ± 1.9 (12)	9.8 ± NA (1)	16.9 ± 3.3 (8)	16.1 ± 4.0 (5)
> 75% leaflet length thickened	0.5% (1/184)	0.6% (1/162)	3.1% (5/160)	1.5% (2/134)
Mean gradient (mmHg)	10.2 ± NA (1)	16.8 ± NA (1)	20.2 ± 6.2 (5)	9.0 ± 4.2 (2)
Number of leaflets with HALT	6.7% (37/552)	2.3% (11/486)	12.7% (61/480)	8.2% (33/402)
0 leaflets thickening	156	155	119	110
1 leaflet thickening	21	4	26	15
2 leaflets thickening	5	2	10	9
3 leaflets thickening	2	1	5	0
Leaflet mobility [‡]				
Unrestricted	85.3% (145/170)	96.8% (149/154)	77.6% (118/152)	83.0% (108/129)
Mean gradient (mmHg)	12.2 ± 0.3 (145)	10.7 ± 0.3 (148)	13.3 ± 0.4 (114)	12.0 ± 0.5 (105)
Partially restricted, restriction limited to base	5.3% (9/170)	1.3% (2/154)	11.8% (18/152)	8.5% (11/129)
Mean gradient (mmHg)	11.4 ± 0.9 (9)	14.6 ± 1.9 (2)	12.5 ± 0.6 (18)	9.9 ± 1.6 (9)
Partially restricted (< 50%)	5.3% (9/170)	1.3% (2/154)	3.9% (6/152)	3.1% (4/129)
Mean gradient (mmHg)	15.5 ± 2.4 (9)	10.3 ± 0.5 (2)	14.0 ± 2.8 (6)	15.6 ± 3.0 (4)
Partially restricted (50%-75%)	3.5% (6/170)	0.0% (0/154)	4.6% (7/152)	3.9% (5/129)
Mean gradient (mmHg)	12.8 ± 1.7 (6)	NA	21.8 ± 3.9 (7)	11.3 ± 3.6 (5)
Largely immobile	0.6% (1/170)	0.6% (1/154)	2.0% (3/152)	0.8% (1/129)
Mean gradient (mmHg)	13.3 ± NA (1)	16.8 ± NA (1)	19.5 ± 8.1 (3)	13.1 ± NA (1)
Number of leaflets partially restricted or largely immobile				
0 leaflet	145	149	118	108
1 leaflet	21	2	22	13
2 leaflets	4	2	8	8
3 leaflets	0	1	4	0

*Continuous measures - mean ± SE (n); categorical measures - % (no./Total no.). The analysis population included all the patients enrolled in the CT substudy and had at least one adequate CT for leaflet assessments.

[†]HALT was defined as: the presence of any hypoattenuated leaflet thickening in any singular leaflet as identified by an independent CT core laboratory. The extent of the hypoattenuated leaflet thickening was graded with regards to the entire leaflet as: None, < 25%, 25-50%, 50-75%, or > 75%. If more than one leaflet had the appearance of HALT, the thickening measure of the most impacted leaflet was used. Presence of any degree of HALT on any one leaflet rendered a finding.

[‡]Leaflet mobility was determined by an independent CT core laboratory and included: unrestricted, partially restricted mobility limited to the base of a leaflet, partially restricted mobility involving more than the base of the leaflet but less than 50% of the leaflet, partially restricted mobility involving more than 50% of the leaflet but less than 75% of the leaflet, and/or a largely immobile leaflet. Presence of any degree of restriction or immobility on any one leaflet rendered a finding.

Table 55: Mean Aortic Gradient at 1 Year Stratified by HALT at 30 Days

	Summary Statistics*			
	HALT at 30 Days		No HALT at 30 Days	
	TAVR (N = 28)	SAVR (N = 7)	TAVR (N = 156)	SAVR (N = 155)
Mean gradient	13.6 ± 1.2 (24)	13.7 ± 2.7 (5)	13.6 ± 0.4 (137)	11.8 ± 0.4 (125)

*Mean ± SE (n). The analysis population included all the patients enrolled in the CT substudy and had an adequate CT for leaflet assessments at 30 days.

Table 56: Mean Aortic Gradient at 1 Year Stratified by Leaflet Mobility at 30 Days

	Summary Statistics*			
	Reduced Leaflet Mobility at 30 Days		Unrestricted at 30 Days	
	TAVR (N = 25)	SAVR (N = 5)	TAVR (N = 145)	SAVR (N = 149)
Mean gradient	13.7 ± 1.28 (23)	14.2 ± 3.48 (4)	13.3 ± 0.4 (124)	11.7 ± 0.4 (119)

*Mean ± SE (n). The analysis population included all the patients enrolled in the CT substudy and had an adequate CT for leaflet assessments at 30 days.

Table 57: All-Cause Mortality, All Stroke or TIA at 1 Year Stratified by HALT at 30 Days

1-Year Endpoint	Kaplan-Meier Rate*			
	HALT at 30 Days		No HALT at 30 Days	
	TAVR (N = 28)	SAVR (N = 7)	TAVR (N = 156)	SAVR (N = 155)
All-cause mortality	0.0% (0)	0.0% (0)	1.3% (2)	1.4% (2)
All stroke	0.0% (0)	0.0% (0)	0.7% (1)	0.0% (0)
TIA	5.6% (1)	0.0% (0)	1.3% (2)	0.0% (0)
All-cause mortality or all stroke or TIA	5.6% (1)	0.0% (0)	3.3% (5)	1.4% (2)

*Kaplan-Meier rate (no. of patients with event). The analysis population included all the patients enrolled in the CT substudy and had an adequate CT for leaflet assessments at 30 days. The Kaplan-Meier analysis used the CT test date as the start date in determining time to event. Presence of any degree of HALT on any one leaflet rendered a finding and inclusion in the HALT cohort.

Table 58: All-Cause Mortality, All Stroke or TIA at 1 Year Stratified by Leaflet Mobility at 30 Days

1-Year Endpoint	Kaplan-Meier Rate*			
	Reduced Leaflet Mobility at 30 Days		Unrestricted at 30 Days	
	TAVR (N = 25)	SAVR (N = 5)	TAVR (N = 145)	SAVR (N = 149)
All-cause mortality	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)
All stroke	0.0% (0)	0.0% (0)	0.0% (0)	0.0% (0)
TIA	6.3% (1)	0.0% (0)	6.3% (1)	0.0% (0)
All-cause mortality or all stroke or TIA	6.3% (1)	0.0% (0)	3.6% (5)	1.4% (2)

*Kaplan-Meier rate (no. of patients with event). The analysis population included all the patients enrolled in the CT substudy and had an adequate CT for leaflet assessments at 30 days. The Kaplan-Meier analysis used the CT test date as the start date in determining time to event. Reduced leaflet mobility included any of the following assessments: partially restricted limited to base, partially restricted involving more than the base but less than 50% of the leaflet, partially restricted involving more than 50% but less than 75% of the leaflet, and/or largely immobile. Presence of any degree of restriction or immobility on any one leaflet rendered a finding and inclusion in the reduced leaflet mobility cohort.

SAPIEN 3 THV IN BICUSPID AORTIC VALVE FOR PATIENTS AT INTERMEDIATE OR GREATER SURGICAL RISK – STS/ACC TRANSCATHETER VALVE THERAPY REGISTRY (TVTR) ANALYSIS

Patient Accountability

At the time of database extract, of the 545 patients in the bicuspid aortic valve cohort, 527 patients were eligible for the 30-day visit, and 486 (92.2%) patients paid a visit within the 30-day follow-up window defined as the period between 21 and 75 days post-procedure. Of the 465 patients eligible for the 1-year visit, 309 (66.5%) paid a visit within the 1 year follow-up window defined as the period between 305 and 425 days post-procedure. A detailed summary of the patient accountability at 30 days and 1 year is shown in Table 59.

Table 59: Patient Visit Accountability

	30-day Visit	1-year Visit
Total patients	545	545
Non-eligible*	18	80
-Death	14	43
-Withdrawal	2	8
-Lost to follow-up	2	29
-Visit not yet due [‡]	0	0
Eligible	527	465
-Follow-up visit completed	486 (92.2%)	309 (66.5%)
-Missed visit [†]	41 (7.8%)	156 (33.5%)

*This includes all patients who exited the study prior to the end of the follow-up visit window and those who have not had the visit.

[‡]Patients have not reached the end of the visit window and have not completed the follow-up visit yet.

[†]Data extract date has exceeded the end of the visit window and the patients have not reported the visit data.

The "Attempted Implant" population consisted of all patients entered into the registry with a bicuspid aortic valve. The "Valve Implant" population consisted of those patients for whom the valve implant procedure has started and a "No" was indicated for both "procedure aborted" and "conversion to open heart surgery." The "Valve Implant" population consists of 540 patients as 5 patients were converted to open heart surgery and did not receive the SAPIEN 3 transcatheter heart valve.

Patient Demographics and Baseline Characteristics

The demographics and baseline characteristics of bicuspid aortic valve patients, as shown in Table 60, present a multimorbid cohort of patients with a mean STS score of 5.5 ± 4.0.

Table 60: Patient Demographics and Baseline Characteristics - Bicuspid Population (Attempted Implant Population)

Demographics and Baseline Characteristics	Summary Statistics*
Age - years	73.4 ± 11.1 (545)
Male sex	349 / 545
Society of Thoracic Surgeons (STS) score	5.5 ± 4.0 (538)
New York Heart Association (NYHA) class	
I/II	106 / 535 (19.8%)
III/IV	429 / 535 (80.2%)
Previous myocardial infarction	119 / 544 (21.9%)
Previous intervention	
Coronary artery bypass grafting (CABG)	101 / 543 (18.6%)
Percutaneous coronary intervention (PCI)	138 / 545 (25.3%)
Prior aortic valvuloplasty	34 / 545 (6.2%)
Cerebrovascular accident (CVA)	56 / 545 (10.3%)
Peripheral vascular disease	128 / 544 (23.5%)

Demographics and Baseline Characteristics	Summary Statistics*
Atrial fibrillation	183 / 545 (33.6%)
Permanent pacemaker	54 / 545 (9.9%)
Porcelain aorta	12 / 545 (2.2%)
Hostile chest	44 / 545 (8.1%)
Echocardiographic findings (Valve Implant Population)	
Valve area - cm ²	0.7 ± 0.2 (524)
Mean aortic valve gradient - mmHg	44.9 ± 15.5 (535)
Mean left ventricular ejection fraction (LVEF) %	52.9 ± 15.5 (534)
Moderate or severe aortic regurgitation	91 / 536 (17.0%)
Moderate or severe mitral regurgitation	101 / 438 (23.1%)

*Continuous measures - Mean ± SD (Total no.); categorical measures - n. / Total no. (%).

Safety and Effectiveness Results

Safety Endpoints

The mortality rates at discharge, 30 days, one-year and the Kaplan-Meier curve for all-cause mortality are shown in Table 61 and Figure 58, respectively. There were a total of 12 deaths reported at 30 days and 43 deaths reported at one year.

Table 61: Death Rate - Bicuspid Population (Attempted Implant Population)

	Discharge*	30 Days [†]	1 Year [†]
All-cause death [‡]	1.8% (10)	2.3% (12)	10.3% (43)
Cardiac death	1.1% (6)	1.3% (7)	3.0% (13)

*Observed rate - % (n).

[†]Kaplan-Meier estimate - % (n)

[‡]Includes all deaths reported in TVTR and identified through CMS linkage.




Figure 58: All-Cause Death Rate - Bicuspid Population (Attempted Implant Population)

The DCRI adjudicated events, including all strokes, TIAs and aortic valve reinterventions at discharge, 30 days and one year are shown in Table 62.

Table 62: Duke Clinical Research Institute Adjudicated Events - Bicuspid Population (Attempted Implant Population)

Events	Discharge*	30 Days [†]	1 Year [†]
All strokes	1.5% (8, 8)	1.9% (10, 10)	2.7% (13, 13)
Ischemic stroke	1.5% (8, 8)	1.9% (10, 10)	2.7% (13, 13)
Hemorrhagic stroke	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)
Transient ischemic attack (TIA)	0.2% (1, 1)	0.2% (1, 1)	0.2% (1, 1)
Aortic valve reintervention	0.2% (1, 1)	0.2% (1, 1)	0.8% (3, 3)

*Observed rate - % (no. of events, no. of subjects with the event)
[†]Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)

Note: At the time of this extract, there is one stroke and one aortic valve reintervention that are pending adjudication.

Site Reported Adverse Events

The site reported adverse events at discharge, 30 days and one year for the bicuspid population are shown in Table 63.

Table 63: Site Reported Adverse Events - Bicuspid Population (Attempted Implant Population)

Events	Discharge*	30 Days [†]	1 Year [†]
Non-valve related readmission	N/A [‡]	8.7% (50, 45)	26.8% (164, 110)
Conduction/native pacer disturbance req pacer	7.3% (40, 40)	8.6% (46, 46)	9.7% (50, 50)
Minor vascular complication	4.6% (25, 25)	5.0% (28, 27)	5.0% (28, 27)
Unplanned vascular surgery or intervention	3.5% (19, 19)	3.5% (19, 19)	3.8% (20, 20)
Cardiac arrest	2.9% (16, 16)	3.0% (16, 16)	3.0% (16, 16)
Atrial fibrillation	2.4% (13, 13)	2.4% (13, 13)	2.4% (13, 13)
Hematoma at access site	2.2% (12, 12)	2.2% (12, 12)	2.2% (12, 12)
Ischemic Stroke	1.5% (8, 8)	2.0% (11, 11)	2.6% (13, 13)
Other bleed	2.0% (12, 11)	2.0% (12, 11)	2.0% (12, 11)
Unplanned other cardiac surgery or intervention	1.3% (7, 7)	1.9% (10, 10)	3.3% (16, 15)
Bleeding at access site	1.8% (10, 10)	1.8% (10, 10)	1.8% (10, 10)
Major vascular complication	1.1% (6, 6)	1.1% (6, 6)	1.1% (6, 6)
Perforation with or w/o tamponade	1.1% (6, 6)	1.1% (6, 6)	1.1% (6, 6)
Myocardial infarction	0.7% (4, 4)	0.9% (5, 5)	1.5% (7, 7)
Percutaneous coronary intervention (PCI)	0.7% (5, 4)	0.9% (6, 5)	1.2% (7, 6)
Valve Related Readmission	N/A [‡]	0.8% (5, 4)	2.1% (12, 9)
Coronary Compression or Obstruction	0.7% (4, 4)	0.7% (4, 4)	0.7% (4, 4)
New requirement for dialysis	0.4% (2, 2)	0.6% (3, 3)	0.9% (4, 4)
Major Bleeding Event	N/A [‡]	0.4% (2, 2)	1.3% (6, 5)
Conduction/native pacer disturbance requiring implantable cardioverter defibrillator (ICD)	0.2% (1, 1)	0.4% (2, 2)	1.6% (6, 6)
Genitourinary (GU) Bleed	0.4% (2, 2)	0.4% (2, 2)	0.4% (2, 2)
Annular Dissection	0.4% (2, 2)	0.4% (2, 2)	0.4% (2, 2)
Aortic Valve Reintervention	0.2% (1, 1)	0.2% (1, 1)	1.0% (4, 4)
Transient Ischemic Attack	0.2% (1, 1)	0.2% (1, 1)	0.2% (1, 1)
Aortic Dissection	0.2% (1, 1)	0.2% (1, 1)	0.2% (1, 1)
Device recapture or retrieval	0.2% (1, 1)	0.2% (1, 1)	0.2% (1, 1)
Retroperitoneal bleeding	0.2% (1, 1)	0.2% (1, 1)	0.2% (1, 1)
Endocarditis	0.0% (0, 0)	0.0% (0, 0)	0.6% (2, 2)

Events	Discharge*	30 Days†	1 Year†
Device Thrombosis	0.0% (0, 0)	0.0% (0, 0)	0.2% (1, 1)
Undetermined Stroke	0.0% (0, 0)	0.0% (0, 0)	0.3% (1, 1)

*Observed rate - % (no. of events, no. of subjects with the event)
† Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)
‡ N/A = Event not collected on case report form at the time period. % (no. of events, no. of subjects with the event)

Effectiveness Endpoints

Valve Performance

The bicuspid aortic valve echocardiographic performance data are summarized in Figure 59 to Figure 61. The mean gradients 44.9 ± 15.5 mmHg at baseline to 12.0 ± 5.2 mmHg at 30 days and 13.4 ± 9.6 mmHg at one year. Moderate/severe PVL was observed in 4.8% of the patients at 30 days and 5.1% of the patients at one year.




Figure 59: Aortic Mean Gradient - Bicuspid Population (Valve Implant Population)

Note: Line plot with mean and standard deviation. The total number of patients at each time point only counted the patients with valid values.




Figure 60: Aortic Regurgitation - Bicuspid Population (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.




Figure 61: Aortic Paravalvular Leak - Bicuspid Population (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.

NYHA Class

The NYHA class distributions at baseline, 30-day visit and one-year visit and the NYHA class changes from baseline to the 30-day visit and to one-year visit are shown in Figure 62 and Table 64, respectively. The majority (84.0% and 82.5%) of the patients had an improved NYHA class at the 30-day visit and one year visit, respectively.




Figure 62: NYHA Functional Class - Bicuspid Population (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.

Table 64: NYHA Changes - Bicuspid Population (Valve Implant Population)

	NYHA Class Change*		
	Improved	Same	Worsened
Baseline to 30-day visit	340/405 (84.0%)	54/405 (13.3%)	11/405 (2.7%)
Baseline to 1-year visit	208/252 (82.5%)	32/252 (12.7%)	12/252 (4.8%)

*n/Total no. (%); the total no. only counted the patients with valid values.

Five Meter Walk Test

The results of the five-meter walk test are summarized in Table 65.

Table 65: Five-Meter Walk Test - Bicuspid Population (Valve Implant Population)

Visit	Five Meter Walk Time (seconds)*
Baseline	8.0 ± 4.8 (411)
30-day visit	6.7 ± 2.8 (119)
Change from baseline to 30-day visit	-1.2 ± 3.5 (101)
1-year visit	6.2 ± 2.4 (43)
Change from baseline to 1-year visit	-1.6 ± 3.6 (35)

*Mean ± SD (Total no.). The total number of patients at each time point only counted the patients with valid values.

Length of Stay

The mean index hospitalization stay was 4.7 days, which included an average of 1.6 days in the intensive care unit (ICU), as summarized in Table 66.

Table 66: Index Hospitalization Stay - Bicuspid Population (Attempted Implant Population)

	Length (days)*
Index hospitalization duration (day)	4.7 ± 3.8 (545)
Intensive care stay (day)	1.6 ± 2.6 (537)

	Length (days)*
*Mean ± SD (Total no.).	

Quality of Life (QoL)

The QoL at baseline, 30 days and one year as measured by the KCCQ overall summary score, is shown in Figure 63. The mean KCCQ summary score improved from 44.0 at baseline to 77.7 at one year.




Figure 63: KCCQ Overall Summary - Bicuspid Population (Valve Implant Population)

Note: Line plot with mean and standard deviation. The total number of patients at each time point only counted the patients with valid values.

Procedural Information

The procedure information is presented in Table 67. The most common delivery approach for the bicuspid population was the transfemoral approach, which was used in 94.7% (516/545) of cases, followed by the transapical and transaortic in 1.3% (7/545) and 1.3% (7/545) of cases, respectively, and other alternative approaches (subclavian, transcarotid, and other) in 2.8% (15/545). The device was successfully implanted in 539/544 (99.1%) of patients; five patients were converted to open heart surgery 0.9% (5/545) due to ventricular rupture (1 patient), annulus rupture (1 patient), coronary occlusion (1 patient) and other (2 patients). There were no cases of valve embolization. Device implant success is defined as correct positioning of a single prosthetic heart valve in the proper anatomical location.

Table 67: Procedural Data Summary - Bicuspid Population (Attempted Implant Population)

	Summary Statistics*
Operator reason for procedure	
Inoperable/Extreme risk	115/545 (21.1%)
High risk	394/545 (72.3%)
Intermediate risk	24/545 (4.4%)
Low risk	12/545 (2.2%)
Implant approach	
Transfemoral	516/545 (94.7%)
Transapical	7/545 (1.3%)
Transaortic	7/545 (1.3%)
Subclavian/axillary	9/545 (1.7%)
Transcarotid	3/545 (0.6%)

	Summary Statistics*
Other	3/545 (0.6%)
Procedure status	
Elective	497/545 (91.2%)
Urgent	47/545 (8.6%)
Emergency	1/545 (0.2%)
Valve size	
20 mm	16/545 (2.9%)
23 mm	95/545 (17.4%)
26 mm	220/545 (40.4%)
29 mm	214/545 (39.3%)
Primary procedure indication	
Aortic stenosis (Primary)	535/545 (98.2%)
Aortic insufficiency (Primary)	1/545 (0.2%)
Mixed aortic stenosis/aortic insufficiency	9/545 (1.7%)
Cardiopulmonary bypass (CPB)	5/545 (0.9%)
CPB status	
Elective	2/5 (40.0%)
Emergent	3/5 (60.0%)
CPB time (min)	52.2 ± 24.2 (5)
Type of anesthesia	
General anesthesia	389/545 (71.4%)
Moderate sedation	151/545 (27.7%)
Epidural	1/545 (0.2%)
Combination	4/545 (0.7%)
Total procedure time (min)	109.3 ± 49.5 (545)
Fluoroscopy time (min)	19.9 ± 10.4 (528)
Device implanted successfully	539/544 (99.1%)
Procedure aborted	0/545 (0.0%)
Conversion to open heart surgery	5/545 (0.9%)
Ventricular rupture	1/5 (20.0%)
Annulus rupture	1/5 (20.0%)
Coronary occlusion	1/5 (20.0%)
Other	2/5 (40.0%)
Mechanical assist device in place at start of procedure	5/545 (0.9%)
Intra-aortic balloon pump (IABP)	3/5 (60.0%)
Catheter based assist device	2/5 (40.0%)
*Categorical measures – no./Total no. (%); continuous measures - mean ± SD (Total no.). The total no. only counted the patients with valid values at the time point.	

SAPIEN 3 THV IN BICUSPID AORTIC VALVE FOR PATIENTS AT LOW SURGICAL RISK – PARTNER 3 BICUSPID REGISTRY ANALYSIS

A. Patient Accountability

At the time of the database lock, a total of 75 patients were enrolled in the registry.

There were three different analysis populations reflective of the single-arm study design: All Enrolled, Attempted Implant(AI), and Valve Implant (VI), which are defined in Table 68.

Table 68: Analysis Populations

Analysis Population	Definition	Number of Patients
All Enrolled	All patients who were approved by the Case Review Board and enrolled in the study.	75
Attempted Implant (AI)	All enrolled patients in whom the index procedure had begun, whether or not the procedure was completed.	71
Valve Implant (VI)	All AI patients who received and retained the study valve during the index procedure.	71

The overall follow-up compliance of the registry is summarized in Table 69.

Table 69: Bicuspid Registry Study Compliance (AI Population)

Patient Disposition and Visit Status	30-Day Visit	1-Year Visit
Total patients	71	71
Ineligible*	0	1
Death	0	1
Withdrawn	0	0
Lost to follow-up	0	0
Exit with other reason	0	0
Eligible [†]	71	70
Follow-up visit completed	71 (100.0%)	69 (97.2%)

*Ineligible includes patients who exited the study prior to the visit or had a pending visit status
[†]Eligible patients = analysis population N – ineligible patients.

B. Study Population Demographics and Baseline Characteristics

The demographics and baseline characteristics of the study population are typical for a TAVR study performed in the U.S., as shown in Table 70.

Table 70: Patient Demographics and Baseline Characteristics (AI Population)

Demographics and Baseline Characteristics	Summary Statistics* (N = 71)
Age (years)	68.5 ± 6.6
Gender	
Male	69.0% (49/71)
Female	31.0% (22/71)
Society of Thoracic Surgeons (STS) score (%)	1.4 ± 0.59
New York Heart Association (NYHA) class	
I/II	76.1% (54/71)
III/IV	23.9% (17/71)
Previous myocardial infarction	7.0% (5/71)
Previous intervention	
Coronary artery bypass grafting (CABG)	0.0% (0/71)
Percutaneous coronary intervention (PCI)	7.0% (5/71)
Stroke or cerebrovascular accident (CVA)	2.8% (2/71)
Peripheral vascular disease (PVD)	5.6% (4/71)
Atrial fibrillation	4.2% (3/71)
Permanent pacemaker or defibrillator	1.4% (1/71)

Demographics and Baseline Characteristics	Summary Statistics* (N = 71)
Hostile chest	0.0% (0/71)
Echocardiographic findings	
Valve area (cm ²)	0.7 ± 0.2 (65)
Mean gradient (mmHg)	56.1 ± 15.5 (71)
Mean left ventricular ejection fraction (LVEF), %	64.0 ± 11.1 (66)
Moderate or severe aortic regurgitation	5.6% (4/71)
Moderate or severe mitral regurgitation	1.4% (1/70)
Morphology	
Bicuspid Sievers Type 0	14.1% (10/71)
Bicuspid Sievers Type 1	84.5% (60/71)
Bicuspid Sievers Type 2	1.4% (1/71)
* Continuous measures - mean ± SD (Total no.); Categorical measures - % (no./Total no.)	

C. Safety and Effectiveness Results

1. Primary Endpoint

The primary endpoint results are presented in Table 71 and Figure 64. The rate of all-cause death, all stroke, and rehospitalization (valve-related or procedure-related and including heart failure) at 1 year was 8.5%.

Table 71: Primary Endpoint Analysis (AI Population)

Endpoint	Kaplan-Meier Rate* (N = 71)
All-cause death, all stroke, or rehospitalization	8.5% (7, 6)
All-cause death	1.4% (1, 1)
All stroke	2.8% (2, 2)
Rehospitalization	5.6% (4, 4)
*Kaplan-Meier rate (no. of events, no. of patients with the event).	




Figure 64: All-Cause Death, All Stroke, and Rehospitalization through 1 Year (AI Population)

Note: Vertical bar represents 95% confidence interval. The confidence intervals are calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

2. Secondary Endpoints

A summary of the secondary endpoints is shown in Table 72. No formal statistical tests were performed.

Table 72: Secondary Endpoints Analysis (AI Population)

No.	Endpoint	Summary Statistics*
1	New onset atrial fibrillation at 30 days	4.2% (3/71)
2	Length of index hospitalization (days)	2.6 ± 0.1 (71)
3	Death, KCCQ < 45 or KCCQ decrease from baseline ≥ 10 points at 30 days	1.4% (1/71)
4	Death or all stroke at 30 days	2.8% (2/71)
5	All stroke at 30 days	2.8% (2/71)

*Continuous measures - mean ± SE (Total no.); Categorical measures – observed rate, % (no./Total no.)

Valve Performance

The increase in EOA and decrease in gradient were sustained through 1 year. The proportion of patients with total AR ≥ moderate was 1.4% at 30 days and 1.5% at 1 year. The proportion of patients with paravalvular regurgitation ≥ moderate was 1.4% at 30 days and 0.0% at 1 year. Prosthesis patient mismatch (PPM) was measured per VARC-2. The percentages of patients with moderate and severe PPM were 21.5% and 9.2%, respectively. However, no deaths were reported in these patients at 1 year. It should be noted that there are limitations in echocardiogram measurement of PPM in TAVR patients, such as errors in stroke volume calculation, inherent limitations in the Bernoulli equation for deriving post-implant gradients, which will overstate the degree of PPM. Relying on the echocardiogram-derived iEOA alone to estimate PPM without taking into consideration the aortic valve velocity and gradient will lead to conflicting PPM assessment and outcomes post-TAVR^[4]. Accounting for pressure recovery, using

predicted versus measured iEOA, and using more accurate CT-LVOT area measurement may be required in the assessment of PPM and prosthetic valve function. Thus, interpretation of the PPM rate should take into consideration the current methodological limitations in PPM measurement and the totality of the clinical outcomes data of the device.




Figure 65: Effective Orifice Area (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 66: Mean Aortic Gradient (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 67: Total Aortic Regurgitation (VI Population)




Figure 68: Paravalvular Regurgitation (VI Population)

New York Heart Association (NYHA) Class

The NYHA classifications by visit are presented in Figure 69. At baseline, 23.9% of patients were in NYHA class III/IV. At 1 year, all patients (100.0%) were in NYHA class I/II.




Figure 69: NYHA Class by Visit (VI Population)

Six-Minute Walk Test (6MWT)

The results for the 6MWT are presented in Figure 70. The patients showed an increase in mean 6MWT distance from 361.6 meters at baseline to 376.4 meters at 30 days and 422.4 meters at 1 year.




Figure 70: 6MWT Distance (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

Quality of Life

KCCQ

The results for the KCCQ overall summary score are presented in Figure 71. The mean score increased from 71.0 at baseline to 91.8 at 30 days and 90.8 at 1 year.




Figure 71: KCCQ Overall Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

EuroQol (EQ-5D)

The results for the EQ-5D visual analog score (VAS) are presented in Figure 72. The mean score was 78.3 at baseline, 87.7 at 30 days, and 86.1 at 1 year.




Figure 72: EQ-5D Visual Analog Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

Short Form (SF)-36

The results for the SF-36 physical component summary score and mental component summary score are presented in Figure 73 and Figure 74, respectively.




Figure 73: SF-36 Physical Component Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 74: SF-36 Mental Component Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.

3. Adverse Events

The Kaplan-Meier estimates of the CEC-adjudicated adverse events through 1 year are shown in Table 73.

Table 73: CEC-Adjudicated Adverse Events through 1 Year (AI Population)

Event	Kaplan-Meier Rate*	
	30 Days	1 Year
All-cause death	0.0% (0, 0)	1.4% (1, 1)
Cardiovascular death	0.0% (0, 0)	1.4% (1, 1)
All stroke	2.8% (2, 2)	2.8% (2, 2)
Disabling stroke	0.0% (0, 0)	0.0% (0, 0)
Non-disabling stroke	2.8% (2, 2)	2.8% (2, 2)
Death or stroke	2.8% (2, 2)	2.8% (3, 2)
Death or disabling stroke	0.0% (0, 0)	1.4% (1, 1)
Major vascular complications	0.0% (0, 0)	0.0% (0, 0)
Life-threatening / disabling, or major bleeding	1.4% (1, 1)	1.4% (1, 1)
Life-threatening /disabling bleeding	0.0% (0, 0)	0.0% (0, 0)
Major bleeding	1.4% (1, 1)	1.4% (1, 1)
Myocardial infarction	0.0% (0, 0)	0.0% (0, 0)
Requirement for renal replacement [†]	0.0% (0, 0)	0.0% (0, 0)
New permanent pacemaker implantation resulting from new or worsened conduction disturbances [‡]	9.9% (7, 7)	11.3% (8, 8)
Coronary obstruction requiring intervention	0.0% (0, 0)	0.0% (0, 0)
New onset atrial fibrillation	2.9% (2, 2)	4.4% (3, 3)
Rehospitalization [§]	4.2% (3, 3)	5.6% (4, 4)
<p>*Kaplan-Meier rate (no. of events, no. of patients with the event). [†]Requirement for renal replacement was based on the site-reported event. All the other events were based on the CEC-adjudicated results. [‡]Patients with pacemaker or ICD at baseline were not counted as new events. [§]Rehospitalization (valve-related or procedure-related and including heart failure).</p>		

4. Other Study Observations

Procedural Information

The procedural data are summarized in Table 74. The mean procedure time was 57.5 minutes. Conscious sedation was used in the majority of patients (62.0%).

Table 74: Procedure (AI Population)

Variable	Summary Statistics* (N = 71)
Procedure time (min)	57.5 ± 29.4 (71)
Anesthesia type	
General	36.6% (26/71)
Conscious sedation	62.0% (44/71)
Conversion from conscious sedation to general anesthesia during the procedure	1.4% (1/71)
Anesthesia time (min)	140.8 ± 47.0 (71)
Concomitant Procedures	12.7% (9/71)

Variable	Summary Statistics* (N = 71)
Procedure aborted	0.0% (0/71)
Conversion from TAVR to SAVR	0.0% (0/71)
Valve size	
20 mm	0.0% (0/71)
23 mm	28.2% (20/71)
26 mm	43.7% (31/71)
29 mm	28.2% (20/71)
Successful access, delivery and retrieval of the device delivery system	100.0% (71/71)
Specify arterial access	
Left percutaneous	17.4% (12/69)
Right percutaneous	82.6% (57/69)
Total fluoroscopy time (min)	17.0 ± 11.2 (71)
BAV performed	54.9% (39/71)
Post dilatation performed	19.7% (14/71)
Number of post dilatations	
1	92.9% (13/14)
2	7.1% (1/14)
3	0.0% (0/14)
More than one SAPIEN 3 THV implanted	0.0% (0/71)
*Continuous measures - mean ± SD (n); Categorical measures - % (no./Total no.).	

SAPIEN 3 THV VALVE-IN-VALVE – STS/ACC TRANSCATHETER VALVE THERAPY REGISTRY (TVTR) ANALYSIS

Patient Accountability

At the time of database extract, of the 314 patients in the aortic valve-in-valve cohort, 299 patients were eligible for the 30-day visit, and 252 (84.3%) patients paid a visit within the 30-day follow-up window defined as the period between the discharge + 1 day or 21 days post-procedure (whichever occurred first) and 75 days post-procedure; of the 311 patients (SAPIEN XT and SAPIEN 3 valve patients combined) in the mitral valve-in-valve cohort, 290 patients were eligible for the 30-day visit, and 244 (84.1%) patients paid a visit within the 30-day follow-up window. A detailed summary of the patient accountability at 30 days for the two cohorts is shown in Table 75.

Table 75: Patient Accountability at 30-Day Follow-Up Visit

	Aortic Valve-in-Valve	Mitral Valve-in-Valve		
		SAPIEN XT	SAPIEN 3	All
Total patients	314	241	70	311
Non-eligible	15	15	6	21
-Death	11	15	4	19
-Withdrawal	0	0	0	0
-Lost to follow-up	1	0	2	2
-Visit not yet due	3	0	0	0
Eligible	299	226	64	290
-Follow-up visit completed	252 (84.3%)	196 (86.7%)	48 (75.0%)	244 (84.1%)
-Missed Visit	47 (15.7%)	30 (13.3%)	16 (25.0%)	46 (15.9%)

The "Attempted Implant" population consisted of all patients for whom the first vascular access was attempted. The "Valve Implant" population consisted of those patients for whom the valve implant procedure has started and a "No" was indicated for both "procedure aborted" and "conversion to open heart surgery." The number of patients in each analysis population of the aortic valve-in-valve and mitral valve-in-valve cohorts is shown in Table 76.

Table 76: Analysis Populations

Analysis Population	Aortic Valve-in-Valve	Mitral Valve-in-Valve		
		SAPIEN XT	SAPIEN 3	All
All Enrolled population	314	241	70	311
Attempted Implant population	314	241	70	311
Valve Implant population	314	236	69	305

Study Population Demographics and Baseline Characteristics

The demographics and baseline characteristics of both the aortic and mitral valve-in-valve patients, as shown in Table 77 and Table 78, present an elderly, multimorbid cohort of patients, consistent with the high operative risk of the populations.

Table 77: Patient Demographics and Baseline Characteristics - Aortic Valve-in-Valve (Attempted Implant Population)

Demographics and Baseline Characteristics	Summary Statistics*
Age – years	74.3 ± 12.10 (313)
Male sex	188/314
Society of Thoracic Surgeons (STS) score	9.0 ± 8.0 (304)
New York Heart Association (NYHA) class	
I/II	45/312 (14.4%)
III/IV	267/312 (85.6%)
Previous myocardial infarction	62/313 (19.8%)
Previous intervention	
Coronary artery bypass grafting (CABG)	119/314 (37.9%)
Percutaneous coronary intervention (PCI)	56/314 (17.8%)
Prior aortic valvuloplasty	10/306 (3.3%)
Cerebrovascular accident (CVA)	46/313 (14.7%)
Peripheral vascular disease	79/314 (25.2%)
Atrial fibrillation	126/314 (40.1%)
Permanent pacemaker	53/314 (16.9%)
Porcelain aorta	19/314 (6.1%)
Hostile chest	58/314 (18.5%)
Echocardiographic findings (Valve Implant Population)	
Valve area - cm ²	0.8 ± 0.4 (230)
Mean aortic-valve gradient – mmHg	39.3 ± 15.8 (251)
Mean left ventricular ejection fraction (LVEF)%	52.2 ± 13.1 (308)
Moderate or severe aortic regurgitation	168/310 (54.2%)
Moderate or severe mitral regurgitation	126/261 (48.3%)

*Continuous measures - Mean ± SD (Total no.); Categorical measures - n. / Total no. (%)

Table 78: Patient Demographics and Baseline Characteristics - Mitral Valve-in-Valve (Attempted Implant Population)

Demographics and Baseline Characteristics	Summary Statistics*		
	SAPIEN XT	SAPIEN 3	All
Age - years	73.9 ± 12.4 (241)	71.5 ± 15.0 (70)	73.4 ± 13.1 (311)
Male sex	88/241 (36.5%)	32/70 (45.7%)	120/311 (38.6%)
Society of Thoracic Surgeons (STS) score	13.2 ± 9.1 (237)	12.2 ± 8.7 (65)	13.0 ± 8.98 (302)
New York Heart Association (NYHA) class			

Demographics and Baseline Characteristics	Summary Statistics*		
	SAPIEN XT	SAPIEN 3	All
I/II	30/238 (12.6%)	3/70 (4.3%)	33/308 (10.7%)
III/IV	208/238 (87.4%)	67/70 (95.7%)	275/308 (89.3%)
Previous myocardial infarction	47/239 (19.7%)	18/70 (25.7%)	65/309 (21.0%)
Previous intervention			
Coronary artery bypass grafting (CABG)	93/236 (39.4%)	28/69 (40.6%)	121/305 (39.7%)
Percutaneous coronary intervention (PCI)	32/238 (13.4%)	9/69 (13.0%)	41/307 (13.4%)
Cerebrovascular accident (CVA)	45/241 (18.7%)	15/70 (21.4%)	60/311 (19.3%)
Peripheral vascular disease	42/239 (17.6%)	6/70 (8.6%)	48/309 (15.5%)
Atrial fibrillation/flutter	155/241 (64.3%)	50/70 (71.4%)	205/311 (65.9%)
Permanent pacemaker	74/240 (30.8%)	20/69 (29.0%)	94/309 (30.4%)
Porcelain aorta	6/240 (2.5%)	1/69 (1.4%)	7/309 (2.3%)
Hostile chest	41/241 (17.0%)	6/70 (8.6%)	47/311 (15.1%)
Echocardiographic findings (Valve Implant Population)			
Mitral valve area - cm ²	1.5 ± 0.9 (153)	1.4 ± 1.0 (46)	1.5 ± 0.88 (199)
Mean mitral-valve gradient - mmHg	12.7 ± 5.5 (215)	13.7 ± 6.2 (65)	12.9 ± 5.65 (280)
Mean left ventricular ejection fraction (LVEF), %	54.4 ± 11.7 (230)	53.8 ± 13.9 (67)	54.3 ± 12.2 (297)
Moderate or severe aortic regurgitation	35/231 (15.2%)	7/67 (10.5%)	42/298 (14.1%)
Moderate or severe mitral regurgitation	149/233 (63.9%)	39/68 (57.4%)	188/301 (62.5%)
*Continuous measures - Mean ± SD (Total no.); categorical measures - n. / Total no. (%). The total no. only counted the patients with valid values.			

Safety and Effectiveness Results

Aortic Valve-in-Valve

Safety Endpoints

The mortality rates at discharge and 30 days and the Kaplan-Meier curve for all-cause mortality for the aortic valve-in-valve cohort are shown in Table 79 and Figure 75, respectively. There were a total of 12 deaths reported at 30 days.

Table 79: Death Rate - Aortic Valve-in-Valve (Attempted Implant Population)

	Discharge*	30 Days†
All-cause death	2.5% (8)	4.5% (12)
Cardiac death	1.3% (4)	2.2% (6)
*Observed rate - % (n)		
†Kaplan-Meier estimate - % (n)		




Figure 75: All-Cause Death Rate - Aortic Valve-in-Valve (Attempted Implant Population)

The DCRI adjudicated events, including all strokes/TIAs and aortic valve reinterventions at discharge and 30 days for the aortic valve-in-valve cohort, are shown in Table 80.

Table 80: Duke Clinical Research Institute Adjudicated Events - Aortic Valve-in-Valve (Attempted Implant Population)

Events	Discharge*	30 Days†
All stroke	1.0% (3, 3)	1.0% (3, 3)
Ischemic stroke	1.0% (3, 3)	1.0% (3, 3)
Hemorrhagic stroke	0.0% (0, 0)	0.0% (0, 0)
Transient ischemic attack (TIA)	0.0% (0, 0)	0.0% (0, 0)
Aortic valve reintervention	0.3% (1, 1)	0.3% (1, 1)

*Observed rate - % (no. of events, no. of subjects with the event)
†Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)

Site Reported Adverse Events

The site reported adverse events at discharge and 30 days for the aortic valve-in-valve cohort are shown in Table 81.

Table 81: Site Reported Adverse Events - Aortic Valve-in-Valve (Attempted Implant Population)

Events	Discharge*	30 Days†
Non-valve related readmission	N/A	5.9% (15, 15)
Minor vascular complication	3.8% (12, 12)	4.3% (13, 13)
Conduction/native pacer disturbance requiring pacer	2.9% (9, 9)	3.0% (9, 9)
Hematoma at access site	2.9% (9, 9)	2.9% (9, 9)
Atrial fibrillation	2.5% (8, 8)	2.6% (8, 8)
Bleeding at access site	2.5% (8, 8)	2.5% (8, 8)

Events	Discharge*	30 Days†
Cardiac arrest	2.5% (8, 8)	2.5% (8, 8)
Unplanned vascular surgery or intervention	1.6% (5, 5)	2.0% (7, 6)
Percutaneous coronary intervention (PCI)	1.3% (4, 4)	1.7% (5, 5)
Other bleed	1.3% (4, 4)	1.3% (4, 4)
Coronary compression or obstruction	1.0% (3, 3)	1.0% (3, 3)
Hemorrhagic stroke	0.6% (2, 2)	1.1% (3, 3)
Life threatening bleeding	N/A	1.1% (3, 3)
Unplanned other cardiac surgery or intervention	1.0% (3, 3)	1.0% (3, 3)
Major bleeding event	N/A	0.8% (2, 2)
Major vascular complication	0.6% (2, 2)	0.6% (3, 2)
Myocardial infarction	0.3% (1, 1)	0.7% (2, 2)
New requirement for dialysis	0.6% (2, 2)	0.8% (2, 2)
Other device related event	0.6% (2, 2)	0.6% (2, 2)
Aortic valve re-intervention	0.0% (0, 0)	0.4% (1, 1)
Conduction/native pacer disturbance requiring implantable cardioverter defibrillator (ICD)	0.3% (1, 1)	0.3% (1, 1)
Device migration	0.3% (1, 1)	0.3% (1, 1)
Gastrointestinal bleeding (GI) bleed	0.3% (1, 1)	0.3% (1, 1)
Transapical related event	0.3% (1, 1)	0.3% (1, 1)
Valve related readmission	N/A	0.4% (1, 1)
Device thrombosis	0.0% (0, 0)	0.0% (0, 0)
* Observed rate - % (no. of events, no. of subjects with the event)		
† Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)		

Effectiveness Endpoints

Valve Performance

The aortic valve-in-valve echocardiographic performance data are summarized in Figure 76 to Figure 78. The mean gradients improved from 39.3 ± 15.8 mmHg at baseline to 21.5 ± 11.3 mmHg at 30 days. Moderate/severe aortic regurgitation was observed in 54.2% of the patients at baseline, which decreased to 1.5% of the patients at 30 days.




Figure 76: Mean Gradient by Visit - Aortic Valve-in-Valve (Valve Implant Population)

Note: Line plot with mean and standard deviation. The total number of patients at each time point only counted the patients with valid values.




Figure 77: Aortic Regurgitation by Visit - Aortic Valve-in-Valve (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.




Figure 78: Paravalvular Regurgitation by Visit - Aortic Valve-in-Valve (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.

NYHA Class

The NYHA class distributions at baseline and the 30-day visit and the NYHA class changes from baseline to the 30-day visit are shown in Figure 79 and Table 82, respectively. The majority (85.4%) of the patients had an improved NYHA class at the 30-day visit.




Figure 79: NYHA Functional Class - Aortic Valve-in-Valve (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.

Table 82: NYHA Class Change - Aortic Valve-in-Valve (Valve Implant Population)

	NYHA Class Change*		
	Improved	Same	Worsened
Baseline to 30-day visit	193/226 (85.4%)	31/226 (13.7%)	2/226 (0.9%)

*n/Total no. (%); the total no. only counted the patients with valid values.

Five-Meter Walk Test

The results of the five-meter walk test are summarized in Table 83.

Table 83: Five-Meter Walk Test - Aortic Valve-in-Valve (Valve Implant Population)

Visit*	Five Meter Walk Time (seconds)†
Baseline	7.6 ± 3.9 (209)
30-day visit	5.9 ± 2.4 (68)
Change from baseline to 30-day visit	-1.4 ± 2.9 (51)

*There were up to 3 five-meter walk tests for each patient at each visit, and the results were averaged.
†Mean ± SD (Total no.). The total number of patients at each time point only counted the patients with valid values.

Length of Stay


The mean index hospitalization stay was 4.9 days, which included an average of 1.8 days in the intensive care unit (ICU), as summarized in Table 84.

Table 84: Index Hospitalization Stay - Aortic Valve-in-Valve (Attempted Implant Population)

	Length (days)*
Index Hospitalization Stay	4.9 ± 3.9 (314)
Intensive Care Stay	1.8 ± 2.6 (311)
*Mean ± SD (Total no.).	

Quality of Life (QoL)

The QoL at baseline and 30 days as measured by the KCCQ clinical summary score is shown in Figure 80. The mean KCCQ summary score improved from 39.4 at baseline to 75.3 at 30 days.

**Figure 80: KCCQ Overall Summary Score - Aortic Valve-in-Valve (Valve Implant Population)**

Note: Line plot with mean and standard deviation. The total number of patients at each time point only counted the patients with valid values.

Procedural Information

The procedure information is presented in Table 85. The most common delivery approach for the aortic valve-in-valve implantation was the transfemoral approach, which was used in 93.0% (292/314) of cases, followed by the transapical approach in 4.1% (13/314) of cases, and other alternative approaches (transaortic, subclavian, and other) in 2.9% (9/314) of cases. There were no aborted procedures or conversions to open heart surgery. The overall device success rate was 88.9% (272/306), which was defined as the following:

- Successful vascular access, delivery, and deployment of the device and successful retrieval of the delivery system, and
- Correct position of the device in the proper anatomical location, and
- Intended performance of the prosthetic heart valve (aortic valve area > 1.2 cm² and mean aortic valve gradient < 20 mmHg or peak velocity < 3 m/s, without moderate or severe prosthetic valve regurgitation), and
- Only one valve implanted in the proper anatomical location.

Table 85: Procedural Data Summary - Aortic Valve-in-Valve (Attempted Implant Population)

Procedural Data	Summary Statistics*
Operator Reason for Procedure	
Inoperable/extreme risk	80/313 (25.6%)
High risk	219/313 (70.0%)
Intermediate risk	10/313 (3.2%)
Low risk	4/313 (1.3%)
Implant Approach	
Transfemoral	292/314 (93.0%)

Procedural Data	Summary Statistics*
Transapical	13/314 (4.1%)
Transaortic	1/314 (0.3%)
Subclavian/axillary	6/314 (1.9%)
Other [†]	2/314 (0.6%)
Prior Valve Type	
Bioprosthetic stented	159/308 (51.6%)
Bioprosthetic stentless	79/308 (25.6%)
Procedure Status	
Elective	231/314 (73.6%)
Urgent	74/314 (23.6%)
Emergency	8/314 (2.5%)
Salvage	1/314 (0.3%)
Valve Size	
20 mm	83/314 (26.4%)
23 mm	130/314 (41.4%)
26 mm	57/314 (18.2%)
29 mm	44/314 (14.0%)
Primary Procedure Indication	
Aortic stenosis (Primary)	95/313 (30.4%)
Aortic insufficiency (Primary)	19/313 (6.1%)
Mixed aortic stenosis/aortic insufficiency	10/313 (3.2%)
Failed bioprosthetic valve	189/313 (60.4%)
Cardiopulmonary Bypass (CPB)	
CPB status	
Elective	4/5 (80.0%)
Emergent	1/5 (20.0%)
CPB time (min)	90.5 ± 140.9 (4)
Type of Anesthesia	
General anesthesia	240/314 (76.4%)
Moderate sedation	72/314 (22.9%)
Epidural	0/314 (0.0%)
Combination	2/314 (0.6%)
Total procedure time (min)	110.7 ± 63.0 (314)
Fluoroscopy time (min)	21.2 ± 16.1 (304)
Device success	
Procedure aborted	0/314 (0.0%)
Conversion to open heart surgery	0/314 (0.0%)
Mechanical assist device in place at start of procedure	
Intra-aortic balloon pump (IABP)	2/5 (40.0%)
Catheter based assist device	3/5 (60.0%)
*Categorical measures – no./Total no. (%); continuous measures - mean ± SD (Total no.). The total no. only counted the patients with valid values at the time point.	
†The data collection form was changed in February 2013 to specify non-transfemoral (non-TF), non-transapical (non-TA) approaches rather than “other”; hence, “other” likely included the non-TF and non-TA approaches.	

Mitral Valve-in-Valve


Safety Endpoints

The mortality rates at discharge and 30 days and the Kaplan-Meier curve for all-cause mortality for the mitral valve-in-valve cohort are shown in Table 86 and Figure 81, respectively. There were 16 reported deaths in the SAPIEN XT valve patients and 4 in the SAPIEN 3 valve patients at 30 days.

Table 86: Death Rate - Mitral Valve-in-Valve (Attempted Implant Population)

Event	Discharge*			30 Days†		
	SAPIEN XT	SAPIEN 3	All	SAPIEN XT	SAPIEN 3	All
All-cause death	5.0% (12)	5.7% (4)	5.1% (16)	6.9% (16)	6.6% (4)	6.8% (20)
Cardiac death	3.7% (9)	4.3% (3)	3.9% (12)	4.2% (10)	4.9% (3)	4.3% (13)

*Observed rate - % (n)
†Kaplan-Meier estimate - % (n)



No. at Risk				
SAPIEN XT	241	216	209	186
SAPIEN 3	70	56	52	45
ALL	311	272	261	231

Figure 81: All-Cause Death Rate - Mitral Valve-in-Valve (Attempted Implant Population)

The DCRI-adjudicated events, including all strokes/TIAs, heart failure readmissions, and mitral valve reinterventions at discharge and 30 days, for the mitral valve-in-valve cohort are shown in Table 87.

Table 87: Duke Clinical Research Institute Adjudicated Events - Mitral Valve-in-Valve (Attempted Implant Population)

Events	Discharge*			30 Days†		
	SAPIEN XT	SAPIEN 3	All	SAPIEN XT	SAPIEN 3	All
All stroke	0.4% (1, 1)	1.4% (1, 1)	0.6% (2, 2)	0.4% (1, 1)	1.5% (1, 1)	0.7% (2, 2)

Events	Discharge*			30 Days†		
	SAPIEN XT	SAPIEN 3	All	SAPIEN XT	SAPIEN 3	All
Ischemic stroke	0.4% (1, 1)	1.4% (1, 1)	0.6% (2, 2)	0.4% (1, 1)	1.5% (1, 1)	0.7% (2, 2)
Hemorrhagic stroke	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)
Transient ischemic attack (TIA)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)
Readmission - heart failure	N/A	N/A	N/A	1.0% (2, 2)	0.0% (0, 0)	0.8% (2, 2)
Mitral valve reintervention	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.5% (1, 1)	0.0% (0, 0)	0.4% (1, 1)

*Observed rate - % (no. of events, no. of subjects with the event)
†Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)

Site Reported Adverse Events

The site reported adverse events at discharge and 30 days for the mitral valve-in-valve cohort are shown in Table 88.

Table 88: Site Reported Adverse Events - Mitral Valve-in-Valve (Attempted Implant Population)

Events	Discharge*			30 Days†		
	SAPIEN XT	SAPIEN 3	All	SAPIEN XT	SAPIEN 3	All
Other bleed	5.4% (13, 13)	4.3% (3, 3)	5.1% (16, 16)	6.1% (14, 14)	4.4% (3, 3)	5.8% (17, 17)
Readmission - not cardiac	N/A	N/A	N/A	5.8% (12, 12)	0.0% (0, 0)	4.6% (12, 12)
Atrial septal defect closure following transseptal catheterization	4.6% (11, 11)	5.7% (4, 4)	4.8% (15, 15)	4.6% (11, 11)	5.7% (4, 4)	4.9% (15, 15)
Cardiac arrest	4.1% (10, 10)	2.9% (2, 2)	3.9% (12, 12)	4.2% (10, 10)	3.2% (2, 2)	4.0% (12, 12)
Unplanned other cardiac surgery or intervention	3.3% (8, 8)	0.0% (0, 0)	2.6% (8, 8)	3.8% (9, 9)	0.0% (0, 0)	3.0% (9, 9)
Atrial fibrillation	3.3% (8, 8)	1.4% (1, 1)	2.9% (9, 9)	3.4% (8, 8)	1.5% (1, 1)	2.9% (9, 9)
New requirement for dialysis	2.9% (7, 7)	1.4% (1, 1)	2.6% (8, 8)	3.0% (7, 7)	1.6% (1, 1)	2.7% (8, 8)
Bleeding at access site	2.5% (6, 6)	1.4% (1, 1)	2.3% (7, 7)	2.5% (6, 6)	1.4% (1, 1)	2.3% (7, 7)
Unplanned vascular surgery or intervention	2.5% (6, 6)	2.9% (2, 2)	2.6% (8, 8)	2.5% (6, 6)	3.2% (2, 2)	2.6% (8, 8)
Perforation with or w/o tamponade	2.1% (5, 5)	0.0% (0, 0)	1.6% (5, 5)	2.1% (5, 5)	0.0% (0, 0)	1.6% (5, 5)
Hematoma at access site	1.2% (3, 3)	0.0% (0, 0)	1.0% (3, 3)	1.3% (3, 3)	0.0% (0, 0)	1.0% (3, 3)
Minor vascular complication	1.2% (3, 3)	1.4% (1, 1)	1.3% (4, 4)	1.2% (3, 3)	1.7% (1, 1)	1.3% (4, 4)
Transapical related event	1.2% (3, 3)	0.0% (0, 0)	1.0% (3, 3)	1.2% (3, 3)	0.0% (0, 0)	1.0% (3, 3)
Transseptal related event	1.2% (3, 3)	0.0% (0, 0)	1.0% (3, 3)	1.2% (3, 3)	0.0% (0, 0)	1.0% (3, 3)
Gastrointestinal bleed	0.8% (2, 2)	1.4% (1, 1)	1.0% (3, 3)	0.9% (2, 2)	1.4% (1, 1)	1.1% (3, 3)

Events	Discharge*			30 Days†		
	SAPIEN XT	SAPIEN 3	All	SAPIEN XT	SAPIEN 3	All
Major vascular complication	0.8% (2, 2)	0.0% (0, 0)	0.6% (2, 2)	0.8% (2, 2)	0.0% (0, 0)	0.6% (2, 2)
Readmission - cardiac	N/A	N/A	N/A	0.9% (2, 2)	0.0% (0, 0)	0.8% (2, 2)
Device embolization	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.5% (1, 1)	0.0% (0, 0)	0.4% (1, 1)
Device migration	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)	0.5% (1, 1)	1.4% (1, 1)	0.7% (2, 2)
Device recapture or retrieval	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)	0.5% (1, 1)	1.4% (1, 1)	0.7% (2, 2)
Genitourinary bleed	0.4% (1, 1)	0.0% (0, 0)	0.3% (1, 1)	0.4% (1, 1)	0.0% (0, 0)	0.3% (1, 1)
Major bleeding event	N/A	N/A	N/A	0.5% (1, 1)	0.0% (0, 0)	0.4% (1, 1)
Non-valve related readmission	N/A	N/A	N/A	0.5% (1, 1)	0.0% (0, 0)	0.4% (1, 1)
Conduction/native pacer disturbance requiring pacer	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)	0.0% (0, 0)	1.5% (1, 1)	0.3% (1, 1)
Device thrombosis	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)
Endocarditis	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)
Life threatening bleeding	N/A	N/A	N/A	0.0% (0, 0)	0.0% (0, 0)	0.0% (0, 0)
Myocardial infarction	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)
Other device related event	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)	0.0% (0, 0)	1.4% (1, 1)	0.3% (1, 1)
Transient ischemic attack	0.4% (1, 1)	0.0% (0, 0)	0.3% (1, 1)	0.4% (1, 1)	0.0% (0, 0)	0.3% (1, 1)
Ischemic stroke	0.4% (1, 1)	1.4% (1, 1)	0.6% (2, 2)	0.4% (1, 1)	1.5% (1, 1)	0.7% (2, 2)
Readmission - heart failure	N/A	N/A	N/A	1.0% (2, 2)	3.8% (2, 2)	1.6% (4, 4)
*Observed rate - % (no. of events, no. of subjects with the event)						
†Kaplan-Meier estimate - % (no. of events, no. of subjects with the event)						

Effectiveness Endpoints

Valve Performance

The mitral valve-in-valve echocardiographic performance data are summarized in Figure 82 to Figure 84. The mean gradients improved from 12.9 mmHg at baseline to 7.1 mmHg at 30 days. Moderate/severe mitral regurgitation was observed in 62.5% of the patients at baseline, which decreased to 2.2% of the patients at 30 days.




Figure 82: Mean Gradient by Visit - Mitral Valve-in-Valve (Valve Implant Population)

Note: Line plot with mean and standard deviation. The total number of patients at each time point only counted the patients with valid values.




Figure 83: Mitral Regurgitation by Visit - Mitral Valve-in-Valve (Valve Implant Population)

Note: Values that are < 1.0% are not labeled in the bar chart. The total number of patients at each time point only counted the patients with valid values.




Figure 84: Paravalvular Regurgitation by Visit - Mitral Valve-in-Valve (Valve Implant Population)

Note: Values that are < 1.0% are not labeled in the bar chart. The total number of patients at each time point only counted the patients with valid values.

NYHA Class

The NYHA class distributions at baseline and the 30-day visit and the NYHA class changes from baseline to the 30-day visit are shown in Figure 85 and Table 89, respectively. The majority (85.6%) of the patients had an improved NYHA class at the 30-day visit.




Figure 85: NYHA Functional Class - Mitral Valve-in-Valve (Valve Implant Population)

Note: The total number of patients at each time point only counted the patients with valid values.

Table 89: NYHA Class Change - Mitral Valve-in-Valve (Valve Implant Population)

		NYHA Class Change*		
		Improved	Same	Worsened
Baseline to 30-day visit	SAPIEN XT	133/159 (83.6%)	24/159 (15.1%)	2/159 (1.3%)
	SAPIEN 3	40/43 (93.0%)	3/43 (7.0%)	0/43 (0.0%)
	All	173/202 (85.6%)	27/202 (13.4%)	2/202 (1.0%)

*n/Total no. (%); the total no. only counted the patients with valid values.

Six-Minute Walk Test (6MWT)

The results of the 6MWT are summarized in Table 90.

Table 90: Six-Minute Walk Test - Mitral Valve-in-Valve (Valve Implant Population)

Visit	6-Minute Walk Distance (feet)*		
	SAPIEN XT	SAPIEN 3	All
Baseline	240.5 ± 366.2 (77)	375.6 ± 370.4 (32)	280.2 ± 370.9 (109)
30-day visit	768.7 ± 480.6 (34)	977.5 ± 597.4 (8)	808.5 ± 503.7 (42)

Visit	6-Minute Walk Distance (feet)*		
	SAPIEN XT	SAPIEN 3	All
Change from baseline to 30 days	479.0 ± 471.3 (20)	457.6 ± 348.1 (5)	474.7 ± 442.9 (25)

*Mean ± SD (Total no.). The total number of patients at each time point only counted the patients with valid values. The 6-minute walk distance was counted as 0 for the 6-minute walk tests not performed due to cardiac reasons.

Length of Stay

The mean index hospitalization stay was 8.5 days, which included an average of 3.4 days in the intensive care unit (ICU), as summarized in Table 91.

Table 91: Index Hospitalization Stay - Mitral Valve-in-Valve (Attempted Implant Population)

	Length (days)*		
	SAPIEN XT	SAPIEN 3	All
Index hospitalization stay	8.8 ± 7.1 (241)	7.6 ± 7.4 (70)	8.5 ± 7.1 (311)
Intensive care stay	3.3 ± 4.8 (234)	3.7 ± 7.1 (63)	3.4 ± 5.3 (297)

*Mean ± SD (Total no.).

Quality of Life (QoL)

The KCCQ clinical summary scores at baseline and 30 days are shown in Figure 86. The mean KCCQ summary score improved from 31.6 at baseline to 68.2 at 30 days.




Figure 86: KCCQ Overall Summary Score - Mitral Valve-in-Valve (Valve Implant Population)

Note: Line plot with mean and standard deviation. The total number of patients at each time point only counted the patients with valid values.

Procedural Information

The procedure information is presented in Table 92. The most common delivery approach for the mitral valve-in-valve implantation was the transapical approach, which was used in 65.3% (203 of 311) of cases, followed by the transseptal approach in 27.0% (84 of 311) of cases, the transfemoral approach in 6.1% (19/311) of cases, and other alternative approaches in 1.6% (5 of 311) of cases. The procedures were considered elective in 71.0% (220/310) of cases, urgent in 27.7% (86/310) of cases, and emergent or salvage in 1.3% (4/310) of cases. Two (2) cases were aborted and 5 were converted to open heart surgery. Overall, the device was implanted successfully in 97.4% (303/311) of the cases, which was defined as correct positioning of a single prosthetic heart valve in the proper anatomical location.

Table 92: Procedural Data Summary - Mitral Valve-in-Valve (Attempted Implant Population)

Procedural Data	Summary Statistics*		
	SAPIEN XT	SAPIEN 3	All
Operator reason for procedure			
Inoperable/extreme risk	96/241 (39.8%)	11/69 (15.9%)	107/310 (34.5%)
High risk	141/241 (58.5%)	52/69 (75.4%)	193/310 (62.3%)
Intermediate risk	4/241 (1.7%)	5/69 (7.2%)	9/310 (2.9%)
Low risk	0/241 (0.0%)	1/69 (1.4%)	1/310 (0.3%)
Implant approach			
Transapical	192/241 (79.7%)	11/70 (15.7%)	203/311 (65.3%)
Transseptal	43/241 (17.8%)	41/70 (58.6%)	84/311 (27.0%)
Femoral artery	4/241 (1.7%)	15/70 (21.4%)	19/311 (6.1%)
Other	2/241 (0.8%)	3/70 (4.3%)	5/311 (1.6%)
Prior valve type			
Bioprosthetic stented	143/180 (79.4%)	35/41 (85.4%)	178/221 (80.5%)
Bioprosthetic stentless	37/180 (20.6%)	6/41 (14.6%)	43/221 (19.5%)
Procedure status			
Elective	173/241 (71.8%)	47/69 (68.1%)	220/310 (71.0%)
Urgent	64/241 (26.6%)	22/69 (31.9%)	86/310 (27.7%)
Emergency	2/241 (0.8%)	0/69 (0.0%)	2/310 (0.6%)
Salvage	2/241 (0.8%)	0/69 (0.0%)	2/310 (0.6%)
Valve size			
23 mm	22/241 (9.1%)	5/70 (7.1%)	27/311 (8.7%)
26 mm	93/241 (38.6%)	24/70 (34.3%)	117/311 (37.6%)
29 mm	126/241 (52.3%)	41/70 (58.6%)	167/311 (53.7%)
Cardiopulmonary bypass	25/241 (10.4%)	2/69 (2.9%)	27/310 (8.7%)
Status of CP Bypass			
Elective	20/25 (80.0%)	0/2 (0.0%)	20/27 (74.1%)
Emergent	5/25 (20.0%)	2/2 (100.0%)	7/27 (25.9%)
CP Bypass Time (min)	38.3 ± 51.2 (24)	148.0 ± 157.0 (2)	46.7 ± 65.4 (26)
Type of anesthesia			
General anesthesia	240/241 (99.6%)	68/69 (98.6%)	308/310 (99.4%)
Moderate sedation	0/241 (0.0%)	1/69 (1.4%)	1/310 (0.3%)
Epidural	0/241 (0.0%)	0/69 (0.0%)	0/310 (0.0%)
Combination	1/241 (0.4%)	0/69 (0.0%)	1/310 (0.3%)
Total procedure time (min)	143.6 ± 60.4 (240)	157.7 ± 107.2 (69)	146.7 ± 73.5 (309)

Procedural Data	Summary Statistics*		
	SAPIEN XT	SAPIEN 3	All
Fluoroscopy time (min)	23.9 ± 20.7 (223)	36.9 ± 27.3 (63)	26.8 ± 22.9 (286)
Device implanted successfully	234/241 (97.1%)	69/70 (98.6%)	303/311 (97.4%)
Procedure aborted	1/241 (0.4%)	1/70 (1.4%)	2/311 (0.6%)
Procedure aborted reason			
Navigation issue after successful access	1/1 (100.0%)	0/1 (0.0%)	1/2 (50.0%)
Other	0/1 (0.0%)	1/1 (100.0%)	1/2 (50.0%)
Procedure aborted action			
Conversion to open heart surgery	0/1 (0.0%)	1/1 (100.0%)	1/2 (50.0%)
Other	1/1 (100.0%)	0/1 (0.0%)	1/2 (50.0%)
Conversion to open heart surgery	4/241 (1.7%)	1/70 (1.4%)	5/311 (1.6%)
Tamponade/bleeding in the heart	4/4 (100.0%)	0/1 (0.0%)	4/5 (80.0%)
Other	0/4 (0.0%)	1/1 (100.0%)	1/5 (20.0%)
Mechanical assist device in place at start of procedure	9/241 (3.7%)	4/70 (5.7%)	13/311 (4.2%)
IABP	7/9 (77.8%)	3/4 (75.0%)	10/13 (76.9%)
Catheter-based assist device	2/9 (22.2%)	1/4 (25.0%)	3/13 (23.1%)
*Categorical measures – no./Total no. (%); continuous measures - mean ± SD (Total no.). The total no. only counted the patients with valid values at the time point.			

SAPIEN 3 AND SAPIEN 3 ULTRA THV THV-IN-THV - STS/ACC TRANSCATHETER VALVE THERAPY REGISTRY (TVTR) ANALYSIS

Accountability of PMA Cohorts

At the time of database extract, 242 of the 263 patients were eligible for the 30-day visit and 216 (89.3%) completed the visit within the 30-day follow-up window, defined as the period between 21 days post-procedure and 75 days post-procedure. At 1 year, 200 patients were eligible for the 1-year visit and 136 (68.0%) completed the visit within the follow-up window, defined as the period between 305 days post-procedure and 425 days post-procedure. A detailed summary of the patient accountability at 30 days and 1 year is shown in Table 93.

Table 93: Patient Visit Accountability (All Population)

	30-day Visit	1-year Visit
Total patients	263	263
Non-eligible	21	63
Death	19	47
Withdrawal	1	3
Lost to follow-up	1	13
Eligible	242	200
Follow-up visit completed	89.3% (216)	68.0% (136)
Missed visit	10.7% (26)	32.0% (64)

The "Attempted Implant" population consisted of all patients for whom the first vascular access was attempted. The "Valve Implant" population consisted of those patients for whom the valve implant procedure has started and a "No" was indicated for both "procedure aborted" and "conversion to open heart surgery". The number of patients in the analysis population is shown in Table 94.

Table 94: Analysis Populations

Analysis Population	Number of Patients
Attempted Implant Population	263
Valve Implant Population	261

Study Population Demographics and Baseline Characteristics

The demographics and baseline characteristics of the patients, as shown in Table 95 present an elderly, multimorbid cohort of patients, consistent with the high operative risk of the populations.

Table 95: Patient Demographics and Baseline Characteristics (AI Population)

Demographics and Baseline Characteristics	Summary Statistics* (N = 263)
Age - years	78.9 ± 10.5
Male sex	55.1% (145/263)
Society of Thoracic Surgeons (STS) score	10.2 ± 8.6 (242)
New York Heart Association (NYHA) class	
I/II	12.3% (32/261)
III/IV	87.7% (229/261)
Previous myocardial infarction	26.6% (70/263)
Previous intervention	
Coronary artery bypass grafting (CABG)	28.1% (74/263)
Percutaneous coronary intervention (PCI)	34.0% (89/262)
Prior aortic valvuloplasty	13.7% (36/263)
Stroke or Cerebrovascular accident (CVA)	18.3% (48/263)
Peripheral vascular disease (PVD)	32.1% (84/262)
Atrial fibrillation/flutter	48.7% (128/263)
Permanent pacemaker	32.2% (84/261)
Porcelain aorta	8.0% (21/262)
Hostile chest	8.7% (23/263)
Echocardiographic findings (Valve Implant Population)	
Valve area (cm ²)	0 ± 0.5 (115)
Mean gradient (mmHg)	29.4 ± 19.0 (135)
Mean left ventricular ejection fraction (LVEF), %	49.3 ± 15.1 (257)
Moderate or severe aortic regurgitation	79.3% (207/261)
Moderate or severe mitral regurgitation	42.1% (98/233)

*Continuous measures - Mean ± SD (Total no.); Categorical measures - % (no./Total no.)

Safety and Effectiveness Results

Safety Endpoints

The Kaplan-Meier estimates of site-reported adverse events through 1 year are presented in Table 96. The Kaplan-Meier curve for all-cause mortality is shown in Figure 87. The all-cause mortality rate was 5.8% at 30 days and 18.2% at 1 year, including a cardiovascular death rate of 2.7% at 30 days and 5.4% at 1 year. Other relatively more frequent adverse events included conduction/native pacer disturbance requiring pacer (8.2% at 30 days and 10.6% at 1 year) and valve-related readmission (4.3% at 30 days and 8.6% at 1 year).

Table 96: Site Reported Adverse Events (AI Population)

Adverse Event	Kaplan-Meier Rate*	
	30 Days (N = 263)	1 Year (N = 263)
All-cause death	5.8% (15, 15)	18.2% (45, 45)
Cardiovascular death	2.7% (7, 7)	5.4% (13, 13)
All stroke	2.3% (6, 6)	2.8% (7, 7)
Ischemic stroke	1.9% (5, 5)	2.4% (6, 6)
Undetermined stroke	0.4% (1, 1)	0.4% (1, 1)

Adverse Event	Kaplan-Meier Rate*	
	30 Days (N = 263)	1 Year (N = 263)
Transient ischemic attack (TIA)	0.8% (2, 2)	1.9% (4, 4)
Major vascular complication	0.4% (1, 1)	1.0% (2, 2)
Major bleeding	1.2% (3, 3)	2.5% (7, 5)
Myocardial infarction	0.8% (3, 2)	3.1% (6, 5)
New requirement for dialysis	0.4% (1, 1)	1.6% (3, 3)
Conduction/native pacer disturbance requiring pacer	8.2% (21, 21)	10.6% (25, 25)
Conduction/native pacer disturbance requiring implantable cardioverter defibrillator (ICD)	0.4% (1, 1)	2.2% (4, 4)
Aortic valve re-intervention	0.4% (1, 1)	1.2% (2, 2)
Unplanned other cardiac surgery or intervention	2.4% (6, 6)	4.5% (9, 9)
Unplanned vascular surgery or intervention	1.5% (4, 4)	2.6% (6, 6)
Device thrombosis	0.4% (1, 1)	1.0% (2, 2)
Valve-related readmission	4.3% (11, 11)	8.6% (22, 18)

* Kaplan-Meier rate - (no. of events, no. of patients with the event).




Figure 87: All-Cause Mortality through 1 Year (AI Population)

Note: The 95% confidence intervals were calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

Effectiveness Endpoints

Valve Performance

The echocardiographic valve performance results are shown in Figure 88 to Figure 90. The decrease in gradients were sustained through 1 year. The mean aortic gradient decreased from 29.4 mmHg at baseline to 14.4 mmHg at 30 days, which was maintained through 1 year (13.3 mmHg). Moderate or severe total aortic regurgitation was observed in 79.3% of the patients at baseline, which decreased to 4.6% at 30 days and 3.4% at 1 year. The proportion of patients with \geq moderate paravalvular regurgitation was 4.5% at 30 days and 2.6% at 1 year.




Figure 88: Mean Aortic Gradient (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each visit time point only counted the patients with valid values.




Figure 89: Total Aortic Regurgitation (VI Population)




Figure 90: Paravalvular Regurgitation (VI Population)

NYHA Functional Class

The NYHA functional class distributions by visit are presented in Figure 91. At baseline, 87.7% of patients were in NYHA III/IV. At 1 year, the majority (88.4%) of patients were in NYHA I/II.




Figure 91: NYHA Class by Visit (VI Population)

Length of Stay

The mean index hospitalization stay was 4.9 days, which included an average of 1.7 days in the intensive care unit (ICU), as summarized in Table 97.

Table 97: Index Hospitalization (AI Population)

	Length of Stay (days)*
Index hospitalization duration	4.9 ± 0.3 (263)
Intensive care stay	1.7 ± 0.2 (255)
* Mean ± SE (Total no.).	

Quality of Life

The results for the KCCQ overall summary score are presented in Figure 92. The mean score increased from 38.7 at baseline to 65.3 and 73.7 at 30 days and 1 year, respectively.




Figure 92: KCCQ Overall Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each time point only counted patients with valid values.

Other Study Observations

Procedural Information

The procedural information is summarized in Table 98. General anesthesia was used in the majority (70%) of patients. Conversion to open heart surgery occurred in two patients due to ventricular rupture and annulus rupture, respectively.

Table 98: Procedural Data Summary (AI Population)

Procedural Data	Summary Statistics*
Operator reason for procedure	
Inoperable/extreme risk	20.8% (54/259)
High risk	69.1% (179/259)
Intermediate risk	9.7% (25/259)
Low risk	0.4% (1/259)
Implant approach	
Transfemoral	95.8% (252/263)
Transapical	1.1% (3/263)
Transaortic	0.8% (2/263)
Subclavian/axillary	0.8% (2/263)
Transseptal	0.4% (1/263)
Transcarotid	1.1% (3/263)
Valve size	

Procedural Data	Summary Statistics*
20 mm	2.3% (6/263)
23 mm	35.0% (92/263)
26 mm	30.8% (81/263)
29 mm	31.9% (84/263)
Cardiopulmonary bypass	0.4% (1/263)
Cardiopulmonary bypass status	
Emergent	100.0% (1/1)
Cardiopulmonary bypass time, minutes	254.0 ± NA (1)
Type of anesthesia	
General anesthesia	70.0% (184/263)
Moderate sedation	29.7% (78/263)
Combination	0.4% (1/263)
Total procedure time, minutes	108.5 ± 4.3 (263)
Device implanted successfully	98.9% (260/263)
Procedure aborted	0.0% (0/263)
Conversion to open heart surgery	0.8% (2/263)
Ventricular rupture	1
Annulus rupture	1
Mechanical assist device in place at start of procedure	0.4% (1/263)
Catheter-based assist device	100.0% (1/1)
*Continuous measures - mean ± SE (n); categorical measures - % (no./Total no.)	

SAPIEN 3 Valve in Ring - STS/ACC Transcatheter Valve Therapy Registry (TVTR) & Mitral Implantation of Transcatheter Valves (MITRAL) Study Analysis

Patient Accountability

At the time of database extract, 205 of the 236 patients were eligible for the 30-day visit and 178 (86.8%) completed the visit within the 30-day follow-up window, defined as the period between 21 days post-procedure and 75 days post-procedure. At 1 year, 152 patients were eligible for the 1-year visit and 103 (67.8%) completed the visit within the follow-up window, defined as the period between 305 days post-procedure and 425 days post-procedure. A detailed summary of the patient accountability at 30 days and 1 year is shown in Table 99.

Table 99: Patient Visit Accountability

	30-day Visit	1-year Visit
Total patients	236	236
Non-eligible	31	84
Death	24	54
Withdrawal	3	5
Lost to follow-up	4	13
Visit not yet due	0	12
Eligible	205	152
Follow-up visit completed	86.8% (178)	67.8% (103)
Missed visit	13.2% (27)	32.2% (49)

The "Attempted Implant (AI)" population consisted of all patients in the dataset. The "Valve Implant (VI)" population consisted of those patients for whom the valve implant procedure has started and a "No" was indicated for both "procedure aborted" and "conversion to open heart surgery" in the case report form of the TVT Registry (no patients in the MITRAL study had an aborted procedure). The numbers of patients in these two analysis populations are shown in Table 100.

Table 100: Analysis Populations

Analysis Population	Number of Patients
Attempted implant population	236
Valve implant population	232

Study Population Demographics and Baseline Characteristics

Patient demographics and baseline characteristics, as shown in Table 101, present an elderly, multimorbid cohort of patients, consistent with the high operative risk of the populations.

Table 101: Patient Demographics and Baseline Characteristics (AI Population)

Demographics and Baseline Characteristics	Summary Statistics*		
	TVT Registry (N = 206)	MITRAL Study (N = 30)	Overall (N = 236)
Age - years	72.1 ± 10.3 (206)	71.7 ± 8.9 (30)	72.1 ± 10.1 (236)
Male sex	47.6% (98/206)	63.3% (19/30)	49.6% (117/236)
Society of Thoracic Surgeons (STS) score	9.4 ± 6.4 (196)	8.7 ± 4.7 (30)	9.3 ± 6.2 (226)
New York Heart Association (NYHA) class			
I/II	18.7% (38/203)	23.3% (7/30)	19.3% (45/233)
III/IV	81.3% (165/203)	76.7% (23/30)	80.7% (188/233)
Previous myocardial infarction	30.2% (62/205)	22.2% (6/27)	29.3% (68/232)
Stroke	15.0% (31/206)	13.8% (4/29)	14.9% (35/235)
Transient Ischemic Attack	8.8% (18/205)	7.1% (2/28)	8.6% (20/233)
Diabetes	31.1% (64/206)	30.0% (9/30)	30.9% (73/236)
Hypertension	88.3% (182/206)	90.0% (27/30)	88.6% (209/236)
Previous intervention			
Coronary artery bypass grafting (CABG)	46.1% (95/206)	63.3% (19/30)	48.3% (114/236)
Percutaneous coronary Intervention (PCI)	25.2% (52/206)	31.0% (9/29)	26.0% (61/235)
Atrial fibrillation/flutter	66.0% (136/206)	70.0% (21/30)	66.5% (157/236)
Permanent pacemaker	23.3% (48/206)	36.7% (11/30)	25.0% (59/236)
Porcelain aorta	1.9% (4/206)	0.0% (0/30)	1.7% (4/236)
Previous implantable cardioverter defibrillator (ICD)	23.3% (48/206)	23.3% (7/30)	23.3% (55/236)
Echocardiographic findings (Valve Implant Population)			
Mitral valve area (cm ²)	1.9 ± 0.9 (115)	2.7 ± 0.8 (30)	2.1 ± 0.9 (145)
Mitral valve mean gradient (mmHg)	8.0 ± 4.6 (171)	7.5 ± 4.8 (30)	7.9 ± 4.6 (201)
Left ventricular ejection fraction (LVEF)	47.0 ± 14.5 (201)	46.3 ± 14.0 (30)	46.9 ± 14.4% (231)
≥ Moderate mitral regurgitation	80.1% (161/201)	66.7% (20/30)	78.4% (181/231)
Annuloplasty ring type			
Partial ring	16.5% (34/206)	26.7% (8/30)	17.8% (42/236)
Circumferential ring	83.5% (172/206)	73.3% (22/30)	82.2% (194/236)
*Continuous measures - Mean ± SD (Total no.); Categorical measures - % (no./Total no.)			

Safety and Effectiveness Results

Safety Endpoints

The Kaplan-Meier estimates of site-reported adverse events through 1 year are presented in Table 102. The Kaplan-Meier curves for all-cause mortality and cardiovascular mortality are shown in Figure 93. The all-cause mortality rate was 10.9% at 30 days and 28.6% at 1 year, including a cardiovascular mortality rate of 5.6% at 30 days and 9.5% at 1 year. Other relatively more frequent adverse events included new requirement for dialysis (6.8% at both 30 days and 1 year), left ventricular outflow tract (LVOT) obstruction (6.6% at 30 days and 7.2% at 1 year), readmission due to heart failure (6.2% at 30 days and 22.4% at 1 year), and non-cardiac readmission (8.9% at 30 days and 32.6% at 1 year).

Table 102: Site-Reported Adverse Events (All Population)

Adverse Event	Kaplan-Meier Rate*	
	30 Days (N = 236)	1 Year (N = 236)
All-cause death	10.9% (24, 24)	28.6% (53, 53)
Cardiovascular death	5.6% (12, 12)	9.5% (18, 18)
All stroke	1.3% (3, 3)	1.3% (3, 3)
Ischemic stroke	0.9% (2, 2)	0.9% (2, 2)
Hemorrhagic stroke	0.4% (1, 1)	0.4% (1, 1)
Transient Ischemic Attack	0.4% (1, 1)	1.4% (2, 2)
Major vascular complication	2.2% (5, 5)	2.2% (5, 5)
Life threatening/Major bleeding	1.4% (3, 3)	4.8% (8, 7)
Myocardial infarction	0.8% (2, 2)	1.5% (4, 3)
New onset atrial fibrillation	2.9% (6, 6)	6.6% (11, 11)
Conduction/native pacer disturbance requiring pacer	0.9% (2, 2)	3.0% (5, 5)
New requirement for dialysis	6.8% (15, 15)	6.8% (15, 15)
Mitral valve reintervention	4.7% (12, 10)	11.4% (22, 19)
Device thrombosis	1.0% (2, 2)	1.0% (2, 2)
Device embolization	0.8% (2, 2)	0.8% (2, 2)
Device migration	0.5% (1, 1)	0.5% (1, 1)
LVOT Obstruction	6.6% (15, 15)	7.2% (16, 16)
Other device related event	4.9% (12, 11)	5.7% (13, 12)
Endocarditis	0.0% (0, 0)	0.8% (1, 1)
Readmission - heart failure	6.2% (13, 12)	22.4% (48, 34)
Readmission - cardiac	3.1% (6, 6)	14.5% (33, 22)
Readmission - non-cardiac	8.9% (18, 18)	32.6% (62, 49)
Unplanned other cardiac surgery or intervention	10.1% (23, 23)	11.7% (25, 25)

*Kaplan-Meier rate - % (no. of events, no. of patients with the event)





Figure 93: Mortality through 1 Year (AI Population)

Note: The 95% confidence intervals were calculated without multiplicity adjustment. The adjusted confidence intervals could be wider than presented here. As such, confidence intervals are provided to illustrate the variability only and should not be used to draw any statistical conclusion.

Effectiveness Endpoints

Valve Performance

The echocardiographic valve performance results are shown in Figure 94 to Figure 96. The mean mitral gradient was 7.9 mmHg at baseline, which was maintained at 30 days (7.9 mmHg) and through 1 year (8.1 mmHg). Moderate or greater total mitral regurgitation was observed in 78.4% of the patients at baseline, which decreased to 5.3% at 30 days and 5.5% at 1 year. The proportion of patients with moderate or greater paravalvular regurgitation was 3.3% at 30 days and 0.0% at 1 year.



Note: Line plot with mean and standard error. The total number of patients at each time point only counted the patients with valid values

Figure 94: Mean Mitral Gradient (VI Population)




Figure 95: Total Mitral Regurgitation (VI Population)




Figure 96: Paravalvular Regurgitation (VI Population)

NYHA Functional Class

The NYHA functional class distributions by visit are presented in Figure 97. At baseline, 80.4% of patients were in NYHA III/IV. At 1 year, the majority (76.3%) of patients were in NYHA I/II.




Figure 97: NYHA Class by Visit (VI Population)

Length of Stay

The mean index hospitalization stay was 7.4 days, which included an average of 2.8 days in the intensive care unit (ICU), as summarized in Table 103.

Table 103: Index Hospitalization (AI Population)

Index hospitalization	Length in days*
Index hospitalization duration	7.4 ± 0.55 (236)
Intensive care stay	2.8 ± 0.31 (233)
*Mean ± SE (Total no.).	

Quality of Life

The results for the KCCQ overall summary score are presented in Figure 98. The mean score increased from 40.3 at baseline to 58.1 and 65.6 at 30 days and 1 year, respectively.




Figure 98: KCCQ Overall Summary Score (VI Population)

Note: Line plot with mean and standard error. The total number of patients at each time point only counted the patients with valid values.

Other Study Observations

Procedural Information

The procedural information is summarized in Table 104. General anesthesia was used in the majority (97.9%) of patients. Conversion to open heart surgery occurred in two patients due to access related problem/injury and device embolization, respectively.

Table 104: Procedural Data Summary (AI Population)

Procedural Data	Summary Statistics* (N = 236)
Operator reason for procedure	
Inoperable/high risk	87.7% (207/236)
Intermediate risk	12.3% (29/236)
Implant approach	
Transseptal	81.8% (193/236)
Transapical	8.5% (20/236)
Femoral artery	8.9% (21/236)
Direct left atrium	0.4% (1/236)
Other	0.4% (1/236)
Valve size	
20 mm	0.4% (1/236)
23 mm	19.9% (47/236)
26 mm	44.5% (105/236)
29 mm	35.2% (83/236)

Procedural Data	Summary Statistics* (N = 236)
Type of anesthesia	
General anesthesia	97.9% (231/236)
Moderate sedation	2.1% (5/236)
Total procedure time (minute)	151.2 ± 82.9 (235)
Device implanted successfully	92.4% (218/236)
Procedure aborted	0.8% (2/236)
Access related	50.0% (1/2)
System issue	50.0% (1/2)
Conversion to open heart surgery	0.8% (2/236)
Access related problem/injury	50.0% (1/2)
Device embolization	50.0% (1/2)
*Continuous measures - mean ± SD (n); categorical measures - % (no./Total no.)	

Subgroup analysis

The Kaplan-Meier curves for all-cause mortality and cardiovascular mortality are shown in Figure 99 and Figure 100 for patients with a partial and circumferential annuloplasty ring, respectively. In patients with a partial annuloplasty ring, the 1-year all-cause and cardiovascular mortality rates were 25.1% and 10.9%, respectively, as compared to the corresponding rates of 29.4% and 9.2% for those with a circumferential annuloplasty ring.





Figure 99: Mortality Through 1 Year (Partial Ring Patients)



No. at Risk

All	194	147	76
Cardiovascular	194	147	76

Figure 100: Mortality Through 1 Year (Circumferential Ring Patients)

SAPIEN 3 Ultra System

Patient Accountability

At the time of the database extract, all 40 patients enrolled were implanted, discharged and completed 30-day follow-up.

Patient Demographics and Baseline Characteristics

The demographics and baseline characteristics are shown in Table 105.

Table 105: Patient Demographics and Baseline Characteristics

Demographics and Baseline Characteristic	Summary Statistics
Age - years	83.4 ± 5.13 (40)
Male sex	24/40 (60.0%)
Society of Thoracic Surgeons (STS) score	3.4 ± 1.27 (40)
New York Heart Association (NYHA) class	
I/II	20/40 (50.0%)
III/IV	20/40 (50.0%)
Previous myocardial infarction	1/40 (2.5%)
Previous intervention	9/40 (22.5%)
Coronary artery bypass grafting (CABG)	3/40 (7.5%)
Percutaneous bypass intervention (PCI)	8/40 (20.0%)
Prior aortic valvuloplasty	0/40 (0.0%)
Cerebrovascular accident (CVA)	1/40 (2.5%)
Peripheral vascular disease	4/40 (10.0%)

Demographics and Baseline Characteristic	Summary Statistics
Atrial fibrillation	19/40 (47.5%)
Prior pacemaker	5/40 (12.5%)
Porcelain aorta	0/40 (0.0%)
Echocardiographic findings	
Valve area - cm ²	0.7 ± 0.16 (40)
Mean aortic-valve gradient -mmHg	51.0 ± 13.17 (40)
Mean left ventricular ejection fraction (LVEF) %	60.6 ± 7.06 (40)
Moderate or severe aortic regurgitation	1/39 (2.6%)
*Continuous measures—Mean ± SD (Total no.); Categorical measures—n./Total no. (%)	

Safety and Effectiveness Results

Primary Endpoint

The primary endpoint was procedural success, defined as freedom from mortality, conversion to surgery, and moderate or severe PVR at exit from the procedure room which was achieved in all subjects as outlined in Table 106.

Table 106: Primary Endpoint Analysis

Primary Endpoint	Results
Overall procedural success	40/40 (100.0%)
Freedom from mortality at exit from procedure room	40/40 (100.0%)
Freedom from conversion to surgery at exit from procedure room	40/40 (100.0%)
Freedom from moderate or severe paravalvular regurgitation at exit from procedure room	40/40 (100.0%)

Secondary Endpoints

There were no major vascular complications, valve migrations, or embolizations through discharge.

Adverse Events

There were no deaths or strokes through 30-days. The selected adverse events for the treated population are presented in Table 107.

Table 107: Selected Adverse Events

Adverse Event	Discharge*	30 Days†
Major vascular complications	0.0% (0, 0)	0.0% (0, 0)
Acute kidney injury (Stage III)	0.0% (0, 0)	0.0% (0, 0)
Life threatening bleeding	0.0% (0, 0)	0.0% (0, 0)
Major bleeding	5.0% (2,2)	5.0% (2,2)
Hematoma	5.0% (2,2)	5.0% (2,2)
Bleeding at access site	15.0% (6,6)	15.0% (6,6)
Dissection	2.5% (1,1)	2.5% (1,1)
Pseudoaneurysm	2.5% (1,1)	2.5% (1,1)
Aortic-valve reintervention	0.0% (0, 0)	0.0% (0, 0)
Endocarditis	0.0% (0, 0)	0.0% (0, 0)
Device thrombosis	0.0% (0, 0)	0.0% (0, 0)
*Observed rate,% (no. of events, no. of subjects with the event)		
†Kaplan-Meier estimate,% (no. of events, no. of subjects with the event)		

The new conduction abnormalities requiring permanent pacemaker implantation through 30-days for the first 20 subjects and the last 20 subjects are presented in Table 108.

Table 108: New Conduction Abnormalities Requiring Permanent Pacemaker Implantation

Adverse Event	First 20 Subjects		Last 20 Subjects	
	Discharge*	30 Days†	Discharge*	30 Days†
Conduction disturbance requiring permanent pacemaker‡	29.4% (5,5)	29.4% (5,5)	5.6% (1,1)	5.6% (1,1)

‡5 Subjects (3 from First 20 subject cohort and 2 from Last 20 subject cohort) with baseline pacemaker were excluded from the analysis.
 *Observed rate,% (no. of events, no. of subjects with the event)
 †Kaplan-Meier estimate,% (no. of events, no. of subjects with the event)


Other Results

Procedural Information

Overall, the mean procedure time was 56.5 ± 26.8 minutes. Conscious sedation was utilized in 95% of the patients with one patient converted to general anesthesia. The valve was placed in the intended position in all cases, there were no aborted implantation procedures or conversion to open heart surgery. Successful access, delivery and retrieval of the device and delivery system occurred in all cases. The average length of stay was 4.1 ± 2.4 days.

Valve Performance

The measurements of effective orifice area, mean gradient, total aortic regurgitation, aortic paravalvular regurgitation (PVL) are presented in Figure 101 to Figure 104. Mean EOA increased and gradients decreased. PVL was trace or none in 85% of the patients.



*Site reported. **Core lab reported.

Figure 101: Effective Orifice Area




Figure 102: Mean Gradient by Valve Size




Figure 103: Total Aortic Regurgitation




Figure 104: Aortic Paravalvular Regurgitation

NYHA

The NYHA Functional Class summary is shown in Figure 105. At 30-day follow-up, 80.0% of subjects experienced improvement in NYHA Class and all subjects were in Class I/II.




Figure 105: NYHA Class by Visit

13.0 References

- [1] Bapat V, Attia R, Thomas M. Effect of Valve Design on the Stent Internal Diameter of a Bioprosthetic Valve: A Concept of True Internal Diameter and Its Implications for the Valve-in-Valve Procedure. *JACC: Cardiovascular Interventions*. Vol. 7, No. 2 2014: 115-127.
- [2] Kappetein AP, Head SJ, Génèreux P, et al. Updated standardized endpoint definitions for transcatheter aortic valve implantation: the Valve Academic Research Consortium-2 consensus document (VARC-2). *Eur J Cardiothorac Surg* 2012;42: S45-60.
- [3] Imbens G. W. (2004) Nonparametric Estimation of Average Treatment Effects under Exogeneity: A Review. *The Review of Economics and Statistics*, February 2004, 86(1): 4–29.
- [4] Abbas AE, et al. (2020) Hemodynamic principles of prosthetic aortic valve evaluation in the transcatheter aortic valve replacement era. *Echocardiography*, 2020;37:738-757.

These products are manufactured and sold under one or more of the following US patent(s): US Patent No. 7,530,253; 7,780,723; 7,895,876; 8,382,826; 8,591,575; 8,690,936; 8,790,387; 9,061,119; 9,301,840; 9,301,841; 9,339,384; 9,393,110; and corresponding foreign patents.



Edwards

2026-04
10033860009 A
© Copyright 2026, Edwards Lifesciences LLC
All rights reserved.

 **Manufacturer**
Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Telephone +1.949.250.2500
+1.800.424.3278
FAX +1.949.250.2525

Web IFU