Edwards

DIRECTORY

English	1	Português.....	25	Eesti	48
Français	4	Česky.....	27	Lietuvių	51
Deutsch	6	Magyar.....	30	Latviešu	53
Español.....	9	Polski	32	Türkçe.....	56
Italiano.....	12	Slovensky.....	35	Русский	58
Nederlands.....	14	Norsk.....	38	Srpski	61
Dansk	17	Suomi.....	40	中文.....	64
Svenska	19	Български.....	43	繁體中文 (台灣)	66
Ελληνικά.....	22	Română	46	한국어	68

English

TruWave Pressure Monitoring Kit with TruWave Disposable Pressure Transducer

Instructions for Use For Single Use Only

Carefully read all instructions for use, warnings, and precautions before using the device.

These are general instructions for setting up a Pressure Monitoring system with the Edwards Lifesciences TruWave disposable pressure transducer. Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

Caution: The use of lipids with the TruWave disposable pressure transducer may compromise product integrity.

Description

The Edwards Lifesciences Pressure Monitoring Kit with TruWave disposable pressure transducer, is a sterile, single-use kit that monitors pressures. The disposable sterile cable (available in 12-inch/30 cm and 48-inch/120 cm lengths) interfaces exclusively with an Edwards Lifesciences cable that is specifically wired for the monitor being used.

The TruWave disposable pressure transducer has a straight, flow-through design across the pressure sensor, and is available with or without the integral flush device.

Indications

The Pressure Monitoring Kit with TruWave disposable pressure transducer is for use on

patients requiring intravascular, intracranial, or intrauterine pressure monitoring.

Contraindications

There are no absolute contraindications for using the TruWave disposable pressure transducer in patients requiring invasive pressure monitoring.

Flush devices should not be used when monitoring intracranial pressure.

Warnings

- **Do not use the flush device during intracranial pressure monitoring.**
- **To avoid over infusing severely fluid-restricted patients such as neonates and children, use an infusion pump in series with the flush device to accurately regulate the minimum amount of flush solution needed to maintain catheter patency while allowing continuous pressure monitoring. Use an IV set that is indicated for use with the IV pump. High pressures, which may be generated by an infusion pump at certain flow rates, may override the flush device restriction, resulting in fast flushing at the rate set by the pump.**
- **Do not allow air bubbles to enter the setup, especially when monitoring arterial pressures.**
- **Avoid contact with any topical cream or ointment that attacks polymeric materials.**
- **Do not expose electrical connections to fluid contact.**
- **Do not autoclave the reusable cable.**
- **This device is designed, intended, and distributed for single use only. Do not resterilize or reuse this device. There are**

no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

- **Do not use this device if it has been exposed to an environment outside the range of 0 °C to 40 °C and 5 % to 90 % RH, or if it has gone beyond its Use By date.**
- **For patients undergoing an MR examination refer to the MRI Safety Information section for specific conditions.**

Equipment

The Edwards TruWave disposable pressure transducers are provided in pressure monitoring kits. Pressure monitoring kits may comprise of the following components or accessories:

1. Pressure tubing
2. Stopcocks
3. Flush device (3 mL/hr or 30 mL/hr)
4. IV set
5. VAMP (closed blood sampling system)
6. Disposable holder
7. Arm strap
8. Compatible monitor cable

Note: Refer to VAMP closed blood sampling system Instructions for Use for closed blood sampling instructions.

Procedure

1. Ensure that the cable is compatible with the monitor being used. Connect the reusable cable to the monitor. Turn the monitor on to allow the electronics to warm up.

Edwards, Edwards Lifesciences, the stylized E logo, TruWave, and VAMP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Warning: Moisture within the connector may result in the device malfunctioning or in inaccurate pressure readings. If this occurs, replace wet part or parts.

Note: Refer to compatible monitor instructions for proper cable to monitor connection.

- Using aseptic technique, remove the transducer and kit from the sterile packaging.
- Add additional components as needed to complete the monitoring system per hospital policy.
- Ensure that all connections are secure.

Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.

- Connect the transducer cable to the reusable cable.
 - Remove all air from the IV flush solution bag, which is heparinized per hospital policy.
- Caution:** If all air is not removed from the bag, air may be forced into the patient's vascular system when the solution is exhausted.
- Close the roller clamp on the IV set and connect it to the IV flush bag. Hang the bag approximately 2 feet (60 cm) above the patient. This height will provide approximately 45 mmHg of pressure to the prime setup.
 - Fill the drip chamber halfway, or per hospital policy, with flush solution by squeezing the drip chamber. Open the roller clamp.
 - To fill system:
 - For transducers without integral flush device (Snap-Tab), fill system per hospital policy.
 - For transducers with integral flush device, flow is provided by pulling on the Snap-Tab, and discontinued by releasing the Snap-Tab.

- For kits with IV sets attached, open the transducer vent port by turning the stopcock handle towards pressure tubing. Allow gravity to deliver flush solution first through the transducer and out through the vent port, then through the remaining pressure tubing by turning the appropriate stopcocks. Remove all air bubbles.

Caution: Significant distortion of the pressure waveform or air emboli can result from air bubbles in the setup.

- Replace all vented caps on sideports of the stopcocks with non-vented caps.
- Mount the transducer either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder.
- Pressurize the IV flush solution bag. Flow rate will vary with pressure across the flush device. The flow rates with the IV bag pressurized to 300 mmHg are as follows:
 - 3 ± 1 mL/hr (DPT with blue Snap-Tab)
- Connect pressure tubing to the catheter per manufacturer's instructions.

- Flush system per hospital policy.

Caution: After each fast-flush operation, observe the drip chamber to verify that the continuous flush rate is as desired (see **Complications**).

Caution: All steps should be done prior to connecting to the patient's catheter or access site.

Zeroing and Calibration

- Adjust the level of the transducer vent port (the fluid-air interface) to correspond to the chamber where pressure is being measured. For example, in cardiac monitoring zero at level of the right atrium. This is at the phlebostatic axis, determined by the intersection of the midaxillary line and the fourth intercostal space.
- Remove the non-vented cap and open the vent port to the atmosphere.
- Adjust the monitor to read zero mmHg.
- Check monitor calibration using procedure recommended by the monitor manufacturer.
- Close the vent port to the atmosphere and replace the non-vented cap.
- System is ready to begin monitoring pressure.

Testing Dynamic Response

The assembly may be tested for dynamic response by observing the pressure waveform on a monitor. Bedside determination of the dynamic response of the catheter, monitor, kit, and transducer system is done after the system is flushed, attached to the patient, zeroed, and calibrated. A square-wave test may be performed by pulling the Snap-Tab and releasing quickly.

Note: Poor dynamic response can be caused by air bubbles, clotting, excessive lengths of tubing, excessively compliant pressure tubing, small bore tubing, loose connections, or leaks.

Routine Maintenance

Follow hospital policies and procedures for frequency of zeroing the transducer and monitor and for replacing and maintaining pressure monitoring lines. The TruWave disposable pressure transducer is precalibrated and has a negligible drift with time (see **Specifications**).

- Adjust zero pressure reference each time level of the patient is changed.

Caution: When rechecking zero or verifying accuracy, ensure that the non-vented cap is removed before opening the transducer vent port to the atmosphere.

- Periodically check fluid path for air bubbles. Ensure that connecting lines and stopcocks remain tightly fitted.
- Periodically observe the drip chamber to verify that the continuous flush rate is as desired.
- The Centers for Disease Control recommends replacing IV lines and disposable or reusable transducers at 96-hour intervals.

MRI Safety Information

The following device was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment:

TruWave disposable pressure transducer

Non-clinical testing demonstrated that these devices are MR Conditional according to the following conditions:

-Static magnetic field of 3-Tesla or less

-Maximum spatial gradient field of 4000-Gauss/cm (40.0 T/m) or less

-This device and the associated cable are not intended for use inside of the bore of the MR system and should not be in contact with the patient

-This device and the associated cable may be in the MR system room but not in operation or connected to a pressure monitoring system during an MR examination

-The pressure tubing connected to the TruWave DPT is MR safe and can be placed inside the bore during an MR examination

Complications

Sepsis/Infection

Positive cultures can result from contamination of the pressure setup. Increased risks of septicemia and bacteremia have been associated with blood sampling, infusing fluids, and catheter related thrombosis.

Air Emboli

Air can enter the patient through stopcocks that are inadvertently left open, from accidental disconnection of the pressure setup, or from flushing residual air bubbles into the patient.

Clotted Catheter and Bleed-Back

If the flush system is not adequately pressurized relative to the patient's blood pressure, blood bleed-back and catheter clotting may occur.

Overinfusion

Excessive flow rates may result from pressures greater than 300 mmHg. This may result in a potentially harmful increase in blood pressure and fluid overdose.

Abnormal Pressure Readings

Pressure readings can change quickly and dramatically because of loss of proper calibration, loose connection, or air in the system.

Warning: Abnormal pressure readings should correlate with the patient's clinical manifestations.

How Supplied

TruWave disposable pressure transducers are supplied sterile in preconnected monitoring kits (either standard design or special order). This product is for single use only. Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

Storage

Store in a cool, dry place.

Temperature Limitation: 0 °C - 40 °C

Humidity Limitation: 5 % - 90 % RH

Shelf Life

The recommended shelf life is as marked on each package.

Resterilization will not extend the shelf life of this product.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Specifications*

TruWave Disposable Pressure Transducer

Operating Pressure Range	-50 to + 300 mmHg
Operating Temperature Range	15 to 40 °C
Sensitivity	5.0µ V/V/mmHg ±1%
Accuracy: Combined effects of Sensitivity, Repeatability, Nonlinearity and Hysteresis	±1 mmHg plus 1% of reading from -50 to +50 mmHg and ±2.5% of reading from +50 to +300 mmHg
Excitation Frequency	dc to 5,000 Hz
Excitation Impedance	350 ohms ± 10% with typical Edwards Monitor Cable attached
Excitation Impedance (DPT)	1800 ohms - 3300 ohms
Phase Shift	<5°
Signal Impedance	300 ohms ± 5 %
Zero Offset	≤±25 mmHg
Zero Thermal Drift	≤±0.3 mmHg/ °C (for products sterilized by EO) ≤ ± 0.5 mmHg/ °C (for products irradiated by e-beam)
Zero Output Drift	±1 mmHg per 8 hours after 20 seconds warm-up
Sensitivity Thermal Drift	≤ 0.4 % / °C
Defibrillator Challenge	withstands 5 repeated discharges of 360 Joules within 5 minutes delivered into a 50 ohm load
Leakage Current	<2µ amps at 120V RMS 60 Hz or at 264V RMS 50 Hz
Overpressure Tolerance	-500 to +5000 mmHg
Shock Resistance	withstands 3 drops from 1 meter
Light Sensitivity	< 1 mmHg at 6 volts excitation when exposed to a 3400 °K tungsten light source at 3000 foot candles (32,293 lm/m ²)
Volumetric Displacement	≤ 0.03 mm ³ /100 mmHg for transducer without flush device
Flow rate across flush device with IV bag pressurized to 300 mmHg	
Blue Snap-Tab	3 ± 1 mL/hr
Yellow Snap-Tab	30 ± 10 mL/hr

* at 6.00VDC and 25 °C unless otherwise stated

All specifications meet or exceed the AAMI BP22 Standard for blood pressure transducers.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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©Copyright 2017, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences AG
Edwards Lifesciences LLC