DIRECTORY

English	1	Português	14	Eesti	27
Français	2	Česky	15	Lietuvių	28
Deutsch	4	Magyar	17	Latviešu	30
Español	5	Polski	18	Türkçe	31
Italiano	7	Slovensky	20	Русский	33
Nederlands	8	Norsk	21	Srpski	34
Dansk	10	Suomi	22	中文	36
Svenska	11	Български	24	繁體中文 (台灣)	37
Ελληνικά	12	Română	25	한국어	38

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Thru-Lumen Embolectomy Catheter

For Single Use Only

Indications

The Fogarty Thru-Lumen Embolectomy Catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels in the arterial system. Additionally, it may also be used for temporary occlusion of blood vessels, infusion of fluids, and blood sampling.

The device is a double lumen catheter with a latex balloon at the distal end. One lumen is used for balloon inflation and is accessed through the gate valve at the proximal end. The thru-lumen is used for infusion of fluids, sampling of blood or guidewire passage through the distal tip. A removable stainless steel stylet(s) is provided.

To remove fibrous or adherent material, alternative devices such as the Fogarty corkscrew catheter and Fogarty graft thrombectomy catheter are recommended.

Contraindications

- The thru-lumen embolectomy catheter should not be used outside the arterial system or for endarterectomy procedures.
- The thru-lumen embolectomy catheter is not recommended for the removal of fibrous, adherent, or calcified material (e.g., chronic clot, atherosclerotic plaque). The catheter is not

designed to withstand the additional pull force needed to remove these materials.

- The thru-lumen embolectomy catheter is not intended for use as a vessel dilator.
- The thru-lumen embolectomy catheter should not be used in the venous system.

Warnings

- Balloon rupture and catheter separation as a result of excessive pull force applied to remove adherent material are the most frequent causes of reported failures. The possibility of balloon rupture must be taken into account when considering the risks involved in any embolectomy procedure.**
- To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume and pull force for each size catheter. See Specification Table.**
- Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.**
- To avoid air embolus in case of balloon rupture, air should not be used for balloon inflation. Carbon dioxide is the only recommended gas. See Instructions.**

Instructions

Precaution: To minimize lateral wall pressures and shear forces on the intimal lining of the vessel wall, the smallest diameter balloon catheter that will remove the embolus should be used.

Prepare

Remove all stylets from the catheter.

Purge

Inflate the balloon with sterile fluid or gas to the maximum recommended volume. Pull a vacuum on the syringe. Repeat until all air is removed. **Note:** For all inflations, use the smallest syringe capable of holding the stated maximum fluid capacity. Attach another syringe filled with sterile, heparinized saline or other sterile solution to the thru-lumen port and purge air from the thru-lumen.

Inspect

The catheter should be inspected with the balloon inflated during purging. A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

Caution: The amount of fluid in the syringe should be checked before each inflation. Do not exceed the recommended maximum inflation volume. See Specification Table.

Place

With the balloon deflated, insert the catheter into the vessel and beyond the obstructing material either directly by cutdown or by percutaneous

Specification Table

Catheter French Size (mm)	3 1.00	4 1.33	5.5 1.83	6 2.0	7 2.3
Maximum Liquid Capacity (ml)	0.15	0.5	0.9	1.4	1.6
Maximum Gas Capacity (ml)	0.4	1.2	1.7	2.3	2.8
Maximum Recommended Pull Force on Inflated Balloon (lbs.)	0.7	1.5	2.0	2.5	3.5
Maximum Recommended Pull Force on Inflated Balloon (kg)	0.32	0.68	0.91	1.13	1.59
Diameter of Inflated Balloon (mm)	5	9	11	13	14
Maximum French Size of Deflated Balloon	4.5	5	6	7	8
Maximum Size of Deflated Balloon (mm)	1.50	1.67	2.0	2.3	2.7
Compatible Guidewire Diameter (in.)	0.018	0.025	0.035	0.035	0.038
Compatible Guidewire Diameter (mm)	0.46	0.64	0.89	0.89	0.97
* Compatible Introducer Size (F)	4.5	5	6	7	8
* Compatible Introducer Size (mm)	1.50	1.67	2.0	2.3	2.7

* Ease of introducing the catheter through hemostasis valve may vary among manufacturers.

technique through a suitable introducer. Heparinized saline or other solutions may be flushed through the catheter thru-lumen as needed during the procedure. Do not advance the catheter against resistance.

Inflate

With the catheter suitably positioned, inflate the balloon with sterile fluid or gas. Balloon inflation should be stopped when the balloon can be felt to engage the vessel wall. **Note:** CO₂ is the preferred medium in cases where gas is used for inflation. It should be noted that CO₂ diffuses through the balloon wall and corrections for gas loss must be made. Use of liquid inflation media will result in slower inflation/deflation rates.

Caution: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest inflated balloon diameter that will remove the obstructing material.

Withdraw

Remove the occlusive material by gently withdrawing the catheter. During withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by controlling the inflation volume.

Caution: Do not exceed the Maximum Recommended Pull Force. See Specification Table.

Caution: Inflation of the balloon is associated with a feeling of resistance to the fluid or gas injection. When no resistance is encountered it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

MRI Information

This product has not been tested for MRI compatibility.

Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection, local hematomas, intimal disruption, arterial dissection, perforation and vessel rupture, hemorrhage, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.
Temperature Limitation: 0 ° - 40 °C,
Humidity Limitation: 5% - 90% RH

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. Do not re-sterilize or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC