DIRECTORY

English	1	Português	12	Eesti	24
Français	2	Česky	14	Lietuvių	25
Deutsch	4	Magyar	15	Latviešu	26
Español	5	Polski	16	Türkçe	28
Italiano	6	Slovensky	18	Русский	29
Nederlands	7	Norsk	19	Srpski	30
Dansk	9	Suomi	20	中文	31
Svenska	10	Български	21	繁體中文 (台灣)	32
Ελληνικά	11	Română	23	한국어	33

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Arterial Embolectomy Catheter with Gate Valve

For Single Use Only

Indications

The Fogarty arterial embolectomy catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels in the arterial system.

To remove fibrous or adherent material, alternative devices such as the Fogarty adherent clot and graft thrombectomy catheter are recommended.

Contraindications

- The arterial embolectomy catheter should not be used outside the arterial system.
- The arterial embolectomy catheter is not recommended for the removal of fibrous, adherent, or calcified material (e.g. chronic clot, atherosclerotic plaque). The catheter is not designed to withstand the additional pull force needed to remove these materials.
- The arterial embolectomy catheter is not intended for use as a vessel dilator.
- The arterial embolectomy catheter should not be used in the venous system.

Warnings

- Balloon rupture and catheter separation as a result of excessive pull force applied to remove adherent material are the most frequent causes

of reported failures. The possibility of balloon rupture must be taken into account when considering the risks involved in any embolectomy procedure.

- To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume and pull force (see **Specification Table**).
- Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.
- To avoid air embolus in case of balloon rupture, air should not be used for balloon inflation. Carbon dioxide is the only recommended gas. See Instructions.

Instructions

Prepare

Fill the inflation syringe with 0.2 ml of carbon dioxide gas. Tightly connect the syringe to the catheter gate valve and set the gate valve to its "Open" position.

Note: For all inflations use the smallest syringe capable of holding the stated maximum gas capacity.

Purge

Inflate the balloon with carbon dioxide gas to the maximum recommended volume. Pull a vacuum on the syringe. Repeat until all air is removed.

Inspect

The catheter should be inspected with the balloon inflated during purging. A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

Caution: The amount of gas in the syringe should be checked before each inflation. Do not exceed the recommended maximum inflation volume. (See **Specification Table**).

Place

With the balloon deflated, insert the catheter into the vessel and beyond the obstructing material.

Inflate

With the catheter suitably positioned, inflate the balloon with carbon dioxide gas. Balloon inflation should be stopped when the balloon can be felt to engage the arterial wall. Gas seepage through the intact balloon will require frequent adjustment of the inflation volume.

Caution: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest inflated balloon diameter that will remove the obstructing material.

Withdraw

Remove the occlusive material by gently withdrawing the catheter. During withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by controlling the inflation volume. **Caution:** Do not exceed the Maximum

Recommended Pull Force. (See **Specification Table**).

Caution: Inflation of the balloon is associated with a feeling of resistance to the gas injection. When no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

MRI Information

This product has not been tested for MRI compatibility.

Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection, local hematomas, intimal disruption, arterial dissection, perforation and vessel rupture, hemorrhage, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.
Warehouse Storage: 0° - 40 °C, 5% - 90% RH.

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration. **Note:** Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: .. 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Sterilized Using Ethylene Oxide

Specification Table

Model Number	12GV0402F
Maximum Gas Capacity (ml)	0.2
Maximum Recommended Pull Force on Inflated Balloon (lbs.)	0.5 (0.2 kg)
Diameter of Inflated Balloon (mm)	4
Stylet	No
Maximum French Size of Deflated Balloon	3.1F (1.03 mm)
Catheter French Size	2F (0.67 mm)
Length (cm)	45

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC