DIRECTORY

English	1	Português	11	Eesti	20
Français	2	Česky	12	Lietuvių	21
Deutsch	3	Magyar	13	Latviešu	22
Español	5	Polski	14	Türkçe	23
Italiano	6	Slovensky	15	Русский	24
Nederlands	7	Norsk	16	Srpski	25
Dansk	8	Suomi	17	中文	26
Svenska	9	Български	18	繁體中文 (台灣)	27
Ελληνικά	10	Română	19	한국어	28

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Corkscrew Catheter

For Single Use Only

Description

The device is a flexible catheter consisting of a polyvinylchloride outer catheter body surrounding a movable, radiopaque stainless steel spring body. A second radiopaque, spiral-shaped stainless steel cable covered by a latex membrane is attached to the distal end of the catheter. A handle to modulate the spiral cable diameter is attached to the inner spring body at the proximal end of the catheter.

Indications

The Fogarty Corkscrew catheter is indicated for the removal of emboli and thrombi from either native vessels or synthetic grafts in the arterial system.

Contraindications

- The Fogarty Corkscrew catheter should not be used outside the arterial system.
- The Fogarty Corkscrew catheter is not recommended for endarterectomy procedures.
- The Fogarty Corkscrew catheter is not intended for use as a vessel dilator.
- The Fogarty Corkscrew catheter should not be used in the venous system.

Warnings

- **Exposure to the atmosphere, handling during insertion, and plaque and other deposits within the blood vessel may cause latex membrane degradation.**
- **To minimize the risk of vessel or membrane damage, do not exceed the maximum recommended pull force (see Catheter Specifications).**

Instructions

Technique

1. Inspect the catheter.
2. With the membrane extended, with the control button in the locked position, place the catheter in the artery or graft so that the membrane is past the clot to be removed.
3. Contract the membrane by pushing down the control button and pulling it backwards. This causes the membrane to expand into a helical configuration so that it grips or engages the clot.
4. Withdraw the entire catheter from the vessel while the membrane maintains slight contact with luminal wall, removing the occluding material as the catheter is withdrawn.

Note: A sensitive touch in controlling the contracted membrane during catheter withdrawal is important so that the membrane diameter can be adjusted to varying luminal

diameters. Contraction of the membrane is associated with a feeling of resistance.

Note: Because of the complexity and variation in the procedure, the choice of surgical technique, appropriately modified in accordance with the previously described warnings, precautions, and technique, is left to the discretion of the individual surgeon.

MRI Information

This product has not been tested for MRI compatibility.

Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection, local hematomas, intimal disruption, arterial dissection, perforation and vessel rupture, hemorrhage, arterial thrombolysis, distal emboli or blood clots or arteriosclerotic plaque, air emboli, aneurysms, arterial spasms, arteriovenous fistula formation, membrane separation and distal embolization.

How Supplied

Contents are sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Edwards, Edwards Lifesciences, the stylized E logo, Corkscrew, Fogarty and Fogarty Corkscrew are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Specification Table

Model	140806	140808	1408010
French Size	4F	5F	6F
(mm)	1.33	1.67	2.0
Usable Length (cm)	80	80	80
Maximum Diameter of Extended Membrane (mm)	2.7	3.1	3.3
Diameter of Contracted Membrane (mm)	6	8	10
Maximum Recommended Pull Force on Contracted Membrane (lbs.)	0.9 (0.4 kg)	1.3 (0.6 kg)	1.8 (0.8 kg)

Storage

Store in a cool, dry place.

Temperature Limitation: 0° - 40 °C.

Humidity Limitation: 5% - 90% RH.

Shelf Life

The recommended shelf life is marked on each package. Storage beyond the expiration date may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: . 0870 606 2040 - Option 4

In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC