DIRECTORY

English	1	Português	14	Eesti	27
Français	3	Česky	16	Lietuvių	28
Deutsch	4	Magyar	17	Latviešu	29
Español	6	Polski	18	Türkçe	31
Italiano	7	Slovensky	20	Русский	32
Nederlands	9	Norsk	21	Srpski	34
Dansk	10	Suomi	22	中文	35
Svenska	11	Български	24	繁體中文 (台灣)	36
Ελληνικά	13	Română	25	한국어	38

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Biliary Balloon Probe

For Single Use Only

Indications

Indications for use include removal of stones and ductal debris and exploration.

Warnings

The use of highly viscous or particulate contrast media is not recommended for balloon inflation as the lumen may become occluded.

The maximum recommended volume for balloon inflation must not be exceeded as over inflation increases the possibility of balloon rupture.

Precautions

As with all catheterization procedures, complications may occur. These include intrahepatic ductal damage, transient bleeding, postoperative cholangitis, and stricture of a primary hepatic duct. Instances of hepatic duct disruption during retrieval of intrahepatic stones have also been reported.

The balloon must be allowed to deflate as it passes through the ampullary region to avoid damage to the sphincter.

Instances of balloon rupture are minimized if care is taken to inflate the balloon slowly and withdraw the catheter gently.

Instructions

Technique

The probe is introduced with the balloon deflated into the common duct through the choledochotomy. When the desired position has been reached, the balloon is inflated by the introduction of a sterile fluid. It may then be gently withdrawn for the removal of debris or manipulated for further exploration. During the period over which this technique has been developed, the Fogarty Biliary Balloon Probe has been found to be a valuable adjunct to presently available techniques. The Fogarty Biliary Balloon Probe does not, however, have utility in removing stones when the probe cannot pass beyond the calculi or when a calculus is in a sacculus.

The Biliary Balloon Probe aids in easy and rapid removal of amorphous intraductal debris and gravel. Numerous passes may be required to remove all material, but care should be taken not to over-distend the common duct during the period of manipulation. Closed-end irrigation of the common duct may be done by the introduction of the probe into the duodenum, followed by inflation and impingement of the inflated balloon against the sphincter mechanism. The irrigation catheter may then be employed. More effective washing is sometimes accomplished by closure of an open sphincter mechanism.

Atraumatic Calibration of the Sphincter Mechanism

The probe is introduced into the distal common duct and passed through the sphincter into the duodenum. Once in the duodenum, the balloon is inflated and then carefully palpated to make certain that it has passed through the sphincter mechanism. While the surgeon keeps the fingers of one hand placed on the balloon in the duodenum, the balloon probe is withdrawn until it impinges at the area of the sphincter. With light traction on the probe, the surgeon deflates the balloon until it passes through the sphincter with ease. As the balloon is passed through the sphincter, the quantity of fluid contained within the syringe is noted. When the balloon re-enters the common duct, the balloon is inflated until it engages the duct wall, and then gently withdrawn. The same amount of fluid which was noted as the probe passed through the sphincter is reintroduced into the balloon. Calibration of the balloon diameter accurately reflects the size of the sphincter. The common duct can be calibrated in the same manner. Mobilization of the duodenum greatly facilitates this maneuver. This technique aids in identification of the retroduodenal and pancreatic position of the common duct and accurately locates the sphincter of Oddi.

For Removal of Intrahepatic Stones

The probe can be bent with the stylet in place to facilitate introduction into the right or left hepatic duct. The stylet may be removed after the probe enters into the hepatic duct. The probe is then

advanced as far as possible with the balloon deflated. As the balloon is gently inflated, the probe is simultaneously removed. After noting when the probe engages the wall of the ductal system, the inflated balloon is carefully withdrawn through the incision in the common duct. Several passes may be necessary. The final exploratory passage in these radicles is carried out by insertion of probes into both the right and left hepatic ducts. The two probes are then simultaneously removed. This maneuver prevents an isolated stone from falling from one hepatic duct into the other during removal. Undue manipulation, probing, or over-distention of the balloon may cause bleeding and should be avoided.

Identification of the Common Duct

Where there is difficulty in identifying the entire extent of the common duct because of previous surgery or inflammation, the use of the probe allows identification of the duct and facilitates dissection. If the gall bladder remains and the cystic duct has been identified, a small incision in the cystic duct with introduction of the probe into the distal common duct is easily accomplished. Similarly, if the common duct has been identified by needle absorption, a small puncture wound is made in the anterior surface of the common duct and the balloon probe is introduced. After entrance into the common duct, palpation of the inflated balloon in various locations facilitates dissection and identification of the common duct and associated anatomy.

Use of Stylet

The stylet supplied with the Biliary Balloon Probe is of a designed stiffness intended to avoid perforation and yet maintain enough shape to facilitate introduction of the probe into areas of sharp angulation.

When introducing the probe into the common duct, entrance into the duodenum may be accomplished with the stylet in place and angulated at approximately 30°. If difficulty is encountered in passing through the sphincter mechanism, the stylet should be unscrewed until it lies free of the probe body. A gentle probing motion with the stylet in this position will allow for passage into the duodenum if the sphincter mechanism is large enough to accommodate the probe body. When extreme pliability of the distal tip of the instrument is desired, the stylet should be withdrawn 1 cm. A probe possessing a longer length of intermediate pliability is achieved by withdrawing the stylet a greater length.

Balloon Inflation

The balloon should be inflated with a sterile, blood compatible fluid that may include radiopaque solutions, if they are highly dilute and non particulate.

Before use, the probe with the balloon inflated to its maximum recommended volume should be inspected by the surgeon.

The amount of fluid in the syringe should be checked before each inflation. If the amount exceeds the balloon’s stated capacity, remove the syringe from the probe and refill to the proper volume. The amount of fluid remaining in the catheter lumen must be taken into account.

MRI Information

This product has not been tested for MRI compatibility.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.
Temperature Limitation: 0° - 40 °C,

Humidity Limitation: 5% - 90% RH.

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration. **Note:** Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: . 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Sterilized Using Ethylene Oxide

Specification Table

French Size	5F	6F (2.0 mm)
	(1.67 mm)	
Maximum Liquid Capacity (ml)	0.4	1.25
Maximum Air Capacity (ml)	0.8	2.5
Maximum Inflated Balloon Diameter (mm)	8	13

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC