DIRECTORY

English	1	Português	12	Eesti	23
Français	2	Česky	13	Lietuvių	24
Deutsch	4	Magyar	14	Latviešu	25
Español	5	Polski	16	Türkçe	26
Italiano	6	Slovensky	17	Русский	28
Nederlands	7	Norsk	18	Srpski	29
Dansk	9	Suomi	19	中文	30
Svenska	10	Български	20	繁體中文 (台灣)	32
Ελληνικά	11	Română	22	한국어	33

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Venous Thrombectomy Catheter

For Single Use Only

Description

Fogarty venous thrombectomy catheters have a long, soft tip which permits advancement past the venous valves without undue trauma. The flexibility required for the venous procedure is supplied by the spring-wound body. The spring is covered with a braiding for additional strength. The gate valve is designed for single-handed operation and to minimize the possibility of leakage. The gate valve has an arrow to indicate the "open" and "closed" positions. There are three available catheter sizes (the 6 French (2.0 mm), 8 French (2.7 mm), and 8/10 French (2.7/3.3 mm)). Each may be quickly identified by the color-coded body. The size and filling capacity are printed on each catheter.

Indications

Fogarty venous thrombectomy catheters are specifically designed for the performance of venous thrombectomy. The Fogarty venous thrombectomy catheter takes the path of least resistance through the obstructing material. Once properly positioned, the balloon is inflated with a sterile medium. By the mechanism of fluid displacement, the balloon maintains uniform contact with the vessel wall as the catheter and occluding material are withdrawn. Additionally, it may be used for the temporary occlusion of blood vessels.

Fogarty venous thrombectomy catheters are designed for use in the venous system only.

Contraindications

Venous thrombectomy catheters are not designed to withstand the additional strain encountered during endarterectomy procedures.

Warnings

- **Air should not be used for inflation in instances where balloon rupture could produce a dangerous air embolus.**
- **The maximum recommended volume should not be exceeded, as overinflation increases the possibility of balloon rupture. (See specification table).**
- **Undue manipulation or overdistension of the balloon can result in significant damage to the vessel lining, increasing the potential for recurrent thrombosis.**
- **Experience, along with the realization of the limitations of the instrument, have been the most significant factors in reducing the incidence of complications.**

Instructions

Balloon Inflation

The balloon should be inflated with a sterile, blood-compatible fluid which may include radiopaque solutions if they are highly dilute and nonparticulate.

Carbon dioxide should be substituted in cases where gas is used as the inflation fluid. However, it should be recognized that CO₂ will readily diffuse through the balloon wall, and corrections for gas loss must be made.

Use of a highly viscous or highly particulate contrast medium is not recommended for balloon inflation as the inflation passage may become occluded. Prior to use, the catheter, with the balloon inflated to its maximum recommended volume, should be carefully inspected by the surgeon.

The smallest syringe which will hold the stated maximum fluid capacity should be used, as the smaller diameter facilitates control of the withdrawal force because greater displacement of the plunger is required to make a given change in a balloon diameter.

The amount of fluid in the syringe should be checked before each inflation. If it is more than the maximum capacity, the syringe should be removed from the catheter and the barrel refilled to the proper volume, taking into account the amount of fluid remaining in the lumen of the catheter.

Inflation of the balloon is associated with a feeling of resistance. When no resistance is encountered, it should be assumed that the balloon has ruptured. Inflation should be discontinued at once, and the catheter withdrawn.

Specification Table/English

Model	320806F	320808F	32080810F
Maximum Liquid Capacity (ml)	1.5	2.25	4.0
Maximum Gas Capacity (ml)	3.0	4.5	8.0
Diameter of Inflated Balloon (mm)	12.0	13.0	19.0
Maximum French size Over Deflated Balloon	7.7F (2.6 mm)	8.8F (2.9 mm)	10.0F (3.3 mm)
Catheter French Size	6.0F (2.0 mm)	8.0F (2.7 mm)	8.0F (2.7 mm)
Length (cm)	80	80	80

MRI Information

This product has not been tested for MRI compatibility.

Complications

As with all catheterization procedures, complications may occur. Incidents such as breakage of the balloon, rupture of the balloon with fragmentation and embolization, separation of the catheter tip with distal embolization, disruption of the intima, and arteriovenous fistula have been reported with the use of thrombectomy catheters.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.

Temperature Limitation: 0° - 40 °C,

Humidity Limitation: 5% - 90% RH.

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: . . . 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC