DIRECTORY

English	1	Português	10	Eesti	19
Français	2	Česky	11	Lietuvių	20
Deutsch	3	Magyar	12	Latviešu	21
Español	4	Polski	13	Türkçe	22
Italiano	5	Slovensky	14	Русский	23
Nederlands	6	Norsk	15	Srpski	24
Dansk	7	Suomi	16	中文	25
Svenska	8	Български	17	繁體中文 (台灣)	26
Ελληνικά	9	Română	18	한국어	27

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Graft Thrombectomy Catheter

For Single Use Only

Description

The Fogarty Graft Thrombectomy catheter is designed to remove organized thrombus and adherent thrombotic material from synthetic bypass grafts, specifically AV fistulas for dialysis access and grafts for iliofemoral bypass. During use, the control button on the device handle is grasped and moved forward and aft to cause the spiral wires at the distal end to extend or retract. The spiral wires are extended during passage of the device through a narrowed region of a vessel. The spiral wires are then retracted to facilitate occlusive material removal.

Indications

General indications for arterial reconstruction of an obstructed graft include disabling claudication, ischemic pain at rest, or ischemic skin lesions or gangrene.

General indications for reconstruction of a synthetic AV fistula internal bridge graft include the inability to properly dialyse the patient.

Contraindications

The Fogarty Graft Thrombectomy catheter is not designed for use in native vessels.

In some cases, an excessive amount of force may be necessary to remove adherent material from within the graft. When this situation occurs, it is often advisable to abandon the thrombectomy approach and consider an alternative method of vascular reconstruction.

Precautions

Do not use excessive force with this device to remove occlusive material from a graft. As with any blind procedure, extreme care and discretion must be exercised at all times.

Instructions

1. Make an appropriate incision to expose the occluded graft.
2. Make an appropriate incision in the graft for introduction of the catheter.
3. Lock the spiral wires of the graft thrombectomy catheter in the extended position and pass the device through the incision.
4. Pass the catheter into the occluded portion of the graft. Retract the wires by pushing down on the control button and pulling it backwards.

Note: The degree of expansion of the wires can be controlled by adjusting the button position of the handle.

5. With the wires retracted, draw the catheter back through the proximal incision and remove any occlusive material lodged around the spiral loops of the tool.

6. Repeat steps 3 through 5 until thrombectomy is complete.
7. If needed, pass balloon embolectomy catheter through the recanalized graft to remove any residual material.
8. Close the proximal and distal (if appropriate) incisions.
9. Employ arteriography to verify patency and distal outflow and then remove clamps or other occlusive devices.

MRI Information

This product has not been tested for MRI compatibility.

Complications

Potential complications include, but are not limited to, perforation, vessel rupture, thrombus formation and embolism.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.

Temperature Limitation: 0° - 40 °C,

Humidity Limitation: 5% - 90% RH.

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: . 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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Reference

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1. Shook, C.D., et al. Another Alternative in Leg Revascularization. *The American Journal of Surgery*, Vol. 148, 224-228, August 1984.
2. Lerwick, E.R. Oscillating Loop Endarterectomy for Peripheral Vascular Reconstruction. *Surgery*, 97:5, 574-584, May 1985.
3. Inahara, T. Endarterectomy for Occlusive Disease of the Aortoiliac and Common Femoral Arteries. *The American Journal of Surgery*, Vol. 124, 235-243, August 1972.
4. LeVeen, H.H. Technical Features in Endarterectomy. *Surgery*, 57:1, 22-27, January 1965.
5. Capdevila, J.M. Aorto-iliac Endarterectomy: A review of 249 Operations. *11th World Congress of the International Cardiovascular Society*, 34-37.
6. Crispin, H.A. Experience with the Vascular Brush. *Journal of Cardiovascular Surgery*, Vol. 28, 45-49, 1987.
7. Darling, R.C., et al. Aorto-iliac Reconstruction. *Surgical Clinics of North America*, Vol. 59:4, 565-579, August 1979.
8. Bengoechea, E., et al. Extensive Endarterectomy of the Aorta, Common and External Iliac Arteries, and Common Femoral Arteries by a Modified LeVeen Method. *Surgery*, 99:1, 537-548, May 1986.

**Symbol Legend • Légende des symboles • Zeichenerklärung • Significado de los símbolos
• Legenda dei simboli • Lijst van symbolen • Symbolforklaring • Förklaring till symbolen • Λεζάντα Συμβόλου**

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V27SL40x16_6_A4.1

Edwards

Edwards Lifesciences Services GmbH

Edisonstrasse 6
85716 Unterschleissheim
Germany

5/18
©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC

One Edwards Way
Irvine, CA 92614 USA

Made in USA

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

WEB IFU
10002946002 A
DOC-0088388 A