DIRECTORY

English	1	Português	10	Eesti	19
Français	2	Česky	11	Lietuvių	20
Deutsch	3	Magyar	12	Latviešu	21
Español	4	Polski	13	Türkçe	22
Italiano	5	Slovensky	14	Русский	23
Nederlands	6	Norsk	15	Srpski	24
Dansk	7	Suomi	16	中文	25
Svenska	8	Български	17	繁體中文 (台灣)	26
Ελληνικά	9	Română	18	한국어	27

English

Cholangiography Catheter For Single Use Only

Indications

The Edwards Lifesciences cholangiography catheter is indicated for use in specific diagnostic procedures providing a conduit for contrast media or saline to be introduced into the cystic and common bile ducts for cholangiographic assessment or irrigation.

The 80 cm length catheter is compatible with the laparoscopic cholecystectomy technique. It may be introduced via a percutaneous introducer sheath through an abdominal puncture or via a reducer through the trocar.

Surgical Benefits

1. The bevel on the catheter tip facilitates introduction into the duct.
2. The distal tip of the catheter is indicated by a black marking that extends 1.5 cm proximally from the tip.
3. The tip is designed to reduce accidental withdrawal during the operative procedure.
4. Syringes attach easily to a 3-way luer-lock hub. Dual ports connect two syringes: one for irrigating, the other for injecting dye.
5. Air bubbles are easily visualized through the clear vinyl body of the catheter.

6. 80 cm catheter length facilitates catheter placement through a trocar or an introducer to access cystic and common bile ducts.

Precautions

Some false-positive and false-negative cholangiograms will be encountered. A false-positive film is usually caused by an air bubble in the duct; thus care must be taken to remove all air bubbles from the system. If a bubble is inadvertently introduced, it should be aspirated and recovered within the syringe.

A false-negative cholangiogram may occur because a small stone is hidden as a result of too large a concentration of contrast medium in the duct system.

Instructions

Prepare the catheter by placing two 20 ml luer-lock syringes on the 3-way stopcock. To avoid entrapped air bubbles, draw saline and contrast media through the cholangiography catheter into their respective syringes. Position the catheter stopcock key exactly midway between the two ports to prevent air bubbles and fluid leakage from the catheter. To facilitate handling of the catheter and securing it in place, remove the syringes.

Isolate the cystic duct at its junction with the common duct. For surgical procedures place a right-angled gall bladder clamp on the cystic duct 1.5 to 2.0 cm from the common duct. Place a ligature of 2.0 silk around the cystic duct between the clamp and common duct and make an incision in the cystic duct. For surgical and laparoscopic procedures introduce the cholangiography catheter into the

cystic duct so that the catheter's terminal port is in close proximity to the common duct; this determination is facilitated by the black mark extending 1.5 cm from the distal tip. For surgical procedures securely tie the ligature around the cystic duct and catheter. The tie should be between the incision and the conical projection, which is near the beveled tip of the catheter.

For laparoscopic procedures secure the catheter in the cystic duct.

Reattach the syringes after the catheter is inserted and secured. Rule out the possibility of air bubbles by gently irrigating with saline. Bubbles can be visualized through the clear vinyl catheter body. The preferred method of operative cholangiography can now be carried out.

MRI Information

This product has not been tested for MRI compatibility.

Complications

Potential complications are those associated with cholangiography procedures and use of catheters in the biliary tract. These include laceration or perforation of the biliary tree, hemorrhage, infection, endothelial damage, scar tissue formation leading to strictures or obstruction, pancreatitis, cholangitis, transient hyperamylasemia, and transient low-grade fever.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.
Temperature Limitation: 0° - 40 °C,
Humidity Limitation: 5% - 90% RH.

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.
Note: Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: . 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize**

or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.
Price, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC