Edwards

DIRECTORY

English	1	Português.....	24	Eesti	45
Français	3	Česky.....	26	Lietuvių	47
Deutsch	6	Magyar.....	28	Latviešu	50
Español.....	9	Polski	31	Türkçe.....	52
Italiano.....	11	Slovensky.....	33	Русский	54
Nederlands.....	14	Norsk.....	36	Srpski	57
Dansk	16	Suomi	38	中文.....	59
Svenska	19	Български.....	40	繁體中文 (台灣)	61
Ελληνικά.....	21	Română	43	한국어	63

English

VAMP

Closed Blood Sampling System

Instructions for Use For Single Use Only

For figures 1 through 8 please refer to pages 66 and 71.

Carefully read all instructions for use, warnings, and precautions before using the VAMP closed blood sampling system.

Description

The Edwards Lifesciences VAMP closed blood sampling system when attached to pressure monitoring lines, provides a safe and convenient method for the withdrawal of blood samples. The blood sampling system is designed for use with all disposable and reusable pressure transducers and for connection to central line catheters and arterial catheters. The VAMP closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Indications

To be used only for blood withdrawal.

The blood sampling system is indicated for use on patients requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters

and central venous catheters, which are attached to pressure monitoring lines.

Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

Warnings

This device is designed, intended, and distributed for single use only. Do not resterilize or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Instructions for Use

Caution: The use of lipids with the VAMP closed blood sampling system may compromise product integrity.

Equipment

1. Flush device or flow controlling device (maximum 4 mL/hr flow rate)
2. Disposable or reusable pressure transducer, if desired
3. VAMP closed blood sampling system that contains one 5 mL VAMP reservoir with integral shut-off valve and one VAMP needleless sampling site

Setup

1. Using aseptic technique, remove the VAMP kit from the sterile package.
2. Attach the distal end with the female luer-lock connector to the transducer stopcock on the disposable transducer or reusable transducer dome (see page 66, Fig. 1).

3. All connections should be secure.

Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.

4. Make sure that the reservoir plunger is down in the closed and locked position.
5. Deliver flush solution first through the transducer and out through the transducer vent port according to its manufacturer's instructions.
6. Replace all vented caps on the sideports of the stopcocks with non-vented caps.
7. To fill the VAMP kit, ensure that the shut-off valve is in the open position indicated by the handle being parallel to the tubing. Hold the reservoir by the plunger flexures. Orient the kit such that the reservoir is in the vertical position and the sampling site is above the reservoir at approximately 45°. In this position, gently squeeze and hold the plunger flexures together to slightly raise the seal to enlarge the fluid path and deliver flush solution slowly through the remainder of the kit. Once the kit is filled, release plunger flexures.

Caution: Remove all air bubbles to reduce the risk of air emboli.

8. Mount the transducer either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder.
9. Mount the reservoir either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder (see page 67, Fig. 2).
10. Pressurize the IV solution bag. Flow rate will vary with pressure across the flush device.

Edwards, Edwards Lifesciences, the stylized E logo, and VAMP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

11. Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter (see page 66, Fig. 1).
12. Zero and calibrate the transducer according to its manufacturer's instructions.

Drawing Blood Samples

Note: Use aseptic technique.

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

Although a variety of techniques can be used for drawing samples, the following guidelines are provided as an aid to the clinician:

Two methods may be used to draw blood samples using the VAMP blood sampling system. Method one uses a sampling syringe with the VAMP needleless cannula and a blood transfer unit (BTU). Method two, direct-line sampling, uses a VAMP direct-draw blood transfer unit (BTU) with a VAMP needleless cannula.

Drawing Blood Samples Using Method One (Syringe and Cannula)

1. Gently squeeze the reservoir plunger flexures together. Pull up on the reservoir plunger until it stops and the reservoir has reached its 5 mL volume capacity (see page 68, Fig. 3A).

Note: If difficulties are experienced in drawing the clearing sample, check the catheter for possible occlusions or restrictions (e.g., positional lines).

Important: Recommended time to pull reservoir plunger up to the fully open position is 3 to 5 seconds.

2. Once the clearing sample has been drawn, close the reservoir shut-off valve by turning the handle perpendicular to the tubing. This will further ensure that the sample drawn is from the patient and not the reservoir.
3. Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy. **Note:** Do not use acetone.
4. To draw a blood sample, use either a VAMP needleless cannula (packaged separately) and syringe, or a packaged VAMP needleless cannula preassembled onto a syringe.

Caution: Do not use a needle through the sampling site.

- a. If an individually packaged VAMP needleless cannula is used:
 - I. Using aseptic technique, peel open the cannula pouch.
 - II. Pick up the cannula by the protective shield (see page 68, Fig. 4).
 - III. Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure.
 - IV. Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
 - V. Push the cannula into the VAMP needleless sampling site and then

draw the required volume of blood into the syringe for sample (see page 69, Fig. 5).

Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.

- VI. Remove the syringe and cannula from the sampling site by **pulling straight out**.
- Caution:** Do not twist the syringe out of the sampling site.
- b. If a VAMP needleless cannula preassembled on a syringe is used:
 - I. Using aseptic technique, peel open the pouch.
 - II. Ensure that the cannula is securely tightened to the syringe.
 - III. Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
 - IV. Push the cannula into the VAMP needleless sampling site and then draw the required volume of blood into the syringe for sample (see page 69, Fig. 5).
- Note:** If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.
- V. Remove the syringe and cannula from the sampling site by **pulling straight out**.
- Caution:** Do not twist the syringe out of the sampling site.
5. Once the sample has been drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing.
6. Push down on the plunger until the flexures lock in place in the fully closed position and all fluid has been reinfused into the line (see page 68, Fig. 3B).

Important: Recommended time to push the reservoir plunger to the fully closed position is 3 to 5 seconds.

Important: The clearing sample should not remain in the reservoir for longer than 3 minutes.

7. Flush the line clear and swab the sampling site, ensuring removal of any excess blood left on the sampling port.

Note: Do not use acetone.

Warning: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

8. To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU) (see page 70, Fig. 6).
 - a. Using aseptic technique, peel open the pouch.
 - b. Ensure that all connections are tight.

- c. Hold the VAMP BTU in one hand and push the cannula on the filled sample syringe through the VAMP BTU needleless injection site.
- d. Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk on the vacuum tube.
- e. Fill the vacuum tube to the desired volume.
- f. Repeat steps (d) and (e) according to the requirements for the patient's blood study.
- g. Discard the VAMP BTU after transferring the blood sample from the syringe into the vacuum tubes.
9. Discard all syringes and cannulae after use according to hospital policy.

Drawing Blood Samples Using Method Two (Direct-Draw)

1. Gently squeeze the reservoir plunger flexures together. Pull up on the reservoir plunger until the plunger stops and the reservoir has reached its 5 mL volume capacity (see page 68, Fig. 3A).

Caution: If difficulties are experienced in drawing the clearing sample, check the catheter for possible occlusions or restrictions (e.g., positional lines).

Important: Recommended time to pull reservoir plunger up to the fully open position is 3 to 5 seconds.

2. Once the clearing sample has been drawn, close the reservoir shut-off valve by turning the handle 90° so that it is perpendicular to the tubing. This will further ensure that the sample drawn is from the patient and not the reservoir.
3. Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. **Note:** Do not use acetone.
4. To draw a blood sample, use the VAMP direct-draw blood transfer unit.

Caution: Do not use a needle through the sampling site.

- a. Using aseptic technique, peel open the pouch.
- b. Ensure that the cannula is securely tightened to the BTU housing.
- c. Position the patient's arm such that the sample site faces upward (see page 71, Fig. 7).

Note: For VAMP closed blood sampling systems that contain small bore tubing the sample site may not be located on the patient's arm (see page 71, Fig. 8).

- d. Push the cannula of the direct-draw BTU into the sampling site.
- e. Insert the selected vacuum tube into the open end of the direct-draw BTU and push until the internal needle of the BTU has punctured the rubber disk on the vacuum tube.

Caution: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity.

- f. Fill the vacuum tube to the desired volume.

Caution: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.

- g. Repeat steps (e) and (f) according to the requirements of the patient's blood study.
- h. When the last sample has been drawn, remove the vacuum tube first and then grasp the VAMP direct-draw BTU by the cannula and **pull straight out**.

Caution: Do not twist the BTU housing or remove it with the vacuum tube still attached.

5. Discard the VAMP direct-draw BTU after use according to hospital policy.
6. Once the sample has been drawn, open the shut-off valve by turning the handle 90° so that the handle is parallel to the tubing.
7. Push down on the plunger until the flexures lock in place in the fully closed position and all fluid has been reinfused into the line (see page 68, Fig. 3B).

Important: Recommended time to push the reservoir plunger to the fully closed position is 3 to 5 seconds.

Important: The clearing sample should not remain in the reservoir for longer than 3 minutes.

8. Flush the line clear and swab the sampling site, ensuring removal of any excess blood left

on the sampling port. **Note:** Do not use acetone.

Warning: Laboratory values should correlate with the patient's clinical manifestations. **Verify accuracy of laboratory values before instituting therapy.**

Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

MRI Safety Information

The VAMP closed blood sampling system is MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize. This product is for single use only.

Storage

Store in a cool, dry place.

Temperature Limitation: 0 °C - 40 °C
Humidity Limitation: 5 % - 90 % RH

Shelf Life

The recommended shelf life is as marked on each package.

Resterilization will not extend the shelf life of this product.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences AG
Edwards Lifesciences LLC