DIRECTORY

English	1	Português.....	26	Eesti	51
Français	4	Česky.....	29	Lietuvių	54
Deutsch	7	Magyar.....	32	Latviešu	56
Español.....	10	Polski	34	Türkçe.....	59
Italiano.....	13	Slovensky.....	37	Русский	62
Nederlands.....	15	Norsk.....	40	Srpski	65
Dansk	18	Suomi	43	中文.....	68
Svenska	21	Български.....	45	繁體中文 (台灣)	70
Ελληνικά.....	23	Română	48	한국어	73

English

VAMP Closed Blood Sampling System with LASS

Instructions for Use For Single Use Only

For figures 1 through 8 please refer to pages 76 through 80.

Carefully read all instructions for use, warnings, and precautions before using the needleless VAMP closed blood sampling system.

Description

The Edwards Lifesciences VAMP closed blood sampling system provides a safe and convenient method for the withdrawal of blood samples from pressure monitoring lines. The blood sampling system is designed for use with disposable pressure transducers and for connection to central line, venous, and arterial catheters where the system can be flushed clear after sampling. The VAMP closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Indications

To be used only for blood withdrawal.

The blood sampling system is indicated for use on patients requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

Edwards, Edwards Lifesciences, the stylized E logo, TruWave, VAMP, VAMP Flex, and VAMP Plus are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

Warnings

This device is designed, intended, and distributed for single use only. Do not resterilize or reuse this device. There is no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Instructions for Use

Caution: The use of lipids with the VAMP closed blood sampling system may compromise product integrity.

Equipment

1. Flush device for flow control (nominal 3 mL/hr flow rate)
2. Disposable pressure transducer
3. A VAMP closed blood sampling system that contains one reservoir with shut-off valve, and sampling site(s)

Setup

1. Using aseptic technique, remove the VAMP kit from the sterile package.
2. If the VAMP kit does not include a pre-connected TruWave disposable pressure transducer, remove the protective caps and attach the distal end with the female luer port to the male luer of a transducer, or other fluid flush device.
3. Ensure that all connections are secure and that the sampling site is in the prime/clear position (see page 76, Fig. 1A).

Note: Wet connections promote over-tightening by lubricating the fittings.

Over-tightened connections may result in cracks or leaks.

4. **For the VAMP Flex reservoir**, turn the reservoir stopcock off to the patient end of the tubing (see page 77, Fig. 2).
5. **For the VAMP Adult reservoir**, raise the plunger to facilitate flow of priming solution (see page 78, Fig. 4A).
6. **For the VAMP Plus reservoir**, push the reservoir plunger to the closed and locked position by pressing on the reservoir extension until it latches in the closed position (see page 77, Fig. 3).
7. Connect the kit to an IV fluid source without pressurizing the fluid. Begin to gravity fill the kit first through the transducer and out through the vent port according to the manufacturer's instructions.
8. Replace all vented caps on the side-ports of the stopcocks with non-vented caps.
9. **To fill the VAMP Flex reservoir:**
 - a. Slowly fill the reservoir to the full 10 mL position by pulling back on the plunger.
 - b. Turn the reservoir stopcock off to the transducer/IV bag (see page 79, Fig. 5).
 - c. Ensure that the sample site(s) are in the prime/clear position (see page 76, Fig. 1A).
 - d. Orient the reservoir such that it is in the vertical position and the sampling site is above the reservoir at approximately 45°, and slowly reinfuse the fluid. Repeat as necessary to remove air bubbles.
 - e. Turn the stopcock off to the reservoir (see page 79, Fig. 6).
10. **To fill the VAMP Adult reservoir:**
 - a. Ensure that the shut-off valve is in the open position indicated by the handle

being parallel to the tubing and that the sample site(s) are in the prime/clear position (see page 76, Fig. 1A).

- b. Hold the reservoir by the plunger flexures. Orient the reservoir in the vertical position and the sampling site above the reservoir at approximately 45°.
- c. Provide flow by pulling the Snap-Tab of the Edwards' TruWave disposable pressure transducer. Deliver priming solution to remove air. Lower plunger to closed position (see page 78, Fig. 4B).

11. To fill the VAMP Plus reservoir:

- a. Ensure that the shut-off valve is in the open position indicated by the handle being parallel to the tubing.
- b. Ensure that the sample site(s) are in the prime/clear position (see page 76, Fig. 1A).
- c. Orient the reservoir in the vertical position and the sampling site above the reservoir at approximately 45°.
- d. Provide flow by pulling the Snap-Tab of the Edwards' TruWave disposable pressure transducer. Deliver priming solution to remove air.

Note: The transducer and reservoir with bracket can be mounted to an IV pole using the Edwards' disposable pressure transducer holder.

12. Position the LASS sampling site(s) to the pressure monitoring position (see page 76, Fig. 1B) and pressurize the IV solution bag per hospital policy and transducer manufacturer's instructions. Flow rate will vary with pressure gradient across the flush device.
13. Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter.
14. Connect pressure transducer cable to monitor.
15. Zero the transducer according to its manufacturer's instructions.

Caution: Remove all air bubbles to reduce the risk of air emboli and reduce loss of pressure signal.

Caution: Failure to position LASS sampling site in pressure monitoring position (see page 76, Fig. 1B) will adversely impact dynamic response of the pressure signal.

Clearing IV Fluid in Preparation for Drawing Blood Samples

Note: Use aseptic technique.

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume is required for coagulation studies.

1. Position sampling site into the prime/clear position (see page 76, Fig. 1A).
2. **For VAMP Flex closed blood sampling system**, turn the reservoir stopcock off to the transducer/IV bag (see page 79, Fig. 5). Fill the reservoir by drawing the plunger. Draw the appropriate clearing volume according to hospital policy. (Recommended draw rate is

1 second for each mL.) Once the clearing volume is drawn, turn the reservoir stopcock off to the patient and catheter (see page 77, Fig. 2). This will ensure that the sample drawn is from the patient and not the reservoir.

3. **For VAMP Adult reservoir**, firmly squeeze the flexures and slowly fill the reservoir at its 5 mL volume capacity. (Recommended draw rate is 1 second for each mL.) Close the shut-off valve by turning the handle perpendicular to the tubing.
4. **For VAMP Plus closed blood sampling system**, release plunger latch and draw back on the reservoir plunger until it stops at its 12 mL volume capacity. (Recommended draw rate is 1 second for each mL.) Close the shut-off valve by turning the handle perpendicular to the tubing.

Note: If difficulties are experienced in drawing the clearing sample, check the catheter for possible occlusions or restrictions.

Drawing Blood Samples Using Syringe or Direct-Draw Method

Caution: Do not use a needle or needleless cannula through the sampling site.

Syringe Method

1. Always use ISO-594 or EN 20594-1 compliant syringes to avoid potential collapsing of the septum after syringe insertion and removal.
2. Prepare the VAMP kit for drawing a blood sample per the method outlined above.
3. Ensure that the LASS sampling site is still in the prime/clear position (see page 76, Fig. 1A).
4. Disinfect the sampling site for access per hospital policy.
5. Carefully connect the syringe to the LASS sampling site by slowly pushing the syringe **straight** into the site using a clockwise, rotating motion until it is fully seated. (See page 80, Fig. 7 (A and B) and Fig. 8 (A and B) for proper syringe insertion techniques.)

Note: Do not use acetone.

- **Do Not** try to insert the syringe at an angle.
 - **Do Not** try to pry open the septum's slit in the sampling site. The slit automatically opens when the syringe is connected.
 - **Do Not** use excessive force to connect the syringe.
6. Remove syringe from the LASS sampling site, rotating counter clockwise. The blood sampling site's septum completely closes after each use.

Direct-Draw Method

1. Always use ISO-594 or EN 20594-1 compliant blood collection unit to avoid potential collapsing of the septum after unit insertion and removal.
2. Prepare the VAMP kit for drawing a blood sample per the method outlined above.
3. Ensure that the sampling site is still in the prime/clear position (see page 76, Fig. 1A).

4. Disinfect the sampling site for access per hospital policy.

Note: Do not use acetone.

5. Carefully connect the direct-draw unit to the LASS sampling site by slowly pushing the unit **straight** into the site using a clockwise, rotating motion until it is fully seated.
 - **Do Not** try to insert the direct-draw unit at an angle.
 - **Do Not** try to pry open the septum's slit in the sampling site. The slit automatically opens when the direct-draw unit is connected.
 - **Do Not** use excessive force to connect the direct-draw unit.
6. Insert the vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube. Fill to the desired volume.

Caution: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity.

7. Repeat step 6 as required to fill additional vacuum tubes.

Caution: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.
8. When the last sample is drawn, remove the vacuum tube first.
9. Remove direct-draw unit from sampling site, rotating counter clockwise. The blood sampling site's septum completely closes after each use.
10. Discard the direct-draw unit after use according to hospital policy.

Caution: The clearing sample should not remain in the reservoir for longer than 2 minutes.

Reinfusing Clearing Volume

1. For VAMP Flex sampling system:

- a. Once the sample is drawn, turn the reservoir stopcock off to the transducer/IV bag (see page 79, Fig. 5).
- b. Ensure that the sample site(s) are in the prime/clear position (see page 76, Fig. 1A). Reinfuse the clearing volume by pushing the plunger down.
- c. If blood has been drawn into the VAMP Flex reservoir, extra steps might be necessary to clear the reservoir. To accomplish this, turn the reservoir stopcock off to the patient (see page 77, Fig. 2) and partially fill the reservoir with 1 to 2 mL of fluid. If the IV bag is properly pressurized, activating the flush device will automatically cause the reservoir to fill. Turn the reservoir stopcock off to the transducer/IV bag (see page 79, Fig. 5) and slowly reinfuse into the line. Repeat as necessary until all blood is removed.

Note: For optimal clearing of the reservoir, it is recommended to repeat this

clearing process multiple times with small volume of fluid (i.e., 1 to 2 mL).

- d. Once the reservoir has been cleared, turn the stopcock off to the reservoir (see page 79, Fig. 6).

2. **For VAMP Adult closed blood sampling system:**

- a. Once the sample is drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing. Turn the handle on the sampling site to the prime/clear position (see page 76, Fig. 1A).
- b. Press together the reservoir flexures extension (see page 78, Fig. 4B) until the plunger stops and is in the closed position.

3. **For VAMP Plus closed blood sampling system:**

- a. Once the sample is drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing. Turn the handle on the sampling site to the prime/clear position (see page 76, Fig. 1A).
- b. Close the reservoir by pushing the plunger until it stops and firmly latches in the closed position.

Note: Recommended time to reinfuse the clearing volume is approximately 1 second for each mL.

4. Flush sample site(s) and line clear of residual blood by pulling the Snap-Tab of the Edwards' TruWave disposable pressure transducer.
5. Return the sampling site to the pressure monitoring position (see page 76, Fig. 1B).
6. Clean the sampling site, ensuring removal of any excess blood.

Note: Do not use acetone.

Caution: Failure to position sampling site in pressure monitoring position (see page 76, Fig. 1B) will adversely impact dynamic response of the pressure signal.

Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

MRI Safety Information

VAMP Flex closed blood sampling system and VAMP Adult closed blood sampling system

The VAMP Flex and VAMP Adult closed blood sampling system are MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

VAMP Plus closed blood sampling system

Warnings

For a patient undergoing an MRI examination, please refer to the MRI Safety Information section for specific conditions to ensure patient safety.

Non-clinical testing demonstrated that the VAMP Plus closed blood sampling system is MR Conditional and can be used in the MRI environment according to the following conditions:

<conditions for
MR safety>

- Static magnetic field of 3-Tesla
- Spatial gradient field of 720-Gauss/cm or less
- This device is not an implant. Therefore, image artifact, RF Heating and SAR are not applicable
- This device is not intended for use inside of the bore of the MR system and should not be in contact with the patient

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Plus closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

This product is for single use only.

Storage

Store in a cool, dry place.

Temperature Limitation: 0 °C - 40 °C

Humidity Limitation: 5% - 90% RH

Shelf Life

The recommended shelf life is as marked on each package. Resterilization will not extend the shelf life of this product.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences AG
Edwards Lifesciences LLC