DIRECTORY

English	1	Português	12	Eesti	22
Français	2	Česky	13	Lietuvių	23
Deutsch	4	Magyar	14	Latviešu	24
Español	5	Polski	15	Türkçe	25
Italiano	6	Slovensky	17	Русский	26
Nederlands	7	Norsk	18	Srpski	28
Dansk	8	Suomi	19	中文	29
Svenska	10	Български	20	繁體中文 (台灣)	30
Ελληνικά	11	Română	21	한국어	31

English

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Fogarty Occlusion Catheter

For Single Use Only

Indications

Fogarty Occlusion Catheters are indicated for temporary vessel occlusion.

The Large Occlusion catheters are intended to be used for temporary occlusion in the aorta, vena cava and internal jugular vein.

The Small Occlusion catheters are intended to be used for temporary occlusion in the peripheral vascular system.

Contraindications

Fogarty Occlusion Catheters are not indicated for endarterectomy procedures or vessel dilation.

Warnings

- Balloon rupture and catheter separation are the most frequent causes of failure for vascular catheters. The possibility of balloon rupture must be taken into account when considering the risks involved in any occlusion procedure.
- To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume for each size catheter (see **Specification Table**).
- Exposure to calcified plaque within the vessel may increase the possibility of balloon rupture.

- Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.
- Air should not be used for inflation in instances where balloon rupture could produce a dangerous air embolus. If carbon dioxide is used, gas seepage through the intact balloon will require frequent adjustment of the inflation volume.

Instructions

Inspection

Prior to use, the catheter should be inspected with the balloon inflated to its recommended balloon inflation volume (see **Specification Table**). A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

Preparation

Obtain the smallest syringe that will hold the balloon's stated maximum fluid or gas volume. Fill the syringe to the volume indicated with a sterile, blood compatible medium. Fluid may include a highly diluted, nonparticulate, radiopaque solution. Connect the syringe to the gate valve.

Catheter Placement

Place the catheter into the desired position through an arteriotomy or venotomy.

Balloon Inflation

Ensure that the gate valve is in the open position. Slowly inject the inflation medium until occlusion is obtained. Close the gate valve to maintain occlusion. Fluoroscopic visualization is recommended where appropriate to ensure proper placement and inflation.

Caution: Overinflation of the balloon may cause vessel damage.

The fluid in the syringe should be checked before each inflation. If the amount exceeds the stated capacity, remove the syringe and refill to the proper volume. The amount of fluid remaining in the catheter must be taken into account.

Withdrawing the Catheter

To deflate the balloon, open the gate valve and pull back on the plunger of the syringe. Withdraw the catheter from the insertion point.

MRI Information

This product has not been tested for MRI compatibility.

Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection, local hematoma, intimal disruption, vessel dissection, perforation and rupture, hemorrhage, thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm,

Specification Table

Model Number	620403F	620404F	620405F	62080814F	62080822F
Maximum Liquid Capacity (ml)	0.2	0.75	1.5	10	43
Maximum Gas Capacity (ml)	0.6	1.7	3.0	20	50
Diameter of Inflated Balloon (mm)	5	9	11	28	45
Maximum French Size of Deflated Balloon	4.3F (1.43 mm)	5.0F (1.67 mm)	6.0F (2.0 mm)	14F (4.7 mm)	22F (7.3 mm)
Catheter French Size	3F (1.00 mm)	4F (1.33 mm)	5F (1.67 mm)	8F (2.7 mm)	8F (2.7 mm)
Length (cm)	40	40	40	80	80

arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

Storage

Store in a cool, dry place.

Temperature Limitation: 0° - 40 °C,

Humidity Limitation: 5% - 90% RH.

Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the recommended time may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours):800.822.9837
Outside the U.S. and Canada
(24 hours):949.250.2222
In the UK:0870 606 2040 - Option 4
In Ireland:01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC