
Edwards

DIRECTORY

English	1	Português.....	25	Eesti	47
Français	4	Česky.....	27	Lietuvių	50
Deutsch	7	Magyar.....	30	Latviešu	52
Español.....	10	Polski	32	Türkçe.....	55
Italiano.....	12	Slovensky.....	35	Русский	57
Nederlands.....	15	Norsk	37	Srpski	60
Dansk	17	Suomi	40	中文.....	62
Svenska	20	Български.....	42	繁體中文 (台灣)	64
Ελληνικά.....	22	Română	45	한국어.....	67

English

VAMP Closed Blood Sampling System with Two Needleless Sampling Sites

Instructions for Use For Single Use Only

For figures 2 through 5 please refer to pages 70 through 73.

Carefully read all instructions for use, warnings, and precautions before using the VAMP closed blood sampling system.

Description

The Edwards Lifesciences VAMP closed blood sampling system provides a safe and convenient method for withdrawing blood samples. The VAMP system is designed for use with disposable and reusable pressure transducers and connection to central line and arterial catheters. The VAMP closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or

cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Indications

To be used only for blood withdrawal.

The blood sampling system is indicated for use on patients requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

Warnings

This device is designed, intended, and distributed for single use only. Do not resterilize or reuse this device. There are no

data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Caution: The use of lipids with the VAMP closed blood sampling system may compromise product integrity.

Instructions for Use

Equipment

1. Flush device or flow controlling device (maximum 4 mL/hour flow rate).
2. Disposable or reusable pressure transducer, if desired.
3. VAMP closed blood sampling system with integral shut-off valve and two VAMP needleless sampling sites.

Setup

1. Using aseptic technique, remove the VAMP kit from the sterile package.
2. Attach the distal end with the female luer-lock connector to the transducer stopcock on the disposable transducer or reusable transducer dome, or to an IV line/IV infusion pump. Ensure that all connections are secure.

Edwards, Edwards Lifesciences, the stylized E logo and VAMP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

3. If a transducer is in the line, deliver flush solution first through the transducer and out through the transducer vent port according to the manufacturer's instructions.
4. Replace all vented caps on the sideports of the stopcocks with non-vented caps.
5. Fill the VAMP system with flush solution using gravity and ensure that the shut-off valve by the clearing volume (distal) access site is in the open position by turning the handle parallel to the tubing. Deliver flush solution slowly through the kit.
Caution: Remove all air bubbles to reduce the risk of air emboli.
6. If a transducer is being used, mount the transducer either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder.
7. Pressurize the IV solution bag. Flow rate will vary with pressure across the flush device.
8. Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter (Figure 1).
9. If necessary, zero and calibrate the transducer according to the transducer manufacturer's instructions.

Drawing the Clearing Volume and Blood Samples

Note: Use aseptic technique.

Drawing Clearing Volume

1. To draw clearing volume, use either a VAMP needleless cannula (packaged separately) and syringe, or a packaged VAMP needleless cannula preassembled onto a syringe.
Caution: Do not insert a needle into the clearing volume (distal) access site.
 - a. Using individually packaged VAMP needleless cannula:
 - i. Using aseptic technique, peel open the cannula pouch.
 - ii. Pick up the cannula by the protective shield (see page 70, Fig. 2).
 - iii. Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure.
 - iv. Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
2. Swab the clearing volume (distal) access site with disinfectant such as alcohol or betadine, depending upon hospital policy.
Note: Do not use acetone.
3. Push the cannula into the VAMP needleless clearing volume (distal) access site, hold in place for approximately 1 to 2 seconds, and draw the clearing volume (see page 71, Fig. 3).

Important: Recommended rate to draw the clearing volume into the syringe is 1 cc/second.

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

Note: If difficulties are experienced in drawing the clearing volume, check the catheter and Edwards kit for possible occlusions or restrictions.

4. Once the clearing volume has been drawn, close the shut-off valve at the clearing volume (distal) access site by turning the handle perpendicular to the tubing. This will ensure that the sample is drawn from the patient and not the distal tubing.
5. While holding the cannula, remove the syringe and cannula from the clearing volume access site by **pulling straight out**.
Caution: Do not twist the syringe out of the sampling site.
Important: Discard the clearing volume according to hospital protocol.
6. Swab the access site to ensure removal of any excess blood.
Note: Do not use acetone.

Drawing Blood Samples from the VAMP Needleless Sampling Site

Although a variety of techniques can be used for drawing samples, the following guidelines are provided as an aid to the clinician:

Two methods may be used to draw blood samples from the proximal VAMP needleless sampling site. Method one uses a sampling syringe with the VAMP needleless cannula. Method two, direct line sampling, uses a VAMP direct-draw unit with an integral VAMP needleless cannula.

Drawing Blood Samples Using Method One (Syringe and Cannula)

1. Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy.
Note: Do not use acetone.
2. To draw a blood sample, use a cannula (packaged separately) and syringe.
Caution: Do not use a needle through the sampling site.
3. Using an individually packaged cannula:
 - a. Using aseptic technique, open the cannula package.
 - b. Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure.
4. Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
5. Push the cannula into the proximal VAMP needleless sampling site and hold in place for approximately 1 to 2 seconds (see page 71, Fig. 3).
6. Draw the required volume of blood into the syringe.

Note: If difficulties are experienced in drawing the sample, check the catheter and VAMP kit for possible occlusions or restrictions.

7. While holding the cannula, remove the syringe and cannula from the sampling site by **pulling straight out**.

Caution: Do not twist the syringe out of the sampling site.

8. Once the last sample has been drawn, swab the sampling site to ensure removal of any residual blood remaining on the sampling port.

Note: Do not use acetone.

Transferring Blood Samples

1. To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU) (see page 72, Fig. 4).
 - a. Using aseptic technique, peel open the pouch.
 - b. Ensure that all connections are tight.
 - c. Hold the VAMP blood transfer unit in one hand and push the cannula on the filled sample syringe through the VAMP blood transfer unit needleless injection site.
 - d. Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk on the vacuum tube.
 - e. Fill the vacuum tube to the desired volume.
 - f. Repeat steps (d) and (e) according to the requirements for the patient's blood study.
2. According to hospital policy, discard the VAMP blood transfer unit, syringes, and cannulae after transferring the blood sample from the syringe into the vacuum tubes.

Drawing Blood Samples Using Method Two (Direct-Draw Method)

1. Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy.
Note: Do not use acetone.
2. To draw a blood sample, use the VAMP direct-draw unit.

Caution: Do not use a needle through the sampling site.

- a. Using aseptic technique, peel open the pouch.
- b. Ensure that the cannula is securely tightened to the direct-draw housing.
- c. Position the sampling site so that it faces upward.
- d. Push the cannula of the direct-draw unit into the sampling site (see page 73, Fig. 5).
- e. Insert the selected vacuum tube into the open end of the VAMP direct-draw unit and push until the internal needle of the VAMP direct-draw unit has punctured the rubber disk on the vacuum tube.

Caution: To prevent backflow of vacuum tube contents (including air) from entering

fluid path, remove vacuum tube before reaching maximum fill capacity.

- f. Fill the vacuum tube to the desired volume.

Note: If difficulties are experienced in drawing the sample, check the catheter and VAMP kit for possible occlusions or restrictions.

- g. Repeat steps (e) and (f) according to the requirements of the patient's blood study.

- h. When the last sample has been drawn, remove the vacuum tube first and then grasp the VAMP direct-draw unit by the cannula and **pull straight out**.

Caution: Do not twist the direct-draw unit housing or remove it with the vacuum tube still attached.

- 3. Discard the VAMP direct-draw unit after use according to hospital policy.
- 4. Once the last sample has been drawn, swab the sampling site to ensure removal of any excess blood remaining on the sampling port.

Note: Do not use acetone.

Flushing the Line

- 1. After the sample has been drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing.
- 2. Flush the line clear.
- 3. Discard all syringes and cannulae after use according to hospital policy.

Warning: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

MRI Safety Information

The VAMP closed blood sampling system is MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize. This product is for single use only.

Storage

Store in a cool, dry place.
Temperature Limitation: 0 °C - 40 °C
Humidity Limitation: 5 % - 90 % RH

Shelf Life

The recommended shelf life is as marked on each package. Resterilization will not extend the shelf life of this product.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences AG
Edwards Lifesciences LLC