Edwards

## DIRECTORY

English .....	1	Português.....	23	Eesti .....	45
Français .....	3	Česky.....	25	Lietuvių .....	47
Deutsch .....	6	Magyar.....	28	Latviešu .....	49
Español.....	8	Polski .....	30	Türkçe.....	52
Italiano.....	11	Slovensky.....	33	Русский .....	54
Nederlands.....	13	Norsk .....	35	Srpski .....	57
Dansk .....	16	Suomi .....	37	中文.....	59
Svenska .....	18	Български.....	39	繁體中文 (台灣) .....	61
Ελληνικά.....	20	Română .....	42	한국어.....	63

### English

## VAMP Plus Closed Blood Sampling System

### Instructions for Use

#### For Single Use Only

For figures 1 through 3 please refer to page 66.

Carefully read all instructions for use, warnings, and precautions before using the VAMP Plus closed blood sampling system.

#### Description

The Edwards Lifesciences VAMP Plus closed blood sampling system provides a safe and convenient method for the withdrawal of blood samples from pressure monitoring lines. The blood sampling system is designed for use with disposable and reusable pressure transducers and for connection to central line catheters, venous, and arterial catheters where the system can be flushed clear after sampling. The VAMP Plus closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

#### Indications

To be used only for blood withdrawal.

The blood sampling system is indicated for use on patients requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

Edwards, Edwards Lifesciences, the stylized E logo, TruWave, VAMP, and VAMP Plus are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

#### Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

#### Warnings

- **This device is designed, intended, and distributed for single use only. Do not resterilize or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.**
- **For a patient undergoing an MRI examination, please refer to the MRI Safety Information section for specific conditions to ensure patient safety.**

#### Instructions for Use

**Caution:** The use of lipids with the VAMP Plus closed blood sampling system may compromise product integrity.

#### Equipment

1. Flush device or flow controlling device (maximum 4 mL/hr flow rate)
2. Disposable or reusable pressure transducer, if desired
3. VAMP Plus closed blood sampling system

#### Setup

1. Using aseptic technique, remove the VAMP Plus kit from the sterile package.
2. If the VAMP Plus kit does not include a preconnected TruWave disposable pressure transducer, remove the protective caps and attach the distal end with the female luer-lock connector to the male luer of a transducer, or other fluid flush device (see page 66, Fig. 1).

3. All connections should be secure.

**Note:** Wet connections promote over tightening by lubricating the fittings. Over tightened connections may result in cracks or leaks.

4. Push the reservoir plunger to the closed and locked position by pressing on the reservoir extension until it latches in the closed position (see page 66, Fig. 3).
5. Connect the kit to an IV fluid source without pressurizing the fluid. Begin to gravity fill the kit first through the transducer and out through the vent port according to its manufacturer's instructions.
6. Replace all vented caps on the sideports of the stopcocks with non-vented caps.
7. To fill the VAMP Plus kit, ensure that the shut off valve is in the open position indicated by the handle being parallel to the tubing. Orient the kit such that the reservoir inlet and outlet ports are in an upward vertical position and the sampling sites are above the reservoir at approximately 45°. Slowly fill and de-bubble the reservoir and each sample site in succession.

**Caution:** Remove all air bubbles to reduce the risk of air emboli and reduce loss of pressure signal.

8. The transducer and reservoir with bracket can be mounted to an IV pole using the Edwards disposal pressure transducer holder (see page 66, Fig. 1). For convenience, the sample site nearest the reservoir can be easily mounted to the holder.
9. Pressurize the IV solution bag per hospital policy and transducer manufacturer's instructions. Flow rate will vary with pressure gradient across the flush device.
10. Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter.

11. Zero the transducer according to its manufacturer's instructions.

### Clearing IV Fluid in Preparation for Drawing Blood Samples

Two methods may be used for drawing the clearing volume into the VAMP Plus reservoir.

**Important:** A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

#### One-Hand Pole Mount Method to Draw Clearing Volume

1. With the VAMP Plus reservoir bracket mounted in the Edwards disposable pressure transducer holder, release the plunger latch. Fill the reservoir by drawing the plunger to the open position by pressing together on the plunger and the bracket extension finger hold (see page 66, Fig. 2) until the plunger stops and the reservoir has reached its 12 mL volume capacity.

**Note:** If difficulties are experienced in drawing the clearing volume, check the catheter for possible occlusions or restrictions.

**Note:** Recommended time to fill the reservoir is approximately 1 second for each mL of clearing volume.

2. Once the clearing volume has been drawn, close the reservoir shut off valve by turning the handle perpendicular to the tubing. This will ensure that the sample drawn is from the patient and not the reservoir.
3. Draw the sample according to one of the methods outlined in the **Drawing Blood Samples** section below.
4. Once the sample has been drawn, open the shut off valve by turning the handle 90° so that it is parallel to the tubing.
5. Reinfuse the clearing volume back to the patient by returning the reservoir plunger to the closed position. Press together the reservoir extension finger hold and the plunger pad (see page 66, Fig. 3) until the plunger stops and firmly latches in the closed position.

**Note:** If difficulties are experienced in reinfusing the clearing volume, check the catheter for possible occlusions or restrictions.

**Note:** Recommended time to push the reservoir plunger to the fully closed position is approximately 1 second for each mL of clearing volume.

**Caution:** The clearing sample should not remain in the reservoir for longer than 2 minutes.

6. Flush the reservoir, sample sites, and line clear of residual blood. Swab the sampling sites, ensuring removal of any excess blood.

**Note:** Do not use acetone.

#### Hand Held Syringe Method to Draw Clearing Volume

1. Remove the VAMP Plus reservoir from the bracket. Hold the reservoir in one hand and release the plunger latch. Fill the reservoir by drawing the plunger to the open position by either pushing the plunger open with your thumb, or pulling with your other hand, until

the plunger stops and the reservoir has reached its 12 mL capacity.

**Note:** If difficulties are experienced in drawing the clearing volume, check the catheter for possible occlusions or restrictions.

**Note:** Recommended time to fill the reservoir is approximately 1 second for each mL of clearing volume.

2. Once the clearing volume has been drawn, close the reservoir shut off valve by turning the handle perpendicular to the tubing. This will ensure that the sample drawn is from the patient and not the reservoir.
3. Draw the sample according to one of the methods outlined in the **Drawing Blood Samples** section below.
4. Once the sample has been drawn, open the shut off valve by turning the handle 90° so that it is parallel to the tubing.
5. Reinfuse the clearing volume back to the patient by pushing the plunger until it stops and firmly latches in the closed position.

**Note:** If difficulties are experienced in reinfusing the clearing volume, check the catheter for possible occlusions or restrictions.

**Note:** Recommended time to push the reservoir plunger to the fully closed position is approximately 1 second for each mL of clearing volume.

**Caution:** The clearing sample should not remain in the reservoir for longer than 2 minutes.

6. Flush the reservoir, samples sites, and line clear of residual blood. Swab the sampling site(s), ensuring removal of any excess blood left on the sampling port.

**Note:** Do not use acetone.

#### Drawing Blood Samples

Two methods may be used to draw blood samples using the VAMP Plus closed blood sampling system.

##### Syringe Method for Drawing Blood Samples

1. Prepare the VAMP Plus reservoir for drawing a blood sample by drawing a clearing volume per one of the methods outlined above.
2. Once the line is clear, swab the VAMP Plus needleless sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy.

**Note:** Do not use acetone.

3. To draw a blood sample, use either a preassembled packaged VAMP needleless cannula and syringe or a VAMP needleless cannula (packaged separately) and separate syringe.

**Caution:** Do not use a needle through the sampling site.

4. Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
5. Push the cannula into the VAMP Plus needleless sampling site and then draw the required volume of blood into the syringe.

**Note:** If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.

6. Remove the syringe and cannula from the sampling site by **pulling straight out**.

**Caution:** Do not twist the syringe out of the sampling site.

7. Once the sample has been drawn, reinfuse the clearing volume according to one of the methods described above.

**Warning:** Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

8. To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU).
  - a. Using aseptic technique peel open the pouch.
  - b. Ensure that all connections are tight.
  - c. Hold the VAMP BTU in one hand and push the cannula on the filled sample syringe through the VAMP BTU needleless injection site.
  - d. Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk of the vacuum tube.
  - e. Fill the vacuum tube to the desired volume.
  - f. Repeat steps (d) and (e) as required to fill additional tubes.
  - g. Discard the VAMP BTU after transferring the blood sample from the syringe into the vacuum tubes.
9. Discard all syringes and cannulae after use according to hospital policy.

##### Direct-Draw Method for Drawing Blood Samples

1. Prepare the VAMP Plus reservoir for drawing a blood sample by drawing a clearing volume per one of the methods outlined above.
2. Once the line is clear, swab the VAMP Plus needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy.

**Note:** Do not use acetone.

3. To draw a blood sample, use the VAMP direct-draw unit.

**Caution:** Do not use a needle through the sampling site.

- a. Using aseptic technique, peel open the pouch.

- b. Ensure that the cannula is securely tightened to the direct-draw unit.

- c. Push the cannula of the direct-draw unit into the sampling site.

- d. Insert the selected vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube and fill to the desired volume.

**Caution:** To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity.

**Caution:** If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.

- e. Repeat step (d) as required to fill additional vacuum tubes.
  - f. When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw unit by the cannula and **pull straight out**.
- Caution:** Do not twist the direct-draw unit housing or remove it with the vacuum tube still attached.
- g. Discard the VAMP direct-draw unit after use according to hospital policy.
4. Once the sample has been drawn, reinfuse the clearing volume according to one of the methods described above.

## Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

## MRI Safety Information

MR Conditional

Non-clinical testing demonstrated that the VAMP Plus closed blood sampling system is MR

Conditional and can be used in the MRI environment according to the following conditions:

- Static magnetic field of 3Tesla or less
- This device should remain outside of the bore of the MR system and should not be in contact with the patient during MR imaging
- Device should be secured (for example taped down) prior to MRI

**Precaution:** Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Plus closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

## How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

This product is for single use only.

## Storage

Store in a cool, dry place.

Temperature Limitation: 0 °C - 40 °C

Humidity Limitation: 5 % - 90 % RH

## Shelf Life

The recommended shelf life is as marked on each package. Resterilization will not extend the shelf life of this product.

## Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): ..... 800.822.9837

Outside the U.S. and Canada

(24 hours): ..... 949.250.2222

In the UK: ..... 0870 606 2040 - Option 4

In Ireland: ..... 01 8211012 Option 4

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

## Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

**Refer to the symbol legend at the end of this document.**

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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©Copyright 2018, Edwards Lifesciences LLC
All rights reserved.
**Edwards Lifesciences Services GmbH**
**Edwards Lifesciences AG**
**Edwards Lifesciences LLC**