DIRECTORY

English.....	1	Português.....	26	Română.....	47
Français.....	4	Česky.....	28	Eesti.....	50
Deutsch.....	7	Magyar.....	31	Lietuvių.....	52
Español.....	10	Polski.....	34	Latviešu.....	55
Italiano.....	12	Slovensky.....	37	Türkçe.....	57
Nederlands.....	15	Norsk.....	39	Русский.....	60
Dansk.....	18	Suomi.....	42	Srpski.....	63
Svenska.....	20	Български.....	44	Hrvatski.....	65
Ελληνικά.....	23				

English

TruWave

Pressure Monitoring Kit with TruWave Disposable Pressure Transducer

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

These are general instructions for setting up a Pressure Monitoring system with the Edwards Lifesciences TruWave disposable pressure transducer. Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

CAUTION: The use of lipids with the TruWave disposable pressure transducer may compromise product integrity.

1.0 Description

The Edwards Lifesciences Pressure Monitoring Kit with TruWave disposable pressure transducer, is a sterile, single-use kit that monitors pressures. The disposable sterile cable (available in 12-inch/30 cm and 48-inch/120 cm lengths) interfaces exclusively with an Edwards Lifesciences cable that is specifically wired for the monitor being used.

The TruWave disposable pressure transducer has a straight, flow-through design across the pressure sensor, and is available with or without the integral flush device.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The known benefits of pressure monitoring include providing information regarding the hemodynamic status of the

patient which may lead to improved data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation. In addition, VAMP blood conservation technology reduces unnecessary blood loss.

This device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and clinical usage of arterial line pressure monitoring as part of their respective institutional guidelines.

2.0 Intended Use/Purpose

The TruWave disposable pressure transducer is intended to monitor pressure when attached to a pressure monitoring catheter and compatible pressure monitoring instrument.

3.0 Indications

The Pressure Monitoring Kit with TruWave disposable pressure transducer is for use on patients requiring intravascular, intracranial, or intrauterine pressure monitoring.

4.0 Contraindications

There are no absolute contraindications for using the TruWave disposable pressure transducer in patients requiring invasive pressure monitoring.

Flush devices should not be used when monitoring intracranial pressure.

5.0 Warnings

- Do not use the flush device during intracranial pressure monitoring.
- To avoid over infusing severely fluid-restricted patients such as neonates and children, use an infusion pump in series with the flush device to accurately regulate the minimum amount of flush solution needed to maintain catheter patency while allowing continuous pressure monitoring. Use an IV set that is indicated for use with the IV pump. High pressures, which may be generated by an infusion pump at certain flow rates, may override the flush device restriction, resulting in fast flushing at the rate set by the pump.
- Do not allow air bubbles to enter the setup, especially when monitoring arterial pressures.
- If the flush system is not adequately pressurized relative to the patient's blood pressure, blood bleed-back and catheter clotting may occur.
- Avoid contact with any topical cream or ointment that attacks polymeric materials.
- Do not expose electrical connections to fluid contact.
- Do not autoclave the reusable cable.

- This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.
- Do not use this device if it has been exposed to an environment outside the range of 0 °C to 40 °C and 5 % to 90 % RH, or if it has gone beyond its Use By date.
- For patients undergoing an MR examination refer to the MRI Safety Information section for specific conditions.
- Some models may contain phthalates, specifically DEHP [Bis (2-ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

-  Type CF applied part

6.0 Equipment

The Edwards TruWave disposable pressure transducers are provided in pressure monitoring kits. Pressure monitoring kits may comprise of the following components or accessories:

- Pressure tubing
- Stopcocks
- Flush device (3 mL/hr or 30 mL/hr)
- IV set
- VAMP (closed blood sampling system)
- Disposable holder
- Arm strap
- Compatible monitor cable

Note: Refer to VAMP closed blood sampling system Instructions for Use for closed blood sampling instructions.

7.0 Procedure

Step	Procedure
1	Ensure that the cable is compatible with the monitor being used. Connect the reusable cable to the monitor. Turn the monitor on to allow the electronics to warm up. WARNING: Moisture within the connector may result in the device malfunctioning or in inaccurate pressure readings. If this occurs, replace wet part or parts. Note: Refer to compatible monitor instructions for proper cable to monitor connection.
2	Using aseptic technique, remove the transducer and kit from the sterile packaging.
3	Add additional components as needed to complete the monitoring system per hospital policy.
4	Ensure that all connections are secure. Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.
5	Connect the transducer cable to the reusable cable.
6	Remove all air from the IV flush solution bag, which is heparinized per hospital policy. CAUTION: If all air is not removed from the bag, air may be forced into the patient's vascular system when the solution is exhausted.
7	Close the roller clamp on the IV set and connect it to the IV flush bag. Hang the bag approximately 2 feet (60 cm) above the patient. This height will provide approximately 45 mmHg of pressure to the prime setup.
8	Fill the drip chamber halfway, or per hospital policy, with flush solution by squeezing the drip chamber. Open the roller clamp.
9	To fill system: a) For transducers without integral flush device (Snap-Tab), fill system per hospital policy. b) For transducers with integral flush device, flow is provided by pulling on the Snap-Tab, and discontinued by releasing the Snap-Tab.
10	For kits with IV sets attached, open the transducer vent port by turning the stopcock handle towards pressure tubing. Allow gravity to deliver flush solution first through the transducer and out through the vent port, then through the remaining pressure tubing by turning the appropriate stopcocks. Remove all air bubbles. CAUTION: Significant distortion of the pressure waveform or air emboli can result from air bubbles in the setup.
11	Replace all vented caps on sideports of the stopcocks with non-vented caps.
12	Mount the transducer either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder.

Step	Procedure
13	Pressurize the IV flush solution bag. Flow rate will vary with pressure across the flush device. The flow rates with the IV bag pressurized to 300 mmHg are as follows: • 3 ± 1 mL/hr (DPT with blue Snap-Tab)
14	Connect pressure tubing to the catheter per manufacturer's instructions.
15	Flush system per hospital policy. CAUTION: After each fast-flush operation, observe the drip chamber to verify that the continuous flush rate is as desired (see Complications). CAUTION: All steps should be done prior to connecting to the patient's catheter or access site.

8.0 Zeroing and Calibration

Step	Procedure
1	Adjust the level of the transducer vent port (the fluid-air interface) to correspond to the chamber where pressure is being measured. For example, in cardiac monitoring zero at level of the right atrium. This is at the phlebostatic axis, determined by the intersection of the midaxillary line and the fourth intercostal space.
2	Remove the non-vented cap and open the vent port to the atmosphere.
3	Adjust the monitor to read zero mmHg.
4	Check monitor calibration using procedure recommended by the monitor manufacturer.
5	Close the vent port to the atmosphere and replace the non-vented cap.
6	System is ready to begin monitoring pressure.

9.0 Testing Dynamic Response

The assembly may be tested for dynamic response by observing the pressure waveform on a monitor. Bedside determination of the dynamic response of the catheter, monitor, kit, and transducer system is done after the system is flushed, attached to the patient, zeroed, and calibrated. A square-wave test may be performed by pulling the Snap-Tab and releasing quickly.

Note: Poor dynamic response can be caused by air bubbles, clotting, excessive lengths of tubing, excessively compliant pressure tubing, small bore tubing, loose connections, or leaks.

10.0 Routine Maintenance

Follow hospital policies and procedures for frequency of zeroing the transducer and monitor and for replacing and maintaining pressure monitoring lines. The TruWave disposable pressure transducer is precalibrated and has a negligible drift with time (see **Specifications**).

Step	Procedure
1	Adjust zero pressure reference each time level of the patient is changed. CAUTION: When rechecking zero or verifying accuracy, ensure that the non-vented cap is removed before opening the transducer vent port to the atmosphere.
2	Periodically check fluid path for air bubbles. Ensure that connecting lines and stopcocks remain tightly fitted.
3	Periodically observe the drip chamber to verify that the continuous flush rate is as desired.

Step	Procedure
4	The Centers for Disease Control recommends replacing IV lines and disposable or reusable transducers at 96-hour intervals.

11.0 MRI Safety Information

MR Conditional

The following device was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment:

TruWave disposable pressure transducer

Non-clinical testing demonstrated that these devices are MR Conditional according to the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient field of 4000-Gauss/cm (40.0 T/m) or less
- This device and the associated cable are not intended for use inside of the bore of the MR system and should not be in contact with the patient
- This device and the associated cable may be in the MR system room but not in operation or connected to a pressure monitoring system during an MR examination
- The pressure tubing connected to the TruWave DPT is MR safe and can be placed inside the bore during an MR examination

12.0 Complications

12.1 Sepsis/Infection

Positive cultures can result from contamination of the pressure setup. Increased risks of septicemia and bacteremia have been associated with blood sampling, infusing fluids, and catheter related thrombosis.

12.2 Air Emboli

Air can enter the patient through stopcocks that are inadvertently left open, from accidental disconnection of the pressure setup, or from flushing residual air bubbles into the patient.

12.3 Overinfusion

Excessive flow rates may result from pressures greater than 300 mmHg. This may result in a potentially harmful increase in blood pressure and fluid overdose.

12.4 Abnormal Pressure Readings

Pressure readings can change quickly and dramatically because of loss of proper calibration, loose connection, or air in the system.

WARNING: Abnormal pressure readings should correlate with the patient's clinical manifestations.

Additional risks include blood loss, thrombi, emboli, hypervolemia, adverse reaction to materials, burns or electrical shock, inappropriate/unintended treatment, and/or delay in monitoring and/or treatment.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

13.0 How Supplied

TruWave disposable pressure transducers are supplied sterile in preconnected monitoring kits (either standard design or special order). This product is for single use only. Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of the packaging integrity prior to use. Do not resterilize.

14.0 Storage

Store in a cool, dry place.

15.0 Shelf Life

The shelf life is as marked on each package.

Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada

(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

17.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

STERILE EO

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

STERILE R

has been sterilized using Irradiation.

18.0 Specifications*

TruWave Disposable Pressure Transducer

Operating Pressure Range	-50 to + 300 mmHg
Operating Temperature Range	15 to 40 °C
Sensitivity	5.0μ V/V/ mmHg ± 1%
Accuracy: Combined effects of Sensitivity, Repeatability, Nonlinearity and Hysteresis	± 1 mmHg plus 1% of reading from -50 to + 50 mmHg and ± 2.5% of reading from +50 to +300 mmHg
Excitation Frequency	dc to 5,000 Hz
Excitation Impedance	350 ohms ± 10% with typical Edwards Monitor Cable attached
Excitation Impedance (DPT)	1800 ohms - 3300 ohms
Phase Shift	< 5°
Signal Impedance	300 ohms ± 5 %
Zero Offset	≤ ± 25 mmHg
Zero Thermal Drift	≤ ± 0.3 mmHg/ °C (for products sterilized by EO) ≤ ± 0.5 mmHg/ °C (for products irradiated by e-beam)
Zero Output Drift	± 1 mmHg per 8 hours after 20 seconds warm-up
Sensitivity Thermal Drift	≤ 0.4 % / °C
Defibrillator Challenge	withstands 5 repeated discharges of 360 Joules within 5 minutes delivered into a 50 ohm load
Leakage Current	< 2μ amps at 120V RMS 60 Hz or at 264V RMS 50 Hz
Overpressure Tolerance	-500 to + 5000 mmHg
Shock Resistance	withstands 3 drops from 1 meter
Light Sensitivity	< 1 mmHg at 6 volts excitation when exposed to a 3400 °K tungsten light source at 3000 foot candles (32,293 lm/m²)
Volumetric Displacement	≤ 0.03 mm³/100 mmHg for transducer without flush device
Flow rate across flush device with IV bag pressurized to 300 mmHg	
Blue Snap-Tab	3 ± 1 mL/hr
Yellow Snap-Tab	30 ± 10 mL/hr

* at 6.00VDC and 25 °C unless otherwise stated

All specifications meet or exceed the AAMI BP22 Standard for blood pressure transducers.

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
	English	Français	Deutsch	Español	Italiano
	Contains or presence of Phthalates: bis (2- ethylhexyl) phthalate (DEHP)	Traces ou présence de phthalates : phtalate de bis (2- éthylhexyle) (DEHP)	Enthält Phthalate oder kann Spuren von Phthalaten enthalten: Bis(2- ethylhexyl phthalat (DEHP)	Contiene o hay presencia de ftalatos: Di[2-etilhexil] ftalato (DEHP)	Contiene o sono presenti ftalati: di (2-etilesil) ftalato (DEHP)
	Single use	Ne pas réutiliser	Zum einmaligen Gebrauch	Un solo uso	Monouso
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CEKennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/Unión Europea	Rappresentante autorizzato nella Comunità Europea/Unione Europea
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Consultare le istruzioni per l'uso sul sito Web
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Consultare le istruzioni per l'uso
	Non-DEHP	Sans DEHP	Nicht DEHP-haltig	Sin-DEHP	Non contiene DEHP
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro

	English	Français	Deutsch	Español	Italiano
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Importer	Importateur	Importeur	Importador	Importatore
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarriersystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile
	Quantity	Quantité	Menge	Cantidad	Quantità
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil, manipular con cuidado	Fragile, maneggiare con cura
	This way up	Flèches indiquant le haut du carton	Oben	Este lado arriba	Alto
Note: Not all symbols may be included in the labeling of this product. ■ Remarque : il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ Hinweis: Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ Nota: Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ Nota: alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto.					

Lijst metsymbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl) ftalaat (DEHP)	Indeholder ftalater: bis (2-ethylhexyl) ftalat (DEHP)	Innehåller eller har spår av ftalater: bis(2-ethylhexyl) ftalat (DEHP)	Περιέχει φθαλκικές ενώσεις ή ίχνη αυτών: φθαλκικός δι (2-αιθυλοεξυλο) εστέρας (DEHP)	Contém ou estão presentes ftalatos: ftalato de bis (2-etil-hexilo) (DEHP)
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Mantener seco
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CEmarkering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμαν-ση CE)	Conformité Européenne (Marcação CE)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas
	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen	Ακολουθείτε τις οδη- γίες χρήσης που ανα- φέρονται στην τοπο- θεσία web	Seguir as instruções de utilização no site
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις οδηγίες χρήσης	Seguir as instruções de utilização
	Niet-DEHP	Ikke-DEHP	DEHP-fri	Δεν περιέχει DEHP	Sem DEHP

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Vervaldatum	Dato for sidste anvendelse	Sista förbrukningsdag	Ημερομηνία λήξης	Data de vencimento
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Ikke-pyrogen	Μη πυρετογόνο	Não pirogénico
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Utilização condicionada em ambiente de RM
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barriersystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Breekbaar, voorzichtigheid geboden	Skrøbelig, skal håndteres forsigtigt	Ömtåligt, hanteras försiktigt	Εύθραστο: χειριστείτε με προσοχή	Frágil, manuseie com cuidado
	Deze kant boven	Denne side opad	Denna sida upp	Αυτή η πλευρά προς τα πάνω	Este lado para cima
Opmerking: Het label van dit product bevat misschien niet alle symbolen. ■ Bemærk: Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ Obs! Alla symboler kanske inte används för märkning av denna produkt. ■ Σημείωση: Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ Nota: poderão não estar incluídos todos os símbolos na rotulagem deste produto.					

Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring

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
	Suomi	Български	Română	Eesti	Lietuvių
	Sisältää tai saattaa sisältää ftalaatteja: bis(2-etyyliheksyyli)ftalaatti (DEHP)	Съдържание или наличие на фталати: бис(2-етилхексил) фталат (DEHP)	Conține sau prezintă ftalați: di (2-etilhexil) ftalat (DEHP)	Sisaldab ftalaate: bis(2-etiülheksüül)ftalaat (DEHP)	Sudėtyje yra ftalatų: bis-(2-etilheksil) ftalatas (DEHP)
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
	Pidä kuivana	Да се пази сухо	A se păstra uscat	Hoida kuivalt	Laikykite sausoje vietoje
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CEmerkintä)	Conformité Européenne (CE мар-кировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhiseiga	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	espectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovaukītės naudojimo instrukcijomis, pateiktomis svetainėje
	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis
	Ei sisällä DEHP:tä	Не съдържа DEHP	Nu conține DEHP	Mitte-DEHP	Ne DEHP

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	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Vērojama ftalātu klātbūtne vai satur ftalātus: bis(2-etilheksil) ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE işareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na webstranici
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Ievērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	DEHP nesaturošs	DEHP değildir	Не содержит диэтилгексилфталата (ДЭГФ)	Bez DEHP	Ne sadrži DEHP

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilenoksidom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Daudzums	Miktar	Количество	Količina	Količina
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Trausls, rīkotos piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore
Piezīme. Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ Not: Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ Примечание. На этикетках данного изделия могут указываться не все символы. ■ Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ Napomena: na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.					

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Made in Dominican Republic
CCT Critical Care Technologies S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

08/24
10014979003 B / DOC-0102465 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH