Edwards

DIRECTORY

English	1	Português.....	28	Eesti.....	55
Français.....	4	Česky.....	31	Lietuvių.....	58
Deutsch.....	7	Magyar	34	Latviešu	61
Español.....	10	Polski.....	37	Türkçe.....	64
Italiano	13	Slovensky	40	Русский	67
Nederlands.....	16	Norsk.....	43	Srpski	71
Dansk	19	Suomi	46	简体中文.....	74
Svenska.....	22	Български	49	한국어	77
Ελληνικά.....	25	Română	52	日本語.....	80

English

HemoSphere ClearSight Module

Carefully read these instructions for use and all contained warnings and cautions before using this product.

Description

The HemoSphere ClearSight module is a reusable device which in conjunction with the HemoSphere advanced monitor and connected compatible pressure controller, heart reference sensor, and finger cuff(s) constitutes the HemoSphere noninvasive system. The HemoSphere noninvasive system measures patient brachial artery blood pressure and provides continuous calculation of cardiac output (CO), stroke volume (SV), pulse pressure variation (PPV), and stroke volume variation (SVV). Systemic vascular resistance (SVR) can also be calculated when a central venous pressure (CVP) value is available. When connected to an Acumen IQ finger cuff additional parameters are available that provide the clinician with information on the likelihood of a hypotensive event.

Edwards, Edwards Lifesciences, the stylized E logo, Acumen, Acumen HPI, Acumen IQ, ClearSight, FloTrac, HemoSphere, HPI, Swan, Swan-Ganz, and TruWave are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.


Figure 1. HemoSphere ClearSight Module (Model HEMCSM10)

Compatibility

The HemoSphere ClearSight module is intended to be used in combination with a compatible monitor, pressure controller, heart reference sensor, and Edwards finger cuff(s) to noninvasively measure blood pressure and additional hemodynamic parameters.

Indications for Use

The HemoSphere advanced monitor when used with the HemoSphere ClearSight module, pressure controller and an Edwards finger cuff are indicated for patients over 18 years of age in which the balance between cardiac function, fluid status and vascular resistance needs continuous assessment. It may be used for monitoring hemodynamic parameters in conjunction with a perioperative goal directed therapy protocol in a hospital

environment. In addition, the noninvasive system is indicated for use in patients with co-morbidities for which hemodynamic optimization is desired and invasive measurements are difficult. The HemoSphere advanced monitor and compatible Edwards finger cuff(s) noninvasively measures blood pressure and associated hemodynamic parameters. Refer to the ClearSight finger cuff and Acumen IQ finger cuff indications for use statements for information on target patient population specific to the finger cuff being used.

The Edwards Lifesciences Acumen Hypotension Prediction Index (HPI) feature provides the clinician with physiological insight into a patient's likelihood of future hypotensive events (defined as mean arterial pressure < 65 mmHg for at least one minute in duration) and the associated hemodynamics. The Acumen HPI feature is intended for use in operating room (OR) patients receiving advanced hemodynamic monitoring. The Acumen HPI feature is considered to be additional quantitative information regarding the patient's physiological condition for reference only and no therapeutic decisions should be made based solely on the Acumen Hypotension Prediction Index (HPI) parameter.

Intended Use

The HemoSphere advanced monitoring platform is intended for use with compatible Edwards Swan-Ganz and oximetry catheters, FloTrac, Acumen IQ, and TruWave DPT sensors, and Edwards finger cuff(s).

A comprehensive list of parameters available while monitoring with the HemoSphere advanced monitor and a connected HemoSphere ClearSight module or connected HemoSphere ClearSight module and oximetry cable are listed in tables 2 and 3.

Contraindications

In some patients with extreme contraction of the smooth muscle in the arteries and arterioles in the lower arm and hand, such as may be present in patients with Raynaud's disease, blood pressure measurement can become impossible.

Warnings

Compliance to IEC 60601-1 is only maintained when the HemoSphere ClearSight module (applied part connection) is connected to a compatible monitoring platform. Connecting external equipment or configuring the system in a way not described in these instructions will not meet this standard. Failure to use the device as instructed may increase the risk of electrical shock to the patient/operator.

Do not modify, service or alter the product in any way. Servicing, alteration or modification may

affect patient/operator safety and/or product performance.

Precautions

The HemoSphere ClearSight module is electrostatic discharge (ESD) sensitive. Do not attempt to open the module housing or use if the housing has been damaged.

Do not force the module into the HemoSphere advanced monitor slot. Apply even pressure to slide and click the module into place.

Do not use if the package or module is damaged. Damage may include cracks, scratches, dents or any signs that the housing may be compromised.

Module Setup

1. Align and insert the HemoSphere ClearSight module into an large technology (L-Tech) module slot. This is located on the top left panel of the HemoSphere advanced monitor. The module will click in place when properly engaged.

Note: Once inserted, the module does not need to be removed unless maintenance or cleaning is necessary.

2. Setup and initiate the compatible monitor.
3. Connect the pressure controller to the HemoSphere ClearSight module.
4. Wrap the pressure controller band around the patient's wrist and attach the pressure controller to the band.
5. Follow detailed instructions contained in the finger cuff IFU for proper cuff size selection and placement.
6. Apply the compatible Edwards finger cuff to the patient's finger and connect cuff to the pressure controller.
7. Follow on-screen prompts for heart reference sensor (HRS) connection and proper zeroing.
8. Refer to the HemoSphere advanced monitor operator's manual, available at eifu.edwards.com, for comprehensive monitoring and troubleshooting instructions.

MRI Information

MR Unsafe

This device is MR-unsafe and poses hazards in the MRI environment. This device contains metallic components, which can experience RF-induced heating in the MRI environment.

Specifications

Physical Specifications

Weight: approximately 0.9 kg (2 lb)
Height: 13 cm (5.1 in)
Width: 14 cm (5.6 in)
Depth: 10 cm (3.9 in)

Operating conditions

Temperature: 10 to 32.5 °C
Humidity range: 20 to 85% non-condensing
Altitude (Atmospheric Pressure): 0 m/0 ft (1,013 hPa) to 3,000 m/9,483 ft (700 hPa)
IP Rating: IPX1

Table 1. Parameter Measurement Specifications

Parameter	Specification	
arterial blood pressure	display range	0 to 300 mmHg
	accuracy ¹	1% of full scale (max 3 mmHg)
cardiac output (CO)	display range	1.0 to 20.0 L/min
	accuracy	Bias $\leq \pm 0.6$ L/min or $\leq 10\%$ (whichever is greater) Precision (1 σ) $\leq \pm 20\%$ over the range of cardiac output from 2 to 20 L/min
	reproducibility ²	$\pm 6\%$
	update rate	20 seconds
¹ Accuracy tested under laboratory conditions compared to a calibrated pressure gauge		
² Coefficient of variation - measured using electronically generated data		

Table 2. HemoSphere ClearSight Module Available Parameters

Abbreviation	Definition	Sub-system technology used	Patient population	Hospital environment
CO	continuous cardiac output	HemoSphere ClearSight module	adult only	operating room and intensive care unit
CI	continuous cardiac index			
DI _{NI}	noninvasive diastolic blood pressure			
dP/dt*	systolic slope			
E _{dyn} *	dynamic arterial elastance			
MAP _{NI}	noninvasive mean arterial blood pressure			
PPV	pulse pressure variation			
PR _{NI}	noninvasive pulse rate			
SV	stroke volume			
SVI	stroke volume index			
SVR	systemic vascular resistance			
SVRI	systemic vascular resistance index			
SVV	stroke volume variation			
SYS _{NI}	noninvasive systolic blood pressure	operating room only		
HPI*	Acumen hypotension prediction index			
*HPI parameters are available when using an Acumen IQ finger cuff, heart reference sensor (HRS), and if the HPI feature is activated. Activation is only available in certain areas. Please contact your local Edwards representative for more information on enabling this advanced feature.				

Table 3. HemoSphere ClearSight Module with Oximetry Cable Available Parameters

Abbreviation	Description	Sub-system technology used	Patient population	Hospital environment
DO ₂	oxygen delivery	HemoSphere ClearSight module and HemoSphere oximetry cable	adult only	operating room and intensive care unit
DO ₂ I	oxygen delivery index			
VO ₂	oxygen consumption			
VO ₂ I	oxygen consumption index			
VO ₂ e	estimated oxygen consumption when ScvO ₂ is being monitored			
VO ₂ le	estimated oxygen consumption index when ScvO ₂ is being monitored			

Maintenance

Do not pull on the pressure controller cable when unplugging it from the HemoSphere ClearSight module. If it is necessary to remove the module from the HemoSphere advanced monitor, press the release button to unlatch and slide module out.

It is recommended to perform a yearly pressure calibration check and send the HemoSphere ClearSight module to a qualified Edwards service center for routine service and preventive maintenance every two years. The HemoSphere ClearSight module is a reusable device and can be surface cleaned using a cloth

dampened with 70% isopropyl alcohol, 2% glutaraldehyde, 10% bleach solution, or quaternary ammonium solution. After cleaning, allow the module to air dry before using.

Cautions:

Do not use any other cleaning agents or spray/pour cleaning solution directly on the module.

Do not steam, radiate, or EtO sterilize the module.

Do not immerse.

Note: The expected useful life of the HemoSphere ClearSight module is 5 years from date of manufacture, at which point the module should be replaced and returned to Edwards Lifesciences. Please contact Technical Support or your local Edwards representative for further assistance.

Storage

The HemoSphere ClearSight module can be stored in the HemoSphere advanced monitor using the environmental limits defined below. If it is necessary to remove the module, store the module in a cool, dry place in original packaging to prevent damage.

Temperature: -18 to 45 °C
Relative humidity: 20 to 90% non-condensing at 45 °C
Altitude: 0 to 19,685 feet (6,000 m)

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Prices, specifications and model availability are subject to change without notice.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

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	Български	Română	Eesti	Lietuvių
	Каталожен номер	Număr de catalog	Kataloogi number	Katalogo numeris
 eifu.edwards.com + 1 888 570 4016	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige veebisaidil toodud kasutusjuhendit	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje
	Сериен номер	Număr de serie	Seerianumber	Serijos numeris
	Производител	Producător	Tootja	Gamintojas
	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Упълномощен представител в Европейската общност	Reprezentant autorizat în Comunitatea Europeană	Volitatud esindaja Euroopa Ühenduses	Igaliotasis atstovas Europos Bendrijoje
	Небезопасно при MR	Incompatibil IRM	MR-ohhtlik	MR nesaugus
IPX1	Защитено от капки оборудване	Echipament protejat împotriva picăturilor de apă	Tilkumiskindlad seadmed	Įranga su apsauga nuo pūslių
IPX4	Предоставя защита срещу пръски вода във всички посоки според стандарт IPX4	Asigură protecție la standardul IPX4 împotriva apei împrăscate în orice direcție	Pakub kaitset igast suunast tulevate veepritsmete eest standardi IPX4 järgi	Apsauga nuo aptaškymo vandeniu bet kuria kryptimi pagal IPX4 standartą
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Свързка на приложената част от тип BF	Racord piesă aplicată tip BF	BF-tüüpi kontaktosa ühendus	BF tipo su pacientu besiliečiančios dalies jungtis
	Не разглобявайте	A nu se demonta	Mitte demonteerida	Neardyti
	Не съдържа олово	Fără plumb	Pliivaba	Be švino
	Това устройство съдържа трансмитер на нейонизиращо лъчение, който може да причини РЧ интерференция с други устройства в близост до това устройство.	Acest dispozitiv conține un emițător de radiații neionizante, care poate cauza interferențe RF cu alte aparate din apropierea acestui dispozitiv.	See seade sisaldab mitteioniseeriva kiirguse saatjat, mis võib tekitada lähedal asuvates teistes seadmetes raadiosageduslikke häiringuid.	Šiame prietaise yra nejonizuojančiosios spinduliuotės siųstuvas, kuris netoliese esantiems prietaisams gali sukelti RD trukdžius.
	Да се пази сухо	A se păstra uscat	Hoida kuivana	Laikykite sausoje vietoje
	Чупливо, да се борава внимателно	Fragil, a se manipula cu atenție	Õrn, käsitseda ettevaatlikult	Trapus, elgtis atsargiai
	Съвместимост с Директивата за ограничаване на употребата на определени опасни вещества (RoHS) – само за Китай	Conformitatea cu restricțiile de utilizare a substanțelor periculoase (RoHS) - numai China	Vastavus ohtlike ainete kasutamise piirangu (RoHS) direktiivile – ainult Hiina	Atitinka pavojingų medžiagų apribojimus (RoHS) – tik Kinijoje

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SENT27SL5.5x7.25.8

Edwards

Edwards Lifesciences Services GmbH

Edisonstrasse 6
85716 Unterschleissheim
Germany

10/19
© Copyright 2019, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC

One Edwards Way
Irvine, CA 92614 USA
Made in USA

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

WEB IFU
10031912001 A
DOC-0117289 A