
Edwards

Heart Reference Sensor • Capteur de référence cardiaque • Sensor de referencia del corazón
• Sensore di riferimento cuore • Hartreferentiesensor

DIRECTORY		
English.....	2	Español..... 8
Français	5	Italiano 11
		Nederlands 14

Heart Reference Sensor (HRS)

REF EVHRS

For figures 1 through 6 please refer to page 18.

Description

The heart reference sensor (also referred to as 'HRS') can be used with the EV1000 clinical platform NI and HemoSphere advanced monitoring platform to noninvasively measure blood pressure and key hemodynamic parameters using a compatible Edwards finger cuff(s) wrapped around the finger. When using this method it is necessary to use a heart reference sensor, which compensates the hydrostatic pressure changes due to difference in height between the finger and heart. The heart reference sensor is a non-serviceable part.

Intended Use

The heart reference sensor is intended for use with a compatible Edwards noninvasive monitoring system - composed of compatible monitor, pressure source (pump), heart reference sensor, compatible Edwards finger cuff(s) and pressure controller - for noninvasive measurement of blood pressure and associated hemodynamic parameters. Refer to the operator's manual of the compatible Edwards monitor being used for specific information on the intended use environment and patient population.

Instructions for Use

Unfold the cover of this Instruction For Use to show the figures corresponding to the instructions.

CAUTION: Do not use a damaged heart reference sensor. This may result in inaccurate measurements or may damage the EV1000 clinical platform NI or HemoSphere advanced monitoring platform.

Connect the Heart Reference Sensor to the Pressure Controller

1. Connect the heart reference sensor to the pressure controller (Fig 2).

Zero the Heart Reference Sensor

Prior to every measurement using the heart reference sensor, it needs to be zeroed.

2. Navigate to the *Zero & Waveform* screen on the EV1000 or HemoSphere monitor. Refer to monitor operator's manual for instructions.
3. Vertically align the two ends of the heart reference sensor (Fig 3) and touch the Zero button - **0** - (Fig 4).
NOTE: Keep both ends of the heart reference sensor vertically aligned until the Zeroing procedure has ended.

Attach the Heart Reference Sensor to the Patient

4. Attach the finger end of the heart reference sensor to a single compatible Edwards finger cuff (Fig 5).
5. Attach the heart end of the heart reference sensor to the patient at the phlebostatic axis (Fig 6), level by using an HRS clip.
NOTE: Changing body position during the measurement may require reapplication of the heart end.

Edwards, Edwards Lifesciences, the stylized E logo, EV1000, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

CAUTION: Improper zeroing of the heart reference sensor may adversely affect the accuracy of the readings.

CAUTION: Make sure that the heart reference sensor is correctly applied so that it can be leveled to the phlebostatic axis.

Cleaning Instructions

The heart reference sensor can be cleaned using the following disinfectants:

- 70% isopropyl alcohol solution
 - 10% sodium hypochlorite water solution
1. Moisten a clean cloth with disinfectant and wipe the surfaces.
 2. Dry the surface with a clean, dry cloth.

WARNING: Do not disinfect the heart reference sensor by autoclave or gas sterilization.

CAUTION: Do not immerse the heart reference sensor.

CAUTION: Do not immerse any cable connectors in fluid.

CAUTION: Conduct periodic inspections of all cables for defects. Do not coil cables tightly when storing.

CAUTION: Clean and store the heart reference sensor after each use.

Preventive Maintenance

The finger component of the heart reference sensor may be damaged if subjected to moderate to significant surface impact. Although the likelihood of damage is small, the resulting displayed values would be biased by the difference in height from the heart to the finger cuff. Even though this damage cannot be seen by looking at the heart reference sensor, it is possible to confirm whether the damage has occurred by following the below procedure prior to each use:

1. Connect the heart reference sensor to the port on compatible Edwards noninvasive monitoring system and go to the zeroing screen.
2. As instructed in the Instructions for Use section, bring the two ends of the heart reference sensor level with each other.
3. Observe the value shown on the zeroing screen below the zeroing button.
4. Raise one end of the heart reference sensor 6 inches (15 cm) above the other end.
5. Observe that the value shown has changed by at least 5 mmHg.
6. Reverse the ends such that the other end is now 6 inches (15 cm) above the first end.
7. Observe the value shown changed in the opposite direction by at least 5 mmHg from the original value.

If the value does not change as described, then the heart reference sensor may have been damaged. Contact your local Technical Support office as indicated in the Technical Assistance section. A replacement unit shall be provided. If the value does change, the heart reference sensor is functioning normally and can be used for hemodynamic monitoring.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837

Outside the U.S. and Canada (24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 Option 4

This document is subject to change without prior notice. Refer to the latest version of the EV1000 clinical platform NI or HemoSphere advanced monitoring platform operator's manual for additional information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Symbol Legend

Symbols	Description
	Catalogue Number
	Serial number
	Quantity
	Use by
	Conformité Européenne (CE Mark)
	Authorized representative in the European Community
	Manufacturer
	Follow instructions for use
	Follow instructions for use on the website
	Separate collection for electrical and electronic equipment in accordance with EC Directive 2012/19/EU
	Caution
Rx only	Caution: Federal (USA) law restricts this device to sale by, or on the order of a physician
	MR Unsafe
	Keep dry
	Fragile, handle with care
	Keep away from direct sunlight
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Medical Device
Note: Not all symbols may be included in the labeling of this product.	

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© Copyright 2019, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC