Edwards

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	7
Español (es).....	10
Italiano (it).....	13
Nederlands (nl).....	16
Dansk (da).....	19
Svenska (sv).....	22
Ελληνικά (el).....	25
Português (pt).....	28
Česky (cs).....	31
Magyar (hu).....	34
Polski (pl).....	37
Slovensky (sk).....	40
Norsk (no).....	43
Suomi (fi).....	46
Български (bg).....	49
Română (ro).....	52
Eesti (et).....	55
Lietuvių (lt).....	58
Latviešu (lv).....	61
Türkçe (tr).....	64
Русский (ru).....	67
Srpski (sr).....	70
Hrvatski (hr).....	73
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike	76
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπομνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sūmbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	79

English

Acumen IQ Fluid Meter

Edwards, Edwards Lifesciences, the stylized E logo, Acumen, Acumen AFM, Acumen IQ, AFM, HemoSphere, TruWave, and VAMP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device. Refer to the HemoSphere advanced monitor operator's manual, available at eifu.edwards.com, for comprehensive monitoring procedures.

For Single Use Only

1.0 Description

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use.

The Acumen IQ fluid meter (reference Figure 1) is a sterile, single use device that measures the flow of fluid delivered to a patient through the intravenous line to which it is connected. When connected to a compatible monitor via the reusable Acumen AFM cable (reference Figure 2), the Acumen IQ fluid meter automatically tracks the fluid being administered to the patient and provides that data to the Assisted Fluid Management (AFM) software feature.

The Acumen IQ fluid meter only operates when used with an Acumen IQ sensor and is part of the AFM software feature. The AFM software feature is intended to work with a single Acumen IQ fluid meter at a time. Refer to the Acumen IQ sensor instructions for use (IFU) for more details on the sensor; refer to the HemoSphere advanced monitor for details on the AFM software feature.

2.0 Intended Use/Purpose

The Acumen IQ fluid meter is a sterile single use device that is intended to be used with the Acumen AFM cable and AFM software feature to inform the user of the rate of flow. The device is intended to be used by qualified personnel or clinicians in a clinical setting for up to 24 hours.

3.0 Indications

The Acumen IQ fluid meter is indicated for surgical patients over 18 years of age to track the fluid being administered to the patient, when used with a compatible hemodynamic monitoring platform.

4.0 Contraindications

There are no contraindications for using the Acumen IQ fluid meter.

5.0 Warnings

- The Acumen IQ fluid meter, when used with Acumen AFM cable and Acumen AFM software feature, should not be used exclusively to treat the patient. A review of the patient's hemodynamics is recommended prior to initiating treatment.
- The Acumen IQ fluid meter is MR Conditional. For a patient undergoing an MRI examination, please refer to the MRI Safety Information section for specific conditions to ensure patient safety.
- Do not modify, service or alter the product in any way. Servicing, alteration or modification may affect patient/operator safety and/or product performance.
- This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.
- If all air is not removed from the intravenous flush solution bag, air may be forced into the patient's vascular system when the solution is exhausted.
- If there are air bubbles in the intravenous line, there might be an impact to the fluid meter accuracy.
- Do not allow air bubbles to enter the setup as air bubbles may lead to air emboli.
- The Acumen IQ fluid meter should not be used for more than 24 hours. Use of the Acumen IQ fluid meter for more than 24 hours may compromise accuracy of the fluid meter measurement.
- This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0. Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- The Acumen IQ fluid meter accuracy is affected by the fluid type and the rate of fluid delivery. The use of the Acumen IQ fluid meter could result in delivery of more than the intended amount per the performance specification provided in Table 1.

6.0 Precautions

The Acumen IQ fluid meter is not compatible with blood products. The Acumen IQ fluid

meter is only compatible with the following IV solutions:

- Sodium Chloride Injection 0.9% (NaCl 0.9%)
- Ringer's lactate solution (RL), also known as sodium lactate solution and Hartmann's solution
- PlasmaLyte
- Dextran 40
- Albumin 5%
- Hetastarch 6%

Use of the Acumen IQ fluid meter with blood products or any other solution or medication may compromise accuracy of the fluid meter measurement.

7.0 Complications

There are some risks related to fluid management, which include, but are not limited to, the following:

- Pulmonary edema
- Pleural effusion
- Bladder rupture
- Myocardial Injury after Noncardiac Surgery (MINS)
- Renal-replacement therapy
- Acute kidney injury
- Anastomic leakage (for bowel surgery)
- Stoma necrosis
- Surgical site infection

Additional risks associated with the use of the device when used with a connected compatible monitor:

- Hypervolemia
- Hypovolemia
- Sepsis/Infection
- Air emboli
- Patient or clinician burns or electrical shock
- Tissue damage

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

8.0 Compatible Devices

The Acumen IQ fluid meter is intended to be used in conjunction with a compatible Edwards HemoSphere advanced monitor, the Acumen AFM cable, and Edwards approved sensors/ closed blood sampling system as part of the Assisted Fluid Management (AFM) solution.

Note: Refer to the HemoSphere advanced monitor operator's manual, available at eifu.edwards.com, for comprehensive monitoring procedures and information on use of the Acumen AFM cable and Acumen AFM software feature.

9.0 Equipment

The Edwards Acumen IQ fluid meter can be provided as part of a kit.

9.1 Kit Components

The kits may be comprised of the following components or accessories:

- Tubing
- Stopcocks
- Caps
- Acumen IQ sensor
- TruWave disposable pressure transducer
- VAMP closed blood sampling system

Note: Refer to the Acumen IQ sensor, TruWave disposable pressure transducer and VAMP closed blood sampling system Instructions for Use for instructions on how to use these devices.

9.2 Additional Components Required

- Intravenous (IV) set
- Venous catheter

- IV solution

9.3 Compatible IV Solutions

- Sodium Chloride Injection 0.9% (NaCl 0.9%)
- Ringer's lactate solution (RL), also known as sodium lactate solution and Hartmann's solution
- PlasmaLyte
- Dextran 40
- Albumin 5%
- Hetastarch 6%

10.0 Procedure

These are general instructions for setting up an Edwards approved sensor and/or Edwards compatible hardware with the Acumen IQ fluid meter. Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

Step	Procedure
1	Remove the Acumen IQ fluid meter from the sterile packaging.
2	Place Acumen IQ fluid meter in line with the intravenous set ensuring connections as illustrated in Figure 3 on page 77. Ensure connections are secure but not overtightened. Precaution: The Acumen IQ fluid meter is not compatible with blood products. The Acumen IQ fluid meter is only compatible with the following IV solutions: • Sodium Chloride Injection 0.9% (NaCl 0.9%) • Ringer's lactate solution (RL), also known as sodium lactate solution and Hartmann's solution • PlasmaLyte • Dextran 40 • Albumin 5% • Hetastarch 6% Use of the Acumen IQ fluid meter with blood products or any other solution or medication may compromise accuracy of the fluid meter measurement.
3	Remove all air from the intravenous flush solution bag. WARNING: If all air is not removed from the intravenous flush solution bag, air may be forced into the patient's vascular system when the solution is exhausted. WARNING: If there are air bubbles in the intravenous line, there might be an impact to the fluid meter accuracy.
4	Close the roller clamp on the intravenous set and connect the intravenous set to the intravenous flush bag. Hang the bag approximately 2 feet (60 cm) above the patient. This height will provide approximately 45 mmHg of pressure to prime the setup.
5	Fill the drip chamber halfway with flush solution by squeezing the drip chamber. Open roller clamp.
6	Prime system using gravity only to decrease fluid turbulence and mitigation of bubbles.
7	Mount the Acumen IQ fluid meter on an IV pole in a vertical position using the appropriate clamp and holder. Slide into place in holder.
8	Connect the extension tubing to venous catheter per manufacturer instructions.

Step	Procedure
9	Flush system per hospital policy. Precaution: If flushing system with a non-compatible IV solution, flush from a port that is distal to the fluid meter. Flushing with a non-compatible IV solution through the fluid meter may compromise accuracy of the fluid meter measurement.
10	Connect the Acumen AFM cable to a compatible HemoSphere advanced monitor at the end indicated by ① in Figure 2 on page 76. Note: Refer to the HemoSphere advanced monitor operator's manual, available at eifu.edwards.com, for comprehensive monitoring procedures and information on use of the Acumen AFM software feature.
11	Setup and initiate the compatible HemoSphere advanced monitor.
12	Connect the Acumen IQ fluid meter to the end of the Acumen AFM cable indicated by ② in Figure 2 on page 76.
13	When unplugging the Acumen IQ fluid meter from the Acumen AFM cable, always pull at the connection site. Do not pull from cable or use tools to disconnect.

11.0 MRI Safety Information

MR Conditional

WARNING: The Acumen IQ fluid meter shall not be connected to the Acumen AFM cable and HemoSphere Advanced Monitor during an MRI examination. Failure to follow this guideline may result in a serious patient injury.

The following device was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International Designation: F2503, Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

11.1 Acumen IQ fluid meter

Non-clinical testing demonstrated that the Acumen IQ fluid meter is MR Conditional under the following conditions:

- Static magnetic field of 3.0 T or less.
- Maximum spatial gradient magnetic field of 3,000 gauss/cm (30 T/m).
- The Acumen IQ fluid meter is not intended for use inside the bore of the MR system and should not be in direct contact with the patient.
- The device may be in the MR system room but not in operation or connected to a monitoring system during an MRI examination.
- All parts of the Acumen IQ fluid meter should be at least 4 cm from the imaging region.

Metal in the Acumen IQ fluid meter is non ferro-magnetic. Using ASTM F2052 and ASTM F2213, there was minimal magnetic force and torque in a 3T MR system.

Image artifact was tested using sequences in ASTM F2119. Artifact was discernible at distances of up to 4 cm from the Acumen IQ fluid meter.

Precaution: Follow the conditions for safe scanning for any accessory devices that are connected to the Acumen IQ fluid meter. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

12.0 Maintenance

Periodically check fluid path for air bubbles. Ensure that connecting lines and stopcocks remain tightly fitted.

WARNING: Do not allow air bubbles to enter the setup as air bubbles may lead to air emboli.

Periodically observe the drip chamber to verify that the continuous flush rate is as desired.

Follow hospital policies and procedures for replacing and maintaining IV/pressure monitoring/fluid lines.

WARNING: The Acumen IQ fluid meter should not be used for more than 24 hours. Use of the Acumen IQ fluid meter for more than 24 hours may compromise accuracy of the fluid meter measurement.

Do not steam, radiate, or EtO sterilize the device. Do not immerse.

13.0 How Supplied

The Acumen IQ fluid meter is supplied sterile and the fluid path is nonpyrogenic if package is

undamaged or unopened. Visually inspect for breaches of packaging integrity prior to use. Do not use if package is opened or damaged. Do not re-sterilize.

14.0 Storage

Store in a cool, dry place.

15.0 Specifications

15.1 Operating Conditions/Use Environment

Temperature Range: 10 to 37 °C

Humidity range: 20 - 80% non-condensing

Altitude (Atmospheric Pressure):
0 m/0 ft (1013 hPa) to 3048 m/10,000 ft (697 hPa)

IP Rating: IPX4

16.0 Shelf Life

The recommended shelf life for the Acumen IQ fluid meter is 24 months from the date of manufacture.

17.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

18.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide

Table 1: Acumen IQ fluid meter Performance Specification

Liquid	Flow Rate (mL/min)		
	Low Flow (0-50)	Medium Flow (50-100)	High Flow (100-150)
Plasmalyte	-1.5 ± 1.4%	-0.6 ± 2.2%	-1.4 ± 1.5%
NaCl 0.9%	1.6 ± 1.8%	2.1 ± 1.6%	1.3 ± 2.4%
Lactated Ringers	-0.0 ± 3.1%	1.0 ± 2.7%	0.2 ± 2.8%
Dextran 40	-6.2 ± 1.8%	-5.1 ± 5.1%	-7.8 ± 2.5%
HESPAN (Hetastarch 6%)	-2.7 ± 2.1%	-3.2 ± 1.8%	-3.4 ± 2.1%
Albumin 5%	-0.9 ± 1.6%	-1.8 ± 2.0%	-2.1 ± 2.8%

*Accuracy tested under laboratory conditions

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

	English	Français	Deutsch	Español	Italiano
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Consultare le istruzioni per l'uso sul sito Web
	Consult instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	MR Conditional	IRM sous conditions	Bedingt MR-sicher	Condicional con respecto a RM	A compatibilità RM condizionata
	Provides protection against water splashing in any direction to IPX4 standard	Assure une protection contre les projections d'eau multidirectionnelles conformément à la norme IPX4	Schutz gegen Spritzwasser aus allen Richtungen gemäß IPX4	Proporciona protección contra las salpicaduras de agua en cualquier dirección hasta el estándar IPX4	Fornisce protezione contro gli spruzzi d'acqua in qualsiasi direzione secondo lo standard IPX4
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	English	Français	Deutsch	Español	Italiano
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Quantity	Quantité	Menge	Cantidad	Quantità
	Importer	Importateur	Importeur	Importador	Importatore
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarrieresystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile
	Contains hazardous substances	Contient des substances dangereuses	Enthält gefährliche Substanzen	Contiene sustancias peligrosas	Contiene sostanze pericolose
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número de modelo
	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen	Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Siga as instruções de uso no site
	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Se bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης	Consulte as instruções de uso
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabricação
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Condicionado em ambiente de RM
	Biedt bescherming tegen rondspetterend water op IPX4-standaard	Yder beskyttelse mod vand, der sprøjter i alle retninger til IPX4-standard	Skyddar mot stänkande vatten från samtliga riktningar enligt standarden IPX4	Παρέχει προστασία από την εκτίναξη νερού προς κάθε κατεύθυνση σύμφωνα με το πρότυπο IPX4	Oferece proteção contra respingos de água em qualquer direção, de acordo com a norma IPX4
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Aviso: a lei federal (nos EUA) restringe este aparelho à venda por um médico ou sob sua solicitação.
	Niet hergebruiken	Må ikke genanvendes	Får inte återanvändas	Μην επαναχρησιμοποιείτε	Não reutilizar
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado utilizando óxido de etileno
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não use se a embalagem estiver danificada e consulte as instruções de uso
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marca CE)
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Vervaldatum	Anvendes inden	Sista förbrukningsdag	Ημερομηνία λήξης	Data de vencimento
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Niet-pyrogeen	Ikke-pyrogen	Ikke-pyrogen	Μη πυρετογόνο	Não pirogênico
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não reesterilizar

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrieresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Bevat gevaarlijke stoffen	Indeholder farlige stoffer	Innehåller farliga ämnen	Περιέχει επικίνδυνες ουσίες	Contém substâncias perigosas
	Unieke instrumentidentificatie	Unik udstysidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único do dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Armazenar em um local fresco e seco

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
	Suomi	Български	Română	Eesti	Lietuvių
	Mallinúmero	Номер на модела	Număr model	Mudeli number	Modelio numeris
	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovaukites naudojimo instrukcijomis, pateiktomis svetainėje
	Katso käyttöohjeet	Направете справка с инструкциите за употреба	Consultați instrucțiunile de utilizare	Palun lugege kasutusjuhiseid!	Žr. naudojimo instrukcijas
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Įgaliojasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Ehdollisesti turvallinen magneettikuvauksessa	Безопасно при MR при определени условия	Condiționat RM	Ohutu magnetresonants-tomograafias teatud tingimustel	Sąlyginis MR
	Kotelointiluokka IPX4 (suojaus roiskavalta vedeltä)	Предоставя защита срещу пръски вода във всички посоки според стандарт IPX4	Asigură protecție la standardul IPX4 împotriva apei împrăscate în orice direcție	Kaitseb mis tahes suunas langevate veepritsmete eest standardi IPX4 järgi.	Apsauga nuo aptaškymo vandeniu bet kuria kryptimi pagal IPX4 standartą
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Ei uudelleenkäytettävä	Да не се използва повторно	A nu se reutiliza	Mitte korduskasutada	Nenaudoti pakartotinai
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Naudoti iki nurodytos datos

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Sisältää vaarallisia aineita	Съдържа опасни вещества	Conține substanțe periculoase	Sisaldab ohtlikke aineid	Sudėtyje yra pavojingų medžiagų
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusiis priemonės identifikatorius
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	levērojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Skatīt lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu	Pogledajte upute za upotrebu
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašćeni predstavnik u Europskoj zajednici / Europskoj uniji

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Nodrošina aizsardzību pret ūdens šļakstīšanos jebkurā virzienā atbilstoši IPX4 standartam.	IPX4 standardına uygun olarak herhangi bir yönden sıçrayan suya karşı koruma sağlar	Обеспечивает защиту от брызг воды (в любом направлении) по стандарту IPX4	Obezbeđuje zaštitu od prskanja vode u bilo kom pravcu prema IPX4 standardu	Osigurava zaštitu od prskanja vode u bilo kojem smjeru u skladu sa standardom IPX4
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Oprez: Savezni zakoni (SAD) zakon ograničava prodaju ovog uređaja na prodaju od strane ili po nalogu lekara.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Daudzums	Miktar	Количество	Količina	Količina
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda dos símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Satur bīstamas vielas	Tehlikeli maddeler içerir	Содержит вредные вещества	Sadrži opasne supstance	Sadrži opasne tvari
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Observação:** Nem todos os símbolos podem ser incluídos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiety niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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Made in Dominican Republic
Edwards Lifesciences AG
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican
Republic

11/23
10035604001 A / DOC-0131376 A
© Copyright 2023,
Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH