DIRECTORY

English (en)	1
Français (fr)	3
Deutsch (de)	4
Español (es)	7
Italiano (it)	8
Nederlands (nl)	10
Dansk (da)	12
Svenska (sv)	14
Ελληνικά (el)	16
Português (pt)	18
Česky (cs)	20
Magyar (hu)	21

Polski (pl)	23
Slovensky (sk)	25
Norsk (no)	27
Suomi (fi)	29
Български (bg)	30
Română (ro)	32
Eesti (et)	34
Lietuvių (lt)	36
Latviešu (lv)	37
Türkçe (tr)	39
Русский (ru)	41
Srpski (sr)	43

Hrvatski (hr)	45
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola	47-55

English

Pressure Monitoring Kit

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

These are general instructions for setting up a typical pressure monitoring system. Since kit configurations and setup procedures vary according to hospital preferences, it is the responsibility of the hospital staff to determine exact policies and procedures.

1.0 Description

The Edwards Lifesciences Pressure Monitoring Kit is a sterile, single-use kit that monitors pressure when attached to a pressure monitoring catheter and a compatible pressure transducer.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The known benefits of pressure monitoring include providing information regarding the hemodynamic status of the patient which may lead to improved data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation.

This device is intended to be used by medical professionals who have been trained in the safe use

of hemodynamic technologies and clinical usage of arterial line pressure monitoring as part of their respective institutional guidelines.

2.0 Intended Use/Purpose

The Pressure Monitoring Set is intended to monitor pressure when attached to a pressure monitoring catheter and compatible pressure monitoring instrument.

3.0 Indications

The Pressure Monitoring Kit is for use on patients requiring intravascular, intracranial, or intrauterine pressure monitoring.

4.0 Contraindications

There are no absolute contraindications for using the Pressure Monitoring Kit in patients requiring invasive pressure monitoring. Flush devices should not be used when monitoring intracranial pressures.

5.0 Warnings

- Do not use the flush device during intracranial pressure monitoring.
- For severely fluid-restricted patients such as children and neonates, use an infusion pump in series with the flush device to accurately regulate the minimum amount of flush solution needed to maintain catheter patency, while allowing continuous pressure monitoring. High pressures, which may be generated by an infusion pump at certain flow rates, may override the flush device restriction, resulting in fast flushing at the rate set by the pump.
- Do not allow air bubbles to enter the pressure monitoring system, especially when monitoring atrial pressures.
- If the flush system is not adequately pressurized relative to the patient's blood pressure, blood bleed-back and catheter clotting may occur.

- Avoid contact with any topical cream or ointment that attacks polymeric materials.
- Do not expose electrical connections to fluid contact.
- Do not autoclave the reusable cable.
- This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.
- Some models may contain phthalates, specifically DEHP [Bis (2 ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

6.0 Procedure

Step	Procedure
1	Using aseptic technique, remove the pressure monitoring kit from the sterile package.
2	Add additional components as needed to complete the monitoring system using procedures recommended by your hospital staff.
3	Ensure that all connections are secure. Note: Wet connections promote overtightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.
4	Remove all air from the IV flush solution bag, which is heparinized per hospital policy. CAUTION: If all air is not removed from the bag, air may be forced into the patient's vascular system when the solution is exhausted.

Step	Procedure
5	Close the roller clamp on the IV set and connect it to the IV flush bag. Hang the bag approximately 2 feet (60 cm) above the patient. This height will provide approximately 45 mmHg of pressure to prime the setup.
6	Fill the drip chamber halfway, or per hospital policy, with flush solution by squeezing the drip chamber. Open the roller clamp.
7	To fill system: a) For transducers without integral flush device (Snap-Tab device), fill system per hospital policy. b) For transducers with integral flush device, flow is provided by pulling on the Snap-Tab device and discontinued by releasing the Snap-Tab device.
8	Allow gravity to deliver flush solution, first through the transducer then through the remaining kit by turning the appropriate stopcocks. Remove all bubbles. CAUTION: Significant distortion of the pressure waveform or air emboli can result from air bubbles in the setup.
9	Replace vented caps on sideports of stopcocks with non-vented caps. CAUTION: Overtightening may damage the transducer.
10	Pressurize the IV flush solution bag. Flow rate will vary with pressure across the flush device, 3 ± 1 mL/hr or 30 ± 10 mL/hr at 300 mmHg. CAUTION: Avoid pressurization above 300 mmHg. Overpressure relief (increased fluid flow rate within the fluid path) may occur above 500 mmHg.
11	Connect pressure tubing to the catheter.
12	Open stopcock to the patient and flush system to fill the catheter. CAUTION: After each fast-flush operation, observe the drip chamber to verify that the continuous flush rate is as desired (see Complications).

7.0 Zeroing and Calibration

Step	Procedure
1	Adjust level of the transducer vent port (the fluid-air interface) to correspond to the chamber where pressure is being measured. (Example: In cardiac monitoring, zero at level of right atrium. This is at the phlebostatic axis, determined by the intersection of the midaxillary line and the fourth intercostal space.)
2	Remove the non-vented cap and open the vent port to the atmosphere.
3	Adjust the monitor to read zero mmHg.

Step	Procedure
4	Check monitor calibration using procedure recommended by the manufacturer.
5	Check transducer performance by applying a known amount of pressure to the transducer through the vent port using either a water or mercury manometer. Make sure calibration is within operating range. For pulmonary artery monitoring, calibrate at 20 mmHg. For arterial pressure monitoring, calibrate at 100 mmHg.
6	Close the vent port to the atmosphere and replace the non-vented cap.
7	Open the transducer to the catheter and begin monitoring pressure.

8.0 Routine Maintenance

Follow hospital policies and procedures for zeroing and calibrating transducer, monitor, and when replacing pressure monitoring lines.

Step	Procedure
1	Adjust zero pressure reference each time level of the patient is changed. CAUTION: When rechecking zero and calibration, turn off the stopcock to the patient before opening the transducer vent port to the atmosphere.
2	Periodically check lines and transducer domes for air bubbles. Ensure that connecting lines and stopcocks remain tightly fitted.
3	Periodically observe the drip chamber to verify that the continuous flush rate is as desired.

9.0 Complications

9.1 Sepsis/Infection

Positive cultures can result from contamination of the pressure setup. Increased risks of septicemia and bacteremia have been associated with blood sampling, infusing fluids, and catheter related thrombosis.

9.2 Over-infusion

Excessive flow rates may result from pressures greater than 300 mmHg. This may result in a potentially harmful increase in blood pressure and fluid overdose.

9.3 Air Emboli

Air can enter the patient through stopcocks that are inadvertently left open, from accidental disconnection of the pressure setup, or from flushing residual air bubbles into the patient.

9.4 Abnormal Pressure Readings

Pressure readings can change quickly and dramatically because of loss of proper calibration, loose connections, or air in the system.

WARNING: Abnormal pressure readings should correlate with the patient's clinical manifestations. Verify transducer function with a known amount of pressure before instituting therapy.

Additional risks include blood loss, thrombi, emboli, hypervolemia, adverse reaction to materials, burns or electrical shock, inappropriate/unintended treatment, and/or delay in monitoring and/or treatment.

10.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use. Do not resterilize.

This product is for single use only.

11.0 Storage

Store in a cool, dry place.

12.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

13.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

14.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Because of continuing product improvements, prices, specifications, and model availability are subject to change without notice.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

STERILE EO

**has been sterilized using Ethylene Oxide.
Alternatively, product bearing the symbol:**

STERILE R

has been sterilized using Irradiation.

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
	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl)ftalaat (DEHP)	Indeholder ftalater: bis (2-ethylhexyl) ftalat (DEHP)	Innehåller eller har spår av ftalater: bis(2-etylhexyl) ftalat (DEHP)	Περιέχει φθαλικές ενώσεις ή ίχνη αυτών: φθαλικός δι (2-αιθυλοεξυλο) εστέρας (DEHP)	Contém ou estão presentes ftalatos: ftalato de bis (2-etil-hexilo) (DEHP)
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό	Manter seco
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Conformité Européenne (CE-markering) van TÜV SÜD Product Service GmbH (aangemelde instantie)	Conformité Européenne (CE-mærke) fra TÜV SÜD Product Service GmbH (bemyndiget organ)	Conformité Européenne (CE-märke) för TÜV SÜD Product Service GmbH (anmält organ)	Conformité Européenne (σήμανση CE) της TÜV SÜD Product Service GmbH (κοινοποιημένος οργανισμός)	Conformité Européenne (Marcação CE) da TÜV SÜD Product Service GmbH (organismo notificado)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing op de website eifu.edwards.com	Se brugsanvisningen på hjemmesiden eifu.edwards.com	Se bruksanvisningen på webbplatsen eifu.edwards.com	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web eifu.edwards.com	Consultar as instruções de utilização no website eifu.edwards.com
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrónicas
	Volg de gebruiksaanwijzing op de website eifu.edwards.com	Følg brugsanvisningen på hjemmesiden eifu.edwards.com	Följ bruksanvisningen på webbplatsen eifu.edwards.com	Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web eifu.edwards.com	Seguir as instruções de utilização no website eifu.edwards.com
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις οδηγίες χρήσης	Seguir as instruções de utilização
	Niet-DEHP	Ikke-DEHP	DEHP-fri	Δεν περιέχει DEHP	Sem DEHP
	Gesteriliseerd met behulp van ethylenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Vervaldatum	Dato for sidste anvendelse	Sista förbrukningsdag	Ημερομηνία λήξης	Data de vencimento
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogénico
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Utilização condicionada em ambiente de RM

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
Opmerking: Het label van dit product bevat misschien niet alle symbolen. ■ Bemærk: Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ Obs! Alla symboler kanske inte används för märkning av denna produkt. ■ Σημείωση: Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ Nota: poderão não estar incluídos todos os símbolos na rotulagem deste produto.					

Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring

	Česky	Magyar	Polski	Slovensky	Norsk
	Obsahuje ftaláty: di-(2-ethylhexyl) ftalát (DEHP)	Ftalátokat tartalmaz, illetve ftalátok lehetnek jelen: di(2-etilhexil)-ftalát (DEHP)	Zawartość lub obecność ftalanów: ftalan bis(2-etyloheksylu) (DEHP)	Obsah alebo prítomnosť ftalátov: di(2-etylhexyl) ftalát (DEHP)	Inneholder eller kan inneholde spor av ftalater: bis(2-etylheksyl)ftalat (DEHP)
	K jednorázovému použití	Egyszer használatos	Do jednorazowego użytku	Na jednorazové použitie	Til engangsbruk
	Chraňte před vlhkem	Tartsa szárazon	Chronić przed wilgocią	Uchovávať v suchu	Oppbevares tørt
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (značka CE) TÜV SÜD Product Service GmbH (notifikovaný subjekt)	A TÜV SÜD Product Service GmbH (bejelentett szervezet) Conformité Européenne (CE-jelölés) jelölése	Conformité Européenne (oznakowanie CE) nadane przez firmę TÜV SÜD Product Service GmbH (jednostkę notyfikowaną)	Conformité Européenne (značka CE) spoločnosti TÜV SÜD Product Service GmbH (notifikovaná osoba)	Conformité Européenne (CE-merke) fra TÜV SÜD Product Service GmbH (teknisk kontrollorgan)
	Autorizovaný zástupce v Evropském společenství / Evropské unii	Hivatalos képviselő az Európai Közösségekben/ Európai Unióban	Autoryzowany przedstawiciel we Wspólnocie Europejskiej/ Unii Europejskiej	Autorizovaný zástupca v krajinách Európskeho spoločenstva/Európskej únie	Autorisert representant i Det europeiske fellesskap / EU
	Výrobce	Gyártó	Producent	Výrobca	Produsent
	Datum výroby	Gyártás ideje	Data produkcji	Dátum výroby	Produksjonsdato
	Číslo šarže	Tételszám	Numer serii	Číslo šarže	Lotnummer

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

	Suomi	Български	Română	Eesti	Lietuvių
	Sisältää tai saattaa sisältää ftalaatteja: bis(2-etyyliheksyyli)ftalaatti (DEHP)	Съдържание или наличие на фталати: бис(2-етилхексил) фталат (DEHP)	Conține sau prezintă ftalați: di (2-etilhexil) ftalat (DEHP)	Sisaldab ftalaate: bis(2-etiülheksüül) ftalaat (DEHP)	Sudėtyje yra ftalatų: bis-(2-etilheksil) ftalatas (DEHP)
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
	Pidä kuivana	Да се пази сухо	A se păstra uscat	Hoida kuivalta	Laikykite sausoje vietoje
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	TÜV SÜD Product Service GmbH:n (ilmoitettu laitos) myöntämä Conformité Européenne (CE) -merkintä	Conformité Européenne (CE) маркировка) на TÜV SÜD Product Service GmbH (нотифициран орган)	Conformité Européenne (marcă CE) al TÜV SÜD Product Service GmbH (organism notificat)	TÜV SÜD Product Service GmbH (teavitatud asutus) Conformité Européenne (CE-märk)	TÜV SÜD Product Service GmbH (notifikuotoji įstaiga) Conformité Européenne (CE žymė)
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/ Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotas atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Katso käyttöohjeet verkkosivustolta eifu.edwards.com	Направете справка с инструкциите за употреба на уебсайта eifu.edwards.com	Consultați instrucțiunile de utilizare de pe site-ul web eifu.edwards.com	Palun lugege kasutusjuhiseid, mille leiame veebisaidilt eifu.edwards.com	Žr. naudojimo instrukcijas, pateikiamas svetainėje eifu.edwards.com
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhiseiga	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Noudata käyttöohjeita, jotka ovat verkkosivustolla eifu.edwards.com.	Спазвайте инструкциите за употреба, дадени на уебсайта eifu.edwards.com	Respectați instrucțiunile de utilizare disponibile pe site-ul web eifu.edwards.com	Järgige kasutusjuhendit, mille leiame veebisaidilt eifu.edwards.com	Vadovaukites naudojimo instrukcijomis, pateikiamomis svetainėje eifu.edwards.com
	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis
	Ei sisällä DEHP:tä	Не съдържа DEHP	Nu conține DEHP	Mitte-DEHP	Ne DEHP

	Suomi	Български	Română	Eesti	Lietuvių
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant
	Viimeinen käyttöpäivä	Срок на годност	Data expirării	Kõlblik kuni	Naudoti iki nurodytos datos
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Ehdollisesti turvallinen magneettikuvauksessa	Безопасно при MR при определени условия	Condiționat RM	Ohutu magnetresonants-tomograafias teatud tingimustel	Sąlyginis MR
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Maahantuojaja	Вносител	Importator	Importija	Importuotojas
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje


Huomautus: kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus:** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simbolai.

Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Vērojama ftalātu klātbūtne vai satur ftalātus: bis(2-etilheksil) ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Conformité Européenne (CE marķējums), ko piešķir TÜV SÜD Product Service GmbH (paziņotā struktūra)	TÜV SÜD Product Service GmbH (onaylı kuruluş) Conformité Européenne (CE İşareti)	Соответствие требованиям директив Европейского Союза (маркировка CE) подтверждено компанией TÜV SÜD Product Service GmbH (нотифицированный орган)	Conformité Européenne (CE znak) kompanije TÜV SÜD Product Service GmbH (ovlaščeno telo)	Conformité Européenne (oznaka CE) tijela TÜV SÜD Product Service GmbH (prijavljeno tijelo)
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlaščeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Evropskoj zajednici / Evropskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatīt lietošanas instrukciju, kas sniegta timeklā vietnē eifu.edwards.com	Web sitesindeki Kullanım talimatlarına başvurun eifu.edwards.com	Следуйте инструкциям по применению на веб-сайте eifu.edwards.com	Pogledajte uputstvo za upotrebu na internet stranici eifu.edwards.com	Pogledajte upute za upotrebu na internetskoj stranici eifu.edwards.com
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Ievērojiet lietošanas instrukciju, kas pieejama timeklā vietnē eifu.edwards.com	Web sitesindeki Kullanım Talimatlarını takip edin eifu.edwards.com	Следуйте инструкциям по применению на веб-сайте eifu.edwards.com	Pratite uputstva za upotrebu na internet stranici eifu.edwards.com	Pridržavajte se uputa za upotrebu na internetskoj stranici eifu.edwards.com
	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	DEHP nesaturošs	DEHP değildir	Не содержит диэтилгексилфталата (ДЭГФ)	Bez DEHP	Ne sadrži DEHP
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена, и см. инструкции по применению	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Daudzums	Miktar	Количество	Količina	Količina
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
Piezīme. Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ Not: Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ Примечание. На этикетках данного изделия могут указываться не все символы. ■ Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ Napomena: na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.					

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Made in Dominican Republic
CCT Critical Care Technologies S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

08/24
10014984002 B / DOC-0131639 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH