Edwards


Fogarty Arterial Embolectomy Catheter with Gate Valve

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	6
Español (es).....	8
Italiano (it).....	10
Nederlands (nl).....	12
Dansk (da).....	14
Svenska (sv).....	16
Ελληνικά (el).....	18
Português (pt).....	20
Česky (cs).....	22
Magyar (hu).....	24
Polski (pl).....	26
Slovensky (sk).....	28
Norsk (no).....	30
Suomi (fi).....	32
Български (bg).....	34
Română (ro).....	36
Eesti (et).....	38
Lietuvių (lt).....	40
Latviešu (lv).....	42
Türkçe (tr).....	44
Русский (ru).....	46
Srpski (sr).....	48
Hrvatski (hr).....	50
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	
	52

English

Fogarty Arterial Embolectomy Catheter with Gate Valve

FG1003

Fogarty Arterial Embolectomy Catheter with Gate Valve

1. Inflation Hub with Gate Valve
2. Catheter Shaft
3. Latex Balloon

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For single use only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

Edwards, Edwards Lifesciences, the stylized E logo, and Fogarty are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

1.0 Description

The Fogarty embolectomy catheters provide a means of clearing emboli and thrombi from the arterial circulatory system through inflation of the balloon and engagement with the arterial wall beyond the vascular obstruction followed by gentle withdrawal of the catheter thereby removing the obstruction from its position.

To remove fibrous or adherent material, alternative devices such as the Fogarty adherent clot and graft thrombectomy catheter are recommended.

The device is intended to be used by medical professionals (surgeons, interventional clinicians such as radiologists), who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

This device aids in clot management and facilitates the removal of emboli and thrombi in the arterial vasculature of the non-central circulatory system.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

2.0 Intended Use/Purpose

The Fogarty arterial embolectomy catheters are intended for the removal of fresh, soft emboli and thrombi from vessels in the arteries of the noncentral circulatory system.

3.0 Indications

The Fogarty embolectomy catheters are indicated for use in adult patients with acute limb ischemia or those requiring therapeutic removal of acute occlusions (such as emboli or thrombi) in the arteries of the non-central circulatory system.

4.0 Contraindications

- The arterial embolectomy catheter should not be used outside the arterial system.
- The arterial embolectomy catheter is not recommended for the removal of fibrous, adherent, or calcified material (e.g. chronic clot, atherosclerotic plaque). The catheter is not designed to withstand the additional pull force needed to remove these materials.
- The arterial embolectomy catheter is not intended for use as a vessel dilator.
- The arterial embolectomy catheter should not be used in the venous system.

5.0 Warnings

- **Balloon rupture and catheter separation as a result of excessive pull force applied to remove adherent material are the most frequent causes of reported failures. The possibility of balloon rupture must be taken into account when considering the risks involved in any embolectomy procedure.**
- **To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume and pull force (see Specification Table).**
- **Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.**
- **To avoid air embolus in case of balloon rupture, air should not be used for balloon inflation. Carbon dioxide is the only recommended gas. See Instructions.**

6.0 Instructions

6.1 Prepare

Fill the inflation syringe with 0.2 ml of carbon dioxide gas. Tightly connect the syringe to the catheter gate valve and set the gate valve to its "Open" position.

Note: For all inflations use the smallest syringe capable of holding the stated maximum gas capacity.

6.2 Purge

Inflate the balloon with carbon dioxide gas to the maximum recommended volume. Pull a vacuum on the syringe. Repeat until all air is removed.

6.3 Inspect

The catheter should be inspected with the balloon inflated during purging. A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

CAUTION: The amount of gas in the syringe should be checked before each inflation. Do not exceed the recommended maximum inflation volume. (See Specification Table).

6.4 Place

With the balloon deflated, insert the catheter into the vessel and beyond the obstructing material.

6.5 Inflate

With the catheter suitably positioned, inflate the balloon with carbon dioxide gas. Balloon inflation should be stopped when the balloon can be felt to engage the arterial wall. Gas seepage through the intact balloon will require frequent adjustment of the inflation volume.

CAUTION: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest inflated balloon diameter that will remove the obstructing material.

6.6 Withdraw

Remove the occlusive material by gently withdrawing the catheter. During withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by controlling the inflation volume.

CAUTION: Do not exceed the Maximum Recommended Pull Force. (See Specification Table).

CAUTION: Inflation of the balloon is associated with a feeling of resistance to the gas injection. When no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

7.0 MRI Information

This product has not been tested for MRI compatibility.

8.0 Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection/sepsis, adverse reaction to device materials, blood loss, local hematomas, hemorrhage, intimal disruption, arterial dissection, perforation and vessel rupture, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

9.0 How Supplied

Contents sterile and nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: . . . 0870 606 2040 - Option 4
In Ireland: . . . 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

14.0 Specification Table

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Model Number	12GV0402F
Maximum Gas Capacity (ml)	0.2
Maximum Recommended Pull Force on Inflated Balloon (lbs.)	0.5 (0.2 kg)
Diameter of Inflated Balloon (mm)	4
Stylet	No
Maximum French Size of Deflated Balloon	3.1F (1.03 mm)
Catheter French Size	2F (0.67 mm)
Length (cm)	45

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




	Nederlands	Dansk	Svenska	Ελληνικά	Português
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site
	Importeur	Importør	Importör	Εισαγωγέας	Importador
SZ	Afmetingen	Størrelse	Storlek	Μέγεθος	Tamanho
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrónicas
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστεριώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogénico
CE 0123	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
UDI	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo

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


	Suomi	Български	Română	Eesti	Lietuvių
	Ei uudelleenkäytettävä	Да не се използва повторно	A nu se reutiliza	Mitte korduskasutada	Nenaudoti pakartotinai
QTY	Määrä	Количество	Cantitate	Kogus	Kiekis
LOT	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
MD	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
EC REP	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/ Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Älä käyttää, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Naudoti iki nurodytos datos
STERILEO	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
BC	Pallon tilavuus	Вместимост на балона	Capacitatea balonului	Balloonni mahutavus	Balionėlio talpa
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
LATEX	Sisältää luonnonkumilateksia tai sen jäämiä	Съдържа или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural	Sisaldab naturaalsed kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
 eifu.edwards.com + 1 888 570 4016	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
SZ	Koko	Размер	Mărime	Suurus	Dydis
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhiseiga	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Älä steriloitu uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepyrogeenne	Nepirogeniškas

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola


	Suomi	Български	Română	Eesti	Lietuvių
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusis priemonės identifikatorius

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Daudzums	Miktar	Количество	Količina	Količina
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom
	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Skatīt lietošanas instrukciju timeklā vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internetskoj stranici
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmērs	Boyut	Размер	Veličina	Veličina
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiecie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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02/24
10000199003 B / DOC-0133125 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH