Fogarty Biliary Balloon Probe

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	7
Español (es).....	10
Italiano (it).....	13
Nederlands (nl).....	16
Dansk (da).....	19
Svenska (sv).....	21
Ελληνικά (el).....	24
Português (pt).....	27
Česky (cs).....	30
Magyar (hu).....	32
Polski (pl).....	35
Slovensky (sk).....	38
Norsk (no).....	41
Suomi (fi).....	43
Български (bg).....	46
Română (ro).....	49
Eesti (et).....	52
Lietuvių (lt).....	54
Latviešu (lv).....	56
Türkçe (tr).....	58
Русский (ru).....	60
Srpski (sr).....	63
Hrvatski (hr).....	65
Figure ■ Figure ■ Abbildung ■ Figura ■ Figura ■ Afbeeldingen ■ Figur ■ Figur ■ Εικόνα ■ Figura ■ Obrázek ■ Ábra ■ Rysunek ■ Obrázok ■ Figurer ■ Kuva ■ Фигура ■ Figura ■ Joonis ■ Paveikslas ■ Attēls ■ Şekil ■ Рисунок ■ Slika ■ Slika.....	67
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	68

English

Fogarty Biliary Balloon Probe

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Fogarty biliary balloon probe is a sterile, single use probe. The device consists of a radiopaque flexible probe shaft with an integrated elastomeric latex balloon at the distal end. A hub at the proximal end is used for balloon inflation. The device facilitates operative procedures in the diagnosis and treatment of symptomatic cholelithiasis.

The device is intended to be used by medical professionals who are trained in the safe use and clinical usage of surgical biliary technologies as part of their respective institutional guidelines.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

2.0 Intended Use/Purpose

The biliary balloon probe is intended for the removal of hepatic calculi and exploration.

3.0 Contraindications

There are no known contraindications.

4.0 Indications

Indications for use include removal of stones and ductal debris and exploration in adult surgical patients.

5.0 Warnings

The use of highly viscous or particulate contrast media is not recommended for balloon inflation as the lumen may become occluded.

The maximum recommended volume for balloon inflation must not be exceeded as over inflation increases the possibility of balloon rupture.

6.0 Precautions

As with all catheterization procedures, complications may occur. These include blood loss/bleeding, peripheral air embolism, infection, adverse reaction to device materials, tissue/vascular trauma/perforation, and pancreatitis.

The balloon must be allowed to deflate as it passes through the ampullary region to avoid damage to the sphincter.

Instances of balloon rupture are minimized if care is taken to inflate the balloon slowly and withdraw the catheter gently.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

7.0 Instructions

7.1 Technique

The probe is introduced with the balloon deflated into the common duct through the choledochotomy. When the desired position has been reached, the balloon is inflated by the introduction of a sterile fluid. It may then be gently withdrawn for the removal of debris or manipulated for further exploration. During the period over which this technique has been developed, the Fogarty biliary balloon probe has been found to be a valuable adjunct to presently available techniques. The Fogarty biliary balloon probe does not, however, have utility in removing stones when the probe cannot pass beyond the calculi or when a calculus is in a sacculus.

The biliary balloon probe aids in easy and rapid removal of amorphous intraductal debris and gravel. Numerous passes may be required to remove all material, but care should be taken not to over-distend the common duct during the period of manipulation. Closed-end irrigation of the common duct may be done by the introduction of the probe into the duodenum, followed by inflation and impingement of the inflated balloon against the sphincter mechanism. The irrigation catheter may then be employed. More effective washing is sometimes accomplished by closure of an open sphincter mechanism.

7.2 Atraumatic Calibration of the Sphincter Mechanism

The probe is introduced into the distal common duct and passed through the sphincter into the duodenum. Once in the duodenum, the balloon is inflated and then carefully palpated to make certain that it has passed through the sphincter mechanism. While the surgeon keeps the fingers of one hand placed on the balloon in the duodenum, the balloon probe is withdrawn until it impinges at the area of the sphincter. With light traction on the probe, the surgeon deflates the balloon until it passes through the sphincter with ease. As the balloon is passed through the sphincter, the quantity of fluid contained within the syringe is noted. When the balloon re-enters the common duct, the balloon is inflated until it engages the duct wall, and then gently withdrawn. The same amount of fluid which was noted as the probe passed through the sphincter is reintroduced into the balloon. Calibration of the balloon diameter accurately reflects the size of the sphincter. The common duct can be calibrated in the same manner. Mobilization of the duodenum greatly facilitates this maneuver. This technique aids in identification of the retroduodenal and pancreatic position of the common duct and accurately locates the sphincter of Oddi.

7.3 For Removal of Intrahepatic Stones

The probe can be bent with the stylet in place to facilitate introduction into the right or left hepatic duct. The stylet may be removed after the probe enters into the hepatic duct. The probe is then advanced as far as possible with the balloon deflated. As the balloon is gently inflated, the probe is simultaneously removed. After noting when the probe engages the wall of the ductal system, the inflated balloon is carefully withdrawn through the incision in the common duct. Several passes may be necessary. The final exploratory passage in these radicles is carried out by insertion of probes into both the right and left hepatic ducts. The two probes are then simultaneously removed. This maneuver prevents an isolated stone from falling from one hepatic duct into the other during removal. Undue manipulation, probing, or over-distention of the balloon may cause bleeding and should be avoided.

7.4 Identification of the Common Duct

Where there is difficulty in identifying the entire extent of the common duct because of previous surgery or inflammation, the use of the probe allows identification of the duct and facilitates dissection. If the gall bladder remains and the cystic duct has been identified, a small incision in the cystic duct with introduction of the probe into the distal common duct is easily accomplished. Similarly, if the common duct has been identified by needle aspiration, a small puncture wound is made in the anterior surface of the common duct and the balloon probe is introduced. After entrance into the common duct, palpation of the inflated balloon in various locations facilitates dissection and identification of the common duct and associated anatomy.

7.5 Use of Stylet

The stylet supplied with the biliary balloon probe is of a designed stiffness intended to avoid perforation and yet maintain enough shape to facilitate introduction of the probe into areas of sharp angulation.

When introducing the probe into the common duct, entrance into the duodenum may be accomplished with the stylet in place and angulated at approximately 30°. If difficulty is encountered in passing through the sphincter mechanism, the stylet should be unscrewed until it lies free of the probe body. A gentle probing motion with the stylet in this position will allow for passage into the duodenum if the sphincter mechanism is large enough to accommodate the probe body. When extreme pliability of the distal tip of the instrument is desired, the stylet should be withdrawn 1 cm. A probe possessing a longer length of intermediate pliability is achieved by withdrawing the stylet a greater length.

8.0 Balloon Inflation

The balloon should be inflated with a sterile, blood compatible fluid that may include radiopaque solutions, if they are highly dilute and non particulate.

Before use, the probe with the balloon inflated to its maximum recommended volume should be inspected by the surgeon.

The amount of fluid in the syringe should be checked before each inflation. If the amount exceeds the balloon's stated capacity, remove the syringe from the probe and refill to the proper volume. The amount of fluid remaining in the catheter lumen must be taken into account.

9.0 MRI Information

This product has not been tested for MRI compatibility.

10.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

11.0 Storage

Store in a cool, dry place.

12.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

13.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

14.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

15.0 Specification Table

French Size	5F (1.67 mm)	6F (2.0 mm)
Maximum Liquid Capacity (ml)	0.4	1.25
Maximum Air Capacity (ml)	0.8	2.5
Maximum Inflated Balloon Diameter (mm)	8	13

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




	Suomi	Български	Română	Eesti	Lietuvių
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
MD	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
SZ	Koko	Размер	Mărime	Suurus	Dydis
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinaï nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
CE 0123	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
UDI	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusiis priemonės identifikatorius
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
#	Mallinumero	Номер на модела	Număr model	Mudeli number	Modelio numeris

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
QTY	Daudzums	Miktar	Количество	Količina	Količina
LOT	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Izmērs	Boyut	Размер	Veličina	Veličina
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik


Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen
 ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli
 ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus
 ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámká:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiety niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámká:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almayış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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05/23
10000132003 B / DOC-0133127 B
© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH