HemoSphere
Advanced Monitor

Quick start guide

Edwards

Edwards

DIRECTORY

English	1	Português	19	Eesti	37
Français	3	Česky	21	Lietuvių	39
Deutsch	5	Magyar	23	Latviešu	41
Español	7	Polski	25	Türkçe	43
Italiano	9	Slovensky	27	Русский	45
Nederlands	11	Norsk	29	Srpski	47
Dansk	13	Suomi	31	中文	49
Svenska	15	Български	33	한국어	51
Ελληνικά	17	Română	35	日本語	53

English

HemoSphere Advanced Monitor Quick Start Guide


Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

This quick start guide provides brief installation instructions for the HemoSphere advanced monitor. Refer to the HemoSphere advanced monitor operator's manual, included on the enclosed USB memory stick, or at the website eifu.edwards.com, for detailed monitoring instructions. The electronic version of the operator's manual is in a portable document format (PDF). Adobe Acrobat Reader is needed for viewing the file (<https://get.adobe.com/reader/>).

Description

The HemoSphere advanced monitor provides Edwards hemodynamic monitoring technologies through modular connections. Refer to the operator's manual for currently available connections and associated monitored parameters.

Edwards, Edwards Lifesciences, the stylized E logo, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

The HemoSphere advanced monitoring platform allows you to see and interact with patient hemodynamic parameters. In conjunction with the compatible sensors and predictive decision support software, the modular HemoSphere platform facilitates proactive clinical decision-making and insight for individualized patient care.

Device performance, including functional characteristics, has been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

Indications for Use

See the HemoSphere advanced monitor operator's manual for details on available modules/cables and their specific indications for use and intended use statements.

Intended Use

The HemoSphere advanced monitor is intended to be used in combination with a compatible Edwards hemodynamic monitoring technology module and/or cable and accompanying Edwards accessories and/or disposables.

Warnings

Compliance to IEC 60601-1 is only maintained when the HemoSphere advanced monitor is connected as described. Connecting external equipment or configuring the system in a way not described here and in the HemoSphere advanced monitor operator's manual will not meet this standard. Failure to use the device as instructed may increase the risk of electrical shock to the patient/operator.

Do not modify, service or alter the product in any way. Servicing, alteration or modification may affect patient/operator safety and/or product performance.

Read the operator's manual included on the enclosed USB stick or at eifu.edwards.com carefully before attempting to use the HemoSphere advanced monitor.

This product contains metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

Precautions

Do not use if the packaging or contents are damaged. Damage may include cracks, scratches, dents or any signs that the monitor, modules, cables, cable housing or system accessories may have been compromised during shipping or unpacking.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

HemoSphere Advanced Monitor Setup

Refer to the HemoSphere advanced monitor operator's manual included on the enclosed USB memory stick or located on the Edwards website (eifu.edwards.com) for comprehensive monitoring procedures, cautions and warnings.

1. Securely position the HemoSphere advanced monitor or mount on a compatible stand.

2. Open the battery door on the lower right panel of the HemoSphere advanced monitor ① and insert the battery.
3. Insert desired hemodynamic technology module ② and/or cable ③. Refer to module or cable instructions for proper insertion/connection procedures.
4. Connect the power cord to the rear panel of the monitor ④ and plug into a hospital grade outlet.

Note: After inserting the battery for the first time, the battery must be fully charged and conditioned. Refer to operator's manual for additional information.

5. Press the power button on the front of the monitor ⑤.

Note: Upon initial startup the display language can be selected which affects the time/date format and units of measurement.

6. Follow instructions contained in the HemoSphere advanced monitor operator's manual for desired hemodynamic technology module or cable monitoring.

MRI Information

MR Unsafe

Do not use the HemoSphere advanced monitor in an MR environment. The monitor is MR unsafe since the device contains metallic components, which can experience RF-induced heating in the MRI environment.

Operating conditions

Temperature: 10 to 32.5 °C
Humidity range: 20 to 90% non-condensing
Altitude (Atmospheric Pressure):
0 ft/0 m (1013 hPa) to
3048 m/10,000 ft (697 hPa)
IP Rating: IPX1

Maintenance

The HemoSphere advanced monitor is a multi-use device and can be surface cleaned using a cloth dampened with cleaning agents containing 70% isopropyl alcohol, 2% glutaraldehyde, 10% bleach solution, or quaternary ammonium solution.

Do not use any other cleaning agents, spray, or pour cleaning solution directly on the monitor.

Storage

Temperature: -18 to 45 °C
Relative humidity: 90% non-condensing at 45 °C
Altitude: 0 to 20,000 feet (6096 m)

Note: For long term storage (>3 months), a maximum storage temperature of 35 °C is recommended. Battery capacity starts to degrade with extended exposure above 35 °C.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: ... 0870 606 2040 - Option 4
In Ireland: 01 821 1012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

This product is manufactured and sold under one or more of the following US patent(s): US Patent No. 6,371,923; and corresponding foreign patents. Additional patents pending.

Prices, specifications and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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© Copyright 2019, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences LLC