Edwards

DIRECTORY						
English (en)	1	Português (pt)	18	Română (ro)	33	
Français (fr)	3	Česky (cs)	20	Eesti (et)	35	
Deutsch (de)	5	Magyar (hu)	22	Lietuvių (lt)	37	
Español (es)	7	Polski (pl)	24	Latviešu (lv)	39	
Italiano (it)	9	Slovensky (sk)	26	Türkçe (tr)	41	
Nederlands (nl)	11	Norsk (no)	28	Русский (ru)	42	
Dansk (da)	13	Suomi (fi)	29	Srpski (sr)	44	
Svenska (sv)	14	Български (bg)	31	Hrvatski (hr)	46	
Ελληνικά (el)	16					
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figuren ■ Figuren ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figuren ■ Kuvat ■ Φωτογραφίες ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike						48
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sumbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola						55

English

Chandler Transluminal V-Pacing Probe D98100 and Flex-Tip Transluminal A-Pacing Probe D98500

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific country. D98500 is not available in EU.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

For figures, see Figure 1 on page 48 through Figure 12 on page 54.

Chandler transluminal V-pacing probe D98100 (Chandler probe):

For use with Swan-Ganz Pacerport (Model 931F75) or A-V Pacerport (Model 991F8) catheters for ventricular pacing only.

Developed in collaboration with John P. Chandler, M.D., Assistant Clinical Professor of Medicine, Yale School of Medicine, New Haven, CT.

Flex-Tip transluminal A-pacing probe D98500 (Flex-Tip probe):

Use only with Swan-Ganz thermomodulation A-V Pacerport catheter (Model 991F8). For atrial pacing only.

1.0 Description

The probe can be inserted with or without the aid of fluoroscopy.

The Chandler transluminal V-pacing probe is recommended for use *in situ* for up to 72 hours. After the catheter is inserted and floated into the pulmonary artery with the right ventricular (RV) port (19 cm from the distal tip) properly placed 1 to 2 cm distal to the tricuspid valve, the Chandler pacing probe is inserted into the Pacerport or A-V Pacerport catheter's RV lumen and advanced into the ventricle for endocardial pacing or the A-pacing probe is inserted into the A-V Pacerport catheter's RA lumen and advanced into the atrium for endocardial pacing.

The probe/J-tip wire is a bipolar, coaxial, wire construction composed of a stainless steel round wire and a PTFE coated, coiled flat wire.

As part of the insertion procedure this product is used for ECG detection during placement, but is not intended for ECG monitoring.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The Swan-Ganz Chandler transluminal V-pacing probe and Flex-Tip transluminal A-pacing probes provide temporary electrical cardiac stimulation to facilitate return to and/or continuance of normal cardiac rhythm.

The intravascular catheter is inserted through the central vein to connect to the right side of the heart and advance towards the pulmonary artery. Route of insertion can be internal jugular vein, femoral, antecubital, or subclavian vein. The cardiac structures in contact are the right atria, right ventricle and pulmonary artery as well as the circulatory system.

This device has not yet been tested in the pediatric population or in pregnant or lactating women.

2.0 Intended Use/Purpose

The Swan-Ganz Chandler transluminal V-pacing probes are intended to serve as therapeutic tools for transvenous ventricular temporary pacing when used with a Swan-Ganz Pacerport catheter. The probe can also be used for intraventricular and electrocardiographic detection during catheter placement.

3.0 Indications

The Chandler transluminal V-pacing probe is a transluminal bipolar pacing probe, which must be used with a Swan-Ganz Pacerport catheter is indicated for use in adult surgical and/or critically ill cardiac patients. This includes, but is not limited to, patients with arrhythmias, acute myocardial infarction, cardiac surgery and coronary angiography.

4.0 Contraindications

Although there are no absolute contraindications to the use of temporary endocardial pacing electrodes, relative contraindications may include patients with recurrent sepsis or with a

hypercoagulable state where the electrode could serve as a focus for septic or bland thrombus formation.

Use of the probe is contraindicated in patients with small hearts where the RV port of the Swan-Ganz Pacerport or A-V Pacerport catheter cannot be placed into the right ventricle without spontaneously wedging the catheter in the pulmonary artery with the balloon deflated. Also, the Chandler probe is not intended for use with any catheter except the Swan-Ganz Pacerport or A-V Pacerport catheter. The Flex-Tip atrial probe is not intended for use with any catheter except the Swan-Ganz thermomodulation A-V Pacerport catheter.

These products contain metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

5.0 Warnings

D98500 may contain phthalates, specifically DEHP (Bis (2-ethylhexyl) phthalate), which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

These devices may contain the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0. Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

This device is designed, intended, and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity, and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Do not modify or alter the product in any way. Alteration or modification may affect patient/operator safety or product performance.

As part of the insertion procedure this product is used for ECG detection during placement, but is not intended for ECG monitoring.

6.0 Cautions

Clinicians using the device should be familiar with the device and understand its applications prior to use.

When handling indwelling leads, the terminal pins or exposed metal (on the product) are not to be touched nor be allowed to contact electrically conductive or wet surfaces, to avoid electrical shock to patient or clinician.

7.0 Insertion

Visually inspect for breaches of packaging integrity prior to use.

The Pacerport (Model 931F75) or A-V Pacerport (Model 991F8) catheter can be inserted at the patient's bedside, usually without the aid of fluoroscopy, by continuous pressure monitoring from the distal and right ventricular lumens. A unipolar electrocardiogram can be recorded from the distal tip electrode by connection to a V lead of a properly isolated electrocardiograph.

Use of the Edwards introducer will facilitate optimal insertion of the probes up to 72 hours after catheter insertion. Use of other manufacturer's introducers will not guarantee probe passage through the catheter within this 72-hour period.

8.0 Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (Type CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

- Swan-Ganz Pacerport (Model 931F75) or A-V Pacerport catheter (Model 991F8)
- Compatible percutaneous sheath hemostasis valve introducer tray, kit or single assembly
- Chandler transluminal V-pacing probe, Model D98100
- Flex-Tip transluminal A-pacing probe, Model D98500
- Ventricular demand external pulse generator, demand external pulse generator (atrial or A-V sequential as appropriate for pacing needs)
- External pulse generator cable adapters
- ECG recorder
- Sterile flush system and pressure transducers
- Bedside ECG and pressure monitor system

In addition, the following items should be immediately available if complications arise during catheter or probe insertion: antiarrhythmic drugs, defibrillator, and respiratory assist equipment.

9.0 Catheter Insertion and Placement

The catheter can be inserted using percutaneous technique cutdown through a jugular, subclavian, femoral (Flex-Tip probe), or antecubital vein. A protective catheter sheath is recommended to aid in maintaining sterility when repositioning of the catheter is necessary.

To facilitate subsequent insertion, insertion of the catheter is best accomplished by using two pressure transducers; one transducer is connected to the distal (PA) lumen, the other to the right ventricular (RV) lumen, which terminates 19 cm from the tip. The ideal placement of the RV port for probe placement is 1 to 2 cm distal to the tricuspid valve. A radiopaque marker is provided at the RV port to confirm port placement by X-ray or fluoroscopy. Refer to the Pacerport or A-V Pacerport catheter package insert for detailed instructions on insertion.

Flex-Tip probe: It is recommended that femoral insertion of the A-V Pacerport catheter (Model 991F8) and Flex-Tip and Chandler probes be done under fluoroscopy.

Also, femoral insertion should be used only when short-term pacing is required (e.g., cath lab procedures), because of the possible placement of the Chandler probe in the RV outflow tract. Please note that with the femoral approach, the RA pacing port may be positioned in the inferior vena cava. This position may prohibit proper atrial probe placement or may require an extended length of probe insertion to obtain atrial capture. Refer to the Pacerport (Models 991F8, 931F75) catheter package insert for detailed instructions on insertion.

A protective catheter shield is recommended to aid in maintaining sterility when repositioning of the catheter is necessary.

Advance the catheter into the pulmonary artery while continuously monitoring both PA and RV lumen pressures. When the catheter tip is at the wedge position, RV port location may vary according to heart size.

- Normal Size Hearts:** At the wedge position, the RV lumen shows RV tracing. Deflate the balloon. Pull the catheter back until the RV port is in the right atrium. Then readvance the catheter until the port is 1 to 2 cm distal to the tricuspid valve.
- Small Hearts:** The RV lumen shows a right atrial (RA) pressure tracing at the wedge position. Deflate the balloon. Advance the catheter slowly while closely monitoring the PA and RV lumen pressures until an RV pressure tracing is first obtained from the RV lumen. Continue to advance the catheter 1 to 2 cm distal to the tricuspid valve for optimal placement of the RV port.

If RV port placement in the right ventricle results in spontaneous wedging, catheter repositioning is required. To insert the probe into the right ventricle, one can advance the probe a centimeter at a time after withdrawing the catheter a centimeter at a time until a pulmonary artery pressure tracing is seen continuously from the distal lumen.

WARNING: In some patients, the catheter might spontaneously wedge (with the balloon deflated) before positioning of the RV port in the right ventricle. Discontinue advancing the catheter. This pacing system is not suitable for use in these patients; however, the catheter can still be used for pressure monitoring, blood sampling, fluid infusion and cardiac output determinations. Do not attempt to insert the probe if the RV port is in the RA. Damage to the tricuspid valve may result. Always make certain that the RV port is inside the ventricle before inserting the probe.

- Enlarged Hearts:** Wedge position is not yet reached when the RV lumen shows an RV pressure tracing. Continue to advance the catheter to obtain a wedge pressure recording. Note the distance the catheter is advanced between the first RV pressure tracing from the RV lumen and wedge position. Deflate the balloon. Withdraw the catheter until an RA pressure is obtained from the RV lumen, then readvance the catheter until the RV port is 1 to 2 cm distal to the tricuspid valve. The catheter tip should be in the pulmonary artery. In these patients it may not be possible to obtain capture and wedge pressure measurements simultaneously.

WARNING: If the RV port is too distal, then the probe may exit the RV port pointed toward the RV outflow tract. This may result in poor thresholds, unstable pacing, and potential damage to the outflow tract and pulmonic valve.

10.0 Pacing Probe Insertion and Placement

WARNING: Handle the probe using sterile technique. Chandler probe: Be certain that the probe is inserted only into the RV lumen (clear extension tube with orange Luer-lock hub). Do not insert the probe into either the proximal (RA) or distal (PA) lumen.

Flex-Tip probe: Be certain that the atrial probe is inserted only into the RA/probe lumen (clear extension tube with yellow Luer-lock hub). Do not insert the probe into the proximal injectate, RV/probe, or distal (PA) lumens.

Before inserting the probe make sure that the catheter portion outside the patient is not coiled, as this will make probe insertion difficult.

Edwards, Edwards Lifesciences, the stylized E logo, Flex-Tip, Pacerport, Swan, and Swan-Ganz are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Step	Procedure
1	Verify lumen patency.
2	Connect the catheter's RV lumen hub to a pressure transducer and verify proper placement of the RV port (1 to 2 cm distal to the tricuspid valve) (See Figure 3 on page 50). For Flex-Tip probe, this will position the RA pacing/infusing port in the right atrium (See Figure 3 on page 50). To prevent catheter movement, secure the catheter at the insertion site. When using a Flex-Tip probe, a Chandler V-pacing probe may be inserted into the RV lumen at this time using the instructions provided for the Chandler probe.
3	Chandler probe: Open the probe package and retract the probe tip into the Tuohy-Borst (T-B) adapter by turning the carousel clockwise . Flex-Tip probe: Open the probe package, gently remove the protective blue tube from the probe tip, and retract the probe "J"-tip into the Tuohy-Borst (T-B) adapter by turning the carousel clockwise .
4	Connect the T-B adapter to the orange RV lumen hub (Chandler probe) or the yellow RA lumen hub (Flex-Tip probe). Be careful not to damage the probe tip (See Figure 4 on page 50).
5	Advance the probe until its depth reference mark (black band) is placed at the zero mark on the clear extension tube of the Chandler RV lumen or the Flex-Tip RA lumen probe (See Figure 5 on page 50). Because of manufacturing tolerances, the tip of the probe is now between the RV/RA port and a point 2 cm proximal to the port. The probe is ready to be advanced into the RV/RA. Note: There may be some resistance as the probe passes through the hemostasis valve of the introducer, curves in the catheter at the subclavian-SVC junction and at the RV/RA port. Resistance at any other point may indicate that the catheter is kinked. Do not force the probe if resistance is met. Note: The J-tip atrial probe encounters more friction than the straight V- pacing probe as it is passed into the catheter lumen. Insertion difficulty may be related to crimping of the catheter tubing at the hemostasis valve of the introducer. If this occurs, slight manipulation of the catheter relative to the introducer will usually allow probe passage. Precaution: The PTFE coating on the probe is a lubricating agent, not an electrical insulator. Do not allow the probe surface to come into contact with any line power equipment because of potential current leakage due to faulty ground, which can cause ventricular fibrillation. When not connected to the external pulse generator, the electrode pin connectors should remain protected.
6	Attach the distal end of the probe contamination sheath to the T-B adapter. Remove the probe and T-B adapter from the dispenser, and, to help maintain sterility of the probe, attach the other end of the sheath to the proximal end of the probe (See Figure 6 on page 51). a) Chandler Probe Only: Connect the distal electrode to a V lead of a properly isolated electrocardiograph (See Figure 7 on page 51 and Figure 8 on page 52). Under continuous ECG monitoring, advance the probe several centimeters until ST segment elevation of the ECG indicates contact with the endocardium. Note: ST elevation is usually seen with the probe out 4 to 5 cm. If the probe is out more than 10 cm, the probe may be in the RV outflow tract. Pull the probe back to 4 to 5 cm and reposition in the RV apex.
7	Connect the distal and proximal electrodes to the negative and positive pulse generator terminals, respectively (See Figure 9 on page 52). For Chandler probe: Determine the pacing threshold. A threshold of 1.0 to 2.0 mA is generally an indication of proper electrode placement. An initial threshold greater than 5 mA indicates poor probe placement, possibly in the RV outflow tract. Withdraw the probe several centimeters and reposition the probe in the RV apex. Best pacing thresholds are obtained with the probe approximately 5 cm out of the RV port. Stable pacing cannot usually be obtained with the probe out less than 3 cm. Note: If transient multifocal PVCs or V tach persist during (or after) probe insertion, pull the catheter back 1 to 2 cm and advance the probe to the RV apex. (See PVCs During/After Insertion.) For Flex-Tip probe: Set the desired mA and rate for atrial capture (See Figure 8 on page 52). Note: To facilitate a connection between the probe and pulse generator, a cable adapter may be required. a) Flex-Tip probe Only: Advance the probe several centimeters until atrial capture is achieved, and determine pacing threshold (See Figure 8 on page 52). The mean initial threshold in clinical studies was 4.4 mA (range 1.1 to 12.0 mA). In clinical studies, the best initial pacing position was with the probe an average of 5.7 cm (range 2 to 11 cm) out of the RA port. Note: The atrial appendage and coronary sinus are often used pacing sites in the atrium. However, the probe does not necessarily go to the atrial appendage - it may go to the lateral atrial wall. It may not be necessary to deploy the full "J" into the atrium. Note: To enhance long-term pacing stability, advance the probe 2 cm extra beyond the best initial pacing position. This creates pressure against the atrial wall and may help maintain pacing in the event of excessive patient movement, vigorous respirations, and changes in patient's position.
8	Firmly tighten the compression nut to secure the probe in place (See Figure 10 on page 53). Aspirate any air from the side port. Institute a continuous or intermittent heparinized flush through the side- port fitting of the T-B adapter.
9	Using a 3-way Luer-Lock stopcock, connect the side port of the T-B adapter to a continuous saline flush device. Using the 5 ml syringe provided, aspirate any air from the sideport, then flush the lumen (See Figure 11 on page 53). Note: When the Chandler probe is in the RV lumen, do not infuse solutions at a rate greater than 30 ml/hr because the solution may back up into the probe contamination sheath.
10	Obtain a chest x-ray film as soon as possible after insertion to document the initial placement.

11.0 MRI Information

MR Unsafe

The Chandler and Flex-Tip probe are MR unsafe as the result of the device containing metallic components, which experience RF-induced heating in the MRI environment; therefore the device poses hazards in all MRI environments.

12.0 Complications

12.1 Loss of Capture

Loss of capture may occur due to inadvertent pulling of the probe off the endocardium, poor initial placement (probe in RV outflow tract), myocardial perforation, vigorous respirations, or patient movement. If the probe is pulled off the endocardium, is in the RV outflow tract, or has perforated the myocardium (see **Ventricular Perforation**), reposition the probe. Transient loss of capture after patient movement is corrected by placing the patient in a supine position and, if necessary, increasing the threshold or repositioning the probe.

12.2 Inability to Wedge the Balloon

If placement required that the catheter be pulled back from initial wedge position, then it may not be possible to obtain wedge pressure upon balloon inflation. Monitor the PA diastolic pressure instead of wedge whenever possible. Intermittent pacing may occur following balloon inflation for wedging.

However, capture is usually regained upon balloon deflation without increasing the pacing threshold. If wedging is required and pacing is no longer needed, advance the catheter to wedge position after turning the pulse generator off and withdrawing the probe completely into the catheter.

12.3 PVCs During/After Probe Insertion (Chandler probe)

Transient multifocal PVCs or V tach may occur due to irritation of the endocardium by the probe tip. Additional probe advancement or catheter manipulation usually resolves the PVCs. If PVCs or V tach persist, pull the catheter back 1 to 2 cm and advance the probe toward the RV apex.

12.4 Inadvertent Atrial Pacing (Chandler probe)

Atrial pacing may result if the RV port is in the atrium rather than in the ventricle. In addition, atrial pacing may occur due to catheter or probe movement into the right atrium. Withdraw the probe completely into the catheter and advance the catheter's RV port into the ventricle. Readvance the probe into the ventricle.

12.5 Inadequate Sensing (Chandler probe)

Inadequate sensing of the demand pulse generator may occur if the probe is partially in the atrium. Reposition the probe in the RV apex to improve sensing after withdrawing the probe into the catheter and advancing the catheter 1 to 2 cm distal to the tricuspid valve.

12.6 Ventricular Perforation (Chandler probe)

Ventricular perforation has been reported with temporary transvenous pulse generator electrodes and usually results in intermittent or failed cardiac pacing. Treatment of ventricular perforation is withdrawal of the electrode back into the ventricle. Perforation can be diagnosed by connecting the distal electrode to the V lead of a battery-powered electrocardiograph. As the electrode is slowly withdrawn, a ventricular ectopic beat occurs when the tip is in the myocardium. The ST segment is markedly elevated and the T wave deeply inverted, producing an endocardial "current of injury" pattern. In rare instances, cardiac tamponade may result. The probe tip is designed to be very soft to minimize injury to the ventricular endocardium. However, to prevent potential damage to the endocardium during open heart surgery, withdraw the probe into the catheter before manipulating the heart.

12.7 Diaphragmatic Stimulation (Flex-Tip probe)

Diaphragmatic stimulation may occur due to phrenic nerve stimulation. Either reposition the probe or catheter or both and/or reduce the mA setting on the pulse generator to minimize or eliminate stimulation of the diaphragm.

12.8 Atrial Perforation (Flex-Tip probe)

The incidence of atrial perforation from temporary atrial pacing is unknown. In one study of 100 patients who had a J-tip electrode, no cases of atrial perforation were reported. In rare instances, cardiac tamponade may result. The probe tip is designed to be very soft to minimize injury to the endocardium. However, to prevent potential damage to the endocardium during open heart surgery, withdraw the probe into the catheter before manipulating the heart.

12.9 Other Complications

Pulmonary artery catheters have also been associated with cardiac arrhythmias, blood loss, cardiac structure/wall injury or damage, tricuspid and pulmonic valve damage, hematoma, perforation of the pulmonary artery, embolism, thrombosis, thrombophlebitis, anaphylaxis, cardiac tissue/artery burn, sepsis/infection, pulmonary infarction, inadvertent stimulation/pacing, knotting, right bundle branch block, complete heart block, and pneumothorax.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

13.0 How Supplied

Chandler transluminal V-pacing and Flex-Tip transluminal A-pacing probes are supplied sterile (unless otherwise stated) and nonpyrogenic. Do not use if package has been previously opened or damaged.

Note: Probes are for single use only. Do not clean or resterilize a used probe.

14.0 Packaging

The probe is provided preloaded in a packaging dispenser designed to aid in probe insertion and to help maintain sterility during insertion. It is therefore recommended that the probe remain inside the package until use.

15.0 Storage

Store in a cool, dry place.

16.0 Operating Conditions

Intended to operate under physiological conditions of the human body.

17.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the recommended time may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

18.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

19.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the latest version of the monitoring system operator's manual for more information.

Refer to the symbol legend at the end of this document.

20.0 Specifications:

Chandler Transluminal V-Pacing Probe (Model D98100)	
For use with a Swan-Ganz Pacerport or A-V Pacerport catheter for ventricular pacing only.	
Usable length (cm)	
Total	135
In ventricle	15
Body Diameter (F)	2.4 (0.80 mm)
Electrodes	Stainless steel with pin connectors (0.08 inch diameter, or 0.2 cm) at proximal end
Distal:	
Length (cm)	1.3
Proximal:	
Length (cm)	15
Contents	
Chandler transluminal V-pacing probe	
ECG adapter (see Figure 12 on page 54)	
Syringe, 5 ml Luer-Lock	
Sterile drape, folded, 18" x 26" (45.72 cm x 66.04 cm)	
Contamination sheath	
Dispenser carousel	

All specifications given are nominal values.

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

Figure 9: Connection of leads to pulse generator

- Figure 9: Connexion des électrodes au générateur d'impulsions
- Abbildung 9: Anschluss der Elektroden an den Impulsgeber
 - Figura 9: Conexión de los electrodos a un generador de pulso
 - Figura 9: Collegamento degli elettrodi al generatore di impulsi
 - Afbeelding 9: Aansluiting van de leads op de pulsgenerator
 - Figur 9: Tilslutning af ledninger til pulsgeneratoren
 - Figur 9: Anslutning av ledare till pulsgenerator
 - Εικόνα 9: Σύνδεση των ηλεκτροδίων στη γεννήτρια ερεθισμάτων
 - Figura 9: Ligação dos eletrodos ao gerador de impulsos
 - Obrázek 9: Připojení elektrod ke generátoru pulsů
 - 9. ábra: A vezetékek csatlakoztatása az impulzusgenerátorhoz
 - Rysunek 9: Podłączanie odprowadzeń do generatora impulsów
 - Obrázok 9: Pripojenie zvodov ku generátoru impulzov
 - Figur 9: Kobling av ledninger til pulsgenerator
 - Kuva 9: Johtojen kytkeminen pulssigeneraattoriin
 - Фигура 9: Свързване на отведения към импулсен генератор
 - Figura 9: Conexiunea derivațiilor la generatorul de impulsuri
 - Joonis 9. Elektroodide ühendamine impulsigeneraatoriga
 - 9 pav. Laidų prijungimas prie impulsų generatoriaus
 - 9. attēls. Pievadu savienošana ar impulsu generatoru
 - Şekil 9: Elektrotların puls üreticine bağlanması
 - Рис. 9. Соединение электродов и кардиостимулятора
 - Slika 9: Povezivanje elektroda sa generatorom pulsa
 - Slika 9: Povezivanje elektroda na generator pulsa


Figure 10: Securing probe to T-B adapter ■ Figure 10 : Fixation de la sonde à l'adaptateur T-B ■ Abbildung 10: Fixieren der Sonde am T-B-Adapter ■ Figura 10: Fijación de la sonda al adaptador T-B ■ Figura 10: Fissaggio della sonda all'adattatore T-B ■ Afbeelding 10: De sonde vastzetten op de T-B -adapter ■ Figur 10: Fastgørelse af probe til T-B adapter ■ Figur 10: Anslutning av sond till T-B adapter ■ Εικόνα 10: Στερέωση μήλης σε προσαρμογέα T-B ■ Figura 10: Fixação da sonda ao adaptador T-B ■ Obrázek 10: Připevnění sondy k adaptéru T-B ■ 10. ábra: A sonda rögzítése a T-B-adapterhez ■ Rysunek 10: Zabezpieczanie sondy do adaptera typu T-B ■ Obrázok 10: Zabezpečenie sondy k adaptéru T-B ■ Figur 10: Festing av proben til en T-B-adapter ■ Kuva 10: Sondin kiinnittäminen T-B-sovittimeen ■ Фигура 10: Фиксация на сондата към адаптера T-B ■ Figura 10: Fixarea sondei de adaptorul T-B ■ Joonis 10. Anduri kinnitamine T-B adapteri külge ■ 10 pav. Zondo pritvirtinimas prie T-B adapterio ■ 10. attēls. Zondes nostiprināšana pie T-B adaptera ■ Şekil 10: Proben T-B adaptörüne bağlanması ■ Рис. 10. Закрепление соединения зонда и переходника T-B ■ Slika 10: Pričvrščanje sonde za T-B adapter ■ Slika 10: Pričvrščanje sonde na adapter T-B

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

Figure 12: PACING PRODUCTS KITS (Models D98100, D98500)

- Figure 12: KITS DE PRODUITS DE STIMULATION (modèles D98100, D98500)
- Abbildung 12: PRODUKTSETS FÜR DIE STIMULATIONSTHERAPIE (Modelle D98100, D98500)
- Figura 12: KITS DE PRODUCTOS DE ESTIMULACIÓN (Modelos D98100 y D98500)
- Figura 12: KIT DI PRODOTTI PER STIMOLAZIONE (modelli D98100, D98500)
- Afbeelding 12: SETS VOOR STIMULATIEPRODUCTEN (modellen D98100, D98500)
- Figur 12: PACING-PRODUKTSÆT (model D98100, D98500)
- Figur 12: STIMULERINGSPRODUKTSATSER (modellerna D98100 och D98500)
- Εικόνα 12: KIT ΠΡΟΪΟΝΤΩΝ ΒΗΜΑΤΟΔΟΤΗΣΗΣ (μοντέλα D98100, D98500)
- Figura 12: KITS DE PRODUTOS DE ESTIMULAÇÃO (modelos D98100, D98500)
- Obrázek 12: SOUPRAVY STIMULAČNÍCH PRODUKTŮ (modely D98100, D98500)
- 12. ábra: INGERLŐTERMÉK-KÉSZLETEK (D98100, D98500 típusok)
- Rysunek 12: ZESTAWY PRODUKTÓW DO STYMULACJI (modele D98100, D98500)
- Obrázok 12: SÚPRAVY STIMULAČNÝCH VÝROBKOV (Modely D98100, D98500)
- Figur 12: PACINGPRODUKTSETT (modellene D98100, D98500)
- Kuva 12: TAHDISTUSTUOTTEIDEN PAKKAUKSET (Mallit D98100, D98500)
- Фигура 12: КОМПЛЕКТИ ПРОДУКТИ ЗА ПЕЙСИРАНЕ (модели D98100, D98500)
- Figura 12: SETURI DE PRODUSE DE STIMULARE (modelele D98100, D98500)
- Joonis 12. STIMULEERIMISTOODETE KOMPLEKTID (mudelid D98100, D98500)
- 12 pav. STIMULIAVIMO GAMINIŲ RINKINIAI (modeliai D98100, D98500)
- 12. attēls. KARDIOSTIMULĀCIJAS IZSTRĀDĀJUMU KOMPLEKTI (modeļi D98100, D98500)
- Şekil 12: PACING ÜRÜNÜ KİTİLERİ (D98100, D98500 Modelleri)
- Рис. 12. КОМПЛЕКТ ИЗДЕЛИЙ ДЛЯ КАРДИОСТИМУЛЯЦИИ (модели D98100, D98500)
- Slika 12: KOMPLETI PROIZVODA ZA PEJSING (modeli D98100, D98500)
- Slika 12: KOMPLETI PROIZVODA ZA ELEKTROSTIMULACIJU (modeli D98100, D98500)

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen

	English	Français	Deutsch	Español	Italiano	Nederlands
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße der Einführhilfe	Tamaño de introductor mínimo	Misura minima dell'introduttore	Minimum introductergrootte
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig
	Size	Taille	Größe	Tamaño	Misura	Afmetingen
	Store in a cool, dry place	Conservé dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarriersystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile	Enkelvoudig steriel barrièresysteem
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website
	Importer	Importateur	Importeur	Importador	Importatore	Importeur
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM	Non compatibile con RM	MRI-onveilig
	Consult instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso	Raadpleeg de gebruiksaanwijzing
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Seguire le istruzioni per l'uso	Volg de gebruiksaanwijzing
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen
	Type B applied part	Pièce appliquée de type B	Anwendungsteil vom Typ B	Pieza aplicada tipo B	Parte applicata di tipo B	Toegepast onderdeel van type B
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union	Representante autorizado en la Comunidad Europea/Unión Europea	Rappresentante autorizzato nella Comunità Europea/Unione Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/Europese Unie
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)	Conformité Européenne (CE-markering)
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo	Unieke instrumentidentificatie
	Contains or presence of Phthalates: bis (2-ethylhexyl) phthalate (DEHP)	Traces ou présence de phthalates : phthalate de bis (2-éthylhexyle) (DEHP)	Enthält Phthalate oder kann Spuren von Phthalaten enthalten: Bis(2-ethylhexyl)phthalat (DEHP)	Contiene o hay presencia de ftalatos: Di(2-etilhexil) ftalato (DEHP)	Contiene o sono presenti ftalati: di (2-etilil) ftalato (DEHP)	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl)ftalaat (DEHP)
	Defibrillation Proof Type CF applied part	Partie appliquée de type CF résistante à la défibrillation	Defibrillationsgeschütztes Anwendungsteil vom Typ CF	Pieza aplicada de tipo CF a prueba de desfibrilación	Parte applicata di tipo CF a prova di defibrillazione	Defibrillatiebestendig toegepast onderdeel van type CF
	Contains hazardous substances	Contient des substances dangereuses	Enthält gefährliche Substanzen	Contiene sustancias peligrosas	Contiene sostanze pericolose	Bevat gevaarlijke stoffen

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen.

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Välisläbimõõt	Išorinis skersmuo	Ārējais diametrs	Dış çap	Внешний диаметр	Spoljašnji prečnik	Vanjski promjer
	Kasutatav pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebljiva dužina	Upotrebljiva dužina
	Sisești väikseim suurus	Minimalus įvediklio dydis	Minimālais ievadītāja izmērs	Minimum introduced boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	Minimalna veličina uvodnog instrumenta
	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	Količina
	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Ballooni mahutavus	Balionių talpa	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Sterilizēritud etiļeenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etiļēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilena-oksidadom	Sterilizirano etilen-oksidom
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Tootmiskuupäev	Pagaminimo data	Izgatašanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Sisaldab naturaalselt kummileaksite	Turi savo sudėtyje ir yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Suurus	Dydis	Izmērs	Boyut	Размер	Velicina	Velicina
	Säilitäda jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Palun lugege veebisaidlil olevaid kasutusjuhiseid	Skatit lietošanas instrukciju timekla vietne	Skatit lietošanas instrukciju timekla vietne	Web sitesindeki kullanim talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na web-sajtu	Pogledajte upute za upotrebu na web-stranici
	Importija	Importuotojas	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Ohtlik magnetresonantstomogra afias	MR nesaugus	Nedrikst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR	Nije sigurno kod pregleda MR-om
	Vaadake kasutusjuhendit	Žr. naudojimo instrukcijas	Skatit lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu	Pogledajte upute za upotrebu
	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis	Ievērojiet lietošanas instrukciju	Kullanım talimatlarını takip edin	Следуйте инструкциям по применению	Pridržavajte se uputstva za upotrebu	Pridržavajte se uputa za upotrebu
	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Mittepürogeenne	Nepirogeniškas	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	B-tüüpi kohaldatav osa	B tipo liečiamoji dalis	B tipa lietotā daļa	B Tipi uygulanan parça	Рабочая часть типа B	Primenjeni deo tipa B	Primjenjen dio tipa B
	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igalotasis atstovas Europos Bendrijoje / Europos Sąjungoje	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Evropskoj zajednici / Evropskoj uniji
	Conformitē Eiropēenne (CE-märk)	Conformitē Eiropēenne (CE žymė)	Conformitē Eiropēenne (CE marķējums)	Conformitē Eiropēenne (CE İşareti)	Conformitē Eiropēenne (маркировка CE)	Conformitē Eiropēenne (CE znak)	Conformitē Eiropēenne (oznaka CE)
	Meditsiiniseade	Medicinos priemonė	Medicīniska ierice	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Ärge kasutage, kui pakend on kahjustatud, ja lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas	Nelietot, ja iepakojums ir bojāts, un skatit lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje	Ievērojiet lietošanas instrukciju, kas pieejama timekla vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa web-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici
	Mudeli number	Modelio numeris	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius	Unikāls ierices identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Sisaldab ftalaate: bis(2-etiülheksüül)ftalaat (DEHP)	Sudėtyje yra ftalatų: bis(2-etilheksil) ftalatas (DEHP)	Vērojama ftalātu klātbūtne vai satur ftalātus: bis(2-etilheksil)ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Defibrillatsioonikindel CF-tüüpi rakendusosa	Defibriliavimo poveikiui atspari CF tipo darbinė dalis	Pret defibrilāciju noturīga CF tipa daļa, kas saskaras ar pacientu	Defibrilasyon Korumalı CF Tipi hastayla temas eden parça	Контактирующий с пациентом элемент типа CF с защитой от разряда дефибриллятора	Primenjeni deo tipa CF otporan na defibrilaciju	Primjenjen dio tipa CF otporan na defibrilaciju
	Sisaldab ohtlikke aineid	Sudėtyje yra pavojingų medžiagų	Satur bīstamas vielas	Tehlikeli maddeler içerir	Содержит вредные вещества	Sadrži opasne supstance	Sadrži opasne tvari

Märkus. Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba:** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezime:** Si izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

01/24
10005783004 A / DOC-0133324 A
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

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Edwards Lifesciences GmbH