Edwards Oximetry Central Venous Catheter/PediaSat Oximetry Catheter For Use with Edwards Oximetry Monitors or Compatible Bedside Modules

Directory

English (en)	1	Português (pt)	28	Română (ro)	52
Français (fr)	4	Česky (cs)	31	Eesti (et)	55
Deutsch (de)	7	Magyar (hu)	34	Lietuvių (lt)	58
Español (es)	10	Polski (pl)	37	Latviešu (lv)	61
Italiano (it)	13	Slovensky (sk)	40	Türkçe (tr)	64
Nederlands (nl)	16	Norsk (no)	43	Русский (ru)	67
Dansk (da)	19	Suomi (fi)	46	Srpski (sr)	70
Svenska (sv)	22	Български (bg)	49	Hrvatski (hr)	73
Ελληνικά (el)	25				
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike					76
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkleien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola					78

English

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

This product does not contain dry natural rubber.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

For Single Use Only

For the Edwards oximetry central venous catheter, please refer to Figure 1 on page 76 through Figure 2 on page 76 and Figure 5 on page 77 through Figure 8 on page 77.

For the PediaSat oximetry catheter, please refer to Figure 3 on page 76 through Figure 8 on page 77.

1.0 Description

The Edwards oximetry central venous catheter and PediaSat oximetry catheter are intended to provide the means for infusion of solutions, measuring pressure, and taking blood samples through the distal lumen (terminates at the catheter tip), the proximal lumen (terminates 7 cm proximal to the tip), and the medial lumen (terminates 5 cm proximal to the tip). The Edwards oximetry central venous catheter and PediaSat oximetry catheter are also intended to provide the means for continuously monitoring oxygen saturation measurement using an Edwards Lifesciences oximetry monitor or compatible bedside module.

These devices are for central venous access. Route of insertion can be internal jugular vein and sub-clavian vein.

Oxygen saturation is monitored by fiberoptic reflectance spectrophotometry. The amount of light absorbed, refracted, and reflected depends on the relative amounts of oxygenated and deoxygenated hemoglobin in the blood.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

When used with a compatible monitoring platform, oximetry catheters offer a means for pressure monitoring and continuously monitoring oxygen saturation measurement to help Intended Users assess cardiovascular function and guide treatment decisions. Intended Users include medical professionals trained in safe use of invasive hemodynamic technologies that are trained in clinical usage of central venous oximetry catheters.

The patient population for Edwards oximetry central venous catheter is adult critically ill or surgical patients.

The patient population for PediaSat oximetry catheter is adult and/or pediatric critically ill or surgical patients.

2.0 Indications

Edwards Oximetry Central Venous Catheters (EOCVC):

Edwards oximetry central venous catheters are indicated for hemodynamic monitoring through blood sampling, pressure monitoring, and oxygen saturation measurement when used with a compatible monitoring platform.

PediaSat Oximetry Catheter

PediaSat oximetry catheter is indicated for hemodynamic monitoring through blood sampling, pressure monitoring, and oxygen saturation measurements when used with a compatible monitoring platform in adults and/or pediatric patients.

3.0 Contraindications

Although there are no absolute contraindications to the use of the oximetry catheters, relative contraindications may include patients with recurrent sepsis or a hyper-coagulable state where the catheter could serve as a focus for septic or bland thrombus formation.

4.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RSTERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

5.0 Recommended Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (Type CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

WARNING: Do not modify or alter the product in any way. Alteration or modification may affect patient/operator safety or product performance.

- Oximetry catheters
- Edwards oximetry monitor or compatible bedside module
- OM2 optical module
- Sterile flush system
- Bedside ECG and pressure monitor system
- Pressure transducers and amplifiers

In addition, the following items should be immediately available: antiarrhythmic drugs, defibrillator, and respiratory assist equipment.

6.0 Monitor Set-up and Calibration for Oxygen Saturation Monitoring

Precaution: Do an *in vitro* calibration before preparing the catheter. Do not get the catheter tip or calibration cup wet before calibration.

Edwards oximetry monitors can be calibrated prior to catheter insertion by performing an *in vitro* calibration. When performing an *in vitro* calibration, do so before preparing the catheter (i.e. flushing the lumens). The catheter tip or calibration cup must not get wet before an *in vitro* calibration is performed. An *in vivo* calibration is required if an *in vitro* calibration is not done. *In vivo* calibration may be used to periodically recalibrate the monitor. Refer to the monitor operator's manual and the catheter's specification sheet for detailed calibration instructions.

Step	Procedure
1	Connect the optical module to the Edwards oximetry monitor or compatible bedside module.
2	Turn on monitor power switch.
3	Follow the calibration procedure for oximetry monitoring.

7.0 Suggested Insertion Procedure

Use sterile technique.

Do not use acetone or isopropyl alcohol on the catheter.

The following items should be immediately available: antiarrhythmic drugs, defibrillator, and respiratory assist equipment.

Any or all of these components may be included in your product configuration.

Step	Procedure
1	Prep the anticipated vessel puncture site as required.
2	Place fenestrated drape over the field.
3	Perform a skin wheal at the insertion site using a 25 gauge (0.5 mm) needle and 3 ml syringe. Used needles should be inserted into a sharps receptacle for disposal per hospital policy.
4	Flush catheter lumens with a sterile solution to assure patency and to avoid introduction of air into the circulation.

Step	Procedure
5	The oximetry catheters can be introduced by cutdown or percutaneous technique over a guidewire, with or without fluoroscopy. Continuous pressure monitoring is suggested during the insertion procedure. a) After antiseptic skin preparation and infiltration with local anesthetic, enter the vessel with locating needle and syringe. Upon aspiration of venous blood, remove the needle and syringe. Note: For EOCVC, suggested: 22 gauge [0.7 mm] Note: For PediaSat catheter, suggested: 21 gauge [0.8 mm] for 4.5F [1.50 mm] catheter, and 19 gauge [1.1 mm] for 5.5F [1.83 mm] catheter b) Attach 5 ml syringe to the catheter over needle. Insert the needle and relocate the vein. c) Upon aspiration of venous blood, remove the needle and syringe, leaving the catheter in place. d) Insert the desired tip (straight or "J") of a guidewire into the placed catheter or, if used, the thin-wall needle (For EOCVC, see Figure 1 on page 76, Figure 2 on page 76; for PediaSat, see Figure 3 on page 76, Figure 4 on page 76). Gentle manipulation may be needed to insert the guidewire. The guidewire should never be forced. If difficulty is met during insertion of the guidewire, completely withdraw the guidewire and reattempt insertion. CAUTION: To avoid possible severing of the guidewire, do not withdraw against the needle bevel. This can result in the need for an additional percutaneous intervention to retrieve the guidewire. CAUTION: Do not withdraw a PTFE-coated guidewire through a metal needle as this may damage the guidewire coating, potentially resulting in an embolism. e) Remove the catheter or thin-wall needle, leaving the guidewire in place. f) Enlarge the insertion site by threading a dilator over the guidewire (the puncture site may also be enlarged with a small scalpel). g) Leaving the guidewire in place, remove the dilator and thread the oximetry catheter over the guidewire. Important: Insertion depths will vary according to the insertion site and the size of the patient. h) Remove the guidewire and assure that blood can be freely aspirated through the distal lumen. For continuous infusion, attach infusion set luer connector to desired lumen hub and infuse or flush/cap per hospital protocol. Note: If placing the oximetry catheter into the internal jugular or subclavian vein, be sure to stop above the junction of the right atrium and the superior vena cava. WARNING: Positioning the distal tip of the catheter in the right atrium or ventricle is NOT recommended (see Complications: Cardiac Perforation).
6	Once in position after the guidewire has been removed, secure the catheter by suturing the integral suture wings to the skin.
7	If desired, the optional suture loop/box clamp can be placed on the catheter and sutured to the skin. a) Place the optional suture loop onto the catheter by spreading the suture loop wings and pressing onto the catheter (See Figure 5 on page 77). b) Snap the box clamp over the optional suture loop to secure both components to the catheter (See Figure 6 on page 77). c) Suture the optional suture loop and box clamp together to the patient to prevent catheter migration (See Figure 7 on page 77). Precaution: The box clamp must be removed from the catheter before attempting guidewire passage prior to catheter exchange to avoid guidewire damage.
8	If placing the catheter tip in the superior vena cava, verify catheter tip position in the superior vena cava by chest X-ray film immediately after insertion and periodically.

8.0 Maintenance and Use *in Situ*

Step	Procedure
1	Adequate maintenance is required to avoid catheter occlusion. Keep pressure monitoring and infusion lumens patent by intermittent flush, continuous, slow infusion with heparinized saline solution or use of a heparin lock (using individual injection caps in conjunction with heparinized saline solution). Precaution: Always disinfect injection caps before entry with syringe needle (see Complications: Sepsis/Infection). Use only a small bore needle (22 gauge (0.7 mm) or smaller) to puncture and inject through the injection caps. To use Interlink Injection sites: a) Ensure that the injection sites are securely connected to the lumen hubs. b) Grasp finger flange to stabilize injection site (See Figure 8 on page 77). c) Swab septum with preferred antiseptic. d) Insert Interlink cannula, attached to an appropriate device, directly through the center of the septum. CAUTION: If a conventional needle must be used, insert a small gauge needle into the perimeter of the septum to avoid fluid leakage while the needle is in place. e) Engage locking features if applicable.
2	For blood sampling, attach blood sampling device to desired lumen hub and draw blood sample per hospital protocol.
3	Periodically check lines and transducers for air bubbles. Ensure that the connecting lines and stopcocks remain tightly fitted.
4	The catheter should remain indwelling as determined by hospital protocol. WARNING: The incidence of complications increases significantly with indwelling periods greater than 72 hours. Prophylactic anticoagulation and infection control measures should be considered with the use of central venous catheters.

9.0 MRI Information

MR	MR Safe
-----------	---------

The oximetry catheters are made from nonmetallic, non-conducting, and nonmagnetic materials. Therefore, the EOCVC and PediaSat catheters are MR-safe, which are items that pose no known hazards in all MR environments.

Precaution: The cables which connect the oximetry catheters to monitors do contain metals and must be disconnected prior to performing the MRI procedure. Failure to do so may cause unintentional removal of catheter from patient.

10.0 Complications

Possible complications include:

10.1 Thrombosis

Thrombi have been shown to form on the surface of catheters after their insertion into the central circulation. Complications associated with thrombosis may include pulmonary emboli and infarction, and septic phlebitis.

10.2 Sepsis/Infection

Positive catheter-tip cultures resulting from contamination and colonization have been reported, as well as incidences of septic and aseptic vegetation in the right heart. Increased risks of septicemia and bacteremia have been associated with blood sampling, infusion of fluids, and catheter-related thrombosis. Preventive measures should be taken to guard against infection (e.g., use of sterile technique, application of topical antibiotic ointment, changing of sterile dressings as indicated by institutional policy, and disinfecting the injection caps before entry with syringe needle) as well as the frequent assessment of the continued need for hemodynamic monitoring.

10.3 Cardiac Perforation

Atrial perforation and subsequent pericardial tamponade have been reported. Preventive measures should include verification of catheter tip position by chest X-ray film and noting insertion depth immediately following insertion. Ideally, the catheter tip should be positioned parallel to the vessel wall and no farther than the junction of the superior vena cava and right atrium.

10.4 Vessel Perforation

Venous perforation and necrosis of the wall of the vein which can lead to perforation, due to a malpositioned catheter, have been reported. Preventive measures should include verification of the catheter tip position by chest X-ray film, noting insertion depth immediately following insertion. Ideally, the catheter tip should be positioned parallel to the vessel wall and no farther than the junction of the superior vena cava and right atrium.

WARNING: If there is any doubt that the catheter tip may not be intravascular, further steps should be taken to identify the exact location of the catheter tip, see Complications for Cardiac Perforation and Vessel Perforation.

10.5 Other Complications

Pressure catheters have also been associated with pneumothorax, embolism, catheter embolism, hemomediastinum/hydropneumothorax, arrhythmia, bleeding, hematoma, tissue damage, vessel perforation, hemothorax, chylothorax, adverse reaction to material, nerve injury, and delay in treatment leading to hemodynamic instability in the event of device malfunction.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance (SSCP). After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

11.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use. Do not resterilize.

The packaging is designed to protect the catheter from damage. Handling outside the package may result in crushing of the catheter body. It is therefore recommended that the catheter remain inside the package until use.

12.0 Storage

Store in a cool, dry place.

13.0 Operating Conditions

Intended to operate under physiological conditions of the human body.

14.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in catheter deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

15.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
 (24 hours): ... 800.822.9837
 Outside the U.S. and Canada
 (24 hours): ... 949.250.2222
 In the UK: ... 0870 606 2040 - Option 4
 In Ireland: ... 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

16.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices subject to change without notice.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

17.0 Specifications

	Oximetry Central Venous Catheter X3816*** 8.5F, 3 Lumen, 16 cm**	Oximetry Central Venous Catheter X3820*** 8.5F, 3 Lumen, 20 cm	PediaSat Oximetry Catheter XT245*** 4.5F (1.50 mm), 5 cm	PediaSat Oximetry Catheter XT248*** 4.5F (1.50 mm), 8 cm	PediaSat Oximetry Catheter XT358*** 5.5F (1.83 mm), 8 cm	PediaSat Oximetry Catheter XT3515*** 5.5F (1.83 mm), 15 cm
Usable Length (cm)	16	20	5	8	8	15
Catheter Body French Size	8.5F (2.8 mm)	8.5F (2.8 mm)	4.5F (1.50 mm)	4.5F (1.50 mm)	5.5F (1.83 mm)	5.5F (1.83 mm)
Recommended Dilator Size	10.5F (3.5 mm)	10.5F (3.5 mm)	5.5F (1.83 mm)	5.5F (1.83 mm)	6.5F (2.2 mm)	6.5F (2.2 mm)
Distance from tip (cm)						
Proximal Port	7	7	1	1	2	2
Medial Port	5	5	NA	NA	1	1
Lumen Volume						
Distal Lumen	0.6 ml	0.7 ml	0.2 cc	0.3 cc	0.2 cc	0.3 cc
Medial Lumen	0.3 ml	0.3 ml	NA	NA	0.2 cc	0.2 cc
Proximal Lumen	0.3 ml	0.3 ml	0.1 cc	0.2 cc	0.2 cc	0.2 cc
Infusion Rate* (ml/min)						
Normal Saline						
Distal Lumen	75	76	16	11	29	21
Medial Lumen	29	26	NA	NA	12	7
Proximal Lumen	31	28	13	8	13	8
Maximum Guidewire Diameter	0.032" (0.81 mm)	0.032" (0.81 mm)	0.018" (0.46 mm)	0.018" (0.46 mm)	0.025" (0.64 mm)	0.025" (0.64 mm)
Distal Lumen						

All oximetry catheters, except those with black oximetry connectors, are compatible with Edwards optical modules. Models ending with the letter "P" are also compatible with Philips optical modules. The Catheter Factor required for *in vivo* calibration with Philips monitors is located on the top of the optical connector.

Catheters are constructed from multi-lumen proprietary polymer blend extrusions.

All specifications given are nominal values.

*Using normal saline at room temperature, 1 meter above insertion site, gravity drip.

**The 8.5F, 3 Lumen, 16 cm catheter is not available on the European market.

***Model number suffix, other than "P" as defined above, denotes packaging and component configuration. Contact Technical Assistance for additional model number information.

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







	English	Français	Deutsch	Español	Italiano
	Number of lumens	Nombre de lumières	Anzahl der Lumen	Número de luces	Numero di lumi
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Single use	Ne pas réutiliser	Zum einmaligen Gebrauch	Un solo uso	Monouso
	Quantity	Quantité	Menge	Cantidad	Quantità
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union	Representante autorizado en la Comunidad Europea/Unión Europea	Rappresentante autorizzato nella Comunità Europea/Unione Europea
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)
	MR Safe	Aucun risque en milieu RM	MR-sicher	Seguro para RM	Compatibile con RM
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung lesen	Siga las instrucciones de uso	Seguire le istruzioni per l'uso
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
	Type CF applied part	Pièce appliquée de type CF	Anwendungsteil vom Typ CF	Pieza aplicada tipo CF	Parte applicata di tipo CF
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione

	English	Français	Deutsch	Español	Italiano
	Importer	Importateur	Importeur	Importador	Importatore
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Double sterile barrier system	Système de barrière stérile double	Doppeltes Sterilbarrieresystem	Un sistema doble de barrera estéril	Sistema a doppia barriera sterile
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarrieresystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web
	Not made with natural rubber latex	Non fabriqué à partir de latex de caoutchouc naturel	Nicht aus Latex hergestellt	No contiene látex de caucho natural.	Non realizzato in lattice di gomma naturale

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Aantal lumina	Antal lumens	Antal lumen	Αριθμός αυλών	Números de lúmenes
	Buitendiameter	Udvendig diameter	Ytterdiameter	Εξωτερική διάμετρος	Diâmetro exterior
	Bruikbare lengte	Anvendelig længde	Brukbar längd	Ωφέλιμο μήκος	Comprimento útil
	Aanbevolen maat voerdraad	Anbefalet guidewirestørrelse	Rekommenderad ledarstorlek	Συνιστώμενο μέγεθος οδηγού σύρματος	Tamanho recomendado do fio-guia
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Lotnummer	Partnummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Vervaldatum	Anvendes inden	Sista förbrukningsdag	Ημερομηνία λήξης	Data de vencimento
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθYLEνίου	Esterilizado com óxido de etileno
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντι- πρόσωπος στην Ευρωπαϊκή Κοινότητα/Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
	MRI-veilig	MR-sikker	MR-säker	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)	Utilização segura em ambiente de RM
	Volg de gebruiksaanwijzing	Følg brugsanvisningen	Följ bruksanvisningen	Ακολουθείτε τις Οδηγίες χρήσης	Seguir as instruções de utilização
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogénico
	Toegepast onderdeel van type CF	Anvendt del af typen CF	Patientansluten del av CF-typ	Εφαρμοζόμενο εξάρτημα τύπου CF	Peça aplicada de tipo CF
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Op een koele en droge plaats bewaren place	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Dubbel steriel barrièresysteem	Dobbelt sterilt barriersystem	Dubbelt sterilt barriärsystem	Σύστημα διπλού αποστεί-ρωμένου φραγμού	Sistema de barreira dupla esterilizada
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barriersystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστεί-ρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωρι-στικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden	Följ bruksanvisningen på webbplatsen	Ακολουθείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Seguir as instruções de utilização no site
	Bevat geen natuurlijk rubberlatex	Ikke fremstillet med naturligt gummilatex	Innehåller inte naturgummilatex	Δεν κατασκευάζεται από φυσικό λάτεξ	Não fabricado com látex de borracha natural

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	Suomi	Български	Română	Eesti	Lietuvių
	Luumenien määrä	Брой лумени	Număr de Lumenuri	Valendike arv	Spindžių skaičius
	Ulkoalkaisija	Външен диаметър	Diametru exterior	Välisläbimõõt	Išorinis skersmuo
	Käyttöpituus	Използваема дължина	Lungimea utilă	Kasutatav pikkus	Naudingasis ilgis
	Ohjainlangan suositeltu koko	Препоръчителен размер телен водач	Mărimea recomandată a firului de ghidaj	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis
	Mallinumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks.	Vienkartinio naudojimo
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de	Kõlblik kuni	Naudoti iki nurodytos datos
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен предста- вител в Европейската общност/Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Sopii magneettikuvaukseen	Безопасно при MR	Sigur în utilizarea la RM	Ohutu magnetresonantsomograafias	MR saugus
	Noudata käyttöohjeita	Спазвайте инструкциите за употреба	Respectați instrucțiunile de utilizare	Järgige kasutusjuhendit	Vadovautis naudojimo instrukcijomis
	Älä steriloï uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuoti
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	CF-tyypin liityntäosa	Приложна част тип CF	Componentă aplicată de tip CF	CF-tüüpi kohaldatav osa	CF tipo liečiamoji dalis
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Основные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Kaksinkertainen steriili sulkujärjestelmä	Двойна стерилна бариерна система	Sistem cu barieră sterilă dublă	Kahekordne steriilne kattesüsteem	Dvigubo sterilaus barjero sistema
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalusis priemonės identifikatorius
	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje
	Ei sisällä luonnonkumilateksia	Не е произведено с естествен каучуков латекс.	Nu este fabricat cu latex din cauciuc natural	Pole valmistatud looduslikust kummilateksist.	Pagaminta be natūralios gumos latekso

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Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

10/24
10009300003 B / DOC-0133326B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU