Edwards

Fogarty Surgical Clamp Inserts

Directory	
English (en).....	1
Français (fr).....	3
Deutsch (de).....	5
Español (es).....	7
Italiano (it).....	9
Nederlands (nl).....	11
Dansk (da).....	13
Svenska (sv).....	15
Ελληνικά (el).....	17
Português (pt).....	19
Česky (cs).....	21
Magyar (hu).....	23
Polski (pl).....	25
Slovensky (sk).....	27
Norsk (no).....	29
Suomi (fi).....	31
Български (bg).....	33
Română (ro).....	35
Eesti (et).....	37
Lietuvių (lt).....	39
Latviešu (lv).....	41
Türkçe (tr).....	43
Русский (ru).....	45
Srpski (sr).....	47
Hrvatski (hr).....	49
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázás ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sēmbol Aizkļamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	51

English

Fogarty Surgical Clamp Inserts

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.


Figure 1

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

Fogarty clamp inserts are intended to be used by medical professionals trained in safe use of invasive hemodynamic technologies who have been trained in clinical usage of temporary atraumatic vascular occlusion technologies as part of their respective institutional guidelines. A unique grasping principle accomplished through the use of replaceable jaw inserts can be incorporated into compatible clamps of the appropriate length to evenly distribute applied pressure while providing adequate traction.

For a combination of atraumatic traction, the EverGrip insert, a compliant insert with a traction pattern, may be used.

Clinical benefits of clamp inserts include providing control of blood flow with minimization of trauma to the tissue.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the

safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

1.1 Replaceable Inserts

1. Multipurpose.
2. Easy to replace.
3. Adaptable to specific procedures.
4. Radiopaque.

1.2 Design Features

Compliant inserts are available in four models: Hydrjaw, Softjaw, Safejaw and EverGrip inserts.

The Hydrjaw insert is designed to be very compliant and is composed of soft synthetic rubber.

The Softjaw insert is designed to be moderately compliant and is composed of soft synthetic rubber.

The Safejaw and EverGrip inserts are constructed of the same material as the Softjaw insert with the added feature of a traction pattern. Soft pyramids of synthetic rubber gently grip the vessel and conform to its profile to provide secure occlusion.

The Hydrjaw, Softjaw, Safejaw and EverGrip inserts contain a metal spine for attachment to a clamp; as a result, they are radiopaque and can be visualized by fluoroscopy.

The Tractionjaw insert is made of rigid polymer which provides both vertical and horizontal security. It has two rows of compound angled teeth; each tooth is uniform with precisely controlled angles, depth, and contour for maximum traction with minimal penetration. Barium sulfate is dispersed within the polymer to make it radiopaque.

All inserts of a given size are interchangeable; hence, they may be used in combinations other than described. For example, for procedures requiring maximum traction at the sacrifice of some gentleness, two traction inserts may be used.

When extreme gentleness is required, a combination of two compliant inserts may be used. This combination also permits the passage of large intraluminal instruments. The use of two compliant inserts on any of the available clamps produces a noncrushing clamp.

See Specification Table for model information.

2.0 Intended Use/Purpose

The Fogarty clamp inserts provide a means of occlusion pressure to control the flow of circulating blood. The Fogarty clamp inserts are intended for atraumatic vessel occlusion.

3.0 Indications

These devices are indicated for vessel occlusion in adult critically ill or surgical patient populations.

Vessel Occlusion

1. Suitable for both veins and arteries.
2. Provides occlusion of atherosclerotic vessels without excessive clamping forces.
3. Minimizes intimal damage and fragmentation of atherosclerotic material.
4. Performs efficiently on either diseased or normal vessels.

CAUTION: Care should be taken to ensure a firm grasp of the vessel when two compliant inserts are used, since traction is reduced and the possibility of slippage is increased.

4.0 Contraindications

There are no known contraindications.

5.0 Complications

This device has been associated with blood loss or bleeding, vascular damage, and adverse reaction to device materials.

6.0 Precaution

It is possible to inadvertently pass a needle through the compliant inserts, thereby attaching the insert to the suture line. If the Hydrjaw insert is damaged it might cause silicone to leak out and reduce the clamping force, so it should be replaced.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

7.0 Instructions

7.1 To Place Inserts in Clamp

Step	Procedure
1	Open the clamp and hold the handle very firmly. Place the pin located on the wide end of the insert in the mounting hole closest to the box lock.
2	Snap the other pin in its mounting hole by pushing forward and down on the insert with sufficient pressure to secure it. Do not attempt to push down directly on this pin; instead, slide it forward with the thumb until it seats in place (see Figure 1).

7.2 To Remove Inserts from Clamp

Step	Procedure
1	Open the clamp and hold by the handle to which the insert is attached.
2	Hold the insert to be removed between the thumb and forefinger of the right hand.
3	Pull outward. The insert will snap out. Do <u>not</u> try to pry off the insert from the middle between the two pins, since this will cause it to bend or break.

8.0 MRI Information

This product has not been tested for MRI compatibility.

14.0 Specification Table

Models	Product Description	Insert Length (mm)
HYDRA33	Hydrajaw/Tractionjaw	33
HYDRA61	Hydrajaw/Tractionjaw	61
HYDRA86	Hydrajaw/Tractionjaw	86
SOFT33	Softjaw/Tractionjaw	33
SOFT61	Softjaw/Tractionjaw	61
SOFT86	Softjaw/Tractionjaw	86
SAFE33	Safejaw/Softjaw	33
SAFE61	Safejaw/Softjaw	61
SAFE86	Safejaw/Softjaw	86
EVERGRIP33	EverGrip/EverGrip	33
EVERGRIP61	EverGrip/EverGrip	61
EVERGRIP86	EverGrip/EverGrip	86
DSAFE33	Safejaw/Safejaw	33
DSAFE61	Safejaw/Safejaw	61
DSAFE86	Safejaw/Safejaw	86
DSOFT33	Softjaw/Softjaw	33
DSOFT61	Softjaw/Softjaw	61
DSOFT86	Softjaw/Softjaw	86
DHYDRA33	Hydrajaw/Hydrajaw	33
DHYDRA61	Hydrajaw/Hydrajaw	61
DHYDRA86	Hydrajaw/Hydrajaw	86
SAFEHY33	Safejaw/Hydrajaw	33
SAFEHY61	Safejaw/Hydrajaw	61
SAFEHY86	Safejaw/Hydrajaw	86

9.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is as marked on each package.

Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

For toll free numbers, visit
<https://eifu.edwards.com/eifu/pages/eifu-contact>.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

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STERILE EO

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Modele	Descrierea produsului	Lungimea inserției (mm)
HYDRA33	Hydrajaw/Tractionjaw	33
HYDRA61	Hydrajaw/Tractionjaw	61
HYDRA86	Hydrajaw/Tractionjaw	86
SOFT33	Softjaw/Tractionjaw	33
SOFT61	Softjaw/Tractionjaw	61
SOFT86	Softjaw/Tractionjaw	86
SAFE33	Safejaw/Softjaw	33
SAFE61	Safejaw/Softjaw	61
SAFE86	Safejaw/Softjaw	86
EVERGRIP33	EverGrip/EverGrip	33
EVERGRIP61	EverGrip/EverGrip	61
EVERGRIP86	EverGrip/EverGrip	86
DSAFE33	Safejaw/Safejaw	33
DSAFE61	Safejaw/Safejaw	61
DSAFE86	Safejaw/Safejaw	86
DSOFT33	Softjaw/Softjaw	33
DSOFT61	Softjaw/Softjaw	61
DSOFT86	Softjaw/Softjaw	86
DHYDRA33	Hydrajaw/Hydrajaw	33
DHYDRA61	Hydrajaw/Hydrajaw	61
DHYDRA86	Hydrajaw/Hydrajaw	86
SAFEHY33	Safejaw/Hydrajaw	33
SAFEHY61	Safejaw/Hydrajaw	61
SAFEHY86	Safejaw/Hydrajaw	86

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STERILE EO

14.0 Tabela sa specifikacijama

Modeli	Opis proizvoda	Dužina uloška (mm)
HYDRA33	Hydrajaw/Tractionjaw	33
HYDRA61	Hydrajaw/Tractionjaw	61
HYDRA86	Hydrajaw/Tractionjaw	86
SOFT33	Softjaw/Tractionjaw	33
SOFT61	Softjaw/Tractionjaw	61
SOFT86	Softjaw/Tractionjaw	86
SAFE33	Safejaw/Softjaw	33
SAFE61	Safejaw/Softjaw	61
SAFE86	Safejaw/Softjaw	86
EVERGRIP33	EverGrip/EverGrip	33
EVERGRIP61	EverGrip/EverGrip	61
EVERGRIP86	EverGrip/EverGrip	86
DSAFE33	Safejaw/Safejaw	33
DSAFE61	Safejaw/Safejaw	61
DSAFE86	Safejaw/Safejaw	86
DSOFT33	Softjaw/Softjaw	33
DSOFT61	Softjaw/Softjaw	61
DSOFT86	Softjaw/Softjaw	86
DHYDRA33	Hydrajaw/Hydrajaw	33
DHYDRA61	Hydrajaw/Hydrajaw	61
DHYDRA86	Hydrajaw/Hydrajaw	86
SAFEHY33	Safejaw/Hydrajaw	33
SAFEHY61	Safejaw/Hydrajaw	61
SAFEHY86	Safejaw/Hydrajaw	86

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
	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Catalogusnummer	Katalognummer	Katalognummer	Αριθμός καταλόγου	Número de catálogo
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Te gebruiken tot	Anvendes inden	Använd före	Ημερομηνία λήξης	Validade
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (Σήμανση CE)	Conformité Européenne (Marcação CE)
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Let op: Federal (USA) law restricts this device to sale by or on the order of a physician.	Forsigtig: Federal (USA) law restricts this device to sale by or on the order of a physician.	Var försiktig: Federal (USA) law restricts this device to sale by or on the order of a physician.	Προσοχή: Federal (USA) law restricts this device to sale by or on the order of a physician.	Aviso: Federal (USA) law restricts this device to sale by or on the order of a physician.

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Bevat geen natuurlijk rubberlatex	Ikke fremstillet med naturligt gummlatex	Innehåller inte naturgummlatex	Δεν κατασκευάζεται από φυσικό λάτεξ	Não fabricado com látex de borracha natural
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrónicas
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Μοναδικό αναγνωριστικό συσκευής	Identificador único de dispositivo
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogénico
	Bruikbare lengte van klemstukken in mm	Anvendelig længde af klemmeindføringsanordninger i mm	Brukbar längd av kläminlägg i mm	Ωφέλιμο μήκος ένθετων σφιγκτήρα σε mm	Comprimento útil dos mordentes de grampo em mm
	Openingswijdte van klemveer in mm.	Størrelse af fjederclipsens kæbeåbning i mm.	Storlek i mm på fjäderklämmans öppna käft.	Μέγεθος ανοίγματος σιαγόνας αγκτήρα με ελατήριο σε mm.	Tamanho de abertura da garra do clipe de mola em mm.
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico

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



	Suomi	Български	Română	Eesti	Lietuvių
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Käytettävä ennen	Срок на годност	Valabil până la	Kõlblik kuni	Naudoti iki
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcaj CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Älä käyttää, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Steriloitu etyleenioksidiilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Tärkeä huomautus: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Atenție: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Perspėjimas: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Ei sisällä luonnonkumilateksia	Не е произведено с естествен каучуков латекс.	Nu este fabricat cu latex din cauciuc natural	Pole valmistatud looduslikust kummilateksist.	Pagaminta be natūralios gumos latekso
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Katso käyttöohjeet verkkosivustolta + 1 888 570 4016	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Älä steriloitu uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotina nesterilizuokite
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sümbole seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Suomi	Български	Română	Eesti	Lietuvių
	Puristinlisäosien käyttöpietus (mm)	Използваема дължина на вложките за клампа в mm	Lungimea utilă a inserțiilor pensei clamp, exprimată în mm	Klambri vahetükkide kasutatav pikkus (mm)	Įterpiamų spaustukų naudingasis ilgis (mm)
	Kiinnityspihtien leukojen aukon koko (mm).	Размер на отваряне на челюстите на пружинната скоба в mm.	Mărimea deschiderii brațelor clipului de tensiune, exprimată în mm.	Vedruklemmiga haaratsi ava suurus (mm)	Spyruoklinio spaustuvo griebtuvo angos dydis (mm).
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Kataloga numurs	Katalog numarası	Номер по каталогу	Kataloški broj	Kataloški broj
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
	Daudzums	Miktar	Количество	Količina	Količina
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Upotrebljivo do	Upotrijebiti do
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Izgatavošanā nav izmantots dabiskā kaučuka latekss	Doğal kauçuk lateksten yapılmamıştır	Не содержит натурального каучукового латекса	Ne sadrži lateks od prirodne gume	Nije izrađeno od prirodnog gumenog lateksa
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internetskoj stranici


Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarazás ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символах ■ Legenda de simboluri ■ Sümböölite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Spaiļu ieliktņu izmantojamais garums (mm)	Klemp insertlerinin kullanılabilir uzunluğu (mm)	Рабочая длина вставок зажима в мм	Upotrebljiva dužina uložka kleme u mm	Upotrebljiva duljina umetaka hvataljke u mm
	Atsperspauļu atvēruma izmērs (mm).	Yaylı klips çene açıklığı boyutu (mm).	Раскрытие браншей пружинного зажима в мм.	Veličina otvora kraka štipaljke sa oprugom u mm.	Veličina otvora čeljusti opružne stezaljke u mm.
	Medicīniska ierīce	Tibbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Stítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Előkelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykietie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notā:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Ši izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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06/24
159439007 C / DOC-133371 C
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH