Edwards

Fogarty Corkscrew Catheter

Directory	
English (en).....	1
Français (fr).....	4
Deutsch (de).....	6
Español (es).....	8
Italiano (it).....	10
Nederlands (nl).....	12
Dansk (da).....	14
Svenska (sv).....	16
Ελληνικά (el).....	18
Português (pt).....	20
Česky (cs).....	22
Magyar (hu).....	24
Polski (pl).....	26
Slovensky (sk).....	28
Norsk (no).....	30
Suomi (fi).....	32
Български (bg).....	34
Română (ro).....	36
Eesti (et).....	38
Lietuvių (lt).....	40
Latviešu (lv).....	42
Türkçe (tr).....	44
Русский (ru).....	46
Srpski (sr).....	48
Hrvatski (hr).....	50
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike.....	52
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola.....	53

English

Fogarty Corkscrew Catheter

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The device is a flexible catheter consisting of a polyvinylchloride outer catheter body surrounding a movable, radiopaque stainless steel spring body. A second radiopaque, spiral-shaped stainless steel cable covered by a latex membrane is attached to the distal end of the catheter. A handle to modulate the spiral cable diameter is attached to the inner spring body at the proximal end of the catheter. The spiral cable design of the Fogarty Corkscrew catheter increases the surface area for entrapping fibrous material.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

This device aids in clot management and is used for restoration of circulatory blood flow to distal arterial system and downstream tissues.

The device is to be used by medical professionals who are trained in the safe use and clinical usage of endovascular catheter technologies as part of their respective institutional guidelines.

2.0 Intended Use/Purpose

The Fogarty Corkscrew catheter is intended to be used for the removal of emboli and thrombi from either native vessels or synthetic grafts in the arterial system.

3.0 Indications

The Fogarty Corkscrew catheter is indicated for use in adult patients with adherent emboli/thrombi in the peripheral arterial vasculature in either native vessels or synthetic grafts.

4.0 Contraindications

- The Fogarty Corkscrew catheter must not be used outside the arterial system.
- The Fogarty Corkscrew catheter is not recommended for endarterectomy procedures.
- The Fogarty Corkscrew catheter is not intended for use as a vessel dilator.
- The Fogarty Corkscrew catheter must not be used in the venous system.

5.0 Warnings

- **Exposure to the atmosphere, handling during insertion, and plaque and other deposits within the blood vessel may cause latex membrane degradation.**
- **To minimize the risk of vessel or membrane damage, do not exceed the maximum recommended pull force (see Catheter Specifications).**

6.0 Instructions

6.1 Technique

Step	Procedure
1	Inspect the catheter.
2	With the membrane extended, with the control button in the locked position, place the catheter in the artery or graft so that the membrane is past the clot to be removed.
3	Contract the membrane by pushing down the control button and pulling it backwards. This causes the membrane to expand into a helical configuration so that it grips or engages the clot.
4	Withdraw the entire catheter from the vessel while the membrane maintains slight contact with luminal wall, removing the occluding material as the catheter is withdrawn. Note: A sensitive touch in controlling the contracted membrane during catheter withdrawal is important so that the membrane diameter can be adjusted to varying luminal diameters. Contraction of the membrane is associated with a feeling of resistance. Note: Because of the complexity and variation in the procedure, the choice of surgical technique, appropriately modified in accordance with the previously described warnings, precautions, and technique, is left to the discretion of the individual surgeon.

7.0 MRI Information

This product has not been tested for MRI compatibility. This device is not intended to be used in an MRI environment and should not be used in this environment as it contains metal.

8.0 Complications

Potential complications include, but are not limited to, vascular damage, including perforation, vessel rupture, dissection and intimal and/or membrane disruption, bleeding/hemorrhage, local hematomas, distal embolization, arterial embolization, thrombosis, air embolization, local or systemic infection, adverse reaction to device materials, arterial spasm, device breakage requiring intervention to remove, aneurysms and arteriovenous fistula formation.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

9.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 800.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

14.0 Specification Table

Model	140806	140808	1408010
French Size	4F	5F	6F
(mm)	1.33	1.67	2.0
Usable Length (cm)	80	80	80
Maximum Diameter of Extended Membrane (mm)	2.7	3.1	3.3
Diameter of Contracted Membrane (mm)	6	8	10
Maximum Recommended Pull Force on Contracted Membrane (lbs.)	0.9 (0.4 kg)	1.3 (0.6 kg)	1.8 (0.8 kg)

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© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.