Directory

English (en)	1	Português (pt)	23	Română (ro)	41
Français (fr)	3	Česky (cs)	25	Eesti (et)	43
Deutsch (de)	6	Magyar (hu)	27	Lietuvių (lt)	46
Español (es)	8	Polski (pl)	29	Latviešu (lv)	48
Italiano (it)	11	Slovensky (sk)	32	Türkçe (tr)	50
Nederlands (nl)	13	Norsk (no)	34	Русский (ru)	52
Dansk (da)	15	Suomi (fi)	36	Srpski (sr)	55
Svenska (sv)	18	Български (bg)	38	Hrvatski (hr)	57
Ελληνικά (el)	20				
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike 60					
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de simbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola 61					

English

VAMP Plus Closed Blood Sampling System

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use

For Single Use Only

For figures, refer to Figure 1 on page 60 through Figure 3 on page 60.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Edwards Lifesciences VAMP Plus closed blood sampling system provides a safe and convenient method for the withdrawal of blood samples from pressure monitoring lines. The blood sampling system is designed for use with disposable and reusable pressure transducers and for connection to central line catheters, venous, and arterial catheters where the system can be flushed clear after sampling. The VAMP Plus closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use of blood sampling technologies as part of their respective institutional guidelines.

Edwards, Edwards Lifesciences, the stylized E logo, TruWave, VAMP, and VAMP Plus are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

3.0 Indications for Use

For adult patients with medical conditions requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

4.0 Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

5.0 Warnings

- **This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.**
- Some models may contain phthalates, specifically DEHP [Bis (2-ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.
- Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

6.0 Instructions for Use

CAUTION: The use of lipids with the VAMP Plus closed blood sampling system may compromise product integrity.

6.1 Equipment

- Flush device or flow controlling device (maximum 4 mL/hr flow rate)
- Disposable or reusable pressure transducer, if desired
- VAMP Plus closed blood sampling system

6.2 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP Plus kit from the sterile package.
2	If the VAMP Plus kit does not include a preconnected TruWave disposable pressure transducer, remove the protective caps and attach the distal end with the female luer-lock connector to the male luer of a transducer, or other fluid flush device (see Figure 1 on page 60).
3	All connections should be secure. Note: Wet connections promote over tightening by lubricating the fittings. Over tightened connections may result in cracks or leaks.
4	Push the reservoir plunger to the closed and locked position by pressing on the reservoir extension until it latches in the closed position (see Figure 3 on page 60).
5	Connect the kit to an IV fluid source without pressurizing the fluid. Begin to gravity fill the kit first through the transducer and out through the vent port according to its manufacturer's instructions.
6	Replace all vented caps on the sideports of the stopcocks with non-vented caps.
7	To fill the VAMP Plus kit, ensure that the shut off valve is in the open position indicated by the handle being parallel to the tubing. Orient the kit such that the reservoir inlet and outlet ports are in an upward vertical position and the sampling sites are above the reservoir at approximately 45°. Slowly fill and de-bubble the reservoir and each sample site in succession. CAUTION: Remove all air bubbles to reduce the risk of air emboli and reduce loss of pressure signal.

Step	Procedure
8	The transducer and reservoir with bracket can be mounted to an IV pole using the Edwards disposal pressure transducer holder (see Figure 1 on page 60). For convenience, the sample site nearest the reservoir can be easily mounted to the holder.
9	Pressurize the IV solution bag per hospital policy and transducer manufacturer's instructions. Flow rate will vary with pressure gradient across the flush device.
10	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter.
11	Zero the transducer according to its manufacturer's instructions.

6.3 Clearing IV Fluid in Preparation for Drawing Blood Samples

Two methods may be used for drawing the clearing volume into the VAMP Plus reservoir.

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

6.4 One-Hand Pole Mount Method to Draw Clearing Volume

Step	Procedure
1	With the VAMP Plus reservoir bracket mounted in the Edwards disposable pressure transducer holder, release the plunger latch. Fill the reservoir by drawing the plunger to the open position by pressing together on the plunger and the bracket extension finger hold (see Figure 2 on page 60) until the plunger stops and the reservoir has reached its 12 mL volume capacity. Note: If difficulties are experienced in drawing the clearing volume, check the catheter for possible occlusions or restrictions. Note: Recommended time to fill the reservoir is approximately 1 second for each mL of clearing volume.
2	Once the clearing volume has been drawn, close the reservoir shut off valve by turning the handle perpendicular to the tubing. This will ensure that the sample drawn is from the patient and not the reservoir.
3	Draw the sample according to one of the methods outlined in the Drawing Blood Samples section below.
4	Once the sample has been drawn, open the shut off valve by turning the handle 90° so that it is parallel to the tubing.
5	Reinfuse the clearing volume back to the patient by returning the reservoir plunger to the closed position. Press together the reservoir extension finger hold and the plunger pad (see Figure 3 on page 60) until the plunger stops and firmly latches in the closed position. Note: If difficulties are experienced in reinfusing the clearing volume, check the catheter for possible occlusions or restrictions.

Step	Procedure
	Note: Recommended time to push the reservoir plunger to the fully closed position is approximately 1 second for each mL of clearing volume. CAUTION: The clearing sample should not remain in the reservoir for longer than 2 minutes.
6	Flush the reservoir, sample sites, and line clear of residual blood. Swab the sampling sites, ensuring removal of any excess blood. Note: Do not use acetone.

6.5 Hand Held Syringe Method to Draw Clearing Volume

Step	Procedure
1	Remove the VAMP Plus reservoir from the bracket. Hold the reservoir in one hand and release the plunger latch. Fill the reservoir by drawing the plunger to the open position by either pushing the plunger open with your thumb, or pulling with your other hand, until the plunger stops and the reservoir has reached its 12 mL capacity. Note: If difficulties are experienced in drawing the clearing volume, check the catheter for possible occlusions or restrictions. Note: Recommended time to fill the reservoir is approximately 1 second for each mL of clearing volume.
2	Once the clearing volume has been drawn, close the reservoir shut off valve by turning the handle perpendicular to the tubing. This will ensure that the sample drawn is from the patient and not the reservoir.
3	Draw the sample according to one of the methods outlined in the Drawing Blood Samples section below.
4	Once the sample has been drawn, open the shut off valve by turning the handle 90° so that it is parallel to the tubing.
5	Reinfuse the clearing volume back to the patient by pushing the plunger until it stops and firmly latches in the closed position. Note: If difficulties are experienced in reinfusing the clearing volume, check the catheter for possible occlusions or restrictions. Note: Recommended time to push the reservoir plunger to the fully closed position is approximately 1 second for each mL of clearing volume. CAUTION: The clearing sample should not remain in the reservoir for longer than 2 minutes.
6	Flush the reservoir, samples sites, and line clear of residual blood. Swab the sampling site(s), ensuring removal of any excess blood left on the sampling port. Note: Do not use acetone.

6.6 Drawing Blood Samples

Two methods may be used to draw blood samples using the VAMP Plus closed blood sampling system.

6.7 Syringe Method for Drawing Blood Samples

Step	Procedure
1	Prepare the VAMP Plus reservoir for drawing a blood sample by drawing a clearing volume per one of the methods outlined above.
2	Once the line is clear, swab the VAMP Plus needleless sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy. Note: Do not use acetone.
3	To draw a blood sample, use either a preassembled packaged VAMP needleless cannula and syringe or a VAMP needleless cannula (packaged separately) and separate syringe. CAUTION: Do not use a needle through the sampling site.
4	Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
5	Push the cannula into the VAMP Plus needleless sampling site and then draw the required volume of blood into the syringe. Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.
6	Remove the syringe and cannula from the sampling site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site.
7	Once the sample has been drawn, reinfuse the clearing volume according to one of the methods described above. WARNING: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.
8	To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU). a) Using aseptic technique peel open the pouch. b) Ensure that all connections are tight. c) Hold the BTU in one hand and push the cannula on the filled sample syringe through the BTU needleless injection site. d) Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk of the vacuum tube. e) Fill the vacuum tube to the desired volume. f) Repeat steps (d) and (e) as required to fill additional tubes. g) Discard the BTU after transferring the blood sample from the syringe into the vacuum tubes.
9	Discard all syringes and cannulae after use according to hospital policy.

6.8 Direct-Draw Method for Drawing Blood Samples

Step	Procedure
1	Prepare the VAMP Plus reservoir for drawing a blood sample by drawing a clearing volume per one of the methods outlined above.
2	Once the line is clear, swab the VAMP Plus needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. Note: Do not use acetone.
3	To draw a blood sample, use the direct-draw unit. CAUTION: Do not use a needle through the sampling site. a) Using aseptic technique, peel open the pouch. b) Ensure that the cannula is securely tightened to the direct-draw unit. c) Push the cannula of the direct-draw unit into the sampling site. d) Insert the selected vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube and fill to the desired volume. CAUTION: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity. CAUTION: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions. e) Repeat step (d) as required to fill additional vacuum tubes. f) When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw unit by the cannula and pull straight out . CAUTION: Do not twist the direct-draw unit housing or remove it with the vacuum tube still attached. g) Discard the direct-draw unit after use according to hospital policy.
4	Once the sample has been drawn, reinfuse the clearing volume according to one of the methods described above.

7.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

8.0 MRI Safety Information

MR Conditional

Non-clinical testing demonstrated that the VAMP Plus closed blood sampling system is MR Conditional and can be used in the MRI environment according to the following conditions:

- Static magnetic field of 3Tesla or less
- This device should remain outside of the bore of the MR system and should not be in contact with the patient during MR imaging

- Device should be secured (for example taped down) prior to MRI

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Plus closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

9.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

This product is for single use only. Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (σήμανση CE)	Conformité Européenne (Marcação CE)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/ Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/União Europeia
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl)ftalaat (DEHP)	Indeholder ftalater: bis (2-ethylhexyl) ftalat (DEHP)	Innehåller eller har spår av ftalater: bis(2-etylhexyl) ftalat (DEHP)	Περιέχει φθαλικές ενώσεις ή ίχνη φθαλικών ενώσεων: φθαλικός δι-(2-αιθυλεξυλ) εστέρας (DEHP)	Contém ou estão presentes ftalatos: ftalato de bis (2-etil-hexilo) (DEHP)
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης στον ιστότοπο	Consultar as instruções de utilização no site
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστεριώνετε	Não voltar a esterilizar
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Te gebruiken tot	Anvendes inden	Använd före	Ημερομηνία λήξης	Validade
	Niet-pyrogeen	Ikke-pyrogen	Ikke-pyrogen	Μη πυρετογόνο	Não pirogénico
	MRI-veilig onder bepaalde voorwaarden	MR-betinget	MR-villkorlig	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR) υπό προϋποθέσεις	Utilização condicionada em ambiente de RM
	Breekbaar, voorzichtigheid geboden	Skrøbelig, skal håndteres forsigtigt	Ömtåligt, hanteras försiktigt	Εύθραυστο: χειριστείτε με προσοχή	Frágil, manusear com cuidado
	Deze kant boven	Denne side opad	Denna sida upp	Αυτή η πλευρά προς τα πάνω	Este lado para cima
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo

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	Suomi	Български	Română	Eesti	Lietuvių
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
EC REP	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност/Европейския съюз	Reprezentant autorizat în Comunitatea Europeană/ Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Sisältää tai saattaa sisältää ftalaatteja: bis (2-etyyliheksyyli)ftalaatti (DEHP)	Съдържа или наличие на фталати: бис (2-етилхексил) фталат (DEHP)	Conține sau prezintă ftalați: di (2-etilhexil) ftalat (DEHP)	Sisaldab ftalaate: bis (2-etiülheksüül)ftalaat (DEHP)	Sudėtyje yra ftalatų: bis-(2-etilheksil) ftalatas (DEHP)
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
LOT	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
QTY	Määrä	Количество	Cantitate	Kogus	Kiekis
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
eifu.edwards.com + 1 888 570 4016	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas

	Suomi	Български	Română	Eesti	Lietuvių
	Älä steriloi uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Käytettävä ennen	Срок на годност	Valabil până la	Kõlblik kuni	Naudoti iki
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Ehdollisesti turvallinen magneettikuvauksessa	Безопасно при MR при определени условия	Condiționat RM	Ohutu magnetresonant-stomograafias teatud tingimustel	Sąlyginis MR
	Särkyvä, käsiteltävä varoen	Чупливо, да се борава внимателно	Fragil, a se manipula cu atenție	Õrn, käsitseda ettevaatlikult	Trapus, elgtis atsargiai
	Tämä puoli ylöspäin	С тази страна нагоре	Cu partea aceasta în sus	See pool üles	Šia puse į viršų
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Mallinumero	Номер на модела	Număr model	Mudeli number	Modelio numeris

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Vienreizējai lietošanai	Tek kullanımlık	Только для однократного использования	Jednokratna upotreba	Za jednokratnu upotrebu
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
EC REP	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Importētājs	İthalatçı	Импортёр	Uvoznik	Uvoznik
	Vērojama ftalātu klātbūtne vai satur ftalātus: bis (2-etilheksil)ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
LOT	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
QTY	Daudzums	Miktar	Количество	Količina	Količina
	Uzmanību!	Dikkat	Предостережение	Opres	Opres
	Skatīt lietošanas instrukciju tiešajā vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internetskoj stranici
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilenoksidom	Sterilizirano etilenoksidom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Derīguma termiņš	Son kullanma tarihi:	Использовать до	Upotrebljivo do	Upotrijebiti do
	Nepiroģēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Izmantojams MR vidē, ievērojot noteiktus nosacījumus	MR Koşullu	Условно безопасно при проведении МРТ	Uslovno bezbedno za MR	Uvjetno siguran kod pregleda MR-om
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela

Note: Not all symbols may be included in the labeling of this product. ■ **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. ■ **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. ■ **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ **Obs!** Alla symboler kanske inte används för märkning av denna produkt. ■ **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ **Nota:** poderão não estar incluídos todos os símbolos na rotulagem deste produto. ■ **Poznámka:** Štítky tohoto výrobku nemusí obsahovat všechny symboly. ■ **Megjegyzés:** Elképzelhető, hogy a termék címkéjén nem szerepel az összes szimbólum. ■ **Uwaga:** Na etykiecie niniejszego produktu mogą nie znajdować się wszystkie symbole. ■ **Poznámka:** Je možné, že na štítku pre tento produkt nie sú uvedené všetky symboly. ■ **Merk:** Alle symboler er kanskje ikke inkludert på produktmerkingen. ■ **Huomautus:** kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забелешка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notā:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus.** Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. ■ **Piezīme.** Ši izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ **Not:** Bu ürünün etiketinde tüm semboller yer alması olabilir. ■ **Примечание.** На этикетках данного изделия могут указываться не все символы. ■ **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ **Napomena:** na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.

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Made in Dominican Republic
CCT Critical Care Technologies S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

01/2025
10014980003 B / DOC-0135112 B
© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH