Edwards

Directory

English (en)	1	Português (pt)	28	Română (ro)	53
Français (fr)	4	Česky (cs)	32	Eesti (et)	56
Deutsch (de)	7	Magyar (hu)	35	Lietuvių (lt)	58
Español (es)	10	Polski (pl)	38	Latviešu (lv)	61
Italiano (it)	13	Slovensky (sk)	41	Türkçe (tr)	64
Nederlands (nl)	16	Norsk (no)	44	Русский (ru)	67
Dansk (da)	20	Suomi (fi)	46	Srpski (sr)	71
Svenska (sv)	22	Български (bg)	49	Hrvatski (hr)	74
Ελληνικά (el)	25				
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Φιγυρι ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike77					
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola80					

English

VAMP Closed Blood Sampling System with LASS

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use

For Single Use Only

For figures, please see Figure 1 on page 77 through Figure 10 on page 79.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Edwards Lifesciences VAMP closed blood sampling system provides a safe and convenient method for the withdrawal of blood samples from pressure monitoring lines. The blood sampling system is designed for use with disposable pressure transducers and for connection to central line, venous, and arterial catheters where the system can be flushed clear after sampling. The VAMP closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use of blood sampling technologies as part of their respective institutional guidelines.

VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

3.0 Indications for Use

For adult patients with medical conditions requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

4.0 Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

5.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Some models may contain phthalates, specifically DEHP [Bis (2- ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

6.0 Instructions for Use

CAUTION: The use of lipids with the VAMP closed blood sampling system may compromise product integrity.

6.1 Equipment

- Flush device for flow control (nominal 3 ml/hr flow rate)
- Disposable pressure transducer
- A VAMP closed blood sampling system that contains one reservoir with shut-off valve, and sampling site(s)

6.2 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP kit from the sterile package.
2	If the VAMP kit does not include a pre-connected TruWave disposable pressure transducer, remove the protective caps and attach the distal end with the female luer port to the male luer of a transducer, or other fluid flush device.

Step	Procedure
3	Ensure that all connections are secure and that the sampling site is in the prime/clear position (see Figure 1 on page 77). Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.
4	For the VAMP Flex reservoir , turn the reservoir stopcock off to the patient end of the tubing (see Figure 3 on page 77).
5	For the VAMP Adult reservoir , raise the plunger to facilitate flow of priming solution (see Figure 5 on page 78).
6	For the VAMP Plus reservoir , push the reservoir plunger to the closed and locked position by pressing on the reservoir extension until it latches in the closed position (see Figure 4 on page 77).
7	Connect the kit to an IV fluid source without pressurizing the fluid. Begin to gravity fill the kit first through the transducer and out through the vent port according to the manufacturer's instructions.
8	Replace all vented caps on the side-ports of the stopcocks with non-vented caps.
9	To fill the VAMP Flex reservoir: a) Slowly fill the reservoir to the full 10 ml position by pulling back on the plunger. b) Turn the reservoir stopcock off to the transducer/IV bag (see Figure 7 on page 78). c) Ensure that the sample site(s) are in the prime/clear position (see Figure 1 on page 77). d) Orient the reservoir such that it is in the vertical position and the sampling site is above the reservoir at approximately 45°, and slowly reinfuse the fluid. Repeat as necessary to remove air bubbles. e) Turn the stopcock off to the reservoir (see Figure 8 on page 78).
10	To fill the VAMP Adult reservoir: a) Ensure that the shut-off valve is in the open position indicated by the handle being parallel to the tubing and that the sample site(s) are in the prime/clear position (see Figure 1 on page 77). b) Hold the reservoir by the plunger flexures. Orient the reservoir in the vertical position and the sampling site above the reservoir at approximately 45°. c) Provide flow by pulling the Snap-Tab of the Edwards' TruWave disposable pressure transducer. Deliver priming solution to remove air. Lower plunger to closed position (see Figure 6 on page 78).
11	To fill the VAMP Plus reservoir: a) Ensure that the shut-off valve is in the open position indicated by the handle being parallel to the tubing. b) Ensure that the sample site(s) are in the prime/clear position (see Figure 1 on page 77). c) Orient the reservoir in the vertical position and the sampling site above the reservoir at approximately 45°. d) Provide flow by pulling the Snap-Tab of the Edwards' TruWave disposable pressure transducer. Deliver priming solution to remove air. Note: The transducer and reservoir with bracket can be mounted to an IV pole using the Edwards' disposable pressure transducer holder.
12	Position the LASS sampling site(s) to the pressure monitoring position (see Figure 2 on page 77) and pressurize the IV solution bag per hospital policy and transducer manufacturer's instructions. Flow rate will vary with pressure gradient across the flush device.
13	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter.
14	Connect pressure transducer cable to monitor.
15	Zero the transducer according to its manufacturer's instructions. CAUTION: Remove all air bubbles to reduce the risk of air emboli and reduce loss of pressure signal.

CAUTION: Failure to position LASS sampling site in pressure monitoring position (see Figure 2 on page 77) will adversely impact dynamic response of the pressure signal.

6.3 Clearing IV Fluid in Preparation for Drawing Blood Samples

Note: Use aseptic technique.

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume is required for coagulation studies.

Step	Procedure
1	Position sampling site into the prime/clear position (see Figure 1 on page 77).
2	For VAMP Flex closed blood sampling system , turn the reservoir stopcock off to the transducer/IV bag (see Figure 7 on page 78). Fill the reservoir by drawing the plunger. Draw the appropriate clearing volume according to hospital policy. (Recommended draw rate is 1 second for each ml.) Once the clearing volume is drawn, turn the reservoir stopcock off to the patient and catheter (see Figure 3 on page 77). This will ensure that the sample drawn is from the patient and not the reservoir.
3	For VAMP Adult reservoir , firmly squeeze the flexures and slowly fill the reservoir at its 5 ml volume capacity. (Recommended draw rate is 1 second for each ml.) Close the shut-off valve by turning the handle perpendicular to the tubing.
4	For VAMP Plus closed blood sampling system , release plunger latch and draw back on the reservoir plunger until it stops at its 12 ml volume capacity. (Recommended draw rate is 1 second for each ml.) Close the shut-off valve by turning the handle perpendicular to the tubing. Note: If difficulties are experienced in drawing the clearing sample, check the catheter for possible occlusions or restrictions.

6.4 Drawing Blood Samples Using Syringe or Direct-Draw Method

CAUTION: Do not use a needle or needleless cannula through the sampling site.

6.4.1 Syringe Method

Step	Procedure
1	Always use ISO-594 or EN 20594-1 compliant syringes to avoid potential collapsing of the septum after syringe insertion and removal.
2	Prepare the VAMP kit for drawing a blood sample per the method outlined above.
3	Ensure that the LASS sampling site is still in the prime/clear position (see Figure 1 on page 77).
4	Disinfect the sampling site for access per hospital policy. Note: Do not use acetone.
5	Carefully connect the syringe to the LASS sampling site by slowly pushing the syringe straight into the site using a clockwise, rotating motion until it is fully seated. (See Figure 9 on page 79 through Figure 10 on page 79 for proper syringe insertion techniques.) Do Not try to insert the syringe at an angle. Do Not try to pry open the septum's slit in the sampling site. The slit automatically opens when the syringe is connected. Do Not use excessive force to connect the syringe.
6	Remove syringe from the LASS sampling site, rotating counter clockwise. The blood sampling site's septum completely closes after each use.

6.4.2 Direct-Draw Method

Step	Procedure
1	Always use ISO-594 or EN 20594-1 compliant blood collection unit to avoid potential collapsing of the septum after unit insertion and removal.
2	Prepare the VAMP kit for drawing a blood sample per the method outlined above.
3	Ensure that the sampling site is still in the prime/clear position (see Figure 1 on page 77).

Step	Procedure
4	Disinfect the sampling site for access per hospital policy. Note: Do not use acetone.
5	Carefully connect the direct-draw unit to the LASS sampling site by slowly pushing the unit straight into the site using a clockwise, rotating motion until it is fully seated. • Do Not try to insert the direct-draw unit at an angle. • Do Not try to pry open the septum's slit in the sampling site. The slit automatically opens when the direct-draw unit is connected. • Do Not use excessive force to connect the direct-draw unit.
6	Insert the vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube. Fill to the desired volume. CAUTION: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity.
7	Repeat step 6 as required to fill additional vacuum tubes. CAUTION: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.
8	When the last sample is drawn, remove the vacuum tube first.
9	Remove direct-draw unit from sampling site, rotating counter clockwise. The blood sampling site's septum completely closes after each use.
10	Discard the direct-draw unit after use according to hospital policy.

CAUTION: The clearing sample should not remain in the reservoir for longer than 2 minutes.

6.5 Reinfusing Clearing Volume

Step	Procedure
1	For VAMP Flex sampling system: - Once the sample is drawn, turn the reservoir stopcock off to the transducer/IV bag (see Figure 7 on page 78). - Ensure that the sample site(s) are in the prime/clear position (see Figure 1 on page 77). Reinfuse the clearing volume by pushing the plunger down. - If blood has been drawn into the VAMP Flex reservoir, extra steps might be necessary to clear the reservoir. To accomplish this, turn the reservoir stopcock off to the patient (see Figure 3 on page 77) and partially fill the reservoir with 1 to 2 ml of fluid. If the IV bag is properly pressurized, activating the flush device will automatically cause the reservoir to fill. Turn the reservoir stopcock off to the transducer/IV bag (see Figure 7 on page 78) and slowly reinfuse into the line. Repeat as necessary until all blood is removed. Note: For optimal clearing of the reservoir, it is recommended to repeat this clearing process multiple times with small volume of fluid (i.e., 1 to 2 ml). - Once the reservoir has been cleared, turn the stopcock off to the reservoir (see Figure 8 on page 78).
2	For VAMP Adult closed blood sampling system: - Once the sample is drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing. Turn the handle on the sampling site to the prime/clear position (see Figure 1 on page 77). - Press together the reservoir flexures extension (see Figure 6 on page 78) until the plunger stops and is in the closed position.
3	For VAMP Plus closed blood sampling system: - Once the sample is drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing. Turn the handle on the sampling site to the prime/clear position (see Figure 1 on page 77). - Close the reservoir by pushing the plunger until it stops and firmly latches in the closed position. Note: Recommended time to reinfuse the clearing volume is approximately 1 second for each ml.

Step	Procedure
4	Flush sample site(s) and line clear of residual blood by pulling the Snap-Tab of the Edwards' TruWave disposable pressure transducer.
5	Return the sampling site to the pressure monitoring position (see Figure 2 on page 77).
6	Clean the sampling site, ensuring removal of any excess blood. Note: Do not use acetone.

CAUTION: Failure to position sampling site in pressure monitoring position (see Figure 2 on page 77) will adversely impact dynamic response of the pressure signal.

7.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

8.0 MRI Safety Information

8.1 VAMP Flex closed blood sampling system and VAMP Adult closed blood sampling system

MR Safe

The VAMP Flex and VAMP Adult closed blood sampling system are MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

8.2 VAMP Plus closed blood sampling system

8.2.1 Warnings

For a patient undergoing an MRI examination, please refer to the MRI Safety Information section for specific conditions to ensure patient safety.

MR Conditional

Non-clinical testing demonstrated that the VAMP Plus closed blood sampling system is MR Conditional and can be used in the MRI environment according to the following conditions:

<conditions for
MR safety>

- Static magnetic field of 3- T
- Spatial gradient field of 720- gauss/cm or less
- This device is not an implant. Therefore, image artifact, RF Heating and SAR are not applicable
- This device is not intended for use inside of the bore of the MR system and should not be in contact with the patient

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Plus closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

9.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

This product is for single use only.

Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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© Copyright 2025, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.