Edwards

Directory

English (en)	1	Português (pt)	27	Română (ro)	50
Français (fr)	4	Česky (cs)	30	Eesti (et)	53
Deutsch (de)	7	Magyar (hu)	33	Lietuvių (lt)	56
Español (es)	10	Polski (pl)	36	Latviešu (lv)	59
Italiano (it)	13	Slovensky (sk)	39	Türkçe (tr)	61
Nederlands (nl)	16	Norsk (no)	42	Русский (ru)	64
Dansk (da)	19	Suomi (fi)	44	Srpski (sr)	67
Svenska (sv)	21	Български (bg)	47	Hrvatski (hr)	70
Ελληνικά (el)	24				
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike 73-75					
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola 76-83					

English

VAMP Jr. Closed Blood Sampling System

For Neonatal and Pediatric Applications

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use For Single Use Only

For figures please refer to Figure 1 on page 73 through Figure 6 on page 75.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Edwards Lifesciences VAMP Jr. blood sampling system provides a safe and convenient method for withdrawing blood samples. The blood sampling system is designed for connection to umbilical, arterial and central line catheters, and for use with disposable and reusable pressure transducers, if desired.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use of blood sampling technologies as part of their respective institutional guidelines.

VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

3.0 Indications for Use

For pediatric patients (including neonates) with medical conditions requiring periodic withdrawal of blood samples from umbilical, arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

4.0 Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial or umbilical applications.

There are no absolute contraindications when used for venous applications.

5.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Some models may contain phthalates, specifically DEHP [Bis (2-ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

6.0 Instructions for Use

CAUTION: The use of lipids with the VAMP Jr. closed blood sampling system may compromise product integrity.

7.0 Equipment

7.1 Pediatric Applications

- Flush device or flow controlling device
- Disposable or reusable pressure transducer, if desired
- VAMP Jr. blood sampling system with 3 ml VAMP Jr. reservoir, one shut-off valve, and VAMP needleless sampling site

7.2 Pediatric and Neonatal Applications

- Flush device or flow controlling device
- Disposable or reusable pressure transducer, if desired
- VAMP Jr. blood sampling system with one 3 ml VAMP Jr. reservoir, two shut-off valves, one VAMP needleless sampling site, and one VAMP needleless access/flush site

8.0 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP Jr. kit from the sterile package.
2	Attach the distal end with the female luer-lock connector to the transducer stopcock on the disposable or reusable transducer dome, or to an IV line/IV pump per hospital protocol (see Figure 1 on page 73).
3	All connections should be secure. Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.
4	Pull reservoir plunger to approximately 1/2 ml.

Edwards, Edwards Lifesciences, the stylized E logo, VAMP, and VAMP Jr. are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Step	Procedure
5	If a transducer is in the line, deliver flush solution first through the transducer and out through its vent port according to your transducer manufacturer's instructions. Note: Replace all vented caps on the sideports of the stopcocks with non-vented caps.
6	To fill the VAMP Jr. kit, ensure that the shut-off valve is in the open position indicated by the handle being parallel to the tubing.
7	Hold the VAMP Jr. kit vertically and deliver flush solution slowly to fill the reservoir.
8	Push reservoir plunger down and fill the remainder of the kit. CAUTION: Remove all air bubbles to reduce the risk of air emboli.
9	If a transducer is being used, mount the transducer per hospital procedure.
10	Pressurize the IV solution bag and set flow rate per hospital policy.
11	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter. Important: Ensure the VAMP Jr. kit and the patient's catheter is a fluid-to-fluid connection.
12	If necessary, zero and calibrate the transducer according to the transducer manufacturer's instructions.

9.0 Drawing the Clearing Volume

Note: Use aseptic technique

Step	Procedure
1	For VAMP Jr. systems that have a single shut-off valve, close transducer or IV line stopcock to patient.
2	For VAMP Jr. systems that have two shut-off valves, close distal shut-off valve by turning the handle perpendicular to the tubing.
3	Pull up on the reservoir plunger until the plunger stops and the reservoir has reached its 3 ml capacity (see Figure 2 on page 74). Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies. Important: Recommended rate to pull the reservoir plunger up to the fully open position is approximately 1 ml every 10 to 15 s (approximately 30 to 45 s to fill reservoir to capacity). Note: If difficulties are experienced or bubbles are observed when drawing the clearing volume, check catheter and kit for possible loose connections, occlusions, or restrictions (e.g. positional lines).

Step	Procedure
4	To remove bubbles in line between the patient connection and the proximal sampling site while drawing clearing volume: a) Close proximal shut-off valve by turning handle perpendicular to the line. b) Insert syringe and cannula into the proximal sampling site. c) Invert kit. d) Aspirate bubbles from sampling site into syringe. e) Turn kit upright. f) Open proximal shut-off valve by turning handle parallel to the line. g) Continue drawing clearing volume.
5	Once the clearing volume has been drawn, close the shut-off valve by turning the handle perpendicular to the tubing. This will ensure that the sample is drawn from the patient and not the reservoir. Note: To remove bubbles from the VAMP Jr. reservoir, see section Reinfusing Clearing Volume from VAMP Jr. Reservoir.

10.0 Drawing Blood Samples from the VAMP Needleless Sampling Site

Although a variety of techniques can be used for drawing samples, the following guidelines are provided as an aid to the clinician:

Once clearing volume has been drawn, two methods may be used to draw blood samples from the proximal VAMP needleless sampling site. Method one uses a sampling syringe with the VAMP needleless cannula. Method two, direct line sampling, uses a direct-draw unit with an integral VAMP needleless cannula.

10.1 Drawing Blood Samples Using Method One (Syringe and Cannula)

Step	Procedure
1	Swab the VAMP Jr. needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. Note: Do not use acetone.
2	Using a cannula (packaged separately; see Figure 3 on page 74) and syringe, draw blood sample in the following manner: CAUTION: Do not use a needle through the sampling site. a) Using aseptic technique, open the cannula package. b) Attach the cannula onto a selected luer-tip syringe (size determined by hospital protocol) by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure. c) Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.

Step	Procedure
2	d) Push the cannula into the proximal VAMP Jr. needleless sampling site and hold in place for approximately 1 to 2 s (see Figure 4 on page 74). e) Draw the required volume of blood into the syringe. Note: If difficulties are experienced in drawing the sample, check the catheter and VAMP Jr. kit for possible occlusions or restrictions. f) While holding the cannula, remove the syringe and cannula from the sampling site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site. g) Once the last sample has been drawn, swab the sampling site to ensure removal of any residual blood remaining on the sampling port. Note: Do not use acetone.

10.2 Transferring Blood Samples

Step	Procedure
1	To transfer the blood sample from the syringe to vacuum tubes, use the Blood Transfer Unit (BTU) (see Figure 4 on page 74). a) Using aseptic technique, peel open the pouch. b) Ensure that all connections are tight. c) Hold the VAMP Jr. blood transfer unit in one hand and push the cannula on the filled sample syringe through the VAMP blood transfer unit needleless injection site. d) Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk on the vacuum tube. e) Fill the vacuum tube to the desired volume. f) Repeat steps (d) and (e) according to the requirements for the patient's blood study.
2	According to hospital policy, discard the VAMP blood transfer unit, syringes, and cannulae after transferring the blood sample from the syringe into the vacuum tubes.

10.3 Drawing Blood Samples Using Method Two (Direct-Draw Method; Recommended for Pediatric Applications Only)

Step	Procedure
1	Swab the sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy. Note: Do not use acetone.
2	To draw a blood sample, use the direct-draw unit. CAUTION: Do not insert a needle into the sampling site. a) Using aseptic technique, peel open the pouch.

Step	Procedure
2	b) Ensure that the cannula is securely tightened to the direct-draw housing. c) Position the sampling site so that it faces upward. d) Push the cannula of the direct-draw unit into the sampling site (see Figure 5 on page 74). e) Insert the selected vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube. f) Begin to fill the vacuum tube. Remove the vacuum tube prior to it reaching its full capacity. Note: If difficulties are experienced in drawing the sample, check the catheter and VAMP Jr. kit for possible occlusions or restrictions. g) Repeat steps (e) and (f) according to the requirements of the patient's blood study. h) When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw unit by the cannula and pull straight out. CAUTION: Do not twist the direct-draw unit housing or remove it with the vacuum tube still attached. i) Discard the direct-draw unit after use according to hospital policy. j) Once the last sample has been drawn, swab the sampling site to ensure removal of any excess blood remaining on the sampling port. Note: Do not use acetone.

11.0 Reinfusing Clearing Volume from the VAMP Jr. Reservoir

Important: The clearing volume should not remain in the reservoir for longer than 3 minutes. It is recommended that the VAMP Jr. kit be changed or evacuated per hospital protocol if the clearing volume remains in the reservoir for longer than 3 minutes.

Step	Procedure
1	Open the proximal shut-off valve by turning the handle parallel to the tubing.
2	Slowly push down on the plunger until it is in the fully closed position and all fluid is re-infused into the line. Important: Recommended rate to push the reservoir plunger to the closed position is approximately 1 ml every 10 to 15 s (approximately 30 to 45 s to empty a completely filled VAMP Jr. reservoir).
3	For a VAMP Jr. system with one shut-off valve, if bubbles are observed in the reservoir, as a result of loose catheter connections, catheter occlusion, or restrictions: a) Open reservoir to approximately 0.5 ml.
3	b) Invert kit. c) Flush the line (see Flushing the Line).

Step	Procedure
4	For a VAMP Jr. system with two shut-off valves, if bubbles are observed in the reservoir, as a result of loose catheter connections, catheter occlusion, or restrictions: a) Open reservoir to approximately 0.5 ml. b) Invert kit. c) Insert syringe and cannula into distal access flush site. d) Open distal shut-off valve. e) Aspirate bubbles from distal access flush site. f) Close distal shut-off valve and close reservoir. g) Flush the line (see Flushing the Line).

12.0 Flushing the Line

Step	Procedure
1	If a flush device is used to clear the line, flush per hospital procedure.
2	If flushing the line with a syringe and a VAMP Jr. system that has two shut-off valves, do the following: a) Fill syringe with flush solution. b) Attach a VAMP needleless cannula to the prefilled syringe by aligning the cannula luer-lock to the luer-tip on the syringe and twist until secure. c) Ensure that the syringe and cannula are free of bubbles and air. d) Insert syringe and cannula into the distal access/flush site (see Figure 6 on page 75). e) Flush line per hospital procedure. Note: The approximate volume required to flush the line is 1.0 to 1.5 ml. f) Open transducer or IV line stopcock to patient.

12.1 Transferring Blood Samples to Small Collection Tubes

When patient size and clinical condition do not permit large blood samples to be drawn:

Step	Procedure
1	Transfer the blood sample from the syringe to blood sample containers per hospital procedure. Note: If a VAMP needleless cannula was used, the cannula may need to be removed prior to transferring the blood sample from the syringe to the blood sample container.
2	Discard syringes and cannulae according to hospital policy after transferring the blood sample from the syringe into the sample containers.

WARNING: Laboratory values should correlate with the patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

13.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

14.0 MRI Safety Information

The VAMP Jr. closed blood sampling system is MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Jr. closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

15.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

Visually inspect for breaches of packaging integrity prior to use.

16.0 Storage

Store in a cool, dry place.

17.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

18.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

19.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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


	English	Français	Deutsch	Español	Italiano
	Contains or presence of Phthalates: bis (2-ethylhexyl) phthalate (DEHP)	Traces ou présence de phtalates : phtalate de bis (2-éthylhexyle) (DEHP)	Enthält Phthalate oder kann Spuren von Phthalaten enthalten: Bis(2-ethylhexyl)phthalat (DEHP)	Contiene o hay presencia de ftalatos: Di[2-etilhexil] ftalato (DEHP)	Contiene o sono presenti ftalati: di (2-etilesil) ftalato (DEHP)
	Single use	Ne pas réutiliser	Zum einmaligen Gebrauch	Un solo uso	Monouso
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)
	Authorized representative in the European Community/European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Quantity	Quantité	Menge	Cantidad	Quantità
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare
	Use by	Utiliser avant	Verwendbar bis	Caducidad	Utilizzare entro
	Importer	Importateur	Importeur	Importador	Importatore
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
	MR Safe	Aucun risque en milieu RM	MR-sicher	Seguro para RM	Compatibile con RM

	English	Français	Deutsch	Español	Italiano
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil; manipular con cuidado	Fragile, maneggiare con cura
	This way up	Flèches indiquant le haut du carton	Oben	Este lado arriba	Alto
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarrieresystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
Note: Not all symbols may be included in the labeling of this product. ■ Remarque : il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ Hinweis: Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ Nota: Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ Nota: alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto.					

Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl)ftalaat (DEHP)	Indeholder ftalater: bis (2-ethylhexyl) ftalat (DEHP)	Innehåller eller har spår av ftalater: bis(2-etylhexyl) ftalat (DEHP)	Περιέχει φθαλικές ενώσεις ή ίχνη αυτών: φθαλικός δι (2-αιθυλοεξυλο) εστέρας (DEHP)	Contém ou estão presentes ftalatos: ftalato de bis (2-etil-hexilo) (DEHP)
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (σήμανση CE)	Conformité Européenne (Marcação CE)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/ Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso

	Nederlands	Dansk	Svenska	Ελληνικά	Português
 eifu.edwards.com + 1 888 570 4016	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστεριώνετε	Não voltar a esterilizar
	Te gebruiken tot	Anvendes inden	Använd före	Ημερομηνία λήξης	Validade
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Niet-pyrogeen	Ikke-pyrogen	Ikke-pyrogen	Μη πυρετογόνο	Não pirogénico
	MRI-veilig	MR-sikker	MR-säker	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)	Utilização segura em ambiente de RM
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Breekbaar, voorzichtigheid geboden	Skrøbelig, skal håndteres forsigtigt	Ömtåligt, hanteras försiktigt	Εύθραστο: χειριστείτε με προσοχή	Frágil, manusear com cuidado
	Deze kant boven	Denne side opad	Denna sida upp	Αυτή η πλευρά προς τα πάνω	Este lado para cima
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barriersystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifiering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
Opmerking: Het label van dit product bevat misschien niet alle symbolen. ■ Bemærk: Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ Obs! Alla symboler kanske inte används för märkning av denna produkt. ■ Σημείωση: Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ Nota: poderão não estar incluídos todos os símbolos na rotulagem deste produto.					

Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring

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	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Vērojama ftalātu klātbūtne vai satur ftalātus: bis(2-etilheksil)ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Daudzums	Miktar	Количество	Količina	Količina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Derīguma termiņš	Son kullanma tarihi	Использовать до	Upotrebljivo do	Upotrijebiti do
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	Nepirogeno
	Drošs lietošanai MR vidē	MR için Güvenli	Безопасно при проведении МРТ	Bezbedno za MR	Sigurno kod pregleda MR-om

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
Piezīme. Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ Not: Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ Примечание. На этикетках данного изделия могут указываться не все символы. ■ Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ Napomena: na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.					

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Made in Dominican Republic
CCT Critical Care Technologies S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

08/24
10014981003 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3274
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH