Edwards


VAMP Jr. Closed Blood Sampling System

Directory		
English (en)	1	
Français (fr)	4	
Deutsch (de)	7	
Español (es)	10	
Italiano (it)	13	
Nederlands (nl)	16	
Dansk (da)	19	
Svenska (sv)	22	
Ελληνικά (el)	25	
Português (pt)	28	Română (ro) 53
Česky (cs)	31	Eesti (et) 56
Magyar (hu)	34	Lietuvių (lt) 59
Polski (pl)	37	Latviešu (lv) 62
Slovensky (sk)	40	Türkçe (tr) 65
Norsk (no)	43	Русский (ru) 68
Suomi (fi)	46	Srpski (sr) 72
Български (bg)	49	Hrvatski (hr) 75
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Φιγυρι ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike 78		
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola 83		

English

VAMP Jr. Closed Blood Sampling System For neonatal and pediatric applications

V1047

Instructions for Use For Single Use Only

For figures please refer to Figure 2 on page 78 through Figure 5 on page 81.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Edwards Lifesciences VAMP Jr. closed blood sampling system provides a safe and convenient method for withdrawing blood samples. The

VAMP Jr. system is designed for use with disposable pressure transducers and for connection to central line and arterial catheters.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use of blood sampling technologies as part of their respective institutional guidelines.

- VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

3.0 Indications for Use

For pediatric patients (including neonates) with medical conditions requiring periodic withdrawal of blood samples from umbilical, arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

4.0 Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial or umbilical applications.

There are no absolute contraindications when used for venous applications.

Edwards, Edwards Lifesciences, the stylized E logo, VAMP, and VAMP Jr. are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

5.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Some models may contain phthalates, specifically DEHP [Bis (2-ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

6.0 Instructions for Use

CAUTION: The use of lipids with the VAMP Jr. closed blood sampling system may compromise product integrity.

6.1 Equipment

- Flush device or flow controlling device (maximum 4 ml/hour flow rate)
- Disposable pressure transducer, if desired
- VAMP Jr. closed blood sampling system with integral shut-off valve and needleless sampling site

6.2 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP Jr. kit from the sterile package.
2	Attach the distal end with the female luer-lock connector to the transducer stopcock on the disposable transducer or reusable transducer dome, or to an IV line/IV infusion pump. Ensure that all connections are secure.
3	If a transducer is in the line, deliver flush solution first through the transducer and out through the transducer vent port according to the manufacturer's instructions.
4	Replace all vented caps on the sideports of the stopcocks with non-vented caps.
5	Fill the VAMP Jr. kit with flush solution using gravity and ensure that the shut-off valve by the clearing volume (distal) access site is in the open position by turning the handle parallel to the tubing. Deliver flush solution slowly through the kit. CAUTION: Remove all air bubbles to reduce the risk of air emboli.
6	If a transducer is being used, mount the transducer either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder.
7	Pressurize the IV solution bag. Flow rate will vary with pressure across the flush device.
8	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter (Figure 1 on page 78). Important: Ensure the VAMP Jr. kit and the patient's catheter is a fluid-to-fluid connection.

Step	Procedure
9	If necessary, zero and calibrate the transducer according to the transducer manufacturer's instructions.

7.0 Drawing Clearing Volume and Blood Samples

Note: Use aseptic technique.

7.1 Drawing Clearing Volume

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

Step	Procedure
1	To draw clearing volume, use a VAMP Jr. needleless cannula (packaged separately) and syringe. Do not insert a needle into the clearing volume (distal) access site. Using VAMP Jr. needleless cannula: a) Using aseptic technique, peel open the cannula pouch. b) Pick up the cannula by the protective shield (see Figure 2 on page 78). c) Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure. d) Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
2	Swab the clearing volume (distal) access site with disinfectant such as alcohol or betadine, depending upon hospital policy. Note: Do not use acetone.
3	Push the cannula into the VAMP Jr. needleless clearing volume (distal) access site, hold in place for approximately 1 to 2 s, and draw the clearing volume (see Figure 3 on page 79). Important: Recommended rate to draw the clearing volume into the syringe is 1 ml every 10 to 15 s. Note: If difficulties are experienced in drawing the clearing volume, check the catheter and VAMP Jr. kit for possible occlusions or restrictions.
4	Once the clearing volume has been drawn, close the 3-way shut-off valve at the clearing volume (distal) access site by turning the handle perpendicular to the tubing. This will ensure that the sample is drawn from the patient and not the distal tubing.
5	Close 1-way shut-off valve.
6	While holding the cannula, remove the syringe and cannula from the clearing volume access site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site. Important: Discard the clearing volume according to hospital protocol.
7	Swab the access site to ensure removal of any excess blood. Note: Do not use acetone.

Step	Procedure
8	Once the syringe and cannula is removed from the clearing volume access site, open 1-way shut-off valve by turning the handle parallel to tubing.

8.0 Drawing Blood Samples from the VAMP Jr. Needleless Sampling Site

Although a variety of techniques can be used for drawing samples, the following guidelines are provided as an aid to the clinician:

Two methods may be used to draw blood samples from the proximal VAMP Jr. needleless sampling site. Method one uses a sampling syringe with the VAMP Jr. needleless cannula. Method two, direct line sampling, uses a direct-draw unit with an integral VAMP Jr. needleless cannula.

8.1 Drawing Blood Samples Using Method One (Syringe and Cannula)

Step	Procedure
1	Swab the VAMP Jr. needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. Note: Do not use acetone.
2	To draw a blood sample, use a cannula (packaged separately) and syringe. CAUTION: Do not use a needle through the sampling site.
3	Using an individually packaged cannula: a) Using aseptic technique, open the cannula package. b) Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure.
4	Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
5	Push the cannula into the proximal VAMP Jr. needleless sampling site and hold in place for approximately 1 to 2 s (see Figure 3 on page 79).
6	Draw the required volume of blood into the syringe. Note: If difficulties are experienced in drawing the sample, check the catheter and VAMP Jr. kit for possible occlusions or restrictions.
7	Close 1-way shut-off valve.
8	While holding the cannula, remove the syringe and cannula from the sampling site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site.
9	Once the last sample has been drawn, swab the sampling site to ensure removal of any residual blood remaining on the sampling port. Note: Do not use acetone.

8.2 Transferring Blood Samples

Step	Procedure
1	To transfer the blood sample from the syringe to vacuum tubes, use the Blood Transfer Unit (BTU) (see Figure 4 on page 80). a) Using aseptic technique, peel open the pouch.

Step	Procedure
	b) Ensure that all connections are tight. c) Hold the VAMP Jr. blood transfer unit in one hand and push the cannula on the filled sample syringe through the VAMP blood transfer unit needleless injection site. d) Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk on the vacuum tube. e) Fill the vacuum tube to the desired volume. f) Repeat steps (d) and (e) according to the requirements for the patient's blood study.
2	According to hospital policy, discard the VAMP blood transfer unit, syringes, and cannulae after transferring the blood sample from the syringe into the vacuum tubes.

8.3 Drawing Blood Samples Using Method Two (Direct-Draw Method); Recommended for Pediatric Applications Only

Step	Procedure
1	Swab the sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy. Note: Do not use acetone.
2	To draw a blood sample, use the direct-draw unit. CAUTION: Do not use a needle through the sampling site. a) Using aseptic technique, peel open the pouch. b) Ensure that the cannula is securely tightened to the direct-draw housing. c) Position the sampling site so that it faces upward. d) Push the cannula of the direct-draw unit into the sampling site (see Figure 5 on page 81). e) Insert the selected vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube. CAUTION: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity. f) Fill the vacuum tube to the desired volume. Note: If difficulties are experienced in drawing the sample, check the catheter and VAMP Jr. kit for possible occlusions or restrictions. g) Repeat steps (e) and (f) according to the requirements of the patient's blood study.

Step	Procedure
	h) When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw unit by the cannula and pull straight out. CAUTION: Do not twist the direct-draw unit housing or remove it with the vacuum tube still attached.
3	Discard the direct-draw unit after use according to hospital policy.
4	Once the last sample has been drawn, swab the sampling site to ensure removal of any excess blood remaining on the sampling port. Note: Do not use acetone.

8.4 Flushing the Line

Step	Procedure
1	After the sample has been drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing.
2	Flush the line clear.
3	Discard all syringes and cannulas after use according to hospital policy.

8.5 Transferring Blood Samples to Small Collection Tubes

When patient size and clinical condition do not permit large blood samples to be drawn:

Step	Procedure
1	Transfer the blood sample from the syringe to blood sample containers per hospital procedure. Note: If a VAMP needleless cannula was used, the cannula may need to be removed prior to transferring the blood sample from the syringe to the blood sample container.
2	Discard syringes and cannulae according to hospital policy after transferring the blood sample from the syringe into the sample containers. WARNING: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

9.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

10.0 MRI Safety Information

The VAMP Jr. blood sampling system is MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Jr. blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

11.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize. This product is for single use only.

Visually inspect for breaches of packaging integrity prior to use.

12.0 Storage

Store in a cool, dry place.

13.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

14.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
 (24 hours): 800.822.9837
 Outside the U.S. and Canada
 (24 hours): 949.250.2222
 In the UK: 0870 606 2040 - Option 4
 In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

15.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (σήμανση CE)	Conformité Européenne (Marcação CE)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/ Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl) ftalaat (DEHP)	Indeholder ftalater: bis (2-ethylhexyl) ftalat (DEHP)	Innehåller eller har spår av ftalater: bis(2-etylhexyl) ftalat (DEHP)	Περιέχει φθαλικές ενώσεις ή ίχνη αυτών: φθαλικός δι (2-αιθυλοεξυλο) εστέρας z(DEHP)	Contém ou estão presentes ftalatos: ftalato de bis (2-etil-hexilo) (DEHP)
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Te gebruiken tot	Anvendes inden	Använd före	Ημερομηνία λήξης	Validade
	Niet-pyrogeen	Ikke-pyrogen	Ikke-pyrogen	Μη πυρετογόνο	Não pirogénico
	MRI-veilig	MR-sikker	MR-säker	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)	Utilização segura em ambiente de RM
	Breekbaar, voorzichtigheid geboden	Skrøbelig, skal håndteres forsigtigt	Ömtåligt, hanteras försiktigt	Εύθραστο: χειριστείτε με προσοχή	Frágil, manusear com cuidado
	Deze kant boven	Denne side opad	Denna sida upp	Αυτή η πλευρά προς τα πάνω	Este lado para cima
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara produkten svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo

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

	Suomi	Български	Română	Eesti	Lietuvių
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks	Vienkartinio naudojimo
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijoje / Europos Sąjungoje
	Älä käyttää, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniiseade	Medicinos priemonė
	Sisältää tai saattaa sisältää ftalaatteja: bis(2-etyyliheksyyli)ftalaatti (DEHP)	Съдържа или наличие на фталати: бис(2-етилхексил)фталат (DEHP)	Conține sau prezintă ftalați: di (2-etilhexil) ftalat (DEHP)	Sisaldab ftalaate: bis(2-etiülheksüül)ftalaat (DEHP)	Sudėtyje yra ftalatu: bis-(2-etilheksil) ftalatas (DEHP)
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas interneto svetainėje
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkanjal või elektroonilise kasutusjuhisega	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Älä steriloitu uudelleen	Не стерилизирайте повторно	A nu se resteriliza	Ärge resteriliseerige	Kartotinai nesterilizuokite
	Maahantuojaja	Вносител	Importator	Importija	Importuotojas
	Steriloitu etyleenioksidilla	Стерилизирано с этиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu
	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant
	Käytettävä ennen	Срок на годност	Valabil până la	Kõlblik kuni	Naudoti iki
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepirogeenne	Nepirogeniškas
	Sopii magneettikuvaukseen	Безопасно при MR	Sigur în utilizarea la RM	Ohutu magnetresonantsto mograafias	MR saugus
	Särkyvä, käsiteltävä varoen	Чупливо, да се борава внимателно	Fragil, a se manipula cu atenție	Õrn, käsitseda ettevaatlikult	Trapus, elgtis atsargiai
	Tämä puoli ylöspäin	С тази страна нагоре	Cu partea aceasta în sus	See pool üles	Šia puse į viršų
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitada jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris

Huomautus: kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ **Забележка:** Възможно е не всички символи да са включени в етикета на този продукт. ■ **Notă:** este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ **Märkus:** Kõik sümboolid ei pruugi antud toote etiketil esineda. ■ **Pastaba:** Šio gaminio etiketėje gali būti pateikti ne visi simbolai.

Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Vienreizējai lietošanai	Tek kullanımlık	Только для одноразового использования	Jednokratna upotreba	Za jednokratnu upotrebu
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformitē Europēenne (CE marķējums)	Conformitē Europēenne (CE İşareti)	Conformitē Europēenne (маркировка CE)	Conformitē Europēenne (CE znak)	Conformitē Europēenne (oznaka CE)
	Pilnvarotais pārstāvis Eiropas Kopienā/ Eiropas Savienībā	Avrupa Birliği'nde/ Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/ Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj uniji
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Vērojama ftalātu klātbūtne vai satur ftalātus: bis(2-etilheksil)ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Daudzums	Miktar	Количество	Količina	Količina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Derīguma termiņš	Son kullanma tarihi:	Использовать до	Upotrebljivo do	Upotrijebiti do
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Drošs lietošanai MR vidē	MR için Güvenli	Безопасно при проведении МРТ	Bezbedno za MR	Sigurno kod pregleda MR-om
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore
	Viena sterilitātes aizsarglāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barjerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela

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EC REP

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

Made in Dominican Republic
CCT Critical Care Technologies S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

08/24
10021365002 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU