## Directory

English (en) ..... 1	Português (pt) ..... 24	Română (ro) ..... 44
Français (fr) ..... 4	Česky (cs) ..... 27	Eesti (et) ..... 47
Deutsch (de) ..... 6	Magyar (hu) ..... 29	Lietuvių (lt) ..... 49
Español (es) ..... 9	Polski (pl) ..... 32	Latviešu (lv) ..... 52
Italiano (it) ..... 12	Slovensky (sk) ..... 34	Türkçe (tr) ..... 54
Nederlands (nl) ..... 14	Norsk (no) ..... 37	Русский (ru) ..... 56
Dansk (da) ..... 17	Suomi (fi) ..... 39	Srpski (sr) ..... 59
Svenska (sv) ..... 19	Български (bg) ..... 41	Hrvatski (hr) ..... 61
Ελληνικά (el) ..... 21		


Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike ..... 64

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola ..... 68

## English

# VAMP Closed Blood Sampling System for venous applications without a flow controlling device

V1032

### Instructions for Use

#### For Single Use Only

For figures, see Figure 1 on page 64 through Figure 9 on page 67.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

#### 1.0 Description

The Edwards Lifesciences VAMP closed blood sampling system provides a safe and convenient method for withdrawing blood samples. The VAMP system is designed for connection to central line and arterial catheters. The VAMP closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use. The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use of blood sampling technologies as part of their respective institutional guidelines.

VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

#### 2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

#### 3.0 Indications for Use

For adult patients with medical conditions requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

#### 4.0 Contraindications

There are no absolute contraindications when used for venous applications.

#### 5.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Some models may contain phthalates, specifically DEHP [Bis (2- ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Edwards, Edwards Lifesciences, the stylized E logo, and VAMP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

## 6.0 Instructions for Use

**CAUTION: The use of lipids with the VAMP closed blood sampling system may compromise product integrity.**

### 6.1 Equipment

- Heparin lock or IV pump/IV line
- VAMP closed blood sampling system with one 5 ml VAMP reservoir, two integral shut-off valves, and one VAMP needleless sampling site.

### 6.2 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP kit from the sterile package.
2	Attach the female luer-lock connector at the distal end of the VAMP kit to a heparin lock or an IV source per hospital procedure, depending on the specific application.
3	Ensure that the reservoir plunger is down in the closed and locked position.
4	Fill the VAMP kit with flush solution using gravity and ensure that the shut-off valves are in the open position by turning the handles parallel to the tubing. Hold the reservoir by the plunger flexures. Orient the kit such that the reservoir is in a vertical position and the sampling site is above the reservoir at approximately 45°. In this position, gently squeeze and hold the plunger flexures together to slightly raise the seal to enlarge the fluid path and deliver flush solution through the kit. Once the kit is filled, release the plunger flexures. <b>CAUTION: Remove all air bubbles to reduce the risk of air emboli.</b>
5	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter (Figure 1 on page 64).

## 7.0 Drawing Blood Samples

**Note: Use aseptic technique.**

Although a variety of techniques can be used for drawing samples, the following guidelines are provided as an aid to the clinician:

### 7.1 Drawing the Clearing Volume

**Important: A minimum clearing volume of two times the dead space should be achieved. Additional drawing volume may be required for coagulation studies.**

Step	Procedure
1	Close the distal shut-off valve by turning the handle perpendicular to the tubing. This will ensure that the fluid is drawn from the patient and not the IV source (see Figure 2 on page 65).
2	Gently squeeze the reservoir plunger flexures together. Pull up on the reservoir plunger until it stops and the reservoir has reached its 5 ml volume capacity (see Figure 3 on page 65). <b>Note: If difficulties are experienced when drawing the clearing volume, check the catheter for possible occlusions or restrictions (e.g., positional lines).</b> <b>Important: Recommended rate to pull the reservoir plunger up to the fully open position is 1 ml/second.</b>

Step	Procedure
3	Once the clearing volume has been drawn, close the proximal shut-off valve by turning the handle perpendicular to the tubing (see Figure 4 on page 65). This will further ensure that the sample is drawn from the patient and not the reservoir.

## 8.0 Drawing Blood Samples from the VAMP Needleless Sampling Site

Two methods may be used to draw blood samples from the VAMP needleless sampling site. Method one uses a sampling syringe with the VAMP needleless cannula and a blood transfer unit (BTU). Method two, direct line sampling, uses a direct-draw with an integral VAMP needleless cannula.

### 8.1 Drawing Blood Samples Using Method One (Syringe and Cannula)

Step	Procedure
1	Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. <b>Note: Do not use acetone.</b>
2	To draw a blood sample, use a VAMP needleless cannula (packaged separately) and syringe. Do not use a needle through the sampling site.
3	Using VAMP needleless cannula: a) Using aseptic technique, peel open the cannula pouch. b) Pick up the cannula by the protective shield (see Figure 6 on page 66). c) Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure.
4	Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
5	Push the cannula into the VAMP needleless sampling site and hold in place for approximately 1 to 2 s (see Figure 7 on page 66).
6	Draw the required volume of blood into the syringe. <b>Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.</b>
7	While holding the cannula, remove the syringe and cannula from the sampling site by <b>pulling straight out</b> . <b>CAUTION: Do not twist the syringe out of the sampling site.</b>

## 8.2 Transferring Blood Samples

Step	Procedure
1	To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU) (see Figure 8 on page 67). a) Using aseptic technique, peel open the pouch. b) Ensure that all connections are tight. c) Hold the BTU in one hand and push the cannula on the filled sample syringe through the BTU needleless injection site. d) Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk on the vacuum tube. e) Fill the vacuum tube to the desired volume. f) Repeat steps (d) and (e) according to the requirements for the patient's blood study.
2	According to hospital policy, discard the BTU, syringes, and cannulae after transferring the blood sample from the syringe into the vacuum tubes.

### 8.3 Drawing Blood Samples Using Method Two (Direct-Draw Method)

Step	Procedure
1	Swab the sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. <b>Note: Do not use acetone.</b>
2	To draw a blood sample, use the direct-draw unit. <b>CAUTION: Do not use a needle through the sampling site.</b> a) Using aseptic technique, peel open the pouch. b) Ensure that the cannula is securely tightened to the direct-draw unit housing. c) Position the sampling site so that it faces upward. d) Push the cannula of the direct-draw unit into the sampling site (see Figure 9 on page 67). e) Insert the selected vacuum tube into the open end of the direct-draw unit and push until the internal needle of the BTU has punctured the rubber disk on the vacuum tube. <b>CAUTION: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity.</b> f) Fill the vacuum tube to the desired volume. <b>Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.</b> g) Repeat steps (e) and (f) according to the requirements of the patient's blood study.

Step	Procedure
	h) When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw unit by the needleless cannula and <b>pull straight out</b> .  <b>CAUTION: Do not twist the BTU housing or remove it with the vacuum tube still attached.</b>
3	Discard the direct-draw unit after use according to hospital policy.

## 9.0 Re-infusing the Clearing Volume

Step	Procedure
1	Once the sample has been drawn, swab the sampling site to ensure the removal of any excess blood remaining on the sampling port. <b>Note: Do not use acetone.</b>
2	Open the proximal shut-off valve by turning the handle parallel to the tubing (see Figure 3 on page 65).
3	Push down on the plunger until the flexures lock into place in the fully closed position and all fluid has been re-infused into the line (see Figure 2 on page 65).  <b>Important: Recommended rate to push the reservoir plunger to the closed position is 1 ml/second.</b>  <b>Important: The clearing volume should not remain in the reservoir for longer than 3 minutes.</b>

## 10.0 Flushing the Line

Step	Procedure
1	Close the proximal shut-off valve and open the distal shut-off valve.
2	Open the pre-attached IV line, or attach a flush syringe to the pre-attached heparin lock.

Step	Procedure
3	Gently squeeze the reservoir plunger flexures together. Pull up on the reservoir plunger until the desired amount of flush solution is in the reservoir (see Figure 5 on page 65).
4	Close the distal shut-off valve and open the proximal shut-off valve.
5	Push down on the plunger until the flexures lock into place in the fully closed position and all fluid has been re-infused into the line (see Figure 2 on page 65).  <b>Important: Recommended rate to push the reservoir plunger to the fully closed position is 1 ml/second.</b>
6	Repeat steps 1 to 5 until the line is cleared of residual blood per hospital policy.
7	Open the distal shut-off valve to resume IV fluid supply.

**WARNING: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.**

## 11.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

## 12.0 MRI Safety Information

MR Safe

The VAMP closed blood sampling system is MR Safe.

**Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.**

## 13.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

This product is for single use only.

Visually inspect for breaches of packaging integrity prior to use.

## 14.0 Storage

Store in a cool, dry place.

## 15.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

## 16.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada  
(24 hours): ..... 800.822.9837  
Outside the U.S. and Canada  
(24 hours): ..... 949.250.2222  
In the UK: ..... 0870 606 2040 - Option 4  
In Ireland: ..... 01 821 1012 - Option 4

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

## 17.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

**Refer to the symbol legend at the end of this document.**

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.
**Edwards Lifesciences GmbH**
**Edwards Lifesciences LLC**
**Edwards Lifesciences B.V.**