Directory

English (en)	1	Português (pt)	28	Română (ro)	51
Français (fr)	4	Česky (cs)	31	Eesti (et)	54
Deutsch (de)	7	Magyar (hu)	33	Lietuvių (lt)	56
Español (es)	10	Polski (pl)	36	Latviešu (lv)	59
Italiano (it)	13	Slovensky (sk)	39	Türkçe (tr)	62
Nederlands (nl)	16	Norsk (no)	42	Русский (ru)	65
Dansk (da)	19	Suomi (fi)	45	Srpski (sr)	68
Svenska (sv)	22	Български (bg)	47	Hrvatski (hr)	71
Ελληνικά (el)	24				

Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Φιγυρι ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike	74
--	----

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola	78
---	----

English

VAMP

Closed Blood Sampling System

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use

For Single Use Only

For figures, see Figure 1 on page 74 through Figure 8 on page 77.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Edwards Lifesciences VAMP closed blood sampling system when attached to pressure monitoring lines, provides a safe and convenient method for the withdrawal of blood samples. The blood sampling system is designed for use with all disposable and reusable pressure transducers and for connection to central line catheters and arterial catheters. The VAMP closed blood sampling system is used for the drawing and retention of heparinized blood from the catheter or cannula within the line, allowing undiluted blood samples to be drawn from an in-line sampling site. At the completion of sample drawing, the mixed heparin and blood solution is reinfused into the patient to reduce fluid loss.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the

safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use of blood sampling technologies as part of their respective institutional guidelines.

VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

3.0 Indications for Use

For adult patients with medical conditions requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

4.0 Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

5.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Some models may contain phthalates, specifically DEHP [Bis (2-ethylhexyl) phthalate], which may

pose a risk of reproductive or developmental harm in pediatric patients, pregnant or nursing women.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

6.0 Instructions for Use

CAUTION: The use of lipids with the VAMP closed blood sampling system may compromise product integrity.

6.1 Equipment

- Flush device or flow controlling device (maximum 4 ml/hr flow rate)
- Disposable or reusable pressure transducer, if desired
- VAMP closed blood sampling system that contains one 5 ml VAMP reservoir with integral shut-off valve and one VAMP needleless sampling site

6.2 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP kit from the sterile package.
2	Attach the distal end with the female luer-lock connector to the transducer stopcock on the disposable transducer or reusable transducer dome (see Figure 1 on page 74).
3	All connections should be secure. Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.
4	Make sure that the reservoir plunger is down in the closed and locked position.

Step	Procedure
5	Deliver flush solution first through the transducer and out through the transducer vent port according to its manufacturer's instructions.
6	Replace all vented caps on the sideports of the stopcocks with non-vented caps.
7	To fill the VAMP kit, ensure that the shut-off valve is in the open position indicated by the handle being parallel to the tubing. Hold the reservoir by the plunger flexures. Orient the kit such that the reservoir is in the vertical position and the sampling site is above the reservoir at approximately 45°. In this position, gently squeeze and hold the plunger flexures together to slightly raise the seal to enlarge the fluid path and deliver flush solution slowly through the remainder of the kit. Once the kit is filled, release plunger flexures. CAUTION: Remove all air bubbles to reduce the risk of air emboli.
8	Mount the transducer either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder.
9	Mount the reservoir either on the patient's body per hospital procedure or on an IV pole using the appropriate clamp and holder (see Figure 2 on page 74).
10	Pressurize the IV solution bag. Flow rate will vary with pressure across the flush device.
11	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter (see Figure 1 on page 74).
12	Zero and calibrate the transducer according to its manufacturer's instructions.

7.0 Drawing Blood Samples

Note: Use aseptic technique.

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

Although a variety of techniques can be used for drawing samples, the following guidelines are provided as an aid to the clinician:

Two methods may be used to draw blood samples using the VAMP blood sampling system. Method one uses a sampling syringe with the VAMP needleless cannula and a blood transfer unit (BTU). Method two, direct-line sampling, uses a direct-draw blood transfer unit (BTU) with a VAMP needleless cannula.

7.1 Drawing Blood Samples Using Method One (Syringe and Cannula)

Step	Procedure
1	Gently squeeze the reservoir plunger flexures together. Pull up on the reservoir plunger until it stops and the reservoir has reached its 5 ml volume capacity (see Figure 3 on page 75). Note: If difficulties are experienced in drawing the clearing sample, check the catheter for possible occlusions or restrictions (e.g., positional lines). Important: Recommended time to pull reservoir plunger up to the fully open position is 3 to 5 s.
2	Once the clearing sample has been drawn, close the reservoir shut-off valve by turning the handle perpendicular to the tubing. This will further ensure that the sample drawn is from the patient and not the reservoir.
3	Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy. Note: Do not use acetone.
4	To draw a blood sample, use either a VAMP needleless cannula (packaged separately) and syringe, or a packaged VAMP needleless cannula preassembled onto a syringe. CAUTION: Do not use a needle through the sampling site. If an individually packaged VAMP needleless cannula is used: a) Using aseptic technique, peel open the cannula pouch. b) Pick up the cannula by the protective shield (see Figure 4 on page 75). c) Attach the cannula onto a selected luer-tip syringe by aligning the cannula luer-lock to the luer-tip on the syringe, and twisting until secure. d) Ensure that the syringe plunger is depressed to the bottom of the syringe barrel. e) Push the cannula into the VAMP needleless sampling site and then draw the required volume of blood into the syringe for sample (see Figure 5 on page 76). Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions. f) Remove the syringe and cannula from the sampling site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site.

Step	Procedure
	If a VAMP needleless cannula preassembled on a syringe is used: a) Using aseptic technique, peel open the pouch. b) Ensure that the cannula is securely tightened to the syringe. c) Ensure that the syringe plunger is depressed to the bottom of the syringe barrel. d) Push the cannula into the VAMP needleless sampling site and then draw the required volume of blood into the syringe for sample (see Figure 5 on page 76). Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions. e) Remove the syringe and cannula from the sampling site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site.
5	Once the sample has been drawn, open the shut-off valve by turning the handle 90° so that it is parallel to the tubing.
6	Push down on the plunger until the flexures lock in place in the fully closed position and all fluid has been reinfused into the line (see Figure 3 on page 75). Important: Recommended time to push the reservoir plunger to the fully closed position is 3 to 5 s. Important: The clearing sample should not remain in the reservoir for longer than 3 minutes.
7	Flush the line clear and swab the sampling site, ensuring removal of any excess blood left on the sampling port. Note: Do not use acetone. WARNING: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

Step	Procedure
8	To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU) (see Figure 6 on page 76). - Using aseptic technique, peel open the pouch. - Ensure that all connections are tight. - Hold the BTU in one hand and push the cannula on the filled sample syringe through the BTU needleless injection site. - Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk on the vacuum tube. - Fill the vacuum tube to the desired volume. - Repeat steps (d) and (e) according to the requirements for the patient's blood study. - Discard the BTU after transferring the blood sample from the syringe into the vacuum tubes.
9	Discard all syringes and cannulae after use according to hospital policy.

7.2 Drawing Blood Samples Using Method Two (Direct-Draw)

Step	Procedure
1	Gently squeeze the reservoir plunger flexures together. Pull up on the reservoir plunger until the plunger stops and the reservoir has reached its 5 ml volume capacity (see Figure 3 on page 75). CAUTION: If difficulties are experienced in drawing the clearing sample, check the catheter for possible occlusions or restrictions (e.g., positional lines). Important: Recommended time to pull reservoir plunger up to the fully open position is 3 to 5 s.
2	Once the clearing sample has been drawn, close the reservoir shut-off valve by turning the handle 90° so that it is perpendicular to the tubing. This will further ensure that the sample drawn is from the patient and not the reservoir.
3	Swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. Note: Do not use acetone.

Step	Procedure
4	To draw a blood sample, use the direct-draw blood transfer unit. CAUTION: Do not use a needle through the sampling site. - Using aseptic technique, peel open the pouch. - Ensure that the cannula is securely tightened to the BTU housing. - Position the patient's arm such that the sample site faces upward (see Figure 7 on page 77). Note: For VAMP closed blood sampling systems that contain small bore tubing the sample site may not be located on the patient's arm (see Figure 8 on page 77). - Push the cannula of the direct-draw BTU into the sampling site. - Insert the selected vacuum tube into the open end of the direct-draw BTU and push until the internal needle of the BTU has punctured the rubber disk on the vacuum tube. CAUTION: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity. - Fill the vacuum tube to the desired volume. CAUTION: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions. - Repeat steps (e) and (f) according to the requirements of the patient's blood study. - When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw BTU by the cannula and pull straight out. CAUTION: Do not twist the BTU housing or remove it with the vacuum tube still attached.
5	Discard the direct-draw BTU after use according to hospital policy.
6	Once the sample has been drawn, open the shut-off valve by turning the handle 90° so that the handle is parallel to the tubing.
7	Push down on the plunger until the flexures lock in place in the fully closed position and all fluid has been reinfused into the line (see Figure 3 on page 75). Important: Recommended time to push the reservoir plunger to the fully closed position is 3 to 5 s. Important: The clearing sample should not remain in the reservoir for longer than 3 minutes.

Step	Procedure
8	Flush the line clear and swab the sampling site, ensuring removal of any excess blood left on the sampling port. Note: Do not use acetone. WARNING: Laboratory values should correlate with the patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.

8.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

9.0 MRI Safety Information

MR MR Safe

The VAMP closed blood sampling system is MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

10.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize. This product is for single use only.

Visually inspect for breaches of packaging integrity prior to use.

11.0 Storage

Store in a cool, dry place.

12.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

13.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
 (24 hours): 800.822.9837
 Outside the U.S. and Canada
 (24 hours): 949.250.2222
 In the UK: 0870 606 2040 - Option 4
 In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

14.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences GmbH
Edwards Lifesciences LLC
Edwards Lifesciences B.V.