Directory		
English (en)	1	
Français (fr)	3	
Deutsch (de)	6	
Español (es)	9	
Italiano (it)	11	
Nederlands (nl)	14	
Dansk (da)	17	
Svenska (sv)	19	
Ελληνικά (el)	22	
Português (pt)	25	
Česky (cs)	27	
Magyar (hu)	30	
Polski (pl)	32	
Slovensky (sk)	35	
Norsk (no)	38	
Suomi (fi)	40	
Български (bg)	43	
Română (ro)	45	
Eesti (et)	48	
Lietuvių (lt)	50	
Latviešu (lv)	53	
Türkçe (tr)	55	
Русский (ru)	58	
Srpski (sr)	61	
Hrvatski (hr)	63	
Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνας ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунки ■ Slike ■ Slike 66		
Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarázat ■ Legenda symboli ■ Vysvětlivky k symbolom ■ Symbolforklaring ■ Merkien selitykset ■ Легенда на символите ■ Legendă de simboluri ■ Sümbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Условные обозначения ■ Legenda sa simbolima ■ Legenda simbola 70		

English

VAMP Flex Closed Blood Sampling System

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific region.

Instructions for Use

For Single Use Only

For figures, see Figure 1 through Figure 4 on page 66.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The Edwards Lifesciences VAMP Flex closed blood sampling system provides a safe and convenient method for the withdrawal of blood samples from pressure monitoring lines. The blood sampling system is designed for use with disposable and reusable pressure transducers and for connection to central line catheters, venous, and arterial catheters where the system can be flushed clear after sampling.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended to be used by medical professionals who have been trained in the safe use of hemodynamic technologies and the clinical use

of blood sampling technologies as part of their respective institutional guidelines.

VAMP blood conservation technology reduces unnecessary blood loss and risk of infection. Additional risks include blood loss, blood splatter, embolus, thrombosis, adverse reaction to device materials, tissue trauma/injury, systemic infection, and/or hemolysis.

2.0 Intended Use/Purpose

VAMP closed blood sampling system is intended to be used only for blood withdrawal.

3.0 Indications for Use

For adult patients with medical conditions requiring periodic withdrawal of blood samples from arterial and central line catheters, including peripherally inserted central catheters and central venous catheters, which are attached to pressure monitoring lines.

4.0 Contraindications

Not to be used without an attached flush device or flow controlling device when used for arterial applications.

There are no absolute contraindications when used for venous applications.

5.0 Warnings

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Some models may contain phthalates, specifically DEHP [Bis (2-ethylhexyl) phthalate], which may pose a risk of reproductive or developmental harm in

pediatric patients, pregnant or nursing women.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

6.0 Instructions for Use

CAUTION: The use of lipids with the VAMP Flex closed blood sampling system may compromise product integrity.

6.1 Equipment

- Flush device or flow controlling device (maximum 4 ml/hr flow rate)
- Disposable or reusable pressure transducer, if desired
- A VAMP Flex closed blood sampling system that contains one 10 ml reservoir with shut-off valve and VAMP needleless sampling site(s)

6.2 Setup

Step	Procedure
1	Using aseptic technique, remove the VAMP Flex kit from the sterile package.
2	If the VAMP Flex kit does not include a preconnected TruWave disposable pressure transducer, remove the protective caps and attach the distal end female luer port of the 4-way stopcock to the male luer of a transducer, or other fluid flush device (see Figure 1 on page 66).
3	All connections should be secure. Note: Wet connections promote over-tightening by lubricating the fittings. Over-tightened connections may result in cracks or leaks.

Step	Procedure
4	Connect the kit to an IV fluid source and without pressurizing the fluid, begin to gravity fill the kit first through the transducer and out through the transducer vent port according to transducer manufacturer's instructions.
5	Replace all vented caps on the sideports of the stopcocks with non-vented caps.
6	To fill the VAMP Flex reservoir, first turn the stopcock off to the open or patient end of the tubing (see Figure 2 on page 67). Slowly fill the reservoir to the full 10 ml position. Turn the stopcock off to the transducer/IV bag (see Figure 3 on page 68). Orient the reservoir such that any air bubbles will flow out of the syringe, and slowly reinfuse the fluid. Repeat as necessary to remove air bubbles. CAUTION: Remove all air bubbles to reduce the risk of air emboli and reduce loss of pressure signal.
7	Turn the stopcock off to the VAMP Flex reservoir (see Figure 4 on page 69), and continue filling the line using gravity fill. Orient the line such that the sampling site(s) are above the reservoir at approximately 45°. De-bubble the sample site(s) if needed.
8	Pressurize the IV solution bag per hospital policy and transducer manufacturer's instructions. Flow rate will vary with pressure gradient across the flush device.
9	Securely connect the proximal end of the kit with the male luer-lock connector to the pre-filled catheter.
10	Zero the transducer according to the transducer manufacturer's instructions.

6.3 Clearing IV Fluid in Preparation for Drawing Blood Samples

Important: A minimum clearing volume of two times the dead space should be achieved. Additional clearing volume may be required for coagulation studies.

Step	Procedure
1	Turn the stopcock off to the transducer/IV bag (see Figure 3 on page 68). Fill the reservoir by drawing the plunger. Draw the appropriate clearing volume according to hospital policy. Note: If difficulties are experienced in drawing the clearing volume, check the catheter and stopcock handle OFF position for possible occlusions or restrictions. Note: Recommended time to fill the reservoir is approximately 1 second for each ml of clearing volume.
2	Once the clearing volume has been drawn, turn the stopcock off to the patient and catheter (see Figure 2). This will further ensure that the sample drawn is from the patient and not the reservoir.

Step	Procedure
3	Draw the sample according to one of the methods outlined in the Drawing Blood Samples section below.
4	Once the sample has been drawn, turn the handle off to the transducer/IV bag (see Figure 3 on page 68).
5	Reinfuse the clearing volume back to the patient by returning the reservoir plunger to the closed position. Turn the stopcock handle to the monitoring position (see Figure 4 on page 69). Note: If difficulties are experienced in reinfusing the clearing volume, check the catheter and stopcock handle OFF position for possible occlusions or restrictions. Note: Recommended time to push the reservoir plunger to the fully closed position is approximately 1 second for each ml of clearing volume. CAUTION: The clearing sample should not remain in the reservoir for longer than 2 minutes.
6	Flush the reservoir, sample site(s), and line clear of residual blood.
7	If blood has been drawn into the VAMP Flex reservoir, it will be necessary to clear the reservoir. To accomplish this, turn the stopcock off to the patient/catheter (see Figure 2 on page 67). Partially fill the reservoir with 1 to 2 ml of fluid. If the IV bag is properly pressurized, activating the flush device will automatically cause the reservoir to fill. Turn the stopcock off to the transducer/IV bag (see Figure 3 on page 68) and slowly reinfuse into the line. Repeat as necessary until all blood is removed. Note: For optimal clearing of the reservoir, it is recommended to repeat this clearing process multiple times with small volume of fluid (i.e., 1 to 2 ml). After all blood has been removed, turn the stopcock to the monitoring position (see Figure 4 on page 69).
8	Swab the sampling site(s), ensuring removal of any excess blood left on the sampling port. Note: Do not use acetone.

6.4 Drawing Blood Samples

Two methods may be used to draw blood samples using the VAMP Flex closed blood sampling system.

6.4.1 Syringe Method for Drawing Blood Samples

Step	Procedure
1	Prepare the VAMP Flex system for drawing a blood sample by drawing a clearing volume per the method outlined above.

Step	Procedure
2	Once the line is clear, swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending upon hospital policy. Note: Do not use acetone.
3	To draw a blood sample, use either a preassembled packaged VAMP needleless cannula and syringe or a VAMP needleless cannula (packaged separately) and separate syringe. CAUTION: Do not use a needle through the sampling site.
4	Ensure that the syringe plunger is depressed to the bottom of the syringe barrel.
5	Push the cannula into the VAMP needleless sampling site and then draw the required volume of blood into the syringe. Note: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.
6	Remove the syringe and cannula from the sampling site by pulling straight out . CAUTION: Do not twist the syringe out of the sampling site.
7	Once the sample has been drawn, reinfuse the clearing volume according to the method described above. WARNING: Laboratory values should correlate with patient's clinical manifestations. Verify accuracy of laboratory values before instituting therapy.
8	To transfer the blood sample from the syringe to vacuum tubes, use the blood transfer unit (BTU). a) Using aseptic technique peel open the pouch. b) Ensure that all connections are tight. c) Hold the BTU in one hand and push the cannula on the filled sample syringe through the BTU needleless injection site. d) Insert the selected vacuum tube into the opening of the BTU until the internal needle has punctured the rubber disk of the vacuum tube. e) Fill the vacuum tube to the desired volume. f) Repeat steps (d) and (e) as required to fill additional tubes. g) Discard the BTU after transferring the blood sample from the syringe into the vacuum tubes.
9	Discard all syringes and cannulas after use according to hospital policy.

6.4.2 Direct-Draw Method for Drawing Blood Samples

Step	Procedure
1	Prepare the VAMP Flex reservoir for drawing a blood sample by drawing a clearing volume per the method outlined above.
2	Once the line is clear, swab the VAMP needleless sampling site with disinfectant such as alcohol or betadine, depending on hospital policy. Note: Do not use acetone.
3	To draw a blood sample, use the direct-draw unit. CAUTION: Do not use a needle through the sampling site. a) Using aseptic technique, peel open the pouch. b) Ensure that the cannula is securely tightened to the direct-draw unit. c) Push the cannula of the direct-draw unit into the sampling site. d) Insert the selected vacuum tube into the open end of the direct-draw unit and push until the internal needle of the direct-draw unit has punctured the rubber disk on the vacuum tube and fill to the desired volume. CAUTION: To prevent backflow of vacuum tube contents (including air) from entering fluid path, remove vacuum tube before reaching maximum fill capacity. CAUTION: If difficulties are experienced in drawing the sample, check the catheter for possible occlusions or restrictions.

Step	Procedure
	e) Repeat step (d) as required to fill additional vacuum tubes. f) When the last sample has been drawn, remove the vacuum tube first and then grasp the direct-draw unit by the cannula and pull straight out. CAUTION: Do not twist the direct-draw unit housing or remove it with the vacuum tube still attached. g) Discard the direct-draw unit after use according to hospital policy.
4	Once the sample has been drawn, reinfuse the clearing volume according to the method described above.

7.0 Routine Maintenance

Since kit configurations and procedures vary according to hospital preferences, it is the responsibility of the hospital to determine exact policies and procedures.

8.0 MRI Safety Information

MR Safe

The VAMP Flex closed blood sampling system is MR Safe.

Precaution: Follow the conditions for safe scanning for any accessory devices (e.g., disposable transducers or reusable transducers) that are connected to the VAMP Flex closed blood sampling system. If the MR safety status for the accessory devices is not known, assume they are MR Unsafe and do not allow them to enter the MR environment.

9.0 How Supplied

Contents sterile and fluid path nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize. This product is for single use only.

Visually inspect for breaches of packaging integrity prior to use.

10.0 Storage

Store in a cool, dry place.

11.0 Shelf Life

The shelf life is as marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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
Figure 4: Stopcock positioned off to VAMP Flex (for pressure measurement)

- **Figure 4 : Robinet d'arrêt en position fermée vers le système VAMP Flex (pour mesurer la pression)**
 - **Abbildung 4: Absperrhahn in Richtung VAMP Flex geschlossen (zur Druckmessung)**
 - **Figura 4: Llave de paso en posición cerrada a VAMP Flex (para medición de la presión)**
 - **Figura 4: Rubinetto di arresto chiuso verso VAMP Flex (per la misurazione della pressione)**
- **Afbeelding 4: Het kraantje naar de VAMP Flex staat op de gesloten positie (voor drukmeting)**
 - **Figur 4: Stophane placeret ud for VAMP Flex (til trykmåling)**
 - **Figur 4: Kranen i avstängt läge mot VAMP Flex (för tryckmätning)**
- **Εικόνα 4: Στρόφιγγα στη θέση OFF (Κλειστή) προς το VAMP Flex (για μέτρηση της πίεσης)**
- **Figura 4: Torneira de passagem fechada para VAMP Flex (para a medição de pressão)**
- **Obrázek 4: Umístění uzavíracího kohoutu na jednotku VAMP Flex (k měření tlaku)**
 - **4. ábra: Zárócsap zárt helyzetben a VAMP Flex felé (nyomásméréshez)**
- **Rysunek 4: Zawór odcinający ustawiony w stronę systemu VAMP Flex (do pomiaru ciśnienia)**
- **Obrázok 4: Ventil v polohe uzatvorenia (OFF) smerom k systému VAMP Flex (na meranie tlaku)**
 - **Figur 4: Stoppekran satt til av for VAMP Flex (for trykkmåling)**
 - **Kuva 4: Hana suljettuna VAMP Flex -sarjaan päin (paineen mittausta varten)**
- **Фигура 4: Спирателно кранче, позиционирано до VAMP Flex (за измерване на налягане)**
- **Figura 4: Robinetul de închidere în poziția închis spre VAMP Flex (pentru măsurarea presiunii)**
 - **Joonis 4. kraani käepide VAMP Flex -i asendis (rõhu mõõtmiseks)**
 - **4 pav. Užsuktas „VAMP Flex“ atšakos čiupas (matuojant spaudimą)**
 - **4. attēls. Noslēgkrāns aizgriezts VAMP Flex virzienā (spiediena mērīšanai)**
- **Şekil 4: Musluk VAMP Flex cihazına doğru kapalı şekilde konumlandırılmış (basınç ölçümü için)**
- **Рис. 4. Запорный кран в выключенном положении для VAMP Flex (в состоянии для измерения давления)**
 - **Slika 4: Slavinica zatvorena, okrenuta ka VAMP Flex sistemu (za merenje pritiska)**
 - **Slika 4: Zaporni ventil okrenut prema kompletu VAMP Flex (za mjerenje tlaka)**

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli

	English	Français	Deutsch	Español	Italiano
	Single use	Ne pas réutiliser	Zum einmaligen Gebrauch	Un solo uso	Monouso
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE Mark)	Conformité européenne (marquage CE)	Conformité Européenne (CE-Kennzeichnung)	Conformité Européenne (Marca CE)	Conformité Européenne (Marchio CE)
	Authorized representative in the European Community/ European Union	Représentant autorisé dans la Communauté européenne/l'Union européenne	Bevollmächtigter in der Europäischen Gemeinschaft/ Europäischen Union	Representante autorizado en la Comunidad Europea/ Unión Europea	Rappresentante autorizzato nella Comunità Europea/ Unione Europea
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
	Importer	Importateur	Importeur	Importador	Importatore
	Contains or presence of Phthalates: bis (2-ethylhexyl) phthalate (DEHP)	Traces ou présence de phtalates : phtalate de bis (2-éthylhexyle) (DEHP)	Enthält Phthalate oder kann Spuren von Phthalaten enthalten: Bis(2-ethylhexyl) phthalat (DEHP)	Contiene o hay presencia de ftalatos: Di[2-etilhexil] ftalato (DEHP)	Contiene o sono presenti ftalati: di (2-etilesil) ftalato (DEHP)
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto
	Quantity	Quantité	Menge	Cantidad	Quantità
	Caution	Avertissement	Vorsicht	Aviso	Attenzione
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web
	Consult instructions for use or consult electronic instructions for use	Consultez le mode d'emploi sur papier ou le mode d'emploi au format électronique	Gebrauchsanweisung beachten oder elektronische Gebrauchsanweisung beachten	Consulte las instrucciones de uso o las instrucciones de uso electrónicas	Consultare le istruzioni per l'uso o le istruzioni per l'uso elettroniche
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene

	English	Français	Deutsch	Español	Italiano
	Use by	Utiliser avant	Verwendbar bis	Caducidad	Utilizzare entro
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno
	MR Safe	Aucun risque en milieu RM	MR-sicher	Seguro para RM	Compatibile con RM
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil; manipular con cuidado	Fragile, maneggiare con cura
	This way up	Flèches indiquant le haut du carton	Oben	Este lado arriba	Alto
	Single sterile barrier system	Système de barrière stérile unique	Einfaches Sterilbarrieresystem	Un único sistema de barrera estéril	Sistema a singola barriera sterile
	Unique device identifier	Identifiant unique du dispositif	Einmalige Produktkennung	Identificador único del dispositivo	Identificatore univoco del dispositivo
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello
Note: Not all symbols may be included in the labeling of this product. ■ Remarque : il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. ■ Hinweis: Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. ■ Nota: Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. ■ Nota: alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto.					

Lijst metsymbolen ■ Symbolforklaring ■ Symbolförklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Voor eenmalig gebruik	Til engangsbrug	För engångsbruk	Μίας χρήσης	Uso único
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-markering)	Conformité Européenne (CE-mærke)	Conformité Européenne (CE-märke)	Conformité Européenne (σήμανση CE)	Conformité Européenne (Marcação CE)
	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap/ Europese Unie	Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union	Auktoriserad representant inom Europeiska gemenskapen/Europeiska unionen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/ Ευρωπαϊκή Ένωση	Representante autorizado na Comunidade Europeia/ União Europeia
	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing	Må ikke anvendes, hvis emballagen er beskadiget, og se brugsanvisningen	Använd inte om förpackningen är skadad och se bruksanvisningen	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης	Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
	Importeur	Importør	Importör	Εισαγωγέας	Importador
	Bevat ftalaten of ftalaten zijn aanwezig: bis (2-ethylhexyl)ftalaat (DEHP)	Indeholder ftalater: bis (2-ethylhexyl) ftalat (DEHP)	Innehåller eller har spår av ftalater: bis(2-etylhexyl) ftalat (DEHP)	Περιέχει φθαλικές ενώσεις ή ίχνη αυτών: φθαλικός δι (2-αιθυλοεξυλο) εστέρας (DEHP)	Contém ou estão presentes ftalatos: ftalato de bis (2-etil-hexilo) (DEHP)
	Fabrikant	Producent	Tillverkare	Κατασκευαστής	Fabricante
	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής	Data de fabrico
	Lotnummer	Partinummer	Lotnummer	Αριθμός Παρτίδας	Número de lote
	Hoeveelheid	Mængde	Antal	Ποσότητα	Quantidade
	Let op	Forsigtig	Var försiktig	Προσοχή	Aviso
	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web	Consultar as instruções de utilização no site
	Raadpleeg de gebruiksaanwijzing of de elektronische gebruiksaanwijzing	Se brugsanvisningen, eller se elektronisk brugsanvisning	Se bruksanvisningen eller den elektroniska bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης ή συμβουλευτείτε τις ηλεκτρονικές οδηγίες χρήσης	Consultar as instruções de utilização ou consultar as instruções de utilização eletrônicas
	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε	Não voltar a esterilizar
	Gesteriliseerd met behulp van ethyleenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου	Esterilizado com óxido de etileno

	Nederlands	Dansk	Svenska	Ελληνικά	Português
	Te gebruiken tot	Anvendes inden	Använd före	Ημερομηνία λήξης	Validade
	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο	Não pirogênico
	MRI-veilig	MR-sikker	MR-säker	Ασφαλές σε περιβάλλον μαγνητικής τομογραφίας (MR)	Utilização segura em ambiente de RM
	Medisch hulpmiddel	Medicinsk udstyr	Medicinteknisk produkt	Ιατροτεχνολογικό προϊόν	Dispositivo médico
	Gesteriliseerd met behulp van straling	Steriliseret ved brug af bestråling	Steriliserad med strålning	Αποστειρωμένο με χρήση ακτινοβολίας	Esterilizado por irradiação
	Breekbaar, voorzichtigheid geboden	Skrøbelig, skal håndteres forsigtigt	Ömtåligt, hanteras försiktigt	Εύθραστο: χειριστείτε με προσοχή	Frágil, manusear com cuidado
	Deze kant boven	Denne side opad	Denna sida upp	Αυτή η πλευρά προς τα πάνω	Este lado para cima
	Enkelvoudig steriel barrièresysteem	Enkelt sterilt barrièresystem	Enkelt sterilt barriärsystem	Σύστημα μονού αποστειρωμένου φραγμού	Sistema de barreira única esterilizada
	Unieke instrumentidentificatie	Unik udstyrsidentifikation	Unik produktidentifisering	Αποκλειστικό αναγνωριστικό ιατροτεχνολογικού προϊόντος	Identificador único de dispositivo
	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο	Guardar num local fresco e seco
	Modelnummer	Modelnummer	Modellnummer	Αριθμός μοντέλου	Número do modelo
Opmerking: Het label van dit product bevat misschien niet alle symbolen. ■ Bemærk: Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. ■ Obs! Alla symboler kanske inte används för märkning av denna produkt. ■ Σημείωση: Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος. ■ Nota: poderão não estar incluídos todos os símbolos na rotulagem deste produto.					

Legenda se symboly ■ Jelmagyarázat ■ Legenda simboli ■ Vysvetlivky k symbolom ■ Symbolforklaring

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	Suomi	Български	Română	Eesti	Lietuvių
	Kertakäyttöinen	За еднократна употреба	De unică folosință	Ühekordseks kasutamiseks.	Vienkartinio naudojimo
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE-merkintä)	Conformité Européenne (CE маркировка)	Conformité Européenne (marcă CE)	Conformité Européenne (CE-märk)	Conformité Européenne (CE žymė)
	Valtuutettu edustaja Euroopan yhteisössä / Euroopan unionissa	Упълномощен представител в Европейската общност / Европейския съюз	Reprezentant autorizat în Comunitatea Europeană / Uniunea Europeană	Volitatud esindaja Euroopa Ühenduses / Euroopa Liidus	Igaliotasis atstovas Europos Bendrijose / Europos Sąjungoje
	Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare	Ärge kasutage, kui pakend on kahjustatud, ja palun lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas
	Maahantuoja	Вносител	Importator	Importija	Importuotojas
	Sisältää tai sisältää ftalaatteja: bis(2-etyyliheksyyli)ftalaatti (DEHP)	Съдържание или наличие на фталати: бис(2-етилхексил) фталат (DEHP)	Conține sau prezintă ftalați: di (2-etilhexil) ftalat (DEHP)	Sisaldab ftlaate: bis(2-etyülheksüül)ftalaat (DEHP)	Sudėtyje yra ftalatų: bis(2-etilheksil) ftalatas (DEHP)
	Valmistaja	Производител	Producător	Tootja	Gamintojas
	Valmistuspäivämäärä	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Eränumero	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Määrä	Количество	Cantitate	Kogus	Kiekis
	Tärkeä huomautus	Внимание	Atenție	Ettevaatust	Perspėjimas
	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web	Palun lugege veebisaidil olevaid kasutusjuhiseid	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje
	Katso käyttöohjeet tai sähköiset käyttöohjeet	Направете справка с инструкциите за употреба или електронните инструкции за употреба	Consultați instrucțiunile de utilizare sau consultați instrucțiunile de utilizare în format electronic	Tutvuge paberkandjal või elektroonilise kasutusjuhiseiga	Žr. naudojimo instrukcijas arba elektronines naudojimo instrukcijas
	Älä steriloitu uudelleen	Не стерилизирайте повторно	A nu se reesteriliza	Ärge reesteriliseerige	Kartotinai nesterilizuokite
	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu

	Suomi	Български	Română	Eesti	Lietuvių
	Käytettävä ennen	Срок на годност	Valabil până la	Kõlblikkusaeg	Naudoti iki
	Pyrogeeniton	Непирогенно	Non pirogen	Mittepürogeenne	Nepirogeniškas
	Sopii magneettikuvaukseen	Безопасно при MR	Sigur în utilizarea la RM	Ohutu magnetresonantstomograafias	MR saugus
	Lääkinnällinen laite	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos priemonė
	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere	Steriliseeritud kiiritamist kasutades	Sterilizuota švitinant
	Särkyvä, käsiteltävä varoen	Чупливо, да се бори внимателно	Fragil, a se manipula cu atenție	Õrn, käsitseda ettevaatlikult	Trapus, elgtis atsargiai
	Tämä puoli ylöspäin	С тази страна нагоре	Cu partea aceasta în sus	See pool üles	Šia puse į viršų
	Yksinkertainen steriili sulkujärjestelmä	Единична стерилна бариерна система	Sistem cu barieră sterilă individuală	Ühekordne steriilne kattesüsteem	Vieno sterilaus barjero sistema
	Laitteen yksilöivä tunnus	Уникален идентификатор на изделието	Identificator unic al dispozitivului	Seadme kordumatu identifitseerimistunnus	Unikalasis priemonės identifikatorius
	Säilytettävä kuivassa ja viileässä	Да се съхранява на хладно, сухо място	Stocați într-un loc rece și uscat	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje
	Mallinnumero	Номер на модела	Număr model	Mudeli number	Modelio numeris
Huomautus: kaikkia symboleja ei välttämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. ■ Забележка: Възможно е не всички символи да са включени в етикета на този продукт. ■ Notă: este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. ■ Märkus. Kõik sümbolid ei pruugi antud toote etiketil esineda. ■ Pastaba. Šio gaminio etiketėje gali būti pateikti ne visi simboliai.					

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Vienreizējai lietošanai	Tek kullanımlık	Только для однократного использования	Jednokratna upotreba	Za jednokratnu upotrebu
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)	Conformité Européenne (oznaka CE)
	Pilnvarotais pārstāvis Eiropas Kopienā/Eiropas Savienībā	Avrupa Birliği'nde/Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе/Европейском союзе	Ovlašćeni predstavnik u Evropskoj zajednici / Evropskoj Uniji	Ovlašteni predstavnik u Europskoj zajednici / Europskoj uniji
	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristite ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Importētājs	İthalatçı	Импортер	Uvoznik	Uvoznik
	Vērojama ftalātu klātbūtne vai satur ftalātus: bis(2-etilheksil) ftalāts (DEHP)	Ftalat içerir veya Ftalat mevcuttur: bis (2-etilheksil) ftalat (DEHP)	Содержит фталаты или изготовлено из фталатов: ДЭГФ (диэтилгексилфталат)	Sadrži ili ima prisustva ftalata: bis (2-etilheksil) ftalat (DEHP)	Sadrži ftalate ili su prisutni ftalati: bis(2-etilheksil)-ftalat (DEHP)
	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Daudzums	Miktar	Количество	Količina	Količina
	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na web-stranici
	Skatiet lietošanas instrukcijas vai elektroniskās lietošanas instrukcijas	Kullanım talimatlarına veya elektronik kullanım talimatlarına başvurun	См. инструкции по применению либо электронные инструкции по применению	Pogledajte uputstvo za upotrebu ili elektronsko uputstvo za upotrebu	Pogledajte upute za upotrebu ili elektroničke upute za upotrebu
	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	Nemojte ponovno sterilizirati
	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidadom	Sterilizirano etilen-oksidadom

	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Upotrebljivo do	Upotrijebiti do
	Nepirogēns	Nonpirojeniktir	Апирогенно	Apirogeno	Nepirogeno
	Drošs lietošanai MR vidē	MR için Güvenli	Безопасно при проведении МРТ	Bezbedno za MR	Sigurno kod pregleda MR-om
	Medicīniska ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	Medicinski uređaj
	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom	Lomljivo, pažljivo rukovati
	Ar šo pusi uz augšu	Bu taraf yukarı	Верхняя сторона	Ova strana mora da bude prema gore	Ova strana prema gore
	Viena sterilitātes aizsargslāņa sistēma	Tek steril bariyer sistemi	Одинарная барьерная система для стерилизации	Sistem sa jednom sterilnom barijerom	Sustav jedne sterilne barijere
	Unikāls ierīces identifikators	Benzersiz cihaz tanımlayıcısı	Уникальный идентификатор устройства	Jedinstveni identifikator sredstva	Jedinstveni identifikator proizvoda
	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	Pohranite na hladnom, suhom mjestu
	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
Piezīme. Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. ■ Not: Bu ürünün etiketinde tüm semboller yer almamış olabilir. ■ Примечание. На этикетках данного изделия могут указываться не все символы. ■ Napomena: Na ambalaži proizvoda se možda ne nalaze svi simboli. ■ Napomena: na oznakama ovog proizvoda ne moraju se nalaziti svi simboli.					

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Made in Dominican Republic
CCT Critical Care Technologies S.R.L.
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

08/24
10014983003 B
© Copyright 2024, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences GmbH