Edwards

Directory

English (en)	1	Português (pt)	54	Română (ro)	101
Français (fr)	7	Česky (cs)	60	Eesti (et)	107
Deutsch (de)	13	Magyar (hu)	66	Lietuvių (lt)	112
Español (es)	19	Polski (pl)	72	Latviešu (lv)	118
Italiano (it)	25	Slovensky (sk)	78	Türkçe (tr)	124
Nederlands (nl)	31	Norsk (no)	84	Русский (ru)	130
Dansk (da)	37	Suomi (fi)	89	Srpski (sr)	136
Svenska (sv)	42	Български (bg)	95	Hrvatski (hr)	142
Ελληνικά (el)	48				

Figures ■ Figures ■ Abbildungen ■ Figuras ■ Figure ■ Afbeeldingen ■ Figurer ■ Figurer ■ Εικόνες ■ Figuras ■ Obrázky ■ Ábrák ■ Rysunki ■ Obrázky ■ Figurer ■ Kuvat ■ Фигури ■ Figuri ■ Joonised ■ Paveikslėliai ■ Attēli ■ Şekiller ■ Рисунок ■ Slike ■ Slike

Symbol Legend ■ Légende des symboles ■ Zeichenerklärung ■ Significado de los símbolos ■ Legenda dei simboli ■ Lijst met symbolen ■ Symbolforklaring ■ Symbolforklaring ■ Υπόμνημα συμβόλων ■ Legenda de símbolos ■ Legenda se symboly ■ Jelmagyarozat ■ Legenda symboli ■ Vysvetlivky k symbolom ■ Symbolforklaring ■ Merkkien selitykset ■ Легенда на символите ■ Legenda de simboluri ■ Sūbolite seletus ■ Simbolių paaiškinimas ■ Simbolu skaidrojums ■ Sembol Açıklamaları ■ Основные обозначения ■ Legenda sa simbolima ■ Legenda simbola

English

Swan-Ganz

Thermodilution Pacerport Catheter: 931F75
Thermodilution A-V Pacerport Catheter: 991F8
Oximetry TD Catheter: 631F55N
Oximetry Pacerport TD Catheter: 780F75M
VIP Oximetry TD Catheter: 782F75M

The devices described herein may not all be licensed in accordance with Canadian law or approved for sale in your specific country.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

Model 931F75 is illustrated in Figure 1 on page 148 Model 780F75M is illustrated in Figure 2 on page 149. The models listed above contain some but not all of the features shown.

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

TD (thermodilution) Pacerport catheter. Please refer to Figure 2 on page 149 and Figure 3 on page 150. Figure 1 on page 148, Figure 3 on page 150 and Figure 4 on page 151. Oximetry TD: Please refer to Figure 2 on page 149 and Figure 4 on page 151.

1.0 Description

The device is used by medical professionals who have been trained in safe use of invasive hemodynamic technologies and clinical usage of pulmonary artery catheters as part of their respective institutional guidelines.

The Swan-Ganz catheters are flow-directed pulmonary artery catheters used to monitor hemodynamic pressures. The distal (pulmonary artery) port on the Swan-Ganz catheters allow sampling of mixed venous blood for the assessment of oxygen transport balance and the calculation of derived parameters such as oxygen consumption, oxygen utilization coefficient, and intrapulmonary shunt fraction when used with compatible hemodynamic monitor and accessories.

The intravascular catheter is inserted through the central vein into the right side of the heart and is advanced towards the pulmonary artery. Route of insertion can be internal jugular, femoral, antecubital and brachial veins. The body parts in contact are atrium, ventricles, pulmonary artery and circulatory system.

Edwards, Edwards Lifesciences, the stylized E logo, Chandler, CO-Set, CO-Set+, Flex-Tip, Hi-Shore, Intro-Flex, Pacerport, Swan, Swan-Ganz, and VIP are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

The device is intended for use in adult critically ill or surgical patient populations. The device has not yet been tested in the pediatric population or in pregnant or lactating women.

1.1 Model 931F75

The Swan-Ganz thermodilution Pacerport catheter (Model 931F75) allows insertion and placement of the Chandler transluminal V-pacing probe (Model D98100). The Swan-Ganz thermodilution Pacerport catheters (Model 931F75) are intended for use in patients who require hemodynamic monitoring when the need for temporary transvenous pacing is anticipated.

The Pacerport catheter (Model 931F75) right ventricular (RV) lumen terminates 19 cm from the tip, and is used for the insertion of Model D98100 Chandler transluminal V-pacing probes into the right ventricle when the tip of the catheter is in the pulmonary artery. When the pacing probe is not inserted (in Model 931F75) the RV lumen may be used for infusing solutions. The proximal lumen may be used for pressure monitoring and bolus injection for cardiac output determinations by the thermodilution technique.

1.2 Model 991F8

The Swan-Ganz thermodilution A-V Pacerport catheter (Model 991F8) allows insertion and placement of the Chandler transluminal V-pacing probe (Model D98100) and Flex-Tip A-pacing probe (Model D98500).

The thermodilution A-V Pacerport catheter (Model 991F8) right ventricular (RV) lumen which terminates 19 cm from the tip is provided for the insertion of a Model D98100 Chandler transluminal V-pacing probe into the right ventricle when the tip of the catheter is in the pulmonary artery. A right atrial (RA) lumen that terminates 30 cm from the tip is provided for the insertion of a Model D98500 Flex-Tip transluminal A-pacing probe into the right atrium when the tip of the catheter is in the pulmonary artery. When the pacing probes are not inserted (in Model 991F8), the ventricular and atrial lumens may be used for right ventricular or atrial pressure monitoring, or infusing solutions.

1.3 Models 631F55N, 780F75M, 782F75M

Model 631F55N Swan-Ganz oximetry TD catheters enable monitoring of cardiac output, and continuous mixed venous oxygen. Models 782F75M Swan-Ganz VIP oximetry TD catheters enable monitoring of cardiac output, and provide an additional VIP lumen that allows for continuous infusion. The VIP lumen terminates at a port located 31 cm from the distal tip. The VIP lumen provide direct access to the right atrium or vena cava and allows for continuous infusion of solutions, pressure monitoring or blood sampling.

Model 780F75M Swan-Ganz oximetry Pacerport TD catheter is intended for use in patients who require hemodynamic monitoring when temporary transvenous pacing is anticipated. The oximetry Pacerport catheter's right ventricular (RV) lumen terminates 19 cm from the tip and is used for the insertion of a model D98100 Chandler transluminal V-pacing probe into the right ventricle when the tip of the catheter is in the pulmonary artery. When the pacing probe is not inserted, the right ventricular lumen may be used for right ventricular pressure monitoring or infusing solutions.

Mixed venous oxygen saturation is monitored by fiberoptic reflectance spectrophotometry. The amount of light absorbed, refracted, and reflected depends on the relative amounts of oxygenated and deoxygenated hemoglobin in the blood. The oximetry lumen terminates at the distal tip. This lumen contains the fibers that transmit the light to the pulmonary artery for measurement of mixed venous oxygen saturation. The proximal injectate lumen terminates at a port located 30 cm from the distal tip for models 780F75M and 782F75M; it terminates at 15 cm for model 631F55N. When the distal tip is located in the pulmonary artery, the proximal injectate port will reside in the right atrium or vena cava, allowing for bolus cardiac output injections, right atrial pressure monitoring, blood sampling, or infusion of solutions.

2.0 Intended Use / Purpose

Swan-Ganz catheters are pulmonary artery catheters intended for short-term use in the central circulatory system for patients who require intracardiac hemodynamic monitoring, blood sampling and infusing solution. Pacerport catheters can provide temporary transvenous pacing when used with a pacing probe. When used with a compatible monitoring platform and accessories, Swan-Ganz catheters offer a comprehensive hemodynamic profile to help clinicians assess cardiovascular function and guide treatment decisions.

3.0 Indications

Swan-Ganz catheters are diagnostic and monitoring tools used for hemodynamic monitoring of adult critically ill patients including but not limited to post major surgical recovery, trauma, sepsis, burns, pulmonary disease, pulmonary failure, cardiac disease including heart failure.

4.0 Contraindications

These products contain metallic components. Do NOT use in a Magnetic Resonance (MR) environment.

Patients with either recurrent sepsis, or with hypercoagulopathy, in which the catheter could serve as a focal point for septic or bland thrombus formation, should not be considered candidates for a balloon flotation catheter.

5.0 Warnings

No absolute contraindications to the use of flow-directed pulmonary artery catheters exist. However, a patient with a left bundle branch block may develop a right bundle branch block during catheter insertion, resulting in complete heart block. In such patients, temporary pacing modes should be immediately available.

Electrocardiographic monitoring during catheter passage is encouraged and is particularly important in the presence of either of the following conditions:

- Complete left bundle branch block, in which the risk of complete heart block is somewhat increased.
- Wolff-Parkinson-White syndrome and Ebstein's malformation, in which the risk of tachyarrhythmias is present.

Do not modify or alter the product in any way. Alteration or modification may affect product performance.

Air should never be used for balloon inflation in any situation where it may enter the arterial circulation, e.g., in all pediatric patients and in adults with suspected right to left intracardiac or intrapulmonary shunts. Bacteria-filtered carbon dioxide is the recommended inflation medium because of its rapid absorption into the blood in the event of balloon rupture within the circulation. Carbon dioxide diffuses through the latex balloon, diminishing the balloon's flow-directed capability after 2 to 3 minutes of inflation.

Do not leave the catheter in a permanent wedge position. Furthermore, avoid lengthy balloon inflation while the catheter is in a wedge position; this occlusive maneuver may result in pulmonary infarction.

This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Cleaning and resterilization will damage the integrity of the latex balloon. Damage may not be obvious during routine inspection.

6.0 Precautions

Failure of a balloon flotation catheter to enter the right ventricle or pulmonary artery is rare, but may occur in patients with an enlarged right atrium or ventricle particularly if the cardiac output is low or in the presence of tricuspid or pulmonic incompetence or pulmonary hypertension. Deep inspiration by the patient during advancement may also facilitate passage.

Clinicians using the device should be familiar with the device and understand its applications prior to use.

7.0 Recommended Equipment

WARNING: Compliance to IEC 60601-1 is only maintained when the catheter or probe (Type CF applied part, defibrillation proof) is connected to a patient monitor or equipment that has a Type CF defibrillation proof rated input connector. If attempting to use a third party monitor or equipment, check with the monitor or equipment manufacturer to ensure IEC 60601-1 compliance and compatibility with the catheter or probe. Failure to ensure monitor or equipment compliance to IEC 60601-1 and catheter or probe compatibility may increase the risk of electrical shock to the patient/operator.

TD Paceport catheter:

- Swan-Ganz thermodilution A-V Paceport catheter (Model 991F8) or Swan-Ganz thermodilution Paceport catheter (Model 931F75)
- Edwards Intro-Flex percutaneous sheath hemostasis valve introducer tray or kit, or single assembly.
- Chandler transluminal V-pacing probe (Model D98100)
- Flex-Tip transluminal A-pacing probe (Model D98500 - when using Model 991F8)
- A-V sequential or ventricular demand external pacemaker
- Any compatible cardiac output computer, compatible injectate probe, and connecting cable
- Sterile flush system and pressure transducers
- Bedside ECG and pressure monitor system
- Percutaneous sheath introducer and contamination shield

Oximetry TD Catheters:

- Swan-Ganz oximetry TD catheter or oximetry Paceport TD catheter or VIP oximetry catheter
- Percutaneous sheath introducer and contamination shield
- An Edwards Lifesciences oximetry monitor or compatible bedside module system (or any compatible cardiac output computer for measuring cardiac output by the bolus thermodilution method)
- Injectate temperature sensing probe
- Connecting cables
- Model OM-2 optical module
- Chandler transluminal V-pacing probe (model D98100) (for use with oximetry Paceport catheters model 780F75M only)
- Ventricular demand external pacemaker (for use with oximetry Paceport catheter model 780F75M and Chandler probe model D98100)
- Sterile flush system and pressure transducers
- Bedside ECG and pressure monitor system

In addition, the following items should be immediately available if complications arise during catheter insertion: antiarrhythmic drugs, defibrillator, respiratory assist equipment and a means for temporary pacing.

All Edwards oximetry catheters, except those with black oximetry connectors, are compatible with Edwards optical modules. Models ending with the letter

'M' are also compatible with Philips optical modules. The Catheter Factor required for *in vivo* calibration with Philips monitors is located on the top of the optical connector.

Monitor Set Up and Calibration for Mixed Venous Oxygen Saturation Monitoring

Precaution: *In vitro* calibration cannot be performed with catheter model 631F55N. For proper calibration, the catheter must be inserted into the patient and an *in vivo* calibration performed (see appropriate operations manual for *in vivo* calibration procedures). The compatible cardiac output computers can be calibrated prior to catheter insertion by performing an *in vitro* calibration. When performing an *in vitro* calibration, do so before preparing the catheter (i.e. flushing the lumens). The catheter tip must not get wet before an *in vitro* calibration is performed. An *in vivo* calibration is required if an *in vitro* calibration is not done.

Note: To avoid damaging the balloon, do not pull the balloon through the silicone gripper.

8.0 Catheter Preparation

Use aseptic technique.

Visually inspect for breaches of packaging integrity prior to use.

Note: Use of a protective catheter sheath is recommended.

Precaution: Avoid forceful wiping or stretching of the catheter during testing and cleaning as not to break the optical fibers and/or thermistor wire circuitry if present.

Step	Procedure
1	Flush catheter lumens with a sterile solution to ensure patency and to remove air.
2	Check balloon integrity by inflating it to the recommended volume. Check for major asymmetry and for leaks by submerging in sterile saline or water. Deflate balloon before insertion.
3	Connect the catheter's injectate and pressure monitoring lumens to the flush system and pressure transducers. Ensure that the lines and transducers are free of air.
4	Test the thermistor's electrical continuity before insertion (refer to computer operations manual for detailed information).

9.0 Insertion Procedure

WARNING: In some patients, the catheter may spontaneously wedge (with the balloon deflated) before positioning of the RV port in the right ventricle. Discontinue advancing the catheter. This pacing system is not suitable for use in these patients; however, the catheter can still be used for pressure monitoring, blood sampling, fluid infusion and cardiac output determinations. Do not attempt to insert the probe if the RV port is in the RA. Damage to the tricuspid valve may result. Always make certain that the RV port is inside the ventricle before inserting the probe.

Swan-Ganz catheters can be inserted at the patient's bedside, without the aid of fluoroscopy, guided by continuous pressure monitoring. For proper RV port positioning for placement of the pacing probes, simultaneous pressure monitoring from the distal and RV lumens is recommended.

Simultaneous pressure monitoring from the distal lumen is recommended. Fluoroscopy is recommended for femoral vein insertion.

Note: The V-pacing and A-pacing probes (in Model 991F8) and the model D98100 Chandler transluminal V-pacing probe should be placed prophylactically into their respective lumens immediately after placement of the catheter. Severe difficulty in probe passage may be encountered if probe insertion is delayed.

Note: Should the catheter require stiffening during insertion, slowly perfuse the catheter with 5 ml to 10 ml of cold sterile saline or 5% dextrose as the catheter is advanced through a peripheral vessel.

Note: The catheter should pass easily through the right ventricle and pulmonary artery and into a wedge position in less than a minute.

Although a variety of techniques can be used for insertion, the following guidelines are provided as an aid to the physician:

It is recommended that a contamination shield be used on the catheter because of the potential need to manipulate the catheter after initial insertion.

Step	Procedure
1	Introduce the catheter into the vein through a sheath introducer using percutaneous insertion using modified Seldinger technique. Note: The catheter can be inserted through a jugular, subclavian vein, or antecubital fossa. Femoral insertion is not recommended. Note: Femoral insertion should be used only when short-term pacing is recommended (e.g., catheterization lab procedures), because of the possible placement of the Chandler probe in the RV outflow tract.
2	Under continuous pressure monitoring, with or without fluoroscopy, gently advance the catheter into the right atrium. Entry of the catheter tip into the thorax is signaled by an increased respiratory fluctuation in pressure. Figure 4 on page 151 shows the characteristic intracardiac and pulmonary pressure waveforms. Note: When the catheter is near the junction of the right atrium and the superior or inferior vena cava of the typical adult patient, the tip has been advanced approximately 40 cm from the right or 50 cm from the left antecubital fossa, 15 to 20 cm from the jugular vein, or 10 to 15 cm from the subclavian vein, or 30 cm from the femoral vein.
3	Using the syringe provided, inflate the balloon with CO₂ or air to the recommended volume maximum. Do not use liquid. Note that an offset arrow on the gate valve indicates the "closed" position. Note: Inflation is usually associated with a feeling of resistance. On release, the syringe plunger should usually spring back. If no resistance to inflation is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may continue to be used for hemodynamic monitoring. However, be sure to take precautions to prevent infusion of air or liquid into the balloon lumen. WARNING: Pulmonary complications may result from improper inflation technique. To avoid damage to the pulmonary artery and possible balloon rupture, do not inflate above the recommended volume.
4	Advance the catheter until pulmonary artery occlusion pressure (PAOP) is obtained, then passively deflate the balloon by removing the syringe from the gate valve. Do not forcefully aspirate as this may damage the balloon. After deflation, re-attach the syringe. Note: Avoid prolonged maneuvers to obtain wedge pressure. If difficulties are encountered, give up the "wedge". Note: Before reinflation with CO₂ or air, completely deflate the balloon by removing the syringe and opening the gate valve. Precaution: It is recommended that the provided syringe be re-attached to the gate valve after balloon deflation to prevent inadvertent injection of liquids into the balloon lumen. Precaution: If a right ventricular pressure tracing is still observed after advancing the catheter several cm beyond the point where the initial right ventricular pressure was observed, the catheter may be looping in the ventricle, which can result in kinking or knotting of the catheter (see Complications). Deflate the balloon and withdraw the catheter into the right atrium. Reinflate the balloon and readvance the catheter to a pulmonary artery wedge position, then deflate the balloon. Precaution: Catheter looping may occur when excessive length has been inserted, which could result in kinking or knotting (see Complications). If the right ventricle is not entered after advancing the catheter 15 cm beyond entry into the right atrium, the catheter may be looping, or the tip may be engaged in a neck vein with only the proximal shaft advancing into the heart. Deflate the balloon and withdraw the catheter until the 20 cm mark is visible. Reinflate the balloon and advance the catheter.

Step	Procedure
5	Reduce or remove any excessive length or loop in the right atrium or ventricle by slowly pulling the catheter back approximately 2 to 3 cm. Precaution: Do not pull the catheter across the pulmonic valve while the balloon is inflated to avoid damage to the valve.
6	Reinflate the balloon to determine the minimum inflation volume necessary to obtain a wedge tracing. If a wedge is obtained with less than the maximum recommended volume (see specifications table for balloon inflation capacity), the catheter must be withdrawn to a position where full inflation volume produces a wedge tracing. WARNING: Do not advance the catheter beyond the wedge position to avoid pulmonary artery rupture. Precaution: Overtightening the proximal Tuohy-Borst adapter of the contamination shield may impair catheter function by potentially compressing and occluding the lumens.
7	Confirm final catheter tip position with chest X-ray.

Note: If using a contamination shield, extend the distal end towards the introducer valve. Extend the proximal end of the catheter contamination shield to desired length, and secure.

Note: After deflation, the catheter tip may tend to recoil towards the pulmonic valve and slip back into the right ventricle, requiring that the catheter be repositioned.

10.0 Positioning the Paceport and A-V Paceport Catheters when Used for Pacing

For precautions and a detailed description of the insertion procedure for the pacing probe, consult the package insert provided with each probe.

Note: Ideal placement of the RV port of the Paceport or A-V Paceport catheter when the Chandler probe is to be inserted is 1 to 2 cm distal to the tricuspid valve.

Step	Procedure
1	During catheter insertion, it is recommended that the distal and RV lumen pressures be simultaneously monitored.
2	Advance the catheter into the pulmonary artery wedge position. Deflate the balloon.
3	Refer to the pacing probe (Models D98100 and D98500) package insert for detailed instructions on insertion.

WARNING: If the RV port is too distal, then the probe may exit the RV port pointed toward the RV outflow tract. This may result in poor thresholds, unstable pacing, and potential damage to the outflow tract and pulmonic valve.

Note: A radiopaque marker is provided at the RV port to aid in port placement and identification by chest X-ray film or fluoroscopy.

Note: Diaphragmatic pacing may occasionally occur; it can usually be alleviated by advancing the catheter 0.5 to 1 cm.

11.0 Guidelines for Femoral Insertion

Precaution: Femoral insertion should be used only when short-term pacing is recommended (e.g., catheterization lab procedures), because of the possible placement of the Chandler probe (Model D98100) in the RV outflow tract.

Precaution: Femoral insertion may lead to excessive catheter length in the right atrium and difficulties in obtaining a pulmonary artery wedge (occlusion) position.

Precaution: With femoral insertion, it is possible to transfix the femoral artery in some situations during percutaneous entry into the vein. Proper femoral vein puncture technique should be followed, including removal of the innermost occluding stylet when the insertion set needle is advanced toward the vein.

Insertion under fluoroscopy is recommended when using the femoral approach.

- When advancing the catheter into the inferior vena cava, the catheter may slip into the opposite iliac vein. Pull the catheter back into the ipsilateral iliac vein, inflate the balloon, and let the bloodstream carry the balloon into the inferior vena cava.
- If the catheter does not pass from the right atrium into the right ventricle, it may be necessary to change the orientation of the tip. Gently rotate the catheter and simultaneously withdraw it several centimeters. Care must be exercised so that the catheter is not kinked as it is rotated.
- If difficulty is encountered when positioning the catheter, a suitably sized guidewire may be inserted to stiffen the catheter.

Precaution: To avoid damage to intracardiac structures, do not advance the guidewire beyond the catheter tip. The tendency for thrombus formation will increase with prolonged guidewire use. Keep the period of time that the guidewire is used to a minimum. Aspirate 2 to 3 ml from the catheter lumen and flush twice after guidewire removal.

- The Paceport catheter, (Model 931F75) and A-V Paceport catheter (Model 991F8), Chandler probes (Model D98100), and A-pacing probes (Model D98500) can be inserted successfully under fluoroscopy using the right femoral vein. However, because of the characteristic shorter catheter loop in the right ventricle, the catheter RV port and Chandler probe (Model D98100) may become oriented towards the RV outflow tract (pulmonary artery) rather than the apex. This orientation may adversely affect stability of long term pacing. In addition, the shorter catheter loop may require advancing the catheter tip into the peripheral pulmonary artery in order to position the catheter RV port distal to the tricuspid valve, possibly resulting in permanent wedge or difficulties in measuring wedge pressure.
- When advancing the catheter into the inferior vena cava, the catheter may slip into the opposite iliac vein. Pull the catheter back into the ipsilateral iliac vein, inflate the balloon, and let the bloodstream carry the balloon into the inferior vena cava.
- If the catheter does not pass from the right atrium into the right ventricle, it may be necessary to change the orientation of the tip. Gently rotate the catheter and simultaneously withdraw it several centimeters. Care must be exercised so that the catheter is not kinked as it is rotated.
- If difficulty is encountered in positioning the catheter, a suitable sized guidewire may be inserted to stiffen the catheter.

Precaution: To avoid damage to intracardiac structures, do not advance the guidewire beyond the catheter tip. The tendency for thrombus formation will increase with the duration of guidewire use. Keep the period of time that the guidewire is used to a minimum; aspirate 2 to 3 ml from the catheter lumen and flush twice after guidewire removal.

Note: Femoral insertion should be used only when short-term pacing is required (e.g., cath lab procedures) because of the possible placement of the Chandler probe in the RV outflow tract.

12.0 Maintenance and Use *in situ*

The catheter should remain indwelling only as long as is required by the patient's condition.

Precaution: The incidence of complications increases significantly with indwelling periods longer than 72 hours.

12.1 Catheter Tip Position

Keep catheter tip centrally located in a main branch of the pulmonary artery near the hilum of the lungs. Do not advance tip too far peripherally. Tip should be kept where full or near full inflation volume is required to produce a wedge tracing. The tip migrates toward periphery during balloon inflation.

12.2 Catheter Tip Migration

Anticipate spontaneous catheter tip migration towards periphery of pulmonary bed. Continuously monitor distal lumen pressure to verify tip position. If wedge tracing is observed when balloon is deflated, pull catheter back. Damage may be caused by prolonged occlusion or over-distension of vessel upon re-inflation of the balloon.

Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal (3 to 5 cm) just before bypass should be considered, as it may help reduce distal migration and prevent permanent catheter wedging postbypass. After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing before inflating the balloon.

Precaution: Over a period of time, the catheter tip may migrate towards the periphery of the pulmonary bed and lodge in a small vessel. Damage may occur either by prolonged occlusion or by over-distention of the vessel upon reinflation of the balloon (see Complications).

PA pressures should be continuously monitored with the alarm parameter set to detect physiologic changes as well as spontaneous wedge.

12.3 Balloon Inflation and Wedge Pressure Measurement

Re-inflation of the balloon should be performed gradually while monitoring pressures. Inflation is usually associated with a feeling of resistance. If no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation at once. The catheter may still be used for hemodynamic monitoring, however, take precautions against infusion of air or liquids into the balloon lumen. During normal catheter use, keep inflation syringe attached to gate valve to prevent inadvertent injection of liquid into the balloon inflation lumen.

Measure wedge pressure only when necessary and only when tip is properly positioned (see above). Avoid prolonged maneuvers to obtain wedge pressure and keep wedge time to a minimum (two respiratory cycles or 10 - 15 seconds), especially in patients with pulmonary hypertension. If

difficulties are encountered, discontinue wedge measurements. In some patients, pulmonary arterial end-diastolic pressure can often be substituted for pulmonary artery wedge pressure if the pressures are nearly identical, obviating the need for repeated balloon inflation.

12.4 Spontaneous Tip Wedging

The catheter may migrate into the distal pulmonary artery and spontaneous tip wedging may occur. To avoid this complication, pulmonary artery pressure should be monitored continuously with a pressure transducer and display monitor.

Forward advancement should never be forced if resistance is encountered.

12.5 Patency

All pressure monitoring lumens should be filled with a sterile, heparinized saline solution (e.g. 500 I.U. heparin in 500 ml saline) and flushed at least once each half hour or by continuous slow infusion. If loss of patency occurs and cannot be corrected by flushing, the catheter should be removed.

12.6 General

Keep pressure monitoring lumens patent by intermittent flush or continuous slow infusion with heparinized saline solution or use of a heparin lock using the provided injection caps with heparinized saline solution. Infusion of viscous solutions (e.g., whole blood or albumin) is not recommended, as they flow too slowly and may occlude the catheter lumen.

To use injection caps:

Step	Procedure
1	Disinfect injection caps before entry with syringe needle (see Complications).
2	Use a small bore needle (22 gauge (0.7 mm) or smaller) to puncture and inject through the injection caps.

WARNING: To avoid pulmonary artery rupture, never flush the catheter when the balloon is wedged in the pulmonary artery.

Periodically check IV lines, pressure lines, and transducers to keep them free of air. Also ensure that connecting lines and stopcocks remain tightly fitted.

13.0 Cardiac Output Determination

To determine cardiac output by thermodilution, a known amount of sterile solution of known temperature is injected into the right atrium or vena cava, and the resultant change in blood temperature is measured in the pulmonary artery by the catheter thermistor. Cardiac output is inversely proportional to the integrated area under the resulting curve. This method has been shown to provide good correlation with the direct Fick method and dye dilution technique for cardiac output determination.

Refer to the appropriate cardiac output computer manual for specific instructions in the use of thermodilution catheters for cardiac output determination. Correction factors or computation constants needed to correct for indicator heat transfer are given in the specifications.

Edwards cardiac output computers require that a computation constant be used to correct for injectate temperature rise as it passes through the catheter. The computation constant is a function of injectate volume, temperature, and catheter dimensions. The computation constants listed in the specifications have been determined *in vitro*.

14.0 MRI Information

The Swan-Ganz device is MR unsafe as the result of the device containing metallic components, which experience RF-induced heating in the MRI environment; therefore the device poses hazards in all MRI environments.

15.0 Complications

Invasive procedures involve some patient risks. Although serious complications are relatively uncommon, the physician is advised before deciding to use the catheter, to consider the potential benefits in relation to the possible complications.

The techniques for insertion, methods of using the catheter to obtain patient data information, and the occurrence of complications is well described in the literature.

Strict adherence to these instructions and awareness of risks reduces the incidence of complications.

Several known complications include:

15.1 Perforation of the Pulmonary Artery

Factors which are associated with the development of fatal pulmonary artery rupture are pulmonary hypertension, advanced age, cardiac surgery with hypothermia and anticoagulation, distal catheter tip migration, and arteriovenous fistula formation and other vascular traumas.

Extreme care should therefore be exercised during the measurement of pulmonary artery wedge pressure in patients with pulmonary hypertension.

In all patients, balloon inflation should be limited to two respiratory cycles, or 10 to 15 seconds.

A central location of the catheter tip near the hilum of the lung may prevent pulmonary artery perforation.

15.2 Pulmonary Infarction

Tip migration with spontaneous wedging, air embolism, and thromboembolism can lead to pulmonary artery infarction.

15.3 Cardiac Arrhythmias

Although usually transient and self-limiting, arrhythmias may occur during insertion, withdrawal, and repositioning of the tip from the pulmonary artery into the right ventricle. Whereas premature ventricular contractions are the most commonly observed arrhythmias, ventricular tachycardia and atrial and ventricular fibrillation have also been reported. ECG monitoring and the immediate availability of antiarrhythmic drugs and defibrillating equipment are recommended. Use of prophylactic lidocaine should be considered to decrease the incidence of ventricular arrhythmias during catheterization.

15.4 Knotting

Flexible catheters have been reported to form knots, most often as a result of looping in the right ventricle. Sometimes the knot can be resolved by insertion of a suitable guidewire and manipulation of the catheter under fluoroscopy. If the knot does not include any intracardiac structures, the knot may be gently tightened and the catheter withdrawn through the site of entry.

15.5 Sepsis/Infection

Positive catheter-tip cultures resulting from contamination and colonization have been reported, as well as incidences of septic and aseptic vegetation in the right heart. Increased risks of septicemia and bacteremia have been associated with blood sampling, the infusing of fluids, and catheter-related thrombosis. Preventive measures should be taken to guard against infection.

15.6 Other Complications

Other complications include right bundle branch block and complete heart block, tricuspid and pulmonic valve damage, pneumothorax, thrombosis, blood loss, cardiac structure/wall injury or damage, hematoma, embolism, anaphylaxis, cardiac tissue/artery burn.

In addition, allergic reactions to latex have been reported. Physicians should identify latex sensitive patients and be prepared to treat allergic reactions promptly.

Refer to <https://meddeviceinfo.edwards.com/> for a Summary of Safety and Clinical Performance for this medical device. After the launch of the European Database on Medical Devices/Eudamed, refer to <https://ec.europa.eu/tools/eudamed> for a SSCP for this medical device.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

15.7 Long-Term Monitoring

The duration of catheterization should be the minimum required by the patient's clinical state since the risk of thromboembolic and infection complications increases with time. The incidence of complications increases significantly with indwelling periods longer than 72 hours. Prophylactic systemic anticoagulation and antibiotic protection should be considered when long-term catheterization (i.e. over 48 hours) is required, as well as in cases involving increased risk of clotting or infection.

16.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged. Do not resterilize.

The packaging is designed to avoid crushing of the catheter and to protect the balloon from exposure to the atmosphere. It is therefore recommended that the catheter remain inside the package until use.

17.0 Storage Conditions

Store in a cool, dry place.

18.0 Operating Conditions/Use Environment

Intended to operate under physiological conditions of the human body in a controlled clinical environment.

19.0 Shelf Life

The shelf life is marked on each package. Storage or usage beyond the expiration date may result in product deterioration and could lead to illness or an adverse event as the device may not function as originally intended.

20.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

21.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Because of continuing product improvements, prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

22.0 Summary Guidelines for Safe Use of Balloon-tipped Pulmonary Artery Catheters

1. Keep catheter tip centrally located in a main branch of the pulmonary artery:

- During insertion, inflate the balloon to the full recommended volume (1.5 ml) and advance the catheter to a pulmonary artery wedge position. Deflate the balloon.
- To reduce or remove any redundant length or loop in the right atrium or ventricle, slowly pull the catheter back 2 to 3 cm.
- Do not advance the catheter tip too far peripherally. Ideally, the catheter tip should be located near the hilum of the lungs. Remember, the tip migrates towards the periphery of the lungs during balloon inflation. Therefore, a central location before inflation is important.
- Keep the tip at all times in a position where a full or near-full (1.0 to 1.5 ml) inflation volume is necessary to produce a “wedge” tracing.

2. Anticipate spontaneous catheter tip migration toward the periphery of the pulmonary bed:

- Reduce any redundant length or loop in the right atrium or ventricle **at the time of insertion** to prevent subsequent peripheral migration (see No. 1).
- Monitor the distal tip pressure continuously to ensure that the catheter is not inadvertently wedged with the balloon deflated (this may induce pulmonary infarction).
- Check catheter position daily by chest X-ray film to detect peripheral placement. If migration has occurred, pull the catheter back to a central pulmonary artery position, carefully avoiding contamination of the insertion site.
- Spontaneous catheter tip migration towards the periphery of the lung occurs during cardiopulmonary bypass. Partial catheter withdrawal (3 to 5 cm) just before bypass should be considered, as withdrawal may help reduce the amount of distal migration and may prevent permanent catheter wedging in the postbypass period. After termination of bypass, the catheter may require repositioning. Check the distal pulmonary artery tracing **before** inflating the balloon.

3. Exercise caution when inflating the balloon:

- If “wedge” is obtained at volumes less than 1.0 ml, pull the catheter back to a position where the full or near-full inflation volume (1.0 to 1.5 ml) produces a wedge pressure tracing.
- Check the distal pressure waveform before inflating the balloon. If the waveform appears dampened or distorted, do not inflate the balloon. The catheter may be wedged with the balloon deflated. Check catheter position.
- When the balloon is reinflated to record wedge pressure, add the inflation medium (CO₂ or air) **slowly** under continuous monitoring of the pulmonary artery pressure waveform. Stop inflating **immediately** when the pulmonary artery tracing is seen to change to pulmonary artery wedge pressure. Remove the syringe to allow rapid balloon deflation, then reattach the syringe to the balloon lumen. Air should never be used for balloon inflation in any situation where air may enter the arterial circulation (see **Insertion Procedure**).
- Never over-inflate the balloon beyond the maximum volume printed on the catheter shaft (1.5 ml). Use the volume-limited syringe provided with the catheter.
- Do not use liquids for balloon inflation; they may be irretrievable and may prevent balloon deflation.
- Keep the syringe attached to the balloon lumen of the catheter to prevent accidental injection of liquids into the balloon.

4. Obtain a pulmonary artery occlusion “wedge” pressure only when necessary:

- If the pulmonary artery diastolic (PAD) and the wedge (PAW) pressures are nearly identical, wedging the balloon may not be necessary: measure PAD pressure instead of PAW as long as the patient's heart rate, blood pressure, cardiac output and clinical state remain stable. However, in states of changing pulmonary arterial and pulmonary venous tone (i.e., sepsis, acute respiratory failure, shock), the relationship between PAD and “wedge” may change with the patient's clinical condition. PAW measurement may be necessary.
- Keep “wedge” time to a minimum (two respiratory cycles or 10 to 15 seconds), especially in patients with pulmonary hypertension.
- Avoid prolonged maneuvers to obtain wedge pressure. If difficulties are encountered, give up the “wedge”.
- Never flush the catheter when the balloon is wedged in the pulmonary artery.

5. Patients at highest risk of pulmonary artery rupture or perforation are elderly patients with pulmonary hypertension. These are usually elderly patients who are undergoing cardiac surgery with anticoagulation and hypothermia. Proximal catheter tip location near the hilum of the lungs may reduce the incidence of pulmonary artery perforation.

Model	Indications	Probe Placement	Pacing Probe
931F75	Temporary Ventricular Pacing Fluid Administration	Ventricle	D98100
991F8	Temporary	Atrium and Ventricle	D98500
	Atrial - Ventricular Pacing		D98100
	Temporary Atrial Pacing only	Atrium	D98500
	Temporary Ventricular Pacing only	Ventricle	D98100
	Fluid Administration		

Computation Constants

Model		931F75	991F8	780F75M 782F75M	631F55N
Injectate Temp. (°C)	Volume (ml)	Computation Constants (CC)***			
0 - 5	10	0.564	0.568	0.564	-
	5	0.262	0.268	0.257	0.265
	3	0.139	0.147	0.143	0.148
	1	-	-	-	0.032
19 - 22	10	0.588	-	0.582	-
	5	0.283	-	0.277	0.294
	3	0.158	-	0.156	0.172
	1	-	-	-	0.049
23 - 25	10	0.612	0.616	0.607	-
	5	0.301	0.302	0.294	0.308
	3	0.177	0.176	0.170	0.183
	1	-	-	-	0.057
Computation Constants* for C0-Set+ Injectate Delivery System					
Cold Injectate					
6 °C - 12 °C	10	0.578	0.553	0.574	-
8 °C - 12 °C	5	-	0.277	-	-
8 °C - 16 °C	5	0.272	-	0.287	0.284
8 °C - 16 °C	3	-	-	-	0.169
Room Temperature Injectate					
18 °C - 25 °C	10	0.592	0.607	0.595	-
18 °C - 25 °C	5	0.290	0.295	0.298	0.306
	3	-	-	-	0.182

*3 ml injectate is not recommended.

***CC = (1.08)C_T(60)/(V_I)

Note: Computation constants for oximetry catheters are the same as for the Swan-Ganz VIP Catheter.

Specifications

Swan-Ganz Thermodilution catheter	Paceport catheter (Model 931F75)	A-V Paceport catheter (Model 991F8)	Oximetry catheter 631F55N**	Oximetry Paceport catheter 780F75M	VIP Oximetry catheter 782F75M
Usable Length (cm)	110	110	75	110	110
Catheter Body French Size	7.5 F (2.5 mm)	8 F (2.7 mm)	5.5 F (1.8 mm)	7.5 F (2.5 mm)	7.5 F (2.5 mm)
Body Color			White	Yellow	Yellow
Balloon Inflation Capacity (ml- CO_2)	1.5	1.5	0.7	1.5	1.5
Diameter of Inflated Balloon (mm)	13	13	8	13	13
Minimum Recommended Introducer Size	8.5 F (2.8 mm)	9 F (3.0 mm)	6.5 F (2.2 mm)	8.5 F (2.8 mm)	8.5 F (2.8 mm)
Distance from Distal Tip (cm)					
Thermoresistor	-	-	1.5	4	4
RA Port	-	27	-	-	-
RV Port	19	19	-	19	-
Proximal Injectate Port	30	30	15	30	30
VIP Infusion Port	-	-	-	-	31
Lumen Volume (ml)					
Proximal Injectate	0.89	0.70	0.51	-	-
Distal	0.88	0.93	0.49	-	-
RA Lumen	-	1.07	-	-	-
RV Lumen	-	1.13	-	-	-
without Probe	1.10	-	-	-	-
with T-B adapter					
attached and probe in lumen	0.88	-	-	-	-
Infusion Rates* (ml/min)					
Distal	5	6	3	6	6
Proximal Injectate	12	8	4	9	8
RA Pacing/Infusion					
with Probe	-	1	-	-	-
without Probe	-	15	-	-	-
RV Pacing/Infusion					
with Probe	0.5	1	-	-	-
without Probe	11	12	-	13	-
VIP Lumen	-	-	-	-	14
Infusion Rate* (ml/min) using D_{50}W					
RV Lumen					
without Probe	3	-	-	-	-
with T-B adapter					
attached and probe in lumen	0.2	-	-	-	-
Radiopaque Marker	distal edge of RV Port	at RA and RV Port	-	-	-
Compatible Guidewire Diameter					
Distal Lumen	0.025 in (0.64 mm)	0.018 in (0.46 mm)	0.012 in (0.30 mm)	0.025 in (0.64 mm)	0.025 in (0.64 mm)
Proximal Injectate Lumen	0.035 in (0.89 mm)	-	-	-	-
RV Lumen	0.035 in (0.89 mm)	0.028 in (0.71 mm)	-	-	-
RA Lumen	-	0.028 in (0.71 mm)	-	-	-
Frequency Response					
Distortion at 10Hz	<3 dB	<3 dB	<3 dB	<3 dB	<3 dB
Distal Lumen					
Total Catheter Sensing Function	-	-	-	0.5%	0.5%

All specifications given are nominal values.

* Using normal saline at room temperature, 1 m above insertion site, gravity drip. Rates represent average values. Infusing blood products or hyperalimentation solutions through the RV lumen is not recommended if the RV lumen will be used for placement of the Chandler probe. All Edwards oximetry catheters, except those with black oximetry connectors, are compatible with Edwards optical modules. Models ending with the letter "M" are also compatible with Philips optical modules. The Catheter Factor required for *in vivo* calibration with Philips monitors is located on the top of the optical connector.

** In situations where the insertion site or patient physiology requires a larger insertion distance, a longer catheter model or larger introducer should be selected.

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	English	Français	Deutsch	Español	Italiano	Nederlands
	Number of lumens	Nombre de lumières	Anzahl der Lumen	Número de luces	Numero di lumi	Aantal lumina
	Exterior diameter	Diamètre externe	Außendurchmesser	Diámetro exterior	Diametro esterno	Buitendiameter
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad
	Lumen size	Taille de la lumière	Lumendurchmesser	Tamaño de la luz	Dimensioni del lume	Lumengrootte
	Model Number	Référence	Modellnummer	Número de modelo	Numero di modello	Modelnummer
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße der Einführhilfe	Tamaño de introductor mínimo	Misura minima dell'introdotore	Minimum introduceergrootte
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum
	Inner diameter	Diamètre interne	Innendurchmesser	Diámetro interior	Diametro interno	Binnendiameter
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto	Droog houden
	Do not use if package is damaged and consult instructions for use	Ne pas utiliser si le conditionnement est endommagé et consulter le mode d'emploi	Bei beschädigter Verpackung nicht verwenden und Gebrauchsanweisung beachten	No lo utilice si el envase está dañado y consulte las instrucciones de uso	Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso	Niet gebruiken als de verpakking is beschadigd en raadpleeg de gebruiksaanwijzing
	Sterile	Stérile	Steril	Estéril	Sterile	Steriel
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethyleenoxide
	Sterilized using irradiation	Stérilisé par irradiation	Durch Bestrahlung sterilisiert	Esterilizado con radiación	Sterilizzato mediante radiazioni	Gesteriliseerd met behulp van straling
	Sterilized using steam or dry heat	Stérilisé à la vapeur d'eau ou à la chaleur sèche	Mit Wasserdampf oder trockener Hitze sterilisiert	Esterilizado con vapor o calor seco	Sterilizzato con vapore o calore secco	Gesteriliseerd met behulp van stoom of droge warmte
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig
	Size	Taille	Größe	Tamaño	Misura	Afmetingen
	Consult instructions for use on the website	Consulter le mode d'emploi sur le site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Consultare le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website

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	Polski	Slovensky	Norsk	Suomi	Български	Română
	Liczba kanałów	Počet lúmenov	Antall lumen	Luumenien määrä	Брой лумени	Număr de Lumenuri
	Średnica zewnętrzna	Vonkajší priemer	Utvendig diameter	Ulkoläpimitta	Външен диаметър	Diametru exterior
	Długość użytkowa	Použitelná dĺžka	Anvendelig lengde	Käyttöpituus	Използваема дължина	Lungimea utilă
	Zalecany rozmiar przewadnika	Odporúčaná veľkosť vodiaceho drôtu	Anbefalt ledveaierstørrelse	Ohjainlangan suositeltu koko	Препоръчителен размер телен водач	Mărimea recomandată a firului de ghidaj
	Rozmiar kanału	Veľkosť lúmenu	Lumenstørrelse	Luumenin koko	Размер на лумена	Mărimea lumenului
	Numer modelu	Číslo modelu	Modellnummer	Mallinnumero	Номер на модела	Număr model
	Minimalny rozmiar introduktora	Minimálna veľkosť zavádzača	Minste tillatte innførerstørrelse	Sisäänviejän vähimmäiskoko	Минимален размер на интродюсера	Mărimea minimă a dispozitivului de introducere
	Przeostroga	Upozornenie	Forsiktig	Tärkeä huomautus	Внимание	Atenție
	Nie używać ponownie	Nepoužívajte opakovane	Må ikke gjenbrukes	Ei uudelleenkäytettävä	Да не се използва повторно	A nu se reutiliza
	Ilość	Množstvo	Antall	Määrä	Количество	Cantitate
	Numer partii	Číslo šarže	Lotnummer	Eränumero	Партиден номер	Număr de lot
	Data przydatności do użycia („Zużyć do”)	Dátum spotreby	Utløpsdato	Viimeinen käyttöpäivä	Срок на годност	A se utiliza până la data de
	Średnica wewnętrzna	Vnútný priemer	Innvendig diameter	Sisäläpimitta	Вътрешен диаметър	Diametru interior
	Pojemność balonu	Objem balónika	Ballongkapasitet	Pallon tilavuus	Вместимост на балона	Capacitatea balonului
	Chronić przed wilgocią	Uchovávať v suchu	Oppbevares tørt	Pidä kuivana	Да се пази сухо	A se păstra uscat
	Nie stosować, jeśli opakowanie jest uszkodzone oraz zapoznać się z instrukcją użycia	Nepoužívajte, ak je obal poškodený, a prečítajte si návod na použitie	Skal ikke brukes hvis pakningen er skadet. Se bruksanvisningen	Älä käyttää, jos pakkaus on vahingoittunut, ja katso käyttöohjeet	Да не се използва, ако опаковката е повредена, и направете справка с инструкциите за употреба	A nu se utiliza dacă ambalajul este deteriorat. Consultați instrucțiunile de utilizare
	Jalowy	Sterilný	Steril	Sterili	Стерилно	Steril
	Wysterylizowano przy użyciu tlenku etylenu	Sterilizované pomocou etylénoxidu	Sterilisert med etylenoksid	Steriloitu etyleenioksidilla	Стерилизирано с етиленов оксид	Sterilizat cu oxid de etilenă
	Wysterylizowano przez napromienienie	Sterilizované pomocou ožarovania	Sterilisert med stråling	Steriloitu säteilyttämällä	Стерилизирано с облъчване	Sterilizat prin iradiere
	Wysterylizowano parą lub suchym, gorącym powietrzem	Sterilizované parou alebo suchým teplom	Sterilisert med stråling eller tørr varme	Lämpö- tai höyrysteriloitu	Стерилизирано чрез автоклавиране или сух стерилизатор	Sterilizat cu abur sau prin căldură uscată
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Producent	Výrobca	Produsent	Valmistaja	Производител	Producător
	Data produkcji	Dátum výroby	Produksjonsdato	Valmistuspäivämäärä	Дата на производство	Data fabricației
	Zawiera lateks pochodzący z naturalnego kauczuku lub wykazuje jego obecność	Obsahuje alebo je prítomný prírodný kaučukový latex	Inneholder eller tilstedeværelse av naturgummilateks	Sisältää luonnokumilateksia tai sen jäämiä	Съдържание или наличие на естествен каучуков латекс	Conține sau prezintă latex din cauciuc natural
	Rozmiar	Veľkosť	Størrelse	Koko	Размер	Mărime
	Zapoznaj się z instrukcją użycia na stronie internetowej	Prečítajte si návod na použitie na webovej stránke	Se bruksanvisningen på nettsiden	Katso käyttöohjeet verkkosivustolta	Направете справка с инструкциите за употреба, намиращи се на уеб сайта	Consultați instrucțiunile de utilizare de pe site-ul web
	Postępować zgodnie z instrukcją stosowania na stronie internetowej	Postupujte podľa návodu na použitie na webovej stránke	Følg bruksanvisningen på nettsiden	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	Hrvatski
	Valendike arv	Spindžių skaičius	Lūmenu skaits	Lümen sayısı	Количество просветов	Broj lumena	Broj lumena
	Vāline diameeter	Išorinis skersmuo	Ārējais diametrs	Dış çap	Внешний диаметр	Vanjski prečnik	Vanjski promjer
	Kasutata pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjliva dužina	Upotrebjliva dužina
	Juhtetraadi soovitusslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Kılavuz Tel İçin Önerilen Boyut	Рекомендованный размер проводника	Preporučena veličina žice vodice	Preporučena veličina žice vodilice
	Valendiku suurus	Spindžio dydis	Lūmena izmērs	Lümen boyutu	Размер просвета	Veličina lumena	Veličina lumena
	Mudeli number	Modelio numeris	Modeļa numurs	Model Numarası	Номер модели	Broj modela	Broj modela
	Sisesti väikeim suurus	Minimalus įvediklio dydis	Minimālais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	Minimalna veličina uvodnog instrumenta
	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	Oprez
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	Nemojte ponovno upotrebljavati
	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	Količina
	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	Broj serije
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	Rok upotrebe
	Sisemine läbimõõt	Vidinis skersmuo	Iekšējais diametrs	İç çap	Внутренний диаметр	Unutrašnji prečnik	Unutrašnji promjer
	Ballooni mahutavus	Balionių talpa	Balona ietilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	Kapacitet balona
	Hoida kuivalt	Laikykite sausoje vietoje	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	Držati suhim
	Ärge kasutage, kui pakend on kahjustatud, ja lugege kasutusjuhiseid	Nenaudoti, jei pakuotė pažeista, ir žr. naudojimo instrukcijas	Nelietot, ja iepakojums ir bojāts, un skatīt lietošanas instrukciju	Ambalaj hasarlıysa kullanmayın ve kullanım talimatlarına başvurun	Не использовать, если упаковка повреждена; использовать в соответствии с инструкцией	Ne koristiti ako je pakovanje oštećeno i pogledajte uputstvo za upotrebu	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Steriilne	Sterilus	Sterils	Steril	Стерильно	Sterilno	Sterilno
	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksīdu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	Sterilizirano etilen-oksidom
	Steriliseeritud kiirgamist kasutades	Sterilizuota švitinant	Sterilizēts apstarojot	Işınlama yoluyla sterilize edilmiştir	Стерилизовано излучением	Sterilizovano zračenjem	Sterilizirano zračenjem
	Steriliseeritud auru või kuiva kuumusega	Sterilizuota naudojant garus arba sausą karštį	Sterilizēts, izmantojot tvaiku vai sausu karstumu	Buhar veya kuru ısı kullanılarak sterilize edilmiştir	Стерилизовано паром или сухим жаром	Sterilizovano korišćenjem pare ili suvom toplotom	Sterilizirano parom ili vrućim, suhim zrakom
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	Proizvođač
	Tootmiskuupeäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	Datum proizvodnje
	Sisaldab naturaalselt kummilateksit	Turi savo sudėtyje ar yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	Sadrži prirodni gumeni lateks ili je prisutan
	Suurus	Dydis	Izmērs	Boyut	Размер	Veličina	Veličina
	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje	Skatīt lietošanas instrukciju tieklā vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstvo za upotrebu na veb-sajtu	Pogledajte upute za upotrebu na internetskoj stranici
	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovaukitės naudojimo instrukcijomis, pateiktomis svetainėje	Ievērojiet lietošanas instrukciju, kas pieejama tieklā vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa veb-sajta	Pridržavajte se uputa za upotrebu na internetskoj stranici

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Edwards

Edwards Lifesciences GmbH
Parking 30
85748 Garching bei München
Germany

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences B.V.
Verlengde Poolseweg 16
4818 CL Breda, Netherlands

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

12/23
10007452006 A / DOC-0138850 A
© Copyright 2023, Edwards Lifesciences LLC
All rights reserved.