## DIRECTORY

|                  |    |                 |    |                |    |
|------------------|----|-----------------|----|----------------|----|
| English .....    | 1  | Português ..... | 28 | Eesti .....    | 55 |
| Français .....   | 4  | Česky .....     | 31 | Lietuvių ..... | 58 |
| Deutsch .....    | 7  | Magyar .....    | 34 | Latviešu ..... | 61 |
| Español .....    | 10 | Polski .....    | 37 | Türkçe .....   | 64 |
| Italiano .....   | 13 | Slovensky ..... | 40 | Русский .....  | 67 |
| Nederlands ..... | 16 | Norsk .....     | 43 | Srpski .....   | 70 |
| Dansk .....      | 19 | Suomi .....     | 46 | 简体中文 .....     | 73 |
| Svenska .....    | 22 | Български ..... | 49 | 한국어 .....      | 76 |
| Ελληνικά .....   | 25 | Română .....    | 52 | 日本語 .....      | 79 |

### English

## Pressure Controller Kit

**REF** PC2K and PCCVR

**Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.**

For figures, see Figure 1 on page 82 through Figure 8 on page 83.

### 1.0 Description

The pressure controller kit consists of a pressure controller (PC2) and accompanying band (PC2B). The pressure controller can be used with the EV1000 clinical platform NI and HemoSphere advanced monitoring platform to noninvasively measure blood pressure and key hemodynamic parameters using a compatible Edwards finger cuff wrapped around the finger. The pressure controller is worn on the patient's wrist and connects to the compatible Edwards noninvasive monitoring system, heart reference sensor and compatible Edwards finger cuff(s) (Figure 1 on page 82).

A pressure controller cover is available as an accessory. The pressure controller cover secures the heart reference sensor into the pressure controller.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use.

## 2.0 Intended Use

The pressure controller is intended for use with an Edwards compatible noninvasive monitoring system - composed of compatible monitor, pressure source (pump), heart reference sensor, compatible Edwards finger cuff(s) and pressure controller - for noninvasive measurement of blood pressure and associated hemodynamic parameters. Refer to the operator's manual of the compatible Edwards monitor being used for specific information on the intended use environment and patient population.

## 3.0 Instructions for Use

### 3.1 Apply the Pressure Controller

| Step | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Wrap the pressure controller band around the patient's wrist (Figure 2 on page 82, step 1).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| 2    | Snap the pressure controller into the plastic sleeve of the band, making sure that the cuff connectors are facing towards the fingers (Figure 2 on page 82, step 2).                                                                                                                                                                                                                                                                                                                                                                                                            |
| 3    | Attach the pressure controller cable to the port on compatible Edwards noninvasive monitoring system (Figure 1 on page 82, 1A).                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 4    | Remove the plastic connector caps in order to connect the finger cuff and heart reference sensor to the pressure controller according to the Instructions for Use supplied with the particular units (Figure 3 on page 82, steps 1-2), or consult the operator's manual of the EV1000 clinical platform NI and HemoSphere advanced monitoring platform.<br><b>Note: It is recommended that the cuff connector caps be kept and used to protect the pressure controller against the ingress of water and dirt when only a single cuff is used (Figure 3 on page 82, step 3).</b> |
| 5    | Once the heart reference sensor (HRS) is connected to the pressure controller, then place the pressure controller cover (PCCVR) on the pressure controller.                                                                                                                                                                                                                                                                                                                                                                                                                     |

**WARNING: Do not overtighten the pressure controller band.**

**WARNING: Do not apply pressure controller band on injured skin as this can cause further injury.**

### 3.2 Remove the Pressure Controller Band

To remove the pressure controller from the pressure controller band:

| Step | Procedure                                                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Pull the sleeve slightly outwards (Figure 4 on page 82, step 1) and tilt the pressure controller to remove it from the sleeve (Figure 4 on page 82, step 2). |

If you have a pressure controller cover (PCCVR), follow the instructions below to apply and remove the PCCVR.

### 3.3 Apply the Pressure Controller Cover

| Step | Procedure                                                                                                                                                                            |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Ensure the heart reference sensor (HRS) is attached prior to attaching the pressure controller cover to the pressure controller.                                                     |
| 2    | Place the pressure controller cover's back mounting notch around the pressure controller cable (Figure 5 on page 83).                                                                |
| 3    | Snap the pressure controller cover over the pressure controller, making sure that the pressure controller cover does not interfere with the heart reference sensor (HRS) connection. |

### 3.4 Remove the Pressure Controller Cover

| Step | Procedure                                                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | Pull the pressure controller cover upwards from the front tab (Figure 6 on page 83) and tilt the pressure controller cover to remove it from the pressure controller (Figure 7 on page 83). |

## 4.0 Cleaning Instructions

The pressure controller and pressure controller cover can be cleaned using the following disinfectants:

- 70% isopropyl alcohol solution
- 10% sodium hypochlorite water solution

| Step | Procedure                                                      |
|------|----------------------------------------------------------------|
| 1    | Moisten a clean cloth with disinfectant and wipe the surfaces. |
| 2    | Dry the surface with a clean dry cloth.                        |

**CAUTION: Do not immerse the pressure controller.**

**CAUTION: Do not disinfect the pressure controller by autoclave or gas sterilization.**

**CAUTION: Do not immerse any cable connectors in fluid.**

**CAUTION: Do not pinch any heart reference sensor tubes or wires under the pressure controller cover during application. Be careful the only wire between the back mounting notch is the pressure controller cable.**

**CAUTION: Do not lift PCCVR from any other point than the front tab (Figure 8 on page 83).**

**CAUTION: Do not use any other cleaning agents aside from those listed.**

The pressure controller band is intended for limited reuse. The operator shall assess whether reuse is appropriate. When reused, the same cleaning instructions can be followed as specified for the pressure controller.

The pressure controller cover is intended for limited reuse. The operator shall assess whether reuse is appropriate. When reused, follow the cleaning instruction above. Replace if damaged.

## 5.0 Accessories

| Model Number | Accessory Description                                      |
|--------------|------------------------------------------------------------|
| PC2K         | Pressure Controller Kit                                    |
| PC2B         | Pressure Controller Band (Multi Pack of 5)                 |
| PC2CCC       | Pressure Controller Cuff Connector Caps (Multi Pack of 25) |
| PCCVR        | Pressure Controller Cover                                  |

## 6.0 Warranty

The pressure controller is warranted for a period of one (1) year from the date of purchase when used in accordance with the directions for use. This warranty does not include the pressure controller band. For a full description of the warranty conditions refer to the operator's manual of the EV1000 clinical platform NI and HemoSphere advanced monitoring platform. Both the pressure controller and the pressure controller band are non-serviceable parts.

## 7.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): ..... 830.822.9837

Outside the U.S. and Canada

(24 hours): ..... 949.250.2222

In the UK: ..... 0870 606 2040 - Option 4

In Ireland: ..... 01 8211012 - Option 4

This document is subject to change without prior notice. Refer to the latest version of the EV1000 clinical platform NI and HemoSphere advanced monitoring platform operator's manual for additional information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

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|                                                                                     | English                                                                                                | Français                                                                                            | Deutsch                                                                                     | Español                                                                                           |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|
| <b>REF</b>                                                                          | Catalogue Number                                                                                       | Référence catalogue                                                                                 | Artikelnummer                                                                               | Número de catálogo                                                                                |
|    | Separate collection for electrical and electronic equipment in accordance with EC Directive 2012/19/EU | Tri sélectif des équipements électriques et électroniques conformément à la Directive CE 2012/19/UE | Getrennte Sammlung von Elektro- und Elektronikgeräten gemäß EU-Richtlinie 2012/19/EU        | Recogida por separado de equipos eléctricos y electrónicos según la directiva de la CE 2012/19/UE |
| <b>MD</b>                                                                           | Medical device                                                                                         | Dispositif médical                                                                                  | Medizinprodukt                                                                              | Producto sanitario                                                                                |
| <b>SN</b>                                                                           | Serial number                                                                                          | Numéro de série                                                                                     | Seriennummer                                                                                | Número de serie                                                                                   |
| <b>LOT</b>                                                                          | Lot number                                                                                             | N° de lot                                                                                           | Chargenbezeichnung                                                                          | Número de lote                                                                                    |
| <b>#</b>                                                                            | Quantity                                                                                               | Quantité                                                                                            | Menge                                                                                       | Cantidad                                                                                          |
|     | Conformité Européenne (CE Mark) of TÜV SÜD Product Service GmbH (notified body)                        | Conformité européenne (marquage CE) de TÜV SÜD Product Service GmbH (organisme notifié)             | Conformité Européenne (CE-Kennzeichnung) der TÜV SÜD Product Service GmbH (Benannte Stelle) | Conformité Européenne (Marca CE) de TÜV SÜD Product Service GmbH (organismo notificado)           |
| <b>EC REP</b>                                                                       | Authorized representative in the European Community                                                    | Représentant autorisé dans la Communauté européenne                                                 | Bevollmächtigter in der Europäischen Gemeinschaft                                           | Representante autorizado en la Comunidad Europea                                                  |
|   | Date of manufacture                                                                                    | Date de fabrication                                                                                 | Herstellungsdatum                                                                           | Fecha de fabricación                                                                              |
|  | Manufacturer                                                                                           | Fabricant                                                                                           | Hersteller                                                                                  | Fabricante                                                                                        |
|  | Follow instructions for use                                                                            | Consulter le mode d'emploi                                                                          | Gebrauchsanweisung befolgen                                                                 | Siga las instrucciones de uso                                                                     |
|   | Follow instructions for use on the website                                                             | Consulter le mode d'emploi sur le site Web                                                          | Siehe Gebrauchsanweisung auf Website                                                        | Siga las instrucciones de uso del sitio web                                                       |
|  | Caution                                                                                                | Avertissement                                                                                       | Vorsicht                                                                                    | Aviso                                                                                             |
|  | MR Unsafe                                                                                              | Risques en milieu RM                                                                                | MR-unsicher                                                                                 | No seguro para RM                                                                                 |
| <b>Rx only</b>                                                                      | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.            | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.         | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.       |
|  | No serviceable parts inside                                                                            | Pas d'accessoires réparables à l'intérieur                                                          | Das Gerät enthält keine Teile, die gewartet werden können                                   | No contiene piezas que requieran mantenimiento                                                    |
|  | Keep dry                                                                                               | Tenir au sec                                                                                        | Vor Nässe schützen                                                                          | Mantener en un lugar seco                                                                         |
|  | Fragile, handle with care                                                                              | Fragile, manipuler avec précaution                                                                  | Vorsicht! Zerbrechlich!                                                                     | Frágil; manipular con cuidado                                                                     |
|  | Keep away from sunlight                                                                                | Tenir à l'abri de la lumière du soleil                                                              | Vor Sonneneinstrahlung schützen                                                             | Mantener alejado de la luz del sol                                                                |
| <b>IP4X</b>                                                                         | Extent of protection against the ingress of objects                                                    | Extension de la protection contre l'entrée d'objets                                                 | Erweiterter Schutz gegen feste Fremdkörper                                                  | Grado de protección contra la entrada de objetos                                                  |
|  | Temperature limit                                                                                      | Limite de température                                                                               | Temperaturgrenze                                                                            | Límite de temperatura                                                                             |

## Symbol Legend • Légende des symboles • Zeichenerklärung • Significado de los símbolos

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|                                                                                     | <b>Български</b>                                                                              | <b>Română</b>                                                                                             | <b>Eesti</b>                                                                                | <b>Lietuvių</b>                                                                             |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|    | Каталожен номер                                                                               | Număr de catalog                                                                                          | Kataloogi number                                                                            | Katalogo numeris                                                                            |
|    | Разделно събиране на електрическо и електронно оборудване според Директивата на ЕС 2012/19/ЕС | Colectarea separată a echipamentelor electrice și electronice, în conformitate cu directiva CE 2012/19/UE | Elektri- ja elektroonikaseadmete eraldi kogumine kooskõlas EÜ direktiiviga 2012/19/EL       | Atskiras elektros ir elektroninės įrangos atliekų surinkimas pagal EB direktyvą 2012/19/ES  |
|    | Медицинско изделие                                                                            | Dispozitiv medical                                                                                        | Meditsiiniseade                                                                             | Medicinos prietaisas                                                                        |
|    | Сериен номер                                                                                  | Număr de serie                                                                                            | Seerianumber                                                                                | Serijos numeris                                                                             |
|    | Партиден номер                                                                                | Număr de lot                                                                                              | Partii number                                                                               | Partijos numeris                                                                            |
|    | Количество                                                                                    | Cantitate                                                                                                 | Kogus                                                                                       | Kiekis                                                                                      |
|     | Conformité Européenne (CE маркировка) на TÜV SÜD Product Service GmbH (нотифициран орган)     | Conformité Européenne (marcă CE) al TÜV SÜD Product Service GmbH (organism notificat)                     | TÜV SÜD Product Service GmbH (teavitatud asutus) Conformité Européenne (CE-märk)            | TÜV SÜD Product Service GmbH (notifikuotoji įstaiga) Conformité Européenne (CE žymė)        |
|     | Упълномощен представител в Европейската общност                                               | Reprezentant autorizat în Comunitatea Europeană                                                           | Volitatud esindaja Euroopa Ühenduses                                                        | Igaliotasis atstovas Europos Bendrijoje                                                     |
|    | Дата на производство                                                                          | Data fabricației                                                                                          | Tootmiskuupäev                                                                              | Pagaminimo data                                                                             |
|   | Производител                                                                                  | Producător                                                                                                | Tootja                                                                                      | Gamintojas                                                                                  |
|  | Следвайте инструкциите за употреба                                                            | Urmați instrucțiunile de utilizare                                                                        | Järgige kasutusjuhendit!                                                                    | Laikykitės naudojimo instrukcijų                                                            |
|   | Спазвайте инструкциите за употреба, дадени на уебсайта                                        | Respectați instrucțiunile de utilizare disponibile pe site-ul web                                         | Järgige kasutusjuhendit, mille leiata veebisaidilt                                          | Vadovaukitės naudojimo instrukcijomis, pateiktomis svetainėje                               |
|  | Внимание                                                                                      | Atenție                                                                                                   | Ettevaatust                                                                                 | Perspėjimas                                                                                 |
|  | Небезопасно при MR                                                                            | Incompatibil RM                                                                                           | Ohtlik magnetresonantstomograafias                                                          | MR nesaugus                                                                                 |
|   | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.   | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.               | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|  | Във вътрешността няма части, подлежащи на сервизиране                                         | Fără componente care pot fi reparate în interior                                                          | Ei sisalda hooldatavaid osasid                                                              | Viduje nėra keičiamų dalių                                                                  |
|  | Да се съхранява сухо                                                                          | A se păstra uscat                                                                                         | Hoidke kuivana                                                                              | Laikykite sausoje vietoje                                                                   |
|  | Чупливо, да се борава внимателно                                                              | Fragil, a se manipula cu atenție                                                                          | Õrn, käsitseda ettevaatlikult                                                               | Trapus, elgtis atsargiai                                                                    |
|  | Да се пази от слънчева светлина                                                               | A se feri de lumina solară                                                                                | Hoida eemal päikesevalgusest                                                                | Saugoti nuo saulės šviesos                                                                  |
|   | Степен на защита срещу проникване на предмети                                                 | Grad de protecție împotriva pătrunderii obiectelor                                                        | Objektide sissetungi eest kaitsmise määrad                                                  | Objektų patekimo į prietaisą apsaugos lygis                                                 |
|  | Температурно ограничение                                                                      | Limită de temperatură                                                                                     | Temperatuuriipiirang                                                                        | Temperatūros apribojimas                                                                    |

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|                                                                                     | Latviešu                                                                                           | Türkçe                                                                                         | Русский                                                                                                                                           | Srpski                                                                                      |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
| <b>REF</b>                                                                          | Kataloga numurs                                                                                    | Katalog numarası                                                                               | Номер по каталогу                                                                                                                                 | Kataloški broj                                                                              |
|    | Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES direktīvu 2012/19/ES             | 2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır | Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU                                              | Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU   |
| <b>MD</b>                                                                           | Medicīnas ierīce                                                                                   | Tıbbi cihaz                                                                                    | Медицинское устройство                                                                                                                            | Medicinsko sredstvo                                                                         |
| <b>SN</b>                                                                           | Sērijas numurs                                                                                     | Seri numarası                                                                                  | Серийный номер                                                                                                                                    | Serijski broj                                                                               |
| <b>LOT</b>                                                                          | Partijas numurs                                                                                    | Lot numarası                                                                                   | Номер партии                                                                                                                                      | Broj partije                                                                                |
| <b>#</b>                                                                            | Daudzums                                                                                           | Miktar                                                                                         | Количество                                                                                                                                        | Količina                                                                                    |
|     | Conformité Européenne (CE marķējums), ko piešķir TÜV SÜD Product Service GmbH (paziņotā struktūra) | TÜV SÜD Product Service GmbH (onaylı kuruluş) Conformité Européenne (CE İşareti)               | Соответствие требованиям директив Европейского Союза (маркировка CE) подтверждено компанией TÜV SÜD Product Service GmbH (нотифицированный орган) | Conformité Européenne (CE znak) kompanije TÜV SÜD Product Service GmbH (ovlašteno telo)     |
| <b>EC REP</b>                                                                       | Pilnvarotais pārstāvis Eiropas Kopienā                                                             | Avrupa Topluluğu'nda Yetkili Temsilci                                                          | Официальный представитель в Европейском сообществе                                                                                                | Ovlašćeni predstavnik u Evropskoj zajednici                                                 |
|   | Izgatavošanas datums                                                                               | Üretim tarihi                                                                                  | Дата производства                                                                                                                                 | Datum proizvodnje                                                                           |
|  | Ražotājs                                                                                           | Üretici                                                                                        | Производитель                                                                                                                                     | Proizvođač                                                                                  |
|  | Levēroiet lietošanas instrukciju                                                                   | Kullanım Talimatlarına uyun                                                                    | См. инструкции по применению                                                                                                                      | Pratite uputstva za upotrebu                                                                |
|   | levēroiet lietošanas instrukciju, kas pieejama tīmekļa vietnē                                      | Web sitesinde bulunan kullanım talimatlarını takip edin                                        | Следуйте инструкциям по применению на веб-сайте                                                                                                   | Pridržavajte se uputstva za upotrebu sa veb-sajta                                           |
|  | Uzmanību!                                                                                          | Dikkat                                                                                         | Предостережение                                                                                                                                   | Oprez                                                                                       |
|  | Nedrīkst lietot MR vidē                                                                            | MR için Güvenli Değildir                                                                       | Опасно при проведении МРТ                                                                                                                         | Nije bezbedno za MR                                                                         |
| <b>Rx only</b>                                                                      | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.        | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.    | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.                                                       | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. |
|  | Nav lietotāja apkopjamu detaļu                                                                     | İçinde bakım yapılabilir parça bulunmamaktadır                                                 | Не содержит обслуживаемых деталей                                                                                                                 | Ne sadrži delove koji se mogu servisirati                                                   |
|  | Uzglabāt sausu                                                                                     | Kuru halde tutun                                                                               | Хранить сухим                                                                                                                                     | Držite suvo                                                                                 |
|  | Trausls, rīkoties piesardzīgi                                                                      | Kırılabilir, dikkatli taşıyın                                                                  | Хрупкое изделие, обращаться с осторожностью                                                                                                       | Lomljivo, rukovati sa oprezom                                                               |
|  | Sargāt no saules gaismas                                                                           | Güneş ışığından uzak tutun                                                                     | Беречь от солнечных лучей                                                                                                                         | Držati zaklonjeno od sunčeve svetlosti                                                      |
| <b>IP4X</b>                                                                         | Aizsardzības līmenis pret priekšmetu iekļūšanu                                                     | Nesne girişine karşı koruma derecesi                                                           | Степень пылезащиты                                                                                                                                | Stepen zaštite od prodora objekata                                                          |
|  | Temperatūras ierobežojums                                                                          | Sıcaklık sınırı                                                                                | Ограничение по температуре                                                                                                                        | Ograničenje temperature                                                                     |

## Simbolu skaidrojums • Sembol Açıklamaları • Условные обозначения • Legenda sa simbolima

|                                                                                                                                                                                                                                                                                                                | Latviešu                                          | Türkçe                                             | Русский                                               | Srpski                                                                |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------|----------------------------------------------------|-------------------------------------------------------|-----------------------------------------------------------------------|
|                                                                                                                                                                                                                               | Mitruma ierobežojums                              | Nem sınırlaması                                    | Ограничение по влажности                              | Ograničenje vlažnosti                                                 |
|                                                                                                                                                                                                                               | Atmosfēras spiediena ierobežojums                 | Atmosfer basıncı sınırlaması                       | Ограничения по атмосферному давлению                  | Ograničenje atmosferskog pritiska                                     |
|                                                                                                                                                                                                                               | Conformité Européenne (CE marķējums)              | Conformité Européenne (CE İşareti)                 | Conformité Européenne (маркировка CE)                 | Conformité Européenne (CE znak)                                       |
|  eifu.edwards.com<br>+1 888 570 4016                                                                                                                                                                                           | Skatīt lietošanas instrukciju tīmekļa vietnē      | Web sitesindeki kullanım talimatlarına başvurun    | См. инструкции по применению на веб-сайте             | Pogledajte uputstva za upotrebu na veb-sajtu                          |
|                                                                                                                                                                                                                               | Importētājs                                       | İthalatçı                                          | Импортер                                              | Uvoznik                                                               |
|                                                                                                                                                                                                                               | Nelietot, ja iepakojums ir bojāts                 | Ambalaj hasarlıysa kullanmayın                     | Не использовать, если упаковка повреждена             | Ne koristite ako je pakovanje oštećeno                                |
|                                                                                                                                                                                                                               | Nesterils                                         | Steril Değildir                                    | Нестерильно                                           | Nesterilno                                                            |
|                                                                                                                                                                                                                               | Skatīt lietošanas instrukciju                     | Kullanım talimatlarına başvurun                    | См. инструкции по применению                          | Pogledajte uputstvo za upotrebu                                       |
|                                                                                                                                                                                                                             | Noņemiet spiediena kontrollera vāku no šī gala.   | Basınç kontrolörü kapağını bu taraftan çıkarın.    | Снимайте крышку регулятора давления с этого конца.    | Uklonite poklopac kontrolera pritiska počevši od ovog kraja.          |
|                                                                                                                                                                                                                             | Nenoņemiet spiediena kontrollera vāku no šī gala. | Basınç kontrolörü kapağını bu taraftan çıkarmayın. | Не снимайте крышку регулятора давления с этого конца. | Nemojte uklanjati poklopac kontrolera pritiska počevši od ovog kraja. |
| <b>Piezīme!</b> Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. • <b>Not:</b> Bu ürünün etiketinde tüm semboller yer almamış olabilir. • <b>Примечание.</b> На этикетках данного изделия могут указываться не все символы. • <b>Napomena:</b> Na ambalaži proizvoda se možda ne nalaze svi simboli. |                                                   |                                                    |                                                       |                                                                       |

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**Edwards Lifesciences Services GmbH**

Edisonstrasse 6  
85716 Unterschleissheim  
Germany

10/20  
10031915003 A / DOC-0142670  
© Copyright 2020, Edwards Lifesciences LLC  
All rights reserved

**Edwards Lifesciences LLC**

One Edwards Way  
Irvine, CA 92614 USA  
Made in USA

**Edwards Lifesciences SA**

Route de l'Etraz 70  
1260 Nyon, Switzerland

Telephone 949.250.2500  
800.424.3278  
FAX 949.250.2525

Web IFU