
Edwards

DIRECTORY

English	1	Português	28	Eesti	55
Français	4	Česky	31	Lietuvių	58
Deutsch	7	Magyar	34	Latviešu	61
Español	10	Polski	37	Türkçe	64
Italiano	13	Slovensky	40	Русский	67
Nederlands	16	Norsk	43	Srpski	70
Dansk	19	Suomi	46	简体中文	73
Svenska	22	Български	49	한국어	76
Ελληνικά	25	Română	52	日本語	79

English

Pressure Controller Kit

REF PC2K and PCCVR

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

For figures, see Figure 1 on page 82 through Figure 8 on page 83.

1.0 Description

The pressure controller kit consists of a pressure controller (PC2) and accompanying band (PC2B). The pressure controller can be used with the EV1000 clinical platform NI and HemoSphere advanced monitoring platform to noninvasively measure blood pressure and key hemodynamic parameters using a compatible Edwards finger cuff wrapped around the finger. The pressure controller is worn on the patient's wrist and connects to the compatible Edwards noninvasive monitoring system, heart reference sensor and compatible Edwards finger cuff(s) (Figure 1 on page 82).

A pressure controller cover is available as an accessory. The pressure controller cover secures the heart reference sensor into the pressure controller.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use.

Edwards, Edwards Lifesciences, the stylized E logo, EV1000, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

2.0 Intended Use

The pressure controller is intended for use with an Edwards compatible noninvasive monitoring system - composed of compatible monitor, pressure source (pump), heart reference sensor, compatible Edwards finger cuff(s) and pressure controller - for noninvasive measurement of blood pressure and associated hemodynamic parameters. Refer to the operator's manual of the compatible Edwards monitor being used for specific information on the intended use environment and patient population.

3.0 Instructions for Use

3.1 Apply the Pressure Controller

Step	Procedure
1	Wrap the pressure controller band around the patient's wrist (Figure 2 on page 82, step 1).
2	Snap the pressure controller into the plastic sleeve of the band, making sure that the cuff connectors are facing towards the fingers (Figure 2 on page 82, step 2).
3	Attach the pressure controller cable to the port on compatible Edwards noninvasive monitoring system (Figure 1 on page 82, 1A).
4	Remove the plastic connector caps in order to connect the finger cuff and heart reference sensor to the pressure controller according to the Instructions for Use supplied with the particular units (Figure 3 on page 82, steps 1-2), or consult the operator's manual of the EV1000 clinical platform NI and HemoSphere advanced monitoring platform. Note: It is recommended that the cuff connector caps be kept and used to protect the pressure controller against the ingress of water and dirt when only a single cuff is used (Figure 3 on page 82, step 3).
5	Once the heart reference sensor (HRS) is connected to the pressure controller, then place the pressure controller cover (PCCVR) on the pressure controller.

WARNING: Do not overtighten the pressure controller band.

WARNING: Do not apply pressure controller band on injured skin as this can cause further injury.

3.2 Remove the Pressure Controller Band

To remove the pressure controller from the pressure controller band:

Step	Procedure
1	Pull the sleeve slightly outwards (Figure 4 on page 82, step 1) and tilt the pressure controller to remove it from the sleeve (Figure 4 on page 82, step 2).

If you have a pressure controller cover (PCCVR), follow the instructions below to apply and remove the PCCVR.

3.3 Apply the Pressure Controller Cover

Step	Procedure
1	Ensure the heart reference sensor (HRS) is attached prior to attaching the pressure controller cover to the pressure controller.
2	Place the pressure controller cover's back mounting notch around the pressure controller cable (Figure 5 on page 83).
3	Snap the pressure controller cover over the pressure controller, making sure that the pressure controller cover does not interfere with the heart reference sensor (HRS) connection.

3.4 Remove the Pressure Controller Cover

Step	Procedure
1	Pull the pressure controller cover upwards from the front tab (Figure 6 on page 83) and tilt the pressure controller cover to remove it from the pressure controller (Figure 7 on page 83).

4.0 Cleaning Instructions

The pressure controller and pressure controller cover can be cleaned using the following disinfectants:

- 70% isopropyl alcohol solution
- 10% sodium hypochlorite water solution

Step	Procedure
1	Moisten a clean cloth with disinfectant and wipe the surfaces.
2	Dry the surface with a clean dry cloth.

CAUTION: Do not immerse the pressure controller.

CAUTION: Do not disinfect the pressure controller by autoclave or gas sterilization.

CAUTION: Do not immerse any cable connectors in fluid.

CAUTION: Do not pinch any heart reference sensor tubes or wires under the pressure controller cover during application. Be careful the only wire between the back mounting notch is the pressure controller cable.

CAUTION: Do not lift PCCVR from any other point than the front tab (Figure 8 on page 83).

CAUTION: Do not use any other cleaning agents aside from those listed.

The pressure controller band is intended for limited reuse. The operator shall assess whether reuse is appropriate. When reused, the same cleaning instructions can be followed as specified for the pressure controller.

The pressure controller cover is intended for limited reuse. The operator shall assess whether reuse is appropriate. When reused, follow the cleaning instruction above. Replace if damaged.

5.0 Accessories

Model Number	Accessory Description
PC2K	Pressure Controller Kit
PC2B	Pressure Controller Band (Multi Pack of 5)
PC2CCC	Pressure Controller Cuff Connector Caps (Multi Pack of 25)
PCCVR	Pressure Controller Cover

6.0 Warranty

The pressure controller is warranted for a period of one (1) year from the date of purchase when used in accordance with the directions for use. This warranty does not include the pressure controller band. For a full description of the warranty conditions refer to the operator's manual of the EV1000 clinical platform NI and HemoSphere advanced monitoring platform. Both the pressure controller and the pressure controller band are non-serviceable parts.

7.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada

(24 hours): 830.822.9837

Outside the U.S. and Canada

(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4

In Ireland: 01 8211012 - Option 4

This document is subject to change without prior notice. Refer to the latest version of the EV1000 clinical platform NI and HemoSphere advanced monitoring platform operator's manual for additional information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

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	English	Français	Deutsch	Español
	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo
	Separate collection for electrical and electronic equipment in accordance with EC Directive 2012/19/EU	Tri sélectif des équipements électriques et électroniques conformément à la Directive CE 2012/19/UE	Getrennte Sammlung von Elektro- und Elektronikgeräten gemäß EU-Richtlinie 2012/19/EU	Recogida por separado de equipos eléctricos y electrónicos según la directiva de la CE 2012/19/UE
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario
	Serial number	Numéro de série	Seriennummer	Número de serie
	Lot number	N° de lot	Chargenbezeichnung	Número de lote
	Quantity	Quantité	Menge	Cantidad
	Conformité Européenne (CE Mark) of TÜV SÜD Product Service GmbH (notified body)	Conformité européenne (marquage CE) de TÜV SÜD Product Service GmbH (organisme notifié)	Conformité Européenne (CE-Kennzeichnung) der TÜV SÜD Product Service GmbH (Benannte Stelle)	Conformité Européenne (Marca CE) de TÜV SÜD Product Service GmbH (organismo notificado)
	Authorized representative in the European Community	Représentant autorisé dans la Communauté européenne	Bevollmächtigter in der Europäischen Gemeinschaft	Representante autorizado en la Comunidad Europea
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación
	Manufacturer	Fabricant	Hersteller	Fabricante
	Follow instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung befolgen	Siga las instrucciones de uso
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web
	Caution	Avertissement	Vorsicht	Aviso
	MR Unsafe	Risques en milieu RM	MR-unsicher	No seguro para RM
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	No serviceable parts inside	Pas d'accessoires réparables à l'intérieur	Das Gerät enthält keine Teile, die gewartet werden können	No contiene piezas que requieran mantenimiento
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener en un lugar seco
	Fragile, handle with care	Fragile, manipuler avec précaution	Vorsicht! Zerbrechlich!	Frágil; manipular con cuidado
	Keep away from sunlight	Tenir à l'abri de la lumière du soleil	Vor Sonneneinstrahlung schützen	Mantener alejado de la luz del sol
IP4X	Extent of protection against the ingress of objects	Extension de la protection contre l'entrée d'objets	Erweiterter Schutz gegen feste Fremdkörper	Grado de protección contra la entrada de objetos
	Temperature limit	Limite de température	Temperaturgrenze	Límite de temperatura

Symbol Legend • Légende des symboles • Zeichenerklärung • Significado de los símbolos

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
	Български	Română	Eesti	Lietuvių
	Каталожен номер	Număr de catalog	Kataloogi number	Katalogo numeris
	Разделно събиране на електрическо и електронно оборудване според Директивата на ЕС 2012/19/ЕС	Colectarea separată a echipamentelor electrice și electronice, în conformitate cu directiva CE 2012/19/UE	Elektri- ja elektroonikaseadmete eraldi kogumine kooskõlas EÜ direktiiviga 2012/19/EL	Atskiras elektros ir elektroninės įrangos atliekų surinkimas pagal EB direktyvą 2012/19/ES
	Медицинско изделие	Dispozitiv medical	Meditsiiniseade	Medicinos prietaisas
	Сериен номер	Număr de serie	Seerianumber	Serijos numeris
	Партиден номер	Număr de lot	Partii number	Partijos numeris
	Количество	Cantitate	Kogus	Kiekis
	Conformité Européenne (CE маркировка) на TÜV SÜD Product Service GmbH (нотифициран орган)	Conformité Européenne (marcaj CE) al TÜV SÜD Product Service GmbH (organism notificat)	TÜV SÜD Product Service GmbH (teavitatud asutus) Conformité Européenne (CE-märk)	TÜV SÜD Product Service GmbH (notifikuotoji įstaiga) Conformité Européenne (CE žymė)
	Упълномощен представител в Европейската общност	Reprezentant autorizat în Comunitatea Europeană	Volitatud esindaja Euroopa Ühenduses	Įgaliotasis atstovas Europos Bendrijoje
	Дата на производство	Data fabricației	Tootmiskuupäev	Pagaminimo data
	Производител	Producător	Tootja	Gamintojas
	Следвайте инструкциите за употреба	Urmați instrucțiunile de utilizare	Järgige kasutusjuhendit!	Laikykitės naudojimo instrukcijų
	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiate veebisaidilt	Vadovaukitės naudojimo instrukcijomis, pateiktomis svetainėje
	Внимание	Atenție	Ettevaatust	Perspėjimas
	Небезопасно при MR	Incompatibil RM	Ohtlik magnetresonantstomograafias	MR nesaugus
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Във вътрешността няма части, подлежащи на сервизиране	Fără componente care pot fi reparate în interior	Ei sisalda hooldatavaid osasid	Viduje nėra keičiamų dalių
	Да се съхранява сухо	A se păstra uscat	Hoidke kuivana	Laikykite sausoje vietoje
	Чупливо, да се борава внимателно	Fragil, a se manipula cu atenție	Õrn, käsitseda ettevaatlikult	Trapus, elgtis atsargiai
	Да се пази от слънчева светлина	A se feri de lumina solară	Hoida eemal päikesevalgusest	Saugoti nuo saulės šviesos
	Степен на защита срещу проникване на предмети	Grad de protecție împotriva pătrunderii obiectelor	Objektide sissetungi eest kaitsmise määr	Objektų patekimo į prietaisą apsaugos lygis
	Температурно ограничение	Limită de temperatură	Temperatuuri piirang	Temperatūros apribojimas

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
	Latviešu	Türkçe	Русский	Srpski
	Kataloga numurs	Katalog numarası	Номер по каталогу	Kataloški broj
	Dalīta elektrisko un elektronisko iekārtu savākšana saskaņā ar ES direktīvu 2012/19/ES	2012/19/EU sayılı EC direktifi uyarınca elektrikli ve elektronik ekipmanlar ayrı toplanmalıdır	Раздельный сбор электрического и электронного оборудования в соответствии с директивой ЕС 2012/19/EU	Odvojeno sakupljanje električne i elektronske opreme u skladu sa direktivom EZ 2012/19/EU
	Medicīnas ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo
	Sērijas numurs	Seri numarası	Серийный номер	Serijski broj
	Partijas numurs	Lot numarası	Номер партии	Broj partije
	Daudzums	Miktar	Количество	Količina
	Conformité Européenne (CE marķējums), ko piešķir TÜV SÜD Product Service GmbH (paziņotā struktūra)	TÜV SÜD Product Service GmbH (onaylı kuruluş) Conformité Européenne (CE İşareti)	Соответствие требованиям директив Европейского Союза (маркировка CE) подтверждено компанией TÜV SÜD Product Service GmbH (нотифицированный орган)	Conformité Européenne (CE znak) kompanije TÜV SÜD Product Service GmbH (ovlašćeno telo)
	Pilnvarotais pārstāvis Eiropas Kopienā	Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе	Ovlašćeni predstavnik u Evropskoj zajednici
	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje
	Ražotājs	Üretici	Производитель	Proizvođač
	Levėrojiet lietošanas instrukciju	Kullanım Talimatlarına uyun	См. инструкции по применению	Pratite uputstva za upotrebu
	levėrojiet lietošanas instrukciju, kas pieejama tīmekļa vietnē	Web sitesinde bulunan kullanım talimatlarını takip edin	Следуйте инструкциям по применению на веб-сайте	Pridržavajte se uputstva za upotrebu sa web-sajta
	Uzmanību!	Dikkat	Предостережение	Oprez
	Nedrikst lietot MR vidē	MR için Güvenli Değildir	Опасно при проведении МРТ	Nije bezbedno za MR
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Nav lietotāja apkopjamu detaļu	İçinde bakım yapılabilir parça bulunmamaktadır	Не содержит обслуживаемых деталей	Ne sadrži delove koji se mogu servisirati
	Uzglabāt sausu	Kuru halde tutun	Хранить сухим	Držite suvo
	Trausls, rīkoties piesardzīgi	Kırılabilir, dikkatli taşıyın	Хрупкое изделие, обращаться с осторожностью	Lomljivo, rukovati sa oprezom
	Sargāt no saules gaismas	Güneş ışığından uzak tutun	Беречь от солнечных лучей	Držati zaklonjeno od sunčeve svetlosti
IP4X	Aizsardzības līmenis pret priekšmetu iekļūšanu	Nesne girişine karşı koruma derecesi	Степень пылезащиты	Stepen zaštite od prodora objekata
	Temperatūras ierobežojums	Sıcaklık sınırı	Ограничение по температуре	Ograničenje temperature

Simbolu skaidrojums • Sembol Açıklamaları • Условные обозначения • Legenda sa simbolima

	Latviešu	Türkçe	Русский	Srpski
	Mitruma ierobežojums	Nem sınırlaması	Ограничение по влажности	Ograničenje vlažnosti
	Atmosfēras spiediena ierobežojums	Atmosfer basıncı sınırlaması	Ограничения по атмосферному давлению	Ograničenje atmosferskog pritiska
	Conformité Européenne (CE marķējums)	Conformité Européenne (CE İşareti)	Conformité Européenne (маркировка CE)	Conformité Européenne (CE znak)
 eifu.edwards.com + 1 888 570 4016	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu
	Importētājs	İthalatçı	Импортер	Uvoznik
	Nelietot, ja iepakojums ir bojāts	Ambalaj hasarlıysa kullanmayın	Не использовать, если упаковка повреждена	Ne koristite ako je pakovanje oštećeno
	Nesterils	Steril Değildir	Нестерильно	Nesterilno
	Skatīt lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu
	Noņemiet spiediena kontrollera vāku no šī gala.	Basınç kontrolörü kapağını bu taraftan çıkarın.	Снимайте крышку регулятора давления с этого конца.	Uklonite poklopac kontrolera pritiska počevši od ovog kraja.
	Nenoņemiet spiediena kontrollera vāku no šī gala.	Basınç kontrolörü kapağını bu taraftan çıkarmayın.	Не снимайте крышку регулятора давления с этого конца.	Nemojte uklanjati poklopac kontrolera pritiska počevši od ovog kraja.
<p>Piezīme! Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. • Not: Bu ürünün etiketinde tüm semboller yer almamış olabilir. • Примечание. На этикетках данного изделия могут указываться не все символы. • Напомена: Na ambalaži proizvoda se možda ne nalaze svi simboli.</p>				

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EC REP

Edwards Lifesciences Services GmbH

Edisonstrasse 6
85716 Unterschleissheim
Germany

01/21
10031915003 B / DOC-0142670 B
© Copyright 2021, Edwards Lifesciences LLC
All rights reserved

Edwards Lifesciences LLC

One Edwards Way
Irvine, CA 92614 USA

Made in USA

Edwards Lifesciences SA

Route de l'Etraz 70
1260 Nyon, Switzerland

Telephone 949.250.2500

800.424.3278

FAX

949.250.2525

Web IFU