
DIRECTORY

English	1	Português	20	Eesti	38
Français	3	Česky	22	Lietuvių	40
Deutsch	5	Magyar	24	Latviešu	42
Español	7	Polski	26	Türkçe	44
Italiano	9	Slovensky	28	Русский	46
Nederlands	11	Norsk	30	Srpski	48
Dansk	13	Suomi	32	简体中文	50
Svenska	15	Български	34	한국어	52
Ελληνικά	17	Română	36	日本語	55

English

Acumen IQ Finger Cuff ClearSight Finger Cuff

REF AIQCS, AIQCM, AIQCL, CSC2S, CSC2M, CSC2L

Acumen IQ and ClearSight Finger Cuff

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

Description

The Acumen IQ and ClearSight finger cuffs, when used with an appropriate Edwards monitoring system, provide continuous, noninvasive hemodynamic monitoring. The Acumen IQ and ClearSight finger cuffs utilize the volume-clamp method to measure blood pressure with an inflatable bladder wrapped around the middle phalanx of the finger.

Device performance, including functional characteristics, have been verified

in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use

The ClearSight system provides continuous blood pressure measurement in a non-invasive manner.

Indications

The Acumen IQ and ClearSight finger cuffs are indicated for patients over 18 years of age to noninvasively measure blood pressure and associated hemodynamic parameters when used with EV1000 clinical platform NI or HemoSphere advanced monitoring platform.

Contraindications

In some patients with extreme contraction of the smooth muscle in the arteries and arterioles in the lower arm and hand, such as may be present in patients with Raynaud's disease, blood pressure measurement can become impossible.

Single Patient Use

The Acumen IQ and ClearSight finger cuffs are designed for single patient use. Upon starting measurement, the finger cuff can be used and re-applied for up to 72 hours on one patient. After 8 hours of continuous monitoring on a single finger, the finger cuff should be re-applied to another finger.

Do not attempt to clean and reuse the finger cuff on more than a single patient.

Warnings

Improper placement, sizing, or alignment of the finger cuff can lead to inaccurate monitoring.

Do not apply the finger cuff on injured skin as this may cause further injury.

To reduce the risk of skin irritation and tissue damage, do not monitor longer than 8 hours continuously on a single finger. To continue to monitor beyond 8 hours, use an additional

Edwards, Edwards Lifesciences, the stylized E logo, Acumen, Acumen IQ, ClearSight, EV1000, and HemoSphere are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

finger cuff on another finger or move the cuff in use to another finger.

Do not use two finger cuffs simultaneously on the same finger.

Measurement on one finger in contradiction with the instructions for use may affect patient comfort and/or lead to minor injuries.

Do not apply finger cuff(s) on a hand or finger when external constriction (that may prevent circulation to the hand or finger) is present.

Do not apply finger cuff(s) on a hand or finger when a second blood pressure measurement device is actively monitoring on the same arm (or hand or finger).

Do not use the finger cuff with magnetic resonance imaging.

Cautions

Do not use a damaged finger cuff. This may result in inaccurate measurements or may damage the Edwards monitoring system.

Never bend a finger cuff to a flat shape as it will damage the finger cuff and affect measurement accuracy.

Excessive ambient lighting may interfere with the finger cuff measurement. Avoid using the finger cuff under close, direct lighting.

Always disconnect the finger cuff when it is not applied to a finger, to prevent damage by accidental overinflation.

The effectiveness of the finger cuffs have not been established in pre-eclamptic patients.

Instructions for Use

Refer to pages 55-57 for figures corresponding to the steps below.

Size the Finger

1. Use the finger cuff sizing aid by wrapping it around the middle phalanx of the index, middle, or ring finger of the patient and pull the smaller end of the color-coded aid through the slot until snug (Figure 1 and Figure 2).
2. Based on the size indicated on the sizing aid, select the appropriate size finger cuff.

Apply the Finger Cuff

3. Gently open the finger cuff and place the finger cuff on the middle phalanx of the sized finger. Ensure the finger cuff is centered between the second and third knuckles, and the two green lines on the inside of the finger cuff (Figure 3).
NOTE: Do not apply the finger cuff on the thumb, small finger, or previously fractured fingers.
4. Align the distal end of the finger with the center line of the alignment tab on the finger cuff (Figure 3).

5. Allow the finger cuff to close around the finger (Figure 4). Ensure that the finger cuff remains aligned and does not rotate.
6. Check that the correct finger cuff size was selected by ensuring that the sizing window is entirely green (Figure 4-B).
NOTE: If there is red in the sizing window, the finger cuff is not the correct size for the selected finger. Repeat Steps 1 through 5 on the selected or an alternate finger to determine the correct size and apply the appropriate size finger cuff to the middle phalanx of that finger.
7. Remove the adhesive backing from the finger cuff and secure around the finger (Figure 4-A). Remove the adhesive backing from the distal end of the alignment tab and secure to the top of the finger cuff (Figure 4-C).
8. Lead the finger cuff cable between two fingers to the back side of the hand (Figure 5).

Connect the Finger Cuff to the System

9. Plug the finger cuff connector into the pressure controller (Figure 5). Refer to the pressure controller instructions for use for more details.
10. If using a Heart Reference Sensor (HRS), attach the finger end of the HRS to the HRS tab on top of the finger cuff (Figure 5). Refer to the HRS instructions for use for more details.
11. Initiate monitoring with the connected Edwards noninvasive monitoring system. To enable advanced parameters, refer to the monitoring system operator's manual for more details.
12. If using double cuff monitoring, repeat Steps 1 through 11 to apply the second finger cuff.

Recommendations

The Edwards noninvasive monitoring system will pause blood pressure measurements periodically per the monitoring system configuration settings. At this time, the fingertip should be checked for perfusion and potential tissue damage. The fingertip should also be checked periodically according to hospital protocol.

MRI Safety

This device is MR-unsafe and poses hazards in the MRI environment. This device contains metallic components, which can experience RF-induced heating in the MRI environment.

Storage Conditions

Store in a cool, dry place away from direct sunlight.
Temperature: -18 ° to 40 °C
Humidity: 10% to 90% RH

Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

Models

Model Number *	Model Description
AIQCS	Acumen IQ Finger Cuff Small
AIQCM	Acumen IQ Finger Cuff Medium
AIQCL	Acumen IQ Finger Cuff Large
CSC2S	ClearSight Finger Cuff Small
CSC2M	ClearSight Finger Cuff Medium
CSC2L	ClearSight Finger Cuff Large

* All models may not be available in all regions.

Warranty

The Acumen IQ and ClearSight finger cuffs are for single patient use only. The Acumen IQ and ClearSight finger cuffs are warrantied at time of delivery to the end user only. The Acumen IQ and ClearSight finger cuffs are a non-serviceable part.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
Inside the U.S. and Canada
(24 hours):800.822.9837
Outside the U.S. and Canada
(24 hours):949.250.2222
In the UK:0870 606 2040 - Option 4
In Ireland:01 8211012 Option 4


Refer to the latest version of the Edwards noninvasive monitoring system operator's manual for more information.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

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© Copyright 2020, Edwards Lifesciences LLC
All rights reserved.
Edwards Lifesciences Services GmbH
Edwards Lifesciences AG
Edwards Lifesciences LLC