
DIRECTORY

English	2	Português.....	16	Eesti	30
Français	3	Česky.....	18	Lietuvių.....	32
Deutsch.....	5	Magyar	19	Latviešu	33
Español.....	7	Polski	21	Türkçe.....	35
Italiano	8	Slovensky	23	Русский	36
Nederlands.....	10	Norsk.....	24	Srpski	38
Dansk	11	Suomi	26	中文	39
Svenska.....	13	Български.....	27	繁體中文 (台灣)	41
Ελληνικά.....	15	Română	29	한국어.....	42

FG1000

Fogarty Biliary Balloon Probe • Sonde à ballonnet pour voies biliaires Fogarty • Fogarty Gallengang-Ballonsonde • Sonda con balón biliar Fogarty • Sonda biliare a palloncino Fogarty • Fogarty -biliaire ballonsonde • Fogarty galdesystemballonprobe • Fogarty biliär ballongsond • Καθετήρας χοληφόρων με μπαλόνι Fogarty • Sonda com balão biliar Fogarty • Biliární balónková sonda Fogarty • Fogarty epeúti ballonos szonda • Sonda do dróg żółciowych Fogarty z balonem • Žilčvodová balóniková sonda Fogarty • Fogarty galleballongprobe • Fogarty -sappitiehytpallosondi • Жлъчна балонна сонда на Fogarty • Sondă biliară cu balon Fogarty • Fogarty sapiteede balloonsond • „Fogarty” balioninis tulžies lataų zondas • Fogarty žultsvadu balonzonde • Fogarty Safra Balon Probu • Билярный баллонный зонд Fogarty • Fogarty bilijarna balon sonda • Fogarty 胆管球囊探头 • Fogarty 膽管氣球探針 • Fogarty 담도 풍선 탐색자

1. Stylet • Stylet • Mandrin • Estilete • Stiletto • Stilet • Stilet • Mandräng • Στελεός • Estilete • Mandrén • Vezetőszonda • Mandryn • Mandrén • Stilett • Mandriini • Стилет • Stilet • Stilett • Vielinis kaištis • Stilets • Stile • Стилет • Mandren • 管心针 • 通管絲 • 스타일렛

2. Inflation Hub • Embout de gonflage • Aufdehnungsanschluss • Conector de inflado • Raccordo di gonfiaggio • Vulhub • Inflationsmuffe • Fyllningsfattning • Ομφαλός διόγκωσης • Extremidade de insuflação • Hrdlo pro plnění • Feltöltőnyílás • Złącze do napełniania balonu • Plniace hrdlo • Fyllingsnav • Täyttöportti • Адаптер за раздуване • Racord pentru umflare • Täitmisjaotur • Išplėtimo movinė jungtis • Uzpildišanas pieslēgvietā • Şişirme Göbeği • Соединитель для раздувания • Čvorište za naduvavanje • 充盈套节 • 充盈接頭 • 팍창 허브

3. Radiopaque Catheter Shaft • Corps de cathéter radio-opaque • Röntgendichter Katheterschaft • Eje de catéter radiopaco • Stelo del catetere radiopaco • Radiopake katheterschacht • Röntgenfast kateterschaft • Röntgentätt kateterschaft • Ακτινοοκιερό στέλεχος καθετήρα • Haste radiopaca do cateter • Rentgenkontrastní díř katétu • Sugárfotó katéterszár • Radiocieniujący trzon cewnika • Röntgenkontrastný driek katétra • Röntgentett kateterschaft • Röntgenpositiivinen katetrin varsi • Рентгеноконтрастен shaft на катетъра • Tijă cateter radioopacă • Röntgenkontrastne kateetri vars • Rentgenokonstrastinis kateterio vamzdelis • Rentgenstarojumu necaurīdīgā katetra ass • Radyoopak Kateter Shaft • Рентгеноконтрастный стержень катетера • Radiološki neprozirno telo katetera • 不透射线导管轴 • 不透射线导管轴 • 방사선 비투과성 카테터 샤프트

4. Latex Balloon • Ballonnet en latex • Latexballon • Balón de látex • Palloncino in lattice • Latexballon • Latexballon • Latexballong • Μπαλόνι από λάτεξ • Balão de látex • Latexový balóněk • Latexballon • Balon lateksowy • Latexový balónik • Lateksballong • Lateksipallo • Латексов балон • Balon din latex • Lateksballoon • Latekso balionelis • Lateksa balons • Lateks Balon • Латексный баллон • Balon od lateksa • 乳胶球囊 • 乳膠球囊 • 라텍스 풍선

Fogarty Biliary Balloon Probe

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

This device assists in the removal of hepatic calculi.

For Single Use Only

1.0 Contraindications

There are no known contraindications.

2.0 Indications

Fogarty biliary balloon probe is indicated for the removal of stones and ductal debris and for exploration.

3.0 Description

The Fogarty biliary balloon probe is a sterile, single-use probe. The device consists of a radiopaque flexible probe shaft with an integrated elastomeric latex balloon at the distal end. A hub at the proximal end is used for balloon inflation.

4.0 Warnings

The use of highly viscous or particulate contrast media is not recommended for balloon inflation as the lumen may become occluded.

The maximum recommended volume for balloon inflation must not be exceeded as over inflation increases the possibility of balloon rupture.

5.0 Precautions

As with all catheterization procedures, complications may occur. These include intrahepatic ductal damage, transient bleeding, postoperative cholangitis, and stricture of a primary hepatic duct. Instances of hepatic duct disruption during retrieval of intrahepatic stones have also been reported.

The balloon must be allowed to deflate as it passes through the ampullary region to avoid damage to the sphincter.

Instances of balloon rupture are minimized if care is taken to inflate the balloon slowly and withdraw the catheter gently.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

6.0 Instructions

6.1 Technique

The probe is introduced with the balloon deflated into the common duct through the choledochotomy. When the desired position has been reached, the balloon is inflated by the introduction of a sterile fluid. It may then be gently withdrawn for the removal of debris or manipulated for further exploration. During the period over which this technique has been developed, the Fogarty biliary balloon probe has been found to be a valuable adjunct to presently available techniques. The Fogarty biliary balloon probe does not, however, have utility in removing stones when the probe cannot pass beyond the calculi or when a calculus is in a sacculus.

The Biliary Balloon Probe aids in easy and rapid removal of amorphous intraductal debris and gravel. Numerous passes may be required to remove all material, but care should be taken not to over-distend the common duct during the period of manipulation. Closed-end irrigation of the common duct may be done by the introduction of the probe into the duodenum, followed by inflation and impingement of the inflated balloon against the sphincter mechanism. The irrigation catheter may then be employed. More effective washing is sometimes accomplished by closure of an open sphincter mechanism.

6.2 Atraumatic Calibration of the Sphincter Mechanism

The probe is introduced into the distal common duct and passed through the sphincter into the duodenum. Once in the duodenum, the balloon is inflated and then carefully palpated to make certain that it has passed through the sphincter mechanism. While the surgeon keeps the fingers of one hand placed on the balloon in the duodenum, the balloon probe is withdrawn until it impinges at the area of the sphincter. With light traction on the probe, the surgeon deflates the balloon until it passes through the sphincter with ease. As the balloon is passed through the sphincter, the quantity of fluid contained within the syringe is noted. When the balloon re-enters the common duct, the balloon is inflated until it engages the duct wall, and then gently withdrawn. The same amount of fluid which was noted as the probe passed through the sphincter is reintroduced into the balloon. Calibration of the balloon diameter accurately reflects the size of the sphincter. The

common duct can be calibrated in the same manner. Mobilization of the duodenum greatly facilitates this maneuver. This technique aids in identification of the retroduodenal and pancreatic position of the common duct and accurately locates the sphincter of Oddi.

6.3 For Removal of Intrahepatic Stones

The probe can be bent with the stylet in place to facilitate introduction into the right or left hepatic duct. The stylet may be removed after the probe enters into the hepatic duct. The probe is then advanced as far as possible with the balloon deflated. As the balloon is gently inflated, the probe is simultaneously removed. After noting when the probe engages the wall of the ductal system, the inflated balloon is carefully withdrawn through the incision in the common duct. Several passes may be necessary. The final exploratory passage in these radicles is carried out by insertion of probes into both the right and left hepatic ducts. The two probes are then simultaneously removed. This maneuver prevents an isolated stone from falling from one hepatic duct into the other during removal. Undue manipulation, probing, or over-distention of the balloon may cause bleeding and should be avoided.

6.4 Identification of the Common Duct

Where there is difficulty in identifying the entire extent of the common duct because of previous surgery or inflammation, the use of the probe allows identification of the duct and facilitates dissection. If the gall bladder remains and the cystic duct has been identified, a small incision in the cystic duct with introduction of the probe into the distal common duct is easily accomplished. Similarly, if the common duct has been identified by needle aspiration, a small puncture wound is made in the anterior surface of the common duct and the balloon probe is introduced. After entrance into the common duct, palpation of the inflated balloon in various locations facilitates dissection and identification of the common duct and associated anatomy.

6.5 Use of Stylet

The stylet supplied with the Biliary Balloon Probe is of a designed stiffness intended to avoid perforation and yet maintain enough shape to facilitate introduction of the probe into areas of sharp angulation.

When introducing the probe into the common duct, entrance into the duodenum may be accomplished with the stylet in place and angulated at approximately 30°. If difficulty is encountered in passing through the sphincter mechanism, the stylet should be unscrewed until it lies free of the probe body. A gentle probing motion with the stylet in this position will allow for passage into the duodenum if the sphincter mechanism is large enough to accommodate the probe body. When extreme pliability of the distal tip of the instrument is desired, the stylet should be withdrawn 1 cm. A probe possessing a longer

length of intermediate pliability is achieved by withdrawing the stylet a greater length.

7.0 Balloon Inflation

The balloon should be inflated with a sterile, blood compatible fluid that may include radiopaque solutions, if they are highly dilute and non particulate.

Before use, the probe with the balloon inflated to its maximum recommended volume should be inspected by the surgeon.

The amount of fluid in the syringe should be checked before each inflation. If the amount exceeds the balloon's stated capacity, remove the syringe from the probe and refill to the proper volume. The amount of fluid remaining in the catheter lumen must be taken into account.

8.0 MRI Information

This product has not been tested for MRI compatibility.

9.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

10.0 Storage

Store in a cool, dry place.
Temperature Limitation: 0° - 40 °C
Humidity Limitation: 5% - 90% RH.

11.0 Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

12.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
Inside the U.S. and Canada
(24 hours): 800.822.9837
Outside the U.S. and Canada
(24 hours): 949.250.2222

In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

13.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

STERILE EO

Specification Table:

French Size	5F (1.67 mm)	6F (2.0 mm)
Maximum Liquid Capacity (ml)	0.4	1.25
Maximum Air Capacity (ml)	0.8	2.5
Maximum Inflated Balloon Diameter (mm)	8	13

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	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk	Svenska	Ελληνικά
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte	Anvendelig længde	Brukbar längd	Ωφέλιμο μήκος
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad	Anbefalet guidewirestørrelse	Rekommenderad ledarstorlek	Συνιστώμενο μέγεθος οδήγου σύρματος
	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo	Numero di catalogo	Catalogusnummer	Katalognummer	Katalognummer	Αριθμός καταλόγου
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße des Einführbestecks	Tamaño de introductor mínimo	Misura minima dell'introduttore	Minimum introducergrootte	Minimum introducerstørrelse	Minsta storlek på införare	Ελάχιστο μέγεθος εισαγωγέα
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op	Forsigtig	Var försiktig	Προσοχή
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken	Må ikke genanvendes	Får inte återanvändas	Μην επαναχρησιμοποιείτε
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid	Mængde	Antal	Ποσότητα
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer	Partnummer	Lotnummer	Αριθμός Παρτίδας
	Authorized representative in the European Community	Représentant autorisé dans la Communauté européenne	Bevollmächtigter in der Europäischen Gemeinschaft	Representante autorizado en la Comunidad Europea	Rappresentante autorizzato nella Comunità Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap	Autoriseret repræsentant i Det Europæiske Fællesskab	Auktoriserad representant inom Europeiska gemenskapen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση
	Do not use if package is damaged	Ne pas utiliser si le conditionnement est endommagé	Bei beschädigter Verpackung nicht verwenden	No lo utilice si el envase está dañado	Non utilizzare se la confezione è danneggiata	Niet gebruiken als de verpakking is beschadigd	Må ikke anvendes, hvis emballagen er beskadiget	Använd inte om förpackningen är skadad	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum	Anvendes inden	Sista förbrukningsdag	Ημερομηνία λήξης
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethylenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit	Ballonkapacitet	Ballongkapacitet	Χωρητικότητα μπαλονιού
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant	Producent	Tillverkare	Κατασκευαστής
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig	Indeholder eller har spor af naturlig gummilatex	Innehåller eller har spår av naturgummilatex	Περιέχει λάτεξ φυσικού καουτσούκ
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr	Medicinsk utrustning	Ιατροτεχνολογικό προϊόν
	Size	Taille	Größe	Tamaño	Misura	Afmetingen	Størrelse	Storlek	Μέγεθος

**Symbol Legend • Légende des symboles • Zeichenerklärung • Significado de los símbolos • Legenda dei simboli
• Lijst met symbolen • Symbolforklaring • Symbolförklaring • Υπόμνημα συμβόλων**

	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk	Svenska	Ελληνικά
	Consult instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Se bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο
	Temperature limit	Limite de température	Temperaturgrenze	Limite de temperatura	Limite di temperatura	Temperatuurlimiet	Temperaturgrænse	Temperaturgräns	Όριο θερμοκρασίας
	Humidity limitation	Limites d'humidité	Zulässiger Luftfeuchtigkeitsbereich	Limitación de humedad	Limiti di umidità	Vochtigheidsbeperkingen	Luftfugtighedsbegrænsning	Luftfuktighetsbegränsning	Αποδεκτό εύρος τιμών υγρασίας
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο
	Conformité Européenne (CE Mark) of TÜV SÜD Product Service GmbH (notified body)	Conformité européenne (marquage CE) de TÜV SÜD Product Service GmbH (organisme notifié)	Conformité Européenne (CE-Kennzeichnung) der TÜV SÜD Product Service GmbH (Benannte Stelle)	Conformité Européenne (Marca CE) de TÜV SÜD Product Service GmbH (organismo notificado)	Conformité Européenne (Marchio CE) di TÜV SÜD Product Service GmbH (ente notificato)	Conformité Européenne (CE-markering) van TÜV SÜD Product Service GmbH (aangemelde instantie)	Conformité Européenne (CE-mærke) fra TÜV SÜD Product Service GmbH (bemyndiget organ)	Conformité Européenne (CE-märke) för TÜV SÜD Product Service GmbH (anmält organ)	Conformité Européenne (σήμανση CE) της TÜV SÜD Product Service GmbH (κοινοποιημένος οργανισμός)
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό
	Importer	Importateur	Importeur	Importador	Importatore	Importeur	Importør	Importör	Εισαγωγέας

Note: Not all symbols may be included in the labeling of this product. • **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. • **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. • **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. • **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. • **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. • **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. • **Obs!** Alla symboler kanske inte används för märkning av denna produkt. • **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος.

Legenda de símbolos • Legenda se Symboly • Jelmagyarázat • Legenda simboli • Vysvetlivky k symbolom
• Symbolforklaring • Merkkien selitykset • Легенда на символите • Legendă de simboluri

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	中文	繁體中文	한국어
	Kasutatav pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjljiva dužina	可用长度	可用長度	사용 가능한 길이
	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Tavsiye edilen kılavuz tel boyutu	Рекомендованный размер проводника	Preporučena veličina žice vodice	推荐的导丝尺寸	建議的導絲大小	가이드와이어 권장 크기
	Kataloogi number	Katalogo numeris	Kataloga numurs	Katalog numarası	Номер по каталогу	Kataloški broj	目录编号	型錄編號	카탈로그 번호
	Sisesti väikseim suurus	Minimalus įvediklio dydis	Minimalais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	最小插管器尺寸	最小導引器尺寸	삽입기 최소 크기
	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	小心	注意	주의 사항
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	不得重复使用	單次使用	재사용 금지
	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	数量	數量	수량
	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	批号	批號	로트 번호
	Volitatud esindaja Euroopa Ühenduses	Įgaliotasis atstovas Europos Bendrijoje	Pilnvarotais pārstāvis Eiropas Kopienā	Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе	Ovlašćeni predstavnik u Evropskoj zajednici	欧洲共同体内授权代表	歐洲共同體授權代表	유럽 공동체 공인 대리점
	Ärge kasutage, kui pakend on kahjustatud.	Nenaudoti, jei pakuotė pažeista	Nelietot, ja iepakojums ir bojāts	Ambalaj hasarlıysa kullanmayın	Не использовать, если упаковка повреждена	Ne koristiti ako je pakovanje oštećeno	如果包装损坏, 请勿使用	包裝有損壞現象時，請勿使用	포장이 손상되었을 경우 사용하지 마십시오.
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	有效期	使用截止日期	사용 기한
	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	采用环氧乙烷气体灭菌	使用環氧乙烷 (Ethylene Oxide) 進行消毒	에틸렌옥사이드를 사용하여 멸균됨
	Ballooni mahutavus	Balioniėlio Talpa	Balona letilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	球囊容量球	球體容量	풍선 용량
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	制造商	製造商	제조업체
	Tootmiskuupäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	制造日期	製造日期	제조일자
	Sisaldab naturaalsed kummilekseit	Turi savo sudėtyje ar yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	含有或存在天然乳胶	含天然橡膠乳膠成分	천연고무 라텍스가 포함되었거나 들어 있습니다.
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	查阅网站上的使用说明	請詳閱網站上的使用說明	웹사이트에 있는 사용 지침을 참조하십시오.
	Meditsiiniseade	Medicinos prietaisas	Medicīnas ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	医疗设备	醫療器械	의료 장치
	Suurus	Dydis	Izmērs	Boyut	Размер	Veličina	尺寸	尺寸	크기

Sümbolite seletus • Simbolių paaiškinimas • Simbolu skaidrojums • Sembol Açıklamaları • Основные обозначения
• Legenda sa simbolima • 符号图例 • 符號圖例 • 기호 범례

	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	中文	繁體中文	한국어
	Vaadake kasutusjuhendit	Žr. naudojimo instrukcijas	Skatīt lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu	查阅使用说明	請詳閱使用說明	사용 지침 참조
	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	存放于阴凉干燥处	存放於陰涼乾燥處	서늘하고 건조한 장소에 보관하십시오.
	Temperatuuripiirang	Temperatūros apribojimas	Temperatūras ierobežojums	Sıcaklık sınırı	Ограничение по температуре	Ograničenje temperature	温度限制	溫度限制	온도 한계
	Niiskuspriirang	Drėgnio apribojimas	Mitruma ierobežojums	Nem sınırlaması	Ограничение по влажности	Ograničenje vlažnosti	湿度限制	濕度限制	습도 제한
	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	请勿重复灭菌	請勿重複滅菌。	재멸균하지 마십시오.
	Mittepüroogenne	Nepirogeniškas	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	无致热源	非熱原性	비발열
	TÜV SÜD Product Service GmbH (teavitatud asutus) Conformité Européenne (CE-märk)	TÜV SÜD Product Service GmbH (notifikuotoji įstaiga) Conformité Européenne (CE žymė)	Conformité Européenne (CE marķējums), ko piešķir TÜV SÜD Product Service GmbH (paziņotā struktūra)	TÜV SÜD Product Service GmbH (onaylı kuruluş) Conformité Européenne (CE İşareti)	Соответствие требованиям директив Европейского Союза (маркировка CE) подтверждено компанией TÜV SÜD Product Service GmbH (нотифицированный орган)	Conformité Européenne (CE znak) kompanije TÜV SÜD Product Service GmbH (ovlaščeno telo)	TÜV SÜD Product Service GmbH (认证机构) 的Conformité Européenne (CE 标志)	TÜV SÜD Product Service GmbH (認證機構) 的 Conformité Européenne (CE 標誌)	TÜV SÜD Product Service GmbH (인증기관)의 Conformité Européenne (CE 마크)
	Hoida kuivalt	Laikykite sausoje vietoje	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	保持干燥	保持乾燥	건조한 상태로보관하십시오.
	Importija	Importuotojas	Importētājs	İthalatçı	Импортёр	Uvoznik	进口商	進口商	수입업체

Märkus. Kõik sümbolid ei pruugi antud toote etiketil esineda. • **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. • **Piezīme:** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. • **Not:** Bu ürünün etiketinde tüm semboller yer almayış olabilir. • **Примечание.** На этикетках данного изделия могут указываться не все символы. • **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. • **注释:** 此产品标签并非包含所有符号。• **註:** 本產品標示中未必包含所有符號。• **참고:** 본 제품의 라벨에 모든 기호가 표시된 것은 아닙니다.

Edwards

EC REP

Edwards Lifesciences Services GmbH

Edisonstrasse 6
85716 Unterschleissheim
Germany

09/20
10000132005 A/DOC-0145537 A
© Copyright 2020, Edwards Lifesciences LLC
All rights reserved

Edwards Lifesciences LLC

One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences SA

Route de l'Etraz 70
1260 Nyon, Switzerland

Telephone

949.250.2500

FAX

800.424.3278

949.250.2525

Web IFU