
DIRECTORY

English	1	Português.....	15	Eesti	28
Français	3	Česky.....	16	Lietuvių.....	30
Deutsch.....	4	Magyar	18	Latviešu	31
Español.....	6	Polski	19	Türkçe	33
Italiano.....	7	Slovensky	21	Русский	34
Nederlands.....	9	Norsk	22	Srpski	36
Dansk	10	Suomi	24	中文	37
Svenska.....	12	Български.....	25	繁體中文（台灣）.....	38
Ελληνικά.....	13	Română	27	한국어.....	40


FG1002

Fogarty Arterial Embolectomy Catheter · Cathéter Fogarty pour embolctomie artérielle · Fogarty Katheter zur arteriellen Embolctomie · Catéter para embolctomia arterial Fogarty · Catetere per embolctomia arteriosa Fogarty · Fogarty · katheter voor arteriële embolctomie · Fogarty arteriealembolctomikateeter · Fogarty arteriell embolctomikateeter · Каѳтеризатор артериальной эмболизатор · Fogarty · Cateter de embolctomia arterial Fogarty · Arteriální embolctomický katétr Fogarty · Fogarty artériás embolctomias kateéter · Cevmik Fogarty do embolctomii tetniczej · Arteriálny embolctomický katéter Fogarty · Fogarty arteriell embolctomikateeter · Valtimossa käytettävä Fogarty · embolctomikatetri · Катетер за артериална емболиктомия на Fogarty · Cateter de embolctomia arterială Fogarty · Fogarty arteriaalembolctomia kateeter · „Fogarty” arteriell embolctomijos kateteris · Fogarty artériálas embolctomijás katets · Fogarty Arteriell embolctomikateeter · Артериальный катетер Fogarty для эмболизатор · Fogarty kateter za arterijsku embolctomiju · Fogarty 动脉取栓导管 · Fogarty 動脈取栓術導管 · Fogarty 동맥 색전 제거용 카테터

1. Styet (applicable models only) - Stylet (modèles applicables uniquement) - Mandrin (nur zutreffende Modelle) - Estilete (solo en los modelos correspondientes) - Stiletto (solo nei modelli applicabili) - Stilet (alleen toepasselijke modellen) - Stilet (kun geldende modellen) - Mandriäng (endast relevanta modeller) - Στελέκος (μόνο για τα πουθενά που ισχύει) - Estilete (apenas modelos aplicáveis) - Mandrin (pouze příslušné modely) - Vezetőszonda (csak megfelelő típusok esetén) - Mandryn (tylko w niektórych modelach) - Mandrin (iba príslušné modely) - Stilet (bare aktuelle modeller) - Mandrini (vain yhteensopivissa mallissa) - Гинет (само за приложимыми моделями) - Stilet (numai pentru modelele aplicabile) - Stilet (aniżeli uchwytliwy modely) - Vezetőszonda (csak megfelelő típusok esetén) - Stilet (tikai piemērojamiem modeļiem) - Stile (sadece geçerli modellerle) - Гинет (только для применимых моделей) - Mandren (samo za važeće modele) - 管心针 (仅适用型号) - 通管針 (僅能適用型號) - 스티일렛 (해당 모델에만 사용됨)
2. Inflation hub - Embout de gonflage - Aufblähungsnuss - Conector de inflado - Raccordo di gonfiaggio - Vuliub - Inflationsmuffe - Fyllningsfattning - Ομφαλός διόγκωσης - Extremité de insuflação - Hrdlo pro plnění - Feltöltőnyílás - Złącze do napełniania balonu - Plniace hrdlo - Fyllningsnavn - Täyrytpörrti - Äännetar za razduvanje - Racord pentru umflare - Tűltetésigép - Isplétimo monolit jungtis - Uzpildināšanas pieslēgvieta - Szjime Göbeği - Соединитель для раздувания - Чворісте за наддування - 充盈套管 - 팍창 허브
3. Catheter Shaft - Corps du cathéter - Katheterschaft - Eje del catéter - Stelo del catetere - Katheterschaft - Katheterschaft - Καθετήρας κοίτηςτρο - Haste do cateter - Dřik katétru - Kátetárszár - Trzon cewnika - Driek katetra - Katheterschaft - Katetrin varsi - Шафт на катетра - Tija cateter - Käteetri varsi - Käteetorio vamizdelis - Katetra ass - Kateter shaft - Крепежный катетра - Telo katetera - 导管轴 - 導管軸 - 카테터 샤프트
4. Latex ballon - Ballonnet en latex - Latexballon - Balón de látex - Palloncino in lattice - Latexballon - Μολών από λάτεξ - Balão de látex - Latexový balónek - Latexballon - Balon lateksowy - Latexový balónek - Lateksballong - Lateksipallo - Латексов ганон - Balon din latex - Lateksballoon - Lateksoball - Lateksballon - Латеко балон - Balon od lateksa - 乳膠球囊 - 乳膠球囊 - 라텍스 풍선

FG1003

Fogarty Arterial Embolectomy Catheter with Gate Valve - Cathéter Fogarty pour embolctomie artérielle avec robinet-vanne - Fogarty Katheter zur arteriellen Embolctomie mit Schieberventil - Catéter para embolctomia arterial Fogarty con válvula de compuerta - Catetere per embolctomia artetiosa Fogarty con valvula di chiusura - Fogarty -katheteri voor arteriële embolctomie met poortklep - Fogarty arterieembolctomiekatheter met schijfventiel - Fogarty arterial embolctomiekatheter met doorventiel - Катетер артериальной эмболектомии с затвором заслончатого типа - Катетер de embolctomia arterial Fogarty con Válvula de Corredor - Arteriální embolctomický katétér Fogarty s uzavíracím ventilem - Fogarty artériás embolctomális katéter kapuzsázeleppel - Cewnik Fogarty do embolctomie katetizacji z zaworem zasuwowym - Arteriálny embolctomický katéter Fogarty so vstupným ventilom - Fogarty arteriell embolctomikatheter med sluseventil - Valtimoissa käytettävä Fogarty -embolctomiakateetri ja luistiventili - Катетер за артериальна емболектомия на Fogarty с шибер - Catetere de embolctomie arterială Fogarty cu valvă obturator - Fogarty arteriaalembolctomia kateteer (silbiriga) - „Fogarty” arterien embolctomijos kateteris su sklende - Fogarty artériás embolctomijos kateteris su aroslegivartu - Giriş Kapaklı Fogarty Arteriyel Embolctomik Kateteri - Артериальный катетер Fogarty для эмболектомии с запорным клапаном - Fogarty kateter za arterijsku embolctomiju sa ulaznim ventilom - 帶閘閥的 Fogarty 動脈取栓技術導管 - 帶閘閥的 Fogarty 動脈取栓技術導管 - 5.1.7.1. 벨브가 있는 Fogarty 동맥 색전제거용 카테터

- | |
|--|
| 1. Inflation bulb with gate Valve - Embout de gonflage avec robinet-vanne - Aufblasungsanschluss mit Schieberventil - Conector de inflado con válvula de compuerta - Racordo di gonfiaggio con valvola di chiusura - Vutub mit poortklep - Inflationsmuffe mit skydeventil - Fyllningsfätrning med portventil - Ουφλοκό βούλκυ με βαλβίδα πύλης - Extremite de soufflage avec vanne - Hrdlo pro plnění s uzavíracím ventilem - Feltöltőnyílás kapuscsatlékol - Złącze do napełniania balonu z zaworem zasuwowym - Płnící hrdlo so vstupním ventilem - Fyllningsfätrning med skruventil - Таўтоўртніа балуны з вентілем - Арматура за раздування с шпіб-р - Raccord pentru umplere cu valvă obturator - Taitmijastuur sibiriga - Ispletimo movine plinosti sa sklence - Uzdplidšanas pieslēgvietā ar noslēgvārstu - Giriş Kapaklı Şişme Gübgesi - Соєдинителъ для раздування с запорным клапаном - Чворіте за надување са улазним ventilom - 帶閘閘的充盈套节 - |
| 2. Catheter Shaft - Corps du cathéter - Katheterschaft - Eje del cateter - Stelo del catetere - Katheterschacht - Kateterschaft - Катетерський коштирпо - Hastе del cateter - Driik katetru - Katerschaz - Trzon cewnika - Driek katetru - Kateterschaft - Katetrin varsi - Шафт на катетра - Tijā cateter - Katetari var - Kateterio vamzdis - Katetra ass - Kateter shaft - Срежняа катетра - Telo katetara - 号管轴 - 号管軸 - 카테터 샤프트 |
| 3. Latex Balloon - Ballonet en latex - Latexballon - Balón de látex - Palloncino in lattice - Latexballon - Latexballon - Млэжовою аю катэц - Balón de látex - Latexový balónek - Latexballon - Balon lateksowy - Latexový balónek - Latexballon - Latexspallo - Латековс балон - Balon din latex - Latexsballoon - Lateksó balonides - Latexs balon - Латековсый балон - Balon од латека - 乳膠球囊 - 乳膠球囊 - 라텍스 풍선 |

Fogarty Arterial Embolectomy Catheter

CAUTION: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For Single Use Only

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

This device aids in clot management.

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions For Use.

1.0 Indications

The Fogarty arterial embolectomy catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels in the arteries of the non-central circulatory system.

To remove fibrous or adherent material, alternative devices such as the Fogarty adherent clot and graft thrombectomy catheter are recommended.

2.0 Description

The Fogarty arterial embolectomy catheter is a sterile, single-use catheter. The device consists of a catheter shaft with an integrated elastomeric latex balloon and an atraumatic tip that is inserted surgically into arterial vessels of the non-central circulatory system. A hub at the proximal end is used for balloon inflation.

3.0 Contraindications

- The arterial embolectomy catheter should not be used outside the arterial system.
- The arterial embolectomy catheter is not recommended for the removal of fibrous, adherent, or calcified material (e.g. chronic clot, atherosclerotic plaque). The catheter is not designed to withstand the additional pull force needed to remove these materials.
- The arterial embolectomy catheter is not intended for use as a vessel dilator.
- The arterial embolectomy catheter should not be used in the venous system.

4.0 Warnings

- Balloon rupture and catheter separation as a result of excessive pull force applied to remove adherent material are the most frequent causes of reported failures. The possibility of balloon rupture must be taken into account when considering the risks involved in any embolectomy procedure.
- To minimize the risk of vessel damage, balloon rupture, or tip detachment, do not exceed the maximum recommended inflation volume and pull force for each size catheter (see **Specification Table**).
- Use of a highly viscous or particulate contrast medium is not recommended for balloon inflation because the inflation lumen may become occluded.
- To avoid air embolus in case of balloon rupture, air should not be used for balloon inflation. Carbon dioxide is the only recommended gas. See Instructions.

5.0 Instructions

5.1 Prepare

Using sterile technique, remove the cap or cap and stiffening stylet from the catheter hub.

5.2 Purge

Inflate the balloon with sterile fluid or gas to the maximum recommended volume. (It is recommended that the 2F (0.67 mm) catheter be inflated with carbon dioxide gas.) Pull a vacuum on the syringe. Repeat until all air is removed.

Note: For all inflations use the smallest syringe capable of holding the stated maximum fluid capacity.

5.3 Inspect

The catheter should be inspected with the balloon inflated during purging. A balloon that does not inflate, leaks, or inflates in a grossly asymmetric (eccentric) manner should not be used.

CAUTION: The amount of fluid in the syringe should be checked before each inflation. Do not exceed the recommended maximum inflation volume. See the specification table.

5.4 Place

With the balloon deflated, insert the catheter into the vessel and beyond the obstructing material.

5.5 Inflate

With the catheter suitably positioned, inflate the balloon with sterile fluid or gas. Balloon inflation should be stopped when the balloon can be felt to engage the arterial wall.

Note: The 2F (0.67 mm) catheter should be inflated with carbon dioxide gas. Gas seepage through the intact balloon will require frequent adjustment of the inflation volume.

CAUTION: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest inflated balloon diameter that will remove the obstructing material.

CAUTION: Air should not be used for inflation in instances where balloon rupture could produce a dangerous air embolus.

5.6 Withdraw

Remove the occlusive material by gently withdrawing the catheter. During withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by controlling the inflation volume.

CAUTION: Do not exceed the Maximum Recommended Pull Force. (See Specification Table).

CAUTION: Inflation of the balloon is associated with a feeling of resistance to the fluid or gas injection. When no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

6.0 MRI Information

This product has not been tested for MRI compatibility.

7.0 Complications

As with all catheterization procedures, complications may occur. These may include local or systemic infection, local hematomas, intimal disruption, arterial dissection, perforation and vessel rupture, hemorrhage, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation, and balloon rupture with fragmentation, tip separation and distal embolization.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

8.0 How Supplied

Contents sterile and nonpyrogenic if package is unopened and undamaged. Do not use if package is opened or damaged.

9.0 Storage

Store in a cool, dry place. Temperature Limitation: 0° - 40 °C, Humidity Limitation: 5% - 90% RH.

10.0 Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.

Note: Reprocessing or resterilization will not extend the indicated shelf life.

11.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone

numbers:
 Inside the U.S. and Canada
 (24 hours): 800.822.9837
 Outside the U.S. and Canada
 (24 hours): 949.250.2222
 In the UK: 0870 606 2040 - Option 4
 In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

12.0 Disposal

After patient contact, treat the device as biohazardous waste. Dispose of in accordance with hospital policy and local regulations.

WARNING: This device is designed, intended and distributed for SINGLE USE ONLY. DO NOT RE-STERILIZE OR REUSE this device. There are no data to support the sterility, non-pyrogenicity and functionality of the device after

reprocessing. Such action could lead to illness or an adverse event as the device may not function as originally intended.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Specification Table:

Catheter French Size (mm)	2F	3F	4F	5F	6F	7F
	0.67	1.00	1.33	1.67	2.0	2.3
Maximum Liquid Capacity (ml)	N/A	0.2	0.75	1.5	2.0	2.5
Maximum Gas Capacity (ml)	0.2	0.6	1.7	3.0	4.5	5.0
Maximum Recommended Pull Force on Inflated Balloon (lbs.)	0.5	0.7	1.5	2.0	2.5	3.5
(kg)	0.23	0.3	0.68	0.9	1.13	1.59
Diameter of Inflated Balloon (mm)	4	5	9	11	13	14
Maximum French Size of Deflated Balloon (mm)	3.4F	4.3F	5.0F	6.0F	7.0F	8.0F
	1.13	1.43	1.67	2.0	2.3	2.7

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	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk	Svenska	Ελληνικά
	Usable length	Longueur utile	Nutzlänge	Longitud útil	Lunghezza utile	Bruikbare lengte	Anvendelig længde	Brukbar längd	Ωφέλιμο μήκος
	Recommended guidewire size	Taille de fil-guide recommandée	Empfohlene Führungsdrahtgröße	Tamaño de guía recomendado	Misura filo guida consigliata	Aanbevolen maat voerdraad	Anbefalet guidewirestørrelse	Rekommenderad ledarstorlek	Συνιστώμενο μέγεθος οδήγου σύρματος
	Catalogue Number	Référence catalogue	Artikelnummer	Número de catálogo	Numero di catalogo	Catalogusnummer	Katalognummer	Katalognummer	Αριθμός καταλόγου
	Minimum introducer size	Taille minimale de l'introducteur	Mindestgröße des Einführbestecks	Tamaño de introductor mínimo	Misura minima dell'introduttore	Minimum introducergrootte	Minimum introducerstørrelse	Minsta storlek på införare	Ελάχιστο μέγεθος εισαγωγέα
	Caution	Avertissement	Vorsicht	Aviso	Attenzione	Let op	Forsigtig	Var försiktig	Προσοχή
	Do not re-use	Ne pas réutiliser	Nicht wiederverwenden	No reutilizar	Non riutilizzare	Niet hergebruiken	Må ikke genanvendes	Får inte återanvändas	Μην επαναχρησιμοποιείτε
	Quantity	Quantité	Menge	Cantidad	Quantità	Hoeveelheid	Mængde	Antal	Ποσότητα
	Lot Number	N° du lot	Chargenbezeichnung	Número de lote	Numero di lotto	Lotnummer	Partnummer	Lotnummer	Αριθμός Παρτίδας
	Authorized representative in the European Community	Représentant autorisé dans la Communauté européenne	Bevollmächtigter in der Europäischen Gemeinschaft	Representante autorizado en la Comunidad Europea	Rappresentante autorizzato nella Comunità Europea	Geautoriseerde vertegenwoordiger in de Europese Gemeenschap	Autoriseret repræsentant i Det Europæiske Fællesskab	Auktoriserad representant inom Europeiska gemenskapen	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση
	Do not use if package is damaged	Ne pas utiliser si le conditionnement est endommagé	Bei beschädigter Verpackung nicht verwenden	No lo utilice si el envase está dañado	Non utilizzare se la confezione è danneggiata	Niet gebruiken als de verpakking is beschadigd	Må ikke anvendes, hvis emballagen er beskadiget	Använd inte om förpackningen är skadad	Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά
	Use-by date	Date d'expiration	Verwendbar bis	Fecha de caducidad	Utilizzare entro	Vervaldatum	Anvendes inden	Sista förbrukningsdag	Ημερομηνία λήξης
	Sterilized using ethylene oxide	Stérilisé à l'oxyde d'éthylène	Mit Ethylenoxid sterilisiert	Esterilizado con óxido de etileno	Sterilizzato con ossido di etilene	Gesteriliseerd met behulp van ethylenoxide	Steriliseret ved brug af ethylenoxid	Steriliserad med etylenoxid	Αποστειρωμένο με χρήση οξειδίου του αιθυλενίου
	Balloon capacity	Capacité du ballonnet	Ballonkapazität	Capacidad del balón	Capacità del palloncino	Balloncapaciteit	Ballonkapacitet	Ballongkapacitet	Χωρητικότητα μπαλονιού
	Manufacturer	Fabricant	Hersteller	Fabricante	Produttore	Fabrikant	Producent	Tillverkare	Κατασκευαστής
	Date of manufacture	Date de fabrication	Herstellungsdatum	Fecha de fabricación	Data di produzione	Fabricagedatum	Fremstillingsdato	Tillverkningsdatum	Ημερομηνία κατασκευής
	Contains or presence of natural rubber latex	Contient ou présence de latex de caoutchouc naturel	Enthält Naturlatex	Contiene o hay presencia de látex de caucho natural	Contiene o è presente lattice di gomma naturale	Bevat natuurlijk rubberlatex of natuurlijk rubberlatex is aanwezig	Indeholder eller har spor af naturlig gummilatex	Innehåller eller har spår av naturgummilatex	Περιέχει λάτεξ φυσικού καουτσούκ
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult instructions for use on the website	Consulter le mode d'emploi sur notre site Web	Gebrauchsanweisung auf der Website beachten	Consulte las instrucciones de uso en el sitio web	Consultare le istruzioni per l'uso sul sito Web	Raadpleeg de gebruiksaanwijzing op de website	Se brugsanvisningen på webstedet	Se bruksanvisningen på webbplatsen	Συμβουλευτείτε τις οδηγίες χρήσης που αναφέρονται στην τοποθεσία web
	Medical device	Dispositif médical	Medizinprodukt	Producto sanitario	Dispositivo medico	Medisch hulpmiddel	Medicinsk udstyr	Medicinsk utrustning	Ιατροτεχνολογικό προϊόν
	Size	Taille	Größe	Tamaño	Misura	Afmetingen	Størrelse	Storlek	Μέγεθος

**Symbol Legend • Légende des symboles • Zeichenerklärung • Significado de los símbolos • Legenda dei simboli
• Lijst met symbolen • Symbolforklaring • Symbolförklaring • Υπόμνημα συμβόλων**

	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk	Svenska	Ελληνικά
	Consult instructions for use	Consulter le mode d'emploi	Gebrauchsanweisung beachten	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Se bruksanvisningen	Συμβουλευτείτε τις οδηγίες χρήσης
	Store in a cool, dry place	Conserver dans un endroit frais et sec	Das Produkt an einem kühlen, trockenen Ort lagern	Guárdese en un lugar fresco y seco	Conservare in un luogo fresco e asciutto	Op een koele en droge plaats bewaren	Skal opbevares køligt og tørt	Förvara svalt och torrt.	Φυλάσσετε σε δροσερό και ξηρό χώρο
	Temperature limit	Limite de température	Temperaturgrenze	Limite de temperatura	Limite di temperatura	Temperatuurlimiet	Temperaturgrænse	Temperaturgräns	Όριο θερμοκρασίας
	Humidity limitation	Limites d'humidité	Zulässiger Luftfeuchtigkeitsbereich	Limitación de humedad	Limiti di umidità	Vochtigheidsbeperkingen	Luftfugtighedsbegrænsning	Luftfuktighetsbegränsning	Αποδεκτό εύρος τιμών υγρασίας
	Do not resterilize	Ne pas restériliser	Nicht resterilisieren	No volver a esterilizar	Non risterilizzare	Niet opnieuw steriliseren	Må ikke resteriliseres	Får inte omsteriliseras	Μην επαναποστειρώνετε
	Non-pyrogenic	Apyrogène	Nicht pyrogen	No pirogénico	Apirogeno	Niet-pyrogeen	Ikke-pyrogen	Icke-pyrogen	Μη πυρετογόνο
	Conformité Européenne (CE Mark) of TÜV SÜD Product Service GmbH (notified body)	Conformité européenne (marquage CE) de TÜV SÜD Product Service GmbH (organisme notifié)	Conformité Européenne (CE-Kennzeichnung) der TÜV SÜD Product Service GmbH (Benannte Stelle)	Conformité Européenne (Marca CE) de TÜV SÜD Product Service GmbH (organismo notificado)	Conformité Européenne (Marchio CE) di TÜV SÜD Product Service GmbH (ente notificato)	Conformité Européenne (CE-markering) van TÜV SÜD Product Service GmbH (aangemelde instantie)	Conformité Européenne (CE-mærke) fra TÜV SÜD Product Service GmbH (bemyndiget organ)	Conformité Européenne (CE-märke) för TÜV SÜD Product Service GmbH (anmält organ)	Conformité Européenne (σήμανση CE) της TÜV SÜD Product Service GmbH (κοινοποιημένος οργανισμός)
	Keep dry	Tenir au sec	Vor Nässe schützen	Mantener seco	Mantenere asciutto	Droog houden	Opbevares tørt	Förvaras torrt	Διατηρείτε στεγνό
	Importer	Importateur	Importeur	Importador	Importatore	Importeur	Importør	Importör	Εισαγωγέας

Note: Not all symbols may be included in the labeling of this product. • **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. • **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. • **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. • **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. • **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. • **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne. • **Obs!** Alla symboler kanske inte används för märkning av denna produkt. • **Σημείωση:** Ενδέχεται να μην συμπεριλαμβάνονται όλα τα σύμβολα στη σήμανση αυτού του προϊόντος.

**Legenda de símbolos • Legenda se Symboly • Jelmagyarázat • Legenda simboli • Vysvetlivky k symbolom
• Symbolforklaring • Merkkien selitykset • Легенда на символите • Legendă de simboluri**

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	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	中文	繁體中文	한국어
	Kasutatav pikkus	Naudingasis ilgis	Izmantojamais garums	Kullanılabilir uzunluk	Рабочая длина	Upotrebjliva dužina	可用长度	可用長度	사용 가능한 길이
	Juhtetraadi soovituslik suurus	Rekomenduojamas kreipiamosios vielos dydis	Ieteicamais vadītājstīgas izmērs	Tavsiye edilen kılavuz tel boyutu	Рекомендованный размер проводника	Preporučena veličina žice vodice	推荐的导丝尺寸	建議的導絲大小	가이드와이어 권장 크기
	Kataloogi number	Katalogo numeris	Kataloga numurs	Katalog numarası	Номер по каталогу	Kataloški broj	目录编号	型錄編號	카탈로그 번호
	Sisesti väikseim suurus	Minimalus įvediklio dydis	Minimālais ievadītāja izmērs	Minimum introdüser boyutu	Минимальный размер интродьюсера	Minimalna veličina uvodnika	最小插管器尺寸	最小導引器尺寸	삽입기 최소 크기
	Ettevaatust	Perspėjimas	Uzmanību!	Dikkat	Предостережение	Oprez	小心	注意	주의 사항
	Mitte korduskasutada	Nenaudoti pakartotinai	Nelietot atkārtoti	Yeniden kullanmayın	Не подлежит повторному использованию	Ne koristiti ponovo	不得重复使用	單次使用	재사용 금지
	Kogus	Kiekis	Daudzums	Miktar	Количество	Količina	数量	數量	수량
	Partii number	Partijos numeris	Partijas numurs	Lot Numarası	Номер партии	Broj partije	批号	批號	로트 번호
	Volitatud esindaja Euroopa Ühenduses	Įgaliotasis atstovas Europos Bendrijoje	Pilnvarotais pārstāvis Eiropas Kopienā	Avrupa Topluluğu'nda Yetkili Temsilci	Официальный представитель в Европейском сообществе	Ovlašćeni predstavnik u Evropskoj zajednici	欧洲共同体内授权代表	歐洲共同體授權代表	유럽 공동체 공인 대리점
	Ärge kasutage, kui pakend on kahjustatud.	Nenaudoti, jei pakuotė pažeista	Nelietot, ja iepakojums ir bojāts	Ambalaj hasarlıysa kullanmayın	Не использовать, если упаковка повреждена	Ne koristiti ako je pakovanje oštećeno	如果包装损坏, 请勿使用	包裝有損壞現象時，請勿使用	포장이 손상되었을 경우 사용하지 마십시오.
	Kõlblik kuni	Naudoti iki nurodytos datos	Derīguma termiņš	Son Kullanma Tarihi	Использовать до	Rok korišćenja	有效期	使用截止日期	사용 기한
	Steriliseeritud etüleenoksiidi kasutades	Sterilizuota etileno oksidu	Sterilizēts ar etilēna oksidu	Etilen oksit kullanılarak sterilize edilmiştir	Стерилизовано этиленоксидом	Sterilizovano etilen-oksidom	采用环氧乙烷气体灭菌	使用環氧乙烷 (Ethylene Oxide) 進行消毒	에틸렌옥사이드를 사용하여 멸균됨
	Ballooni mahutavus	Balionių Talpa	Balona letilpība	Balon Kapasitesi	Рабочий объем баллона	Kapacitet balona	球囊容量球	球體容量	풍선 용량
	Tootja	Gamintojas	Ražotājs	Üretici	Производитель	Proizvođač	制造商	製造商	제조업체
	Tootmiskuupäev	Pagaminimo data	Izgatavošanas datums	Üretim tarihi	Дата производства	Datum proizvodnje	制造日期	製造日期	제조일자
	Sisaldab naturaalseid kummilekseit	Turi savo sudėtyje ar yra natūralios gumos latekso	Satur dabiskā kaučuka lateksu lielā vai mazā daudzumā	Doğal kauçuk lateks içerir veya mevcuttur	Содержит натуральный латекс или изготовлено из натурального латекса	Sadrži ili ima prisustva prirodne lateks gume	含有或存在天然乳胶	含天然橡膠乳膠成分	천연고무 라텍스가 포함되었거나 들어 있습니다.
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
 eifu.edwards.com + 1 888 579 4016	Palun lugege veebisaidil olevaid kasutusjuhiseid!	Žr. naudojimo instrukcijas, pateiktas interneto svetainėje	Skatīt lietošanas instrukciju tīmekļa vietnē	Web sitesindeki kullanım talimatlarına başvurun	См. инструкции по применению на веб-сайте	Pogledajte uputstva za upotrebu na veb-sajtu	查阅网站上的使用说明	請詳閱網站上的使用說明	웹사이트에 있는 사용 지침을 참조하십시오.
	Meditsiiniseade	Medicinos prietaisas	Medicīnas ierīce	Tıbbi cihaz	Медицинское устройство	Medicinsko sredstvo	医疗设备	醫療器械	의료 장치
	Suurus	Dydis	Izmērs	Boyut	Размер	Veličina	尺寸	尺寸	크기

Sümbolite seletus • Simbolių paaiškinimas • Simbolu skaidrojums • Sembol Açıklamaları • Основные обозначения
• Legenda sa simbolima • 符号图例 • 符號圖例 • 기호 범례

	Eesti	Lietuvių	Latviešu	Türkçe	Русский	Srpski	中文	繁體中文	한국어
	Vaadake kasutusjuhendit	Žr. naudojimo instrukcijas	Skatīt lietošanas instrukciju	Kullanım talimatlarına başvurun	См. инструкции по применению	Pogledajte uputstvo za upotrebu	查阅使用说明	請詳閱使用說明	사용 지침 참조
	Säilitage jahedas ja kuivas kohas.	Laikyti vėsioje, sausoje vietoje	Uzglabāt vēsā, sausā vietā	Serin, kuru yerde saklayın	Хранить в прохладном, сухом месте	Čuvati na hladnom i suvom mestu	存放于阴凉干燥处	存放於陰涼乾燥處	서늘하고 건조한 장소에 보관하십시오.
	Temperatuuripiirang	Temperatūros apribojimas	Temperatūras ierobežojums	Sıcaklık sınırı	Ограничение по температуре	Ograničenje temperature	温度限制	溫度限制	온도 한계
	Niiskuspriirang	Drėgnio apribojimas	Mitruma ierobežojums	Nem sınırlaması	Ограничение по влажности	Ograničenje vlažnosti	湿度限制	濕度限制	습도 제한
	Ärge resteriliseerige	Kartotinai nesterilizuokite	Nesterilizēt atkārtoti	Yeniden sterilize etmeyin	Не стерилизовать повторно	Nemojte ponovo sterilisati	请勿重复灭菌	請勿重複滅菌。	재멸균하지 마십시오.
	Mittepüroogenne	Nepirogeniškas	Nepirogēns	Nonpirogeniktir	Апирогенно	Apirogeno	无致热源	非熱原性	비발열
	TÜV SÜD Product Service GmbH (teavitatud asutus) Conformité Européenne (CE-märk)	TÜV SÜD Product Service GmbH (notifikuotoji įstaiga) Conformité Européenne (CE žymė)	Conformité Européenne (CE marķējums), ko piešķir TÜV SÜD Product Service GmbH (paziņotā struktūra)	TÜV SÜD Product Service GmbH (onaylı kuruluş) Conformité Européenne (CE İşareti)	Соответствие требованиям директив Европейского Союза (маркировка CE) подтверждено компанией TÜV SÜD Product Service GmbH (нотифицированный орган)	Conformité Européenne (CE znak) kompanije TÜV SÜD Product Service GmbH (ovlaščeno telo)	TÜV SÜD Product Service GmbH (认证机构) 的Conformité Européenne (CE 标志)	TÜV SÜD Product Service GmbH (認證機構) 的 Conformité Européenne (CE 標誌)	TÜV SÜD Product Service GmbH (인증기관)의 Conformité Européenne (CE 마크)
	Hoida kuivalt	Laikykite sausoje vietoje	Uzglabāt sausu	Kuru halde tutun	Беречь от влаги	Držite suvo	保持干燥	保持乾燥	건조한 상태로보관하십시오.
	Importija	Importuotojas	Importētājs	İthalatçı	Импортёр	Uvoznik	进口商	進口商	수입업체

Märkus. Kõik sümbolid ei pruugi antud toote etiketil esineda. • **Pastaba.** Šio gaminio etiketėje gali būti pateikti ne visi simboliai. • **Piezīme:** Šī izstrādājuma marķējumā var nebūt iekļauti visi simboli. • **Not:** Bu ürünün etiketinde tüm semboller yer almayış olabilir. • **Примечание.** На этикетках данного изделия могут указываться не все символы. • **Napomena:** Na ambalaži proizvoda se možda ne nalaze svi simboli. • **注释:** 此产品标签并非包含所有符号。• **註:** 本產品標示中未必包含所有符號。• **참고:** 본 제품의 라벨에 모든 기호가 표시된 것은 아닙니다.

Edwards

EC REP

Edwards Lifesciences Services GmbH

Edisonstrasse 6
85716 Unterschleissheim
Germany

09/20
10000186006 A/DOC-0145541 A
© Copyright 2020, Edwards Lifesciences LLC
All rights reserved

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in USA

Edwards Lifesciences SA
Route de l'Etraz 70
1260 Nyon, Switzerland

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU