
Edwards

ForeSight Jr

English

Non Adhesive Small Sensor

REF FSESNS

< 8 kg (cerebral), < 5 kg (non-cerebral) – Pediatric

Carefully read these instructions for use, which address the warnings, precautions, and residual risks for this medical device.

1.0 Description

The sensor, when used in combination with the ForeSight oximeter cable or in combination with the FORE-SIGHT ELITE absolute tissue oximeter, is a single use applied part that measures hemoglobin allowing the clinician to continuously and accurately determine absolute levels of blood oxygenation saturation in the tissue (StO₂).

Device performance, including functional characteristics, have been verified in a comprehensive series of testing to support the safety and performance of the device for its intended use when used in accordance with the established Instructions for Use.

The benefits of hemodynamic monitoring may lead to improved data driven clinical decision making for medically necessary intervention and/or clinical re-evaluation.

2.0 Indications for Use

When used in conjunction with the FORE-SIGHT ELITE absolute tissue oximeter or ForeSight oximeter cable:

The small sensor is indicated for monitoring of absolute regional hemoglobin oxygen saturation of blood under the sensor in individuals at risk for reduced flow or no-flow ischemic states. It is intended for cerebral use on pediatric subjects < 8 kg and non-cerebral use on pediatric subjects < 5 kg.

3.0 Contraindications

The sensor is contraindicated for use on patients:

- With a physical site area too limited for proper sensor placement
- Undergoing an MRI scan because of associate risk of injury

4.0 Warnings

- Assess the sensor site at least every 12 hours, or more often as required by the institution's protocol.
- Remove the sensor if the circulatory condition or skin integrity has deteriorated.
- Do not attach the sensor to damaged or irritated skin.
- Do not lay patient on the sensor or cable.
- Do not attach the sensor to skin with unapproved devices, such as head bands, hats, wraps, etc.
- Do not place the sensor or accessories over eyes, nose, or mouth.
- Do not cut the sensor. Cutting the sensor can result in injury to the patient.
- Do not use in an MRI environment.
- When used in settings with LED lighting, sensors should be covered with a light blocker prior to connection to the preamp cable, as some high intensity systems can interfere with the sensor's near infrared light detection.
- The use of barrier film dressings other than Tegaderm could affect the accuracy of the StO₂ readings.

5.0 Cautions

- Avoid positioning the sensor over hair, air sinus, hematoma, birthmark, or broken skin.
- Avoid attaching the sensor to sites with excess adipose, ascites, or edema.
- The materials used in the manufacture of the sensor were NOT designed for reuse. Reuse can cause the sensor not to perform as intended.

6.0 Instructions

Step	Procedure
1	Remove the sensor from the package. Carefully inspect the sensor for damage. Discard and replace if damage is found.
2	Select sensor location on the monitor.
3	Clean and dry the sensor site (Figure 1).

Step	Procedure
4	Apply the sensor to the patient: a) Cerebral Use (Figure 2a): Select the site on the forehead well above the eyebrow and just below the hairline. b) Non-cerebral Use (Figures 2b, 2c): Select the site that provides ideal access to the desired tissue, for example, latissimus dorsi (flank), external oblique or abdomen.
5	Size and cut the sensor band length to fit the patient (Figure 3). Shorten the sensor band away from the patient. Attach the sensor band to the patient with the print facing out. Do not cut the sensor band while on the patient. Do not cut any other part of the sensor. Note: Do not over tighten the sensor band as pressure can be transferred to the baby. Note: You may use Tegaderm between the sensor and skin in patients with delicate skin or edema.
6	Insert the sensor connector straight into the sensor cable connector until it snaps into place (Figure 4). Use the bedsheet clip to secure the cable and prevent pulling on the sensor.
7	If needed, fold the sensor flat cable to route it in the desired direction (Figure 4).

For use only with ForeSight oximeter cable software version 2.5.7 or above and with FORE-SIGHT ELITE absolute tissue oximeter software version 4.5.6 or above.

7.0 Disposal

Sensors are designed for single-patient use, and are not to be reprocessed. Use a new sensor for each patient and discard it after use. Disposal should follow in accordance with local hospital and institution policies. Re-used sensors present a risk of cross-contamination or infection.

8.0 Storage

Store in a cool, dry place.
Temperature Limitation: -20 ° - 50 °C
Humidity Limitation: 15% - 90% RH

9.0 Shelf Life

The recommended shelf life is as marked on each package. Storage beyond the expiration date may result in product deterioration.

10.0 Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada (24 hours): 800.822.9837
Outside the U.S. and Canada (24 hours): 949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 - Option 4

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Users and/or patients should report any serious incidents to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Refer to the symbol legend at the end of this document.

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	English	Français	Deutsch	Español	Italiano	Nederlands	Dansk
	Temperature limit	Limite de température	Temperaturgrenze	Límite de temperatura	Limiti di temperatura	Temperatuurlimiet	Temperaturgrænse
	Humidity limitation	Limites d'humidité	Zulässiger Luftfeuchtigkeitsbereich	Limitación de humedad	Limiti di umidità	Vochtigheidsbeperkingen	Luftfugtighedsbegrænsning
	Atmospheric pressure limitation	Limites de pression atmosphérique	Zulässiger Luftdruckbereich	Limitación de presión atmosférica	Limiti di pressione atmosferica	Atmosferische drukkbeperkingen	Atmosfærisk trykbegrænsning
	Non-sterile	Non stérile	Unsteril	No estéril	Non sterile	Niet-steriel	Ikke-steril
	Not made with natural rubber latex	Non fabriqué à partir de latex de caoutchouc naturel	Nicht aus Latex hergestellt	No contiene látex de caucho natural.	Non realizzato in lattice di gomma naturale	Bevat geen natuurlijk rubberlatex	Ikke fremstillet med naturligt gummilatex
	Follow instructions for use on the website	Consulter le mode d'emploi sur le site Web	Siehe Gebrauchsanweisung auf Website	Siga las instrucciones de uso del sitio web	Seguire le istruzioni per l'uso sul sito Web	Volg de gebruiksaanwijzing op de website	Følg brugsanvisningen på hjemmesiden
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Importer	Importateur	Importeur	Importador	Importatore	Importeur	Importør

Note: Not all symbols may be included in the labeling of this product. • **Remarque :** il est possible que certains symboles n'apparaissent pas sur les étiquettes de ce produit. • **Hinweis:** Unter Umständen sind nicht alle Zeichen auf dem Produktetikett aufgeführt. • **Nota:** Es posible que no todos los símbolos aparezcan en el etiquetado de este producto. • **Nota:** alcuni simboli potrebbero non essere stati inseriti sull'etichetta del prodotto. • **Opmerking:** Het label van dit product bevat misschien niet alle symbolen. • **Bemærk:** Alle symbolerne er muligvis ikke inkluderet på produktmærkaterne.

**Symbolförklaring • Υπόμνημα συμβόλων • Legenda dos símbolos • Legenda se symboly
• Jelmagyarázat • Legenda simboli • Vysvetlivky k symbolom**

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	Norsk	Suomi	Български	Română	Eesti	Lietuvių	Latviešu
	Temperaturgrense	Lämpötilaraja	Температурно ограничение	Limită de temperatură	Temperatuuripiirang	Temperatūros apribojimas	Temperatūras ierobežojums
	Fuktighetsbegränsning	Kosteusrajoitus	Ограничение за влажността	Limita de umiditate	Niiskuspiirang	Drėgnio apribojimas	Mitruma ierobežojums
	Atmosfærisk trykkgrensning	Ilmanpainerajat	Ограничение за атмосферното налягане	Limite presiune atmosferică	Atmosfäärilise rõhu piirang	Ribinė atmosferos slėgio reikšmė	Atmosfēras spiediena ierobežojums
	Usteril	Epästeriili	Нестерилно	Nesteril	Mittersteriilne	Nesterilus	Nesterils
	Fremstilt uten naturgummilateks	Ei sisällä luonnonkumilateksia	Не е произведено с естествен каучуков латекс.	Nu este fabricat cu latex din cauciuc natural	Pole valmistatud looduslikust kummilateksist.	Pagaminta be natūralios gumos latekso	Izgatavošanā nav izmantots dabiskā kaučuka latekss
	Følg bruksanvisningen på nettsiden	Noudata käyttöohjeita, jotka ovat verkkosivustolla	Спазвайте инструкциите за употреба, дадени на уебсайта	Respectați instrucțiunile de utilizare disponibile pe site-ul web	Järgige kasutusjuhendit, mille leiute veebisaidilt	Vadovautis naudojimo instrukcijomis, pateikiamomis svetainėje	Levėrojietai lietošanas instrukciju, kas pieejama tīmekļa vietnē
Rx only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Importør	Maahantuoja	Вносител	Importator	Importija	Importuotojas	Importētājs
Merk: Alle symboler er kanskje ikke inkludert på produktmerkingen. • Huomautus: kaikkia symboleja ei vältämättä ole käytetty tämän tuotteen pakkausmerkinnöissä. • Забележка: Възможно е не всички символи да са включени в етикета на този продукт. • Notă: Este posibil ca nu toate simbolurile să fie incluse pe eticheta acestui produs. • Mārkus. Kōik sūmbolid ei pruugi antud toote etiketil esineda. • Pastaba. Šio gaminio etiketėje gali būti pateikti ne visi simboliai. • Piezīme. Ši izstrādājuma marķējumā var nebūt iekļauti visi simboli.							

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CE 0123

03/21
10042621001 A / DOC-0147803 A
© Copyright 2021, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in China

Edwards Lifesciences SA
Route de l'Etraz 70
1260 Nyon, Switzerland

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

Web IFU

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Edwards Lifesciences Services GmbH